

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH2 2+) -5.50

logP (XH +) -0.05 ±0.07 (n=50) logP (neutral X) 1.93 ±0.01 (n=50)

18C-01005 Points 1 to 35

M15_octanol concentration factor 0.875
Carbonate 0.2232 mM
Acidity error -0.23667 mM

18C-01005 Points 36 to 71

M15_octanol concentration factor 0.909
Carbonate 0.1837 mM
Acidity error -0.20987 mM

18C-01005 Points 72 to 105

M15_octanol concentration factor 0.947
Carbonate 0.0524 mM
Acidity error -0.05926 mM

Warnings and errors

Errors None

Warnings One or more logP values out of range

Sample logD and percent species

Comment
Stomach pH
<u>'</u>
Blood pH
В



Assay ID:


Sample name: M15_octanol Assay name:

pH-metric high logP 18C-01005

Experiment start time: 3/1/2018 6:16:08 AM

Analyst: Pion Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r






Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-01005 Points 1 to 35

Overall results

RMSD 0.176
Average ionic strength 0.159 M
Average temperature 25.0°C
Partition ratio 0.0208 : 1

Analyte concentration range 3855.0 µM to 4000.4 µM

Total points considered 23 of 35

Warnings and errors

Errors None

Warnings One or more logP values out of range

Four-Plus parameters

à	Alpha	0.130	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r
à	S	0.9970	3/1/2018 6:16:08 AM	C:\Sirius T3\HCl18B27.t3r
à	jΗ	8.0	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r
à	jОН	-0.4	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r


Titrants


0.50 M HCI 0.993513 3/1/2018 6:16:08 AM C:\Sirius_T3\HCl18B27.t3r 0.50 M KOH 0.999845 3/1/2018 6:16:08 AM C:\Sirius_T3\KOH18B27.t3r

Sample

Φ	M15_octanol concentration factor	0.875
	Base pKa 1	2.58
	Base pKa 2	5.30
	logP (XH2 2+)	-5.50
₩	logP (XH +)	-3.48
₩	logP (neutral X)	1.91

Sample graphs








Assay ID:

Sample name: M15_octanol Experiment start time: 3/1/2018 6:16:08 AM Assay name:

pH-metric high logP Analyst: Pion Instrument ID: T312060 18C-01005

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r


Sample graphs (continued)



Sample logD and percent species


рН	M15_octanol logD	M15_octanol M15_octanolH2	M15_octanol M15 octanolH	M15_octanol M15_octanol	M15_octanol M15_octanolH2*	M15_octanol M15_octanolH*	M15_octanol M15_octanol*	
1.000	-3.93	97.44 %	2.56 %	0.00 %	0.00 %	0.00 %	0.00 %	
1.200	-3.56	96.00 %	4.00 %	0.00 %	0.00 %	0.00 %	0.00 %	Stomach pH
2.000	-2.07	79.15 %	20.82 %	0.01 %	0.00 %	0.00 %	0.02 %	
3.000	-0.53	27.28 %	71.75 %	0.36 %	0.00 %	0.00 %	0.61 %	
4.000	0.58	3.24 %	85.24 %	4.27 %	0.00 %	0.00 %	7.24 %	
5.000	1.43	0.16 %	42.47 %	21.28 %	0.00 %	0.00 %	36.08 %	
6.000	1.83	0.00 %	6.89 %	34.54 %	0.00 %	0.00 %	58.56 %	
6.500	1.89	0.00 %	2.29 %	36.25 %	0.00 %	0.00 %	61.46 %	
7.000	1.90	0.00 %	0.73 %	36.83 %	0.00 %	0.00 %	62.44 %	
7.400	1.91	0.00 %	0.29 %	36.99 %	0.00 %	0.00 %	62.71 %	Blood pH
8.000	1.91	0.00 %	0.07 %	37.07 %	0.00 %	0.00 %	62.85 %	·
9.000	1.91	0.00 %	0.01 %	37.10 %	0.00 %	0.00 %	62.89 %	
10.000	1.91	0.00 %	0.00 %	37.10 %	0.00 %	0.00 %	62.90 %	
11.000	1.91	0.00 %	0.00 %	37.10 %	0.00 %	0.00 %	62.90 %	
12.000	1.91	0.00 %	0.00 %	37.10 %	0.00 %	0.00 %	62.90 %	


Carbonate and acidity



0.223 mM Acidity error -0.237 mM

Other graphs








Assay name: pH-metric high logP Analyst: Pion 18C-01005 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-01005 Points 36 to 71

Overall results

RMSD 0.148
Average ionic strength 0.166 M
Average temperature 25.0°C
Partition ratio 0.0576 : 1

Analyte concentration range 3420.0 µM to 3556.1 µM

Total points considered 27 of 36

Warnings and errors

Errors None

Warnings One or more logP values out of range

Four-Plus parameters

à	Alpha	0.130	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r
à	S	0.9970	3/1/2018 6:16:08 AM	C:\Sirius T3\HCl18B27.t3r
à	jΗ	8.0	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r
à	jОН	-0.4	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r


Titrants


0.50 M HCI 0.993513 3/1/2018 6:16:08 AM C:\Sirius_T3\HCl18B27.t3r 0.50 M KOH 0.999845 3/1/2018 6:16:08 AM C:\Sirius_T3\KOH18B27.t3r

Sample

₩	M15_octanol concentration factor	0.909
	Base pKa 1	2.58
	Base pKa 2	5.30
	logP (XH2 2+)	-5.50
₩	logP (XH +)	-4.43
₩	logP (neutral X)	1.93

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	
	logD	M15_octanolH2	M15_octanolH	M15_octanol	M15_octanolH2*	M15_octanolH*	M15_octanol*	
1.000	-3.95	97.44 %	2.56 %	0.00 %	0.00 %	0.00 %	0.00 %	
1.200	-3.56	96.00 %	4.00 %	0.00 %	0.00 %	0.00 %	0.00 %	Stomach pH
2.000	-2.05	79.13 %	20.81 %	0.01 %	0.00 %	0.00 %	0.05 %	•
3.000	-0.51	26.97 %	70.93 %	0.36 %	0.00 %	0.00 %	1.74 %	
4.000	0.59	2.85 %	74.97 %	3.76 %	0.00 %	0.00 %	18.42 %	
5.000	1.45	0.10 %	25.24 %	12.65 %	0.00 %	0.00 %	62.01 %	
6.000	1.85	0.00 %	3.27 %	16.39 %	0.00 %	0.00 %	80.34 %	
6.500	1.90	0.00 %	1.06 %	16.76 %	0.00 %	0.00 %	82.18 %	
7.000	1.92	0.00 %	0.34 %	16.89 %	0.00 %	0.00 %	82.78 %	
7.400	1.93	0.00 %	0.13 %	16.92 %	0.00 %	0.00 %	82.94 %	Blood pH
8.000	1.93	0.00 %	0.03 %	16.94 %	0.00 %	0.00 %	83.03 %	
9.000	1.93	0.00 %	0.00 %	16.94 %	0.00 %	0.00 %	83.05 %	
10.000	1.93	0.00 %	0.00 %	16.94 %	0.00 %	0.00 %	83.06 %	
11.000	1.93	0.00 %	0.00 %	16.94 %	0.00 %	0.00 %	83.06 %	
12.000	1.93	0.00 %	0.00 %	16.94 %	0.00 %	0.00 %	83.06 %	

Carbonate and acidity



Other graphs






pH-metric high logP Assay name: Analyst: Pion 18C-01005 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-01005 Points 72 to 105

Overall results

RMSD 0.401
Average ionic strength 0.172 M
Average temperature 25.0°C
Partition ratio 0.1600 : 1

Analyte concentration range 2881.0 µM to 2978.5 µM

Total points considered 21 of 34

Warnings and errors

Errors None

Warnings One or more logP values out of range

Four-Plus parameters

à	Alpha	0.130	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r
à	S	0.9970	3/1/2018 6:16:08 AM	C:\Sirius T3\HCl18B27.t3r
à	jΗ	8.0	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r
à	jОН	-0.4	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r


Titrants


0.50 M HCI 0.993513 3/1/2018 6:16:08 AM C:\Sirius_T3\HCl18B27.t3r 0.50 M KOH 0.999845 3/1/2018 6:16:08 AM C:\Sirius_T3\KOH18B27.t3r

Sample

$\mathbf{\nabla}$	M15_octanol concentration factor	0.947
	Base pKa 1	2.58
a	Base pKa 2	5.30
	logP (XH2 2+)	-5.50
₩	logP (XH +)	-0.01
₩	logP (neutral X)	1.92

Sample graphs








Assay ID:

Sample name: M15_octanol Experiment start time: 3/1/2018 6:16:08 AM Assay name:

pH-metric high logP Analyst: Pion Instrument ID: 18C-01005 T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Sample graphs (continued)



Sample logD and percent species


рН	M15_octanol logD	M15_octanol M15_octanolH2	M15_octanol M15_octanolH	M15_octanol M15_octanol	M15_octanol M15_octanolH2*	M15_octanol M15_octanolH*	M15_octanol M15_octanol*	
1.000	-1.60	97.05 %	2.55 %	0.00 %	0.00 %	0.40 %	0.00 %	
1.200	-1.41	95.40 %	3.98 %	0.00 %	0.00 %	0.62 %	0.00 %	Stomach pH
2.000	-0.68	76.59 %	20.15 %	0.01 %	0.00 %	3.12 %	0.14 %	r
3.000	0.00	23.66 %	62.22 %	0.31 %	0.00 %	9.63 %	4.18 %	
4.000	0.68	1.99 %	52.22 %	2.62 %	0.00 %	8.08 %	35.10 %	
5.000	1.46	0.05 %	11.93 %	5.98 %	0.00 %	1.85 %	80.20 %	
6.000	1.85	0.00 %	1.36 %	6.83 %	0.00 %	0.21 %	91.60 %	
6.500	1.90	0.00 %	0.44 %	6.90 %	0.00 %	0.07 %	92.59 %	
7.000	1.91	0.00 %	0.14 %	6.93 %	0.00 %	0.02 %	92.91 %	
7.400	1.92	0.00 %	0.06 %	6.93 %	0.00 %	0.01 %	93.00 %	Blood pH
8.000	1.92	0.00 %	0.01 %	6.94 %	0.00 %	0.00 %	93.05 %	·
9.000	1.92	0.00 %	0.00 %	6.94 %	0.00 %	0.00 %	93.06 %	
10.000	1.92	0.00 %	0.00 %	6.94 %	0.00 %	0.00 %	93.06 %	
11.000	1.92	0.00 %	0.00 %	6.94 %	0.00 %	0.00 %	93.06 %	
12.000	1.92	0.00 %	0.00 %	6.94 %	0.00 %	0.00 %	93.06 %	


Carbonate and acidity



Carbonate 0.052 mM Acidity error -0.059 mM

Other graphs






Reported at: 3/2/2018 2:23:43 PM



Assay name: pH-metric high logP Analyst: Pion 18C-01005 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name	M15_octanol	2/27/2018 5:03:03 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.001970 g	2/28/2018 4:25:24 PM	User entered value
Formula weight	209.25 g/mol	2/27/2018 5:03:03 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	209.25	2/27/2018 5:03:03 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	2	2/27/2018 5:03:03 PM	User entered value
Sample is a	Base	2/27/2018 5:03:03 PM	User entered value
pKa 1	2.58	2/27/2018 5:03:03 PM	User entered value
pKa 2	5.30	2/27/2018 5:03:03 PM	User entered value
logP (XH2 2+)	-5.50	2/28/2018 2:10:35 PM	User entered value
logp (XH +)	-4.96	2/28/2018 2:10:28 PM	User entered value
logP (neutral X)	1.92	2/28/2018 2:10:15 PM	User entered value

Events

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
9:23.6	Initial pH = 7.47									
12:24.5	Data point 1		0.09563 mL				0.00623	0.50576	0.00043	10.0 s
13:10.8	Data point 2		0.09563 mL				-0.00953	0.35232	0.00079	10.0 s
13:46.7	Data point 3	2.20002 mL	0.09563 mL	0.04570 mL	0.05000 mL	2.464	-0.00246	0.19027	0.00028	10.0 s
	Data point 4	2.20002 mL	0.09563 mL	0.05576 mL	0.05000 mL	2.670	0.01234	0.53634	0.00083	10.0 s
14:58.0	Data point 5	2.20002 mL	0.09563 mL	0.06286 mL	0.05000 mL	2.900	0.00043	0.01947	0.00015	10.5 s
15:34.0	Data point 6	2.20002 mL	0.09563 mL	0.06776 mL	0.05000 mL	3.097	-0.00497	0.18772	0.00057	10.0 s
16:09.5	Data point 7	2.20002 mL	0.09563 mL	0.07131 mL	0.05000 mL	3.307	-0.00870	0.18868	0.00099	10.0 s
16:45.1	Data point 8	2.20002 mL	0.09563 mL	0.07385 mL	0.05000 mL	3.525	-0.00071	0.08489	0.00012	10.5 s
17:21.0	Data point 9		0.09563 mL				-0.00085	0.13611	0.00011	10.0 s
17:56.5	Data point 10	2.20002 mL	0.09563 mL	0.07723 mL	0.05000 mL	3.901	-0.00101	0.13288	0.00014	10.0 s
18:42.4	Data point 11	2.20002 mL	0.09563 mL	0.07869 mL	0.05000 mL	4.105	-0.00256	0.29892	0.00023	10.0 s
19:17.9	Data point 12	2.20002 mL	0.09563 mL	0.08022 mL	0.05000 mL	4.378	-0.00592	0.25762	0.00058	10.0 s
19:58.5	Data point 13	2.20002 mL	0.09563 mL	0.08154 mL	0.05000 mL	4.569	-0.00172	0.03384	0.00046	10.0 s
20:44.3	Data point 14	2.20002 mL	0.09563 mL	0.08316 mL	0.05000 mL	4.761	-0.00067	0.00810	0.00037	10.0 s
21:30.1	Data point 15	2.20002 mL	0.09563 mL	0.08511 mL	0.05000 mL	4.956	-0.00175	0.25727	0.00017	10.5 s
22:06.1	Data point 16	2.20002 mL	0.09563 mL	0.08732 mL	0.05000 mL	5.158	-0.00440	0.41157	0.00034	10.0 s
	Data point 17	2.20002 mL	0.09563 mL	0.08935 mL	0.05000 mL	5.362	0.00350	0.06911	0.00066	10.5 s
	Data point 18	2.20002 mL	0.09563 mL	0.09099 mL	0.05000 mL	5.592	-0.01026	0.39926	0.00080	10.0 s
23:53.0	Data point 19	2.20002 mL	0.09563 mL	0.09222 mL	0.05000 mL	5.831	-0.00232	0.01329	0.00099	10.0 s
24:33.6	Data point 20		0.09563 mL				-0.01229	0.83805	0.00066	10.5 s
	Data point 21		0.09563 mL				-0.01579	0.76334	0.00089	12.5 s
25:52.6	Data point 22	2.20002 mL	0.09563 mL	0.09412 mL	0.05000 mL	6.573	-0.01546	0.70398	0.00091	14.5 s
26:37.6	Data point 23	2.20002 mL	0.09563 mL	0.09450 mL	0.05000 mL	6.938	-0.01891	0.87947	0.00100	48.5 s
	Data point 24	2.20002 mL	0.09563 mL	0.09476 mL	0.05000 mL	7.316	-0.02494	0.93500	0.00127	Timed out at 59.5 s
29:37.6	Data point 25	2.20002 mL	0.09563 mL	0.09492 mL	0.05000 mL	7.818	-0.04146	0.96788	0.00208	Timed out at 59.5 s
31:13.2	Data point 26	2.20002 mL	0.09563 mL	0.09499 mL	0.05000 mL	8.108	-0.02610	0.95271	0.00132	Timed out at 59.5 s
32:48.8	Data point 27	2.20002 mL	0.09563 mL	0.09506 mL	0.05000 mL	8.299	-0.01949	0.95727	0.00098	45.5 s
	Data point 28	2.20002 mL	0.09563 mL	0.09518 mL	0.05000 mL	8.587	-0.01701	0.74040	0.00098	24.5 s
	Data point 29		0.09563 mL				-0.01958	0.96654	0.00098	19.5 s
	Data point 30		0.09563 mL				-0.01834		0.00094	
	Data point 31		0.09563 mL				-0.01816		0.00094	
	Data point 32		0.09563 mL				-0.01968		0.00099	
	Data point 33		0.09563 mL				-0.01921		0.00096	
	Data point 34		0.09563 mL				-0.00595		0.00056	
	Data point 35		0.09563 mL				-0.01774	-	0.00091	

Reported at: 3/2/2018 2:23:43 PM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Events (continued)

FACII12	(Continueu)									
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	•	dpH/dt time
41:39.6	Data point 36						-0.00850		0.00043	10.0 s
42:26.1	Data point 37						0.00527	0.50918	0.00036	
43:02.0	Data point 38						-0.01330		0.00086	
43:37.7	Data point 39						-0.00134		0.00034	
44:13.3	Data point 40						-0.00056		0.00022	
44:48.9	Data point 41						0.00129	0.22583		10.0 s
45:24.4	Data point 42						0.00052	0.05566	0.00011	10.0 s
45:59.9	Data point 43						-0.00387		0.00047	
46:35.4	Data point 44						-0.00150		0.00020	
47:11.4	Data point 45						0.00050	0.03983	0.00012	
47:46.8	Data point 46						-0.00502			
48:22.3	Data point 47						-0.00521		0.00064	
49:08.2	Data point 48						-0.00615		0.00049	
49:43.7	Data point 49						-0.00269		0.00025	
50:19.1	Data point 50						-0.01206		0.00090	
50:54.6 51:30.6	Data point 51 Data point 52						0.00999	0.53425	0.00068	
							-0.01513		0.00091	
52:11.7	Data point 53						-0.01166		0.00097	
52:48.2 53:25.6	Data point 54						-0.00111 -0.01512	0.00552	0.00073	
53.25.6 54:04.0	Data point 55 Data point 56								0.00092 0.00091	
	Data point 57						-0.01714			
54:43.4 55:36.4	Data point 58						-0.01867		0.00097	
56:28.8	Data point 59						-0.01996 -0.01534		0.00099	
50.26.6 57:33.2	Data point 60						-0.01953		0.00077 0.00099	
57.55.2 58:52.9	Data point 61						-0.01933			Timed out at
30.32.9	Data point of	2.20002 IIIL	0.20009 IIIL	0.19932 IIIL	0.13000 IIIL	1.301	-0.03001	0.97510	0.00104	59.5 s
1:00:28.6	Data point 62	2.20002 mL	0.20089 mL	0.19939 mL	0.15000 mL	7.553	-0.03803	0.97454	0.00190	Timed out at
1.00.20.0	Data point 02	2.200022	0.20000 1112	0.10000 1112	0.100001112	1.000	0.0000	0.07 10 1	0.00.00	59.5 s
1:01:59.1	Data point 63	2.20002 mL	0.20089 mL	0.19944 mL	0.15000 mL	7.763	-0.03211	0.96281	0.00162	Timed out at
										59.5 s
1:03:34.7	Data point 64	2.20002 mL	0.20089 mL	0.19951 mL	0.15000 mL	8.039	-0.02677	0.95360	0.00135	Timed out at
4.05.40.0	Data maint OF	0.000001	0.000001	0.400501	0.45000	0.000	0.04005	0.74455	0.00005	59.5 s
	Data point 65						-0.01665		0.00095	
	Data point 66						-0.01458		0.00087	
	Data point 67						-0.01709		0.00091	
	Data point 68 Data point 69						-0.01489		0.00088	
							-0.01158		0.00094 0.00088	
	Data point 70 Data point 71								0.00086	
	Data point 72						0.00067	0.00237	0.00098	
	Data point 73						0.00007	0.85670	0.00008	
	Data point 74						0.01730	0.80261	0.00093	
	Data point 75						0.01203	0.64955	0.00071	
	Data point 76						0.00730	0.38392	0.00043	
	Data point 77						0.00165	0.05640	0.00013	
	Data point 78						-0.00318		0.00014	
	Data point 79						0.00558	0.58627	0.00036	
	Data point 80						-0.00357	0.40444	0.00038	
	Data point 81						0.00000	0.00000	0.00020	
	Data point 82						-0.00271		0.00030	
	Data point 83						-0.00694		0.00073	
	Data point 84						0.00192	0.03274	0.00070	
	Data point 85						-0.00572		0.00060	
	Data point 86						-0.00051		0.00052	
	Data point 87						-0.00660		0.00089	
	Data point 88						-0.00192		0.00087	
	•									

Reported at: 3/2/2018 2:23:43 PM

Assay Events



Sample name: M15_octanol Experiment start time: 3/1/2018 6:16:08 AM

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Events (continued)

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
1:22:33.9	Data point 89	2.20002 mL	0.30889 mL	0.30562 mL	0.45000 mL	5.547	-0.00385	0.05953	0.00078	11.5 s
1:23:10.9	Data point 90	2.20002 mL	0.30889 mL	0.30600 mL	0.45000 mL	5.773	-0.01234	0.38543	0.00098	18.0 s
1:23:54.3	Data point 91	2.20002 mL	0.30889 mL	0.30626 mL	0.45000 mL	5.990	-0.01294	0.62761	0.00081	22.0 s
1:24:41.7	Data point 92	2.20002 mL			0.45000 mL				0.00085	30.5 s
1:25:42.9	Data point 93	2.20002 mL	0.30889 mL	0.30661 mL	0.45000 mL	6.538	-0.02270	0.93887	0.00116	Timed out at 59.5 s
1:27:13.4	Data point 94	2.20002 mL	0.30889 mL	0.30670 mL	0.45000 mL	6.787	-0.03247	0.95665	0.00164	Timed out at 59.5 s
1:28:38.8	Data point 95	2.20002 mL	0.30889 mL	0.30677 mL	0.45000 mL	6.955	-0.02853	0.95075	0.00145	Timed out at 59.5 s
1:30:04.2	Data point 96	2.20002 mL	0.30889 mL	0.30682 mL	0.45000 mL	6.926	-0.00738	0.28183	0.00069	10.5 s
1:30:45.2	Data point 97	2.20002 mL	0.30889 mL	0.30699 mL	0.45000 mL	8.028	-0.06462	0.99131	0.00321	Timed out at 59.5 s
1:32:20.9	Data point 98	2.20002 mL	0.30889 mL	0.30710 mL	0.45000 mL	8.410	-0.01600	0.62423	0.00100	23.5 s
1:33:25.2	Data point 99	2.20002 mL	0.30889 mL	0.30724 mL	0.45000 mL	8.685	-0.01186	0.41042	0.00092	15.5 s
1:34:16.4		2.20002 mL	0.30889 mL	0.30739 mL	0.45000 mL	8.958	-0.00728	0.18983	0.00082	14.5 s
1:35:06.5	•	2.20002 mL	0.30889 mL	0.30757 mL	0.45000 mL	9.218	-0.01378	0.67757	0.00083	13.0 s
1:35:55.3		2.20002 mL					-0.01183	0.35300	0.00098	11.5 s
1:36:37.4		2.20002 mL	0.30889 mL	0.30821 mL	0.45000 mL	9.725	-0.00923	0.21248	0.00099	11.5 s
1:37:24.7	Data point 104						-0.01830	0.91815	0.00094	11.0 s
1:38:06.3 1:38:25.8	Data point 105 Assay volumes					10.057	-0.01398	0.91248	0.00072	10.5 s



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

i liename. C.iSinus_iSilvien	tap/20100220_exp20	5_10gP_13-2\16\	2-01005_W15_0Ctano	i_pn-metric mgn logi
Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	10.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings				
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium				
ISA water volume	2.20 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.050 mL			
Partition solvent added	Automatic			
After partition addition, stir for	1 seconds			
Sample Sonication				
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	300 seconds			
After sonication stir for	5 seconds			
Sample Dissolution	.,			
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution				
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge				
Perform a carbonate purge	No			
Temperature Control				
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			
Stir speed of	50%			
Titration 1				
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	50%			
Titration 2				
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume				
Additional partition solvent added	Automatic			
After nH adjust stir for	30 seconds			

Reported at: 3/2/2018 2:23:43 PM

30 seconds

15 seconds

55%

After pH adjust stir for

Stir to allow partitioning for

Stirrer speed for partitioning



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3		•	•	•
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.300 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00100 dpH/dt			
Stability timeout after	60 seconds			

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.130	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r
Four-Plus S	0.9970	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r
Four-Plus jH	8.0	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r
Four-Plus jOH	-0.4	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r
Base concentration factor	1.000	3/1/2018 6:16:08 AM	C:\Sirius_T3\KOH18B27.t3r
Acid concentration factor	0.994	3/1/2018 6:16:08 AM	C:\Sirius_T3\HCl18B27.t3r

Instrument Settings

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module Dispenser 0 Syringe volume Firmware version	Water 2.5 mL 1.2.1(r2)	T3DM1200361	3/31/2009 5:24:52 AM 3/31/2009 5:25:05 AM
Titrant Dispenser 2 Syringe volume Firmware version	Water (0.15 M KCI) Acid 0.5 mL 1.2.1(r2)	02-06-2018	2/27/2018 10:05:59 AM 3/31/2009 5:25:11 AM
Titrant Dispenser 1 Syringe volume Firmware version	Acid (0.5 M HCI) Base 0.5 mL 1.2.1(r2)	02-27-2018	2/27/2018 10:27:22 AM 3/31/2009 5:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2)	9/22/2017	2/27/2018 10:21:22 AM 3/31/2009 5:26:24 AM
Distribution valve 5 Firmware version	Distribution Valve 1.1.3		3/31/2009 5:28:19 AM
Port A Port B Dispenser 3 Syringe volume Firmware version	Methanol (80%, 0.15 M KCI) Cyclohexane Buffer 0.5 mL	09-26-17 11-01-17	2/7/2018 9:42:01 AM 2/27/2018 10:37:57 AM 8/3/2010 5:05:16 AM
Titrant Dispenser 6	1.2.1(r2) Dodecane Octanol	2018/01/31	2/28/2018 10:18:04 AM 10/22/2010 10:52:43 AM

Reported at: 3/2/2018 2:23:43 PM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

anon annon comingo (commuca)			
Setting	Value	Batch Id	Install date
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)	04 24 2040	2/27/2019 0:50:25 AM
Titrant Titrator	Octanol	01-31-2018 T3TM1200161	2/27/2018 9:59:35 AM 3/31/2009 5:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2	1311111200101	3/3 1/2009 3.24.17 AW
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0923	1/23/2018 2:01:00 PM
E0 calibration	+3.93 mV		3/1/2018 6:16:36 AM
Filling solution	3M KCI	KCL097	2/27/2018 9:49:43 AM
Liquids			
Wash 1	50% IPA:50% Water		2/28/2018 10:23:32 AM
Wash 2	0.5% Trition X-100 in H20		2/28/2018 10:23:34 AM
Buffer position 1	pH7 Wash		2/28/2018 10:24:06 AM
Buffer position 2 Storage position	pH 7		2/28/2018 10:24:08 AM 2/28/2018 10:21:14 AM
Wash water	8.4e+003 mL	02-27-2018	2/27/2018 9:54:39 AM
Waste	7.1e+003 mL	02-27-2010	11/28/2017 10:36:29 AM
Temperature controller	7.10.000 IIIE		8/5/2010 6:35:13 AM
Turbidity detector			3/31/2009 5:24:45 AM
Spectrometer		074811	11/23/2010 11:22:28 AM
Dip probe		10196	
Wavelength coefficient A0	183.333		
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	112:08:55		11/23/2010 11:22:28 AM
Calibrated on	2/27/2018 10:40:38 AM		
Integration time	40 10		
Scans averaged Autoloader	10	T3AL1200345	11/10/2015 9:34:13 AM
Left-right axis firmware version	1.17 Al1Dl2DO2 Stepper 2	13AL12003 1 3	11/10/2013 9.54.13 AW
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Configuration	_		
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s) 1.3 mL		
Titrant tube volume Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s 30%		
E0 calibration preparation stir speed E0 calibration buffer wash stir duration	5 S		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		
	•		



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Value Batch Id Install date alibration stir duration 5 s alibration wash pump volume 20.0 mL alibration wash stir duration 5 s alibration wash stir speed 30% nse height 10000

Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00

Experiment Log

- [2:49] Air gap created for Water (0.15 M KCI)
- [2:49] Air gap created for Acid (0.5 M HCI)
- [2:50] Air gap created for Base (0.5 M KOH)
- [2:50] Air gap released for Water (0.15 M KCI)
- [2:54] Titrator arm moved over Titration position
- [2:54] Titration 1 of 3
- [2:54] Adding initial titrants
- [2:54] Automatically add 2.20000 mL of water
- [3:31] Dispensed 2.200024 mL of Water (0.15 M KCI)
- [3:35] Titrator arm moved over Drain
- [9:16] Titrator arm moved to Titration position
- [9:16] Argon flow rate set to 100
- [9:16] Stirrer speed set to 10
- [9:21] Automatically add 0.05000 mL of Octanol
- [9:23] Dispensed 0.050000 mL of Octanol
- [9:24] Initial pH = 7.47
- [9:24] Iterative adjust 7.47 -> 2.00
- [9:24] pH 7.47 -> 2.00
- [9:27] Air gap released for Acid (0.5 M HCl)
- [9:27] Dispensed 0.095626 mL of Acid (0.5 M HCI)
- [9:32] Holding pH 2.00
- [11:33] Stirrer speed set to 0
- [11:33] Stirrer speed set to 50
- [11:33] Iterative adjust 1.95 -> 2.00
- [11:33] pH 1.95 -> 2.00
- [11:33] Air gap released for Base (0.5 M KOH)
- [11:34] Dispensed 0.007785 mL of Base (0.5 M KOH)
- [12:25] Stirrer speed set to 0
- [12:35] Datapoint id 1 collected
- [12:35] Stirrer speed set to 50
- [12:40] pH 2.03 -> 2.23
- [12:40] Using cautious pH adjust
- [12:40] Dispensed 0.011783 mL of Base (0.5 M KOH)
- [12:45] Stepping pH = 2.11
- [12:46] Dispensed 0.009196 mL of Base (0.5 M KOH)
- [12:51] Stepping pH = 2.21
- [12:51] Dispensed 0.001740 mL of Base (0.5 M KOH)
- [12:56] Stepping pH = 2.23
- [13:11] Stirrer speed set to 0
- [13:21] Datapoint id 2 collected
- [13:21] Charge balance equation is out by 3.7%
- [13:21] Stirrer speed set to 50



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [13:26] pH 2.24 -> 2.44
- [13:26] Using charge balance adjust
- [13:27] Dispensed 0.015193 mL of Base (0.5 M KOH)
- [13:47] Stirrer speed set to 0
- [13:57] Datapoint id 3 collected
- [13:57] Charge balance equation is out by 9.8%
- [13:57] Stirrer speed set to 50
- [14:02] pH 2.47 -> 2.67
- [14:02] Using charge balance adjust
- [14:02] Dispensed 0.010066 mL of Base (0.5 M KOH)
- [14:22] Stirrer speed set to 0
- [14:32] Datapoint id 4 collected
- [14:32] Charge balance equation is out by 0.9%
- [14:32] Stirrer speed set to 50
- [14:38] pH 2.68 -> 2.88
- [14:38] Using charge balance adjust
- [14:38] Dispensed 0.007103 mL of Base (0.5 M KOH)
- [14:58] Stirrer speed set to 0
- [15:09] Datapoint id 5 collected
- [15:09] Charge balance equation is out by 10.7%
- [15:09] Stirrer speed set to 50
- [15:14] pH 2.91 -> 3.11
- [15:14] Using charge balance adjust
- [15:14] Dispensed 0.004892 mL of Base (0.5 M KOH)
- [15:34] Stirrer speed set to 0
- [15:44] Datapoint id 6 collected
- [15:44] Charge balance equation is out by -6.0%
- [15:44] Stirrer speed set to 50
- [15:49] pH 3.11 -> 3.31
- [15:49] Using charge balance adjust
- [15:49] Dispensed 0.003551 mL of Base (0.5 M KOH)
- [16:10] Stirrer speed set to 0
- [16:20] Datapoint id 7 collected
- [16:20] Charge balance equation is out by 0.0%
- [16:20] Stirrer speed set to 50
- [16:25] pH 3.31 -> 3.51
- [16:25] Using charge balance adjust
- [16:25] Dispensed 0.002540 mL of Base (0.5 M KOH)
- [16:45] Stirrer speed set to 0
- [16:56] Datapoint id 8 collected
- [16:56] Charge balance equation is out by 4.9%
- [16:56] Stirrer speed set to 50
- [17:01] pH 3.53 -> 3.73
- [17:01] Using charge balance adjust
- [17:01] Dispensed 0.001858 mL of Base (0.5 M KOH)
- [17:21] Stirrer speed set to 0
- [17:31] Datapoint id 9 collected
- [17:31] Charge balance equation is out by -2.1%
- [17:31] Stirrer speed set to 50
- [17:36] pH 3.74 -> 3.94
- [17:36] Using charge balance adjust
- [17:36] Dispensed 0.001529 mL of Base (0.5 M KOH)
- [17:57] Stirrer speed set to 0
- [18:07] Datapoint id 10 collected
- [18:07] Charge balance equation is out by -17.7%
- [18:07] Stirrer speed set to 50
- [18:12] pH 3.91 -> 4.11
- [18:12] Using cautious pH adjust
- [18:12] Dispensed 0.000729 mL of Base (0.5 M KOH)



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [18:17] Stepping pH = 3.99
- [18:17] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [18:22] Stepping pH = 4.09
- [18:22] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [18:27] Stepping pH = 4.11
- [18:43] Stirrer speed set to 0
- [18:53] Datapoint id 11 collected
- [18:53] Charge balance equation is out by -0.4%
- [18:53] Stirrer speed set to 50
- [18:58] pH 4.11 -> 4.31
- [18:58] Using charge balance adjust
- [18:58] Dispensed 0.001529 mL of Base (0.5 M KOH)
- [19:18] Stirrer speed set to 0
- [19:28] Datapoint id 12 collected
- [19:28] Charge balance equation is out by 33.2%
- [19:28] Stirrer speed set to 50
- [19:33] pH 4.38 -> 4.58
- [19:33] Using cautious pH adjust
- [19:33] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [19:38] Stepping pH = 4.51
- [19:38] Dispensed 0.000400 mL of Base (0.5 M KOH)
- [19:44] Stepping pH = 4.57
- [19:59] Stirrer speed set to 0
- [20:09] Datapoint id 13 collected
- [20:09] Charge balance equation is out by 28.2%
- [20:09] Stirrer speed set to 50
- [20:14] pH 4.57 -> 4.77
- [20:14] Using cautious pH adjust
- [20:14] Dispensed 0.001035 mL of Base (0.5 M KOH)
- [20:19] Stepping pH = 4.70
- [20:19] Dispensed 0.000494 mL of Base (0.5 M KOH)
- [20:24] Stepping pH = 4.76
- [20:24] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [20:29] Stepping pH = 4.77
- [20:44] Stirrer speed set to 0
- [20:54] Datapoint id 14 collected
- [20:54] Charge balance equation is out by 22.3%
- [20:54] Stirrer speed set to 50
- [21:00] pH 4.77 -> 4.97
- [21:00] Using cautious pH adjust
- [21:00] Dispensed 0.001129 mL of Base (0.5 M KOH)
- [21:05] Stepping pH = 4.88
- [21:05] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [21:10] Stepping pH = 4.95
- [21:10] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [21:15] Stepping pH = 4.96
- [21:30] Stirrer speed set to 0
- [21:41] Datapoint id 15 collected
- [21:41] Charge balance equation is out by 12.8%
- [21:41] Stirrer speed set to 50
- [21:46] pH 4.96 -> 5.16
- [21:46] Using charge balance adjust
- [21:46] Dispensed 0.002211 mL of Base (0.5 M KOH)
- [22:06] Stirrer speed set to 0
- [22:16] Datapoint id 16 collected
- [22:16] Charge balance equation is out by -2.8%
- [22:16] Stirrer speed set to 50
- [22:21] pH 5.16 -> 5.36
- [22:21] Using charge balance adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [22:22] Dispensed 0.002023 mL of Base (0.5 M KOH)
- [22:42] Stirrer speed set to 0
- [22:52] Datapoint id 17 collected
- [22:52] Charge balance equation is out by -0.1%
- [22:52] Stirrer speed set to 50
- [22:57] pH 5.37 -> 5.57
- [22:57] Using charge balance adjust
- [22:57] Dispensed 0.001646 mL of Base (0.5 M KOH)
- [23:18] Stirrer speed set to 0
- [23:28] Datapoint id 18 collected
- [23:28] Charge balance equation is out by 12.1%
- [23:28] Stirrer speed set to 50
- [23:33] pH 5.60 -> 5.80
- [23:33] Using charge balance adjust
- [23:33] Dispensed 0.001223 mL of Base (0.5 M KOH)
- [23:53] Stirrer speed set to 0
- [24:03] Datapoint id 19 collected
- [24:03] Charge balance equation is out by 16.5%
- [24:03] Stirrer speed set to 50
- [24:08] pH 5.84 -> 6.04
- [24:08] Using cautious pH adjust
- [24:08] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [24:13] Stepping pH = 5.92
- [24:14] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [24:19] Stepping pH = 6.03
- [24:34] Stirrer speed set to 0
- [24:44] Datapoint id 20 collected
- [24:44] Charge balance equation is out by 3.8%
- [24:44] Stirrer speed set to 50
- [24:49] pH 6.05 -> 6.25
- [24:49] Using charge balance adjust
- [24:49] Dispensed 0.000564 mL of Base (0.5 M KOH)
- [25:10] Stirrer speed set to 0
- [25:22] Datapoint id 21 collected
- [25:22] Charge balance equation is out by 16.9%
- [25:22] Stirrer speed set to 50
- [25:27] pH 6.29 -> 6.49
- [25:27] Using cautious pH adjust
- [25:27] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [25:32] Stepping pH = 6.33
- [25:33] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [25:38] Stepping pH = 6.52
- [25:53] Stirrer speed set to 0
- [26:07] Datapoint id 22 collected
- [26:07] Charge balance equation is out by -49.2%
- [26:07] Stirrer speed set to 50
- [26:12] pH 6.58 -> 6.78
- [26:12] Using cautious pH adjust
- [26:12] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [26:17] Stepping pH = 6.60
- [26:18] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [26:23] Stepping pH = 6.84
- [26:38] Stirrer speed set to 0
- [27:26] Datapoint id 23 collected
- [27:26] Charge balance equation is out by -83.4%
- [27:26] Stirrer speed set to 50
- [27:31] pH 6.94 -> 7.14
- [27:31] Using cautious pH adjust
- [27:32] Dispensed 0.000047 mL of Base (0.5 M KOH)



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [27:37] Stepping pH = 6.96
- [27:37] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [27:42] Stepping pH = 7.04
- [27:42] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [27:47] Stepping pH = 7.24
- [28:02] Stirrer speed set to 0
- [29:02] Datapoint id 24 collected
- [29:02] Charge balance equation is out by -171.1%
- [29:02] Stirrer speed set to 50
- [29:07] pH 7.30 -> 7.50
- [29:07] Using cautious pH adjust
- [29:07] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [29:12] Stepping pH = 7.32
- [29:12] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [29:17] Stepping pH = 7.36
- [29:17] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [29:23] Stepping pH = 7.64
- [29:38] Stirrer speed set to 0
- [30:38] Datapoint id 25 collected
- [30:38] Charge balance equation is out by -247.2%
- [30:38] Stirrer speed set to 50
- [30:43] pH 7.80 -> 8.00
- [30:43] Using cautious pH adjust
- [30:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [30:48] Stepping pH = 7.88
- [30:48] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [30:53] Stepping pH = 7.97
- [30:53] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [30:58] Stepping pH = 8.09
- [31:13] Stirrer speed set to 0
- [32:13] Datapoint id 26 collected
- [32:13] Charge balance equation is out by -210.8%
- [32:13] Stirrer speed set to 50
- [32:18] pH 8.06 -> 8.26
- [32:18] Using cautious pH adjust
- [32:19] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [32:24] Stepping pH = 8.11
- [32:24] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [32:29] Stepping pH = 8.16
- [32:29] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [32:34] Stepping pH = 8.26
- [32:49] Stirrer speed set to 0
- [33:35] Datapoint id 27 collected
- [33:35] Charge balance equation is out by -276.2%
- [33:35] Stirrer speed set to 50
- [33:40] pH 8.28 -> 8.48
- [33:40] Using cautious pH adjust
- [33:40] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [33:45] Stepping pH = 8.29
- [33:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [33:50] Stepping pH = 8.33
- [33:50] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [33:55] Stepping pH = 8.43
- [33:55] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [34:00] Stepping pH = 8.55
- [34:15] Stirrer speed set to 0
- [34:40] Datapoint id 28 collected
- [34:40] Charge balance equation is out by -522.0%
- [34:40] Stirrer speed set to 50



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [34:45] pH 8.59 -> 8.79
- [34:45] Using cautious pH adjust
- [34:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [34:50] Stepping pH = 8.63
- [34:50] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [34:55] Stepping pH = 8.69
- [34:55] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [35:00] Stepping pH = 8.75
- [35:00] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [35:05] Stepping pH = 8.81
- [35:21] Stirrer speed set to 0
- [35:40] Datapoint id 29 collected
- [35:40] Charge balance equation is out by -241.0%
- [35:40] Stirrer speed set to 50
- [35:45] pH 8.84 -> 9.04
- [35:45] Using cautious pH adjust
- [35:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [35:50] Stepping pH = 8.86
- [35:50] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [35:55] Stepping pH = 8.93
- [35:56] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [36:01] Stepping pH = 9.03
- [36:16] Stirrer speed set to 0
- [36:44] Datapoint id 30 collected
- [36:44] Charge balance equation is out by -213.1%
- [36:44] Stirrer speed set to 50
- [36:49] pH 9.07 -> 9.27
- [36:49] Using cautious pH adjust
- [36:49] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [36:54] Stepping pH = 9.11
- [36:54] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [36:59] Stepping pH = 9.21
- [36:59] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [37:04] Stepping pH = 9.27
- [37:19] Stirrer speed set to 0
- [37:34] Datapoint id 31 collected [37:34] Charge balance equation is out by -90.2%
- [37:34] Stirrer speed set to 50
- [37:39] pH 9.28 -> 9.48
- [37:39] Using cautious pH adjust
- [37:39] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [37:44] Stepping pH = 9.31
- [37:44] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [37:49] Stepping pH = 9.46
- [37:49] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [37:55] Stepping pH = 9.51
- [38:10] Stirrer speed set to 0
- [38:24] Datapoint id 32 collected
- [38:24] Charge balance equation is out by -88.1%
- [38:24] Stirrer speed set to 50
- [38:29] pH 9.52 -> 9.72
- [38:29] Using cautious pH adjust
- [38:29] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [38:34] Stepping pH = 9.57
- [38:34] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [38:39] Stepping pH = 9.70
- [38:39] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [38:44] Stepping pH = 9.73
- [38:59] Stirrer speed set to 0



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [39:12] Datapoint id 33 collected
- [39:12] Charge balance equation is out by -50.1%
- [39:12] Stirrer speed set to 50
- [39:17] pH 9.74 -> 9.94
- [39:17] Using cautious pH adjust
- [39:17] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [39:22] Stepping pH = 9.82
- [39:22] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [39:27] Stepping pH = 9.91
- [39:27] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [39:32] Stepping pH = 9.94
- [39:48] Stirrer speed set to 0 [39:58] Datapoint id 34 collected
- [39:58] Charge balance equation is out by -18.1%
- [39:58] Stirrer speed set to 50
- [40:03] pH 9.95 -> 10.05
- [40:03] Using cautious pH adjust
- [40:03] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [40:08] Stepping pH = 9.98
- [40:08] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [40:13] Stepping pH = 10.04
- [40:28] Stirrer speed set to 0
- [40:38] Datapoint id 35 collected
- [40:38] Charge balance equation is out by -23.6%
- [40:38] Titration 2 of 3
- [40:38] Adding initial titrants
- [40:38] Automatically add 0.10000 mL of Octanol
- [40:41] Dispensed 0.100000 mL of Octanol
- [40:41] Stirrer speed set to 10
- [40:42] Stirrer speed set to 55
- [40:42] Iterative adjust 10.05 -> 2.00
- [40:42] pH 10.05 -> 2.00
- [40:44] Dispensed 0.100000 mL of Acid (0.5 M HCI)
- [40:49] pH 2.03 -> 2.00
- [40:49] Dispensed 0.005268 mL of Acid (0.5 M HCI)
- [41:40] Stirrer speed set to 0
- [41:50] Datapoint id 36 collected
- [41:50] Stirrer speed set to 55
- [41:55] pH 1.96 -> 2.16
- [41:55] Using cautious pH adjust
- [41:55] Dispensed 0.014722 mL of Base (0.5 M KOH)
- [42:00] Stepping pH = 2.05
- [42:01] Dispensed 0.011054 mL of Base (0.5 M KOH)
- [42:06] Stepping pH = 2.14
- [42:06] Dispensed 0.001999 mL of Base (0.5 M KOH)
- [42:11] Stepping pH = 2.17
- [42:26] Stirrer speed set to 0
- [42:36] Datapoint id 37 collected
- [42:36] Charge balance equation is out by 5.8%
- [42:36] Stirrer speed set to 55
- [42:41] pH 2.17 -> 2.37
- [42:41] Using charge balance adjust
- [42:42] Dispensed 0.018791 mL of Base (0.5 M KOH)
- [43:02] Stirrer speed set to 0
- [43:12] Datapoint id 38 collected
- [43:12] Charge balance equation is out by 6.9%
- [43:12] Stirrer speed set to 55
- [43:17] pH 2.39 -> 2.59
- [43:17] Using charge balance adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [43:18] Dispensed 0.012371 mL of Base (0.5 M KOH)
- [43:38] Stirrer speed set to 0
- [43:48] Datapoint id 39 collected
- [43:48] Charge balance equation is out by 4.3%
- [43:48] Stirrer speed set to 55
- [43:53] pH 2.61 -> 2.81
- [43:53] Using charge balance adjust
- [43:53] Dispensed 0.008514 mL of Base (0.5 M KOH)
- [44:13] Stirrer speed set to 0
- [44:23] Datapoint id 40 collected
- [44:23] Charge balance equation is out by 6.2%
- [44:23] Stirrer speed set to 55
- [44:29] pH 2.83 -> 3.03
- [44:29] Using charge balance adjust
- [44:29] Dispensed 0.005974 mL of Base (0.5 M KOH)
- [44:49] Stirrer speed set to 0
- [44:59] Datapoint id 41 collected
- [44:59] Charge balance equation is out by 9.4%
- [44:59] Stirrer speed set to 55
- [45:04] pH 3.05 -> 3.25
- [45:04] Using charge balance adjust
- [45:04] Dispensed 0.004210 mL of Base (0.5 M KOH)
- [45:25] Stirrer speed set to 0
- [45:35] Datapoint id 42 collected
- [45:35] Charge balance equation is out by -1.8%
- [45:35] Stirrer speed set to 55
- [45:40] pH 3.25 -> 3.45
- [45:40] Using charge balance adjust
- [45:40] Dispensed 0.003151 mL of Base (0.5 M KOH)
- [46:00] Stirrer speed set to 0
- [46:10] Datapoint id 43 collected
- [46:10] Charge balance equation is out by 8.8%
- [46:10] Stirrer speed set to 55
- [46:15] pH 3.47 -> 3.67
- [46:15] Using charge balance adjust
- [46:15] Dispensed 0.002399 mL of Base (0.5 M KOH)
- [46:36] Stirrer speed set to 0
- [46:46] Datapoint id 44 collected
- [46:46] Charge balance equation is out by -5.2%
- [46:46] Stirrer speed set to 55
- [46:51] pH 3.67 -> 3.87
- [46:51] Using charge balance adjust
- [46:51] Dispensed 0.002070 mL of Base (0.5 M KOH)
- [47:11] Stirrer speed set to 0
- [47:22] Datapoint id 45 collected
- [47:22] Charge balance equation is out by -7.7%
- [47:22] Stirrer speed set to 55
- [47:27] pH 3.86 -> 4.06
- [47:27] Using charge balance adjust
- [47:27] Dispensed 0.001999 mL of Base (0.5 M KOH)
- [47:47] Stirrer speed set to 0
- [47:57] Datapoint id 46 collected
- [47:57] Charge balance equation is out by 5.7%
- [47:57] Stirrer speed set to 55
- [48:02] pH 4.07 -> 4.27
- [48:02] Using charge balance adjust
- [48:02] Dispensed 0.002117 mL of Base (0.5 M KOH)
- [48:22] Stirrer speed set to 0
- [48:32] Datapoint id 47 collected



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [48:32] Charge balance equation is out by 24.4%
- [48:32] Stirrer speed set to 55
- [48:38] pH 4.31 -> 4.51
- [48:38] Using cautious pH adjust
- [48:38] Dispensed 0.001152 mL of Base (0.5 M KOH)
- [48:43] Stepping pH = 4.41
- [48:43] Dispensed 0.000870 mL of Base (0.5 M KOH)
- [48:48] Stepping pH = 4.50
- [48:48] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [48:53] Stepping pH = 4.51
- [49:08] Stirrer speed set to 0
- [49:18] Datapoint id 48 collected
- [49:18] Charge balance equation is out by 5.4%
- [49:18] Stirrer speed set to 55
- [49:23] pH 4.51 -> 4.71
- [49:23] Using charge balance adjust
- [49:24] Dispensed 0.002352 mL of Base (0.5 M KOH)
- [49:44] Stirrer speed set to 0
- [49:54] Datapoint id 49 collected
- [49:54] Charge balance equation is out by -0.6%
- [49:54] Stirrer speed set to 55
- [49:59] pH 4.71 -> 4.91
- [49:59] Using charge balance adjust
- [49:59] Dispensed 0.002211 mL of Base (0.5 M KOH)
- [50:19] Stirrer speed set to 0
- [50:29] Datapoint id 50 collected
- [50:29] Charge balance equation is out by 5.6%
- [50:29] Stirrer speed set to 55
- [50:34] pH 4.93 -> 5.13
- [50:34] Using charge balance adjust
- [50:35] Dispensed 0.001881 mL of Base (0.5 M KOH)
- [50:55] Stirrer speed set to 0
- [51:05] Datapoint id 51 collected
- [51:05] Charge balance equation is out by 13.1%
- [51:05] Stirrer speed set to 55
- [51:10] pH 5.16 -> 5.36
- [51:10] Using charge balance adjust
- [51:11] Dispensed 0.001435 mL of Base (0.5 M KOH)
- [51:31] Stirrer speed set to 0
- [51:41] Datapoint id 52 collected
- [51:41] Charge balance equation is out by 15.8%
- [51:41] Stirrer speed set to 55
- [51:46] pH 5.39 -> 5.59
- [51:46] Using cautious pH adjust
- [51:46] Dispensed 0.000517 mL of Base (0.5 M KOH)
- [51:52] Stepping pH = 5.45
- [51:52] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [51:57] Stepping pH = 5.60
- [52:12] Stirrer speed set to 0
- [52:23] Datapoint id 53 collected
- [52:23] Charge balance equation is out by -12.6%
- [52:23] Stirrer speed set to 55
- [52:28] pH 5.65 -> 5.85
- [52:28] Using charge balance adjust
- [52:28] Dispensed 0.000659 mL of Base (0.5 M KOH)
- [52:48] Stirrer speed set to 0
- [53:00] Datapoint id 54 collected
- [53:00] Charge balance equation is out by 0.5%
- [53:00] Stirrer speed set to 55



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [53:05] pH 5.86 -> 6.06
- [53:05] Using charge balance adjust
- [53:06] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [53:26] Stirrer speed set to 0
- [53:39] Datapoint id 55 collected
- [53:39] Charge balance equation is out by -8.3%
- [53:39] Stirrer speed set to 55
- [53:44] pH 6.05 -> 6.25
- [53:44] Using charge balance adjust
- [53:44] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [54:04] Stirrer speed set to 0
- [54:18] Datapoint id 56 collected
- [54:18] Charge balance equation is out by -4.8%
- [54:18] Stirrer speed set to 55
- [54:23] pH 6.24 -> 6.44
- [54:23] Using charge balance adjust
- [54:23] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [54:44] Stirrer speed set to 0
- [55:11] Datapoint id 57 collected
- [55:11] Charge balance equation is out by -4.5%
- [55:11] Stirrer speed set to 55
- [55:16] pH 6.43 -> 6.63
- [55:16] Using charge balance adjust
- [55:16] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [55:37] Stirrer speed set to 0
- [56:04] Datapoint id 58 collected
- [56:04] Charge balance equation is out by -2.3%
- [56:04] Stirrer speed set to 55
- [56:09] pH 6.63 -> 6.83
- [56:09] Using charge balance adjust
- [56:09] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [56:29] Stirrer speed set to 0
- [57:08] Datapoint id 59 collected
- [57:08] Charge balance equation is out by -7.5%
- [57:08] Stirrer speed set to 55
- [57:13] pH 6.81 -> 7.01
- [57:13] Using charge balance adjust
- [57:13] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [57:33] Stirrer speed set to 0
- [58:22] Datapoint id 60 collected
- [58:22] Charge balance equation is out by -25.4%
- [58:22] Stirrer speed set to 55
- [58:28] pH 6.94 -> 7.14
- [58:28] Using cautious pH adjust
- [58:28] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [58:33] Stepping pH = 6.96
- [58:33] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [58:38] Stepping pH = 7.14
- [58:53] Stirrer speed set to 0
- [59:53] Datapoint id 61 collected
- [59:53] Charge balance equation is out by -82.5%
- [59:53] Stirrer speed set to 55
- [59:58] pH 7.28 -> 7.48
- [59:58] Using cautious pH adjust
- [59:58] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:00:03] Stepping pH = 7.34
- [1:00:03] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:00:08] Stepping pH = 7.41
- [1:00:09] Dispensed 0.000024 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [1:00:14] Stepping pH = 7.51
- [1:00:29] Stirrer speed set to 0
- [1:01:29] Datapoint id 62 collected
- [1:01:29] Charge balance equation is out by -78.6%
- [1:01:29] Stirrer speed set to 55
- [1:01:34] pH 7.53 -> 7.73
- [1:01:34] Using cautious pH adjust
- [1:01:34] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:01:39] Stepping pH = 7.61
- [1:01:39] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:01:44] Stepping pH = 7.72
- [1:01:59] Stirrer speed set to 0
- 1:02:59 Datapoint id 63 collected
- [1:02:59] Charge balance equation is out by -84.1%
- [1:02:59] Stirrer speed set to 55
- [1:03:04] pH 7.73 -> 7.93
- [1:03:04] Using cautious pH adjust
- [1:03:04] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:03:09] Stepping pH = 7.79
- [1:03:10] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:03:15] Stepping pH = 7.89
- [1:03:15] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:03:20] Stepping pH = 7.99
- [1:03:35] Stirrer speed set to 0
- [1:04:35] Datapoint id 64 collected
- [1:04:35] Charge balance equation is out by -269.2%
- [1:04:35] Stirrer speed set to 55
- [1:04:40] pH 8.03 -> 8.23
- [1:04:40] Using cautious pH adjust
- [1:04:40] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:04:45] Stepping pH = 8.10
- [1:04:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:04:50] Stepping pH = 8.18
- [1:04:50] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:04:55] Stepping pH = 8.27
- [1:05:10] Stirrer speed set to 0
- [1:05:35] Datapoint id 65 collected
- [1:05:35] Charge balance equation is out by -326.7%
- [1:05:35] Stirrer speed set to 55
- [1:05:40] pH 8.29 -> 8.49
- [1:05:40] Using cautious pH adjust
- [1:05:40] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:05:45] Stepping pH = 8.32
- [1:05:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:05:50] Stepping pH = 8.36
- [1:05:51] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:05:56] Stepping pH = 8.37
- [1:05:56] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:06:01] Stepping pH = 8.53
- [1:06:16] Stirrer speed set to 0
- [1:06:45] Datapoint id 66 collected
- [1:06:45] Charge balance equation is out by -901.8%
- [1:06:45] Stirrer speed set to 55
- [1:06:50] pH 8.59 -> 8.79
- [1:06:50] Using cautious pH adjust
- [1:06:50] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:06:55] Stepping pH = 8.62
- [1:06:55] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:07:00] Stepping pH = 8.66



Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [1:07:00] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:07:05] Stepping pH = 8.73
- [1:07:05] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:07:11] Stepping pH = 8.78
- [1:07:26] Stirrer speed set to 0
- [1:07:42] Datapoint id 67 collected
- [1:07:42] Charge balance equation is out by -378.7%
- [1:07:42] Stirrer speed set to 55
- [1:07:47] pH 8.81 -> 9.01
- [1:07:47] Using cautious pH adjust
- [1:07:47] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:07:52] Stepping pH = 8.81
- [1:07:53] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [1:07:58] Stepping pH = 8.81
- [1:07:58] Dispensed 0.000564 mL of Base (0.5 M KOH)
- [1:08:03] Stepping pH = 9.06
- [1:08:18] Stirrer speed set to 0
- [1:08:32] Datapoint id 68 collected
- [1:08:32] Charge balance equation is out by -1,458.4%
- [1:08:32] Stirrer speed set to 55
- [1:08:37] pH 9.20 -> 9.40
- [1:08:37] Using cautious pH adjust
- [1:08:37] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:08:42] Stepping pH = 9.21
- [1:08:42] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:08:47] Stepping pH = 9.28
- [1:08:47] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [1:08:53] Stepping pH = 9.44
- [1:09:08] Stirrer speed set to 0
- [1:09:26] Datapoint id 69 collected
- [1:09:26] Charge balance equation is out by -266.5%
- [1:09:26] Stirrer speed set to 55
- [1:09:31] pH 9.51 -> 9.71
- [1:09:31] Using cautious pH adjust
- [1:09:31] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:09:36] Stepping pH = 9.51
- [1:09:36] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [1:09:41] Stepping pH = 9.70
- [1:09:56] Stirrer speed set to 0
- [1:10:09] Datapoint id 70 collected
- [1:10:09] Charge balance equation is out by -98.1%
- [1:10:09] Stirrer speed set to 55
- [1:10:14] pH 9.79 -> 9.99
- [1:10:14] Using cautious pH adjust
- [1:10:14] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [1:10:19] Stepping pH = 9.80
- [1:10:19] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [1:10:24] Stepping pH = 10.04
- [1:10:39] Stirrer speed set to 0
- [1:10:49] Datapoint id 71 collected
- [1:10:49] Charge balance equation is out by -90.6%
- [1:10:49] Titration 3 of 3
- [1:10:49] Adding initial titrants
- [1:10:49] Automatically add 0.30000 mL of Octanol
- [1:10:56] Dispensed 0.300000 mL of Octanol
- [1:10:56] Stirrer speed set to 10
- [1:10:58] Stirrer speed set to 60
- [1:10:58] Iterative adjust 10.07 -> 2.00
- [1:10:58] pH 10.07 -> 2.00



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [1:11:00] Dispensed 0.100000 mL of Acid (0.5 M HCI)
- [1:11:05] pH 2.05 -> 2.00
- [1:11:05] Dispensed 0.007996 mL of Acid (0.5 M HCI)
- [1:11:56] Stirrer speed set to 0
- [1:12:06] Datapoint id 72 collected
- [1:12:06] Stirrer speed set to 60
- [1:12:11] pH 1.98 -> 2.18
- [1:12:11] Using cautious pH adjust
- [1:12:11] Dispensed 0.015334 mL of Base (0.5 M KOH)
- [1:12:16] Stepping pH = 2.06
- [1:12:17] Dispensed 0.011900 mL of Base (0.5 M KOH)
- [1:12:22] Stepping pH = 2.16
- [1:12:22] Dispensed 0.001834 mL of Base (0.5 M KOH)
- [1:12:27] Stepping pH = 2.18
- [1:12:42] Stirrer speed set to 0
- [1:12:52] Datapoint id 73 collected
- [1:12:52] Charge balance equation is out by 5.2%
- [1:12:52] Stirrer speed set to 60
- [1:12:57] pH 2.19 -> 2.39
- [1:12:57] Using charge balance adjust
- [1:12:58] Dispensed 0.019450 mL of Base (0.5 M KOH)
- [1:13:18] Stirrer speed set to 0
- [1:13:28] Datapoint id 74 collected
- [1:13:28] Charge balance equation is out by 10.4%
- [1:13:28] Stirrer speed set to 60
- [1:13:33] pH 2.42 -> 2.62
- [1:13:33] Using charge balance adjust
- [1:13:34] Dispensed 0.012676 mL of Base (0.5 M KOH)
- [1:13:54] Stirrer speed set to 0
- [1:14:04] Datapoint id 75 collected
- [1:14:04] Charge balance equation is out by 3.9%
- [1:14:04] Stirrer speed set to 60
- [1:14:09] pH 2.63 -> 2.83
- [1:14:09] Using charge balance adjust
- [1:14:09] Dispensed 0.008843 mL of Base (0.5 M KOH)
- [1:14:29] Stirrer speed set to 0
- [1:14:39] Datapoint id 76 collected
- [1:14:39] Charge balance equation is out by 9.8%
- [1:14:39] Stirrer speed set to 60
- [1:14:45] pH 2.86 -> 3.06
- [1:14:45] Using charge balance adjust
- [1:14:45] Dispensed 0.006256 mL of Base (0.5 M KOH)
- [1:15:05] Stirrer speed set to 0
- [1:15:15] Datapoint id 77 collected
- [1:15:15] Charge balance equation is out by -0.7%
- [1:15:15] Stirrer speed set to 60
- [1:15:20] pH 3.06 -> 3.26
- [1:15:20] Using charge balance adjust
- [1:15:20] Dispensed 0.004704 mL of Base (0.5 M KOH)
- [1:15:41] Stirrer speed set to 0
- [1:15:51] Datapoint id 78 collected
- [1:15:51] Charge balance equation is out by 11.5%
- [1:15:51] Stirrer speed set to 60
- [1:15:56] pH 3.28 -> 3.48
- [1:15:56] Using charge balance adjust
- [1:15:56] Dispensed 0.003622 mL of Base (0.5 M KOH)
- [1:16:16] Stirrer speed set to 0
- [1:16:26] Datapoint id 79 collected
- [1:16:26] Charge balance equation is out by 4.3%
- Reported at: 3/2/2018 2:23:43 PM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [1:16:26] Stirrer speed set to 60
- [1:16:31] pH 3.50 -> 3.70
- [1:16:31] Using charge balance adjust
- [1:16:31] Dispensed 0.003057 mL of Base (0.5 M KOH)
- [1:16:52] Stirrer speed set to 0
- [1:17:02] Datapoint id 80 collected
- [1:17:02] Charge balance equation is out by 17.0%
- [1:17:02] Stirrer speed set to 60
- [1:17:07] pH 3.73 -> 3.93
- [1:17:07] Using cautious pH adjust
- [1:17:07] Dispensed 0.001411 mL of Base (0.5 M KOH)
- [1:17:12] Stepping pH = 3.83
- [1:17:12] Dispensed 0.001011 mL of Base (0.5 M KOH)
- [1:17:17] Stepping pH = 3.91
- [1:17:17] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:17:22] Stepping pH = 3.93
- [1:17:37] Stirrer speed set to 0
- [1:17:48] Datapoint id 81 collected
- [1:17:48] Charge balance equation is out by 7.0%
- [1:17:48] Stirrer speed set to 60
- [1:17:53] pH 3.92 -> 4.12
- [1:17:53] Using charge balance adjust
- [1:17:53] Dispensed 0.002752 mL of Base (0.5 M KOH)
- [1:18:13] Stirrer speed set to 0
- [1:18:23] Datapoint id 82 collected
- [1:18:23] Charge balance equation is out by 8.4%
- [1:18:23] Stirrer speed set to 60
- [1:18:29] pH 4.14 -> 4.34
- [1:18:29] Using charge balance adjust
- [1:18:29] Dispensed 0.002611 mL of Base (0.5 M KOH)
- [1:18:49] Stirrer speed set to 0
- [1:18:59] Datapoint id 83 collected
- [1:18:59] Charge balance equation is out by 20.8%
- [1:18:59] Stirrer speed set to 60
- [1:19:04] pH 4.38 -> 4.58
- [1:19:04] Using cautious pH adjust
- [1:19:04] Dispensed 0.001152 mL of Base (0.5 M KOH)
- [1:19:09] Stepping pH = 4.47
- [1:19:09] Dispensed 0.000847 mL of Base (0.5 M KOH)
- [1:19:15] Stepping pH = 4.56
- [1:19:15] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:19:20] Stepping pH = 4.58
- [1:19:35] Stirrer speed set to 0
- [1:19:45] Datapoint id 84 collected
- [1:19:45] Charge balance equation is out by 7.2%
- [1:19:45] Stirrer speed set to 60
- [1:19:50] pH 4.57 -> 4.77
- [1:19:50] Using charge balance adjust
- [1:19:50] Dispensed 0.001929 mL of Base (0.5 M KOH)
- [1:20:10] Stirrer speed set to 0
- [1:20:20] Datapoint id 85 collected
- [1:20:20] Charge balance equation is out by -8.5%
- [1:20:20] Stirrer speed set to 60
- [1:20:25] pH 4.75 -> 4.95
- [1:20:25] Using charge balance adjust
- [1:20:26] Dispensed 0.001505 mL of Base (0.5 M KOH)
- [1:20:46] Stirrer speed set to 0
- [1:20:56] Datapoint id 86 collected
- [1:20:56] Charge balance equation is out by -4.5%



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [1:20:56] Stirrer speed set to 60
- [1:21:01] pH 4.95 -> 5.15
- [1:21:01] Using charge balance adjust
- [1:21:01] Dispensed 0.001129 mL of Base (0.5 M KOH)
- [1:21:22] Stirrer speed set to 0
- [1:21:32] Datapoint id 87 collected
- [1:21:32] Charge balance equation is out by -5.7%
- [1:21:32] Stirrer speed set to 60
- [1:21:37] pH 5.14 -> 5.34
- [1:21:37] Using charge balance adjust
- [1:21:37] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [1:21:58] Stirrer speed set to 0
- [1:22:09] Datapoint id 88 collected
- [1:22:09] Charge balance equation is out by -0.8%
- [1:22:09] Stirrer speed set to 60
- [1:22:14] pH 5.34 -> 5.54
- [1:22:14] Using charge balance adjust
- [1:22:14] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [1:22:34] Stirrer speed set to 0
- [1:22:46] Datapoint id 89 collected
- [1:22:46] Charge balance equation is out by 1.2%
- [1:22:46] Stirrer speed set to 60
- [1:22:51] pH 5.55 -> 5.75
- [1:22:51] Using charge balance adjust
- [1:22:51] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [1:23:11] Stirrer speed set to 0
- [1:23:29] Datapoint id 90 collected
- [1:23:29] Charge balance equation is out by 9.5%
- [1:23:29] Stirrer speed set to 60
- [1:23:34] pH 5.77 -> 5.97
- [1:23:34] Using charge balance adjust
- [1:23:34] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [1:23:54] Stirrer speed set to 0
- [1:24:16] Datapoint id 91 collected
- [1:24:16] Charge balance equation is out by 7.6%
- [1:24:16] Stirrer speed set to 60
- [1:24:22] pH 6.00 -> 6.20
- [1:24:22] Using charge balance adjust
- [1:24:22] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:24:42] Stirrer speed set to 0
- [1:25:12] Datapoint id 92 collected
- [1:25:12] Charge balance equation is out by 18.0%
- [1:25:12] Stirrer speed set to 60
- [1:25:18] pH 6.23 -> 6.43
- [1:25:18] Using cautious pH adjust
- [1:25:18] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:25:23] Stepping pH = 6.29
- [1:25:23] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:25:28] Stepping pH = 6.50
- [1:25:43] Stirrer speed set to 0
- [1:26:43] Datapoint id 93 collected
- [1:26:43] Charge balance equation is out by -15.4%
- [1:26:43] Stirrer speed set to 60
- [1:26:48] pH 6.54 -> 6.74
- [1:26:48] Using cautious pH adjust
- [1:26:48] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:26:53] Stepping pH = 6.63
- [1:26:53] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:26:58] Stepping pH = 6.80



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [1:27:14] Stirrer speed set to 0
- [1:28:14] Datapoint id 94 collected
- [1:28:14] Charge balance equation is out by 12.0%
- [1:28:14] Stirrer speed set to 60
- [1:28:19] pH 6.78 -> 6.98
- [1:28:19] Using charge balance adjust
- [1:28:19] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:28:39] Stirrer speed set to 0
- [1:29:39] Datapoint id 95 collected
- [1:29:39] Charge balance equation is out by -11.1%
- [1:29:39] Stirrer speed set to 60
- [1:29:44] pH 6.97 -> 7.17
- [1:29:44] Using charge balance adjust
- [1:29:44] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:30:04] Stirrer speed set to 0
- [1:30:15] Datapoint id 96 collected
- [1:30:15] Charge balance equation is out by -76.1%
- [1:30:15] Stirrer speed set to 60
- [1:30:20] pH 6.91 -> 7.11
- [1:30:20] Using cautious pH adjust
- [1:30:20] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:30:25] Stepping pH = 6.90
- [1:30:25] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:30:30] Stepping pH = 7.91
- [1:30:45] Stirrer speed set to 0
- [1:31:45] Datapoint id 97 collected
- [1:31:45] Charge balance equation is out by -214.0%
- [1:31:45] Stirrer speed set to 60
- [1:31:50] pH 8.09 -> 8.29
- [1:31:50] Using cautious pH adjust
- [1:31:51] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:31:56] Stepping pH = 8.14
- [1:31:56] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:32:01] Stepping pH = 8.14
- [1:32:01] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:32:06] Stepping pH = 8.36
- [1:32:21] Stirrer speed set to 0
- [1:32:45] Datapoint id 98 collected
- [1:32:45] Charge balance equation is out by -580.1%
- [1:32:45] Stirrer speed set to 60
- [1:32:50] pH 8.43 -> 8.63
- [1:32:50] Using cautious pH adjust
- [1:32:50] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:32:55] Stepping pH = 8.45
- [1:32:55] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:33:00] Stepping pH = 8.48
- [1:33:00] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:33:05] Stepping pH = 8.61
- [1:33:05] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:33:10] Stepping pH = 8.68
- [1:33:25] Stirrer speed set to 0
- [1:33:41] Datapoint id 99 collected
- [1:33:41] Charge balance equation is out by -464.7%
- [1:33:41] Stirrer speed set to 60
- [1:33:46] pH 8.69 -> 8.89
- [1:33:46] Using cautious pH adjust
- [1:33:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:33:51] Stepping pH = 8.71
- [1:33:51] Dispensed 0.000047 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

- [1:33:56] Stepping pH = 8.77
- [1:33:56] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:34:01] Stepping pH = 8.90
- [1:34:17] Stirrer speed set to 0
- [1:34:31] Datapoint id 100 collected
- [1:34:31] Charge balance equation is out by -252.1%
- [1:34:31] Stirrer speed set to 60
- [1:34:36] pH 8.96 -> 9.16
- [1:34:36] Using cautious pH adjust
- [1:34:36] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:34:41] Stepping pH = 8.98
- [1:34:41] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:34:46] Stepping pH = 9.09
- [1:34:46] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:34:52] Stepping pH = 9.18
- [1:35:07] Stirrer speed set to 0
- [1:35:20] Datapoint id 101 collected
- [1:35:20] Charge balance equation is out by -179.3%
- [1:35:20] Stirrer speed set to 60
- [1:35:25] pH 9.23 -> 9.43
- [1:35:25] Using cautious pH adjust
- [1:35:25] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:35:30] Stepping pH = 9.24
- [1:35:30] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [1:35:35] Stepping pH = 9.36
- [1:35:35] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:35:40] Stepping pH = 9.43
- [1:35:55] Stirrer speed set to 0
- [1:36:07] Datapoint id 102 collected
- [1:36:07] Charge balance equation is out by -161.3%
- [1:36:07] Stirrer speed set to 60
- [1:36:12] pH 9.46 -> 9.66
- [1:36:12] Using cautious pH adjust
- [1:36:12] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:36:17] Stepping pH = 9.47
- [1:36:17] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [1:36:22] Stepping pH = 9.67
- [1:36:38] Stirrer speed set to 0
- [1:36:49] Datapoint id 103 collected
- [1:36:49] Charge balance equation is out by -94.3%
- [1:36:49] Stirrer speed set to 60
- [1:36:54] pH 9.73 -> 9.93
- [1:36:54] Using cautious pH adjust
- [1:36:54] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:36:59] Stepping pH = 9.76
- [1:36:59] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [1:37:05] Stepping pH = 9.91
- [1:37:05] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:37:10] Stepping pH = 9.95
- [1:37:25] Stirrer speed set to 0
- [1:37:36] Datapoint id 104 collected
- [1:37:36] Charge balance equation is out by -63.7%
- [1:37:36] Stirrer speed set to 60
- [1:37:41] pH 9.97 -> 10.05
- [1:37:41] Using cautious pH adjust
- [1:37:41] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:37:46] Stepping pH = 9.98
- [1:37:46] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [1:37:51] Stepping pH = 10.04





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01005 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01005_M15_octanol_pH-metric high logP.t3r

Experiment Log (continued)

[1:38:06] Stirrer speed set to 0

[1:38:17] Datapoint id 105 collected

[1:38:17] Charge balance equation is out by -43.3%

[1:38:17] Argon flow rate set to 0

[1:38:21] Titrator arm moved over Titration position

Reported at: 3/2/2018 2:23:43 PM