

Sample name: **M07**  
 Assay name: **UV-metric pKa**  
 Assay ID: **17I-16002**  
 Filename: **C:\Sirius\_T3\17I-16002\_M07\_UV-metric pKa.t3r**

Experiment start time: **9/15/2017 11:47:51 PM**  
 Analyst: **Dorothy Levorse**  
 Instrument ID: **T311053**

## Results

pKa 1 **6.07**  
 RMSD **0.002 0.002**  
 Chi squared **0.0116**  
 PCA calculated number of pKas **3**  
 Average ionic strength **0.158 M**  
 Average temperature **24.9°C**  
 Analyte concentration range **100.6 µM to 91.2 µM**

Number of pKas source **Predicted**  
 Wavelength clipping **230.0 nm to 450.0 nm**  
 pH clipping **1.276 to 12.701**


## Warnings and errors

Errors None  
 Warnings PCA calculation disagrees with predicted number of pKas

## Assay Settings

| Setting                     | Value            | Original Value | Date/Time changed | Imported from |
|-----------------------------|------------------|----------------|-------------------|---------------|
| Buffer in use               | Yes              |                |                   |               |
| Buffer type                 | Phosphate Buffer |                |                   |               |
| <b>Assay Medium</b>         |                  |                |                   |               |
| Volume of buffer introduced | 0.025000 mL      |                |                   |               |
| Add buffer manually         | Manual           |                |                   |               |









## Graphs



Sample name: **M07**  
 Assay name: **UV-metric pKa**  
 Assay ID: **17I-16002**  
 Filename: **C:\Sirius\_T3\17I-16002\_M07\_UV-metric pKa.t3r**

Experiment start time: **9/15/2017 11:47:51 PM**  
 Analyst: **Dorothy Leverse**  
 Instrument ID: **T311053**


## Graphs (continued)



Sample name: **M07**  
 Assay name: **UV-metric pKa**  
 Assay ID: **171-16002**  
 Filename: **C:\Sirius\_T3\171-16002\_M07\_UV-metric pKa.t3r**

Experiment start time: **9/15/2017 11:47:51 PM**  
 Analyst: **Dorothy Levorse**  
 Instrument ID: **T311053**

## Graphs (continued)



## Assay Model

| Settings                          | Value      | Date/Time changed    | Imported from      |
|-----------------------------------|------------|----------------------|--------------------|
| Sample name                       | M07        | 9/15/2017 3:35:09 PM | User entered value |
| Sample by                         | Volume     |                      | Default value      |
| Sample volume                     | 0.0030 mL  | 9/15/2017 3:35:09 PM | User entered value |
| Solvent                           | DMSO       |                      | Default value      |
| Sample concentration              | 0.053600 M | 9/15/2017 3:35:09 PM | User entered value |
| Solubility                        | Unknown    |                      | Default value      |
| Molecular weight                  | 235.28     | 9/15/2017 3:35:16 PM | User entered value |
| Individual pKa ionic environments | No         |                      | Default value      |
| Number of pKas                    | 1          | 9/15/2017 3:35:09 PM | User entered value |
| Sample is a                       | Base       | 9/15/2017 3:35:09 PM | User entered value |
| pKa 1                             | 5.70       | 9/15/2017 3:35:09 PM | User entered value |
| logp (XH +)                       | -10.00     |                      | Default value      |
| logP (neutral X)                  | -10.00     | 9/15/2017 3:35:09 PM | User entered value |

## Events

| Time   | Event                           | Water      | Acid       | Base       | Buffer     | pH    | dpH/dt   | pH R-squared | pH SD   |
|--------|---------------------------------|------------|------------|------------|------------|-------|----------|--------------|---------|
| 3:12.3 | Dark spectrum                   |            |            |            |            |       |          |              |         |
| 3:13.8 | Reference spectrum              |            |            |            |            |       |          |              |         |
| 3:41.4 | Volume reset due to vial change |            |            |            |            |       |          |              |         |
| 5:11.8 | Initial pH = 7.61               |            |            |            |            |       |          |              |         |
| 6:24.8 | Data point 4                    | 1.50000 mL | 0.07084 mL | 0.00000 mL | 0.02500 mL | 1.776 | -0.00744 | 0.73426      | 0.00043 |
| 6:53.6 | Data point 5                    | 1.50000 mL | 0.07084 mL | 0.02547 mL | 0.02500 mL | 1.978 | 0.00458  | 0.17102      | 0.00055 |
| 7:10.8 | Data point 6                    | 1.50000 mL | 0.07084 mL | 0.04219 mL | 0.02500 mL | 2.184 | 0.00516  | 0.10363      | 0.00079 |
| 7:27.7 | Data point 7                    | 1.50000 mL | 0.07084 mL | 0.05252 mL | 0.02500 mL | 2.390 | -0.01054 | 0.48727      | 0.00075 |
| 7:44.5 | Data point 8                    | 1.50000 mL | 0.07084 mL | 0.05894 mL | 0.02500 mL | 2.587 | 0.00719  | 0.47472      | 0.00051 |
| 8:01.3 | Data point 9                    | 1.50000 mL | 0.07084 mL | 0.06301 mL | 0.02500 mL | 2.826 | 0.00143  | 0.07186      | 0.00026 |

Sample name: **M07**  
 Assay name: **UV-metric pKa**  
 Assay ID: **17I-16002**  
 Filename: **C:\Sirius\_T3\17I-16002\_M07\_UV-metric pKa.t3r**

Experiment start time: **9/15/2017 11:47:51 PM**  
 Analyst: **Dorothy Levorse**  
 Instrument ID: **T311053**

## Events (continued)

| Time    | Event         | Water      | Acid       | Base       | Buffer     | pH     | dpH/dt   | pH R-squared | pH SD   | dpH/dt time |
|---------|---------------|------------|------------|------------|------------|--------|----------|--------------|---------|-------------|
| 8:33.4  | Data point 10 | 1.50000 mL | 0.07084 mL | 0.06540 mL | 0.02500 mL | 3.020  | 0.00892  | 0.78278      | 0.00050 | 10.0 s      |
| 8:50.0  | Data point 11 | 1.50000 mL | 0.07084 mL | 0.06691 mL | 0.02500 mL | 3.211  | 0.01230  | 0.89789      | 0.00064 | 10.0 s      |
| 9:06.6  | Data point 12 | 1.50000 mL | 0.07084 mL | 0.06787 mL | 0.02500 mL | 3.313  | 0.00426  | 0.47301      | 0.00031 | 10.0 s      |
| 9:38.6  | Data point 13 | 1.50000 mL | 0.07084 mL | 0.06881 mL | 0.02500 mL | 3.522  | 0.01169  | 0.87789      | 0.00062 | 10.0 s      |
| 10:05.2 | Data point 14 | 1.50000 mL | 0.07084 mL | 0.06931 mL | 0.02500 mL | 3.741  | 0.02154  | 0.93404      | 0.00110 | 10.0 s      |
| 10:21.9 | Data point 15 | 1.50000 mL | 0.07084 mL | 0.06959 mL | 0.02500 mL | 4.021  | 0.03685  | 0.89652      | 0.00192 | 10.0 s      |
| 10:43.4 | Data point 16 | 1.50000 mL | 0.07084 mL | 0.06987 mL | 0.02500 mL | 4.551  | 0.06049  | 0.91743      | 0.00312 | 10.0 s      |
| 11:10.2 | Data point 17 | 1.50000 mL | 0.07084 mL | 0.07001 mL | 0.02500 mL | 4.788  | 0.09184  | 0.95650      | 0.00469 | 10.5 s      |
| 11:32.4 | Data point 18 | 1.50000 mL | 0.07084 mL | 0.07011 mL | 0.02500 mL | 5.072  | 0.09633  | 0.97713      | 0.00481 | 13.5 s      |
| 11:57.6 | Data point 19 | 1.50000 mL | 0.07084 mL | 0.07016 mL | 0.02500 mL | 5.295  | 0.09549  | 0.97449      | 0.00478 | 12.5 s      |
| 12:21.8 | Data point 20 | 1.50000 mL | 0.07084 mL | 0.07023 mL | 0.02500 mL | 5.565  | 0.04991  | 0.96241      | 0.00251 | 10.0 s      |
| 12:43.5 | Data point 21 | 1.50000 mL | 0.07084 mL | 0.07027 mL | 0.02500 mL | 5.775  | 0.02473  | 0.56978      | 0.00164 | 10.0 s      |
| 13:00.1 | Data point 22 | 1.50000 mL | 0.07084 mL | 0.07034 mL | 0.02500 mL | 6.030  | -0.00122 | 0.00203      | 0.00134 | 10.0 s      |
| 13:21.6 | Data point 23 | 1.50000 mL | 0.07084 mL | 0.07041 mL | 0.02500 mL | 6.247  | 0.02532  | 0.48669      | 0.00180 | 10.0 s      |
| 13:43.3 | Data point 24 | 1.50000 mL | 0.07084 mL | 0.07051 mL | 0.02500 mL | 6.499  | 0.01222  | 0.17135      | 0.00146 | 10.0 s      |
| 14:10.2 | Data point 25 | 1.50000 mL | 0.07084 mL | 0.07060 mL | 0.02500 mL | 6.707  | 0.01306  | 0.27836      | 0.00124 | 10.0 s      |
| 14:42.3 | Data point 26 | 1.50000 mL | 0.07084 mL | 0.07072 mL | 0.02500 mL | 6.953  | 0.05092  | 0.71439      | 0.00300 | 10.0 s      |
| 15:09.1 | Data point 27 | 1.50000 mL | 0.07084 mL | 0.07081 mL | 0.02500 mL | 7.161  | 0.08383  | 0.84268      | 0.00451 | 10.0 s      |
| 15:35.7 | Data point 28 | 1.50000 mL | 0.07084 mL | 0.07088 mL | 0.02500 mL | 7.456  | 0.07496  | 0.77394      | 0.00420 | 11.5 s      |
| 16:04.0 | Data point 29 | 1.50000 mL | 0.07084 mL | 0.07095 mL | 0.02500 mL | 7.826  | 0.08509  | 0.82167      | 0.00469 | 13.5 s      |
| 16:34.5 | Data point 30 | 1.50000 mL | 0.07084 mL | 0.07103 mL | 0.02500 mL | 8.239  | 0.07265  | 0.70504      | 0.00432 | 13.5 s      |
| 17:04.8 | Data point 31 | 1.50000 mL | 0.07084 mL | 0.07110 mL | 0.02500 mL | 8.591  | 0.06777  | 0.70972      | 0.00397 | 11.5 s      |
| 17:33.1 | Data point 32 | 1.50000 mL | 0.07084 mL | 0.07117 mL | 0.02500 mL | 8.835  | 0.07881  | 0.81206      | 0.00432 | 11.0 s      |
| 18:06.1 | Data point 33 | 1.50000 mL | 0.07084 mL | 0.07126 mL | 0.02500 mL | 9.058  | 0.06474  | 0.87630      | 0.00342 | 10.0 s      |
| 18:33.0 | Data point 34 | 1.50000 mL | 0.07084 mL | 0.07138 mL | 0.02500 mL | 9.300  | 0.02884  | 0.79122      | 0.00160 | 10.0 s      |
| 19:05.0 | Data point 35 | 1.50000 mL | 0.07084 mL | 0.07152 mL | 0.02500 mL | 9.525  | 0.01333  | 0.62782      | 0.00083 | 10.0 s      |
| 19:26.7 | Data point 36 | 1.50000 mL | 0.07084 mL | 0.07168 mL | 0.02500 mL | 9.725  | 0.00622  | 0.17959      | 0.00073 | 10.0 s      |
| 19:58.5 | Data point 37 | 1.50000 mL | 0.07084 mL | 0.07194 mL | 0.02500 mL | 9.931  | -0.00266 | 0.11481      | 0.00039 | 10.0 s      |
| 20:25.4 | Data point 38 | 1.50000 mL | 0.07084 mL | 0.07222 mL | 0.02500 mL | 10.126 | -0.00793 | 0.58188      | 0.00051 | 10.0 s      |
| 20:42.0 | Data point 39 | 1.50000 mL | 0.07084 mL | 0.07260 mL | 0.02500 mL | 10.316 | -0.01107 | 0.92163      | 0.00057 | 10.0 s      |
| 20:58.7 | Data point 40 | 1.50000 mL | 0.07084 mL | 0.07314 mL | 0.02500 mL | 10.484 | -0.01598 | 0.94608      | 0.00081 | 10.0 s      |
| 21:30.8 | Data point 41 | 1.50000 mL | 0.07084 mL | 0.07422 mL | 0.02500 mL | 10.676 | -0.01260 | 0.88860      | 0.00066 | 10.0 s      |
| 22:07.8 | Data point 42 | 1.50000 mL | 0.07084 mL | 0.07594 mL | 0.02500 mL | 10.871 | -0.01405 | 0.89635      | 0.00074 | 10.0 s      |
| 22:34.9 | Data point 43 | 1.50000 mL | 0.07084 mL | 0.07829 mL | 0.02500 mL | 11.061 | -0.01459 | 0.92807      | 0.00075 | 10.0 s      |
| 23:01.7 | Data point 44 | 1.50000 mL | 0.07084 mL | 0.08116 mL | 0.02500 mL | 11.251 | -0.01209 | 0.89220      | 0.00063 | 10.0 s      |
| 23:18.5 | Data point 45 | 1.50000 mL | 0.07084 mL | 0.08570 mL | 0.02500 mL | 11.417 | -0.00955 | 0.91159      | 0.00049 | 10.0 s      |
| 23:50.9 | Data point 46 | 1.50000 mL | 0.07084 mL | 0.09393 mL | 0.02500 mL | 11.614 | -0.00899 | 0.93119      | 0.00046 | 10.0 s      |
| 24:23.4 | Data point 47 | 1.50000 mL | 0.07084 mL | 0.10745 mL | 0.02500 mL | 11.811 | -0.00933 | 0.84361      | 0.00050 | 10.0 s      |
| 24:50.9 | Data point 48 | 1.50000 mL | 0.07084 mL | 0.12801 mL | 0.02500 mL | 12.002 | -0.00561 | 0.69123      | 0.00033 | 10.0 s      |
| 25:24.0 | Data point 49 | 1.50000 mL | 0.07084 mL | 0.16409 mL | 0.02500 mL | 12.201 | -0.00522 | 0.79337      | 0.00029 | 10.0 s      |
| 27:24.3 | Assay volumes | 1.75000 mL | 0.24443 mL | 0.16409 mL | 0.02500 mL |        |          |              |         |             |

## Assay Settings

| Setting                             | Value              | Original Value | Date/Time changed | Imported from |
|-------------------------------------|--------------------|----------------|-------------------|---------------|
| <b>General Settings</b>             |                    |                |                   |               |
| Analyst name                        | Dorothy Levorse    |                |                   |               |
| Separate reference vial             | Yes                |                |                   |               |
| <b>Standard Experiment Settings</b> |                    |                |                   |               |
| Number of titrations                | 1                  |                |                   |               |
| Minimum pH                          | 1.800              |                |                   |               |
| Maximum pH                          | 12.200             |                |                   |               |
| pH step between points of           | 0.200              |                |                   |               |
| Minimum titrant addition            | 0.00002 mL         |                |                   |               |
| Maximum titrant addition            | 0.10000 mL         |                |                   |               |
| Argon flow rate                     | 100%               |                |                   |               |
| Start titration using               | Cautious pH adjust |                |                   |               |
| <b>Advanced General Settings</b>    |                    |                |                   |               |

Sample name: **M07**  
 Assay name: **UV-metric pKa**  
 Assay ID: **17I-16002**  
 Filename: **C:\Sirius\_T3\17I-16002\_M07\_UV-metric pKa.t3r**

Experiment start time: **9/15/2017 11:47:51 PM**  
 Analyst: **Dorothy Levorse**  
 Instrument ID: **T311053**

## Assay Settings (continued)

| Setting                             | Value            | Original Value | Date/Time changed | Imported from |
|-------------------------------------|------------------|----------------|-------------------|---------------|
| Detect turbidity using              | Spectrometer     |                |                   |               |
| Monitor at a wavelength of          | 500.0 nm         |                |                   |               |
| Absorbance threshold of             | 0.100            |                |                   |               |
| Collect turbidity sensor data       | No               |                |                   |               |
| Stir after titrant addition for     | 5 seconds        |                |                   |               |
| For titrant addition, stir at       | 15%              |                |                   |               |
| <b><i>Titrant Pre-Dose</i></b>      |                  |                |                   |               |
| Titrant pre-dose                    | None             |                |                   |               |
| <b><i>Assay Medium</i></b>          |                  |                |                   |               |
| Cosolvent in use                    | No               |                |                   |               |
| ISA water volume                    | 1.50 mL          |                |                   |               |
| Water added                         | Automatic        |                |                   |               |
| After water addition, stir for      | 5 seconds        |                |                   |               |
| At a speed of                       | 15%              |                |                   |               |
| Buffer in use                       | Yes              |                |                   |               |
| Buffer type                         | Phosphate Buffer |                |                   |               |
| Volume of buffer introduced         | 0.025000 mL      |                |                   |               |
| Add buffer manually                 | Manual           |                |                   |               |
| After medium addition, stir for     | 5 seconds        |                |                   |               |
| <b><i>Sample Sonication</i></b>     |                  |                |                   |               |
| Sonicate                            | No               |                |                   |               |
| <b><i>Sample Dissolution</i></b>    |                  |                |                   |               |
| Perform a dissolution stage         | No               |                |                   |               |
| <b><i>Carbonate purge</i></b>       |                  |                |                   |               |
| Perform a carbonate purge           | No               |                |                   |               |
| <b><i>Temperature Control</i></b>   |                  |                |                   |               |
| Wait for temperature                | Yes              |                |                   |               |
| Required start temperature          | 25.0°C           |                |                   |               |
| Acceptable deviation                | 0.5°C            |                |                   |               |
| Time to wait                        | 60 seconds       |                |                   |               |
| Stir speed of                       | 15%              |                |                   |               |
| <b><i>Titration 1</i></b>           |                  |                |                   |               |
| Titrate from                        | Low to high pH   |                |                   |               |
| Adjust to start pH                  | Yes              |                |                   |               |
| After pH adjust stir for            | 10 seconds       |                |                   |               |
| <b><i>Data Point Stability</i></b>  |                  |                |                   |               |
| Stir during data point collection   | Yes              |                |                   |               |
| For point collection, stir at       | 15%              |                |                   |               |
| Delay before data point collection  | 0 seconds        |                |                   |               |
| Number of points to average         | 20 points        |                |                   |               |
| Time interval between points        | 0.50 seconds     |                |                   |               |
| Required maximum standard deviation | 0.00500 dpH/dt   |                |                   |               |
| Stability timeout after             | 60 seconds       |                |                   |               |
| <b><i>Experiment cleanup</i></b>    |                  |                |                   |               |
| Adjust pH to cleanup                | To start pH      |                |                   |               |
| And then stir for                   | 60 seconds       |                |                   |               |
| For cleaning, stir at               | 20%              |                |                   |               |
| Then add water volume               | 0.25 mL          |                |                   |               |
| And then stir for                   | 30 seconds       |                |                   |               |

## Calibration Settings

| Setting                   | Value  | Date/Time changed     | Imported from             |
|---------------------------|--------|-----------------------|---------------------------|
| Four-Plus alpha           | 0.112  | 9/15/2017 11:47:51 PM | C:\Sirius_T3\HCl17I15.t3r |
| Four-Plus S               | 1.0006 | 9/15/2017 11:47:51 PM | C:\Sirius_T3\HCl17I15.t3r |
| Four-Plus jH              | 0.7    | 9/15/2017 11:47:51 PM | C:\Sirius_T3\HCl17I15.t3r |
| Four-Plus jOH             | -0.6   | 9/15/2017 11:47:51 PM | C:\Sirius_T3\HCl17I15.t3r |
| Base concentration factor | 1.015  | 9/15/2017 11:47:51 PM | C:\Sirius_T3\KOH17I11.t3r |

Sample name: **M07**  
 Assay name: **UV-metric pKa**  
 Assay ID: **17I-16002**  
 Filename: **C:\Sirius\_T3\17I-16002\_M07\_UV-metric pKa.t3r**

Experiment start time: **9/15/2017 11:47:51 PM**  
 Analyst: **Dorothy Levorse**  
 Instrument ID: **T311053**

## Calibration Settings (continued)

| Setting                   | Value | Date/Time changed     | Imported from             |
|---------------------------|-------|-----------------------|---------------------------|
| Acid concentration factor | 1.003 | 9/15/2017 11:47:51 PM | C:\Sirius_T3\HCl17I15.t3r |

## Instrument Settings

| Setting                          | Value                      | Batch Id    | Install date           |
|----------------------------------|----------------------------|-------------|------------------------|
| Instrument owner                 | Merck                      |             |                        |
| Instrument ID                    | T311053                    |             |                        |
| Instrument type                  | T3 Simulator               |             |                        |
| Software version                 | 1.1.3.0                    |             |                        |
| Dispenser module                 |                            | T3DM1100253 | 3/31/2009 5:24:52 AM   |
| Dispenser 0                      | Water                      |             | 3/31/2009 5:25:05 AM   |
| Syringe volume                   | 2.5 mL                     |             |                        |
| Firmware version                 | 1.2.1(r2)                  |             |                        |
| Titration                        | Water (0.15 M KCl)         | 8-18-17     | 9/8/2017 8:22:43 AM    |
| Dispenser 2                      | Acid                       |             | 3/31/2009 5:25:11 AM   |
| Syringe volume                   | 0.5 mL                     |             |                        |
| Firmware version                 | 1.2.1(r2)                  |             |                        |
| Titration                        | Acid (0.5 M HCl)           | 166940      | 9/8/2017 8:21:27 AM    |
| Dispenser 1                      | Base                       |             | 3/31/2009 5:25:21 AM   |
| Syringe volume                   | 0.5 mL                     |             |                        |
| Firmware version                 | 1.2.1(r2)                  |             |                        |
| Titration                        | Base (0.5 M KOH)           | 01/06/17    | 9/8/2017 8:20:03 AM    |
| Dispenser 5                      | Cosolvent                  |             | 3/31/2009 5:26:24 AM   |
| Syringe volume                   | 2.5 mL                     |             |                        |
| Firmware version                 | 1.2.1(r2)                  |             |                        |
| Distribution valve 5             | Distribution Valve         |             | 3/31/2009 5:28:19 AM   |
| Firmware version                 | 1.1.3                      |             |                        |
| Port A                           | Methanol (80%, 0.15 M KCl) | 8-15-17     | 9/13/2017 11:23:11 AM  |
| Dispenser 3                      | Buffer                     |             | 8/3/2010 5:05:16 AM    |
| Syringe volume                   | 0.5 mL                     |             |                        |
| Firmware version                 | 1.2.1(r2)                  |             |                        |
| Titration                        | Phosphate Buffer           |             | 9/12/2017 11:32:29 AM  |
| Dispenser 6                      | Octanol                    |             | 10/22/2010 10:52:43 AM |
| Syringe volume                   | 0.5 mL                     |             |                        |
| Firmware version                 | 1.2.1(r2)                  |             |                        |
| Titration                        | Octanol                    | 9-14-17     | 9/14/2017 9:30:38 AM   |
| Titration                        |                            | T3TM1100153 | 3/31/2009 5:24:17 AM   |
| Horizontal axis firmware version | 1.17 AI1DI2DO2 Stepper 2   |             |                        |
| Vertical axis firmware version   | 1.17 AI1DI2DO2 Stepper 2   |             |                        |
| Chassis I/O firmware version     | 1.11 AI1DI0DO4 Norgren I/O |             |                        |
| Probe I/O firmware version       | 1.1.1                      |             |                        |
| Electrode                        | T3 Electrode               | T3E0769     | 8/15/2017 9:21:54 AM   |
| E0 calibration                   | -8.33 mV                   |             | 9/15/2017 11:48:15 PM  |
| Filling solution                 | 3M KCl                     | KCL095      | 9/13/2017 8:16:19 AM   |
| Liquids                          |                            |             |                        |
| Wash 1                           | 50% IPA:50% Water          |             | 9/15/2017 8:38:18 AM   |
| Wash 2                           | 0.5% Triton X-100 in H2O   |             | 9/15/2017 8:38:22 AM   |
| Buffer position 1                | pH7 Wash                   |             | 9/15/2017 8:38:24 AM   |
| Buffer position 2                | pH 7                       |             | 9/15/2017 8:38:27 AM   |
| Storage position                 |                            |             | 9/15/2017 8:38:55 AM   |
| Wash water                       | 3.6e+003 mL                | 9-11-17     | 9/11/2017 3:28:43 PM   |
| Waste                            | 6.5e+003 mL                |             | 9/11/2017 3:28:49 PM   |
| Temperature controller           |                            |             | 8/5/2010 6:35:13 AM    |
| Turbidity detector               |                            |             | 3/31/2009 5:24:45 AM   |
| Spectrometer                     |                            | 072390      | 11/23/2010 11:22:28 AM |
| Dip probe                        |                            | 11086       |                        |
| Wavelength coefficient A0        | 185.563                    |             |                        |
| Wavelength coefficient A1        | 2.17439                    |             |                        |





Sample name: **M07**  
Assay name: **UV-metric pKa**  
Assay ID: **17I-16002**  
Filename: **C:\Sirius\_T3\17I-16002\_M07\_UV-metric pKa.t3r**

Experiment start time: **9/15/2017 11:47:51 PM**  
Analyst: **Dorothy Levorse**  
Instrument ID: **T311053**

## Instrument Settings (continued)

| Setting                                     | Value                      | Batch Id    | Install date           |
|---|----------------------------|-------------|------------------------|
| Wavelength coefficient A2                   | -0.000285622               |             |                        |
| Total lamp lit time                         | 114:03:31                  |             | 11/23/2010 11:22:28 AM |
| Calibrated on                               | 9/6/2017 8:33:02 AM        |             |                        |
| Integration time                            | 11                         |             |                        |
| Scans averaged                              | 10                         |             |                        |
| Autoloader                                  |                            | T3AL1100237 | 11/10/2015 9:34:13 AM  |
| Left-right axis firmware version            | 1.17 AI1DI2DO2 Stepper 2   |             |                        |
| Front-back axis firmware version            | 1.17 AI1DI2DO2 Stepper 2   |             |                        |
| Vertical axis firmware version              | 1.17 AI1DI2DO2 Stepper 2   |             |                        |
| Chassis I/O firmware version                | 1.11 AI1DI0DO4 Norgren I/O |             |                        |
| Configuration                               |                            |             |                        |
| Alternate titration position                | Titration position         |             |                        |
| Alternate reference position                | Reference position         |             |                        |
| Maximum standard vial volume                | 3.50 mL                    |             |                        |
| Maximum alternate vial volume               | 25.00 mL                   |             |                        |
| Automatic action idle period                | 5 minute(s)                |             |                        |
| Titrant tube volume                         | 1.3 mL                     |             |                        |
| Syringe flush count                         | 3.50                       |             |                        |
| Flowing wash pump volume                    | 20.0 mL                    |             |                        |
| Flowing wash stir duration                  | 5 s                        |             |                        |
| Flowing wash stir speed                     | 30%                        |             |                        |
| Solvent wash stir duration                  | 5 s                        |             |                        |
| Solvent wash stir speed                     | 30%                        |             |                        |
| Surfactant wash stir duration               | 5 s                        |             |                        |
| Surfactant wash stir speed                  | 30%                        |             |                        |
| E0 calibration minimum number of points     | 10                         |             |                        |
| E0 calibration maximum standard deviation   | 0.01500                    |             |                        |
| E0 calibration timeout period               | 60 s                       |             |                        |
| E0 calibration stir duration                | 5 s                        |             |                        |
| E0 calibration preparation stir speed       | 30%                        |             |                        |
| E0 calibration buffer wash stir duration    | 5 s                        |             |                        |
| E0 calibration buffer wash stir speed       | 30%                        |             |                        |
| E0 calibration reading stir speed           | 0%                         |             |                        |
| Spectrometer calibration stir duration      | 5 s                        |             |                        |
| Spectrometer calibration stir speed         | 30%                        |             |                        |
| Spectrometer calibration wash pump volume   | 20.0 mL                    |             |                        |
| Spectrometer calibration wash stir duration | 5 s                        |             |                        |
| Spectrometer calibration wash stir speed    | 30%                        |             |                        |
| Overhead dispense height                    | 10000                      |             |                        |

## Refinement Settings

| Setting                               | Value        | Default value |
|---------------------------------------|--------------|---------------|
| Turbidity detection method            | Spectrometer | Spectrometer  |
| Turbidity wavelength to assess        | 500.0 nm     | 500.0 nm      |
| Turbidity maximum absorbance          | 0.100        | 0.100         |
| Turbidity probe threshold             | 50.00        | 50.00         |
| Exclude turbid points                 | Yes          | Yes           |
| Low intensity warning threshold       | 100          | 100           |
| Minimum absorbance change threshold   | 0.100        | 0.100         |
| Eigenvector autocorrelation threshold | 0.80         | 0.80          |
| Maximum RMSD severe warning           | 0.250        | 0.250         |
| Maximum RMSD warning                  | 0.050        | 0.050         |