

Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

18C-24004 Instrument ID: **T312060**

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH +) 0.77 ±0.03 (n=50) logP (neutral X) 3.95 ±0.02 (n=50)

RMSD 0.422

18C-24004 Points 1 to 23

M04_octanol concentration factor 0.620
Carbonate 0.0570 mM
Acidity error 0.07992 mM

18C-24004 Points 24 to 50

M04_octanol concentration factor 0.641
Carbonate 0.1306 mM
Acidity error -0.26645 mM

18C-24004 Points 51 to 77

M04_octanol concentration factor 0.822
Carbonate 0.1176 mM
Acidity error -0.01548 mM

Warnings and errors

Errors None Warnings None

Sample logD and percent species

рН	M04_octanol	M04_octanol	M04_octanol	M04_octanol	M04_octanol	Comment
	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.78	14.18 %	0.00 %	84.45 %	1.37 %	
1.200	0.79	14.07 %	0.00 %	83.78 %	2.15 %	Stomach pH
2.000	0.84	12.63 %	0.00 %	75.19 %	12.18 %	
3.000	1.19	6.03 %	0.01 %	35.87 %	58.10 %	
4.000	2.01	0.97 %	0.01 %	5.76 %	93.26 %	
5.000	2.94	0.10 %	0.01 %	0.61 %	99.27 %	
6.000	3.67	0.01 %	0.01 %	0.06 %	99.92 %	
6.500	3.84	0.00 %	0.01 %	0.02 %	99.97 %	
7.000	3.92	0.00 %	0.01 %	0.01 %	99.98 %	
7.400	3.94	0.00 %	0.01 %	0.00 %	99.99 %	Blood pH
8.000	3.95	0.00 %	0.01 %	0.00 %	99.99 %	
9.000	3.95	0.00 %	0.01 %	0.00 %	99.99 %	
10.000	3.95	0.00 %	0.01 %	0.00 %	99.99 %	
11.000	3.95	0.00 %	0.01 %	0.00 %	99.99 %	
12.000	3.95	0.00 %	0.01 %	0.00 %	99.99 %	



Sample name: Assay name:

Assay ID:


M04_octanol pH-metric high logP

18C-24004

Experiment start time: 3/24/2018 4:11:23 AM Analyst:

Dorothy Levorse

Instrument ID: T312060






Filename:

Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM Analyst: Assay name: pH-metric high logP **Dorothy Levorse** Assay ID:

18C-24004 Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

Graphs (continued)





Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM Assay name: pH-metric high logP Analyst: **Dorothy Levorse** Assay ID:

18C-24004 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-24004 Points 1 to 23

Overall results

RMSD 0.635 Average ionic strength 0.157 M Average temperature 24.9°C Partition ratio 0.0185:1

Analyte concentration range 3292.8 µM to 3384.5 µM

Total points considered 19 of 23

Warnings and errors

None Errors

Warnings Sample concentration factor out of range

Four-Plus parameters

Alpha	0.119	3/24/2018 4:11:23 AM	C:\Sirius_T3\HCl18C23.t3r
S	0.9972	3/24/2018 4:11:23 AM	C:\Sirius_T3\HCl18C23.t3r
jΗ	0.9	3/24/2018 4:11:23 AM	C:\Sirius_T3\HCl18C23.t3r
jОН	-0.3	3/24/2018 4:11:23 AM	C:\Sirius_T3\HCl18C23.t3r


Titrants


0.50 M HCI 0.997124 3/24/2018 4:11:23 AM C:\Sirius_T3\HCl18C23.t3r 0.50 M KOH 1.003190 3/24/2018 4:11:23 AM C:\Sirius_T3\KOH18C23.t3r

Sample

M04 octanol concentration factor 0.620 Base pKa 1 5.97 logP (XH+) 0.19 logP (neutral X) 3.88

Sample graphs







Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP

18C-24004


C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

Experiment start time: 3/24/2018 4:11:23 AM Analyst:

Dorothy Levorse

Instrument ID: T312060

Sample graphs (continued)



Sample logD and percent species


рН	M04_octanol	M04_octanol	M04_octanol	M04_octanol	M04_octanol	Comment
	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.21	97.07 %	0.00 %	2.78 %	0.15 %	
1.200	0.22	96.99 %	0.00 %	2.78 %	0.23 %	Stomach pH
2.000	0.37	95.80 %	0.01 %	2.74 %	1.44 %	
3.000	0.99	84.72 %	0.09 %	2.43 %	12.76 %	
4.000	1.91	39.28 %	0.42 %	1.13 %	59.18 %	
5.000	2.87	6.17 %	0.66 %	0.18 %	92.99 %	
6.000	3.59	0.65 %	0.70 %	0.02 %	98.63 %	
6.500	3.77	0.21 %	0.70 %	0.01 %	99.08 %	
7.000	3.84	0.07 %	0.71 %	0.00 %	99.23 %	
7.400	3.87	0.03 %	0.71 %	0.00 %	99.27 %	Blood pH
8.000	3.88	0.01 %	0.71 %	0.00 %	99.29 %	
9.000	3.88	0.00 %	0.71 %	0.00 %	99.29 %	
10.000	3.88	0.00 %	0.71 %	0.00 %	99.29 %	
11.000	3.88	0.00 %	0.71 %	0.00 %	99.29 %	
12.000	3.88	0.00 %	0.71 %	0.00 %	99.29 %	

Carbonate and acidity



Carbonate 0.057 mM Acidity error 0.080 mM

Other graphs





Assay ID: Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP

18C-24004


Experiment start time: 3/24/2018 4:11:23 AM

Analyst: **Dorothy Levorse**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

Other graphs (continued)





Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM Assay name: pH-metric high logP Analyst: **Dorothy Levorse** Assay ID:

T312060 18C-24004 Instrument ID:

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-24004 Points 24 to 50

Overall results

RMSD 0.240 Average ionic strength 0.163 M Average temperature 25.0°C Partition ratio 0.0642:1

Analyte concentration range 2954.9 µM to 3043.2 µM

Total points considered 22 of 27

Warnings and errors

Errors None

Warnings Sample concentration factor out of range

Four-Plus parameters

Alpha 0.119 3/24/2018 4:11:23 AM C:\Sirius_T3\HCl18C23.t3r 0.9972 3/24/2018 4:11:23 AM C:\Sirius_T3\HCl18C23.t3r S jΗ 0.9 3/24/2018 4:11:23 AM C:\Sirius_T3\HCl18C23.t3r jOH -0.33/24/2018 4:11:23 AM C:\Sirius_T3\HCl18C23.t3r


Titrants


0.50 M HCI 0.997124 3/24/2018 4:11:23 AM C:\Sirius T3\HCl18C23.t3r 0.50 M KOH 1.003190 3/24/2018 4:11:23 AM C:\Sirius_T3\KOH18C23.t3r

Sample

M04_octanol concentration factor 0.641 Base pKa 1 5.97 logP(XH +)0.19 logP (neutral X) 3.92

Sample graphs







Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP

18C-24004


Experiment start time: 3/24/2018 4:11:23 AM Analyst:

Dorothy Levorse

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M04_octanol	M04_octanol	M04_octanol	M04_octanol	M04_octanol	Comment
	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.21	90.48 %	0.00 %	9.00 %	0.52 %	
1.200	0.23	90.21 %	0.00 %	8.97 %	0.82 %	Stomach pH
2.000	0.39	86.44 %	0.01 %	8.60 %	4.96 %	
3.000	1.02	59.74 %	0.06 %	5.94 %	34.25 %	
4.000	1.95	14.61 %	0.16 %	1.45 %	83.78 %	
5.000	2.91	1.71 %	0.18 %	0.17 %	97.94 %	
6.000	3.63	0.17 %	0.19 %	0.02 %	99.62 %	
6.500	3.81	0.06 %	0.19 %	0.01 %	99.75 %	
7.000	3.88	0.02 %	0.19 %	0.00 %	99.79 %	
7.400	3.90	0.01 %	0.19 %	0.00 %	99.81 %	Blood pH
8.000	3.92	0.00 %	0.19 %	0.00 %	99.81 %	
9.000	3.92	0.00 %	0.19 %	0.00 %	99.81 %	
10.000	3.92	0.00 %	0.19 %	0.00 %	99.81 %	
11.000	3.92	0.00 %	0.19 %	0.00 %	99.81 %	
12.000	3.92	0.00 %	0.19 %	0.00 %	99.81 %	

Carbonate and acidity



Other graphs





Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP

18C-24004


Experiment start time: 3/24/2018 4:11:23 AM

Analyst: **Dorothy Levorse**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

Other graphs (continued)





Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM Assay name: pH-metric high logP Analyst: **Dorothy Levorse** Assay ID:

T312060 18C-24004 Instrument ID:

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-24004 Points 51 to 77

Overall results

RMSD 0.309 Average ionic strength 0.168 M 25.0°C Average temperature Partition ratio 0.2800:1

Analyte concentration range 2304.6 µM to 2362.2 µM

Total points considered 22 of 27

Warnings and errors

Errors None Warnings None

Four-Plus parameters

Alpha 0.119 3/24/2018 4:11:23 AM C:\Sirius_T3\HCl18C23.t3r 0.9972 3/24/2018 4:11:23 AM C:\Sirius_T3\HCl18C23.t3r S jΗ 0.9 3/24/2018 4:11:23 AM C:\Sirius_T3\HCl18C23.t3r jOH -0.33/24/2018 4:11:23 AM C:\Sirius_T3\HCl18C23.t3r


Titrants


0.50 M HCI 0.997124 3/24/2018 4:11:23 AM C:\Sirius T3\HCl18C23.t3r 0.50 M KOH 1.003190 3/24/2018 4:11:23 AM C:\Sirius_T3\KOH18C23.t3r

Sample

M04_octanol concentration factor 0.822 Base pKa 1 5.97 logP(XH +)0.19 logP (neutral X) 3.67

Sample graphs







Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP

18C-24004


C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

Experiment start time: 3/24/2018 4:11:23 AM Analyst:

Dorothy Levorse

Instrument ID: T312060

Sample graphs (continued)



Sample logD and percent species


рН	M04_octanol	M04_octanol	M04_octanol	M04_octanol	M04_octanol	Comment
	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.20	69.08 %	0.00 %	29.96 %	0.97 %	
1.200	0.21	68.69 %	0.00 %	29.79 %	1.52 %	Stomach pH
2.000	0.31	63.55 %	0.01 %	27.56 %	8.88 %	
3.000	0.81	35.31 %	0.04 %	15.31 %	49.33 %	
4.000	1.71	6.49 %	0.07 %	2.81 %	90.63 %	
5.000	2.66	0.71 %	0.08 %	0.31 %	98.91 %	
6.000	3.38	0.07 %	0.08 %	0.03 %	99.82 %	
6.500	3.56	0.02 %	0.08 %	0.01 %	99.89 %	
7.000	3.63	0.01 %	0.08 %	0.00 %	99.91 %	
7.400	3.65	0.00 %	0.08 %	0.00 %	99.92 %	Blood pH
8.000	3.66	0.00 %	0.08 %	0.00 %	99.92 %	
9.000	3.67	0.00 %	0.08 %	0.00 %	99.92 %	
10.000	3.67	0.00 %	0.08 %	0.00 %	99.92 %	
11.000	3.67	0.00 %	0.08 %	0.00 %	99.92 %	
12.000	3.67	0.00 %	0.08 %	0.00 %	99.92 %	

Carbonate and acidity



Carbonate 0.118 mM Acidity error -0.015 mM

Other graphs





Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP


18C-24004


C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

Experiment start time: 3/24/2018 4:11:23 AM Analyst:

Dorothy Levorse

Instrument ID: T312060







Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

18C-24004 Instrument ID: **T312060**

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

Assay Model

Assay ID:

Settings	Value	Date/Time changed	Imported from
Sample name	M04_octanol	3/9/2018 4:34:21 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.001450 g	3/23/2018 5:00:45 PM	User entered value
Formula weight	269.73 g/mol	3/9/2018 4:34:21 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	269.73	3/9/2018 4:34:21 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	1	3/9/2018 4:34:21 PM	User entered value
Sample is a	Base	3/9/2018 4:34:21 PM	User entered value
pKa 1	5.97	3/9/2018 4:34:21 PM	User entered value
logp (XH +)	0.19	3/9/2018 4:34:33 PM	User entered value
logP (neutral X)	3.50	3/23/2018 2:32:46 PM	User entered value

Events

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
5:00.3	Initial pH = 8.35	4 50000!	0.05445!	0.00747!	0.00004!	0.040	0.00540	0.05044	0.00050	10.00
7:60.0	Data point 1		0.05115 mL						0.00050	
8:46.2	Data point 2		0.05115 mL					0.41042	0.00018	
9:21.8	Data point 3		0.05115 mL						0.00095	
9:57.4	Data point 4		0.05115 mL						0.00047	
10:32.9			0.05115 mL						0.00022	
11:08.3	•		0.05115 mL					0.00064	0.00011	
11:43.9	Data point 7		0.05115 mL						0.00010	
12:19.7			0.05115 mL						0.00013	
	Data point 9		0.05115 mL						0.00039	
13:31.1			0.05115 mL						0.00050	
14:06.6			0.05115 mL						0.00068	
14:47.7			0.05115 mL						0.00050	
	Data point 13		0.05115 mL						0.00084	
	Data point 14		0.05115 mL						0.00085	
	Data point 15		0.05115 mL						0.00088	
	Data point 16		0.05115 mL						0.00099	
18:12.1	Data point 17		0.05115 mL						0.00091	
	Data point 18		0.05115 mL						0.00092	
20:28.4	Data point 19	1.50000 mL	0.05115 mL	0.05122 mL	0.03001 mL	7.844	-0.08873	0.99672	0.00439	Timed out at
										59.5 s
21:58.9	Data point 20	1.50000 mL	0.05115 mL	0.05127 mL	0.03001 mL	8.410	-0.04014	0.97727	0.00201	Timed out at
										59.5 s
23:29.4	Data point 21	1.50000 mL	0.05115 mL	0.05132 mL	0.03001 mL	8.719	-0.02194	0.97381	0.00110	Timed out at
										59.5 s
25:05.0	Data point 22		0.05115 mL						0.00098	
26:05.0	Data point 23		0.05115 mL						0.00097	42.0 s
27:41.8	Data point 24		0.10503 mL					0.14587	0.00044	10.0 s
28:28.0	Data point 25	1.50000 mL	0.10503 mL	0.06816 mL	0.11002 mL	2.160	0.00273	0.60068	0.00017	10.5 s
29:04.3	Data point 26	1.50000 mL	0.10503 mL	0.07947 mL	0.11002 mL	2.371	0.00238	0.28453	0.00022	10.0 s
29:39.9	Data point 27	1.50000 mL	0.10503 mL	0.08652 mL	0.11002 mL	2.601	0.00077	0.08864	0.00013	10.5 s
30:15.9	Data point 28	1.50000 mL	0.10503 mL	0.09087 mL	0.11002 mL	2.792	-0.01433	0.72766	0.00083	10.5 s
30:51.8	Data point 29	1.50000 mL	0.10503 mL	0.09393 mL	0.11002 mL	2.985	0.01516	0.73138	0.00088	10.0 s
31:27.3	Data point 30	1.50000 mL	0.10503 mL	0.09626 mL	0.11002 mL	3.161	-0.01511	0.64124	0.00093	11.0 s
32:14.1	Data point 31	1.50000 mL	0.10503 mL	0.09828 mL	0.11002 mL	3.353	-0.00936	0.31606	0.00082	10.0 s
32:49.5			0.10503 mL						0.00016	10.5 s
	Data point 33		0.10503 mL						0.00035	
34:06.0	•		0.10503 mL						0.00029	
34:41.4	•		0.10503 mL						0.00032	
	Data point 36		0.10503 mL						0.00039	
05.50.7		4 50000	0.40500	0.40040	0.44000	4 740	0.00000	0.00007	0.0000	10.5

1.50000 mL 0.10503 mL 0.10310 mL 0.11002 mL 4.712 0.00092 0.00227

35:52.7 Data point 37

0.00095 10.5 s



Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM pH-metric high logP Analyst: Assay name: **Dorothy Levorse** Assay ID:

18C-24004 Instrument ID: T312060

Filename: $C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric\ high\ logP.t3r$

Events (continued)

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
36:28.6	Data point 38	1 50000 ml	0 10503 ml	0 10332 ml	0.11002 mL	4 945	-0 00605	0.33921	0.00051	
37:05.0	Data point 39				0.11002 mL				0.00086	
37:47.6	Data point 40				0.11002 mL			0.09101	0.00076	
38:31.1	Data point 41				0.11002 mL			0.74805	0.00083	
39:17.2	Data point 42				0.11002 mL			0.88071	0.00092	
40:25.8	Data point 43				0.11002 mL			0.88671	0.00096	
41:50.7	Data point 44		0.10503 mL			6.931	-0.02021	0.98765		Timed out
11.00.7	Data point 44	1.00000 1112	0.10000 IIIL	0.10001 IIIL	0.11002 IIIL	0.001	0.02021	0.00100	0.00100	at 59.5 s
43:21.2	Data point 45	1.50000 mL	0.10503 mL	0.10386 mL	0.11002 mL	7.358	-0.05882	0.99495	0.00291	Timed out
44:51.6	Data point 46	1.50000 mL	0.10503 mL	0.10390 mL	0.11002 mL	7.860	-0.06638	0.98979	0.00330	
	- · · · · · ·	4 = 0000	0.40=00	0.40000	0.44000					at 59.5 s
46:27.2	Data point 47	1.50000 mL	0.10503 mL	0.10398 mL	0.11002 mL	8.405	-0.03574	0.99078	0.00177	Timed out at 59.5 s
47:57.6	Data point 48	1.50000 mL	0.10503 mL	0.10402 mL	0.11002 mL	8.751	-0.01947	0.94555	0.00099	58.0 s
49:26.2	Data point 49				0.11002 mL			0.97546	0.00079	
50:38.7	Data point 50				0.11002 mL			0.92626	0.00097	28.5 s
52:25.2	Data point 51	1.50000 mL	0.16148 mL	0.10421 mL	0.51002 mL	1.964	-0.00413	0.66600	0.00025	10.5 s
53:11.9	Data point 52	1.50000 mL	0.16148 mL	0.12244 mL	0.51002 mL	2.173	0.00880	0.67097	0.00053	10.0 s
53:47.6	Data point 53	1.50000 mL	0.16148 mL	0.13434 mL	0.51002 mL	2.382	0.01052	0.53800	0.00071	10.0 s
54:23.2	Data point 54				0.51002 mL			0.05969	0.00031	10.0 s
54:58.7	Data point 55	1.50000 mL	0.16148 mL	0.14713 mL	0.51002 mL	2.812	-0.00306	0.69301	0.00018	10.5 s
55:34.7	Data point 56				0.51002 mL			0.20419	0.00082	10.0 s
56:10.2	Data point 57				0.51002 mL			0.02683	0.00037	
56:46.2	Data point 58				0.51002 mL			0.19314	0.00092	10.0 s
57:21.7	Data point 59				0.51002 mL			0.14616	0.00059	10.0 s
58:07.4	Data point 60				0.51002 mL			0.00501	0.00089	
58:53.1	Data point 61				0.51002 mL			0.00430	0.00093	
59:38.9	Data point 62				0.51002 mL			0.38768	0.00096	
1:00:28.1	Data point 63				0.51002 mL			0.00611	0.00072	
1:01:18.9					0.51002 mL			0.00902	0.00061	10.5 s
	Data point 65				0.51002 mL			0.92089	0.00096	
1:02:36.9					0.51002 mL			0.92187	0.00096	
1:03:32.5					0.51002 mL			0.67266	0.00092	
1:04:21.6					0.51002 mL			0.96692	0.00096	
1:05:43.7	Data point 69	1.50000 mL	0.16148 mL	0.16070 mL	0.51002 mL	6.858	-0.05249	0.96745	0.00264	Timed out at 59.5 s
1:07:14.2	Data point 70	1.50000 mL	0.16148 mL	0.16075 mL	0.51002 mL	7.233	-0.05610	0.97962	0.00280	Timed out at 59.5 s
1:08:44.7	Data point 71	1.50000 mL	0.16148 mL	0.16079 mL	0.51002 mL	7.571	-0.05662	0.99069	0.00281	
1:10:15.2	Data point 72	1.50000 mL	0.16148 mL	0.16084 mL	0.51002 mL	7.982	-0.05547	0.99203	0.00275	Timed out
1:11:45.7	Data point 73	1.50000 mL	0.16148 mL	0.16089 mL	0.51002 mL	8.235	-0.05197	0.96437	0.00261	at 59.5 s Timed out
1:13:16.2	Data point 74	1.50000 mL	0.16148 mL	0.16094 mL	0.51002 mL	8.444	-0.02913	0.98679	0.00145	at 59.5 s Timed out
4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4	D	4 50000	0.40445	0.40005	0.54005	0.000	0.646==	0.00700	0.00055	at 59.5 s
	Data point 75				0.51002 mL				0.00098	
	Data point 76				0.51002 mL				0.00091	
	Data point 77 Assay volumes		0.16148 mL 0.16148 mL		0.51002 mL 0.51002 mL	8.988	-0.01977	0.96525	0.00099	28.5 S



Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM Assay name: pH-metric high logP **Dorothy Levorse**

Assay ID:	18C-24004	alatan/004000000	_	strument ID:	T312060			
Filename:	C:\Sirius_i 3\lvi	entap\20180323_exp3	3_logP_1 3-2\180	5-24004_M04_octar	ol_pH-metric high log			
Assay Set	Assay Settings							
Setting		Value	Original Value	Date/Time change	d Imported from			
General Sett	tings		_		-			
Analyst name)	Dorothy Levorse						
Standard Ex	periment Settings	•						
Number of tit		3						
Minimum pH		2.000						
Maximum pH		9.000						
pH step betw	een points of	0.200						
Minimum titra	ant addition	0.00002 mL						
Maximum titr	ant addition	0.10000 mL						
Argon flow ra	te	100%						
Start titration		Cautious pH adjust						
Advanced G	eneral Settings	, ,						
Detect turbidi	ty using	None						
	ity sensor data	No						
Collect UV sp		No						
Stir after titra	nt addition for	5 seconds						
For titrant add	dition, stir at	10%						
Titrant Pre-D								
Titrant pre-do	se	None						
Assay Mediu								
ISA water vol		1.50 mL						
Water added		Automatic						
Partition solve	ent type	Octanol						
Partition volu	me .	0.030 mL						
Partition solve	ent added	Automatic						
After partition	addition, stir for	1 seconds						
Sample Son								
Sonicate		Yes						
Adjust nH for	conjugation	No						

Adjust pH for sonication No

Sonicate for 60 seconds After sonication stir for 5 seconds

Sample Dissolution

Perform a dissolution stage Yes

Adjust and hold pH for dissolution To start pH Stir to dissolve for 120 seconds

For dissolution, stir at 10%

Carbonate purge

Perform a carbonate purge No

Temperature Control

Wait for temperature Yes 25.0°C Required start temperature Acceptable deviation 0.5°C Time to wait 60 seconds

Stir speed of 50%

Titration 1

Titrate from Low to high pH

Adjust to start pH Yes After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds

Stirrer speed for partitioning Titration 2

Stirrer speed for partitioning

Titrate from Low to high pH Add additional water 0.00 mL

50%

55%

Additional partition solvent volume 0.080 mL Additional partition solvent added Automatic After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds

Report by: Dorothy Levorse 3/26/2018 12:02:53 PM



Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM pH-metric high logP Analyst: Assay name: **Dorothy Levorse** Assay ID:

18C-24004 Instrument ID: T312060

Filename: $C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric\ high\ logP.t3r$

Assay Settings (continued)

Value	Original Value	Date/Time changed	Imported from
	_	_	•
Low to high pH			
0.00 mL			
0.400 mL			
Automatic			
30 seconds			
15 seconds			
60%			
No			
0 seconds			
20 points			
0.50 seconds			
0.00100 dpH/dt			
60 seconds			
	Low to high pH 0.00 mL 0.400 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt	Low to high pH 0.00 mL 0.400 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt	Low to high pH 0.00 mL 0.400 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.119	3/24/2018 4:11:23 AM	C:\Sirius_T3\HCl18C23.t3r
Four-Plus S	0.9972	3/24/2018 4:11:23 AM	C:\Sirius_T3\HCl18C23.t3r
Four-Plus jH	0.9	3/24/2018 4:11:23 AM	C:\Sirius_T3\HCl18C23.t3r
Four-Plus jOH	-0.3	3/24/2018 4:11:23 AM	C:\Sirius_T3\HCl18C23.t3r
Base concentration factor	1.003	3/24/2018 4:11:23 AM	C:\Sirius_T3\KOH18C23.t3r
Acid concentration factor	0.997	3/24/2018 4:11:23 AM	C:\Sirius_T3\HCl18C23.t3r

Instrument Settings

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module Dispenser 0 Syringe volume Firmware version	Water 2.5 mL 1.2.1(r2)	T3DM1200361	3/31/2009 6:24:52 AM 3/31/2009 6:25:05 AM
Titrant Dispenser 2 Syringe volume Firmware version	Water (0.15 M KCI) Acid 0.5 mL 1.2.1(r2)	02-06-2018	3/16/2018 11:09:18 AM 3/31/2009 6:25:11 AM
Titrant Dispenser 1 Syringe volume Firmware version	Acid (0.5 M HCI) Base 0.5 mL 1.2.1(r2)	03-16-2018	3/16/2018 10:56:23 AM 3/31/2009 6:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2)	3/22/2018	3/23/2018 9:34:17 AM 3/31/2009 6:26:24 AM
Distribution valve 5 Firmware version Port A Port B	Distribution Valve 1.1.3 Methanol (80%, 0.15 M KCI) Cyclohexane	02-08-2018 11-01-17	3/31/2009 6:28:19 AM 3/6/2018 10:28:59 AM 2/27/2018 11:37:57 AM
Dispenser 3 Syringe volume Firmware version Titrant Dispenser 6	Buffer 0.5 mL 1.2.1(r2) Dodecane Octanol	2018/01/31	8/3/2010 6:05:16 AM 2/28/2018 11:18:04 AM 10/22/2010 11:52:43 AM




Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

Assay ID: 18C-24004 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting Syringe volume	Value 0.5 mL	Batch Id	Install date
Firmware version Titrant	1.2.1(r2) Octanol	01-31-2018	2/27/2018 10:59:35 AM
Titrator Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2	T3TM1200161	3/31/2009 6:24:17 AM
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1DI0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0923	1/23/2018 3:01:00 PM
E0 calibration	+4.67 mV	KCL097	3/24/2018 4:11:51 AM
Filling solution Liquids	3M KCI	KCL097	3/23/2018 9:29:07 AM
Wash 1	50% IPA:50% Water		3/23/2018 9:29:12 AM
Wash 2	0.5% Trition X-100 in H20		3/23/2018 9:29:15 AM
Buffer position 1	pH7 Wash		3/23/2018 9:29:19 AM
Buffer position 2	pH 7		3/23/2018 9:29:21 AM
Storage position			3/23/2018 9:30:23 AM
Wash water	7.6e+003 mL	03-12-2018	3/12/2018 9:25:04 AM
Waste	2.6e+003 mL		3/12/2018 9:24:49 AM
Temperature controller			8/5/2010 7:35:13 AM 3/31/2009 6:24:45 AM
Turbidity detector Spectrometer		074811	11/23/2010 12:22:28 PM
Dip probe		10196	11/23/2010 12.22.20 FW
Wavelength coefficient A0	183.333	10100	
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	162:53:01		11/23/2010 12:22:28 PM
Calibrated on	2/27/2018 11:40:38 AM		
Integration time	40		
Scans averaged	10	T2AL 1200245	11/10/2015 10:34:13 AM
Autoloader Left-right axis firmware version	1.17 Al1Dl2DO2 Stepper 2	T3AL1200345	11/10/2015 10.34.13 AW
Front-back axis firmware version	1.17 Al1DI2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1DI0DO4 Norgren I/O		
Configuration	_		
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume Automatic action idle period	25.00 mL 5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir speed	5 s 30%		
Surfactant wash stir speed E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		





Sample name: M04_octanol Experiment start time: 3/24/2018 4:11:23 AM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

Assay ID: 18C-24004 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24004_M04_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting	Value	Batch Id	Install date		
Spectrometer calibration stir duration	5 s				
Spectrometer calibration stir speed	30%				
Spectrometer calibration wash pump volume	20.0 mL				
Spectrometer calibration wash stir duration	5 s				
Spectrometer calibration wash stir speed	30%				
Overhead dispense height	10000				
, ,					

Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00