

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

## pH-metric Result

logP (XH +) 0.91 ±0.07 (n=50) logP (neutral X) 3.05 ±0.02 (n=50)

#### 18C-06007 Points 1 to 23

M09\_octanol concentration factor 1.046
Carbonate 0.0318 mM
Acidity error 0.09711 mM

#### 18C-06007 Points 24 to 50

M09\_octanol concentration factor 0.952
Carbonate 0.0262 mM
Acidity error 0.24086 mM

#### 18C-06007 Points 51 to 74

M09\_octanol concentration factor 0.937
Carbonate 0.0860 mM
Acidity error 0.30528 mM

## Warnings and errors

Errors None Warnings None

рΗ

# Sample logD and percent species

M09\_octanol M09\_octanol M09\_octanol

-	_					
	logD	M09_octanolH	M09_octanol	M09_octanolH*	M09_octanol*	
1.000	0.92	10.79 %	0.00 %	88.69 %	0.52 %	
1.200	0.92	10.76 %	0.00 %	88.42 %	0.82 %	Stomach pH
2.000	0.94	10.31 %	0.00 %	84.72 %	4.96 %	•
3.000	1.11	7.13 %	0.03 %	58.56 %	34.29 %	
4.000	1.73	1.74 %	0.07 %	14.32 %	83.86 %	
5.000	2.54	0.20 %	0.09 %	1.67 %	98.03 %	
6.000	2.96	0.02 %	0.09 %	0.17 %	99.72 %	
6.500	3.02	0.01 %	0.09 %	0.05 %	99.85 %	
7.000	3.04	0.00 %	0.09 %	0.02 %	99.89 %	
7.400	3.05	0.00 %	0.09 %	0.01 %	99.90 %	Blood pH
8.000	3.05	0.00 %	0.09 %	0.00 %	99.91 %	•
9.000	3.05	0.00 %	0.09 %	0.00 %	99.91 %	
10.000	3.05	0.00 %	0.09 %	0.00 %	99.91 %	
11.000	3.05	0.00 %	0.09 %	0.00 %	99.91 %	
12.000	3.05	0.00 %	0.09 %	0.00 %	99.91 %	

M09\_octanol

M09\_octanol Comment



Sample name: M09\_octanol Assay name:


pH-metric high logP

18C-06007 Assay ID: Filename:

Experiment start time: 3/6/2018 3:40:58 PM

**Pion** Analyst: Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r



Log R


pH (Concentration scale)



pH-metric high logP Analyst: Pion Assay name: Assay ID: 18C-06007 Instrument ID: T312060 Filename:

C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# **Graphs** (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 1 of 3 18C-06007 Points 1 to 23

#### Overall results

RMSD 0.369
Average ionic strength 0.158 M
Average temperature 24.9°C
Partition ratio 0.0185 : 1

Analyte concentration range 3130.9 µM to 3230.2 µM

Total points considered 19 of 23

## Warnings and errors

Errors None Warnings None

## Four-Plus parameters

 Alpha
 0.124
 3/6/2018 3:40:58 PM
 C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r

 S
 0.9973
 3/6/2018 3:40:58 PM
 C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r

 jH
 0.9
 3/6/2018 3:40:58 PM
 C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r


 jOH
 -0.7
 3/6/2018 3:40:58 PM
 C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r


#### Titrants

#### Sample

M09\_octanol concentration factor
 M09\_octanol stoichiometry
 Chloride stoichiometry
 Base pKa 1
 logP (XH +)
 logP (neutral X)
 1.046
 1.000
 5.37
 0.75
 3.05

#### Sample graphs








Assay ID:

Sample name: M09\_octanol Experiment start time: 3/6/2018 3:40:58 PM

Assay name: pH-metric high logP Analyst: Pion Instrument ID: T312060 18C-06007

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r


# Sample graphs (continued)




## Sample logD and percent species

рН	M09_octanol	M09_octanol	M09_octanol	_	M09_octanol	Comment
	logD	M09_octanolH		M09_octanolH*	M09_octanol*	
1.000	0.75	90.46 %	0.00 %	9.46 %	0.08 %	
1.200	0.76	90.41 %	0.01 %	9.45 %	0.13 %	Stomach pH
2.000	0.79	89.77 %	0.04 %	9.39 %	0.80 %	•
3.000	1.02	83.46 %	0.36 %	8.72 %	7.46 %	
4.000	1.71	48.98 %	2.09 %	5.12 %	43.81 %	
5.000	2.53	9.55 %	4.07 %	1.00 %	85.38 %	
6.000	2.96	1.05 %	4.50 %	0.11 %	94.34 %	
6.500	3.02	0.34 %	4.54 %	0.04 %	95.09 %	
7.000	3.04	0.11 %	4.55 %	0.01 %	95.34 %	
7.400	3.05	0.04 %	4.55 %	0.00 %	95.40 %	Blood pH
8.000	3.05	0.01 %	4.55 %	0.00 %	95.44 %	·
9.000	3.05	0.00 %	4.55 %	0.00 %	95.45 %	
10.000	3.05	0.00 %	4.55 %	0.00 %	95.45 %	
11.000	3.05	0.00 %	4.55 %	0.00 %	95.45 %	
12.000	3.05	0.00 %	4.55 %	0.00 %	95.45 %	

# Carbonate and acidity



# Other graphs






Pion pH-metric high logP Assay name: Analyst: 18C-06007 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 2 of 3 18C-06007 Points 24 to 50

#### Overall results

RMSD 0.910
Average ionic strength 0.164 M
Average temperature 25.0°C
Partition ratio 0.0407 : 1

Analyte concentration range 2862.7 µM to 2955.4 µM

Total points considered 21 of 27

## Warnings and errors

Errors None Warnings None

## Four-Plus parameters


#### Titrants


0.50 M KOH 0.999845 3/6/2018 3:40:58 PM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M09\_octanol concentration factor 0.952
M09\_octanol stoichiometry 1.000
Chloride stoichiometry 1.000
Base pKa 1 5.37
logP (XH +) 0.76
logP (neutral X) 2.97

### Sample graphs








Assay ID:

Sample name: M09\_octanol Experiment start time: 3/6/2018 3:40:58 PM

Assay name: pH-metric high logP Analyst: Pion Instrument ID: T312060 18C-06007

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# Sample graphs (continued)



## Sample logD and percent species


рН	M09_octanol logD	M09_octanol M09_octanolH	M09_octanol	M09_octanol M09_octanolH*	M09_octanol*	Comment
1.000	0.76	8 <del>0</del> .91 %	0.00 %	<del>-</del> 18.96 %	0.13 %	
1.200	0.76	80.84 %	0.01 %	18.94 %	0.21 %	Stomach pH
2.000	0.79	79.94 %	0.03 %	18.73 %	1.30 %	•
3.000	0.99	71.38 %	0.30 %	16.73 %	11.59 %	
4.000	1.64	34.48 %	1.47 %	8.08 %	55.97 %	
5.000	2.45	5.59 %	2.38 %	1.31 %	90.72 %	
6.000	2.88	0.60 %	2.54 %	0.14 %	96.72 %	
6.500	2.94	0.19 %	2.55 %	0.04 %	97.21 %	
7.000	2.96	0.06 %	2.56 %	0.01 %	97.37 %	
7.400	2.97	0.02 %	2.56 %	0.01 %	97.41 %	Blood pH
8.000	2.97	0.01 %	2.56 %	0.00 %	97.43 %	
9.000	2.97	0.00 %	2.56 %	0.00 %	97.44 %	
10.000	2.97	0.00 %	2.56 %	0.00 %	97.44 %	
11.000	2.97	0.00 %	2.56 %	0.00 %	97.44 %	
12.000	2.97	0.00 %	2.56 %	0.00 %	97.44 %	


# **Carbonate and acidity**



Carbonate 0.026 mM Acidity error 0.241 mM

# Other graphs





Experiment start time: 3/6/2018 3:40:58 PM Pion

T312060



Assay name:

Assay ID: Filename:

Sample name: M09\_octanol


pH-metric high logP 18C-06007

Instrument ID:

 $C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric\ high\ logP.t3r$ 

Analyst:

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 3 of 3 18C-06007 Points 51 to 74

#### Overall results

RMSD 0.405
Average ionic strength 0.170 M
Average temperature 25.0°C
Partition ratio 0.0926 : 1

Analyte concentration range 2547.3 µM to 2626.2 µM

Total points considered 19 of 24

## Warnings and errors

Errors None Warnings None

# Four-Plus parameters


Alpha 0.124 3/6/2018 3:40:58 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r S 0.9973 3/6/2018 3:40:58 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r jH 0.9 3/6/2018 3:40:58 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r jOH -0.7 3/6/2018 3:40:58 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r


#### Titrants

#### Sample

M09\_octanol concentration factor
 M09\_octanol stoichiometry
 Chloride stoichiometry
 Base pKa 1
 logP (XH +)
 logP (neutral X)
 0.937
 1.000
 5.37
 0.76
 3.01

#### Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r


# Sample graphs (continued)




## Sample logD and percent species

рН	M09_octanol logD	M09_octanol M09_octanolH	M09_octanol M09_octanol	M09_octanol M09_octanolH*	M09_octanol M09_octanol*	Comment
1.000	0.76	65.23 %	0.00 %	34.50 %	0.27 %	
1.200	0.76	65.13 %	0.00 %	34.45 %	0.42 %	Stomach pH
2.000	0.79	63.69 %	0.03 %	33.69 %	2.60 %	
3.000	1.00	51.51 %	0.22 %	27.25 %	21.02 %	
4.000	1.68	17.69 %	0.75 %	9.36 %	72.19 %	
5.000	2.50	2.34 %	1.00 %	1.24 %	95.43 %	
6.000	2.92	0.24 %	1.03 %	0.13 %	98.60 %	
6.500	2.98	0.08 %	1.03 %	0.04 %	98.85 %	
7.000	3.00	0.02 %	1.03 %	0.01 %	98.93 %	
7.400	3.01	0.01 %	1.03 %	0.01 %	98.95 %	Blood pH
8.000	3.01	0.00 %	1.03 %	0.00 %	98.96 %	
9.000	3.01	0.00 %	1.03 %	0.00 %	98.97 %	
10.000	3.01	0.00 %	1.03 %	0.00 %	98.97 %	
11.000	3.01	0.00 %	1.03 %	0.00 %	98.97 %	
12.000	3.01	0.00 %	1.03 %	0.00 %	98.97 %	

# Carbonate and acidity



# Other graphs





Assay ID: Filename:

Sample name: M09\_octanol Assay name:

pH-metric high logP

18C-06007


Experiment start time: 3/6/2018 3:40:58 PM

Pion Analyst: Instrument ID:

T312060

C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

## Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name	M09_octanol	2/27/2018 4:56:17 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.001470 g	3/6/2018 3:39:12 PM	User entered value
Formula weight	287.74 g/mol	2/27/2018 4:45:45 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	251.28	2/27/2018 4:45:45 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	1	2/27/2018 4:45:45 PM	User entered value
Sample is a	Base	2/27/2018 4:45:45 PM	User entered value
pKa 1	5.37	2/27/2018 4:45:45 PM	User entered value
logp (XH +)	0.76	3/2/2018 3:27:23 PM	User entered value
logP (neutral X)	3.27	3/2/2018 3:27:17 PM	User entered value
Stoichiometry	1.00000		Default value
Aprotic counterion name	Chloride		From standards.xml file
Stoichiometry	1.00		From standards.xml file
Charge per counterion	-1		From standards.xml file

## **Events**

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
6:16.5	Initial pH = $3.92$									
9:16.0	Data point 1		0.04544 mL							
10:02.2			0.04544 mL						0.00032	
10:38.4	Data point 3		0.04544 mL						0.00018	10.5 s
11:14.3	Data point 4		0.04544 mL						0.00091	10.0 s
11:49.8			0.04544 mL						0.00086	10.0 s
12:25.3	Data point 6		0.04544 mL						0.00041	10.0 s
13:00.8	Data point 7	1.50000 mL	0.04544 mL	0.04452 mL	0.03001 mL	3.357	-0.00824	0.86493	0.00044	10.0 s
	Data point 8		0.04544 mL						0.00053	
14:32.2	Data point 9	1.50000 mL	0.04544 mL	0.04694 mL	0.03001 mL	3.727	-0.00887	0.87089	0.00047	10.5 s
15:18.4	Data point 10		0.04544 mL						0.00067	
15:54.0	Data point 11	1.50000 mL	0.04544 mL	0.04969 mL	0.03001 mL	4.093	-0.01852	0.93335	0.00095	11.0 s
	Data point 12		0.04544 mL						0.00090	
17:21.2	Data point 13		0.04544 mL					0.93747	0.00088	
	Data point 14		0.04544 mL						0.00097	
	Data point 15		0.04544 mL						0.00093	
19:23.3	Data point 16		0.04544 mL						0.00094	
20:08.9	Data point 17		0.04544 mL						0.00097	16.0 s
20:55.1	Data point 18		0.04544 mL						0.00093	
21:43.1	Data point 19	1.50000 mL	0.04544 mL	0.05555 mL	0.03001 mL	5.717	-0.01934	0.94249	0.00098	
22:35.7			0.04544 mL						0.00099	
23:39.3	Data point 21	1.50000 mL	0.04544 mL	0.05602 mL	0.03001 mL	7.443	-0.10565	0.99712	0.00522	Timed out at
										59.5 s
	Data point 22		0.04544 mL						0.00098	
26:50.1			0.04544 mL						0.00100	
28:26.7			0.10228 mL							
29:12.7			0.10228 mL							
29:48.4	Data point 26		0.10228 mL						0.00038	10.5 s
	Data point 27		0.10228 mL							
31:00.1			0.10228 mL						0.00069	
	Data point 29		0.10228 mL						0.00072	10.5 s
32:11.6			0.10228 mL						0.00038	10.0 s
32:57.3	•	1.50000 mL	0.10228 mL	0.10285 mL	0.07001 mL	3.527	-0.00844	0.34234	0.00071	10.0 s
33:48.3	•		0.10228 mL						0.00075	
34:24.3	Data point 33	1.50000 mL	0.10228 mL	0.10595 mL	0.07001 mL	3.910	-0.00623	0.66641	0.00038	10.0 s
34:59.8	Data point 34		0.10228 mL						0.00033	10.5 s
35:51.2	Data point 35	1.50000 mL	0.10228 mL	0.10840 mL	0.07001 mL	4.265	-0.01218	0.79615	0.00067	10.0 s



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# **Events (continued)**

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
36:36.9	Data point 36	1.50000 mL	0.10228 mL	0.10927 mL	0.07001 mL	4.453	-0.01574	0.77656	0.00088	10.0 s
37:12.3	Data point 37		0.10228 mL				-0.01114	0.70830	0.00065	
38:03.7	Data point 38		0.10228 mL						0.00081	
38:45.7	Data point 39		0.10228 mL					0.00984	0.00067	
39:30.9	Data point 40		0.10228 mL					0.93748	0.00099	
40:15.0	Data point 41		0.10228 mL					0.93150	0.00094	
41:02.4	Data point 42		0.10228 mL					0.95517	0.00094	
41:58.0	Data point 43		0.10228 mL						0.00098	
43:16.1	Data point 44		0.10228 mL						0.00414	
40.10.1	Data point 44	1.50000 IIIL	0.10220 IIIL	0.11170 IIIL	0.07001 IIIL	7.000	-0.00000	0.00104	0.00414	at 59.5 s
44:46.6	Data point 45	1.50000 mL	0.10228 mL	0.11183 mL	0.07001 mL	7.214	-0.05929	0.99332	0.00294	
40 47 4	5	4.50000	0.40000	0.44400	0.07004	0.404	0.05500	0.00004	0.00075	at 59.5 s
46:17.1	Data point 46	1.50000 mL	0.10228 mL	0.11192 mL	0.07001 mL	8.104	-0.05533	0.99034	0.00275	Timed out at 59.5 s
47:52.7	Data point 47	1.50000 mL	0.10228 mL	0.11199 mL	0.07001 mL	8.361	-0.02231	0.99013	0.00111	Timed out
					0.0.00		0.0220	0.000.0		at 59.5 s
49:28.4	Data point 48		0.10228 mL					0.96662	0.00099	
50:52.3	Data point 49		0.10228 mL					0.95683	0.00096	
51:60.0	Data point 50		0.10228 mL					0.93622	0.00098	27.5 s
53:28.1	Data point 51		0.16298 mL					0.11888	0.00042	10.0 s
54:14.4	Data point 52	1.50000 mL	0.16298 mL	0.13079 mL	0.17001 mL	2.162	-0.00496	0.56826	0.00032	10.0 s
54:50.1	Data point 53	1.50000 mL	0.16298 mL	0.14332 mL	0.17001 mL	2.380	-0.01003	0.32764	0.00087	10.0 s
55:25.7	Data point 54	1.50000 mL	0.16298 mL	0.15113 mL	0.17001 mL	2.605	-0.00142	0.12787	0.00020	10.0 s
56:01.3	Data point 55	1.50000 mL	0.16298 mL	0.15614 mL	0.17001 mL	2.814	-0.00875	0.26635	0.00084	10.0 s
56:36.8	Data point 56		0.16298 mL					0.61153	0.00020	10.0 s
57:12.3	Data point 57	1.50000 mL	0.16298 mL	0.16232 mL	0.17001 mL	3.248	-0.01113	0.63073	0.00069	10.0 s
57:47.8	Data point 58	1.50000 mL	0.16298 mL	0.16449 mL	0.17001 mL	3.462	-0.00502	0.86633	0.00027	10.0 s
58:23.2	Data point 59	1.50000 mL	0.16298 mL	0.16625 mL	0.17001 mL	3.691	-0.01227	0.89636	0.00064	10.0 s
58:58.6	Data point 60	1.50000 mL	0.16298 mL	0.16761 mL	0.17001 mL	3.895	-0.01553	0.75391	0.00088	25.0 s
59:49.0	Data point 61		0.16298 mL					0.79855	0.00092	
1:00:24.5	Data point 62		0.16298 mL					0.66946	0.00027	
1:01:15.9	Data point 63		0.16298 mL					0.01107	0.00098	10.5 s
1:02:07.3	Data point 64		0.16298 mL					0.81862	0.00098	
1:02:58.6	Data point 65		0.16298 mL						0.00099	
1:03:40.2	•	1.50000 mL		0.17161 mL			-0.01888	0.92060	0.00097	
1:04:26.8	Data point 67	1.50000 mL		0.17180 mL				0.91497	0.00099	
1:05:18.4	Data point 68	1.50000 mL		0.17192 mL				0.89914	0.00099	
	Data point 69		0.16298 mL					0.94473	0.00097	
	Data point 70		0.16298 mL							Timed out
1.01.10.0	Data point 10		0.10200 1112	0.112101112	0.170011112		0.10700	0.00100	0.00000	at 59.5 s
1:09:24.1	Data point 71	1.50000 mL	0.16298 mL	0.17225 mL	0.17001 mL	8.041	-0.05145	0.99163	0.00255	Timed out
4.40.50 =	D-4 170	4.50000	0.40000	0.47000	0.47004	0.450	0.00000	0.05400	0.00446	at 59.5 s
1:10:59.7	Data point 72	1.50000 mL	0.16298 mL	0.17232 mL	0.17001 mL	8.452	-0.02336	0.95422	0.00118	Timed out at 59.5 s
1:12:40.5	Data point 73	1.50000 mL	0.16298 mL	0.17241 mL	0.17001 mL	8.754	-0.01965	0.97546	0.00098	
	Data point 74		0.16298 mL						0.00088	
	Access volumes		0.10200 mL			3.5.5	5.5.520	J.J	2.00000	_0.00

1:14:16.2 Assay volumes 1.50000 mL 0.16298 mL 0.17255 mL 0.17001 mL



Assay name: pH-metric high logP Analyst: Pion Assay ID: Instrument ID: 18C-06007 T312060

Filename: C:\Sirius_T3\Meh	tap\20180306_exp3	0_logP_T3-2\180	C-06007_M09_octano	I_pH-metric high logF
Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	9.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings				
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose	1070			
Titrant pre-dose	None			
Assay Medium	140110			
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.030 mL			
Partition solvent added	Automatic			
After partition addition, stir for	1 seconds			
	i seconus			
Sample Sonication	Voo			
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	120 seconds			
After sonication stir for	5 seconds			
Sample Dissolution	V			
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution	To start pH			
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge				
Perform a carbonate purge	No			
Temperature Control				
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			

Time to wait 60 seconds Stir speed of 50%

Titration 1

Titrate from Low to high pH

Adjust to start pH Yes After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds Stirrer speed for partitioning 50%

Titration 2

Titrate from Low to high pH Add additional water 0.00 mL Additional partition solvent volume 0.040 mL Additional partition solvent added Automatic After pH adjust stir for 30 seconds

Stir to allow partitioning for 15 seconds Stirrer speed for partitioning 55%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# Assay Settings (continued)

Value	<b>Original Value</b>	Date/Time changed	Imported from
	•	_	•
Low to high pH			
0.00 mL			
0.100 mL			
Automatic			
30 seconds			
15 seconds			
60%			
No			
0 seconds			
20 points			
0.50 seconds			
0.00100 dpH/dt			
60 seconds			
	Low to high pH 0.00 mL 0.100 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt	Low to high pH 0.00 mL 0.100 mL Automatic 30 seconds 15 seconds 60%  No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt	Low to high pH 0.00 mL 0.100 mL Automatic 30 seconds 15 seconds 60%  No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt

# Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.124	3/6/2018 3:40:58 PM	C:\Sirius_T3\18C-06006_Blank standardisation.t3r
Four-Plus S	0.9973	3/6/2018 3:40:58 PM	C:\Sirius_T3\18C-06006_Blank standardisation.t3r
Four-Plus jH			C:\Sirius_T3\18C-06006_Blank standardisation.t3r
Four-Plus jOH	-0.7	3/6/2018 3:40:58 PM	C:\Sirius_T3\18C-06006_Blank standardisation.t3r
Base concentration factor	1.000	3/6/2018 3:40:58 PM	C:\Sirius_T3\KOH18B27.t3r
Acid concentration factor	0.989	3/6/2018 3:40:58 PM	C:\Sirius_T3\18C-06006_Blank standardisation.t3r

# Instrument Settings

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module Dispenser 0 Syringe volume Firmware version	Water 2.5 mL 1.2.1(r2)	T3DM1200361	3/31/2009 5:24:52 AM 3/31/2009 5:25:05 AM
Titrant Dispenser 2 Syringe volume Firmware version	Water (0.15 M KCI) Acid 0.5 mL 1.2.1(r2)	02-06-2018	2/27/2018 10:05:59 AM 3/31/2009 5:25:11 AM
Titrant Dispenser 1 Syringe volume Firmware version	Acid (0.5 M HCI) Base 0.5 mL 1.2.1(r2)	02-27-2018	2/27/2018 10:27:22 AM 3/31/2009 5:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version Distribution valve 5	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2) Distribution Valve	9/22/2017	2/27/2018 10:21:22 AM 3/31/2009 5:26:24 AM 3/31/2009 5:28:19 AM
Firmware version Port A Port B Dispenser 3 Syringe volume	1.1.3 Methanol (80%, 0.15 M KCl) Cyclohexane Buffer 0.5 mL	02-08-2018 11-01-17	3/6/2018 9:28:59 AM 2/27/2018 10:37:57 AM 8/3/2010 5:05:16 AM
Firmware version Titrant Dispenser 6	1.2.1(r2) Dodecane Octanol	2018/01/31	2/28/2018 10:18:04 AM 10/22/2010 10:52:43 AM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# Instrument Settings (continued)

Setting Syringe volume	<b>Value</b> 0.5 mL	Batch Id	Install date
Firmware version	1.2.1(r2)	04 24 2040	0/07/0040 0.50.05 AM
Titrant Titrator	Octanol	01-31-2018 T3TM1200161	2/27/2018 9:59:35 AM 3/31/2009 5:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2	1011111200101	6/6 // 2000 G.Z 1. 17 / HVI
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 AI1DI0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1	T05000	4/00/0040 0 04 00 DM
Electrode	T3 Electrode	T3E0923	1/23/2018 2:01:00 PM
E0 calibration Filling solution	+6.73 mV 3M KCI	KCL097	3/6/2018 3:41:43 PM 3/6/2018 9:23:20 AM
Liquids	SW NOI	NOLU91	3/0/2010 9.23.20 AIVI
Wash 1	50% IPA:50% Water		3/6/2018 9:24:32 AM
Wash 2	0.5% Trition X-100 in H20		3/6/2018 9:24:35 AM
Buffer position 1	pH7 Wash		3/6/2018 9:24:38 AM
Buffer position 2	pH 7		3/6/2018 9:24:40 AM
Storage position			3/6/2018 9:24:07 AM
Wash water	6.1e+003 mL	02-27-2018	2/27/2018 9:54:39 AM
Waste	9.4e+003 mL		11/28/2017 10:36:29 AM 8/5/2010 6:35:13 AM
Temperature controller Turbidity detector			3/31/2009 5:24:45 AM
Spectrometer		074811	11/23/2010 11:22:28 AM
Dip probe		10196	11/20/2010 11:22:20 / 111
Wavelength coefficient A0	183.333		
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	123:01:40		11/23/2010 11:22:28 AM
Calibrated on	2/27/2018 10:40:38 AM		
Integration time Scans averaged	40 10		
Autoloader	10	T3AL1200345	11/10/2015 9:34:13 AM
Left-right axis firmware version	1.17 Al1Dl2DO2 Stepper 2	10/12/2000/10	11/10/2010 0:01:10/11/1
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Configuration	Titratian position		
Alternate titration position Alternate reference position	Titration position Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL `´		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration Solvent wash stir speed	5 s 30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration E0 calibration buffer wash stir speed	5 s 30%		
E0 calibration reading stir speed	0%		



Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		
· · · · · · · · · · · · · · · · · · ·			

### Refinement Settings

Setting	Value	Default value
•	value	Delauit Value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00
• •		

## Experiment Log

- [1:59] Air gap released for Acid (0.5 M HCI)
- [2:54] Air gap created for Water (0.15 M KCI)
- [2:54] Air gap created for Acid (0.5 M HCI)
- [2:55] Air gap created for Base (0.5 M KOH)
- [2:55] Air gap released for Water (0.15 M KCI)
- [2:59] Titrator arm moved over Titration position
- [2:59] Titration 1 of 3
- [2:59] Adding initial titrants
- [2:59] Automatically add 1.50000 mL of water
- [3:24] Dispensed 1.500000 mL of Water (0.15 M KCI)
- [3:28] Titrator arm moved over Drain
- [6:10] Titrator arm moved to Titration position
- [6:10] Argon flow rate set to 100 [6:10] Stirrer speed set to 10
- [6:15] Automatically add 0.03000 mL of Octanol
- [6:16] Dispensed 0.030009 mL of Octanol
- [6:17] Initial pH = 3.92
- [6:17] Iterative adjust 3.92 -> 2.00
- [6:17] pH 3.92 -> 2.00
- [6:18] Air gap released for Acid (0.5 M HCI)
- [6:19] Dispensed 0.045437 mL of Acid (0.5 M HCI)
- [6:24] Holding pH 2.00
- [8:24] Stirrer speed set to 0
- [8:24] Stirrer speed set to 50
- [8:24] Iterative adjust 1.94 -> 2.00
- [8:24] pH 1.94 -> 2.00
- [8:25] Air gap released for Base (0.5 M KOH)
- [8:26] Dispensed 0.006091 mL of Base (0.5 M KOH)
- [9:16] Stirrer speed set to 0
- [9:26] Datapoint id 1 collected
- [9:26] Stirrer speed set to 50
- [9:31] pH 2.01 -> 2.21
- [9:31] Using cautious pH adjust
- [9:32] Dispensed 0.007714 mL of Base (0.5 M KOH)
- [9:37] Stepping pH = 2.10
- [9:37] Dispensed 0.006185 mL of Base (0.5 M KOH)
- [9:42] Stepping pH = 2.19
- [9:42] Dispensed 0.000847 mL of Base (0.5 M KOH)
- [9:47] Stepping pH = 2.21
- [10:03] Stirrer speed set to 0
- [10:13] Datapoint id 2 collected
- [10:13] Charge balance equation is out by 4.3%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [10:13] Stirrer speed set to 50
- [10:18] pH 2.21 -> 2.41
- [10:18] Using charge balance adjust
- [10:19] Dispensed 0.009690 mL of Base (0.5 M KOH)
- [10:39] Stirrer speed set to 0
- [10:49] Datapoint id 3 collected
- [10:49] Charge balance equation is out by 1.3%
- [10:49] Stirrer speed set to 50
- [10:54] pH 2.42 -> 2.62
- [10:54] Using charge balance adjust
- [10:54] Dispensed 0.006044 mL of Base (0.5 M KOH)
- [11:15] Stirrer speed set to 0
- [11:25] Datapoint id 4 collected
- [11:25] Charge balance equation is out by 13.5%
- [11:25] Stirrer speed set to 50
- [11:30] pH 2.66 -> 2.86
- [11:30] Using charge balance adjust
- [11:30] Dispensed 0.003669 mL of Base (0.5 M KOH)
- [11:50] Stirrer speed set to 0
- [12:00] Datapoint id 5 collected
- [12:00] Charge balance equation is out by 3.7%
- [12:00] Stirrer speed set to 50
- [12:05] pH 2.87 -> 3.07
- [12:05] Using charge balance adjust
- [12:05] Dispensed 0.002469 mL of Base (0.5 M KOH)
- [12:26] Stirrer speed set to 0
- [12:36] Datapoint id 6 collected
- [12:36] Charge balance equation is out by 10.9%
- [12:36] Stirrer speed set to 50
- [12:41] pH 3.10 -> 3.30
- [12:41] Using charge balance adjust
- [12:41] Dispensed 0.001811 mL of Base (0.5 M KOH)
- [13:01] Stirrer speed set to 0
- [13:11] Datapoint id 7 collected
- [13:11] Charge balance equation is out by 30.3%
- [13:11] Stirrer speed set to 50
- [13:16] pH 3.36 -> 3.56
- [13:16] Using cautious pH adjust
- [13:16] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [13:21] Stepping pH = 3.49
- [13:21] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [13:27] Stepping pH = 3.54
- [13:27] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [13:32] Stepping pH = 3.55
- [13:47] Stirrer speed set to 0
- [13:57] Datapoint id 8 collected
- [13:57] Charge balance equation is out by 19.6%
- [13:57] Stirrer speed set to 50
- [14:02] pH 3.55 -> 3.75
- [14:02] Using cautious pH adjust
- [14:02] Dispensed 0.000729 mL of Base (0.5 M KOH)
- [14:07] Stepping pH = 3.68
- [14:07] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [14:12] Stepping pH = 3.72
- [14:12] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [14:17] Stepping pH = 3.74
- [14:33] Stirrer speed set to 0
- [14:43] Datapoint id 9 collected
- [14:43] Charge balance equation is out by 16.9%
- Reported at: 3/9/2018 11:22:53 AM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [14:43] Stirrer speed set to 50
- [14:48] pH 3.73 -> 3.93
- [14:48] Using cautious pH adjust
- [14:48] Dispensed 0.000706 mL of Base (0.5 M KOH)
- [14:53] Stepping pH = 3.84
- [14:53] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [14:58] Stepping pH = 3.89
- [14:59] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [15:04] Stepping pH = 3.92
- [15:19] Stirrer speed set to 0
- [15:29] Datapoint id 10 collected
- [15:29] Charge balance equation is out by 2.0%
- [15:29] Stirrer speed set to 50
- [15:34] pH 3.92 -> 4.12
- [15:34] Using charge balance adjust
- [15:34] Dispensed 0.001364 mL of Base (0.5 M KOH)
- [15:54] Stirrer speed set to 0
- [16:05] Datapoint id 11 collected
- [16:05] Charge balance equation is out by -13.7%
- [16:05] Stirrer speed set to 50
- [16:10] pH 4.10 -> 4.30
- [16:10] Using charge balance adjust
- [16:11] Dispensed 0.001223 mL of Base (0.5 M KOH)
- [16:31] Stirrer speed set to 0
- [16:41] Datapoint id 12 collected
- [16:41] Charge balance equation is out by -26.8%
- [16:41] Stirrer speed set to 50
- [16:46] pH 4.25 -> 4.45
- [16:46] Using cautious pH adjust
- [16:46] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [16:51] Stepping pH = 4.33
- [16:51] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [16:56] Stepping pH = 4.42
- [16:56] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [17:01] Stepping pH = 4.43
- [17:01] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [17:06] Stepping pH = 4.45
- [17:22] Stirrer speed set to 0
- [17:33] Datapoint id 13 collected
- [17:33] Charge balance equation is out by -31.1%
- [17:33] Stirrer speed set to 50
- [17:38] pH 4.44 -> 4.64
- [17:38] Using cautious pH adjust
- [17:38] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [17:43] Stepping pH = 4.55
- [17:43] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [17:49] Stepping pH = 4.59
- [17:49] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [17:54] Stepping pH = 4.64
- [18:09] Stirrer speed set to 0
- [18:21] Datapoint id 14 collected
- [18:21] Charge balance equation is out by -9.3%
- [18:21] Stirrer speed set to 50
- [18:26] pH 4.64 -> 4.84
- [18:26] Using charge balance adjust
- [18:26] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [18:46] Stirrer speed set to 0
- [18:58] Datapoint id 15 collected
- [18:58] Charge balance equation is out by -14.0%



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [18:58] Stirrer speed set to 50
- [19:03] pH 4.83 -> 5.03
- [19:03] Using charge balance adjust
- [19:03] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [19:24] Stirrer speed set to 0
- [19:39] Datapoint id 16 collected
- [19:39] Charge balance equation is out by -38.2%
- [19:39] Stirrer speed set to 50
- [19:44] pH 4.98 -> 5.18
- [19:44] Using cautious pH adjust
- [19:44] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [19:49] Stepping pH = 5.01
- [19:49] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [19:54] Stepping pH = 5.25
- [20:09] Stirrer speed set to 0
- [20:25] Datapoint id 17 collected
- [20:25] Charge balance equation is out by -62.0%
- [20:25] Stirrer speed set to 50
- [20:30] pH 5.23 -> 5.43
- [20:30] Using cautious pH adjust [20:30] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [20:35] Stepping pH = 5.25
- [20:35] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [20:40] Stepping pH = 5.49
- [20:55] Stirrer speed set to 0
- [21:13] Datapoint id 18 collected
- [21:13] Charge balance equation is out by -85.2%
- [21:13] Stirrer speed set to 50
- [21:18] pH 5.49 -> 5.69
- [21:18] Using cautious pH adjust
- [21:18] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [21:23] Stepping pH = 5.50
- [21:23] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [21:28] Stepping pH = 5.74 [21:43] Stirrer speed set to 0
- [22:00] Datapoint id 19 collected [22:00] Charge balance equation is out by -94.3%
- [22:00] Stirrer speed set to 50
- [22:06] pH 5.76 -> 5.96
- [22:06] Using cautious pH adjust
- [22:06] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [22:11] Stepping pH = 5.77
- [22:11] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [22:16] Stepping pH = 5.85[22:16] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [22:21] Stepping pH = 6.19
- [22:36] Stirrer speed set to 0
- [23:09] Datapoint id 20 collected
- [23:09] Charge balance equation is out by -210.5%
- [23:09] Stirrer speed set to 50
- [23:14] pH 6.23 -> 6.43
- [23:14] Using cautious pH adjust
- [23:14] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [23:19] Stepping pH = 6.23
- [23:19] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [23:25] Stepping pH = 7.24
- [23:40] Stirrer speed set to 0
- [24:40] Datapoint id 21 collected
- [24:40] Charge balance equation is out by -211.2%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [24:40] Stirrer speed set to 50
- [24:45] pH 7.39 -> 7.59
- [24:45] Using cautious pH adjust
- [24:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [24:50] Stepping pH = 7.36
- [24:50] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [24:55] Stepping pH = 7.36
- [24:55] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [25:00] Stepping pH = 8.58
- [25:15] Stirrer speed set to 0
- [26:15] Datapoint id 22 collected
- [26:15] Charge balance equation is out by -1,276.8%
- [26:15] Stirrer speed set to 50
- [26:20] pH 8.72 -> 8.92
- [26:20] Using cautious pH adjust
- [26:20] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [26:25] Stepping pH = 8.71
- [26:25] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [26:30] Stepping pH = 8.89
- [26:30] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [26:35] Stepping pH = 9.02
- [26:50] Stirrer speed set to 0
- [27:28] Datapoint id 23 collected
- [27:28] Charge balance equation is out by -331.4%
- [27:28] Titration 2 of 3
- [27:28] Adding initial titrants
- [27:28] Automatically add 0.04000 mL of Octanol
- [27:29] Dispensed 0.040005 mL of Octanol
- [27:29] Stirrer speed set to 10
- [27:30] Stirrer speed set to 55
- [27:30] Iterative adjust 9.02 -> 2.00
- [27:30] pH 9.02 -> 2.00
- [27:31] Dispensed 0.054351 mL of Acid (0.5 M HCl)
- [27:36] pH 2.02 -> 2.00
- [27:37] Dispensed 0.002493 mL of Acid (0.5 M HCI)
- [28:27] Stirrer speed set to 0
- [28:37] Datapoint id 24 collected
- [28:37] Stirrer speed set to 55
- [28:42] pH 1.97 -> 2.17
- [28:42] Using cautious pH adjust
- [28:42] Dispensed 0.009149 mL of Base (0.5 M KOH)
- [28:47] Stepping pH = 2.07
- [28:48] Dispensed 0.005691 mL of Base (0.5 M KOH)
- [28:53] Stepping pH = 2.14
- [28:53] Dispensed 0.002070 mL of Base (0.5 M KOH)
- [28:58] Stepping pH = 2.17
- [29:13] Stirrer speed set to 0
- [29:23] Datapoint id 25 collected
- [29:23] Charge balance equation is out by 7.6%
- [29:23] Stirrer speed set to 55
- [29:28] pH 2.17 -> 2.37
- [29:28] Using charge balance adjust
- [29:29] Dispensed 0.011430 mL of Base (0.5 M KOH)
- [29:49] Stirrer speed set to 0
- [29:59] Datapoint id 26 collected
- [29:59] Charge balance equation is out by 13.7%
- [29:59] Stirrer speed set to 55
- [30:04] pH 2.41 -> 2.61
- [30:04] Using charge balance adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [30:05] Dispensed 0.006797 mL of Base (0.5 M KOH)
- [30:25] Stirrer speed set to 0
- [30:35] Datapoint id 27 collected
- [30:35] Charge balance equation is out by -1.0%
- [30:35] Stirrer speed set to 55
- [30:40] pH 2.61 -> 2.81
- [30:40] Using charge balance adjust
- [30:40] Dispensed 0.004492 mL of Base (0.5 M KOH)
- [31:00] Stirrer speed set to 0
- [31:10] Datapoint id 28 collected
- [31:10] Charge balance equation is out by 4.7%
- [31:10] Stirrer speed set to 55
- [31:16] pH 2.82 -> 3.02
- [31:16] Using charge balance adjust
- [31:16] Dispensed 0.003081 mL of Base (0.5 M KOH)
- [31:36] Stirrer speed set to 0
- [31:46] Datapoint id 29 collected
- [31:46] Charge balance equation is out by 6.2%
- [31:46] Stirrer speed set to 55
- [31:52] pH 3.04 -> 3.24
- [31:52] Using charge balance adjust
- [31:52] Dispensed 0.002305 mL of Base (0.5 M KOH)
- [32:12] Stirrer speed set to 0
- [32:22] Datapoint id 30 collected
- [32:22] Charge balance equation is out by 46.1%
- [32:22] Stirrer speed set to 55
- [32:27] pH 3.34 -> 3.54
- [32:27] Using cautious pH adjust
- [32:27] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [32:32] Stepping pH = 3.46
- [32:32] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [32:37] Stepping pH = 3.51
- [32:37] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [32:43] Stepping pH = 3.53
- [32:58] Stirrer speed set to 0
  - 2:00] Detancint id 21 cellecter
- [33:08] Datapoint id 31 collected
- [33:08] Charge balance equation is out by 15.3%
- [33:08] Stirrer speed set to 55
- [33:13] pH 3.53 -> 3.73
- [33:13] Using cautious pH adjust
- [33:13] Dispensed 0.000823 mL of Base (0.5 M KOH)
- [33:18] Stepping pH = 3.64
- [33:18] Dispensed 0.000494 mL of Base (0.5 M KOH)
- [33:23] Stepping pH = 3.70
- [33:23] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [33:28] Stepping pH = 3.72
- [33:28] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [33:34] Stepping pH = 3.73
- [33:49] Stirrer speed set to 0
- [33:59] Datapoint id 32 collected
- [33:59] Charge balance equation is out by -0.4%
- [33:59] Stirrer speed set to 55
- [34:04] pH 3.72 -> 3.92
- [34:04] Using charge balance adjust
- [34:04] Dispensed 0.001458 mL of Base (0.5 M KOH)
- [34:25] Stirrer speed set to 0
- [34:35] Datapoint id 33 collected
- [34:35] Charge balance equation is out by -7.5%
- [34:35] Stirrer speed set to 55



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [34:40] pH 3.91 -> 4.11
- [34:40] Using charge balance adjust
- [34:40] Dispensed 0.001223 mL of Base (0.5 M KOH)
- [35:00] Stirrer speed set to 0
- [35:11] Datapoint id 34 collected
- [35:11] Charge balance equation is out by -20.6%
- [35:11] Stirrer speed set to 55
- [35:16] pH 4.08 -> 4.28
- [35:16] Using cautious pH adjust
- [35:16] Dispensed 0.000517 mL of Base (0.5 M KOH)
- [35:21] Stepping pH = 4.16
- [35:21] Dispensed 0.000470 mL of Base (0.5 M KOH)
- [35:26] Stepping pH = 4.24
- [35:26] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [35:31] Stepping pH = 4.27
- [35:31] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [35:36] Stepping pH = 4.27
- [35:51] Stirrer speed set to 0
- [36:02] Datapoint id 35 collected [36:02] Charge balance equation is out by -18.3%
- [36:02] Stirrer speed set to 55
- [36:07] pH 4.27 -> 4.47
- [36:07] Using cautious pH adjust
- [36:07] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [36:12] Stepping pH = 4.35
- [36:12] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [36:17] Stepping pH = 4.42
- [36:17] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [36:22] Stepping pH = 4.46
- [36:37] Stirrer speed set to 0
- [36:47] Datapoint id 36 collected
- [36:47] Charge balance equation is out by -13.4%
- [36:47] Stirrer speed set to 55
- [36:52] pH 4.46 -> 4.66
- [36:52] Using charge balance adjust
- [36:52] Dispensed 0.000564 mL of Base (0.5 M KOH)
- [37:13] Stirrer speed set to 0
- [37:23] Datapoint id 37 collected
- [37:23] Charge balance equation is out by -25.1%
- [37:23] Stirrer speed set to 55
- [37:28] pH 4.62 -> 4.82
- [37:28] Using cautious pH adjust
- [37:28] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [37:33] Stepping pH = 4.67
- [37:34] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [37:39] Stepping pH = 4.79
- [37:39] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [37:44] Stepping pH = 4.80
- [37:44] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [37:49] Stepping pH = 4.86
- [38:04] Stirrer speed set to 0
- [38:16] Datapoint id 38 collected
- [38:16] Charge balance equation is out by -66.8%
- [38:16] Stirrer speed set to 55
- [38:21] pH 4.87 -> 5.07
- [38:21] Using cautious pH adjust
- [38:21] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [38:26] Stepping pH = 4.90
- [38:26] Dispensed 0.000259 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [38:31] Stepping pH = 5.07
- [38:46] Stirrer speed set to 0
- [39:01] Datapoint id 39 collected
- [39:01] Charge balance equation is out by -50.1%
- [39:01] Stirrer speed set to 55
- [39:06] pH 5.08 -> 5.28
- [39:06] Using cautious pH adjust
- [39:06] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [39:11] Stepping pH = 5.10
- [39:11] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [39:16] Stepping pH = 5.32
- [39:31] Stirrer speed set to 0
- [39:45] Datapoint id 40 collected
- [39:45] Charge balance equation is out by -73.8%
- [39:45] Stirrer speed set to 55
- [39:50] pH 5.34 -> 5.54
- [39:50] Using cautious pH adjust
- [39:50] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [39:55] Stepping pH = 5.34
- [39:55] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [40:00] Stepping pH = 5.59
- [40:15] Stirrer speed set to 0
- [40:32] Datapoint id 41 collected
- [40:32] Charge balance equation is out by -98.5%
- [40:32] Stirrer speed set to 55
- [40:37] pH 5.61 -> 5.81
- [40:37] Using cautious pH adjust
- [40:37] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [40:43] Stepping pH = 5.62
- [40:43] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [40:48] Stepping pH = 5.97
- [41:03] Stirrer speed set to 0
- [41:28] Datapoint id 42 collected
- [41:28] Charge balance equation is out by -92.4%
- [41:28] Stirrer speed set to 55
- [41:33] pH 5.99 -> 6.19
- [41:33] Using cautious pH adjust
- [41:33] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [41:38] Stepping pH = 6.00
- [41:38] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [41:43] Stepping pH = 6.30
- [41:58] Stirrer speed set to 0
- [42:46] Datapoint id 43 collected
- [42:46] Charge balance equation is out by -90.5%
- [42:46] Stirrer speed set to 55
- [42:51] pH 6.33 -> 6.53
- [42:51] Using cautious pH adjust
- [42:51] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [42:56] Stepping pH = 6.35
- [42:56] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [43:01] Stepping pH = 6.95
- [43:16] Stirrer speed set to 0
- [44:16] Datapoint id 44 collected
- [44:16] Charge balance equation is out by -77.9%
- [44:16] Stirrer speed set to 55
- [44:22] pH 7.16 -> 7.36
- [44:22] Using cautious pH adjust
- [44:22] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [44:27] Stepping pH = 7.21



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [44:27] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [44:32] Stepping pH = 7.36
- [44:47] Stirrer speed set to 0
- [45:47] Datapoint id 45 collected
- [45:47] Charge balance equation is out by -160.1%
- [45:47] Stirrer speed set to 55
- [45:52] pH 7.07 -> 7.27
- [45:52] Using cautious pH adjust
- [45:52] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [45:57] Stepping pH = 7.02
- [45:57] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [46:02] Stepping pH = 7.98
- [46:17] Stirrer speed set to 0
- [47:17] Datapoint id 46 collected
- [47:17] Charge balance equation is out by -347.1%
- [47:17] Stirrer speed set to 55
- [47:23] pH 8.04 -> 8.24
- [47:23] Using cautious pH adjust
- [47:23] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [47:28] Stepping pH = 8.03
- [47:28] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [47:33] Stepping pH = 8.16
- [47:33] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [47:38] Stepping pH = 8.35
- [47:53] Stirrer speed set to 0
- [48:53] Datapoint id 47 collected
- [48:53] Charge balance equation is out by -699.5%
- [48:53] Stirrer speed set to 55
- [48:58] pH 8.43 -> 8.63
- [48:58] Using cautious pH adjust
- [48:58] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [49:03] Stepping pH = 8.46
- [49:03] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [49:08] Stepping pH = 8.57
- [49:09] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [49:14] Stepping pH = 8.69
- [49:29] Stirrer speed set to 0
- [50:12] Datapoint id 48 collected
- [50:12] Charge balance equation is out by -340.0%
- [50:12] Stirrer speed set to 55
- [50:17] pH 8.69 -> 8.89
- [50:17] Using cautious pH adjust
- [50:17] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [50:22] Stepping pH = 8.69
- [50:22] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [50:27] Stepping pH = 8.77
- [50:27] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [50:32] Stepping pH = 8.86
- [50:32] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [50:38] Stepping pH = 8.93
- [50:53] Stirrer speed set to 0
- [51:25] Datapoint id 49 collected
- [51:25] Charge balance equation is out by -371.7%
- [51:25] Stirrer speed set to 55
- [51:30] pH 8.95 -> 9.05
- [51:30] Using cautious pH adjust
- [51:30] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [51:35] Stepping pH = 8.95
- [51:35] Dispensed 0.000024 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [51:40] Stepping pH = 8.96
- [51:40] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [51:45] Stepping pH = 9.13
- [52:00] Stirrer speed set to 0
- [52:28] Datapoint id 50 collected
- [52:28] Charge balance equation is out by -530.3%
- [52:28] Titration 3 of 3
- [52:28] Adding initial titrants
- [52:28] Automatically add 0.10000 mL of Octanol
- [52:30] Dispensed 0.100000 mL of Octanol
- [52:30] Stirrer speed set to 10
- [52:31] Stirrer speed set to 60
- [52:31] Iterative adjust 9.13 -> 2.00
- [52:31] pH 9.13 -> 2.00
- [52:33] Dispensed 0.057314 mL of Acid (0.5 M HCI)
- [52:38] pH 2.03 -> 2.00
- [52:38] Dispensed 0.003387 mL of Acid (0.5 M HCl)
- [53:28] Stirrer speed set to 0
- [53:39] Datapoint id 51 collected
- [53:39] Stirrer speed set to 60
- [53:44] pH 1.97 -> 2.17
- [53:44] Using cautious pH adjust
- [53:44] Dispensed 0.009901 mL of Base (0.5 M KOH)
- [53:49] Stepping pH = 2.06
- [53:49] Dispensed 0.006914 mL of Base (0.5 M KOH)
- [53:54] Stepping pH = 2.14
- [53:55] Dispensed 0.001693 mL of Base (0.5 M KOH)
- [54:00] Stepping pH = 2.17
- [54:15] Stirrer speed set to 0
- [54:25] Datapoint id 52 collected
- [54:25] Charge balance equation is out by 6.5%
- [54:25] Stirrer speed set to 60
- [54:30] pH 2.17 -> 2.37
- [54:30] Using charge balance adjust
- [54:30] Dispensed 0.012535 mL of Base (0.5 M KOH)
- [54:50] Stirrer speed set to 0
- [55:00] Datapoint id 53 collected
- [55:00] Charge balance equation is out by 6.7%
- [55:00] Stirrer speed set to 60
- [55:06] pH 2.38 -> 2.58
- [55:06] Using charge balance adjust
- [55:06] Dispensed 0.007808 mL of Base (0.5 M KOH)
- [55:26] Stirrer speed set to 0
- [55:36] Datapoint id 54 collected
- [55:36] Charge balance equation is out by 10.7%
- [55:36] Stirrer speed set to 60
- [55:41] pH 2.61 -> 2.81
- [55:41] Using charge balance adjust
- [55:41] Dispensed 0.005009 mL of Base (0.5 M KOH)
- [56:02] Stirrer speed set to 0
- [56:12] Datapoint id 55 collected
- [56:12] Charge balance equation is out by 2.5%
- [56:12] Stirrer speed set to 60
- [56:17] pH 2.82 -> 3.02
- [56:17] Using charge balance adjust
- [56:17] Dispensed 0.003551 mL of Base (0.5 M KOH)
- [56:37] Stirrer speed set to 0
- [56:47] Datapoint id 56 collected
- [56:47] Charge balance equation is out by 13.2%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [56:47] Stirrer speed set to 60
- [56:52] pH 3.05 -> 3.25
- [56:52] Using charge balance adjust
- [56:52] Dispensed 0.002634 mL of Base (0.5 M KOH)
- [57:13] Stirrer speed set to 0
- [57:23] Datapoint id 57 collected
- [57:23] Charge balance equation is out by 1.0%
- [57:23] Stirrer speed set to 60
- [57:28] pH 3.25 -> 3.45
- [57:28] Using charge balance adjust
- [57:28] Dispensed 0.002164 mL of Base (0.5 M KOH)
- [57:48] Stirrer speed set to 0
- [57:58] Datapoint id 58 collected
- [57:58] Charge balance equation is out by 5.5%
- [57:58] Stirrer speed set to 60
- [58:03] pH 3.46 -> 3.66
- [58:03] Using charge balance adjust
- [58:03] Dispensed 0.001764 mL of Base (0.5 M KOH)
- [58:24] Stirrer speed set to 0
- 58:34 Datapoint id 59 collected
- [58:34] Charge balance equation is out by 13.2%
- [58:34] Stirrer speed set to 60
- [58:39] pH 3.69 -> 3.89
- [58:39] Using charge balance adjust
- [58:39] Dispensed 0.001364 mL of Base (0.5 M KOH)
- [58:59] Stirrer speed set to 0
- [59:24] Datapoint id 60 collected
- [59:24] Charge balance equation is out by 0.3%
- [59:24] Stirrer speed set to 60
- [59:29] pH 3.90 -> 4.10
- [59:29] Using charge balance adjust
- [59:29] Dispensed 0.001011 mL of Base (0.5 M KOH)
- [59:49] Stirrer speed set to 0
- [59:59] Datapoint id 61 collected
- [59:59] Charge balance equation is out by -14.0%
- [59:59] Stirrer speed set to 60
- [1:00:05] pH 4.08 -> 4.28
- [1:00:05] Using charge balance adjust
- [1:00:05] Dispensed 0.000776 mL of Base (0.5 M KOH)
- [1:00:25] Stirrer speed set to 0
- [1:00:35] Datapoint id 62 collected
- [1:00:35] Charge balance equation is out by -23.3%
- [1:00:35] Stirrer speed set to 60
- [1:00:40] pH 4.24 -> 4.44
- [1:00:40] Using cautious pH adjust
- [1:00:41] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [1:00:46] Stepping pH = 4.28
- [1:00:46] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [1:00:51] Stepping pH = 4.41
- [1:00:51] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:00:56] Stepping pH = 4.42
- [1:00:56] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [1:01:01] Stepping pH = 4.46
- [1:01:16] Stirrer speed set to 0
- [1:01:27] Datapoint id 63 collected
- [1:01:27] Charge balance equation is out by -64.5%
- [1:01:27] Stirrer speed set to 60
- [1:01:32] pH 4.46 -> 4.66
- [1:01:32] Using cautious pH adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [1:01:32] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:01:37] Stepping pH = 4.52
- [1:01:37] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [1:01:42] Stepping pH = 4.62
- [1:01:42] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:01:47] Stepping pH = 4.63
- [1:01:47] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:01:53] Stepping pH = 4.66
- [1:02:08] Stirrer speed set to 0
- [1:02:18] Datapoint id 64 collected
- [1:02:18] Charge balance equation is out by -59.0%
- [1:02:18] Stirrer speed set to 60
- [1:02:23] pH 4.67 -> 4.87
- [1:02:23] Using cautious pH adjust
- [1:02:23] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [1:02:28] Stepping pH = 4.71
- [1:02:28] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:02:34] Stepping pH = 4.84
- [1:02:34] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:02:39] Stepping pH = 4.85
- [1:02:39] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:02:44] Stepping pH = 4.86
- [1:02:59] Stirrer speed set to 0
  - 1:02:00] Othrer speed set to o
- [1:03:10] Datapoint id 65 collected
- [1:03:10] Charge balance equation is out by -59.7%
- [1:03:10] Stirrer speed set to 60
- [1:03:15] pH 4.87 -> 5.07
- [1:03:15] Using cautious pH adjust
- [1:03:15] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:03:20] Stepping pH = 4.88
- [1:03:20] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [1:03:25] Stepping pH = 5.12
- [1:03:41] Stirrer speed set to 0
- [1:03:57] Datapoint id 66 collected
- [1:03:57] Charge balance equation is out by -86.7%
- [1:03:57] Stirrer speed set to 60
- [1:04:02] pH 5.14 -> 5.34
- [1:04:02] Using cautious pH adjust
- [1:04:02] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:04:07] Stepping pH = 5.15
- [1:04:07] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:04:12] Stepping pH = 5.43
- [1:04:27] Stirrer speed set to 0
- [1:04:48] Datapoint id 67 collected
- [1:04:48] Charge balance equation is out by -88.8%
- [1:04:48] Stirrer speed set to 60
- [1:04:53] pH 5.48 -> 5.68
- [1:04:53] Using cautious pH adjust
- [1:04:53] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:04:59] Stepping pH = 5.49
- [1:04:59] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:05:04] Stepping pH = 5.73
- [1:05:19] Stirrer speed set to 0
- [1:05:57] Datapoint id 68 collected
- [1:05:57] Charge balance equation is out by -92.7%
- [1:05:57] Stirrer speed set to 60
- [1:06:02] pH 5.81 -> 6.01
- [1:06:02] Using cautious pH adjust
- [1:06:02] Dispensed 0.000024 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

- [1:06:07] Stepping pH = 5.82
- [1:06:08] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:06:13] Stepping pH = 6.20
- [1:06:28] Stirrer speed set to 0
- [1:07:18] Datapoint id 69 collected
- [1:07:18] Charge balance equation is out by -91.3%
- [1:07:18] Stirrer speed set to 60
- [1:07:23] pH 6.15 -> 6.35
- [1:07:23] Using cautious pH adjust
- [1:07:23] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:07:29] Stepping pH = 6.11
- [1:07:29] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:07:34] Stepping pH = 7.70
- [1:07:49] Stirrer speed set to 0
- [1:08:49] Datapoint id 70 collected
- [1:08:49] Charge balance equation is out by -243.1%
- [1:08:49] Stirrer speed set to 60
- [1:08:54] pH 7.85 -> 8.05
- [1:08:54] Using cautious pH adjust
- [1:08:54] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:08:59] Stepping pH = 7.89
- [1:08:59] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:09:04] Stepping pH = 7.95
- [1:09:04] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:09:09] Stepping pH = 8.07
- [1:09:24] Stirrer speed set to 0
- [1:10:24] Datapoint id 71 collected
- [1:10:24] Charge balance equation is out by -683.0%
- [1:10:24] Stirrer speed set to 60
- [1:10:30] pH 8.14 -> 8.34
- [1:10:30] Using cautious pH adjust
- [1:10:30] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:10:35] Stepping pH = 8.19
- [1:10:35] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:10:40] Stepping pH = 8.32
- [1:10:40] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:10:45] Stepping pH = 8.48
- [1:11:00] Stirrer speed set to 0
- [1:12:00] Datapoint id 72 collected
- [1:12:00] Charge balance equation is out by -550.4%
- [1:12:00] Stirrer speed set to 60
- [1:12:05] pH 8.50 -> 8.70
- [1:12:05] Using cautious pH adjust
- [1:12:05] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:12:10] Stepping pH = 8.51
- [1:12:10] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:12:15] Stepping pH = 8.60
- [1:12:15] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:12:21] Stepping pH = 8.68
- [1:12:21] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:12:26] Stepping pH = 8.77
- [1:12:41] Stirrer speed set to 0
- [1:13:09] Datapoint id 73 collected
- [1:13:09] Charge balance equation is out by -398.4%
- [1:13:09] Stirrer speed set to 60
- [1:13:14] pH 8.76 -> 8.96
- [1:13:14] Using cautious pH adjust
- [1:13:14] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:13:19] Stepping pH = 8.75

## **Experiment Log**



Sample name: M09\_octanol Experiment start time: 3/6/2018 3:40:58 PM

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-06007 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06007\_M09\_octanol\_pH-metric high logP.t3r

# **Experiment Log (continued)**

[1:13:19] Dispensed 0.000071 mL of Base (0.5 M KOH)

[1:13:24] Stepping pH = 8.87

[1:13:24] Dispensed 0.000047 mL of Base (0.5 M KOH)

[1:13:29] Stepping pH = 9.01

[1:13:45] Stirrer speed set to 0

[1:14:08] Datapoint id 74 collected

[1:14:08] Charge balance equation is out by -368.3%

[1:14:08] Argon flow rate set to 0

[1:14:11] Titrator arm moved over Titration position