

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH2 2+) -5.50

logP (XH +) 0.28 ±0.16 (n=50) logP (neutral X) 2.04 ±0.02 (n=50)

18C-03006 Points 1 to 32

M15_octanol concentration factor 0.894
Carbonate 0.2191 mM
Acidity error -0.25037 mM

18C-03006 Points 33 to 65

M15_octanol concentration factor 0.950
Carbonate 0.2177 mM
Acidity error -0.29629 mM

18C-03006 Points 66 to 99

M15_octanol concentration factor 1.015
Carbonate 0.2265 mM
Acidity error -0.47876 mM

Warnings and errors

Errors None

Warnings One or more logP values out of range

Sample logD and percent species

рН	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	Comment
	logD	M15_octanolH2	M15_octanolH	M15_octanol	M15_octanolH2*	M15_octanolH*	M15_octanol*	
1.000	-1.31	92.86 %	2.44 %	0.00 %	0.00 %	4.68 %	0.01 %	
1.200	-1.11	89.13 %	3.72 %	0.00 %	0.00 %	7.12 %	0.03 %	Stomach pH
2.000	-0.39	56.12 %	14.76 %	0.01 %	0.00 %	28.29 %	0.82 %	•
3.000	0.25	9.86 %	25.93 %	0.13 %	0.00 %	49.68 %	14.41 %	
4.000	0.84	0.44 %	11.68 %	0.59 %	0.00 %	22.38 %	64.91 %	
5.000	1.58	0.01 %	1.70 %	0.85 %	0.00 %	3.25 %	94.20 %	
6.000	1.97	0.00 %	0.18 %	0.89 %	0.00 %	0.34 %	98.59 %	
6.500	2.02	0.00 %	0.06 %	0.89 %	0.00 %	0.11 %	98.94 %	
7.000	2.04	0.00 %	0.02 %	0.89 %	0.00 %	0.03 %	99.05 %	
7.400	2.04	0.00 %	0.01 %	0.89 %	0.00 %	0.01 %	99.09 %	Blood pH
8.000	2.04	0.00 %	0.00 %	0.89 %	0.00 %	0.00 %	99.10 %	
9.000	2.04	0.00 %	0.00 %	0.89 %	0.00 %	0.00 %	99.11 %	l
10.000	2.04	0.00 %	0.00 %	0.89 %	0.00 %	0.00 %	99.11 %	ļ
11.000	2.04	0.00 %	0.00 %	0.89 %	0.00 %	0.00 %	99.11 %	
12.000	2.04	0.00 %	0.00 %	0.89 %	0.00 %	0.00 %	99.11 %	



Sample name: M15_octanol Assay name:

pH-metric high logP

18C-03006


Experiment start time: 3/3/2018 8:26:22 AM

Pion Analyst: Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Graphs

Assay ID: Filename:






Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-03006 Points 1 to 32

Overall results

RMSD 0.221
Average ionic strength 0.158 M
Average temperature 25.0°C
Partition ratio 0.0207 : 1

Analyte concentration range 3659.4 µM to 3810.2 µM

Total points considered 23 of 32

Warnings and errors

Errors None

Warnings One or more logP values out of range

Four-Plus parameters


Alpha	0.111	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r
S	0.9988	3/3/2018 8:26:22 AM	C:\Sirius T3\HCl18C02.t3r
jΗ	1.0	3/3/2018 8:26:22 AM	C:\Sirius T3\HCl18C02.t3r
jОН	-0.8	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r


Titrants

Sample

₩	M15_octanol concentration factor	0.894
	Base pKa 1	2.58
	Base pKa 2	5.30
	logP (XH2 2+)	-5.50
₩	logP (XH +)	-2.63
₩	logP (neutral X)	1.84

Sample graphs








Assay ID:

Sample name: M15_octanol Experiment start time: 3/3/2018 8:26:22 AM
Assay name: pH-metric high logP Analyst: Pion

pH-metric high logP Analyst: Pion 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M15_octanol logD	M15_octanol M15_octanolH2	M15_octanol M15_octanolH	M15_octanol M15_octanol		M15_octanol M15_octanolH*	M15_octanol M15_octanol*	
1.000	-3.82	97.44 %	2.56 %	0.00 %	0.00 %	0.00 %	0.00 %	
1.200	-3.50	96.00 %	4.00 %	0.00 %	0.00 %	0.00 %	0.00 %	Stomach pH
2.000	-2.11	79.15 %	20.82 %	0.01 %	0.00 %	0.00 %	0.01 %	•
3.000	-0.60	27.30 %	71.82 %	0.36 %	0.00 %	0.00 %	0.52 %	
4.000	0.50	3.28 %	86.21 %	4.32 %	0.00 %	0.00 %	6.18 %	
5.000	1.36	0.17 %	45.00 %	22.55 %	0.00 %	0.00 %	32.28 %	
6.000	1.76	0.00 %	7.58 %	38.01 %	0.00 %	0.00 %	54.40 %	
6.500	1.81	0.00 %	2.53 %	40.09 %	0.00 %	0.00 %	57.38 %	
7.000	1.83	0.00 %	0.81 %	40.80 %	0.00 %	0.00 %	58.39 %	
7.400	1.84	0.00 %	0.33 %	41.00 %	0.00 %	0.00 %	58.68 %	Blood pH
8.000	1.84	0.00 %	0.08 %	41.10 %	0.00 %	0.00 %	58.82 %	
9.000	1.84	0.00 %	0.01 %	41.13 %	0.00 %	0.00 %	58.86 %	
10.000	1.84	0.00 %	0.00 %	41.13 %	0.00 %	0.00 %	58.87 %	
11.000	1.84	0.00 %	0.00 %	41.13 %	0.00 %	0.00 %	58.87 %	
12.000	1.84	0.00 %	0.00 %	41.13 %	0.00 %	0.00 %	58.87 %	

Carbonate and acidity



Other graphs






Pion pH-metric high logP Assay name: Analyst: 18C-03006 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-03006 Points 33 to 65

Overall results

RMSD 0.481
Average ionic strength 0.167 M
Average temperature 25.0°C
Partition ratio 0.0576 : 1

Analyte concentration range 3242.7 µM to 3373.5 µM

Total points considered 31 of 33

Warnings and errors

Errors None

Warnings One or more logP values out of range

Four-Plus parameters


Alpha	0.111	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r
S	0.9988	3/3/2018 8:26:22 AM	C:\Sirius T3\HCl18C02.t3r
jΗ	1.0	3/3/2018 8:26:22 AM	C:\Sirius T3\HCl18C02.t3r
jОН	-0.8	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r


Titrants

Sample

Φ	M15_octanol concentration factor	0.950
	Base pKa 1	2.58
	Base pKa 2	5.30
	logP (XH2 2+)	-5.50
₩	logP (XH +)	-2.46
₩	logP (neutral X)	1.99

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	
	logD	M15_octanolH2	M15_octanolH	M15_octanol	M15_octanolH2*	M15_octanolH*	M15_octanol*	
1.000	-3.66	97.44 %	2.56 %	0.00 %	0.00 %	0.00 %	0.00 %	
1.200	-3.35	96.00 %	4.00 %	0.00 %	0.00 %	0.00 %	0.00 %	Stomach pH
2.000	-1.97	79.12 %	20.81 %	0.01 %	0.00 %	0.00 %	0.06 %	•
3.000	-0.45	26.90 %	70.76 %	0.35 %	0.00 %	0.01 %	1.97 %	
4.000	0.65	2.78 %	73.13 %	3.67 %	0.00 %	0.01 %	20.41 %	
5.000	1.51	0.09 %	23.28 %	11.67 %	0.00 %	0.00 %	64.96 %	
6.000	1.91	0.00 %	2.95 %	14.77 %	0.00 %	0.00 %	82.28 %	
6.500	1.96	0.00 %	0.95 %	15.08 %	0.00 %	0.00 %	83.97 %	
7.000	1.98	0.00 %	0.30 %	15.18 %	0.00 %	0.00 %	84.52 %	
7.400	1.98	0.00 %	0.12 %	15.21 %	0.00 %	0.00 %	84.67 %	Blood pH
8.000	1.98	0.00 %	0.03 %	15.22 %	0.00 %	0.00 %	84.75 %	
9.000	1.99	0.00 %	0.00 %	15.22 %	0.00 %	0.00 %	84.77 %	
10.000	1.99	0.00 %	0.00 %	15.22 %	0.00 %	0.00 %	84.78 %	
11.000	1.99	0.00 %	0.00 %	15.22 %	0.00 %	0.00 %	84.78 %	
12.000	1.99	0.00 %	0.00 %	15.22 %	0.00 %	0.00 %	84.78 %	

Carbonate and acidity



Other graphs






Pion pH-metric high logP Assay name: Analyst: 18C-03006 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-03006 Points 66 to 99

Overall results

RMSD 0.728
Average ionic strength 0.173 M
Average temperature 25.0°C
Partition ratio 0.1596 : 1

Analyte concentration range 2725.9 µM to 2821.3 µM

Total points considered 33 of 34

Warnings and errors

Errors None

Warnings One or more logP values out of range

Four-Plus parameters


Alpha	0.111	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r
S	0.9988	3/3/2018 8:26:22 AM	C:\Sirius T3\HCl18C02.t3r
jΗ	1.0	3/3/2018 8:26:22 AM	C:\Sirius T3\HCl18C02.t3r
jОН	-0.8	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r


Titrants

Sample

Φ	M15_octanol concentration factor	1.015
	Base pKa 1	2.58
	Base pKa 2	5.30
⊕	logP (XH2 2+)	-5.50
₩	logP (XH +)	0.54
₩	logP (neutral X)	2.20

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	M15_octanol	
	logD	M15_octanolH2	M15_octanolH	M15_octanol	M15_octanolH2*	M15_octanolH*	M15_octanol*	
1.000	-1.05	96.09 %	2.53 %	0.00 %	0.00 %	1.38 %	0.00 %	
1.200	-0.86	93.93 %	3.92 %	0.00 %	0.00 %	2.14 %	0.01 %	Stomach pH
2.000	-0.14	70.90 %	18.65 %	0.01 %	0.00 %	10.21 %	0.23 %	•
3.000	0.48	18.47 %	48.58 %	0.24 %	0.00 %	26.59 %	6.11 %	
4.000	1.02	1.31 %	34.56 %	1.73 %	0.00 %	18.91 %	43.49 %	
5.000	1.74	0.03 %	6.83 %	3.42 %	0.00 %	3.74 %	85.98 %	
6.000	2.12	0.00 %	0.76 %	3.79 %	0.00 %	0.41 %	95.05 %	
6.500	2.17	0.00 %	0.24 %	3.82 %	0.00 %	0.13 %	95.81 %	
7.000	2.19	0.00 %	0.08 %	3.83 %	0.00 %	0.04 %	96.06 %	
7.400	2.19	0.00 %	0.03 %	3.83 %	0.00 %	0.02 %	96.12 %	Blood pH
8.000	2.20	0.00 %	0.01 %	3.83 %	0.00 %	0.00 %	96.16 %	
9.000	2.20	0.00 %	0.00 %	3.83 %	0.00 %	0.00 %	96.17 %	
10.000	2.20	0.00 %	0.00 %	3.83 %	0.00 %	0.00 %	96.17 %	
11.000	2.20	0.00 %	0.00 %	3.83 %	0.00 %	0.00 %	96.17 %	
12.000	2.20	0.00 %	0.00 %	3.83 %	0.00 %	0.00 %	96.17 %	

Carbonate and acidity



Other graphs





Assay name:

Filename:

Sample name: M15_octanol

pH-metric high logP

18C-03006 Assay ID:


Experiment start time: 3/3/2018 8:26:22 AM

Pion Analyst: Instrument ID:

T312060

C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name		2/27/2018 5:03:03 PM	
Sample by	Weight	0/0/00/00 = 00 00 00	Default value
Sample weight		3/2/2018 5:09:20 PM	
Formula weight	209.25 g/mol	2/27/2018 5:03:03 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	209.25	2/27/2018 5:03:03 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	2	2/27/2018 5:03:03 PM	User entered value
Sample is a	Base	2/27/2018 5:03:03 PM	User entered value
pKa 1	2.58	2/27/2018 5:03:03 PM	User entered value
pKa 2	5.30	2/27/2018 5:03:03 PM	User entered value
logP (XH2 2+)	-5.50	2/28/2018 2:10:35 PM	User entered value
logp (XH +)	-4.96	2/28/2018 2:10:28 PM	User entered value
logP (neutral X)	1.92	2/28/2018 2:10:15 PM	User entered value

Events

Nitial pH = 6.94 12:16.5 Data point 1 2.20002 mL 0.09546 mL 0.00000 mL 0.05000 mL 2.195 0.00038 0.01682 0.00030 10.0 s 13:04.2 Data point 2 2.20002 mL 0.09546 mL 0.02439 mL 0.05000 mL 2.195 0.000638 0.59345 0.00041 10.0 s 13:40.0 Data point 3 2.20002 mL 0.09546 mL 0.05160 mL 0.05000 mL 2.395 0.00048 0.03777 0.00012 10.0 s 14:55.8 Data point 4 2.20002 mL 0.09546 mL 0.05160 mL 0.05000 mL 2.392 0.00030 0.05330 0.00098 10.0 s 15:31.3 Data point 5 2.20002 mL 0.09546 mL 0.05506 mL 0.05000 mL 2.777 0.00803 0.16330 0.00018 10.0 s 16:06.8 Data point 7 2.20002 mL 0.09546 mL 0.06928 mL 0.05000 mL 2.394 0.00078 0.00078 0.00078 10.0 s 16:06.8 Data point 7 2.20002 mL 0.09546 mL 0.07272 mL 0.05000 mL 3.157 0.01018 0.35970 0.00089 11.5 s 17:30.2 Data point 9 2.20002 mL 0.09546 mL 0.07272 mL 0.05000 mL 3.574 0.00106 0.00051 10.0 s 18:05.7 Data point 10 2.20002 mL 0.09546 mL 0.07658 mL 0.05000 mL 3.742 0.00016 0.00051 10.0 s 18:05.7 Data point 10 2.20002 mL 0.09546 mL 0.07658 mL 0.05000 mL 3.742 0.00016 0.00051 10.0 s 19:32.3 Data point 12 2.20002 mL 0.09546 mL 0.08281 mL 0.05000 mL 3.742 0.00015 0.02012 0.00052 10.0 s 10.	3	•									
Section Sect	Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	•
13:04.2 Data point 2 2,20002 mL 0,09546 mL 0,02439 mL 0,05000 mL 2,195 0,00638 0,59345 0,00041 10.0 s 13:40.0 Data point 3 2,20002 mL 0,09546 mL 0,04048 mL 0,05000 mL 2,388 0,01047 0,48508 0,00100 14.5 s 14:20.2 Data point 4 2,20002 mL 0,09546 mL 0,05060 mL 0,05000 mL 2,582 0,00048 0,03777 0,00012 10.0 s 16:56.8 Data point 5 2,20002 mL 0,09546 mL 0,06526 mL 0,05000 mL 2,777 0,00080 1,05303 0,00010 10.0 s 16:06.8 Data point 7 2,20002 mL 0,09546 mL 0,06526 mL 0,05000 mL 3,777 0,00080 1,03803 0,00010 10.0 s 16:06.8 Data point 7 2,20002 mL 0,09546 mL 0,06526 mL 0,05000 mL 3,777 0,00080 1,03803 0,00010 10.0 s 16:06.8 Data point 7 2,20002 mL 0,09546 mL 0,07272 mL 0,05000 mL 3,777 0,00080 1,03803 0,00010 10.0 s 17:30 Data point 9 2,20002 mL 0,09546 mL 0,07272 mL 0,05000 mL 3,774 0,00016 0,01066 0,00051 10.0 s 18:05.7 Data point 10 2,20002 mL 0,09546 mL 0,07272 mL 0,05000 mL 3,744 0,00046 0,00051 10.0 s 18:05.7 Data point 11 2,20002 mL 0,09546 mL 0,07843 mL 0,05000 mL 3,744 0,00048 0,00212 0,00052 10.0 s 18:05.8 Data point 12 2,20002 mL 0,09546 mL 0,07843 mL 0,05000 mL 3,744 0,00048 0,00212 0,00052 10.0 s 18:05.8 Data point 12 2,20002 mL 0,09546 mL 0,07843 mL 0,05000 mL 4,715 0,00048 0,00212 0,00052 10.0 s 19:05.8 Data point 12 2,20002 mL 0,09546 mL 0,08848 mL 0,05000 mL 4,715 0,00048 0,00212 0,00052 10.0 s 19:05.2 Data point 12 2,20002 mL 0,09546 mL 0,08848 mL 0,05000 mL 5,749 0,00048 0,00021 0,00052 10.0 s	9:11.9	Initial pH = 6.94									
13:40.0 Data point 3 2.20002 mL 0.09546 mL 0.04048 mL 0.05000 mL 2.388 0.01407 0.48508 0.00100 14.5 s	12:16.5	Data point 1	2.20002 mL	0.09546 mL	0.00000 mL	0.05000 mL	1.982	0.00078	0.01682	0.00030	10.0 s
14:20.2 Data point 4 2.20002 mL 0.09546 mL 0.05160 mL 0.05000 mL 2.582 -0.00048 0.03777 0.00012 10.0 s 15:31.3 Data point 5 2.20002 mL 0.09546 mL 0.050526 mL 0.05000 mL 2.994 -0.00078 0.13803 0.00010 10.0 s 15:31.3 Data point 6 2.20002 mL 0.09546 mL 0.050526 mL 0.05000 mL 2.994 -0.00078 0.13803 0.00010 10.0 s 10:0 s	13:04.2	Data point 2	2.20002 mL	0.09546 mL	0.02439 mL	0.05000 mL	2.195	0.00638	0.59345	0.00041	10.0 s
14:55.8 Data point 5	13:40.0	Data point 3	2.20002 mL	0.09546 mL	0.04048 mL	0.05000 mL	2.388	0.01407	0.48508	0.00100	14.5 s
15:31.3 Data point 6	14:20.2	Data point 4						-0.00048	0.03777	0.00012	10.0 s
16:06.8 Data point 7	14:55.8	Data point 5	2.20002 mL	0.09546 mL	0.05952 mL	0.05000 mL	2.777	-0.00803	0.16330	0.00098	10.0 s
16:54.2 Data point 8 2,20002 mL 0.09546 mL 0.07272 mL 0.05000 mL 3.378 -0.00291 0.66595 0.00018 10.5 s 17:30.2 Data point 10 2,20002 mL 0.09546 mL 0.07658 mL 0.05000 mL 3.742 -0.000515 0.52674 0.00035 10.0 s 18:46.3 Data point 11 2,20002 mL 0.09546 mL 0.07658 mL 0.05000 mL 3.742 -0.00048 0.00212 0.00052 10.0 s 20:18.1 Data point 13 2,20002 mL 0.09546 mL 0.08488 mL 0.05000 mL 4.715 0.01154 0.48986 0.00081 10.0 s 20:18.1 Data point 14 2,20002 mL 0.09546 mL 0.08488 mL 0.05000 mL 4.912 -0.00828 0.71346 0.00081 10.0 s 21:29.0 Data point 15 2,20002 mL 0.09546 mL 0.09064 mL 0.05000 mL 5.499 0.01048 0.79345 0.00099 10.0 s 22:319.5 Data point 17 2,20002 mL 0.09546 mL 0.09546 mL 0.05000 mL			2.20002 mL	0.09546 mL	0.06526 mL	0.05000 mL	2.994	-0.00078	0.13803	0.00010	10.0 s
17:30.2 Data point 9	16:06.8	Data point 7						-0.01081	0.35970	0.00089	11.5 s
18:05.7 Data point 10								-0.00291	0.60595	0.00018	10.5 s
18:46.3 Data point 11 2.20002 mL 0.09546 mL 0.07843 mL 0.05000 mL 0.00500 mL 0.05000 m	17:30.2	Data point 9	2.20002 mL	0.09546 mL	0.07491 mL	0.05000 mL	3.574	0.00106	0.01066	0.00051	10.0 s
19:32.3 Data point 12			2.20002 mL	0.09546 mL	0.07658 mL	0.05000 mL	3.742	-0.00515	0.52674	0.00035	10.0 s
20:18.1 Data point 13	18:46.3	Data point 11								0.00052	10.0 s
20:53.6 Data point 14										0.00081	10.0 s
21:29.0 Data point 15	20:18.1	Data point 13						-0.00828	0.71346	0.00048	10.0 s
22:05.1 Data point 16			2.20002 mL	0.09546 mL	0.08699 mL	0.05000 mL	5.098	-0.00493	0.05980	0.00099	10.0 s
22:43.5 Data point 17 23:19.5 Data point 18 23:19.5 Data point 18 2.20002 mL 0.09546 mL 0.09290 mL 0.05000 mL 5.946 0.01316 0.44366 0.00098 10.0 s 23:55.0 Data point 19 2.20002 mL 0.09546 mL 0.09393 mL 0.05000 mL 6.380 0.00989 0.27105 0.00094 11.5 s 25:07.8 Data point 21 2.20002 mL 0.09546 mL 0.09546 mL 0.09396 mL 0.05000 mL 6.380 0.00989 0.27105 0.00094 11.5 s 25:07.8 Data point 21 2.20002 mL 0.09546 mL 0.09546 mL 0.09424 mL 0.05000 mL 6.380 0.00934 0.03282 0.00095 11.5 s 27:13.3 Data point 22 2.20002 mL 0.09546 mL 0.09546 mL 0.0944 mL 0.05000 mL 6.380 0.01004 0.037219 0.00094 11.5 s 0.00095 11.5 s 0.00100 12.5 c	21:29.0	Data point 15						-0.01048		0.00058	10.5 s
23:19.5 Data point 18			2.20002 mL	0.09546 mL	0.09064 mL	0.05000 mL	5.499			0.00098	13.0 s
23:55.0 Data point 19		•						-0.00818	0.73332	0.00047	10.5 s
24:30.9 Data point 20								-0.01316	0.44366		
25:07.8 Data point 21		•									
25:49.9 Data point 22											
27:13.3 Data point 23		•									
28:31.5 Data point 24											
30:12.4 Data point 25											
at 59.5 s 31:48.0 Data point 26	28:31.5	Data point 24	2.20002 mL	0.09546 mL	0.09508 mL	0.05000 mL	7.813	-0.04110	0.97438	0.00206	
33:23.8 Data point 27	30:12.4	Data point 25	2.20002 mL	0.09546 mL	0.09518 mL	0.05000 mL	8.240	-0.03638	0.97370	0.00182	
34:36.2 Data point 28	31:48.0	Data point 26	2.20002 mL	0.09546 mL	0.09525 mL	0.05000 mL	8.540	-0.01949	0.94537	0.00099	55.0 s
35:42.0 Data point 29	33:23.8	Data point 27						-0.01251	0.96457	0.00063	36.5 s
36:35.3 Data point 30			2.20002 mL	0.09546 mL	0.09551 mL	0.05000 mL	9.132	-0.01942	0.95772	0.00098	30.0 s
37:28.0 Data point 31	35:42.0	Data point 29	2.20002 mL	0.09546 mL	0.09570 mL	0.05000 mL	9.391	-0.01949	0.95986	0.00098	17.5 s
38:13.3 Data point 32	36:35.3	Data point 30	2.20002 mL	0.09546 mL	0.09595 mL	0.05000 mL	9.643	-0.01556	0.97273	0.00078	17.0 s
39:26.5 Data point 33	37:28.0	Data point 31	2.20002 mL	0.09546 mL	0.09628 mL	0.05000 mL	9.887	-0.01964	0.96938	0.00098	14.5 s
40:13.0 Data point 34 2.20002 mL 0.20245 mL 0.12578 mL 0.15000 mL 2.181 -0.00937 0.91245 0.00048 10.5 s 40:49.3 Data point 35 2.20002 mL 0.20245 mL 0.14367 mL 0.15000 mL 2.372 0.01563 0.60347 0.00099 14.0 s	38:13.3	Data point 32	2.20002 mL	0.09546 mL	0.09666 mL	0.05000 mL	10.078	-0.01833	0.96056	0.00092	11.5 s
40:49.3 Data point 35 2.20002 mL 0.20245 mL 0.14367 mL 0.15000 mL 2.372 0.01563 0.60347 0.00099 14.0 s	39:26.5	Data point 33	2.20002 mL	0.20245 mL	0.09666 mL	0.15000 mL	1.977	-0.01121	0.68100	0.00067	10.0 s
·	40:13.0	Data point 34						-0.00937		0.00048	10.5 s
	40:49.3	Data point 35	2.20002 mL	0.20245 mL	0.14367 mL	0.15000 mL	2.372	0.01563	0.60347	0.00099	14.0 s
			2.20002 mL	0.20245 mL	0.15602 mL	0.15000 mL	2.562	-0.00403	0.15580	0.00050	10.0 s

Reported at: 3/6/2018 1:24:22 PM



Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-03006 Instrument ID: T312060

Filename:

 $C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric\ high\ logP.t3r$

Events (continued)

Events (continueu)									
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared		dpH/dt time
42:04.6	Data point 37						-0.00114		0.00010	
42:40.6	Data point 38						0.00091	0.07557	0.00016	
43:16.5	Data point 39						-0.00686		0.00076	
43:52.1	Data point 40						-0.00225		0.00022	
44:27.6	Data point 41						-0.00311		0.00020	
45:03.0	Data point 42						-0.00336		0.00019	
45:38.5	Data point 43						-0.00227		0.00017	
46:24.3	Data point 44						-0.00909		0.00083	
47:03.3	Data point 45						-0.00440		0.00072	
47:38.8	Data point 46						0.00070	0.00237	0.00071	
48:14.2	Data point 47						-0.00945		0.00096	
48:49.7	Data point 48						-0.00439		0.00040	
49:25.6	Data point 49						-0.00627		0.00044	
50:01.6	Data point 50						-0.01045		0.00064	
50:37.5	Data point 51						-0.01202		0.00076	
51:13.9	Data point 52						-0.01959		0.00098	
51:51.4	Data point 53						-0.01017		0.00092	
52:29.4	Data point 54						-0.01877		0.00095	
53:23.6	Data point 55						-0.05573			59.5 s
54:59.2	Data point 56	2.20002 mL	0.20245 mL	0.20158 mL	0.15000 mL	7.661	-0.08712	0.98891	0.00433	Timed out at 59.5 s
56:40.0	Data point 57	2.20002 mL	0.20245 mL	0.20167 mL	0.15000 mL	7.903	-0.06287	0.97801	0.00314	Timed out at 59.5 s
58:20.7	Data point 58	2.20002 mL	0.20245 mL	0.20176 mL	0.15000 mL	8.166	-0.04207	0.98021	0.00210	Timed out at 59.5 s
1:00:11.8	Data point 59	2 20002 ml	0 20245 ml	0 20195 ml	0.15000 ml	8 624	-0.01762	0 98142	0.00088	
	Data point 60						-0.01366		0.00097	
	Data point 61						-0.00856		0.00075	
1:03:24.9	•						-0.01587		0.00092	
	Data point 63						-0.01697		0.00089	
1:04:52.1	Data point 64						-0.01307		0.00095	
	Data point 65						-0.01920		0.00099	
1:06:58.1							-0.01558	0.84604	0.00084	10.0 s
1:07:44.6	Data point 67	2.20002 mL	0.31409 mL	0.23408 mL	0.45000 mL	2.187	0.01489	0.60774	0.00094	12.5 s
1:08:23.0	Data point 68	2.20002 mL	0.31409 mL	0.25320 mL	0.45000 mL	2.385	0.00088	0.09557	0.00014	10.0 s
1:08:58.7	Data point 69	2.20002 mL	0.31409 mL	0.26618 mL	0.45000 mL	2.573	-0.00649	0.12500	0.00091	10.0 s
1:09:34.4	Data point 70	2.20002 mL	0.31409 mL	0.27545 mL	0.45000 mL	2.787	-0.00191	0.06611	0.00037	10.0 s
1:10:09.9	Data point 71	2.20002 mL	0.31409 mL	0.28206 mL	0.45000 mL	2.981	0.00572	0.62627	0.00036	
	Data point 72						-0.01361		0.00096	
	Data point 73						0.00911	0.20975	0.00098	
1:11:57.5							-0.00440		0.00060	
	Data point 75						0.00458	0.65602	0.00028	
	Data point 76						-0.00707		0.00081	
	Data point 77						-0.00994		0.00095	
	Data point 78						-0.00154		0.00059	
	Data point 79						-0.00829		0.00084	
	Data point 80						0.01636	0.66406	0.00099	
	Data point 81						0.00586		0.00081	
	Data point 82						-0.00746		0.00075	
	Data point 83						-0.01384		0.00086	
	Data point 84						0.00783		0.00093	
	Data point 85						-0.01816		0.00094	
	Data point 86						-0.00051		0.00098	
	Data point 87						-0.01711		0.00095	
1.22.22.2	Data point 88	2.20002 mL	U.3 14U9 ML	U.31787 ML	U.45UUU ML	7.035	-0.09484	0.98262	0.004/3	Timed out at

Reported at: 3/6/2018 1:24:22 PM

59.5 s

Assay Events



Sample name: M15_octanol Experiment start time: 3/3/2018 8:26:22 AM

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Events (continued)

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
1:23:57.9	Data point 89	2.20002 mL	0.31409 mL	0.31192 mL	0.45000 mL	7.282	-0.07791	0.98221	0.00388	Timed out
										at 59.5 s
1:25:28.4	Data point 90	2.20002 mL	0.31409 mL	0.31197 mL	0.45000 mL	7.488	-0.07538	0.98886	0.00374	Timed out
										at 59.5 s
1:27:04.1	Data point 91	2.20002 mL	0.31409 mL	0.31204 mL	0.45000 mL	7.764	-0.08258	0.98709	0.00411	Timed out
										at 59.5 s
1:28:39.7	Data point 92	2.20002 mL	0.31409 mL	0.31211 mL	0.45000 mL	8.037	-0.06736	0.98943	0.00334	Timed out
										at 59.5 s
1:30:15.3	Data point 93	2.20002 mL	0.31409 mL	0.31221 mL	0.45000 mL	8.330	-0.04664	0.96334	0.00235	
										at 59.5 s
1:31:51.0		2.20002 mL		0.31242 mL			-0.01725	0.85754	0.00092	21.0 s
1:32:47.7		2.20002 mL	0.31409 mL	0.31268 mL	0.45000 mL	9.073	-0.00613	0.09179	0.00100	27.5 s
1:33:50.9	Data point 96	2.20002 mL	0.31409 mL	0.31305 mL	0.45000 mL	9.365	-0.01654	0.78160	0.00093	14.5 s
1:34:41.2		2.20002 mL	0.31409 mL	0.31345 mL	0.45000 mL	9.616	-0.01302	0.60064	0.00083	12.5 s
1:35:29.4		2.20002 mL	0.31409 mL	0.31388 mL	0.45000 mL	9.835	0.00468	0.05381	0.00100	21.5 s
1:36:26.7	Data point 99	2.20002 mL	0.31409 mL	0.31439 mL	0.45000 mL	10.049	0.00989	0.27168	0.00094	11.5 s
1:36:47.2	Assay volumes	2.20002 mL	0.31409 mL	0.31439 mL	0.45000 mL					



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Filename: C:\Sirius_T3\Meh	tap\20180302_exp2	9_logP_T3-2\180	C-03006_M15_octano	I_pH-metric high log
Assay Settings				
Setting <i>General Settings</i>	Value	Original Value	Date/Time changed	Imported from
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	10.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings				
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No .			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium	0.001			
ISA water volume	2.20 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.050 mL			
Partition solvent added	Automatic			
After partition addition, stir for	1 seconds			
Sample Sonication Sonicate	Voc			
	Yes No			
Adjust pH for sonication Sonicate for	300 seconds			
After sonication stir for	5 seconds			
Sample Dissolution	3 36601103			
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution				
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge	1070			
Perform a carbonate purge	No			
Temperature Control	140			
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			
Stir speed of	50%			
Titration 1				
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	50%			
Titration 2				
Titrate from	Low to high pH			
Add additional water	0.00 mL			
A delitional montition and continues	0.400 mal			

Reported at: 3/6/2018 1:24:22 PM

After pH adjust stir for

Stir to allow partitioning for

Stirrer speed for partitioning

Additional partition solvent volume 0.100 mL Additional partition solvent added Automatic

30 seconds

15 seconds

55%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from	
Titration 3					
Titrate from	Low to high pH				
Add additional water	0.00 mL				
Additional partition solvent volume	0.300 mL				
Additional partition solvent added	Automatic				
After pH adjust stir for	30 seconds				
Stir to allow partitioning for	15 seconds				
Stirrer speed for partitioning	60%				
Data Point Stability					
Stir during data point collection	No				
Delay before data point collection	0 seconds				
Number of points to average	20 points				
Time interval between points	0.50 seconds				
Required maximum standard deviation	0.00100 dpH/dt				
Stability timeout after	60 seconds				

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.111	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r
Four-Plus S	0.9988	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r
Four-Plus jH	1.0	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r
Four-Plus jOH	-0.8	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r
Base concentration factor	1.000	3/3/2018 8:26:22 AM	C:\Sirius_T3\KOH18B27.t3r
Acid concentration factor	0.999	3/3/2018 8:26:22 AM	C:\Sirius_T3\HCl18C02.t3r

Instrument Settings

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module		T3DM1200361	3/31/2009 5:24:52 AM
Dispenser 0	Water		3/31/2009 5:25:05 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)	00 00 0010	0/07/0040 40 05 50 414
Titrant	Water (0.15 M KCI)	02-06-2018	2/27/2018 10:05:59 AM
Dispenser 2 Syringe volume	Acid 0.5 mL		3/31/2009 5:25:11 AM
Firmware version	1.2.1(r2)		
Titrant	Acid (0.5 M HCI)	02-27-2018	2/27/2018 10:27:22 AM
Dispenser 1	Base	02 27 2010	3/31/2009 5:25:21 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Base (0.5 M KOH)	9/22/2017	2/27/2018 10:21:22 AM
Dispenser 5	Cosolvent		3/31/2009 5:26:24 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		0/04/0000 5:00:40 ANA
Distribution valve 5	Distribution Valve		3/31/2009 5:28:19 AM
Firmware version Port A	1.1.3 Methanol (80%, 0.15 M KCl)	09-26-17	2/7/2018 9:42:01 AM
Port B	Cyclohexane	11-01-17	2/27/2018 10:37:57 AM
Dispenser 3	Buffer	11-01-17	8/3/2010 5:05:16 AM
Syringe volume	0.5 mL		0,0,2010 0.00.107
Firmware version	1.2.1(r2)		
Titrant	Dodecane	2018/01/31	2/28/2018 10:18:04 AM
Dispenser 6	Octanol		10/22/2010 10:52:43 AM

Reported at: 3/6/2018 1:24:22 PM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Syringe volume Firmware version	0.5 mL 1.2.1(r2)		
Titrant	Octanol	01-31-2018	2/27/2018 9:59:35 AM
Titrator	Octarior		3/31/2009 5:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2	1011111200101	0,0 1,2000 0.2 1.17 7
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0923	1/23/2018 2:01:00 PM
E0 calibration	+4.91 mV		3/3/2018 8:26:50 AM
Filling solution	3M KCI	KCL097	3/2/2018 9:43:24 AM
Liquids	FOO/ IDA-FOO/ Motor		2/2/2010 0:45:12 AM
Wash 1 Wash 2	50% IPA:50% Water 0.5% Trition X-100 in H20		3/2/2018 9:45:12 AM 3/2/2018 9:45:15 AM
Buffer position 1	pH7 Wash		3/2/2018 9:45:18 AM
Buffer position 2	pH 7		3/2/2018 9:45:21 AM
Storage position	pi i i		3/2/2018 9:44:44 AM
Wash water	6.9e+003 mL	02-27-2018	2/27/2018 9:54:39 AM
Waste	8.5e+003 mL	0 2	11/28/2017 10:36:29 AM
Temperature controller			8/5/2010 6:35:13 AM
Turbidity detector			3/31/2009 5:24:45 AM
Spectrometer		074811	11/23/2010 11:22:28 AM
Dip probe		10196	
Wavelength coefficient A0	183.333		
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		44/00/0040 44:00:00 ABA
Total lamp lit time	120:41:49		11/23/2010 11:22:28 AM
Calibrated on Integration time	2/27/2018 10:40:38 AM 40		
Scans averaged	10		
Autoloader	10	T3AI 1200345	11/10/2015 9:34:13 AM
Left-right axis firmware version	1.17 Al1Dl2DO2 Stepper 2	10/12/2000	11/10/2010 0.01.10744
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Configuration			
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period Titrant tube volume	5 minute(s) 1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation			
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s 30%		
E0 calibration preparation stir speed E0 calibration buffer wash stir duration	30% 5 s		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		

Reported at: 3/6/2018 1:24:22 PM



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

mL					
mL					
00					
	00				

Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00
Experiment Log		

- [2:38] Air gap created for Water (0.15 M KCI)
- [2:38] Air gap created for Acid (0.5 M HCI)
- [2:38] Air gap created for Base (0.5 M KOH)
- [2:39] Air gap released for Water (0.15 M KCI)
- [2:43] Titrator arm moved over Titration position
- [2:43] Titration 1 of 3
- [2:43] Adding initial titrants
- [2:43] Automatically add 2.20000 mL of water
- [3:19] Dispensed 2.200024 mL of Water (0.15 M KCI)
- [3:24] Titrator arm moved over Drain
- [9:05] Titrator arm moved to Titration position
- [9:05] Argon flow rate set to 100
- [9:05] Stirrer speed set to 10
- [9:10] Automatically add 0.05000 mL of Octanol
- [9:11] Dispensed 0.050000 mL of Octanol
- [9:12] Initial pH = 6.94
- [9:12] Iterative adjust 6.94 -> 2.00
- [9:12] pH 6.94 -> 2.00
- [9:15] Air gap released for Acid (0.5 M HCI)
- [9:16] Dispensed 0.091322 mL of Acid (0.5 M HCI)
- [9:21] pH 2.02 -> 2.00
- [9:21] Dispensed 0.002563 mL of Acid (0.5 M HCI)
- [9:26] pH 2.01 -> 2.00
- [9:26] Dispensed 0.001576 mL of Acid (0.5 M HCI)
- [9:31] Holding pH 2.00
- [11:32] Stirrer speed set to 0
- [11:32] Stirrer speed set to 50
- [11:32] Iterative adjust 1.99 -> 2.00
- [12:17] Stirrer speed set to 0
- [12:27] Datapoint id 1 collected
- [12:27] Stirrer speed set to 50 [12:32] pH 1.99 -> 2.19
- [12:32] Using cautious pH adjust
- [12:33] Air gap released for Base (0.5 M KOH)
- [12:34] Dispensed 0.012465 mL of Base (0.5 M KOH)
- [12:39] Stepping pH = 2.08
- [12:39] Dispensed 0.009243 mL of Base (0.5 M KOH)
- [12:44] Stepping pH = 2.16
- [12:44] Dispensed 0.002681 mL of Base (0.5 M KOH)
- [12:49] Stepping pH = 2.19
- [13:05] Stirrer speed set to 0
- [13:15] Datapoint id 2 collected
- Reported at: 3/6/2018 1:24:22 PM



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [13:15] Charge balance equation is out by 2.1%
- [13:15] Stirrer speed set to 50
- [13:20] pH 2.20 -> 2.40
- [13:20] Using charge balance adjust
- [13:20] Dispensed 0.016087 mL of Base (0.5 M KOH)
- [13:40] Stirrer speed set to 0
- [13:55] Datapoint id 3 collected
- [13:55] Charge balance equation is out by -5.3%
- [13:55] Stirrer speed set to 50
- [14:00] pH 2.39 -> 2.59
- [14:00] Using charge balance adjust
- [14:00] Dispensed 0.011124 mL of Base (0.5 M KOH)
- [14:21] Stirrer speed set to 0
- [14:31] Datapoint id 4 collected
- [14:31] Charge balance equation is out by -5.5%
- [14:31] Stirrer speed set to 50
- [14:36] pH 2.59 -> 2.79
- [14:36] Using charge balance adjust
- [14:36] Dispensed 0.007926 mL of Base (0.5 M KOH)
- [14:56] Stirrer speed set to 0
- [15:06] Datapoint id 5 collected
- [15:06] Charge balance equation is out by -5.2%
- [15:06] Stirrer speed set to 50
- [15:11] pH 2.78 -> 2.98
- [15:11] Using charge balance adjust
- [15:12] Dispensed 0.005738 mL of Base (0.5 M KOH)
- [15:32] Stirrer speed set to 0
- [15:42] Datapoint id 6 collected
- [15:42] Charge balance equation is out by 5.1%
- [15:42] Stirrer speed set to 50
- [15:47] pH 3.00 -> 3.20
- [15:47] Using charge balance adjust
- [15:47] Dispensed 0.004022 mL of Base (0.5 M KOH)
- [16:07] Stirrer speed set to 0
- [16:19] Datapoint id 7 collected
- [16:19] Charge balance equation is out by -21.3%
- [16:19] Stirrer speed set to 50
- [16:24] pH 3.16 -> 3.36
- [16:24] Using cautious pH adjust
- [16:24] Dispensed 0.001529 mL of Base (0.5 M KOH)
- [16:29] Stepping pH = 3.23
- [16:29] Dispensed 0.001693 mL of Base (0.5 M KOH)
- [16:34] Stepping pH = 3.35
- [16:34] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [16:39] Stepping pH = 3.38
- [16:55] Stirrer speed set to 0
- [17:05] Datapoint id 8 collected
- [17:05] Charge balance equation is out by -11.9%
- [17:05] Stirrer speed set to 50
- [17:10] pH 3.38 -> 3.58
- [17:10] Using charge balance adjust
- [17:10] Dispensed 0.002187 mL of Base (0.5 M KOH)
- [17:31] Stirrer speed set to 0
- [17:41] Datapoint id 9 collected
- [17:41] Charge balance equation is out by -4.7%
- [17:41] Stirrer speed set to 50
- [17:46] pH 3.58 -> 3.78
- [17:46] Using charge balance adjust
- [17:46] Dispensed 0.001670 mL of Base (0.5 M KOH)



pH-metric high logP **Pion** Assay name: Analyst: Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [18:06] Stirrer speed set to 0
- [18:16] Datapoint id 10 collected
- [18:16] Charge balance equation is out by -19.9%
- [18:16] Stirrer speed set to 50
- [18:21] pH 3.75 -> 3.95
- [18:21] Using cautious pH adjust
- [18:21] Dispensed 0.000729 mL of Base (0.5 M KOH)
- [18:26] Stepping pH = 3.79
- [18:26] Dispensed 0.001129 mL of Base (0.5 M KOH)
- [18:32] Stepping pH = 3.95
- [18:47] Stirrer speed set to 0 [18:57] Datapoint id 11 collected
- [18:57] Charge balance equation is out by -27.1%
- [18:57] Stirrer speed set to 50
- [19:02] pH 3.98 -> 4.18
- [19:02] Using cautious pH adjust
- [19:02] Dispensed 0.000706 mL of Base (0.5 M KOH)
- [19:07] Stepping pH = 4.04
- [19:07] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [19:12] Stepping pH = 4.04
- [19:12] Dispensed 0.002752 mL of Base (0.5 M KOH)
- [19:18] Stepping pH = 4.69
- [19:33] Stirrer speed set to 0
- [19:43] Datapoint id 12 collected
- [19:43] Charge balance equation is out by -214.2%
- [19:43] Stirrer speed set to 50
- [19:48] pH 4.72 -> 4.92
- [19:48] Using cautious pH adjust
- [19:48] Dispensed 0.001058 mL of Base (0.5 M KOH)
- [19:53] Stepping pH = 4.82
- [19:53] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [19:58] Stepping pH = 4.90
- [19:58] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [20:03] Stepping pH = 4.92
- [20:18] Stirrer speed set to 0
- [20:29] Datapoint id 13 collected
- [20:29] Charge balance equation is out by 2.5%
- [20:29] Stirrer speed set to 50
- [20:34] pH 4.92 -> 5.12
- [20:34] Using charge balance adjust
- [20:34] Dispensed 0.002117 mL of Base (0.5 M KOH)
- [20:54] Stirrer speed set to 0
- [21:04] Datapoint id 14 collected
- [21:04] Charge balance equation is out by -10.3%
- [21:04] Stirrer speed set to 50
- [21:09] pH 5.10 -> 5.30
- [21:09] Using charge balance adjust
- [21:09] Dispensed 0.001976 mL of Base (0.5 M KOH)
- [21:29] Stirrer speed set to 0
- [21:40] Datapoint id 15 collected
- [21:40] Charge balance equation is out by -6.8%
- [21:40] Stirrer speed set to 50
- [21:45] pH 5.30 -> 5.50
- [21:45] Using charge balance adjust
- [21:45] Dispensed 0.001670 mL of Base (0.5 M KOH)
- [22:05] Stirrer speed set to 0
- [22:19] Datapoint id 16 collected
- [22:19] Charge balance equation is out by 1.3%
- [22:19] Stirrer speed set to 50



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [22:24] pH 5.50 -> 5.70
- [22:24] Using charge balance adjust
- [22:24] Dispensed 0.001317 mL of Base (0.5 M KOH)
- [22:44] Stirrer speed set to 0
- [22:54] Datapoint id 17 collected
- [22:54] Charge balance equation is out by 10.1%
- [22:54] Stirrer speed set to 50
- [23:00] pH 5.73 -> 5.93
- [23:00] Using charge balance adjust
- [23:00] Dispensed 0.000941 mL of Base (0.5 M KOH)
- [23:20] Stirrer speed set to 0
- [23:30] Datapoint id 18 collected
- [23:30] Charge balance equation is out by 7.6%
- [23:30] Stirrer speed set to 50
- [23:35] pH 5.95 -> 6.15
- [23:35] Using charge balance adjust
- [23:35] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [23:55] Stirrer speed set to 0
- [24:06] Datapoint id 19 collected
- [24:06] Charge balance equation is out by 8.2%
- [24:06] Stirrer speed set to 50
- [24:11] pH 6.18 -> 6.38
- [24:11] Using charge balance adjust
- [24:11] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [24:31] Stirrer speed set to 0
- [24:43] Datapoint id 20 collected
- [24:43] Charge balance equation is out by 0.3%
- [24:43] Stirrer speed set to 50
- [24:48] pH 6.39 -> 6.59
- [24:48] Using charge balance adjust
- [24:48] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [25:08] Stirrer speed set to 0
- [25:20] Datapoint id 21 collected
- [25:20] Charge balance equation is out by -64.1%
- [25:20] Stirrer speed set to 50
- [25:25] pH 6.47 -> 6.67
- [25:25] Using cautious pH adjust
- [25:25] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [25:30] Stepping pH = 6.48
- [25:30] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [25:35] Stepping pH = 6.72
 - 5.50] Otepping pri = 0.72
- [25:50] Stirrer speed set to 0
- [26:38] Datapoint id 22 collected
- [26:38] Charge balance equation is out by -84.6%
- [26:38] Stirrer speed set to 50
- [26:43] pH 6.88 -> 7.08
- [26:43] Using cautious pH adjust
- [26:43] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [26:48] Stepping pH = 6.89
- [26:48] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [26:53] Stepping pH = 6.99
- [26:53] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [26:59] Stepping pH = 7.16
- [27:14] Stirrer speed set to 0
- [27:56] Datapoint id 23 collected
- [27:56] Charge balance equation is out by -162.3%
- [27:56] Stirrer speed set to 50
- [28:01] pH 7.29 -> 7.49
- [28:01] Using cautious pH adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [28:01] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [28:07] Stepping pH = 7.31
- [28:07] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [28:12] Stepping pH = 7.38
- [28:12] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [28:17] Stepping pH = 7.59
- [28:32] Stirrer speed set to 0
- [29:32] Datapoint id 24 collected
- [29:32] Charge balance equation is out by -197.2%
- [29:32] Stirrer speed set to 50
- [29:37] pH 7.82 -> 8.02
- [29:37] Using cautious pH adjust
- [29:37] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [29:42] Stepping pH = 7.84
- [29:42] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [29:47] Stepping pH = 7.88
- [29:47] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [29:53] Stepping pH = 7.99
- [29:53] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [29:58] Stepping pH = 8.15
- [30:13] Stirrer speed set to 0
- [31:13] Datapoint id 25 collected
- [31:13] Charge balance equation is out by -439.5%
- [31:13] Stirrer speed set to 50
- [31:18] pH 8.28 -> 8.48
- [31:18] Using cautious pH adjust
- [31:18] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [31:23] Stepping pH = 8.34
- [31:23] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [31:28] Stepping pH = 8.41
- [31:28] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [31:28] Dispensed 0.00002 [31:33] Stepping pH = 8.51
- [31:48] Stirrer speed set to 0
- [32:44] Datapoint id 26 collected
- [32:44] Charge balance equation is out by -243.3%
- [32:44] Stirrer speed set to 50
- [32:49] pH 8.54 -> 8.74
- [32:49] Using cautious pH adjust
- [32:49] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [32:54] Stepping pH = 8.57
- [32:54] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [32:59] Stepping pH = 8.63
- [32:59] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [33:04] Stepping pH = 8.72
- [33:04] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [33:09] Stepping pH = 8.81
- [33:24] Stirrer speed set to 0
- [34:01] Datapoint id 27 collected
- [34:01] Charge balance equation is out by -320.9%
- [34:01] Stirrer speed set to 50
- [34:06] pH 8.87 -> 9.07
- [34:06] Using cautious pH adjust
- [34:06] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [34:11] Stepping pH = 8.90
- [34:11] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [34:16] Stepping pH = 8.98
- [34:16] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [34:21] Stepping pH = 9.08
- [34:37] Stirrer speed set to 0



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [35:07] Datapoint id 28 collected
- [35:07] Charge balance equation is out by -173.9%
- [35:07] Stirrer speed set to 50
- [35:12] pH 9.15 -> 9.35
- [35:12] Using cautious pH adjust
- [35:12] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [35:17] Stepping pH = 9.18
- [35:17] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [35:22] Stepping pH = 9.28
- [35:22] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [35:27] Stepping pH = 9.36
- [35:42] Stirrer speed set to 0
- [36:00] Datapoint id 29 collected
- [36:00] Charge balance equation is out by -130.4%
- [36:00] Stirrer speed set to 50
- [36:05] pH 9.40 -> 9.60
- [36:05] Using cautious pH adjust
- [36:05] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [36:10] Stepping pH = 9.42
- [36:10] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [36:15] Stepping pH = 9.57
- [36:15] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [36:21] Stepping pH = 9.63
- [36:36] Stirrer speed set to 0
- [36:53] Datapoint id 30 collected
- [36:53] Charge balance equation is out by -89.0%
- [36:53] Stirrer speed set to 50
- [36:58] pH 9.65 -> 9.85
- [36:58] Using cautious pH adjust
- [36:58] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [37:03] Stepping pH = 9.69
- [37:03] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [37:08] Stepping pH = 9.83
- [37:08] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [37:13] Stepping pH = 9.88
- [37:28] Stirrer speed set to 0
- [37:43] Datapoint id 31 collected
- [37:43] Charge balance equation is out by -50.9%
- [37:43] Stirrer speed set to 50
- [37:48] pH 9.89 -> 10.05
- [37:48] Using cautious pH adjust
- [37:48] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [37:53] Stepping pH = 9.93
- [37:53] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [37:59] Stepping pH = 10.04
- [38:14] Stirrer speed set to 0
- [38:25] Datapoint id 32 collected
- [38:25] Charge balance equation is out by -35.0%
- [38:25] Titration 2 of 3
- [38:25] Adding initial titrants
- [38:25] Automatically add 0.10000 mL of Octanol
- [38:28] Dispensed 0.100000 mL of Octanol
- [38:28] Stirrer speed set to 10
- [38:29] Stirrer speed set to 55
- [38:29] Iterative adjust 10.08 -> 2.00
- [38:29] pH 10.08 -> 2.00
- [38:31] Dispensed 0.097789 mL of Acid (0.5 M HCl)
- [38:36] pH 2.06 -> 2.00
- [38:37] Dispensed 0.009196 mL of Acid (0.5 M HCl)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [39:27] Stirrer speed set to 0
- [39:37] Datapoint id 33 collected
- [39:37] Stirrer speed set to 55
- [39:42] pH 1.98 -> 2.18
- [39:42] Using cautious pH adjust
- [39:42] Dispensed 0.013711 mL of Base (0.5 M KOH)
- [39:48] Stepping pH = 2.06
- [39:48] Dispensed 0.011947 mL of Base (0.5 M KOH)
- [39:53] Stepping pH = 2.15
- [39:53] Dispensed 0.003457 mL of Base (0.5 M KOH)
- [39:58] Stepping pH = 2.18
- [40:13] Stirrer speed set to 0
- [40:24] Datapoint id 34 collected
- [40:24] Charge balance equation is out by -6.2%
- [40:24] Stirrer speed set to 55
- [40:29] pH 2.18 -> 2.38
- [40:29] Using charge balance adjust
- [40:30] Dispensed 0.017897 mL of Base (0.5 M KOH)
- [40:50] Stirrer speed set to 0
- [41:04] Datapoint id 35 collected
- [41:04] Charge balance equation is out by -6.3%
- [41:04] Stirrer speed set to 55
- [41:09] pH 2.38 -> 2.58
- [41:09] Using charge balance adjust
- [41:09] Dispensed 0.012347 mL of Base (0.5 M KOH)
- [41:29] Stirrer speed set to 0
- [41:39] Datapoint id 36 collected
- [41:39] Charge balance equation is out by -7.2%
- [41:39] Stirrer speed set to 55
- [41:44] pH 2.57 -> 2.77
- [41:44] Using charge balance adjust
- [41:45] Dispensed 0.008796 mL of Base (0.5 M KOH)
- [42:05] Stirrer speed set to 0
- [42:16] Datapoint id 37 collected
- [42:16] Charge balance equation is out by -0.8%
- [42:16] Stirrer speed set to 55
- [42:21] pH 2.77 -> 2.97
- [42:21] Using charge balance adjust
- [42:21] Dispensed 0.006256 mL of Base (0.5 M KOH)
- [42:41] Stirrer speed set to 0
- [42:51] Datapoint id 38 collected
- [42:51] Charge balance equation is out by -4.0%
- [42:51] Stirrer speed set to 55
- [42:56] pH 2.97 -> 3.17
- [42:56] Using charge balance adjust
- [42:57] Dispensed 0.004563 mL of Base (0.5 M KOH)
- [43:17] Stirrer speed set to 0
- [43:27] Datapoint id 39 collected
- [43:27] Charge balance equation is out by -5.3%
- [43:27] Stirrer speed set to 55
- [43:32] pH 3.16 -> 3.36
- [43:32] Using charge balance adjust
- [43:32] Dispensed 0.003387 mL of Base (0.5 M KOH)
- [43:52] Stirrer speed set to 0
- [44:02] Datapoint id 40 collected
- [44:02] Charge balance equation is out by -3.5%
- [44:03] Stirrer speed set to 55
- [44:08] pH 3.36 -> 3.56
- [44:08] Using charge balance adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [44:08] Dispensed 0.002587 mL of Base (0.5 M KOH)
- [44:28] Stirrer speed set to 0
- [44:38] Datapoint id 41 collected
- [44:38] Charge balance equation is out by -11.2%
- [44:38] Stirrer speed set to 55
- [44:43] pH 3.55 -> 3.75
- [44:43] Using charge balance adjust
- [44:43] Dispensed 0.002140 mL of Base (0.5 M KOH)
- [45:03] Stirrer speed set to 0
- [45:13] Datapoint id 42 collected
- [45:13] Charge balance equation is out by -2.6%
- [45:13] Stirrer speed set to 55
- [45:19] pH 3.75 -> 3.95
- [45:19] Using charge balance adjust
- [45:19] Dispensed 0.001929 mL of Base (0.5 M KOH)
- [45:39] Stirrer speed set to 0
- [45:49] Datapoint id 43 collected
- [45:49] Charge balance equation is out by 18.3%
- [45:49] Stirrer speed set to 55
- [45:54] pH 3.98 -> 4.18
- [45:54] Using cautious pH adjust
- [45:54] Dispensed 0.000988 mL of Base (0.5 M KOH)
- [45:59] Stepping pH = 4.08
- [45:59] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [46:04] Stepping pH = 4.16
- [46:05] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [46:10] Stepping pH = 4.19
- [46:25] Stirrer speed set to 0
- [46:38] Datapoint id 44 collected
- [46:38] Charge balance equation is out by 2.9%
- [46:38] Stirrer speed set to 55
- [46:43] pH 4.19 -> 4.39
- [46:43] Using charge balance adjust
- [46:44] Dispensed 0.002117 mL of Base (0.5 M KOH)
- [47:04] Stirrer speed set to 0
- [47:14] Datapoint id 45 collected
- [47:14] Charge balance equation is out by -2.0%
- [47:14] Stirrer speed set to 55
- [47:19] pH 4.39 -> 4.59
- [47:19] Using charge balance adjust
- [47:19] Dispensed 0.002234 mL of Base (0.5 M KOH)
- [47:39] Stirrer speed set to 0
- [47:49] Datapoint id 46 collected
- [47:49] Charge balance equation is out by -1.0%
- [47:49] Stirrer speed set to 55
- [47:54] pH 4.59 -> 4.79
- [47:54] Using charge balance adjust
- [47:54] Dispensed 0.002187 mL of Base (0.5 M KOH)
- [48:15] Stirrer speed set to 0
- [48:25] Datapoint id 47 collected
- [48:25] Charge balance equation is out by 11.4%
- [48:25] Stirrer speed set to 55
- [48:30] pH 4.81 -> 5.01
- [48:30] Using charge balance adjust
- [48:30] Dispensed 0.001976 mL of Base (0.5 M KOH)
- [48:50] Stirrer speed set to 0
- [49:01] Datapoint id 48 collected
- [49:01] Charge balance equation is out by 14.0%
- [49:01] Stirrer speed set to 55



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [49:06] pH 5.04 -> 5.24
- [49:06] Using charge balance adjust
- [49:06] Dispensed 0.001576 mL of Base (0.5 M KOH)
- [49:26] Stirrer speed set to 0
- [49:37] Datapoint id 49 collected
- [49:37] Charge balance equation is out by 3.3%
- [49:37] Stirrer speed set to 55
- [49:42] pH 5.26 -> 5.46
- [49:42] Using charge balance adjust
- [49:42] Dispensed 0.001176 mL of Base (0.5 M KOH)
- [50:02] Stirrer speed set to 0
- [50:13] Datapoint id 50 collected
- [50:13] Charge balance equation is out by -3.5%
- [50:13] Stirrer speed set to 55
- [50:18] pH 5.46 -> 5.66
- [50:18] Using charge balance adjust
- [50:18] Dispensed 0.000847 mL of Base (0.5 M KOH)
- [50:38] Stirrer speed set to 0
- [50:49] Datapoint id 51 collected
- [50:49] Charge balance equation is out by 7.1%
- [50:49] Stirrer speed set to 55
- [50:54] pH 5.68 -> 5.88
- [50:54] Using charge balance adjust
- [50:54] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [51:14] Stirrer speed set to 0
- [51:26] Datapoint id 52 collected
- [51:26] Charge balance equation is out by -2.4%
- [51:26] Stirrer speed set to 55
- [51:31] pH 5.88 -> 6.08
- [51:31] Using charge balance adjust
- [51:32] Dispensed 0.000400 mL of Base (0.5 M KOH)
- [51:52] Stirrer speed set to 0
- [52:04] Datapoint id 53 collected
- [52:04] Charge balance equation is out by -12.9%
- [52:04] Stirrer speed set to 55
- [52:09] pH 6.06 -> 6.26
- [52:09] Using charge balance adjust
- [52:10] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [52:30] Stirrer speed set to 0
- [52:53] Datapoint id 54 collected
- [52:53] Charge balance equation is out by 17.6%
- [52:53] Stirrer speed set to 55
- [52:58] pH 6.30 -> 6.50
- [52:58] Using cautious pH adjust
- [52:59] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [53:04] Stepping pH = 6.32
- [53:04] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [53:09] Stepping pH = 6.75
- [53:24] Stirrer speed set to 0
- [54:24] Datapoint id 55 collected
- [54:24] Charge balance equation is out by -87.2%
- [54:24] Stirrer speed set to 55
- [54:29] pH 7.05 -> 7.25
- [54:29] Using cautious pH adjust
- [54:29] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [54:34] Stepping pH = 7.07
- [54:34] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [54:39] Stepping pH = 7.11
- [54:39] Dispensed 0.000118 mL of Base (0.5 M KOH)



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [54:45] Stepping pH = 7.46
- [55:00] Stirrer speed set to 0
- [56:00] Datapoint id 56 collected
- [56:00] Charge balance equation is out by -289.0%
- [56:00] Stirrer speed set to 55
- [56:05] pH 7.69 -> 7.89
- [56:05] Using cautious pH adjust
- [56:05] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [56:10] Stepping pH = 7.74
- [56:10] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [56:15] Stepping pH = 7.79
- [56:15] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [56:20] Stepping pH = 7.84
- [56:20] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [56:25] Stepping pH = 7.93 [56:40] Stirrer speed set to 0
- [57:40] Datapoint id 57 collected
- [57:40] Charge balance equation is out by -381.5%
- [57:40] Stirrer speed set to 55
- [57:45] pH 7.90 -> 8.10
- [57:45] Using cautious pH adjust
- [57:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [57:51] Stepping pH = 7.94
- [57:51] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [57:56] Stepping pH = 8.00
- [57:56] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [58:01] Stepping pH = 8.08
- [58:01] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [58:06] Stepping pH = 8.17
- [58:21] Stirrer speed set to 0
- [59:21] Datapoint id 58 collected [59:21] Charge balance equation is out by -472.6%
- [59:21] Stirrer speed set to 55
- [59:26] pH 8.19 -> 8.39
- [59:26] Using cautious pH adjust
- [59:26] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [59:31] Stepping pH = 8.22
- [59:31] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [59:37] Stepping pH = 8.27
- [59:37] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [59:42] Stepping pH = 8.33
- [59:42] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [59:47] Stepping pH = 8.37
- [59:47] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [59:52] Stepping pH = 8.35
- [59:52] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [59:57] Stepping pH = 8.58 [1:00:12] Stirrer speed set to 0
- [1:00:49] Datapoint id 59 collected
- [1:00:49] Charge balance equation is out by -1,013.4%
- [1:00:49] Stirrer speed set to 55
- [1:00:54] pH 8.64 -> 8.84
- [1:00:54] Using cautious pH adjust
- [1:00:54] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:00:59] Stepping pH = 8.66
- [1:00:59] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:01:04] Stepping pH = 8.73
- [1:01:04] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:01:09] Stepping pH = 8.83



Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [1:01:09] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:01:15] Stepping pH = 8.89
- [1:01:30] Stirrer speed set to 0
- [1:01:49] Datapoint id 60 collected
- [1:01:49] Charge balance equation is out by -319.6%
- [1:01:49] Stirrer speed set to 55
- [1:01:54] pH 8.89 -> 9.09
- [1:01:54] Using cautious pH adjust
- [1:01:54] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:01:59] Stepping pH = 8.91
- [1:01:59] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:02:04] Stepping pH = 8.99
- [1:02:05] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:02:10] Stepping pH = 9.11
- [1:02:25] Stirrer speed set to 0
- [1:02:55] Datapoint id 61 collected
- [1:02:55] Charge balance equation is out by -216.7%
- [1:02:55] Stirrer speed set to 55
- [1:03:00] pH 9.15 -> 9.35
- [1:03:00] Using cautious pH adjust
- [1:03:00] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:03:05] Stepping pH = 9.17
- [1:03:05] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:03:10] Stepping pH = 9.37
- [1:03:25] Stirrer speed set to 0
- [1:03:38] Datapoint id 62 collected
- [1:03:38] Charge balance equation is out by -91.0%
- [1:03:38] Stirrer speed set to 55
- [1:03:43] pH 9.43 -> 9.63
- [1:03:43] Using cautious pH adjust
- [1:03:43] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:03:48] Stepping pH = 9.45
- [1:03:48] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [1:03:53] Stepping pH = 9.64
- [1:04:08] Stirrer speed set to 0
- [1:04:22] Datapoint id 63 collected
- [1:04:22] Charge balance equation is out by -92.6%
- [1:04:22] Stirrer speed set to 55
- [1:04:27] pH 9.71 -> 9.91
- [1:04:27] Using cautious pH adjust
- [1:04:27] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:04:32] Stepping pH = 9.75
- [1:04:32] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [1:04:37] Stepping pH = 9.91
- [1:04:53] Stirrer speed set to 0
- [1:05:11] Datapoint id 64 collected
- [1:05:11] Charge balance equation is out by -37.0%
- [1:05:11] Stirrer speed set to 55
- [1:05:16] pH 9.95 -> 10.05
- [1:05:16] Using cautious pH adjust
- [1:05:16] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:05:21] Stepping pH = 9.95
- [1:05:21] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [1:05:26] Stepping pH = 10.07
- [1:05:42] Stirrer speed set to 0
- [1:05:52] Datapoint id 65 collected
- [1:05:52] Charge balance equation is out by -89.5%
- [1:05:52] Titration 3 of 3
- [1:05:52] Adding initial titrants



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [1:05:52] Automatically add 0.30000 mL of Octanol
- [1:05:59] Dispensed 0.300000 mL of Octanol
- [1:05:59] Stirrer speed set to 10
- [1:06:00] Stirrer speed set to 60
- [1:06:00] Iterative adjust 10.11 -> 2.00
- [1:06:00] pH 10.11 -> 2.00
- [1:06:03] Dispensed 0.100000 mL of Acid (0.5 M HCI)
- [1:06:08] pH 2.07 -> 2.00
- [1:06:08] Dispensed 0.011642 mL of Acid (0.5 M HCI)
- [1:06:58] Stirrer speed set to 0
- [1:07:08] Datapoint id 66 collected
- [1:07:08] Stirrer speed set to 60
- [1:07:14] pH 1.98 -> 2.18
- [1:07:14] Using cautious pH adjust
- 1.07.14] Osing Cautious pri aujust
- [1:07:14] Dispensed 0.014887 mL of Base (0.5 M KOH)
- [1:07:19] Stepping pH = 2.06
- [1:07:20] Dispensed 0.012512 mL of Base (0.5 M KOH)
- [1:07:25] Stepping pH = 2.15
- [1:07:25] Dispensed 0.003175 mL of Base (0.5 M KOH)
- [1:07:30] Stepping pH = 2.18
- [1:07:45] Stirrer speed set to 0
- [1:07:58] Datapoint id 67 collected
- [1:07:58] Charge balance equation is out by -2.6%
- [1:07:58] Stirrer speed set to 60
- [1:08:03] pH 2.19 -> 2.39
- [1:08:03] Using charge balance adjust
- [1:08:03] Dispensed 0.019120 mL of Base (0.5 M KOH)
- 1:08:23 Stirrer speed set to 0
- [1:08:33] Datapoint id 68 collected
- [1:08:33] Charge balance equation is out by -2.3%
- [1:08:33] Stirrer speed set to 60
- [1:08:39] pH 2.39 -> 2.59
- [1:08:39] Using charge balance adjust
- [1:08:39] Dispensed 0.012982 mL of Base (0.5 M KOH)
- [1:08:59] Stirrer speed set to 0
- [1:09:09] Datapoint id 69 collected
- [1:09:09] Charge balance equation is out by -8.3%
- [1:09:09] Stirrer speed set to 60
- [1:09:14] pH 2.58 -> 2.78
- [1:09:14] Using charge balance adjust
- [1:09:15] Dispensed 0.009266 mL of Base (0.5 M KOH)
- [1:09:35] Stirrer speed set to 0
- [1:09:45] Datapoint id 70 collected
- [1:09:45] Charge balance equation is out by 2.2%
- [1:09:45] Stirrer speed set to 60
- [1:09:50] pH 2.79 -> 2.99
- [1:09:50] Using charge balance adjust
- [1:09:50] Dispensed 0.006609 mL of Base (0.5 M KOH)
- [1:10:10] Stirrer speed set to 0
- [1:10:21] Datapoint id 71 collected
- [1:10:21] Charge balance equation is out by -5.3%
- [1:10:21] Stirrer speed set to 60
- [1:10:26] pH 2.98 -> 3.18
- [1:10:26] Using charge balance adjust
- [1:10:26] Dispensed 0.004986 mL of Base (0.5 M KOH)
- [1:10:46] Stirrer speed set to 0
- [1:10:57] Datapoint id 72 collected
- [1:10:57] Charge balance equation is out by -2.9%
- [1:10:57] Stirrer speed set to 60



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [1:11:02] pH 3.18 -> 3.38
- [1:11:02] Using charge balance adjust
- [1:11:02] Dispensed 0.003857 mL of Base (0.5 M KOH)
- [1:11:22] Stirrer speed set to 0
- [1:11:32] Datapoint id 73 collected
- [1:11:32] Charge balance equation is out by 9.5%
- [1:11:32] Stirrer speed set to 60
- [1:11:37] pH 3.41 -> 3.61
- [1:11:37] Using charge balance adjust
- [1:11:38] Dispensed 0.003104 mL of Base (0.5 M KOH)
- [1:11:58] Stirrer speed set to 0
- [1:12:08] Datapoint id 74 collected
- [1:12:08] Charge balance equation is out by 0.7%
- [1:12:08] Stirrer speed set to 60
- [1:12:13] pH 3.61 -> 3.81
- [1:12:13] Using charge balance adjust
- [1:12:13] Dispensed 0.002775 mL of Base (0.5 M KOH)
- [1:12:33] Stirrer speed set to 0
- [1:12:44] Datapoint id 75 collected
- [1:12:44] Charge balance equation is out by 12.4%
- [1:12:44] Stirrer speed set to 60
- [1:12:49] pH 3.83 -> 4.03
- [1:12:49] Using charge balance adjust
- [1:12:49] Dispensed 0.002658 mL of Base (0.5 M KOH)
- [1:13:09] Stirrer speed set to 0
- [1:13:27] Datapoint id 76 collected
- [1:13:27] Charge balance equation is out by 13.3%
- [1:13:27] Stirrer speed set to 60
- [1:13:33] pH 4.05 -> 4.25
- [1:13:33] Using charge balance adjust
- [1:13:33] Dispensed 0.002540 mL of Base (0.5 M KOH)
- [1:13:53] Stirrer speed set to 0
- [1:14:03] Datapoint id 77 collected
- [1:14:03] Charge balance equation is out by -4.8%
- [1:14:03] Stirrer speed set to 60
- [1:14:08] pH 4.24 -> 4.44
- [1:14:08] Using charge balance adjust
- [1:14:08] Dispensed 0.002375 mL of Base (0.5 M KOH)
- [1:14:28] Stirrer speed set to 0
- [1:14:38] Datapoint id 78 collected
- [1:14:38] Charge balance equation is out by 0.7%
- [1:14:38] Stirrer speed set to 60
- [1:14:44] pH 4.44 -> 4.64
- [1:14:44] Using charge balance adjust
- [1:14:44] Dispensed 0.002046 mL of Base (0.5 M KOH)
- [1:15:04] Stirrer speed set to 0
- [1:15:14] Datapoint id 79 collected
- [1:15:14] Charge balance equation is out by 0.9%
- [1:15:14] Stirrer speed set to 60
- [1:15:19] pH 4.64 -> 4.84
- [1:15:19] Using charge balance adjust
- [1:15:19] Dispensed 0.001646 mL of Base (0.5 M KOH)
- [1:15:39] Stirrer speed set to 0
- [1:16:00] Datapoint id 80 collected
- [1:16:00] Charge balance equation is out by 5.4%
- [1:16:00] Stirrer speed set to 60
- [1:16:05] pH 4.85 -> 5.05
- [1:16:05] Using charge balance adjust
- [1:16:05] Dispensed 0.001223 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [1:16:25] Stirrer speed set to 0
- [1:16:44] Datapoint id 81 collected
- [1:16:44] Charge balance equation is out by 18.8%
- [1:16:44] Stirrer speed set to 60
- [1:16:49] pH 5.09 -> 5.29
- [1:16:49] Using cautious pH adjust
- [1:16:49] Dispensed 0.000400 mL of Base (0.5 M KOH)
- [1:16:55] Stepping pH = 5.16
- [1:16:55] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [1:17:00] Stepping pH = 5.28
- [1:17:00] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:17:05] Stepping pH = 5.31
- 1.17.00] Otopping pri = 0.01
- [1:17:20] Stirrer speed set to 0
- [1:17:31] Datapoint id 82 collected
- [1:17:31] Charge balance equation is out by -9.2%
- [1:17:31] Stirrer speed set to 60
- [1:17:36] pH 5.32 -> 5.52
- [1:17:36] Using charge balance adjust
- [1:17:36] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [1:17:56] Stirrer speed set to 0
- [1:18:08] Datapoint id 83 collected
- [1:18:08] Charge balance equation is out by 10.6%
- [1:18:08] Stirrer speed set to 60
- [1:18:13] pH 5.55 -> 5.75
- [1:18:13] Using charge balance adjust
- [1:18:13] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [1:18:33] Stirrer speed set to 0
- [1:18:51] Datapoint id 84 collected
- [1:18:51] Charge balance equation is out by 9.2%
- [1:18:51] Stirrer speed set to 60
- [1:18:56] pH 5.77 -> 5.97
- [1:18:56] Using charge balance adjust
- [1:18:57] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [1:19:17] Stirrer speed set to 0
- [1:19:32] Datapoint id 85 collected [1:19:32] Charge balance equation is out by -5.5%
- [1:19:32] Stirrer speed set to 60
- [1:19:37] pH 5.96 -> 6.16
- [1:19:37] Using charge balance adjust
- [1:19:37] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:19:57] Stirrer speed set to 0
- [1:20:27] Datapoint id 86 collected
- [1:20:27] Charge balance equation is out by 5.3%
- [1:20:27] Stirrer speed set to 60
- [1:20:32] pH 6.18 -> 6.38
- [1:20:32] Using charge balance adjust
- [1:20:32] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:20:53] Stirrer speed set to 0
- [1:21:52] Datapoint id 87 collected
- [1:21:52] Charge balance equation is out by 37.0%
- [1:21:52] Stirrer speed set to 60
- [1:21:57] pH 6.45 -> 6.65
- [1:21:57] Using cautious pH adjust
- [1:21:57] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:22:02] Stepping pH = 6.47
- [1:22:02] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:22:08] Stepping pH = 6.91
- [1:22:23] Stirrer speed set to 0
- [1:23:23] Datapoint id 88 collected
- Reported at: 3/6/2018 1:24:22 PM



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [1:23:23] Charge balance equation is out by -88.9%
- [1:23:23] Stirrer speed set to 60
- [1:23:28] pH 7.01 -> 7.21
- [1:23:28] Using cautious pH adjust
- [1:23:28] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:23:33] Stepping pH = 7.02
- [1:23:33] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:23:38] Stepping pH = 7.16
- [1:23:38] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:23:43] Stepping pH = 7.31
- [1:23:58] Stirrer speed set to 0
- [1:24:58] Datapoint id 89 collected
- [1:24:58] Charge balance equation is out by -137.7%
- [1:24:58] Stirrer speed set to 60
- [1:25:03] pH 7.32 -> 7.52
- [1:25:03] Using cautious pH adjust
- [1:25:03] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:25:09] Stepping pH = 7.42
- [1:25:09] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:25:14] Stepping pH = 7.52
- [1:25:29] Stirrer speed set to 0
- [1:26:29] Datapoint id 90 collected
- [1:26:29] Charge balance equation is out by -46.2%
- [1:26:29] Stirrer speed set to 60
- [1:26:34] pH 7.52 -> 7.72
- [1:26:34] Using cautious pH adjust
- [1:26:34] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:26:39] Stepping pH = 7.60
- [1:26:39] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:26:44] Stepping pH = 7.66
- [1:26:44] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:26:49] Stepping pH = 7.76
- [1:27:04] Stirrer speed set to 0
- [1:28:04] Datapoint id 91 collected
- [1:28:04] Charge balance equation is out by -201.7%
- [1:28:04] Stirrer speed set to 60
- [1:28:10] pH 7.83 -> 8.03
- [1:28:10] Using cautious pH adjust
- [1:28:10] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:28:15] Stepping pH = 7.86
- [1:28:15] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:28:20] Stepping pH = 7.99
- [1:28:20] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:28:25] Stepping pH = 8.13
- [1:28:40] Stirrer speed set to 0
- [1:29:40] Datapoint id 92 collected
- [1:29:40] Charge balance equation is out by -321.0%
- [1:29:40] Stirrer speed set to 60
- [1:29:45] pH 8.07 -> 8.27
- [1:29:45] Using cautious pH adjust
- [1:29:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:29:50] Stepping pH = 8.07
- [1:29:50] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:29:55] Stepping pH = 8.24
- [1:29:56] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:30:01] Stepping pH = 8.36
- [1:30:16] Stirrer speed set to 0
- [1:31:16] Datapoint id 93 collected
- [1:31:16] Charge balance equation is out by -434.1%
- Reported at: 3/6/2018 1:24:22 PM



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

- [1:31:16] Stirrer speed set to 60
- [1:31:21] pH 8.37 -> 8.57
- [1:31:21] Using cautious pH adjust
- [1:31:21] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:31:26] Stepping pH = 8.37
- [1:31:26] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:31:31] Stepping pH = 8.37
- [1:31:31] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [1:31:36] Stepping pH = 8.72
- [1:31:51] Stirrer speed set to 0
- [1:32:12] Datapoint id 94 collected
- [1:32:12] Charge balance equation is out by -893.2%
- [1:32:12] Stirrer speed set to 60
- [1:32:18] pH 8.83 -> 9.03
- [1:32:18] Using cautious pH adjust
- [1:32:18] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:32:23] Stepping pH = 8.84
- [1:32:23] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:32:28] Stepping pH = 8.87
- [1:32:28] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [1:32:33] Stepping pH = 9.03
- [1:32:48] Stirrer speed set to 0
- [1:33:16] Datapoint id 95 collected
- [1:33:16] Charge balance equation is out by -410.9%
- [1:33:16] Stirrer speed set to 60
- [1:33:21] pH 9.09 -> 9.29
- [1:33:21] Using cautious pH adjust
- [1:33:21] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:33:26] Stepping pH = 9.10
- [1:33:26] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [1:33:31] Stepping pH = 9.14
- [1:33:31] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [1:33:36] Stepping pH = 9.32
- [1:33:51] Stirrer speed set to 0
- [1:34:06] Datapoint id 96 collected
- [1:34:06] Charge balance equation is out by -341.4%
- [1:34:06] Stirrer speed set to 60
- [1:34:11] pH 9.38 -> 9.58
- [1:34:11] Using cautious pH adjust
- [1:34:11] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:34:16] Stepping pH = 9.39
- [1:34:16] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [1:34:21] Stepping pH = 9.52
- [1:34:21] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:34:26] Stepping pH = 9.60
- [1:34:42] Stirrer speed set to 0
- [1:34:54] Datapoint id 97 collected
- [1:34:54] Charge balance equation is out by -153.1%
- [1:34:54] Stirrer speed set to 60
- [1:34:59] pH 9.62 -> 9.82
- [1:34:59] Using cautious pH adjust
- [1:34:59] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [1:35:04] Stepping pH = 9.66
- [1:35:04] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [1:35:10] Stepping pH = 9.78
- [1:35:10] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:35:15] Stepping pH = 9.83
- [1:35:30] Stirrer speed set to 0
- [1:35:51] Datapoint id 98 collected

Experiment Log



Sample name: M15_octanol Experiment start time: 3/3/2018 8:26:22 AM

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-03006 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180302_exp29_logP_T3-2\18C-03006_M15_octanol_pH-metric high logP.t3r

Experiment Log (continued)

[1:35:51] Charge balance equation is out by -75.4%

[1:35:51] Stirrer speed set to 60

[1:35:57] pH 9.84 -> 10.04

[1:35:57] Using cautious pH adjust

[1:35:57] Dispensed 0.000188 mL of Base (0.5 M KOH)

[1:36:02] Stepping pH = 9.90

[1:36:02] Dispensed 0.000235 mL of Base (0.5 M KOH)

[1:36:07] Stepping pH = 10.00

[1:36:07] Dispensed 0.000094 mL of Base (0.5 M KOH)

[1:36:12] Stepping pH = 10.04

[1:36:27] Stirrer speed set to 0

[1:36:39] Datapoint id 99 collected

[1:36:39] Charge balance equation is out by -35.9%

[1:36:39] Argon flow rate set to 0

[1:36:43] Titrator arm moved over Titration position

Reported at: 3/6/2018 1:24:22 PM