

18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

### pH-metric Result

logP (XH +) 0.59 ±0.04 (n=50) logP (neutral X) 4.03 ±0.01 (n=50)

#### 18C-03012 Points 1 to 26

M02 octanol concentration factor 0.955 Carbonate 0.2235 mM Acidity error -1.42826 mM

#### 18C-03012 Points 27 to 50

M02\_octanol concentration factor 0.892 Carbonate 0.1033 mM -1.59633 mM Acidity error

#### 18C-03012 Points 51 to 75

M02 octanol concentration factor 1.218 Carbonate 0.1516 mM Acidity error -1.31894 mM

### Warnings and errors

Errors None Warnings None

рΗ

# Sample logD and percent species

M02\_octanol M02\_octanol M02\_octanol

•	_	<b>—</b>	<b>-</b>	<b>—</b>	<b>—</b>	
	logD	M02_octanolH				
1.000	0.69	17.03 %	0.00 %	65.96 %	17.01 %	
1.200	0.74	15.49 %	0.00 %	59.99 %	24.51 %	Stomach pH
2.000	1.14	6.73 %	0.01 %	26.06 %	67.20 %	
3.000	2.01	0.95 %	0.01 %	3.70 %	95.34 %	
4.000	2.96	0.10 %	0.01 %	0.39 %	99.51 %	
5.000	3.71	0.01 %	0.01 %	0.04 %	99.94 %	
6.000	3.99	0.00 %	0.01 %	0.00 %	99.99 %	
6.500	4.01	0.00 %	0.01 %	0.00 %	99.99 %	
7.000	4.02	0.00 %	0.01 %	0.00 %	99.99 %	
7.400	4.03	0.00 %	0.01 %	0.00 %	99.99 %	Blood pH
8.000	4.03	0.00 %	0.01 %	0.00 %	99.99 %	
9.000	4.03	0.00 %	0.01 %	0.00 %	99.99 %	
10.000	4.03	0.00 %	0.01 %	0.00 %	99.99 %	
11.000	4.03	0.00 %	0.01 %	0.00 %	99.99 %	
12.000	4.03	0.00 %	0.01 %	0.00 %	99.99 %	

M02\_octanol

M02\_octanol Comment



Filename:

Sample name: M02\_octanol Assay name: Assay ID:

pH-metric high logP


18C-03012

Experiment start time: 3/3/2018 3:08:32 PM

Analyst: **Dorothy Levorse** 

Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r






Sample name: M02\_octanol Experiment start time: 3/3/2018 3:08:32 PM pH-metric high logP Analyst: Assay name: **Dorothy Levorse** Assay ID:

18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

# **Graphs** (continued)





Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 1 of 3 18C-03012 Points 1 to 26

### Overall results

RMSD 0.169
Average ionic strength 0.158 M
Average temperature 24.9°C
Partition ratio 0.0123 : 1

Analyte concentration range 4430.5 µM to 4567.1 µM

Total points considered 19 of 26

### Warnings and errors

Errors None

Warnings Excessive acidity error present

### Four-Plus parameters

 Alpha
 0.111
 3/3/2018 3:08:32 PM
 C:\Sirius\_T3\HCl18C02.t3r

 S
 0.9988
 3/3/2018 3:08:32 PM
 C:\Sirius\_T3\HCl18C02.t3r

 jH
 1.0
 3/3/2018 3:08:32 PM
 C:\Sirius\_T3\HCl18C02.t3r


 jOH
 -0.8
 3/3/2018 3:08:32 PM
 C:\Sirius\_T3\HCl18C02.t3r


### Titrants

#### Sample

M02\_octanol concentration factor 0.955
Base pKa 1 5.03
logP (XH +) 1.12
logP (neutral X) 4.07

### Sample graphs





Analyst:

Experiment start time: 3/3/2018 3:08:32 PM

**Dorothy Levorse** 

T312060




Assay ID:

Sample name: M02\_octanol Assay name: pH-metric high logP

Instrument ID: 18C-03012

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

### Sample graphs (continued)



### Sample logD and percent species


рН	M02_octanol logD	M02_octanol M02_octanolH	M02_octanol M02_octanol	M02_octanol M02_octanolH*	M02_octanol M02_octanol*	
1.000	1.16	84.96 %	0.01 %	13.89 %	1.14 %	
1.200	1.18	84.39 %	0.01 %	13.80 %	1.79 %	Stomach pH
2.000	1.38	77.00 %	0.07 %	12.59 %	10.33 %	•
3.000	2.08	39.77 %	0.37 %	6.50 %	53.36 %	
4.000	3.01	6.81 %	0.64 %	1.11 %	91.44 %	
5.000	3.75	0.73 %	0.68 %	0.12 %	98.46 %	
6.000	4.02	0.07 %	0.69 %	0.01 %	99.22 %	
6.500	4.05	0.02 %	0.69 %	0.00 %	99.28 %	
7.000	4.06	0.01 %	0.69 %	0.00 %	99.30 %	
7.400	4.07	0.00 %	0.69 %	0.00 %	99.31 %	Blood pH
8.000	4.07	0.00 %	0.69 %	0.00 %	99.31 %	•
9.000	4.07	0.00 %	0.69 %	0.00 %	99.31 %	
10.000	4.07	0.00 %	0.69 %	0.00 %	99.31 %	
11.000	4.07	0.00 %	0.69 %	0.00 %	99.31 %	
12.000	4.07	0.00 %	0.69 %	0.00 %	99.31 %	

### Carbonate and acidity



Carbonate 0.223 mM Acidity error -1.428 mM

# Other graphs





Assay ID: Filename:

Sample name: M02\_octanol Assay name:

pH-metric high logP


18C-03012

Experiment start time: 3/3/2018 3:08:32 PM Analyst: **Dorothy Levorse** 

Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Sample name: M02\_octanol Experiment start time: 3/3/2018 3:08:32 PM Analyst: Assay name: pH-metric high logP **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-03012

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 2 of 3 18C-03012 Points 27 to 50

#### Overall results

RMSD 0.235 Average ionic strength 0.163 M 25.0°C Average temperature Partition ratio 0.0290:1

Analyte concentration range 4080.4 µM to 4209.4 µM

Total points considered 18 of 24

### Warnings and errors

Errors None

Warnings Excessive acidity error present

#### Four-Plus parameters

Alpha 0.111 3/3/2018 3:08:32 PM C:\Sirius\_T3\HCl18C02.t3r 0.9988 3/3/2018 3:08:32 PM C:\Sirius\_T3\HCl18C02.t3r S jΗ 1.0 3/3/2018 3:08:32 PM C:\Sirius\_T3\HCl18C02.t3r jOH -0.83/3/2018 3:08:32 PM C:\Sirius\_T3\HCl18C02.t3r


# Titrants


0.50 M HCI 0.999058 3/3/2018 3:08:32 PM C:\Sirius T3\HCI18C02.t3r 0.50 M KOH 0.999845 3/3/2018 3:08:32 PM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M02 octanol concentration factor 0.892 Base pKa 1 5.03 logP(XH +)1.47 logP (neutral X) 4.27

#### Sample graphs







Assay ID:

Filename:

Sample name: M02\_octanol Assay name:


pH-metric high logP 18C-03012

Experiment start time: 3/3/2018 3:08:32 PM Analyst: **Dorothy Levorse** 

Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

### Sample graphs (continued)



### Sample logD and percent species


M02_octanol	M02_octanol	_	_	_	Comment
logD	M02_octanolH	M02_octanol	M02_octanolH*	M02_octanol*	
1.50	52.20 %	0.00 %	45.17 %	2.62 %	
1.51	51.41 %	0.01 %	44.49 %	4.09 %	Stomach pH
1.67	42.22 %	0.04 %	36.54 %	21.21 %	
2.30	14.50 %	0.14 %	12.55 %	72.82 %	
3.21	1.92 %	0.18 %	1.66 %	96.25 %	
3.95	0.20 %	0.18 %	0.17 %	99.45 %	
4.22	0.02 %	0.19 %	0.02 %	99.78 %	
4.25	0.01 %	0.19 %	0.01 %	99.80 %	
4.26	0.00 %	0.19 %	0.00 %	99.81 %	
4.27	0.00 %	0.19 %	0.00 %	99.81 %	Blood pH
4.27	0.00 %	0.19 %	0.00 %	99.81 %	
4.27	0.00 %	0.19 %	0.00 %	99.81 %	
4.27	0.00 %	0.19 %	0.00 %	99.81 %	
4.27	0.00 %	0.19 %	0.00 %	99.81 %	
4.27	0.00 %	0.19 %	0.00 %	99.81 %	
	1.50 1.51 1.67 2.30 3.21 3.95 4.22 4.25 4.26 4.27 4.27 4.27 4.27	logD         M02_octanolH           1.50         52.20 %           1.51         51.41 %           1.67         42.22 %           2.30         14.50 %           3.21         1.92 %           3.95         0.20 %           4.22         0.02 %           4.25         0.01 %           4.26         0.00 %           4.27         0.00 %           4.27         0.00 %           4.27         0.00 %           4.27         0.00 %           4.27         0.00 %           4.27         0.00 %           4.27         0.00 %           4.27         0.00 %	logD         M02_octanolH         M02_octanol           1.50         52.20 %         0.00 %           1.51         51.41 %         0.01 %           1.67         42.22 %         0.04 %           2.30         14.50 %         0.14 %           3.21         1.92 %         0.18 %           3.95         0.20 %         0.18 %           4.22         0.02 %         0.19 %           4.25         0.01 %         0.19 %           4.26         0.00 %         0.19 %           4.27         0.00 %         0.19 %           4.27         0.00 %         0.19 %           4.27         0.00 %         0.19 %           4.27         0.00 %         0.19 %           4.27         0.00 %         0.19 %           4.27         0.00 %         0.19 %	logD         M02_octanolH         M02_octanol         M02_octanolH*           1.50         52.20 %         0.00 %         45.17 %           1.51         51.41 %         0.01 %         44.49 %           1.67         42.22 %         0.04 %         36.54 %           2.30         14.50 %         0.14 %         12.55 %           3.21         1.92 %         0.18 %         0.17 %           3.95         0.20 %         0.18 %         0.17 %           4.22         0.02 %         0.19 %         0.02 %           4.25         0.01 %         0.19 %         0.00 %           4.26         0.00 %         0.19 %         0.00 %           4.27         0.00 %         0.19 %         0.00 %           4.27         0.00 %         0.19 %         0.00 %           4.27         0.00 %         0.19 %         0.00 %           4.27         0.00 %         0.19 %         0.00 %           4.27         0.00 %         0.19 %         0.00 %           4.27         0.00 %         0.19 %         0.00 %           4.27         0.00 %         0.19 %         0.00 %	logD         M02_octanolH         M02_octanol         M02_octanolH*         M02_octanol*           1.50         52.20 %         0.00 %         45.17 %         2.62 %           1.51         51.41 %         0.01 %         44.49 %         4.09 %           1.67         42.22 %         0.04 %         36.54 %         21.21 %           2.30         14.50 %         0.14 %         12.55 %         72.82 %           3.21         1.92 %         0.18 %         1.66 %         96.25 %           3.95         0.20 %         0.18 %         0.17 %         99.45 %           4.22         0.02 %         0.19 %         0.02 %         99.78 %           4.25         0.01 %         0.19 %         0.01 %         99.80 %           4.26         0.00 %         0.19 %         0.00 %         99.81 %           4.27         0.00 %         0.19 %         0.00 %         99.81 %           4.27         0.00 %         0.19 %         0.00 %         99.81 %           4.27         0.00 %         0.19 %         0.00 %         99.81 %           4.27         0.00 %         0.19 %         0.00 %         99.81 %           4.27         0.00 %         0.19 %

### Carbonate and acidity



Carbonate 0.103 mM Acidity error -1.596 mM

### Other graphs





Assay ID: Filename:

Sample name: M02\_octanol Assay name:


pH-metric high logP

18C-03012

Experiment start time: 3/3/2018 3:08:32 PM **Dorothy Levorse** Analyst:

Instrument ID: T312060

# Other graphs (continued)





Sample name: M02\_octanol Experiment start time: 3/3/2018 3:08:32 PM Analyst: Assay name: pH-metric high logP **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-03012

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

### pH-metric high logP Titration 3 of 3 18C-03012 Points 51 to 75

#### Overall results

RMSD 0.321 Average ionic strength 0.169 M 25.0°C Average temperature Partition ratio 0.1643:1

Analyte concentration range 3372.3 µM to 3470.6 µM

Total points considered 20 of 25

### Warnings and errors

Errors None

Warnings Excessive acidity error present

#### Four-Plus parameters

Alpha 0.111 3/3/2018 3:08:32 PM C:\Sirius\_T3\HCl18C02.t3r 0.9988 3/3/2018 3:08:32 PM C:\Sirius\_T3\HCl18C02.t3r S jΗ 1.0 3/3/2018 3:08:32 PM C:\Sirius\_T3\HCl18C02.t3r jOH -0.83/3/2018 3:08:32 PM C:\Sirius\_T3\HCl18C02.t3r


### Titrants


0.50 M HCI 0.999058 3/3/2018 3:08:32 PM C:\Sirius T3\HCI18C02.t3r 0.50 M KOH 0.999845 3/3/2018 3:08:32 PM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M02\_octanol concentration factor 1.218 Base pKa 1 5.03 logP(XH +)1.33 logP (neutral X) 4.47

### Sample graphs








Assay ID:

Sample name: M02\_octanol Experiment start time: 3/3/2018 3:08:32 PM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

**18C-03012** Instrument ID: **T312060** 

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

# Sample graphs (continued)




### Sample logD and percent species

рН	M02_octanol logD	M02_octanol M02_octanolH	M02_octanol M02_octanol	M02_octanol M02_octanolH*	M02_octanol M02_octanol*	Comment
1.000	1.39	20.05 %	0.00 %	<del>7</del> 0.82 %	9.14 %	
1.200	1.41	19.03 %	0.00 %	67.22 %	13.74 %	Stomach pH
2.000	1.69	11.00 %	0.01 %	38.86 %	50.13 %	•
3.000	2.47	2.00 %	0.02 %	7.05 %	90.94 %	
4.000	3.41	0.22 %	0.02 %	0.77 %	99.00 %	
5.000	4.16	0.02 %	0.02 %	0.08 %	99.88 %	
6.000	4.43	0.00 %	0.02 %	0.01 %	99.97 %	
6.500	4.46	0.00 %	0.02 %	0.00 %	99.98 %	
7.000	4.47	0.00 %	0.02 %	0.00 %	99.98 %	
7.400	4.47	0.00 %	0.02 %	0.00 %	99.98 %	Blood pH
8.000	4.47	0.00 %	0.02 %	0.00 %	99.98 %	·
9.000	4.47	0.00 %	0.02 %	0.00 %	99.98 %	
10.000	4.47	0.00 %	0.02 %	0.00 %	99.98 %	
11.000	4.47	0.00 %	0.02 %	0.00 %	99.98 %	
12.000	4.47	0.00 %	0.02 %	0.00 %	99.98 %	

### Carbonate and acidity



# Other graphs





Assay ID: Filename:

Sample name: M02\_octanol Assay name:


pH-metric high logP

18C-03012

Experiment start time: 3/3/2018 3:08:32 PM Analyst: **Dorothy Levorse** 

Instrument ID: T312060

# Other graphs (continued)





Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

### Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name	M02_octanol	12/6/2017 4:20:03 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.002080 g	3/2/2018 5:10:06 PM	User entered value
ormula weight	289.26 g/mol	12/6/2017 4:20:03 PM	User entered value
Solubility	Unknown		Default value
Nolecular weight	289.26	12/6/2017 4:20:03 PM	User entered value
ndividual pKa ionic environments	No		Default value
lumber of pKas	1	12/6/2017 4:20:03 PM	User entered value
Sample is a	Base	12/6/2017 4:20:03 PM	User entered value
Ka 1	5.03	12/6/2017 4:20:03 PM	User entered value
ogp (XH +)	1.32	3/2/2018 3:38:13 PM	User entered value
ogP (neutral X)	4.10	3/2/2018 3:38:07 PM	User entered value

#### **Events**

<b>Time</b> 5:59.7	Event Initial pH = 9.32	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
8:59.2	Data point 1	1.50000 mL	0.05336 mL	0.00113 mL	0.01999 mL	2.071	0.00065	0.01575	0.00025	10.0 s
9:45.5	Data point 2		0.05336 mL					0.13165	0.00060	10.0 s
10:21.2	Data point 3		0.05336 mL					0.41455	0.00027	10.5 s
10:57.3	Data point 4	1.50000 mL	0.05336 mL	0.03057 mL	0.01999 mL	2.669	-0.00471	0.73726	0.00027	10.0 s
11:32.8	Data point 5		0.05336 mL						0.00069	10.0 s
	Data point 6		0.05336 mL						0.00067	
12:43.7			0.05336 mL						0.00062	
13:19.6			0.05336 mL						0.00099	
	Data point 9		0.05336 mL						0.00092	
	Data point 10		0.05336 mL						0.00098	
15:44.0			0.05336 mL						0.00092	
16:24.9	Data point 12		0.05336 mL						0.00095	
	Data point 13		0.05336 mL						0.00098	
18:25.4			0.05336 mL						0.00091	
19:50.9	Data point 15		0.05336 mL							Timed out at 59.5 s
21:21.5	Data point 16	1.50000 mL	0.05336 mL	0.04871 mL	0.01999 mL	4.975	-0.02384	0.95065	0.00121	Timed out at 59.5 s
22:57.2	Data point 17	1.50000 mL	0.05336 mL	0.04892 mL	0.01999 mL	5.443	-0.01935	0.93030	0.00099	56.0 s
24:23.8	Data point 18	1.50000 mL	0.05336 mL	0.04906 mL	0.01999 mL	5.658	-0.01989	0.98952	0.00099	46.5 s
25:40.9	Data point 19	1.50000 mL	0.05336 mL	0.04920 mL	0.01999 mL	6.311	-0.02340	0.95639	0.00118	Timed out at 59.5 s
27:11.5	Data point 20	1.50000 mL	0.05336 mL	0.04934 mL	0.01999 mL	7.143	-0.06220	0.99094	0.00309	Timed out at 59.5 s
28:41.9	Data point 21	1.50000 mL	0.05336 mL	0.04939 mL	0.01999 mL	7.536	-0.06100	0.98752	0.00303	Timed out at 59.5 s
30:07.3	Data point 22	1.50000 mL	0.05336 mL	0.04941 mL	0.01999 mL	7.646	-0.04673	0.98883	0.00232	Timed out at 59.5 s
31:37.8	Data point 23	1.50000 mL	0.05336 mL	0.04946 mL	0.01999 mL	8.129	-0.05480	0.98555	0.00273	Timed out at 59.5 s
33:08.3	Data point 24	1.50000 mL	0.05336 mL	0.04951 mL	0.01999 mL	8.458	-0.03117	0.96668	0.00157	Timed out at 59.5 s
34:43.9	Data point 25	1.50000 mL	0.05336 mL	0.04958 mL	0.01999 mL	8.734	-0.01915	0.91772	0.00099	
36:12.1	Data point 26		0.05336 mL						0.00096	
			0.10861 mL						0.00076	
38:27.6			0.10861 mL						0.00071	
39:03.3	Data point 29		0.10861 mL						0.00042	
39:38.9	Data point 30		0.10861 mL						0.00012	
40:15.0	Data point 31	1.50000 mL	0.10861 mL	0.09087 mL	0.05000 mL	2.809	-0.00504	0.92869	0.00026	10.0 s
40:50.5	Data point 32	1.50000 mL	0.10861 mL	0.09497 mL	0.05000 mL	3.014	-0.00539	0.66005	0.00033	10.0 s

Reported at: 3/6/2018 2:16:39 PM



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

# **Events (continued)**

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
41:26.0	Data point 33	1.50000 mL	0.10861 mL	0.09795 mL	0.05000 mL	3.211	-0.00572	0.74071	0.00033	10.0 s
42:01.6	Data point 34		0.10861 mL					0.48807	0.00096	
42:52.5	Data point 35		0.10861 mL					0.19061	0.00065	
43:28.4	Data point 36		0.10861 mL					0.86567	0.00080	
44:09.9	Data point 37		0.10861 mL					0.92718	0.00100	
44:47.3	Data point 38		0.10861 mL					0.92624	0.00093	
45:25.7	Data point 39		0.10861 mL					0.89982	0.00100	14.0 s
46:05.1	Data point 40		0.10861 mL					0.91342	0.00096	
46:50.7	Data point 41		0.10861 mL					0.98349	0.00099	
47:44.2	Data point 42		0.10861 mL					0.91486	0.00093	
49:02.1	Data point 43		0.10861 mL					0.90215	0.00091	
50:27.2	Data point 44		0.10861 mL					0.99241	0.00033	Timed out
30.27.2	Data point 44	1.50000 IIIL	0.10001 IIIL	0.10327 IIIL	0.03000 IIIL	0.074	-0.03440	0.99241	0.00171	at 59.5 s
51:57.7	Data point 45	1.50000 mL	0.10861 mL	0.10336 mL	0.05000 mL	7.259	-0.06329	0.99069	0.00314	Timed out
53:33.3	Data point 46	1.50000 mL	0.10861 mL	0.10343 mL	0.05000 mL	7.768	-0.05820	0.98124	0.00290	at 59.5 s Timed out
	p									at 59.5 s
55:03.8	Data point 47	1.50000 mL	0.10861 mL	0.10348 mL	0.05000 mL	8.292	-0.03946	0.97572	0.00197	Timed out
	= p									at 59.5 s
56:34.3	Data point 48	1.50000 mL	0.10861 mL	0.10353 mL	0.05000 mL	8.557	-0.01940	0.93969	0.00099	57.5 s
58:07.5	Data point 49		0.10861 mL					0.06822	0.00099	
59:13.7	Data point 50		0.10861 mL					0.91843	0.00094	
1:00:42.8	Data point 51		0.16823 mL					0.10876	0.00043	
1:01:29.1	Data point 52		0.16823 mL					0.64737	0.00099	
	Data point 53		0.16823 mL					0.16357	0.00099	
	Data point 54		0.16823 mL					0.00759	0.00012	
	Data point 55		0.16823 mL					0.75034	0.00055	
1:04:00.1	Data point 56		0.16823 mL					0.31455	0.00100	
1:04:37.1	Data point 57		0.16823 mL					0.56497	0.00093	
1:05:27.9			0.16823 mL					0.15736	0.00029	
	Data point 59		0.16823 mL					0.87321	0.00039	
1:06:44.4			0.16823 mL					0.94298	0.00051	10.5 s
	Data point 61		0.16823 mL					0.75314	0.00070	10.0 s
1:08:00.9			0.16823 mL					0.01015	0.00096	
1:08:39.8			0.16823 mL					0.39352	0.00084	
1:09:22.3	•							0.88825	0.00089	
1:10:06.3				0.16345 mL				0.91841	0.00003	
	Data point 66		0.16823 mL					0.01098	0.00030	
	Data point 67		0.16823 mL							Timed out
	·									at 59.5 s
1:14:03.5	Data point 68	1.50000 mL	0.16823 mL	0.16362 mL	0.30000 mL	6.671	-0.04606	0.97751	0.00230	Timed out at 59.5 s
1:15:39.1	Data point 69	1.50000 mL	0.16823 mL	0.16369 mL	0.30000 mL	7.163	-0.06931	0.98791	0.00344	Timed out at 59.5 s
1:17:09.6	Data point 70	1.50000 mL	0.16823 mL	0.16376 mL	0.30000 mL	7.616	-0.08650	0.98058	0.00432	Timed out at 59.5 s
1:18:40.1	Data point 71	1.50000 mL	0.16823 mL	0.16381 mL	0.30000 mL	8.042	-0.06118	0.94917	0.00310	Timed out
1:20:15.7	Data point 72	1.50000 mL	0.16823 mL	0.16388 mL	0.30000 mL	8.385	-0.04572	0.94304	0.00233	at 59.5 s Timed out
1.21.51 2	Data point 73	1 50000 ml	0.16823 mL	0 16395 ml	0.30000 ml	8 716	_N N1062	n 98797	0.00098	at 59.5 s
	Data point 73		0.16823 mL						0.00098	
	Data point 75		0.16823 mL						0.00071	
	Assay volumes					5.002	0.01243	0.70 <b>2</b> 07	5.00034	10.0 3



Assay ID: Filename:	18C-03012	ehtap\20180302_exp2			「312060 N. pH-metric high logi
		ептарки 100302_ехри	9_l0gF_13-2\10\	5-03012_WI02_OCIANG	n_pri-inetric mgn logi
Assay Setti	ings				
Setting		Value	Original Value	Date/Time changed	Imported from
General Setti	ngs				
Analyst name		Dorothy Levorse			
Standard Exp	eriment Settings				
Number of titra	ations	3			
Minimum pH		2.000			
Maximum pH		9.000			
pH step betwe	en points of	0.200			
Minimum titrar	nt addition	0.00002 mL			
Maximum titra	nt addition	0.10000 mL			
Argon flow rate	е	100%			
Start titration u	ısing	Cautious pH adjust			
Advanced Ge	eneral Settings	. ,			
Detect turbidity	y using	None			
Collect turbidit		No			
Collect UV spe		No			
Stir after titran	t addition for	5 seconds			
For titrant add	ition, stir at	10%			
Titrant Pre-Do	ose				
Titrant pre-dos	se	None			
Assay Mediul	m				
ISA water volu	me	1.50 mL			
Water added		Automatic			
Partition solve	nt type	Octanol			
Partition volum		0.020 mL			
Partition solve		Automatic			
	addition, stir for	1 seconds			
Sample Sonic					

Sonicate Yes Adjust pH for sonication No

Sonicate for 120 seconds

After sonication stir for 5 seconds

Sample Dissolution

Perform a dissolution stage Yes

Adjust and hold pH for dissolution To start pH Stir to dissolve for 120 seconds

For dissolution, stir at 10%

Carbonate purge

Perform a carbonate purge No

Temperature Control

Wait for temperature Yes 25.0°C Required start temperature Acceptable deviation 0.5°C Time to wait 60 seconds

Stir speed of 50%

Titration 1

Titrate from Low to high pH

Adjust to start pH Yes After pH adjust stir for 30 seconds

Stir to allow partitioning for 15 seconds Stirrer speed for partitioning 50%

Titration 2

Titrate from Low to high pH Add additional water 0.00 mL

Additional partition solvent volume 0.030 mL Automatic Additional partition solvent added After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds

Stirrer speed for partitioning 55%

Reported at: 3/6/2018 2:16:39 PM

Page 15 of 30



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

# Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3		_	_	
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.250 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00100 dpH/dt			
Stability timeout after	60 seconds			

### Calibration Settings

Value	Date/Time changed	Imported from
0.111	3/3/2018 3:08:32 PM	C:\Sirius_T3\HCl18C02.t3r
0.9988	3/3/2018 3:08:32 PM	C:\Sirius_T3\HCl18C02.t3r
1.0	3/3/2018 3:08:32 PM	C:\Sirius_T3\HCl18C02.t3r
-0.8	3/3/2018 3:08:32 PM	C:\Sirius_T3\HCl18C02.t3r
1.000	3/3/2018 3:08:32 PM	C:\Sirius_T3\KOH18B27.t3r
0.999	3/3/2018 3:08:32 PM	C:\Sirius_T3\HCl18C02.t3r
	0.111 0.9988 1.0 -0.8 1.000	0.111 3/3/2018 3:08:32 PM 0.9988 3/3/2018 3:08:32 PM 1.0 3/3/2018 3:08:32 PM -0.8 3/3/2018 3:08:32 PM 1.000 3/3/2018 3:08:32 PM

# Instrument Settings

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module Dispenser 0 Syringe volume Firmware version	Water 2.5 mL 1.2.1(r2)	T3DM1200361	3/31/2009 5:24:52 AM 3/31/2009 5:25:05 AM
Titrant Dispenser 2 Syringe volume Firmware version	Water (0.15 M KCI) Acid 0.5 mL 1.2.1(r2)	02-06-2018	2/27/2018 10:05:59 AM 3/31/2009 5:25:11 AM
Titrant Dispenser 1 Syringe volume Firmware version	Acid (0.5 M HCI) Base 0.5 mL 1.2.1(r2)	02-27-2018	2/27/2018 10:27:22 AM 3/31/2009 5:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version Distribution valve 5	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2) Distribution Valve	9/22/2017	2/27/2018 10:21:22 AM 3/31/2009 5:26:24 AM 3/31/2009 5:28:19 AM
Firmware version Port A Port B Dispenser 3 Syringe volume	1.1.3 Methanol (80%, 0.15 M KCI) Cyclohexane Buffer 0.5 mL	09-26-17 11-01-17	2/7/2018 9:42:01 AM 2/27/2018 10:37:57 AM 8/3/2010 5:05:16 AM
Firmware version Titrant Dispenser 6	1.2.1(r2) Dodecane Octanol	2018/01/31	2/28/2018 10:18:04 AM 10/22/2010 10:52:43 AM

Reported at: 3/6/2018 2:16:39 PM



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

# Instrument Settings (continued)

gc (co			
Setting	Value	Batch Id	Install date
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)	04 24 2040	2/27/2019 0:50:25 AM
Titrant Titrator	Octanol	01-31-2018 T3TM1200161	2/27/2018 9:59:35 AM 3/31/2009 5:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2	1311111200101	3/3 1/2009 3.24.17 AW
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0923	1/23/2018 2:01:00 PM
E0 calibration	+6.14 mV		3/3/2018 3:09:00 PM
Filling solution	3M KCI	KCL097	3/2/2018 9:43:24 AM
Liquids			
Wash 1	50% IPA:50% Water		3/2/2018 9:45:12 AM
Wash 2	0.5% Trition X-100 in H20		3/2/2018 9:45:15 AM
Buffer position 1	pH7 Wash		3/2/2018 9:45:18 AM
Buffer position 2	pH 7		3/2/2018 9:45:21 AM
Storage position Wash water	6.7e+003 mL	02-27-2018	3/2/2018 9:44:44 AM 2/27/2018 9:54:39 AM
Waste	8.8e+003 mL	02-21-2010	11/28/2017 10:36:29 AM
Temperature controller	0.00 · 000 IIIE		8/5/2010 6:35:13 AM
Turbidity detector			3/31/2009 5:24:45 AM
Spectrometer		074811	11/23/2010 11:22:28 AM
Dip probe		10196	
Wavelength coefficient A0	183.333		
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	120:41:49		11/23/2010 11:22:28 AM
Calibrated on	2/27/2018 10:40:38 AM		
Integration time	40		
Scans averaged	10	T2AL 1200245	11/10/2015 0:24:12 AM
Autoloader Left-right axis firmware version	1.17 AI1DI2DO2 Stepper 2	T3AL1200345	11/10/2015 9:34:13 AM
Front-back axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Configuration	3 - 3		
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume Flowing wash stir duration	20.0 mL		
Flowing wash stir speed	5 s 30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed E0 calibration reading stir speed	30% 0%		
Lo calibration reading still speed	U / U		



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

### Instrument Settings (continued)

Setting	Value	Batch Id	I Install date
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		
· •			

### Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance		0.100
Turbidity probe threshold	50.00	50.00
· ·		

### **Experiment Log**

[2:37] Air gap created for Water (0.15 M KC	<b>(</b> [(
[2:38] Air gap created for Acid (0.5 M HCI)	,

- [2:38] Air gap created for Base (0.5 M KOH)
- [2:39] Air gap released for Water (0.15 M KCI)
- [2:42] Titrator arm moved over Titration position
- [2:42] Titration 1 of 3
- [2:42] Adding initial titrants
- [2:42] Automatically add 1.50000 mL of water
- [3:08] Dispensed 1.500000 mL of Water (0.15 M KCI)
- [3:12] Titrator arm moved over Drain
- [5:53] Titrator arm moved to Titration position
- [5:53] Argon flow rate set to 100
- [5:53] Stirrer speed set to 10
- 5:58 Automatically add 0.02000 mL of Octanol
- [5:59] Dispensed 0.019991 mL of Octanol
- 6:00 Initial pH = 9.32
- [6:00] Iterative adjust 9.32 -> 2.00
- [6:00] pH 9.32 -> 2.00
- [6:02] Air gap released for Acid (0.5 M HCl)
- [6:02] Dispensed 0.053363 mL of Acid (0.5 M HCI)
- [6:08] Holding pH 2.00
- [8:08] Stirrer speed set to 0
- [8:08] Stirrer speed set to 50
- [8:08] Iterative adjust 1.99 -> 2.00
- [8:08] pH 1.99 -> 2.00
- [8:08] Air gap released for Base (0.5 M KOH)
- [8:09] Dispensed 0.001129 mL of Base (0.5 M KOH)
- [8:59] Stirrer speed set to 0
- [9:09] Datapoint id 1 collected
- [9:09] Stirrer speed set to 50
- [9:15] pH 2.08 -> 2.28
- [9:15] Using cautious pH adjust
- [9:15] Dispensed 0.006726 mL of Base (0.5 M KOH)
- [9:20] Stepping pH = 2.16
- [9:20] Dispensed 0.005738 mL of Base (0.5 M KOH)
- [9:25] Stepping pH = 2.25
- [9:25] Dispensed 0.001834 mL of Base (0.5 M KOH)
- [9:31] Stepping pH = 2.28
- [9:46] Stirrer speed set to 0
- [9:56] Datapoint id 2 collected
- [9:56] Charge balance equation is out by -6.2%
- [9:56] Stirrer speed set to 50



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [10:01] pH 2.29 -> 2.49
- [10:01] Using charge balance adjust
- [10:01] Dispensed 0.008866 mL of Base (0.5 M KOH)
- [10:21] Stirrer speed set to 0
- [10:32] Datapoint id 3 collected
- [10:32] Charge balance equation is out by -5.7%
- [10:32] Stirrer speed set to 50
- [10:37] pH 2.48 -> 2.68
- [10:37] Using charge balance adjust
- [10:37] Dispensed 0.006279 mL of Base (0.5 M KOH)
- [10:57] Stirrer speed set to 0
- [11:07] Datapoint id 4 collected
- [11:07] Charge balance equation is out by -6.6%
- [11:07] Stirrer speed set to 50
- [11:13] pH 2.68 -> 2.88
- [11:13] Using charge balance adjust
- [11:13] Dispensed 0.004704 mL of Base (0.5 M KOH)
- [11:33] Stirrer speed set to 0
- [11:43] Datapoint id 5 collected
- [11:43] Charge balance equation is out by -0.6%
- [11:43] Stirrer speed set to 50
- [11:48] pH 2.88 -> 3.08
- [11:48] Using charge balance adjust
- [11:48] Dispensed 0.003622 mL of Base (0.5 M KOH)
- [12:08] Stirrer speed set to 0
- [12:18] Datapoint id 6 collected
- [12:18] Charge balance equation is out by 0.1%
- [12:18] Stirrer speed set to 50
- [12:24] pH 3.09 -> 3.29
- [12:24] Using charge balance adjust
- [12:24] Dispensed 0.002846 mL of Base (0.5 M KOH)
- [12:44] Stirrer speed set to 0
- [12:54] Datapoint id 7 collected
- [12:54] Charge balance equation is out by -4.2%
- [12:54] Stirrer speed set to 50
- [12:59] pH 3.28 -> 3.48
- [12:59] Using charge balance adjust
- [13:01] Dispensed 0.002211 mL of Base (0.5 M KOH)
- [13:20] Stirrer speed set to 0
- [13:31] Datapoint id 8 collected
- [13:31] Charge balance equation is out by 21.6%
- [13:31] Stirrer speed set to 50
- [13:36] pH 3.53 -> 3.73
- [13:36] Using cautious pH adjust
- [13:36] Dispensed 0.000776 mL of Base (0.5 M KOH)
- [13:41] Stepping pH = 3.67
- [13:41] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [13:46] Stepping pH = 3.70
- [13:46] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [13:51] Stepping pH = 3.72
- [14:07] Stirrer speed set to 0
- [14:18] Datapoint id 9 collected
- [14:18] Charge balance equation is out by 19.9%
- [14:18] Stirrer speed set to 50
- [14:23] pH 3.69 -> 3.89
- [14:23] Using cautious pH adjust
- [14:23] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [14:28] Stepping pH = 3.82
- [14:29] Dispensed 0.000235 mL of Base (0.5 M KOH)



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [14:34] Stepping pH = 3.86
- [14:34] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [14:39] Stepping pH = 3.88
- [14:54] Stirrer speed set to 0
- [15:08] Datapoint id 10 collected
- [15:08] Charge balance equation is out by 17.0%
- [15:08] Stirrer speed set to 50
- [15:14] pH 3.87 -> 4.07
- [15:14] Using cautious pH adjust
- [15:14] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [15:19] Stepping pH = 4.00
- [15:19] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [15:24] Stepping pH = 4.02
- [15:24] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [15:29] Stepping pH = 4.09
- [15:44] Stirrer speed set to 0
- [16:00] Datapoint id 11 collected
- [16:00] Charge balance equation is out by 4.0%
- [16:00] Stirrer speed set to 50
- [16:05] pH 4.05 -> 4.25
- [16:05] Using charge balance adjust
- [16:05] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [16:25] Stirrer speed set to 0
- [16:53] Datapoint id 12 collected
- [16:53] Charge balance equation is out by -3.7%
- [16:53] Stirrer speed set to 50
- [16:58] pH 4.27 -> 4.47
- [16:58] Using charge balance adjust
- [16:58] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [17:19] Stirrer speed set to 0
- [17:55] Datapoint id 13 collected
- [17:55] Charge balance equation is out by -71.4%
- [17:55] Stirrer speed set to 50
- [18:00] pH 4.37 -> 4.57
- [18:00] Using cautious pH adjust
- [18:00] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [18:05] Stepping pH = 4.42
- [18:05] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [18:11] Stepping pH = 4.67
- [18:26] Stirrer speed set to 0
- [19:21] Datapoint id 14 collected
- [19:21] Charge balance equation is out by -17.2%
- [19:21] Stirrer speed set to 50
- [19:26] pH 4.58 -> 4.78
- [19:26] Using cautious pH adjust
- [19:26] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [19:31] Stepping pH = 4.63
- [19:31] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [19:36] Stepping pH = 4.85
- [19:51] Stirrer speed set to 0
- [20:51] Datapoint id 15 collected
- [20:51] Charge balance equation is out by -20.7%
- [20:51] Stirrer speed set to 50
- [20:56] pH 4.78 -> 4.98
- [20:56] Using cautious pH adjust
- [20:56] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [21:01] Stepping pH = 4.82
- [21:01] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [21:07] Stepping pH = 5.08



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [21:22] Stirrer speed set to 0
- [22:22] Datapoint id 16 collected
- [22:22] Charge balance equation is out by -35.6%
- [22:22] Stirrer speed set to 50
- [22:27] pH 5.00 -> 5.20
- [22:27] Using cautious pH adjust
- [22:27] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [22:32] Stepping pH = 5.11
- [22:32] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [22:37] Stepping pH = 5.09
- [22:37] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [22:42] Stepping pH = 5.51
- [22:57] Stirrer speed set to 0
- [23:53] Datapoint id 17 collected
- [23:53] Charge balance equation is out by -140.4%
- [23:53] Stirrer speed set to 50
- [23:59] pH 5.46 -> 5.66
- [23:59] Using cautious pH adjust
- [23:59] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [24:04] Stepping pH = 5.45
- [24:04] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [24:09] Stepping pH = 5.76
- [24:24] Stirrer speed set to 0
- [25:11] Datapoint id 18 collected
- [25:11] Charge balance equation is out by -211.0%
- [25:11] Stirrer speed set to 50
- [25:16] pH 5.83 -> 6.03
- [25:16] Using cautious pH adjust
- [25:16] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [25:21] Stepping pH = 5.82
- [25:21] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [25:26] Stepping pH = 6.32
- [25:41] Stirrer speed set to 0
- [26:41] Datapoint id 19 collected
- [26:41] Charge balance equation is out by -217.8%
- [26:41] Stirrer speed set to 50
- [26:46] pH 6.25 -> 6.45
- [26:46] Using cautious pH adjust
- [26:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [26:51] Stepping pH = 6.23
- [26:51] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [26:57] Stepping pH = 7.03
- [27:12] Stirrer speed set to 0
- [28:12] Datapoint id 20 collected
- [28:12] Charge balance equation is out by -227.6%
- [28:12] Stirrer speed set to 50
- [28:17] pH 7.10 -> 7.30
- [28:17] Using cautious pH adjust
- [28:17] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [28:22] Stepping pH = 7.16
- [28:22] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [28:27] Stepping pH = 7.53
- [28:42] Stirrer speed set to 0
- [29:42] Datapoint id 21 collected
- [29:42] Charge balance equation is out by -180.3%
- [29:42] Stirrer speed set to 50
- [29:47] pH 7.45 -> 7.65
- [29:47] Using cautious pH adjust
- [29:47] Dispensed 0.000024 mL of Base (0.5 M KOH)



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [29:52] Stepping pH = 7.65
- [30:07] Stirrer speed set to 0
- [31:08] Datapoint id 22 collected
- [31:08] Charge balance equation is out by -150.3%
- [31:08] Stirrer speed set to 50
- [31:13] pH 7.60 -> 7.80
- [31:13] Using cautious pH adjust
- [31:13] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [31:18] Stepping pH = 7.75
- [31:18] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [31:23] Stepping pH = 8.16
- [31:38] Stirrer speed set to 0
- [32:38] Datapoint id 23 collected
- [32:38] Charge balance equation is out by -495.8%
- [32:38] Stirrer speed set to 50
- [32:43] pH 8.06 -> 8.26
- [32:43] Using cautious pH adjust
- [32:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [32:48] Stepping pH = 8.02
- [32:48] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [32:53] Stepping pH = 8.38
- [33:08] Stirrer speed set to 0
- [34:09] Datapoint id 24 collected
- [34:09] Charge balance equation is out by -485.1%
- [34:09] Stirrer speed set to 50
- [34:14] pH 8.42 -> 8.62
- [34:14] Using cautious pH adjust
- [34:14] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [34:19] Stepping pH = 8.39
- [34:19] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [34:24] Stepping pH = 8.57
- [34:24] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [34:29] Stepping pH = 8.82
- [34:44] Stirrer speed set to 0
- [35:42] Datapoint id 25 collected
- [35:42] Charge balance equation is out by -501.8%
- [35:42] Stirrer speed set to 50
- [35:47] pH 8.72 -> 8.92
- [35:47] Using cautious pH adjust
- [35:47] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [35:52] Stepping pH = 8.69
- [35:52] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [35:57] Stepping pH = 9.04
- [36:12] Stirrer speed set to 0
- [36:48] Datapoint id 26 collected
- [36:48] Charge balance equation is out by -260.2%
- [36:48] Titration 2 of 3
- [36:48] Adding initial titrants
- [36:48] Automatically add 0.03000 mL of Octanol
- [36:49] Dispensed 0.030009 mL of Octanol
- [36:49] Stirrer speed set to 10
- [36:50] Stirrer speed set to 55
- [36:50] Iterative adjust 9.10 -> 2.00
- [36:50] pH 9.10 -> 2.00
- [36:51] Dispensed 0.055245 mL of Acid (0.5 M HCI)
- [37:41] Stirrer speed set to 0
- [37:52] Datapoint id 27 collected
- [37:52] Stirrer speed set to 55 [37:57] pH 1.99 -> 2.19
- \_\_\_\_\_



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [37:57] Using cautious pH adjust
- [37:57] Dispensed 0.008913 mL of Base (0.5 M KOH)
- [38:02] Stepping pH = 2.08
- [38:02] Dispensed 0.006279 mL of Base (0.5 M KOH)
- [38:07] Stepping pH = 2.17
- [38:08] Dispensed 0.001176 mL of Base (0.5 M KOH)
- [38:13] Stepping pH = 2.19
- [38:28] Stirrer speed set to 0
- [38:38] Datapoint id 28 collected
- [38:38] Charge balance equation is out by 8.2%
- [38:38] Stirrer speed set to 55
- [38:43] pH 2.21 -> 2.41
- [38:43] Using charge balance adjust
- [38:43] Dispensed 0.011406 mL of Base (0.5 M KOH)
- [39:04] Stirrer speed set to 0
- [39:14] Datapoint id 29 collected
- [39:14] Charge balance equation is out by -5.6%
- [39:14] Stirrer speed set to 55
- [39:19] pH 2.40 -> 2.60
- [39:19] Using charge balance adjust
- [39:19] Dispensed 0.007949 mL of Base (0.5 M KOH)
- [39:39] Stirrer speed set to 0
- [39:50] Datapoint id 30 collected
- [39:50] Charge balance equation is out by 8.8%
- [39:50] Stirrer speed set to 55
- [39:55] pH 2.63 -> 2.83
- [39:55] Using charge balance adjust
- [39:55] Dispensed 0.005480 mL of Base (0.5 M KOH)
- [40:15] Stirrer speed set to 0
- [40:25] Datapoint id 31 collected
- [40:25] Charge balance equation is out by -9.5%
- [40:25] Stirrer speed set to 55
- [40:30] pH 2.82 -> 3.02
- [40:30] Using charge balance adjust
- [40:31] Dispensed 0.004092 mL of Base (0.5 M KOH)
- [40:51] Stirrer speed set to 0
- [41:01] Datapoint id 32 collected
- [41:01] Charge balance equation is out by -0.9%
- [41:01] Stirrer speed set to 55
- [41:06] pH 3.02 -> 3.22
- [41:06] Using charge balance adjust
- [41:06] Dispensed 0.002987 mL of Base (0.5 M KOH)
- [41:26] Stirrer speed set to 0
- [41:36] Datapoint id 33 collected
- [41:36] Charge balance equation is out by -4.3%
- [41:36] Stirrer speed set to 55
- [41:41] pH 3.22 -> 3.42
- [41:41] Using charge balance adjust
- [41:42] Dispensed 0.002140 mL of Base (0.5 M KOH)
- [42:02] Stirrer speed set to 0
- [42:12] Datapoint id 34 collected
- [42:12] Charge balance equation is out by 40.8%
- [42:12] Stirrer speed set to 55
- [42:17] pH 3.51 -> 3.71
- [42:17] Using cautious pH adjust
- [42:17] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [42:22] Stepping pH = 3.63
- [42:22] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [42:27] Stepping pH = 3.69



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [42:27] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [42:32] Stepping pH = 3.69
- [42:32] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [42:38] Stepping pH = 3.74
- [42:53] Stirrer speed set to 0
- [43:03] Datapoint id 35 collected
- [43:03] Charge balance equation is out by 12.2%
- [43:03] Stirrer speed set to 55
- [43:08] pH 3.74 -> 3.94
- [43:08] Using charge balance adjust
- [43:08] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [43:29] Stirrer speed set to 0
- [43:40] Datapoint id 36 collected
- [43:40] Charge balance equation is out by 38.4%
- [43:40] Stirrer speed set to 55
- [43:45] pH 4.03 -> 4.23
- [43:45] Using cautious pH adjust
- [43:45] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [43:50] Stepping pH = 4.11
- [43:50] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [43:55] Stepping pH = 4.23
- [44:10] Stirrer speed set to 0
- [44:22] Datapoint id 37 collected
- [44:22] Charge balance equation is out by 7.1%
- [44:22] Stirrer speed set to 55
- [44:27] pH 4.25 -> 4.45
- [44:27] Using charge balance adjust
- [44:27] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [44:48] Stirrer speed set to 0
- [45:01] Datapoint id 38 collected
- [45:01] Charge balance equation is out by -1.3%
- [45:01] Stirrer speed set to 55
- [45:06] pH 4.47 -> 4.67
- [45:06] Using charge balance adjust
- [45:06] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [45:26] Stirrer speed set to 0
- [45:40] Datapoint id 39 collected
- [45:40] Charge balance equation is out by -8.1%
- [45:40] Stirrer speed set to 55
- [45:45] pH 4.69 -> 4.89
- [45:45] Using charge balance adjust
- [45:45] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [46:05] Stirrer speed set to 0
- [46:20] Datapoint id 40 collected
- [46:20] Charge balance equation is out by -45.8%
- [46:20] Stirrer speed set to 55
- [46:25] pH 4.84 -> 5.04
- [46:25] Using cautious pH adjust
- [46:26] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [46:31] Stepping pH = 4.87
- [46:31] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [46:36] Stepping pH = 5.07
- [46:51] Stirrer speed set to 0
- [47:14] Datapoint id 41 collected
- [47:14] Charge balance equation is out by -57.4%
- [47:14] Stirrer speed set to 55
- [47:19] pH 5.15 -> 5.35
- [47:19] Using cautious pH adjust
- [47:19] Dispensed 0.000024 mL of Base (0.5 M KOH)



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [47:24] Stepping pH = 5.16
- [47:24] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [47:29] Stepping pH = 5.51
- [47:44] Stirrer speed set to 0
- [48:27] Datapoint id 42 collected
- [48:27] Charge balance equation is out by -88.9%
- [48:27] Stirrer speed set to 55
- [48:32] pH 5.71 -> 5.91
- [48:32] Using cautious pH adjust
- [48:32] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [48:37] Stepping pH = 5.76
- [48:37] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [48:42] Stepping pH = 5.83
- [48:42] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [48:47] Stepping pH = 6.06
- [49:02] Stirrer speed set to 0
- [49:57] Datapoint id 43 collected
- [49:57] Charge balance equation is out by -92.8%
- [49:57] Stirrer speed set to 55
- [50:02] pH 6.20 -> 6.40
- [50:02] Using cautious pH adjust
- [50:02] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [50:07] Stepping pH = 6.23
- [50:07] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [50:12] Stepping pH = 6.63
- [50:27] Stirrer speed set to 0
- [51:27] Datapoint id 44 collected
- [51:27] Charge balance equation is out by -69.3%
- [51:27] Stirrer speed set to 55
- [51:32] pH 6.63 -> 6.83
- [51:32] Using cautious pH adjust
- [51:33] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [51:38] Stepping pH = 6.63
- [51:38] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [51:43] Stepping pH = 7.07
- [51:58] Stirrer speed set to 0
- [52:58] Datapoint id 45 collected
- [52:58] Charge balance equation is out by -214.8%
- [52:58] Stirrer speed set to 55
- [53:03] pH 7.25 -> 7.45
- [53:03] Using cautious pH adjust
- [53:03] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [53:08] Stepping pH = 7.27
- [53:08] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [53:13] Stepping pH = 7.40
- [53:13] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [53:18] Stepping pH = 7.65
- [53:33] Stirrer speed set to 0
- [54:34] Datapoint id 46 collected
- [54:34] Charge balance equation is out by -403.0%
- [54:34] Stirrer speed set to 55
- [54:39] pH 7.78 -> 7.98
- [54:39] Using cautious pH adjust
- [54:39] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [54:44] Stepping pH = 7.88
- [54:44] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [54:49] Stepping pH = 8.12
- [55:04] Stirrer speed set to 0
- [56:04] Datapoint id 47 collected



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [56:04] Charge balance equation is out by -511.7%
- [56:04] Stirrer speed set to 55
- [56:09] pH 8.34 -> 8.54
- [56:09] Using cautious pH adjust
- [56:09] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [56:14] Stepping pH = 8.38
- [56:14] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [56:19] Stepping pH = 8.53
- [56:35] Stirrer speed set to 0
- [57:32] Datapoint id 48 collected
- [57:32] Charge balance equation is out by -253.0%
- [57:32] Stirrer speed set to 55
- [57:37] pH 8.57 -> 8.77
- [57:37] Using cautious pH adjust
- [57:37] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [57:42] Stepping pH = 8.59
- [57:42] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [57:47] Stepping pH = 8.67
- [57:48] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [57:53] Stepping pH = 8.77
- [58:08] Stirrer speed set to 0
- [58:38] Datapoint id 49 collected
- [58:38] Charge balance equation is out by -278.9%
- [58:38] Stirrer speed set to 55
- [58:43] pH 8.82 -> 9.02
- [58:43] Using cautious pH adjust
- [58:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [58:49] Stepping pH = 8.83
- [58:49] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [58:54] Stepping pH = 8.96
- [58:54] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [58:59] Stepping pH = 9.04
- [59:14] Stirrer speed set to 0
- [59:39] Datapoint id 50 collected
- [59:39] Charge balance equation is out by -176.9%
- [59:39] Titration 3 of 3
- [59:39] Adding initial titrants
- [59:39] Automatically add 0.25000 mL of Octanol
- [59:45] Dispensed 0.250000 mL of Octanol
- [59:45] Stirrer speed set to 10
- [59:46] Stirrer speed set to 60
- [59:46] Iterative adjust 9.02 -> 2.00
- [59:46] pH 9.02 -> 2.00
- [59:47] Dispensed 0.056656 mL of Acid (0.5 M HCI)
- [59:53] pH 2.03 -> 2.00
- [59:53] Dispensed 0.002963 mL of Acid (0.5 M HCI)
- [1:00:43] Stirrer speed set to 0
- [1:00:53] Datapoint id 51 collected
- [1:00:53] Stirrer speed set to 60
- [1:00:58] pH 1.99 -> 2.19
- [1:00:58] Using cautious pH adjust
- [1:00:58] Dispensed 0.009784 mL of Base (0.5 M KOH)
- [1:01:04] Stepping pH = 2.06
- [1:01:04] Dispensed 0.008984 mL of Base (0.5 M KOH)
- [1:01:09] Stepping pH = 2.16
- [1:01:09] Dispensed 0.002117 mL of Base (0.5 M KOH)
- [1:01:14] Stepping pH = 2.19
- [1:01:29] Stirrer speed set to 0
- [1:01:45] Datapoint id 52 collected
- Reported at: 3/6/2018 2:16:39 PM



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [1:01:45] Charge balance equation is out by -6.7%
- [1:01:45] Stirrer speed set to 60
- [1:01:50] pH 2.19 -> 2.39
- [1:01:50] Using charge balance adjust
- [1:01:50] Dispensed 0.013053 mL of Base (0.5 M KOH)
- [1:02:11] Stirrer speed set to 0
- [1:02:24] Datapoint id 53 collected
- [1:02:24] Charge balance equation is out by -6.1%
- [1:02:24] Stirrer speed set to 60
- [1:02:29] pH 2.38 -> 2.58
- [1:02:29] Using charge balance adjust
- [1:02:29] Dispensed 0.008984 mL of Base (0.5 M KOH)
- [1:02:49] Stirrer speed set to 0
- [1:02:59] Datapoint id 54 collected
- [1:02:59] Charge balance equation is out by 4.6%
- [1:02:59] Stirrer speed set to 60
- [1:03:04] pH 2.59 -> 2.79
- [1:03:04] Using charge balance adjust
- [1:03:05] Dispensed 0.006021 mL of Base (0.5 M KOH)
- [1:03:25] Stirrer speed set to 0
- [1:03:35] Datapoint id 55 collected
- [1:03:35] Charge balance equation is out by -0.8%
- [1:03:35] Stirrer speed set to 60
- [1:03:40] pH 2.79 -> 2.99
- [1:03:40] Using charge balance adjust
- [1:03:40] Dispensed 0.004116 mL of Base (0.5 M KOH)
- [1:04:00] Stirrer speed set to 0
- [1:04:12] Datapoint id 56 collected
- [1:04:12] Charge balance equation is out by 9.3%
- [1:04:12] Stirrer speed set to 60
- [1:04:17] pH 3.02 -> 3.22
- [1:04:17] Using charge balance adjust
- [1:04:17] Dispensed 0.002658 mL of Base (0.5 M KOH)
- [1:04:37] Stirrer speed set to 0
- [1:04:47] Datapoint id 57 collected [1:04:47] Charge balance equation is out by 29.0%
- [1:04:47] Charge balance equation is [1:04:47] Stirrer speed set to 60
- [1:04:52] pH 3.29 -> 3.49
- [1:04:52] Using cautious pH adjust
- [1:04:52] Dispensed 0.000776 mL of Base (0.5 M KOH)
- [1:04:57] Stepping pH = 3.39
- [1:04:58] Dispensed 0.000494 mL of Base (0.5 M KOH)
- [1:05:03] Stepping pH = 3.47
- [1:05:03] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:05:08] Stepping pH = 3.48
- [1:05:08] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:05:13] Stepping pH = 3.49
- [1:05:28] Stirrer speed set to 0
- [1:05:39] Datapoint id 58 collected
- [1:05:39] Charge balance equation is out by 4.0%
- [1:05:39] Stirrer speed set to 60
- [1:05:44] pH 3.50 -> 3.70
- [1:05:44] Using charge balance adjust
- [1:05:44] Dispensed 0.000988 mL of Base (0.5 M KOH)
- [1:06:04] Stirrer speed set to 0
- [1:06:14] Datapoint id 59 collected
- [1:06:14] Charge balance equation is out by 16.0%
- [1:06:14] Stirrer speed set to 60
- [1:06:19] pH 3.74 -> 3.94



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [1:06:19] Using cautious pH adjust
- [1:06:19] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [1:06:24] Stepping pH = 3.79
- [1:06:24] Dispensed 0.000400 mL of Base (0.5 M KOH)
- [1:06:29] Stepping pH = 3.99 [1:06:45] Stirrer speed set to 0
- [1:06:55] Datapoint id 60 collected
- [1:06:55] Charge balance equation is out by -17.6%
- [1:06:55] Stirrer speed set to 60
- [1:07:00] pH 4.01 -> 4.21
- [1:07:00] Using cautious pH adjust
- [1:07:00] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [1:07:05] Stepping pH = 4.07
- [1:07:05] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:07:11] Stepping pH = 4.21
- [1:07:26] Stirrer speed set to 0
- [1:07:36] Datapoint id 61 collected
- [1:07:36] Charge balance equation is out by -10.9%
- [1:07:36] Stirrer speed set to 60
- [1:07:41] pH 4.22 -> 4.42
- [1:07:41] Using charge balance adjust
- [1:07:41] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:08:01] Stirrer speed set to 0
- [1:08:15] Datapoint id 62 collected
- [1:08:15] Charge balance equation is out by 8.4%
- [1:08:15] Stirrer speed set to 60
- [1:08:20] pH 4.47 -> 4.67
- [1:08:20] Using charge balance adjust
- [1:08:20] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [1:08:40] Stirrer speed set to 0
- [1:08:52] Datapoint id 63 collected
- [1:08:52] Charge balance equation is out by -24.3%
- [1:08:52] Stirrer speed set to 60
- [1:08:57] pH 4.65 -> 4.85
- [1:08:57] Using cautious pH adjust
- [1:08:57] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:09:02] Stepping pH = 4.69
- [1:09:02] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:09:07] Stepping pH = 4.89
- [1:09:22] Stirrer speed set to 0
- [1:09:36] Datapoint id 64 collected
- [1:09:36] Charge balance equation is out by -31.4%
- [1:09:36] Stirrer speed set to 60
- [1:09:41] pH 4.94 -> 5.14
- [1:09:41] Using cautious pH adjust
- [1:09:41] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:09:46] Stepping pH = 4.94
- [1:09:46] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [1:09:51] Stepping pH = 5.39
- [1:10:07] Stirrer speed set to 0
- [1:10:54] Datapoint id 65 collected
- [1:10:54] Charge balance equation is out by -203.1%
- [1:10:54] Stirrer speed set to 60
- [1:10:59] pH 5.56 -> 5.76
- [1:10:59] Using cautious pH adjust
- [1:10:59] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:11:04] Stepping pH = 5.59
- [1:11:04] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:11:10] Stepping pH = 5.78



Assay ID: 18C-03012 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [1:11:25] Stirrer speed set to 0
- [1:12:08] Datapoint id 66 collected
- [1:12:08] Charge balance equation is out by -75.9%
- [1:12:08] Stirrer speed set to 60
- [1:12:13] pH 6.06 -> 6.26
- [1:12:13] Using cautious pH adjust
- [1:12:13] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:12:18] Stepping pH = 6.18
- [1:12:18] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:12:23] Stepping pH = 6.42
- [1:12:38] Stirrer speed set to 0
- [1:13:38] Datapoint id 67 collected
- [1:13:38] Charge balance equation is out by 2.6%
- [1:13:38] Stirrer speed set to 60
- [1:13:43] pH 6.36 -> 6.56
- [1:13:43] Using charge balance adjust
- [1:13:44] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:14:04] Stirrer speed set to 0
- [1:15:04] Datapoint id 68 collected
- [1:15:04] Charge balance equation is out by 56.6%
- [1:15:04] Stirrer speed set to 60
- [1:15:09] pH 6.79 -> 6.99
- [1:15:09] Using cautious pH adjust
- [1:15:09] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:15:14] Stepping pH = 6.87
- [1:15:14] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:15:19] Stepping pH = 6.91
- [1:15:19] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:15:24] Stepping pH = 7.04
- [1:15:39] Stirrer speed set to 0
- [1:16:39] Datapoint id 69 collected
- [1:16:39] Charge balance equation is out by -89.8%
- [1:16:39] Stirrer speed set to 60
- [1:16:44] pH 7.20 -> 7.40
- [1:16:44] Using cautious pH adjust
- [1:16:44] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:16:50] Stepping pH = 7.19
- [1:16:50] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:16:55] Stepping pH = 7.44
- [1:17:10] Stirrer speed set to 0
- [1:18:10] Datapoint id 70 collected
- [1:18:10] Charge balance equation is out by -272.1%
- [1:18:10] Stirrer speed set to 60
- [1:18:15] pH 7.69 -> 7.89
- [1:18:15] Using cautious pH adjust
- [1:18:15] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:18:20] Stepping pH = 7.81
- [1:18:20] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:18:25] Stepping pH = 8.01
- [1:18:40] Stirrer speed set to 0
- [1:19:40] Datapoint id 71 collected
- [1:19:40] Charge balance equation is out by -351.4%
- [1:19:40] Stirrer speed set to 60
- [1:19:45] pH 8.10 -> 8.30
- [1:19:45] Using cautious pH adjust
- [1:19:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:19:51] Stepping pH = 8.16
- [1:19:51] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:19:56] Stepping pH = 8.26



Assay ID:

Sample name: M02\_octanol Experiment start time: 3/3/2018 3:08:32 PM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

**18C-03012** Instrument ID: **T312060** 

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03012\_M02\_octanol\_pH-metric high logP.t3r

- [1:19:56] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:20:01] Stepping pH = 8.39
- [1:20:16] Stirrer speed set to 0
- [1:21:16] Datapoint id 72 collected
- [1:21:16] Charge balance equation is out by -565.1%
- [1:21:16] Stirrer speed set to 60
- [1:21:21] pH 8.46 -> 8.66
- [1:21:21] Using cautious pH adjust
- [1:21:21] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:21:26] Stepping pH = 8.52
- [1:21:26] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:21:31] Stepping pH = 8.61
- [1:21:31] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:21:36] Stepping pH = 8.71
- [1:21:52] Stirrer speed set to 0
- [1:22:25] Datapoint id 73 collected
- [1:22:25] Charge balance equation is out by -283.1%
- [1:22:25] Stirrer speed set to 60
- [1:22:30] pH 8.76 -> 8.96
- [1:22:30] Using cautious pH adjust
- [1:22:30] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:22:35] Stepping pH = 8.80
- [1:22:35] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:22:41] Stepping pH = 8.85
- [1:22:41] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:22:46] Stepping pH = 8.91
- [1:22:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:22:51] Stepping pH = 8.98
- [1:23:06] Stirrer speed set to 0
- [1:23:30] Datapoint id 74 collected
- [1:23:30] Charge balance equation is out by -271.8%
- [1:23:30] Stirrer speed set to 60
- [1:23:36] pH 8.99 -> 9.05
- [1:23:36] Using cautious pH adjust
- [1:23:36] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:23:41] Stepping pH = 9.01
- [1:23:56] Stirrer speed set to 0
- [1:24:09] Datapoint id 75 collected
- [1:24:09] Charge balance equation is out by -70.9%
- [1:24:09] Argon flow rate set to 0
- [1:24:13] Titrator arm moved over Titration position