

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

### pH-metric Result

logP (XH +) 0.44 ±0.03 (n=50) logP (neutral X) 3.29 ±0.01 (n=50)

#### 18B-28011 Points 1 to 22

M07\_octanol concentration factor 1.039
Carbonate 0.1202 mM
Acidity error 0.21327 mM

#### 18B-28011 Points 23 to 45

M07\_octanol concentration factor 1.041
Carbonate 0.1876 mM
Acidity error 0.19535 mM

#### 18B-28011 Points 46 to 79

M07\_octanol concentration factor 0.988
Carbonate 0.1406 mM
Acidity error 2.18588 mM

## Warnings and errors

Errors None Warnings None

# Sample logD and percent species

рН	M07_octanol logD	M07_octanol M07_octanolH	M07_octanol M07_octanol	M07_octanol M07_octanolH*	M07_octanol M07_octanol*	
1.000	0.45	2 <del>6</del> .38 %	0.00 %	<del>7</del> 3.19 %	0.43 %	
1.200	0.45	26.31 %	0.00 %	73.00 %	0.69 %	Stomach pH
2.000	0.47	25.38 %	0.00 %	70.44 %	4.18 %	
3.000	0.65	18.44 %	0.02 %	51.18 %	30.36 %	
4.000	1.28	4.94 %	0.04 %	13.70 %	81.31 %	
5.000	2.19	0.59 %	0.05 %	1.65 %	97.71 %	
6.000	2.95	0.06 %	0.05 %	0.17 %	99.72 %	
6.500	3.15	0.02 %	0.05 %	0.05 %	99.88 %	
7.000	3.24	0.01 %	0.05 %	0.02 %	99.93 %	
7.400	3.27	0.00 %	0.05 %	0.01 %	99.94 %	Blood pH
8.000	3.28	0.00 %	0.05 %	0.00 %	99.95 %	
9.000	3.29	0.00 %	0.05 %	0.00 %	99.95 %	
10.000	3.29	0.00 %	0.05 %	0.00 %	99.95 %	
11.000	3.29	0.00 %	0.05 %	0.00 %	99.95 %	
12.000	3.29	0.00 %	0.05 %	0.00 %	99.95 %	



Sample name: M07\_octanol Assay name:

pH-metric high logP


Assay ID: 18B-28011 Filename:

Experiment start time: 2/28/2018 4:26:20 PM

Analyst: Pion Instrument ID:

T312060 C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r








pH-metric high logP Analyst: **Pion** Assay name: Assay ID: 18B-28011 Instrument ID: T312060 Filename:

 $C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric\ high\ logP.t3r$ 

# **Graphs** (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

## pH-metric high logP Titration 1 of 3 18B-28011 Points 1 to 22

#### Overall results

RMSD 0.037
Average ionic strength 0.157 M
Average temperature 25.0°C
Partition ratio 0.0123 : 1

Analyte concentration range 2362.2 µM to 2434.4 µM

Total points considered 17 of 22

### Warnings and errors

Errors None Warnings None

## Four-Plus parameters


Alpha 0.130 2/28/2018 4:26:20 PM C:\Sirius\_T3\HCl18B27.t3r S 0.9970 2/28/2018 4:26:20 PM C:\Sirius\_T3\HCl18B27.t3r jH 0.8 2/28/2018 4:26:20 PM C:\Sirius\_T3\HCl18B27.t3r jOH -0.4 2/28/2018 4:26:20 PM C:\Sirius\_T3\HCl18B27.t3r


#### **Titrants**

#### Sample

M07\_octanol concentration factor 1.039
Base pKa 1 6.07
logP (XH +) 0.50
logP (neutral X) 3.27

#### Sample graphs







Assay ID:


Filename:

Sample name: M07\_octanol Experiment start time: 2/28/2018 4:26:20 PM

Assay name: pH-metric high logP Analyst: Pion Instrument ID: T312060 18B-28011

C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r


# Sample graphs (continued)




### Sample logD and percent species


рН	M07_octanol	M07_octanol	M07_octanol	M07_octanol	M07_octanol	Comment
	logD	M07_octanolH	M07_octanol	M07_octanolH*	M07_octanol*	
1.000	0.50	96.24 %	0.00 %	3.74 %	0.02 %	
1.200	0.50	96.23 %	0.00 %	3.74 %	0.03 %	Stomach pH
2.000	0.52	96.07 %	0.01 %	3.73 %	0.19 %	
3.000	0.68	94.41 %	0.08 %	3.67 %	1.84 %	
4.000	1.28	80.48 %	0.69 %	3.13 %	15.71 %	
5.000	2.17	32.51 %	2.77 %	1.26 %	63.46 %	
6.000	2.93	4.67 %	3.98 %	0.18 %	91.17 %	
6.500	3.13	1.53 %	4.11 %	0.06 %	94.30 %	
7.000	3.22	0.49 %	4.16 %	0.02 %	95.33 %	
7.400	3.25	0.20 %	4.17 %	0.01 %	95.63 %	Blood pH
8.000	3.27	0.05 %	4.18 %	0.00 %	95.77 %	
9.000	3.27	0.00 %	4.18 %	0.00 %	95.82 %	
10.000	3.27	0.00 %	4.18 %	0.00 %	95.82 %	
11.000	3.27	0.00 %	4.18 %	0.00 %	95.82 %	
12.000	3.27	0.00 %	4.18 %	0.00 %	95.82 %	

# Carbonate and acidity



## Other graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

## pH-metric high logP Titration 2 of 3 18B-28011 Points 23 to 45

#### Overall results

RMSD 0.084
Average ionic strength 0.162 M
Average temperature 25.0°C
Partition ratio 0.0408 : 1

Analyte concentration range 2150.3 µM to 2218.8 µM

Total points considered 16 of 23

### Warnings and errors

Errors None Warnings None

### Four-Plus parameters


Alpha 0.130 2/28/2018 4:26:20 PM C:\Sirius\_T3\HCl18B27.t3r S 0.9970 2/28/2018 4:26:20 PM C:\Sirius\_T3\HCl18B27.t3r jH 0.8 2/28/2018 4:26:20 PM C:\Sirius\_T3\HCl18B27.t3r jOH -0.4 2/28/2018 4:26:20 PM C:\Sirius\_T3\HCl18B27.t3r


#### Titrants

#### Sample

M07\_octanol concentration factor 1.041
Base pKa 1 6.07
logP (XH +) 0.50
logP (neutral X) 3.31

### Sample graphs








Assay ID:

Sample name: M07\_octanol Experiment start time: 2/28/2018 4:26:20 PM

Assay name: pH-metric high logP Analyst: Pion Instrument ID: T312060 18B-28011

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

# Sample graphs (continued)




## Sample logD and percent species

рН	M07_octanol	M07_octanol	M07_octanol	M07_octanol	M07_octanol	Comment
	logD	M07_octanolH	M07_octanol	M07_octanolH*	M07_octanol*	
1.000	0.50	88.51 %	0.00 %	11.43 %	0.06 %	
1.200	0.50	88.48 %	0.00 %	11.42 %	0.10 %	Stomach pH
2.000	0.52	88.01 %	0.01 %	11.36 %	0.63 %	
3.000	0.69	83.26 %	0.07 %	10.75 %	5.92 %	
4.000	1.31	54.08 %	0.46 %	6.98 %	38.48 %	
5.000	2.21	12.01 %	1.02 %	1.55 %	85.42 %	
6.000	2.97	1.37 %	1.16 %	0.18 %	97.29 %	
6.500	3.17	0.44 %	1.18 %	0.06 %	98.33 %	
7.000	3.26	0.14 %	1.18 %	0.02 %	98.66 %	
7.400	3.29	0.06 %	1.18 %	0.01 %	98.76 %	Blood pH
8.000	3.31	0.01 %	1.18 %	0.00 %	98.80 %	
9.000	3.31	0.00 %	1.18 %	0.00 %	98.82 %	
10.000	3.31	0.00 %	1.18 %	0.00 %	98.82 %	
11.000	3.31	0.00 %	1.18 %	0.00 %	98.82 %	
12.000	3.31	0.00 %	1.18 %	0.00 %	98.82 %	

# **Carbonate and acidity**



# Other graphs






Assay name: pH-metric high logP Analyst: Pion Instrument ID: T312060 Assay ID: 18B-28011 Filename:

C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 3 of 3 18B-28011 Points 46 to 79

### Overall results

RMSD 0.443
Average ionic strength 0.169 M
Average temperature 25.0°C
Partition ratio 0.1746 : 1

Analyte concentration range 1775.6 µM to 1831.9 µM

Total points considered 23 of 34

#### Warnings and errors

Errors None

Warnings Excessive acidity error present

### Four-Plus parameters

Alpha	0.130	2/28/2018 4:26:20 PM	C:\Sirius_T3\HCl18B27.t3r
S	0.9970	2/28/2018 4:26:20 PM	C:\Sirius_T3\HCl18B27.t3r
jΗ	8.0	2/28/2018 4:26:20 PM	C:\Sirius_T3\HCl18B27.t3r
jОН	-0.4	2/28/2018 4:26:20 PM	C:\Sirius_T3\HCl18B27.t3r


#### Titrants


0.50 M HCI 0.993513 2/28/2018 4:26:20 PM C:\Sirius\_T3\HCl18B27.t3r 0.50 M KOH 0.999845 2/28/2018 4:26:20 PM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M07\_octanol concentration factor 0.988
Base pKa 1 6.07
logP (XH +) 0.50
logP (neutral X) 3.30

#### Sample graphs







Assay ID:


Filename:

Sample name: M07\_octanol Experiment start time: 2/28/2018 4:26:20 PM

Assay name: pH-metric high logP Analyst: Pion Instrument ID: T312060 18B-28011

C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r


# Sample graphs (continued)




### Sample logD and percent species


рН	M07_octanol	M07_octanol	M07_octanol	_	M07_octanol	Comment
	logD	M07_octanolH	M07_octanol	M07_octanolH*	M07_octanol*	
1.000	0.50	64.31 %	0.00 %	35.50 %	0.19 %	
1.200	0.50	64.24 %	0.00 %	35.46 %	0.30 %	Stomach pH
2.000	0.52	63.22 %	0.01 %	34.90 %	1.88 %	•
3.000	0.69	54.07 %	0.05 %	29.85 %	16.04 %	
4.000	1.30	22.09 %	0.19 %	12.19 %	65.53 %	
5.000	2.20	3.20 %	0.27 %	1.76 %	94.77 %	
6.000	2.96	0.33 %	0.28 %	0.18 %	99.20 %	
6.500	3.16	0.11 %	0.29 %	0.06 %	99.55 %	
7.000	3.25	0.03 %	0.29 %	0.02 %	99.66 %	
7.400	3.28	0.01 %	0.29 %	0.01 %	99.69 %	Blood pH
8.000	3.30	0.00 %	0.29 %	0.00 %	99.71 %	
9.000	3.30	0.00 %	0.29 %	0.00 %	99.71 %	
10.000	3.30	0.00 %	0.29 %	0.00 %	99.71 %	
11.000	3.30	0.00 %	0.29 %	0.00 %	99.71 %	
12.000	3.30	0.00 %	0.29 %	0.00 %	99.71 %	

# Carbonate and acidity



## Other graphs








Assay name: pH-metric high logP Analyst: Pion Instrument ID: T312060 Assay ID: 18B-28011 Filename:

C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

## Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name	M07_octanol	2/27/2018 4:29:24 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.000900 g	2/28/2018 4:24:06 PM	User entered value
Formula weight	235.28 g/mol	2/27/2018 4:29:24 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	235.28	2/27/2018 4:29:24 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	1	2/27/2018 4:29:24 PM	User entered value
Sample is a	Base	2/27/2018 4:29:24 PM	User entered value
pKa 1	6.07	2/27/2018 4:29:24 PM	User entered value
logp (XH +)	0.50	2/28/2018 1:33:04 PM	User entered value
logP (neutral X)	3.44	2/28/2018 1:33:10 PM	User entered value

#### **Events**

Events	•									
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
8:55.6	Initial pH = 6.32									tillio
11:55.1	Data point 1		0.04887 mL				-0.00759	0.43163	0.00057	10.0 s
12:41.2	Data point 2		0.04887 mL				-0.00650	0.83933	0.00035	10.5 s
13:17.4			0.04887 mL				-0.00986		0.00055	10.0 s
	Data point 4		0.04887 mL				-0.00451		0.00027	
	Data point 5		0.04887 mL				-0.00773		0.00063	
	Data point 6		0.04887 mL				-0.00302		0.00032	
	Data point 7	1.50000 mL	0.04887 mL	0.03911 mL	0.01999 mL	3.259	-0.00406	0.30171	0.00036	10.0 s
16:14.6		1.50000 mL	0.04887 mL	0.04012 mL	0.01999 mL	3.431	-0.00336	0.75731	0.00019	10.0 s
	Data point 9	1.50000 mL	0.04887 mL	0.04099 mL	0.01999 mL	3.667	-0.00529	0.80865	0.00029	10.0 s
	Data point 10		0.04887 mL				-0.01649		0.00086	
	Data point 11		0.04887 mL				-0.00742	0.86888	0.00039	10.0 s
	Data point 12		0.04887 mL				-0.01526		0.00084	
	Data point 13		0.04887 mL				-0.01458		0.00081	
20:19.9	Data point 14		0.04887 mL				-0.01834		0.00099	
21:00.3			0.04887 mL				-0.01816		0.00095	
	Data point 16		0.04887 mL				-0.01774		0.00093	
	Data point 17		0.04887 mL				-0.01937		0.00100	
	Data point 18		0.04887 mL				-0.01860		0.00098	
24:36.8			0.04887 mL				-0.01882		0.00095	
25:52.5	Data point 20	1.50000 mL	0.04887 mL	0.04967 mL	0.01999 mL	8.879	-0.02286	0.13576	0.00307	Timed out
										at 59.5 s
	Data point 21		0.04887 mL				-0.01571			16.5 s
28:15.3			0.04887 mL				-0.01833		0.00095	
29:33.1			0.10350 mL				-0.00524		0.00040	
	Data point 24		0.10350 mL				-0.01414	0.83789	0.00076	
30:56.0			0.10350 mL				0.00314	0.09186	0.00051	
	Data point 26		0.10350 mL				-0.00349		0.00023	
	Data point 27		0.10350 mL				-0.00230		0.00015	10.0 s
32:52.8			0.10350 mL				-0.00226		0.00016	
	Data point 29		0.10350 mL				0.00205	0.02071	0.00070	
	Data point 30		0.10350 mL				-0.00570		0.00036	10.0 s
	Data point 31		0.10350 mL				-0.00228		0.00017	10.5 s
35:30.6			0.10350 mL				-0.00812		0.00077	
	Data point 33		0.10350 mL				-0.00941		0.00084	
36:41.4			0.10350 mL				0.00167	0.00821	0.00091	10.0 s
	Data point 35		0.10350 mL				-0.00558		0.00038	
37:52.8			0.10350 mL				-0.00463		0.00036	10.5 s
38:28.8			0.10350 mL				-0.00336	0.10816	0.00050	10.5 s
39:20.2	Data point 38	1.50000 mL	0.10350 mL	0.10237 mL	0.06999 mL	5.067	-0.00648	0.11307	0.00095	10.5 s
11000	D 1 1100	4 50000 '	0.400=0 .	0 4000 .			0 00 100	0.07000	0 00001	44 =

40:06.4 Data point 39 1.50000 mL 0.10350 mL 0.10285 mL 0.06999 mL 5.275 0.00498 0.07222

Reported at: 3/2/2018 1:05:03 PM

0.00091 11.5 s



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

## **Events (continued)**

	•									
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	-	dpH/dt time
40:48.5	Data point 40	1.50000 mL	0.10350 mL	0.10325 mL	0.06999 mL	5.581	-0.01086	0.46682	0.00078	12.0 s
41:36.2	Data point 41	1.50000 mL	0.10350 mL	0.10388 mL	0.06999 mL	6.795	-0.04882	0.99621	0.00241	Timed out
	-									at 59.5 s
43:27.1	Data point 42	1.50000 mL	0.10350 mL	0.10433 mL	0.06999 mL	9.241	-0.00847	0.26017	0.00082	18.5 s
44:16.0	Data point 43	1.50000 mL	0.10350 mL	0.10459 mL	0.06999 mL	9.587	-0.01617	0.74261	0.00093	21.0 s
45:07.4	Data point 44	1.50000 mL	0.10350 mL	0.10501 mL	0.06999 mL	9.949	-0.01841	0.93681	0.00094	13.5 s
45:51.4	Data point 45	1.50000 mL	0.10350 mL	0.10541 mL	0.06999 mL	10.160	-0.00819	0.56417	0.00054	10.0 s
47:05.4	Data point 46	1.50000 mL	0.16268 mL	0.10541 mL	0.31999 mL	1.946	-0.00597	0.52848	0.00041	10.0 s
47:51.7	Data point 47			0.13290 mL			0.00614	0.38583	0.00049	10.0 s
48:37.7	Data point 48			0.14372 mL			-0.00203		0.00016	
49:13.2	Data point 49			0.15122 mL			0.01311	0.84025	0.00071	
49:48.7	Data point 50	1.50000 mL	0.16268 mL	0.15588 mL	0.31999 mL	2.842	0.00403	0.34986	0.00034	10.0 s
50:24.2	Data point 51			0.15887 mL			-0.00378		0.00022	
51:10.1	Data point 52			0.16108 mL			0.00384	0.15940	0.00048	
51:45.5	Data point 53			0.16279 mL			-0.00301		0.00062	
52:20.9	Data point 54			0.16425 mL			0.00201	0.08022	0.00035	
52:56.4	Data point 55			0.16550 mL			-0.00573		0.00031	
53:31.7	Data point 56			0.16651 mL			-0.00755		0.00040	
54:07.0	Data point 57			0.16731 mL			-0.00644		0.00086	
54:42.4	Data point 58			0.16790 mL			-0.01161		0.00080	
55:33.2	Data point 59			0.16858 mL			-0.00119		0.00080	
56:24.4	Data point 60			0.16900 mL			-0.00586		0.00093	
57:05.5	Data point 61			0.16926 mL			0.00326	0.04495	0.00076	
57:47.1	Data point 62			0.16947 mL			-0.00065		0.00070	
58:34.2	Data point 63			0.16964 mL			0.00764	0.15851	0.00095	
59:22.4	Data point 64			0.16980 mL			-0.01710		0.00088	
1:00:33.7	Data point 65			0.16990 mL			-0.01881		0.00095	
1:02:03.8	Data point 66			0.16997 mL			-0.05982			Timed out
	·									at 59.5 s
1:03:34.3	Data point 67	1.50000 mL	0.16268 mL	0.17001 mL	0.31999 mL	7.040	-0.07343	0.98170	0.00366	Timed out at 59.5 s
1:04:59.6	Data point 68	1.50000 mL	0.16268 mL	0.17004 mL	0.31999 mL	7.096	-0.04310	0.98785	0.00214	
	·									at 59.5 s
1:06:30.2	Data point 69	1.50000 mL	0.16268 mL	0.17008 mL	0.31999 mL	7.517	-0.07608	0.99474	0.00377	Timed out
	·									at 59.5 s
1:08:00.6	Data point 70	1.50000 mL	0.16268 mL	0.17013 mL	0.31999 mL	7.886	-0.06164	0.99333	0.00306	Timed out
	·									at 59.5 s
1:09:31.1	Data point 71	1.50000 mL	0.16268 mL	0.17018 mL	0.31999 mL	8.015	-0.04517	0.98471	0.00225	
1:11:11.9	Data point 72	1.50000 mL	0.16268 mL	0.17034 mL	0.31999 mL	8.445	-0.02344	0.95432	0.00119	at 59.5 s Timed out
										at 59.5 s
1:12:47.5	Data point 73	1.50000 mL	0.16268 mL	0.17051 mL	0.31999 mL	8.787	-0.01486	0.61766	0.00093	
	Data point 74			0.17063 mL			-0.01484		0.00094	
	Data point 75			0.17079 mL			-0.01384		0.00098	
1:15:21.1	Data point 76			0.17096 mL			-0.01653		0.00096	
1:16:07.3				0.17119 mL			-0.01049		0.00077	
	Data point 78			0.17145 mL			-0.00202		0.00062	
	Data point 79			0.17168 mL			-0.01309		0.00095	

1:17:51.5 Assay volumes 1.50000 mL 0.16268 mL 0.17168 mL 0.31999 mL



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

i licitatile. C.\Sirius_13\weii	tap/20100220_exp20	5_10gP_13-2\10E	5-200   1_WIO7_OCIAIIO	_pn-illettic iligii logr
Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	10.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings				
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium				
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.020 mL			
Partition solvent added	Automatic			
After partition addition, stir for	1 seconds			
Sample Sonication				
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	300 seconds			
After sonication stir for	5 seconds			
Sample Dissolution	0 00001100			
Perform a dissolution stage	Yes			
	To start pH			
Stir to dissolve for	120 seconds			
For dissolve for	10%			
The state of the s	10 /0			
Carbonate purge	No			
Perform a carbonate purge	INU			
Temperature Control	Voo			
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			
Stir speed of	50%			
Titration 1				
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	50%			
Titration 2				
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.050 mL			
Additional partition colvent added	Automatic			

Reported at: 3/2/2018 1:05:03 PM

After pH adjust stir for

Stir to allow partitioning for

Stirrer speed for partitioning

Additional partition solvent added Automatic

30 seconds

15 seconds

55%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

# Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3		•	_	•
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.250 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00100 dpH/dt			
Stability timeout after	60 seconds			

# Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.130	2/28/2018 4:26:20 PM	C:\Sirius_T3\HCl18B27.t3r
Four-Plus S	0.9970	2/28/2018 4:26:20 PM	C:\Sirius_T3\HCl18B27.t3r
Four-Plus jH	8.0	2/28/2018 4:26:20 PM	C:\Sirius_T3\HCl18B27.t3r
Four-Plus jOH	-0.4	2/28/2018 4:26:20 PM	C:\Sirius_T3\HCl18B27.t3r
Base concentration factor	1.000	2/28/2018 4:26:20 PM	C:\Sirius_T3\KOH18B27.t3r
Acid concentration factor	0.994	2/28/2018 4:26:20 PM	C:\Sirius_T3\HCl18B27.t3r

# Instrument Settings

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module		T3DM1200361	3/31/2009 5:24:52 AM
Dispenser 0	Water		3/31/2009 5:25:05 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)	00 00 0010	0/07/0040 40 05 50 414
Titrant	Water (0.15 M KCI)	02-06-2018	2/27/2018 10:05:59 AM
Dispenser 2 Syringe volume	Acid 0.5 mL		3/31/2009 5:25:11 AM
Firmware version	1.2.1(r2)		
Titrant	Acid (0.5 M HCI)	02-27-2018	2/27/2018 10:27:22 AM
Dispenser 1	Base	02 27 2010	3/31/2009 5:25:21 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Base (0.5 M KOH)	9/22/2017	2/27/2018 10:21:22 AM
Dispenser 5	Cosolvent		3/31/2009 5:26:24 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		0/04/0000 5:00:40 AM
Distribution valve 5	Distribution Valve		3/31/2009 5:28:19 AM
Firmware version Port A	1.1.3 Methanol (80%, 0.15 M KCl)	09-26-17	2/7/2018 9:42:01 AM
Port B	Cyclohexane	11-01-17	2/27/2018 10:37:57 AM
Dispenser 3	Buffer	11-01-17	8/3/2010 5:05:16 AM
Syringe volume	0.5 mL		0,0,2010 0.00.107
Firmware version	1.2.1(r2)		
Titrant	Dodecane	2018/01/31	2/28/2018 10:18:04 AM
Dispenser 6	Octanol		10/22/2010 10:52:43 AM

Reported at: 3/2/2018 1:05:03 PM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

# Instrument Settings (continued)

Setting Syringe volume	Value 0.5 mL	Batch Id	Install date
Firmware version	1.2.1(r2)		
Titrant	Octanol	01-31-2018	2/27/2018 9:59:35 AM
Titrator		T3TM1200161	3/31/2009 5:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1DI0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1 T3 Electrode	T2E0022	1/22/2010 2:01:00 DM
Electrode E0 calibration	+3.96 mV	T3E0923	1/23/2018 2:01:00 PM 2/28/2018 4:27:04 PM
Filling solution	3M KCI	KCL097	2/27/2018 9:49:43 AM
Liquids		NOLUU1	2/2//2010 3.43.43 AW
Wash 1	50% IPA:50% Water		2/28/2018 10:23:32 AM
Wash 2	0.5% Trition X-100 in H20		2/28/2018 10:23:34 AM
Buffer position 1	pH7 Wash		2/28/2018 10:24:06 AM
Buffer position 2	pH 7		2/28/2018 10:24:08 AM
Storage position			2/28/2018 10:21:14 AM
Wash water	8.8e+003 mL	02-27-2018	2/27/2018 9:54:39 AM
Waste	6.7e+003 mL		11/28/2017 10:36:29 AM
Temperature controller			8/5/2010 6:35:13 AM
Turbidity detector Spectrometer		074811	3/31/2009 5:24:45 AM 11/23/2010 11:22:28 AM
Dip probe		10196	11/23/2010 11.22.28 AW
Wavelength coefficient A0	183.333	10130	
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	112:08:55		11/23/2010 11:22:28 AM
Calibrated on	2/27/2018 10:40:38 AM		
Integration time	40		
Scans averaged	10	T041 400004F	
Autoloader	1 17 AI1DIODOO Ctorner 2	T3AL1200345	11/10/2015 9:34:13 AM
Left-right axis firmware version Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Configuration	,		
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50 20.0 mL		
Flowing wash pump volume Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation			
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s 30%		
E0 calibration preparation stir speed E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		

Reported at: 3/2/2018 1:05:03 PM



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

### Instrument Settings (continued)

Setting	Value	Batch Id	Install date	
Spectrometer calibration stir duration	5 s			
Spectrometer calibration stir speed	30%			
Spectrometer calibration wash pump volume	20.0 mL			
Spectrometer calibration wash stir duration	5 s			
Spectrometer calibration wash stir speed	30%			
Overhead dispense height	10000			

### Refinement Settings

Value	Default value
None	None
500.0 nm	500.0 nm
0.100	0.100
50.00	50.00
	None 500.0 nm 0.100

### Experiment Log

- [1:57] Air gap released for Acid (0.5 M HCI) [1:57] Air gap released for Base (0.5 M KOH)
- [2:33] Air gap created for Water (0.15 M KCI) [2:33] Air gap created for Acid (0.5 M HCI)
- [2:33] Air gap created for Base (0.5 M KOH)
- [2:34] Air gap released for Water (0.15 M KCl)
- [2:38] Titrator arm moved over Titration position
- [2:38] Titration 1 of 3
- [2:38] Adding initial titrants
- [2:38] Automatically add 1.50000 mL of water
- [3:03] Dispensed 1.500000 mL of Water (0.15 M KCI)
- [3:07] Titrator arm moved over Drain
- [8:49] Titrator arm moved to Titration position
- [8:49] Argon flow rate set to 100 [8:49] Stirrer speed set to 10
- [8:54] Automatically add 0.02000 mL of Octanol
- [8:54] Dispensed 0.019991 mL of Octanol
- [8:55] Initial pH = 6.32
- [8:55] Iterative adjust 6.32 -> 2.00 [8:55] pH 6.32 -> 2.00
- [8:57] Air gap released for Acid (0.5 M HCI)
- [8:58] Dispensed 0.048871 mL of Acid (0.5 M HCI)
- [9:03] Holding pH 2.00
- [11:03] Stirrer speed set to 0
- [11:03] Stirrer speed set to 50
- [11:03] Iterative adjust 1.98 -> 2.00
- [11:03] pH 1.98 -> 2.00
- [11:04] Air gap released for Base (0.5 M KOH)
- [11:04] Dispensed 0.002446 mL of Base (0.5 M KOH)
- [11:55] Stirrer speed set to 0
- [12:05] Datapoint id 1 collected
- [12:05] Stirrer speed set to 50
- [12:10] pH 2.02 -> 2.22
- [12:10] Using cautious pH adjust
- [12:10] Dispensed 0.007643 mL of Base (0.5 M KOH)
- [12:15] Stepping pH = 2.11
- [12:16] Dispensed 0.005339 mL of Base (0.5 M KOH)
- [12:21] Stepping pH = 2.19
- [12:21] Dispensed 0.001458 mL of Base (0.5 M KOH)
- [12:26] Stepping pH = 2.22
- [12:41] Stirrer speed set to 0
- [12:52] Datapoint id 2 collected



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [12:52] Charge balance equation is out by 5.5%
- [12:52] Stirrer speed set to 50
- [12:57] pH 2.22 -> 2.42
- [12:57] Using charge balance adjust
- [12:57] Dispensed 0.009525 mL of Base (0.5 M KOH)
- [13:17] Stirrer speed set to 0
- [13:27] Datapoint id 3 collected
- [13:27] Charge balance equation is out by 6.0%
- [13:27] Stirrer speed set to 50
- [13:32] pH 2.44 -> 2.64
- [13:32] Using charge balance adjust
- [13:32] Dispensed 0.005738 mL of Base (0.5 M KOH)
- [13:53] Stirrer speed set to 0
- [14:03] Datapoint id 4 collected
- [14:03] Charge balance equation is out by 5.9%
- [14:03] Stirrer speed set to 50
- [14:08] pH 2.66 -> 2.86
- [14:08] Using charge balance adjust
- [14:08] Dispensed 0.003481 mL of Base (0.5 M KOH)
- [14:28] Stirrer speed set to 0
- [14:38] Datapoint id 5 collected
- [14:38] Charge balance equation is out by 14.5%
- [14:38] Stirrer speed set to 50
- [14:43] pH 2.89 -> 3.09
- [14:43] Using charge balance adjust
- [14:43] Dispensed 0.002070 mL of Base (0.5 M KOH)
- [15:04] Stirrer speed set to 0
- [15:14] Datapoint id 6 collected
- [15:14] Charge balance equation is out by -10.8%
- [15:14] Stirrer speed set to 50
- [15:19] pH 3.08 -> 3.28
- [15:19] Using charge balance adjust
- [15:19] Dispensed 0.001411 mL of Base (0.5 M KOH)
- [15:39] Stirrer speed set to 0
- [15:49] Datapoint id 7 collected
- [15:49] Charge balance equation is out by -9.7%
- [15:49] Stirrer speed set to 50
- [15:54] pH 3.26 -> 3.46
- [15:54] Using charge balance adjust
- [15:54] Dispensed 0.001011 mL of Base (0.5 M KOH)
- [16:14] Stirrer speed set to 0
- [16:24] Datapoint id 8 collected
- [16:24] Charge balance equation is out by -16.8%
- [16:24] Stirrer speed set to 50
- [16:29] pH 3.44 -> 3.64
- [16:29] Using cautious pH adjust
- [16:30] Dispensed 0.000400 mL of Base (0.5 M KOH)
- [16:35] Stepping pH = 3.50
- [16:35] Dispensed 0.000470 mL of Base (0.5 M KOH)
- [16:40] Stepping pH = 3.66
- [16:55] Stirrer speed set to 0
- [17:05] Datapoint id 9 collected
- [17:05] Charge balance equation is out by -8.3%
- [17:05] Stirrer speed set to 50
- [17:10] pH 3.67 -> 3.87
- [17:10] Using charge balance adjust
- [17:10] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [17:30] Stirrer speed set to 0
- [17:41] Datapoint id 10 collected



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [17:41] Charge balance equation is out by 15.6%
- [17:41] Stirrer speed set to 50
- [17:46] pH 3.91 -> 4.11
- [17:46] Using cautious pH adjust
- [17:46] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [17:51] Stepping pH = 4.02
- [17:51] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [17:56] Stepping pH = 4.10
- [17:56] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [18:01] Stepping pH = 4.11
- [18:16] Stirrer speed set to 0
- [18:26] Datapoint id 11 collected
- [18:26] Charge balance equation is out by 11.4%
- [18:26] Stirrer speed set to 50
- [18:31] pH 4.10 -> 4.30
- [18:31] Using charge balance adjust
- [18:32] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [18:52] Stirrer speed set to 0
- [19:02] Datapoint id 12 collected
- [19:02] Charge balance equation is out by 29.2%
- [19:02] Stirrer speed set to 50
- [19:07] pH 4.36 -> 4.56
- [19:07] Using cautious pH adjust
- [19:07] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [19:12] Stepping pH = 4.48
- [19:12] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [19:17] Stepping pH = 4.55
- [19:17] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [19:22] Stepping pH = 4.55
- [19:22] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [19:28] Stepping pH = 4.62
- [19:43] Stirrer speed set to 0
- [19:54] Datapoint id 13 collected
- [19:54] Charge balance equation is out by -12.3%
- [19:54] Stirrer speed set to 50
- [19:59] pH 4.61 -> 4.81
- [19:59] Using charge balance adjust
- [19:59] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [20:20] Stirrer speed set to 0
- [20:35] Datapoint id 14 collected
- [20:35] Charge balance equation is out by -8.1%
- [20:35] Stirrer speed set to 50
- [20:40] pH 4.81 -> 5.01
- [20:40] Using charge balance adjust
- [20:40] Dispensed 0.000847 mL of Base (0.5 M KOH)
- [21:00] Stirrer speed set to 0
- [21:13] Datapoint id 15 collected
- [21:13] Charge balance equation is out by -18.0%
- [21:13] Stirrer speed set to 50
- [21:18] pH 5.00 -> 5.20
- [21:18] Using cautious pH adjust
- [21:18] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [21:23] Stepping pH = 5.07
- [21:23] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [21:28] Stepping pH = 5.17
- [21:28] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [21:33] Stepping pH = 5.14
- [21:33] Dispensed 0.000870 mL of Base (0.5 M KOH)
- [21:39] Stepping pH = 5.45



pH-metric high logP **Pion** Assay name: Analyst: Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [21:54] Stirrer speed set to 0
- [22:09] Datapoint id 16 collected
- [22:09] Charge balance equation is out by -134.1%
- [22:09] Stirrer speed set to 50
- [22:14] pH 5.41 -> 5.61
- [22:14] Using cautious pH adjust
- [22:14] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [22:19] Stepping pH = 5.41
- [22:19] Dispensed 0.001035 mL of Base (0.5 M KOH)
- [22:24] Stepping pH = 6.09
- [22:39] Stirrer speed set to 0
- [23:05] Datapoint id 17 collected
- [23:05] Charge balance equation is out by -201.4%
- [23:05] Stirrer speed set to 50
- [23:10] pH 5.99 -> 6.19
- [23:10] Using cautious pH adjust [23:10] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [23:15] Stepping pH = 5.99
- [23:15] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [23:20] Stepping pH = 6.23
- [23:35] Stirrer speed set to 0
- [24:01] Datapoint id 18 collected
- [24:01] Charge balance equation is out by -98.3%
- [24:01] Stirrer speed set to 50
- [24:06] pH 6.26 -> 6.46
- [24:06] Using cautious pH adjust
- [24:06] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [24:11] Stepping pH = 6.26
- [24:11] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [24:16] Stepping pH = 6.37
- [24:16] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [24:21] Stepping pH = 6.59
- [24:37] Stirrer speed set to 0
- [25:17] Datapoint id 19 collected
- [25:17] Charge balance equation is out by -178.7%
- [25:17] Stirrer speed set to 50
- [25:22] pH 6.71 -> 6.91
- [25:22] Using cautious pH adjust
- [25:22] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [25:27] Stepping pH = 6.72
- [25:27] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [25:32] Stepping pH = 6.75
- [25:32] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [25:37] Stepping pH = 8.75
- [25:52] Stirrer speed set to 0
- [26:52] Datapoint id 20 collected
- [26:52] Charge balance equation is out by -408.8%
- [26:52] Stirrer speed set to 50
- [26:57] pH 8.97 -> 9.17
- [26:57] Using cautious pH adjust
- [26:57] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [27:03] Stepping pH = 8.97
- [27:03] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [27:08] Stepping pH = 8.98
- [27:08] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [27:13] Stepping pH = 9.83
- [27:28] Stirrer speed set to 0
- [27:45] Datapoint id 21 collected
- [27:45] Charge balance equation is out by -942.4%



**Pion** Assay name: pH-metric high logP Analyst: Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [27:45] Stirrer speed set to 50
- [27:50] pH 9.85 -> 10.05
- [27:50] Using cautious pH adjust
- [27:50] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [27:55] Stepping pH = 9.87
- [27:55] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [28:00] Stepping pH = 10.13
- [28:15] Stirrer speed set to 0
- [28:34] Datapoint id 22 collected
- [28:34] Charge balance equation is out by -73.7%
- [28:34] Titration 2 of 3
- [28:34] Adding initial titrants
- [28:34] Automatically add 0.05000 mL of Octanol
- [28:35] Dispensed 0.050000 mL of Octanol
- [28:35] Stirrer speed set to 10
- [28:36] Stirrer speed set to 55
- [28:36] Iterative adjust 10.14 -> 2.00
- [28:36] pH 10.14 -> 2.00
- [28:37] Dispensed 0.052658 mL of Acid (0.5 M HCI)
- [28:42] pH 2.02 -> 2.00
- [28:43] Dispensed 0.001976 mL of Acid (0.5 M HCI)
- [29:33] Stirrer speed set to 0
- [29:43] Datapoint id 23 collected
- [29:43] Stirrer speed set to 55
- [29:48] pH 1.97 -> 2.17
- [29:48] Using cautious pH adjust
- [29:49] Dispensed 0.009196 mL of Base (0.5 M KOH)
- [29:54] Stepping pH = 2.06
- [29:54] Dispensed 0.006515 mL of Base (0.5 M KOH)
- [29:59] Stepping pH = 2.14
- [29:59] Dispensed 0.001505 mL of Base (0.5 M KOH)
- [30:04] Stepping pH = 2.16
- [30:20] Stirrer speed set to 0
- [30:30] Datapoint id 24 collected
- [30:30] Charge balance equation is out by 6.5%
- [30:30] Stirrer speed set to 55
- [30:35] pH 2.16 -> 2.36
- [30:35] Using charge balance adjust
- [30:36] Dispensed 0.011618 mL of Base (0.5 M KOH)
- [30:56] Stirrer speed set to 0
- [31:06] Datapoint id 25 collected
- [31:06] Charge balance equation is out by 20.5%
- [31:06] Stirrer speed set to 55
- [31:11] pH 2.41 -> 2.61
- [31:11] Using cautious pH adjust
- [31:11] Dispensed 0.003293 mL of Base (0.5 M KOH)
- [31:16] Stepping pH = 2.51
- [31:16] Dispensed 0.002023 mL of Base (0.5 M KOH)
- [31:21] Stepping pH = 2.58
- [31:21] Dispensed 0.000823 mL of Base (0.5 M KOH)
- [31:27] Stepping pH = 2.60
- [31:42] Stirrer speed set to 0
- [31:52] Datapoint id 26 collected
- [31:52] Charge balance equation is out by 6.8%
- [31:52] Stirrer speed set to 55
- [31:57] pH 2.61 -> 2.81
- [31:57] Using charge balance adjust
- [31:57] Dispensed 0.004257 mL of Base (0.5 M KOH)
- [32:17] Stirrer speed set to 0



pH-metric high logP **Pion** Assay name: Analyst: Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [32:27] Datapoint id 27 collected
- [32:27] Charge balance equation is out by 11.6%
- [32:27] Stirrer speed set to 55
- [32:32] pH 2.83 -> 3.03
- [32:32] Using charge balance adjust
- [32:32] Dispensed 0.002634 mL of Base (0.5 M KOH)
- [32:53] Stirrer speed set to 0
- [33:03] Datapoint id 28 collected
- [33:03] Charge balance equation is out by 26.9%
- [33:03] Stirrer speed set to 55
- [33:08] pH 3.09 -> 3.29
- [33:08] Using cautious pH adjust
- [33:08] Dispensed 0.000823 mL of Base (0.5 M KOH)
- [33:13] Stepping pH = 3.19
- [33:13] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [33:18] Stepping pH = 3.27
- [33:18] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [33:23] Stepping pH = 3.28
- [33:23] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [33:28] Stepping pH = 3.30
- [33:44] Stirrer speed set to 0
- [33:54] Datapoint id 29 collected
- [33:54] Charge balance equation is out by 1.7%
- [33:54] Stirrer speed set to 55
- [33:59] pH 3.30 -> 3.50
- [33:59] Using charge balance adjust
- [33:59] Dispensed 0.001246 mL of Base (0.5 M KOH)
- [34:19] Stirrer speed set to 0
- [34:29] Datapoint id 30 collected
- [34:29] Charge balance equation is out by -0.4%
- [34:29] Stirrer speed set to 55
- [34:34] pH 3.51 -> 3.71
- [34:34] Using charge balance adjust
- [34:34] Dispensed 0.001105 mL of Base (0.5 M KOH)
- [34:54] Stirrer speed set to 0
- [35:05] Datapoint id 31 collected
- [35:05] Charge balance equation is out by 1.1%
- [35:05] Stirrer speed set to 55
- [35:10] pH 3.71 -> 3.91
- [35:10] Using charge balance adjust
- [35:10] Dispensed 0.001058 mL of Base (0.5 M KOH)
- [35:30] Stirrer speed set to 0
- [35:40] Datapoint id 32 collected
- [35:40] Charge balance equation is out by 2.5%
- [35:40] Stirrer speed set to 55
- [35:45] pH 3.92 -> 4.12
- [35:45] Using charge balance adjust
- [35:46] Dispensed 0.001058 mL of Base (0.5 M KOH)
- [36:06] Stirrer speed set to 0
- [36:16] Datapoint id 33 collected
- [36:16] Charge balance equation is out by 3.6%
- [36:16] Stirrer speed set to 55
- [36:21] pH 4.13 -> 4.33
- [36:21] Using charge balance adjust
- [36:21] Dispensed 0.001011 mL of Base (0.5 M KOH)
- [36:41] Stirrer speed set to 0
- [36:51] Datapoint id 34 collected
- [36:51] Charge balance equation is out by -3.3%
- [36:51] Stirrer speed set to 55



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [36:56] pH 4.34 -> 4.54
- [36:56] Using charge balance adjust
- [36:56] Dispensed 0.000894 mL of Base (0.5 M KOH)
- [37:17] Stirrer speed set to 0
- [37:27] Datapoint id 35 collected
- [37:27] Charge balance equation is out by -9.8%
- [37:27] Stirrer speed set to 55
- [37:32] pH 4.53 -> 4.73
- [37:32] Using charge balance adjust
- [37:32] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [37:53] Stirrer speed set to 0
- [38:03] Datapoint id 36 collected
- [38:03] Charge balance equation is out by -11.7%
- [38:03] Stirrer speed set to 55
- [38:08] pH 4.72 -> 4.92
- [38:08] Using charge balance adjust
- [38:08] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [38:29] Stirrer speed set to 0
- [38:39] Datapoint id 37 collected
- [38:39] Charge balance equation is out by -30.8%
- [38:39] Stirrer speed set to 55
- [38:44] pH 4.87 -> 5.07
- [38:44] Using cautious pH adjust
- [38:44] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [38:49] Stepping pH = 4.93
- [38:49] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [38:54] Stepping pH = 5.01
- [38:55] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [39:00] Stepping pH = 5.05
- [39:00] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [39:05] Stepping pH = 5.07
- [39:20] Stirrer speed set to 0
- [39:30] Datapoint id 38 collected
- [39:30] Charge balance equation is out by -55.1%
- [39:30] Stirrer speed set to 55
- [39:36] pH 5.09 -> 5.29
- [39:36] Using cautious pH adjust
- [39:36] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [39:41] Stepping pH = 5.13
- [39:41] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [39:46] Stepping pH = 5.28
- [39:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [39:51] Stepping pH = 5.29
- [40:06] Stirrer speed set to 0
- [40:18] Datapoint id 39 collected
- [40:18] Charge balance equation is out by -46.4%
- [40:18] Stirrer speed set to 55
- [40:23] pH 5.30 -> 5.50
- [40:23] Using cautious pH adjust
- [40:23] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [40:28] Stepping pH = 5.32
- [40:28] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [40:33] Stepping pH = 5.58
- [40:48] Stirrer speed set to 0
- [41:00] Datapoint id 40 collected
- [41:00] Charge balance equation is out by -83.4%
- [41:00] Stirrer speed set to 55
- [41:05] pH 5.62 -> 5.82
- [41:05] Using cautious pH adjust



pH-metric high logP **Pion** Assay name: Analyst: Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [41:05] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [41:11] Stepping pH = 5.63
- [41:11] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [41:16] Stepping pH = 5.66
- [41:16] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [41:21] Stepping pH = 6.69
- [41:36] Stirrer speed set to 0
- [42:36] Datapoint id 41 collected
- [42:36] Charge balance equation is out by -392.1%
- [42:36] Stirrer speed set to 55
- [42:41] pH 6.97 -> 7.17
- [42:41] Using cautious pH adjust
- [42:41] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [42:46] Stepping pH = 7.02
- [42:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [42:51] Stepping pH = 7.05
- [42:51] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [42:56] Stepping pH = 7.08
- [42:56] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [43:02] Stepping pH = 7.10
- [43:02] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [43:07] Stepping pH = 7.10
- [43:07] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [43:12] Stepping pH = 9.27
- [43:27] Stirrer speed set to 0
- [43:45] Datapoint id 42 collected
- [43:45] Charge balance equation is out by -1,782.6%
- [43:45] Stirrer speed set to 55
- [43:50] pH 9.30 -> 9.50
- [43:50] Using cautious pH adjust
- [43:50] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [43:56] Stepping pH = 9.30
- [43:56] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [44:01] Stepping pH = 9.60
- [44:16] Stirrer speed set to 0
- [44:37] Datapoint id 43 collected
- [44:37] Charge balance equation is out by -200.1%
- [44:37] Stirrer speed set to 55
- [44:42] pH 9.58 -> 9.78
- [44:42] Using cautious pH adjust
- [44:42] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [44:47] Stepping pH = 9.58
- [44:47] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [44:52] Stepping pH = 9.97
- [45:07] Stirrer speed set to 0
- [45:21] Datapoint id 44 collected
- [45:21] Charge balance equation is out by -201.1%
- [45:21] Stirrer speed set to 55
- [45:26] pH 9.95 -> 10.05
- [45:26] Using cautious pH adjust
- [45:26] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [45:31] Stepping pH = 9.94
- [45:31] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [45:36] Stepping pH = 10.17
- [45:51] Stirrer speed set to 0
- [46:01] Datapoint id 45 collected
- [46:01] Charge balance equation is out by -210.8%
- [46:01] Titration 3 of 3
- [46:01] Adding initial titrants



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [46:01] Automatically add 0.25000 mL of Octanol
- [46:07] Dispensed 0.250000 mL of Octanol
- [46:07] Stirrer speed set to 10
- [46:08] Stirrer speed set to 60
- [46:08] Iterative adjust 10.17 -> 2.00
- [46:08] pH 10.17 -> 2.00
- [46:10] Dispensed 0.055715 mL of Acid (0.5 M HCI)
- [46:15] pH 2.03 -> 2.00
- [46:15] Dispensed 0.003457 mL of Acid (0.5 M HCI)
- [47:05] Stirrer speed set to 0
- [47:15] Datapoint id 46 collected
- [47:15] Stirrer speed set to 60
- [47:20] pH 1.96 -> 2.16
- [47:20] Using cautious pH adjust
- [47:21] Dispensed 0.010113 mL of Base (0.5 M KOH)
- [47:26] Stepping pH = 2.02
- [47:26] Dispensed 0.010865 mL of Base (0.5 M KOH)
- [47:31] Stepping pH = 2.09
- [47:31] Dispensed 0.006515 mL of Base (0.5 M KOH)
- [47:36] Stepping pH = 2.18
- [47:51] Stirrer speed set to 0
- [48:01] Datapoint id 47 collected
- [48:01] Charge balance equation is out by -35.9%
- [48:01] Stirrer speed set to 60
- [48:06] pH 2.19 -> 2.39
- [48:06] Using cautious pH adjust
- [48:07] Dispensed 0.005903 mL of Base (0.5 M KOH)
- [48:12] Stepping pH = 2.29
- [48:12] Dispensed 0.003481 mL of Base (0.5 M KOH)
- [48:17] Stepping pH = 2.36
- [48:17] Dispensed 0.001435 mL of Base (0.5 M KOH)
- [48:22] Stepping pH = 2.39
- [48:37] Stirrer speed set to 0
- [48:47] Datapoint id 48 collected
- [48:47] Charge balance equation is out by 8.3%
- [48:47] Stirrer speed set to 60
- [48:52] pH 2.39 -> 2.59
- [48:52] Using charge balance adjust
- [48:53] Dispensed 0.007502 mL of Base (0.5 M KOH)
- [49:13] Stirrer speed set to 0
- [49:23] Datapoint id 49 collected
- [49:23] Charge balance equation is out by 9.3%
- [49:23] Stirrer speed set to 60
- [49:28] pH 2.62 -> 2.82
- [49:28] Using charge balance adjust
- [49:28] Dispensed 0.004657 mL of Base (0.5 M KOH)
- [49:48] Stirrer speed set to 0
- [49:59] Datapoint id 50 collected
- [49:59] Charge balance equation is out by 13.0%
- [49:59] Stirrer speed set to 60
- [50:04] pH 2.85 -> 3.05
- [50:04] Using charge balance adjust
- [50:04] Dispensed 0.002987 mL of Base (0.5 M KOH)
- [50:24] Stirrer speed set to 0
- [50:34] Datapoint id 51 collected
- [50:34] Charge balance equation is out by 15.0%
- [50:34] Stirrer speed set to 60
- [50:39] pH 3.08 -> 3.28
- [50:39] Using cautious pH adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [50:39] Dispensed 0.001058 mL of Base (0.5 M KOH)
- [50:44] Stepping pH = 3.17
- [50:44] Dispensed 0.000823 mL of Base (0.5 M KOH)
- [50:50] Stepping pH = 3.24
- [50:50] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [50:55] Stepping pH = 3.27
- [51:10] Stirrer speed set to 0
- [51:20] Datapoint id 52 collected
- [51:20] Charge balance equation is out by -4.5%
- [51:20] Stirrer speed set to 60
- [51:25] pH 3.28 -> 3.48
- [51:25] Using charge balance adjust
- [51:25] Dispensed 0.001717 mL of Base (0.5 M KOH)
- [51:45] Stirrer speed set to 0
- [51:55] Datapoint id 53 collected
- [51:55] Charge balance equation is out by -5.2%
- [51:55] Stirrer speed set to 60
- [52:00] pH 3.47 -> 3.67
- [52:00] Using charge balance adjust
- [52:00] Dispensed 0.001458 mL of Base (0.5 M KOH)
- [52:21] Stirrer speed set to 0
- [52:31] Datapoint id 54 collected
- [52:31] Charge balance equation is out by -0.7%
- [52:31] Stirrer speed set to 60
- [52:36] pH 3.68 -> 3.88
- [52:36] Using charge balance adjust
- [52:36] Dispensed 0.001246 mL of Base (0.5 M KOH)
- [52:56] Stirrer speed set to 0
- [53:06] Datapoint id 55 collected
- [53:06] Charge balance equation is out by 3.3%
- [53:06] Stirrer speed set to 60
- [53:11] pH 3.89 -> 4.09
- [53:11] Using charge balance adjust
- [53:11] Dispensed 0.001011 mL of Base (0.5 M KOH)
- [53:31] Stirrer speed set to 0
- [53:41] Datapoint id 56 collected
- [53:41] Charge balance equation is out by -1.8%
- [53:41] Stirrer speed set to 60
- [53:47] pH 4.10 -> 4.30
- [53:47] Using charge balance adjust
- [53:47] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [54:07] Stirrer speed set to 0
- [54:17] Datapoint id 57 collected
- [54:17] Charge balance equation is out by -9.4%
- [54:17] Stirrer speed set to 60
- [54:22] pH 4.29 -> 4.49
- [54:22] Using charge balance adjust
- [54:22] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [54:42] Stirrer speed set to 0
- [54:52] Datapoint id 58 collected
- [54:52] Charge balance equation is out by -20.0%
- [54:52] Stirrer speed set to 60
- [54:57] pH 4.45 -> 4.65
- [54:57] Using cautious pH adjust
- [54:57] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [55:02] Stepping pH = 4.52
- [55:02] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [55:08] Stepping pH = 4.61
- [55:08] Dispensed 0.000094 mL of Base (0.5 M KOH)



pH-metric high logP **Pion** Assay name: Analyst: Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [55:13] Stepping pH = 4.62
- [55:13] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [55:18] Stepping pH = 4.67
- [55:33] Stirrer speed set to 0
- [55:43] Datapoint id 59 collected
- [55:43] Charge balance equation is out by -49.9%
- [55:43] Stirrer speed set to 60
- [55:49] pH 4.68 -> 4.88
- [55:49] Using cautious pH adjust
- [55:49] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [55:54] Stepping pH = 4.75
- [55:54] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [55:59] Stepping pH = 4.86
- [55:59] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [56:04] Stepping pH = 4.86
- [56:04] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [56:09] Stepping pH = 4.92
- [56:24] Stirrer speed set to 0
- [56:35] Datapoint id 60 collected
- [56:35] Charge balance equation is out by -46.6%
- [56:35] Stirrer speed set to 60
- [56:40] pH 4.93 -> 5.13
- [56:40] Using cautious pH adjust
- [56:40] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [56:45] Stepping pH = 4.97
- [56:45] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [56:50] Stepping pH = 5.15
- [57:05] Stirrer speed set to 0
- [57:16] Datapoint id 61 collected
- [57:16] Charge balance equation is out by -38.7%
- [57:16] Stirrer speed set to 60
- [57:21] pH 5.16 -> 5.36
- [57:21] Using cautious pH adjust
- [57:21] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [57:27] Stepping pH = 5.17
- [57:27] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [57:32] Stepping pH = 5.35
- [57:47] Stirrer speed set to 0
- [57:58] Datapoint id 62 collected [57:58] Charge balance equation is out by -87.8%
- [57:58] Stirrer speed set to 60
- [58:03] pH 5.38 -> 5.58
- [58:03] Using cautious pH adjust
- [58:03] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [58:09] Stepping pH = 5.41
- [58:09] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [58:14] Stepping pH = 5.54
- [58:14] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [58:19] Stepping pH = 5.57
- [58:34] Stirrer speed set to 0
- [58:47] Datapoint id 63 collected
- [58:47] Charge balance equation is out by -84.7%
- [58:47] Stirrer speed set to 60
- [58:52] pH 5.63 -> 5.83
- [58:52] Using cautious pH adjust
- [58:52] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [58:57] Stepping pH = 5.64
- [58:57] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [59:02] Stepping pH = 5.74



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

```
[59:02] Dispensed 0.000047 mL of Base (0.5 M KOH)
```

- [59:07] Stepping pH = 5.86
- [59:22] Stirrer speed set to 0
- [59:58] Datapoint id 64 collected
- [59:58] Charge balance equation is out by -167.1%
- [59:58] Stirrer speed set to 60
- [1:00:03] pH 5.99 -> 6.19
- [1:00:03] Using cautious pH adjust
- [1:00:03] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:00:08] Stepping pH = 6.02
- [1:00:08] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:00:13] Stepping pH = 6.15
- [1:00:13] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:00:18] Stepping pH = 6.30
- [1:00:33] Stirrer speed set to 0
- [1:01:33] Datapoint id 65 collected
- [1:01:33] Charge balance equation is out by -96.8%
- [1:01:33] Stirrer speed set to 60
- [1:01:38] pH 6.35 -> 6.55
- [1:01:38] Using cautious pH adjust
- [1:01:38] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:01:43] Stepping pH = 6.40
- [1:01:43] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:01:48] Stepping pH = 6.63
- [1:02:04] Stirrer speed set to 0
- [1:03:04] Datapoint id 66 collected
- [1:03:04] Charge balance equation is out by -27.0%
- [1:03:04] Stirrer speed set to 60
- [1:03:09] pH 6.73 -> 6.93
- [1:03:09] Using cautious pH adjust
- [1:03:09] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:03:14] Stepping pH = 6.79
- [1:03:14] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:03:19] Stepping pH = 7.01
- [1:03:34] Stirrer speed set to 0
- [1:04:34] Datapoint id 67 collected
- [1:04:34] Charge balance equation is out by -13.9%
- [1:04:34] Stirrer speed set to 60
- [1:04:39] pH 7.04 -> 7.24
- [1:04:39] Using charge balance adjust
- [1:04:39] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:04:59] Stirrer speed set to 0
- [1:05:59] Datapoint id 68 collected
- [1:05:59] Charge balance equation is out by -70.5%
- [1:05:59] Stirrer speed set to 60
- [1:06:05] pH 7.21 -> 7.41
- [1:06:05] Using cautious pH adjust
- [1:06:05] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:06:10] Stepping pH = 7.21
- [1:06:10] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:06:15] Stepping pH = 7.44
- [1:06:30] Stirrer speed set to 0
- [1:07:30] Datapoint id 69 collected
- [1:07:30] Charge balance equation is out by -155.1%
- [1:07:30] Stirrer speed set to 60
- [1:07:35] pH 7.73 -> 7.93
- [1:07:35] Using cautious pH adjust
- [1:07:35] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:07:40] Stepping pH = 7.83



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [1:07:40] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:07:45] Stepping pH = 8.04
- [1:08:00] Stirrer speed set to 0
- [1:09:00] Datapoint id 70 collected
- [1:09:00] Charge balance equation is out by -346.1%
- [1:09:00] Stirrer speed set to 60
- [1:09:06] pH 7.76 -> 7.96
- [1:09:06] Using cautious pH adjust
- [1:09:06] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:09:11] Stepping pH = 7.79
- [1:09:11] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:09:16] Stepping pH = 8.01
- [1:09:31] Stirrer speed set to 0
- [1:10:31] Datapoint id 71 collected
- [1:10:31] Charge balance equation is out by -354.4%
- [1:10:31] Stirrer speed set to 60
- [1:10:36] pH 8.11 -> 8.31
- [1:10:36] Using cautious pH adjust
- [1:10:36] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:10:41] Stepping pH = 8.17
- [1:10:41] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:10:46] Stepping pH = 8.19
- [1:10:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:10:51] Stepping pH = 8.17
- [1:10:51] Dispensed 0.000094 mL of Base (0.5 M KOH)
- 1.10.01] Disperioca 0.0000
- [1:10:57] Stepping pH = 8.47
- [1:11:12] Stirrer speed set to 0
- [1:12:12] Datapoint id 72 collected
- [1:12:12] Charge balance equation is out by -1,435.6%
- [1:12:12] Stirrer speed set to 60
- [1:12:17] pH 8.43 -> 8.63
- [1:12:17] Using cautious pH adjust
- [1:12:17] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:12:22] Stepping pH = 8.42
- [1:12:22] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:12:27] Stepping pH = 8.45
- [1:12:27] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:12:32] Stepping pH = 8.81
- [1:12:47] Stirrer speed set to 0
- [1:13:08] Datapoint id 73 collected
- [1:13:08] Charge balance equation is out by -948.3%
- [1:13:08] Stirrer speed set to 60
- [1:13:13] pH 8.83 -> 9.03
- [1:13:13] Using cautious pH adjust
- [1:13:13] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:13:18] Stepping pH = 8.84
- [1:13:18] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:13:23] Stepping pH = 8.94
- [1:13:23] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:13:28] Stepping pH = 9.00
- [1:13:28] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:13:33] Stepping pH = 9.04
- [1:13:49] Stirrer speed set to 0
- [1:14:02] Datapoint id 74 collected
- [1:14:02] Charge balance equation is out by -262.5%
- [1:14:02] Stirrer speed set to 60
- [1:14:07] pH 9.05 -> 9.25
- [1:14:07] Using cautious pH adjust
- [1:14:07] Dispensed 0.000024 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18B-28011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18B-28011\_M07\_octanol\_pH-metric high logP.t3r

- [1:14:13] Stepping pH = 9.04
- [1:14:13] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:14:18] Stepping pH = 9.33
- [1:14:33] Stirrer speed set to 0
- [1:14:50] Datapoint id 75 collected
- [1:14:50] Charge balance equation is out by -206.3%
- [1:14:50] Stirrer speed set to 60
- [1:14:55] pH 9.34 -> 9.54
- [1:14:55] Using cautious pH adjust
- [1:14:56] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:15:01] Stepping pH = 9.37
- [1:15:01] Dispensed 0.000118 mL of Base (0.5 M KOH)
- 1:15:06 Stepping pH = 9.53
- [1:15:21] Stirrer speed set to 0
- [1:15:31] Datapoint id 76 collected
- [1:15:31] Charge balance equation is out by -78.6%
- [1:15:31] Stirrer speed set to 60
- [1:15:37] pH 9.54 -> 9.74
- [1:15:37] Using cautious pH adjust
- [1:15:37] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:15:42] Stepping pH = 9.57
- [1:15:42] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:15:47] Stepping pH = 9.72
- [1:15:47] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:15:52] Stepping pH = 9.73
- [1:16:07] Stirrer speed set to 0
- [1:16:18] Datapoint id 77 collected
- [1:16:18] Charge balance equation is out by -62.5%
- [1:16:18] Stirrer speed set to 60
- [1:16:23] pH 9.73 -> 9.93
- [1:16:23] Using cautious pH adjust
- [1:16:23] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:16:28] Stepping pH = 9.78
- [1:16:28] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [1:16:33] Stepping pH = 9.93
- [1:16:48] Stirrer speed set to 0
- [1:16:59] Datapoint id 78 collected
- [1:16:59] Charge balance equation is out by -29.9%
- [1:16:59] Stirrer speed set to 60
- [1:17:04] pH 9.92 -> 10.05
- [1:17:04] Using cautious pH adjust
- [1:17:04] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:17:09] Stepping pH = 9.95
- [1:17:09] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:17:14] Stepping pH = 10.04
- [1:17:29] Stirrer speed set to 0
- [1:17:42] Datapoint id 79 collected
- [1:17:42] Charge balance equation is out by -36.7%
- [1:17:42] Argon flow rate set to 0
- [1:17:46] Titrator arm moved over Titration position