

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

## pH-metric Result

logP (XH +) -7.57 ±1.49 (n=50) logP (neutral X) 2.89 ±0.02 (n=50)

RMSD 3.854

#### 18C-16016 Points 2 to 15

M13\_octanol concentration factor 0.878
Carbonate 0.0003 mM
Acidity error 2.29946 mM

#### 18C-16016 Points 16 to 39

M13\_octanol concentration factor 0.780
Carbonate 0.0000 mM
Acidity error 2.26129 mM

#### 18C-16016 Points 40 to 63

M13\_octanol concentration factor 0.706
Carbonate 0.1381 mM
Acidity error 2.26901 mM

#### Warnings and errors

Errors None

Warnings One or more logP values out of range

#### Sample logD and percent species

рН	M13_octanol logD	M13_octanol M13_octanolH	M13_octanol M13_octanol	M13_octanol M13_octanolH*	M13_octanol M13_octanol*	
1.000	-1.88	9 <del>8</del> .68 %	0.00 %	0.00 %	1.31 %	
1.200	-1.68	97.93 %	0.00 %	0.00 %	2.07 %	Stomach pH
2.000	-0.88	88.23 %	0.01 %	0.00 %	11.76 %	
3.000	0.12	42.84 %	0.07 %	0.00 %	57.08 %	
4.000	1.12	6.97 %	0.12 %	0.00 %	92.91 %	
5.000	2.06	0.74 %	0.13 %	0.00 %	99.13 %	
6.000	2.69	0.07 %	0.13 %	0.00 %	99.80 %	
6.500	2.82	0.02 %	0.13 %	0.00 %	99.85 %	
7.000	2.87	0.01 %	0.13 %	0.00 %	99.87 %	
7.400	2.88	0.00 %	0.13 %	0.00 %	99.87 %	Blood pH
8.000	2.89	0.00 %	0.13 %	0.00 %	99.87 %	
9.000	2.89	0.00 %	0.13 %	0.00 %	99.87 %	
10.000	2.89	0.00 %	0.13 %	0.00 %	99.87 %	
11.000	2.89	0.00 %	0.13 %	0.00 %	99.87 %	
12.000	2.89	0.00 %	0.13 %	0.00 %	99.87 %	



Sample name: M13\_octanol Assay name:


pH-metric high logP

Assay ID: 18C-16016 Filename:

Experiment start time: 3/16/2018 8:34:17 PM

Pion Analyst: Instrument ID: T311053

C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r






Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# **Graphs** (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 1 of 3 18C-16016 Points 2 to 15

#### Overall results

RMSD 0.073
Average ionic strength 0.158 M
Average temperature 24.8°C
Partition ratio 0.0624 : 1

Analyte concentration range 2710.7 µM to 2790.7 µM

Total points considered 8 of 14

#### Warnings and errors

Errors None

Warnings One or more logP values out of range

Excessive acidity error present

## Four-Plus parameters


Alpha 0.167 3/16/2018 8:34:16 PM C:\Sirius\_T3\HCl18C16.t3r S 0.9932 3/16/2018 8:34:16 PM C:\Sirius\_T3\HCl18C16.t3r jH 0.7 3/16/2018 8:34:16 PM C:\Sirius\_T3\HCl18C16.t3r jOH -0.9 3/16/2018 8:34:16 PM C:\Sirius\_T3\HCl18C16.t3r


#### Titrants

#### Sample

M13\_octanol concentration factor 0.878
Base pKa 1 5.77
logP (XH +) -4.56
logP (neutral X) 2.90

#### Sample graphs








Assay ID:

Sample name: M13\_octanol Experiment start time: 3/16/2018 8:34:17 PM

Assay name: pH-metric high logP Analyst: Pion Instrument ID: 18C-16016 T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# Sample graphs (continued)




#### Sample logD and percent species


рН	M13_octanol	M13_octanol	M13_octanol	M13_octanol	M13_octanol	Comment
	logD	M13_octanolH	M13_octanol	M13_octanolH*	M13_octanol*	
1.000	-1.87	99.91 %	0.00 %	0.00 %	0.08 %	
1.200	-1.67	99.86 %	0.00 %	0.00 %	0.13 %	Stomach pH
2.000	-0.87	99.15 %	0.02 %	0.00 %	0.83 %	
3.000	0.13	92.10 %	0.16 %	0.00 %	7.75 %	
4.000	1.12	53.81 %	0.91 %	0.00 %	45.27 %	
5.000	2.06	10.43 %	1.77 %	0.00 %	87.79 %	
6.000	2.70	1.15 %	1.96 %	0.00 %	96.89 %	
6.500	2.83	0.37 %	1.97 %	0.00 %	97.66 %	
7.000	2.87	0.12 %	1.98 %	0.00 %	97.91 %	
7.400	2.89	0.05 %	1.98 %	0.00 %	97.98 %	Blood pH
8.000	2.90	0.01 %	1.98 %	0.00 %	98.01 %	
9.000	2.90	0.00 %	1.98 %	0.00 %	98.02 %	
10.000	2.90	0.00 %	1.98 %	0.00 %	98.02 %	
11.000	2.90	0.00 %	1.98 %	0.00 %	98.02 %	
12.000	2.90	0.00 %	1.98 %	0.00 %	98.02 %	

# Carbonate and acidity



# Other graphs








Assay name: pH-metric high logP Analyst: **Pion** Instrument ID: T311053 Assay ID: 18C-16016 Filename:

C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 2 of 3 18C-16016 Points 16 to 39

#### Overall results

RMSD 6.380
Average ionic strength 0.164 M
Average temperature 24.9°C
Partition ratio 0.1752 : 1

Analyte concentration range 2287.9 µM to 2354.2 µM

Total points considered 13 of 24

#### Warnings and errors

Errors None

Warnings One or more logP values out of range

Excessive acidity error present

## Four-Plus parameters


Alpha 0.167 3/16/2018 8:34:16 PM C:\Sirius\_T3\HCl18C16.t3r S 0.9932 3/16/2018 8:34:16 PM C:\Sirius\_T3\HCl18C16.t3r jH 0.7 3/16/2018 8:34:16 PM C:\Sirius\_T3\HCl18C16.t3r jOH -0.9 3/16/2018 8:34:16 PM C:\Sirius\_T3\HCl18C16.t3r


#### Titrants

#### Sample

M13\_octanol concentration factor 0.780
Base pKa 1 5.77
logP (XH +) -4.56
logP (neutral X) 2.88

#### Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r


# Sample graphs (continued)




#### Sample logD and percent species


рН	M13_octanol logD	M13_octanolH	M13_octanol	M13_octanol M13_octanolH*	M13_octanol M13_octanol*	
1.000	-1.89	99.77 %	0.00 %	0.00 %	0.23 %	
1.200	-1.69	99.64 %	0.00 %	0.00 %	0.36 %	Stomach pH
2.000	-0.89	97.77 %	0.02 %	0.00 %	2.22 %	•
3.000	0.11	81.41 %	0.14 %	0.00 %	18.45 %	
4.000	1.10	30.46 %	0.52 %	0.00 %	69.02 %	
5.000	2.04	4.20 %	0.71 %	0.00 %	95.09 %	
6.000	2.68	0.44 %	0.74 %	0.00 %	98.82 %	
6.500	2.81	0.14 %	0.74 %	0.00 %	99.12 %	
7.000	2.86	0.04 %	0.74 %	0.00 %	99.21 %	
7.400	2.87	0.02 %	0.74 %	0.00 %	99.24 %	Blood pH
8.000	2.88	0.00 %	0.74 %	0.00 %	99.25 %	·
9.000	2.88	0.00 %	0.74 %	0.00 %	99.26 %	
10.000	2.88	0.00 %	0.74 %	0.00 %	99.26 %	
11.000	2.88	0.00 %	0.74 %	0.00 %	99.26 %	
12.000	2.88	0.00 %	0.74 %	0.00 %	99.26 %	

# **Carbonate and acidity**



# Other graphs








Assay name: pH-metric high logP Analyst: **Pion** Instrument ID: T311053 Assay ID: 18C-16016 Filename:

C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 3 of 3 18C-16016 Points 40 to 63

#### Overall results

RMSD 0.265
Average ionic strength 0.169 M
Average temperature 25.0°C
Partition ratio 0.4371 : 1

Analyte concentration range 1749.3 µM to 1790.0 µM

Total points considered 16 of 24

#### Warnings and errors

Errors None

Warnings One or more logP values out of range

Excessive acidity error present

## Four-Plus parameters


Alpha	a 0.167	3/16/2018 8:34:16 PM	C:\Sirius_T3\HCl18C16.t3r
S	0.9932	3/16/2018 8:34:16 PM	C:\Sirius_T3\HCl18C16.t3r
jΗ	0.7	3/16/2018 8:34:16 PM	C:\Sirius_T3\HCl18C16.t3r
jOH	-0.9	3/16/2018 8:34:16 PM	C:\Sirius_T3\HCl18C16.t3r


#### **Titrants**

#### Sample

M13\_octanol concentration factor 0.706
Base pKa 1 5.77
logP (XH +) -4.56
logP (neutral X) 2.92

#### Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# Sample graphs (continued)




#### Sample logD and percent species


рН	M13_octanol	M13_octanol	M13_octanol	M13_octanol	M13_octanol	Comment
	logD	M13_octanolH	M13_octanol	M13_octanolH*	M13_octanol*	
1.000	-1.85	99.39 %	0.00 %	0.00 %	0.61 %	
1.200	-1.65	99.03 %	0.00 %	0.00 %	0.96 %	Stomach pH
2.000	-0.85	94.19 %	0.02 %	0.00 %	5.79 %	
3.000	0.15	61.87 %	0.11 %	0.00 %	38.02 %	
4.000	1.14	13.96 %	0.24 %	0.00 %	85.80 %	
5.000	2.08	1.60 %	0.27 %	0.00 %	98.13 %	
6.000	2.72	0.16 %	0.28 %	0.00 %	99.56 %	
6.500	2.84	0.05 %	0.28 %	0.00 %	99.67 %	
7.000	2.89	0.02 %	0.28 %	0.00 %	99.71 %	
7.400	2.91	0.01 %	0.28 %	0.00 %	99.72 %	Blood pH
8.000	2.92	0.00 %	0.28 %	0.00 %	99.72 %	
9.000	2.92	0.00 %	0.28 %	0.00 %	99.72 %	
10.000	2.92	0.00 %	0.28 %	0.00 %	99.72 %	
11.000	2.92	0.00 %	0.28 %	0.00 %	99.72 %	
12.000	2.92	0.00 %	0.28 %	0.00 %	99.72 %	

# Carbonate and acidity



# Other graphs








Pion Assay name: pH-metric high logP Analyst: Instrument ID: T311053 Assay ID: 18C-16016 Filename:

C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-16016 Instrument ID: T311053

 $C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric\ high\ logP.t3r$ Filename:

# Assay Model

Settings	Value	Date/Time changed	Imported from	
Sample name	M13_octanol	2/27/2018 5:57:49 PM	User entered value	
Sample by	Weight		Default value	
Sample weight	0.001360 g	3/16/2018 5:08:03 PM	User entered value	
Formula weight	295.34 g/mol	2/27/2018 5:57:49 PM	User entered value	
Solubility	Unknown		Default value	
Molecular weight	295.34	2/27/2018 5:57:49 PM	User entered value	
Individual pKa ionic environments	No		Default value	
Number of pKas	1	2/27/2018 5:57:49 PM	User entered value	
Sample is a	Base	2/27/2018 5:57:49 PM	User entered value	
pKa 1	5.77	2/27/2018 5:57:49 PM	User entered value	
logp (XH +)	-4.56	3/2/2018 4:30:48 PM	User entered value	
logP (neutral X)	2.99	3/2/2018 4:30:43 PM	User entered value	
Events				

Events	3									
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/e time
4:09.8	Manual volume addition				0.10000 mL					ume
4:10.9 7:10.4	Initial pH = 3.69 Data point 2	1.50000 mL	0.04605 mL	0.00405 mL	0.10000 mL	2.007	0.00659	0.78015	0.00037	10.0
7:56.4	Data point 3	1.50000 mL	0.04605 mL	0.01919 mL	0.10000 mL	2.211	0.00314	0.24211	0.00032	10.0
8:32.1	Data point 4	1.50000 mL	0.04605 mL	0.02949 mL	0.10000 mL	2.417	-0.01053	0.38713	0.00084	10.0
9:07.7	Data point 5	1.50000 mL	0.04605 mL	0.03591 mL	0.10000 mL	2.660	-0.00938	0.73483	0.00054	10.0
9:53.5	Data point 6	1.50000 mL	0.04605 mL	0.03982 mL	0.10000 mL	2.855	-0.00211	0.30085	0.00019	10.0
10:29.0	Data point 7	1.50000 mL	0.04605 mL	0.04236 mL	0.10000 mL	3.139	-0.00782	0.88102	0.00041	10.5
11:20.4	Data point 8	1.50000 mL	0.04605 mL	0.04393 mL	0.10000 mL	3.394	-0.00736	0.41781	0.00056	10.0
11:55.9	Data point 9	1.50000 mL	0.04605 mL	0.04525 mL	0.10000 mL	3.686	-0.01800	0.93937	0.00092	s 12.0
12:43.6	Data point 10	1.50000 mL	0.04605 mL	0.04657 mL	0.10000 mL	3.982	-0.01971	0.97981	0.00098	14.5
13:23.6	Data point 11	1.50000 mL	0.04605 mL	0.04779 mL	0.10000 mL	4.181	-0.01960	0.95605	0.00099	s 17.5
14:06.5	Data point 12	1.50000 mL	0.04605 mL	0.04892 mL	0.10000 mL	4.317	-0.01919	0.90540	0.00100	s 19.0
15:01.2	Data point 13	1.50000 mL	0.04605 mL	0.05047 mL	0.10000 mL	4.588	-0.01947	0.95881	0.00098	36.5
16:13.4	Data point 14	1.50000 mL	0.04605 mL	0.05153 mL	0.10000 mL	5.212	0.01166	0.44692	0.00086	_
17:16.0	Data point 15	1.50000 mL	0.04605 mL	0.05275 mL	0.10000 mL	9.929	-0.01846	0.92278	0.00095	s 43.5
18:57.2	Data point 16	1.50000 mL	0.10329 mL	0.05275 mL	0.30000 mL	1.971	-0.00110	0.01132	0.00051	10.0
19:43.4	Data point 17	1.50000 mL	0.10329 mL	0.07140 mL	0.30000 mL	2.179	0.00116	0.00555	0.00077	
20:19.1	Data point 18	1.50000 mL	0.10329 mL	0.08333 mL	0.30000 mL	2.393	-0.00552	0.08741	0.00092	s 10.0
20:54.7	Data point 19	1.50000 mL	0.10329 mL	0.09073 mL	0.30000 mL	2.610	-0.00664	0.54896	0.00044	10.0
21:30.3	Data point 20	1.50000 mL	0.10329 mL	0.09541 mL	0.30000 mL	2.793	-0.01670	0.91771	0.00086	s 10.5 s
1										



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# **Events (continued)**

Events (	(Continueu)									
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
22:06.3	Data point 21	1.50000 mL	0.10329 mL	0.09875 mL	0.30000 mL	3.048	-0.00688	0.24255	0.00069	10.0 s
22:52.1	Data point 22		0.10329 mL						0.00072	10.0 s
23:37.8	Data point 23		0.10329 mL						0.00052	
24:28.7	Data point 24		0.10329 mL						0.00097	
25:14.5	Data point 25		0.10329 mL						0.00097	
25:50.4	Data point 26		0.10329 mL						0.00100	
27:04.9	Data point 27		0.10329 mL						0.00098	
27:56.3	Data point 28		0.10329 mL						0.00086	
28:40.9	Data point 29		0.10329 mL						0.00081	
29:18.3	Data point 30		0.10329 mL						0.00091	
30:03.2	Data point 31	1.50000 mL	0.10329 mL	0.10882 mL	0.30000 mL	7.326	-0.14623	0.99588	0.00723	Timed out at 59.5 s
31:33.7	Data point 32	1.50000 mL	0.10329 mL	0.10896 mL	0.30000 mL	7.576	-0.08507	0.99600	0.00421	Timed out
33:04.1	Data point 33	1.50000 mL	0.10329 mL	0.10903 mL	0.30000 mL	7.881	-0.07193	0.99419	0.00357	at 59.5 s Timed out
24.24 5	Data naint 24	1 E0000 l	0.40220	0.40000	0.20000!	7.040	0.04500	0.00227	0.00004	at 59.5 s
34:34.5	Data point 34	1.50000 ML	0.10329 mL	0.10908 ML	0.30000 ML	7.949	-0.04523	0.99327	0.00224	Timed out at 59.5 s
36:05.0	Data point 35	1.50000 mL	0.10329 mL	0.10913 mL	0.30000 mL	8.208	-0.03919	0.98798	0.00195	Timed out at 59.5 s
37:40.6	Data point 36	1.50000 mL	0.10329 mL	0.10920 mL	0.30000 mL	8.410	-0.02748	0.97538	0.00137	Timed out at 59.5 s
39:11.0	Data point 37	1.50000 mL	0.10329 mL	0.10927 mL	0.30000 mL	8.578	-0.01916	0.95293	0.00097	
40:32.1	Data point 38		0.10329 mL						0.00095	
41:42.6	Data point 39		0.10329 mL						0.00098	
43:28.9	Data point 40		0.16319 mL						0.00048	
44:15.2	Data point 41	1.50000 mL	0.16319 mL	0.12869 mL	0.80000 mL	2.174	-0.01641	0.68236	0.00098	10.0 s
44:50.9	Data point 42	1.50000 mL	0.16319 mL	0.14189 mL	0.80000 mL	2.382	-0.00586	0.20619	0.00064	10.0 s
45:26.6	Data point 43		0.16319 mL						0.00044	
46:02.1	Data point 44		0.16319 mL						0.00073	10.0 s
46:37.6	Data point 45		0.16319 mL					0.83294	0.00054	
47:13.6	Data point 46		0.16319 mL						0.00086	
47:54.2	Data point 47		0.16319 mL					0.10317	0.00059	
48:39.8	Data point 48		0.16319 mL						0.00066	
49:30.7	Data point 49		0.16319 mL					0.70214	0.00075	
50:21.6	Data point 50		0.16319 mL					0.00324	0.00044	
51:12.9	Data point 51		0.16319 mL					0.79034	0.00051	
51:48.8	Data point 52		0.16319 mL						0.00084	
52:30.8	Data point 53		0.16319 mL						0.00098	
53:26.7 54:23.9	Data point 54		0.16319 mL						0.00093	
	Data point 55		0.16319 mL						0.00099	
55:21.0 56:22.0	Data point 56		0.16319 mL						0.00097 0.00094	
56.22.0 57:30.5	Data point 57 Data point 58		0.16319 mL 0.16319 mL							Timed out
	·									at 59.5 s
59:00.9	Data point 59	1.50000 mL	0.16319 mL	U.168/2 ML	บ.ชบบบบ mL	7.887	-0.09593	0.99713	0.00475	Timed out at 59.5 s
1:00:36.6	Data point 60	1.50000 mL	0.16319 mL	0.16886 mL	0.80000 mL	8.519	-0.01911	0.89364	0.00100	56.5 s
1:02:08.7	Data point 61	1.50000 mL	0.16319 mL	0.16898 mL	0.80000 mL	8.736	-0.01754	0.87325	0.00093	29.0 s
1:03:13.4	Data point 62	1.50000 mL	0.16319 mL	0.16910 mL	0.80000 mL	8.932	-0.01508	0.70146	0.00089	
	Data point 63 Assay volumes		0.16319 mL 0.16319 ml			9.107	-0.01562	0.77923	0.00087	14.5 s
1000.0	55a, 75iaiii66		555 TO THE	5 O D E 1 111L	5.55555 IIIL					



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	9.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings				
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium	110110			
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.100 mL			
Partition solvent added	Manual in advance			
After partition addition, stir for	1 seconds			
Sample Sonication	1 30001103			
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	60 seconds			
After sonication stir for	5 seconds			
Sample Dissolution	0 00001100			
Perform a dissolution stage	Yes			
	To start pH			
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge	10 70			
Perform a carbonate purge	No			
Temperature Control	INU			
Wait for temperature	Yes			
Required start temperature	25.0°C			
· · · · · · · · · · · · · · · · · · ·	0.5°C			
Acceptable deviation Time to wait				
	60 seconds			
Stir speed of Titration 1	50%			
	Low to high nU			
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	50%			
Titration 2	1 ( 12 1 22			
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume				
Additional partition solvent added	Automatic			

Report by: Dorothy Levorse 3/23/2018 1:20:42 PM

30 seconds

15 seconds

55%

After pH adjust stir for

Stir to allow partitioning for

Stirrer speed for partitioning



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from	
Titration 3					
Titrate from	Low to high pH				
Add additional water	0.00 mL				
Additional partition solvent volume	0.500 mL				
Additional partition solvent added	Automatic				
After pH adjust stir for	30 seconds				
Stir to allow partitioning for	15 seconds				
Stirrer speed for partitioning	60%				
Data Point Stability					
Stir during data point collection	No				
Delay before data point collection	0 seconds				
Number of points to average	20 points				
Time interval between points	0.50 seconds				
Required maximum standard deviation	0.00100 dpH/dt				
Stability timeout after	60 seconds				

# Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.167	3/16/2018 8:34:16 PM	C:\Sirius_T3\HCl18C16.t3r
Four-Plus S	0.9932	3/16/2018 8:34:16 PM	C:\Sirius_T3\HCl18C16.t3r
Four-Plus jH	0.7	3/16/2018 8:34:16 PM	C:\Sirius_T3\HCl18C16.t3r
Four-Plus jOH	-0.9	3/16/2018 8:34:16 PM	C:\Sirius_T3\HCl18C16.t3r
Base concentration factor	1.010	3/16/2018 8:34:17 PM	C:\Sirius_T3\KOH18C16.t3r
Acid concentration factor	0.990	3/16/2018 8:34:16 PM	C:\Sirius_T3\HCl18C16.t3r

# Instrument Settings

Setting	Value	Batch Id	Install date
Instrument owner	Merck		
Instrument ID	T311053		
Instrument type	T3 Simulator		
Software version	1.1.3.0		
Dispenser module		T3DM1100253	3/31/2009 6:24:52 AM
Dispenser 0	Water		3/31/2009 6:25:05 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Water (0.15 M KCI)	2-6-18	3/5/2018 11:53:26 AM
Dispenser 2	Acid		3/31/2009 6:25:11 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Acid (0.5 M HCI)	03-16-2018	3/16/2018 11:26:11 AM
Dispenser 1	Base		3/31/2009 6:25:21 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Base (0.5 M KOH)	09-22-17	3/16/2018 11:19:32 AM
Dispenser 5	Cosolvent		3/31/2009 6:26:24 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Distribution valve 5	Distribution Valve		3/31/2009 6:28:19 AM
Firmware version	1.1.3		
Port A	Methanol (80%, 0.15 M KCl)	03-16-2018	3/16/2018 11:40:37 AM
Port B	Cyclohexane		1/16/2018 1:42:16 PM
Port C	MeCN (50%, 0.15 M KCI)	10-30-17	3/16/2018 11:52:45 AM
Dispenser 3	Buffer		8/3/2010 6:05:16 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Dodecane	1-31-2018	2/28/2018 10:07:18 AM




Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Dispenser 6	Octanol		10/22/2010 11:52:43 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Octanol	1-31-2018	1/31/2018 5:24:59 PM
Titrator		T3TM1100153	3/31/2009 6:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0769	8/15/2017 10:21:54 AM
E0 calibration	-10.42 mV		3/16/2018 8:34:44 PM
Filling solution	3M KCI	KCL095	3/16/2018 11:15:26 AM
Liquids			
Wash 1	50% IPA:50% Water		3/16/2018 11:15:44 AM
Wash 2	0.5% Trition X-100 in H20		3/16/2018 11:15:47 AM
Buffer position 1	pH7 Wash		3/16/2018 11:15:49 AM
Buffer position 2	pH 7		3/16/2018 11:15:51 AM
Storage position	4.0 - + 0.00 1	0.07.40	3/16/2018 11:16:23 AM
Wash water	4.3e+003 mL	2-27-18	2/27/2018 10:29:26 AM
Waste	5.9e+003 mL		2/27/2018 10:29:30 AM
Temperature controller			8/5/2010 7:35:13 AM
Turbidity detector		070000	3/31/2009 6:24:45 AM
Spectrometer		072390	11/23/2010 12:22:28 PM
Dip probe	185.563	11086	
Wavelength coefficient A0 Wavelength coefficient A1	2.17439		
Wavelength coefficient A2	-0.000285622		
Total lamp lit time	804:59:49		11/23/2010 12:22:28 PM
Calibrated on	2/22/2018 5:36:24 PM		11/23/2010 12.22.20 FW
Integration time	17		
Scans averaged	10		
Autoloader	10	T3AL1100237	11/10/2015 10:34:13 AM
Left-right axis firmware version	1.17 Al1Dl2DO2 Stepper 2	10/12/10020/	117 107 20 10 10.0 11.10 7 11.1
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Configuration	ğ		
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration E0 calibration buffer wash stir speed	5 s 30%		
Lo campiation bullet wash still speed	JU /0		





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-16016 Instrument ID: T311053

Filename: C:\Sirius\_T3\Mehtap\20180316\_exp32\_logP\_T3-1\18C-16016\_M13\_octanol\_pH-metric high logP.t3r

# Instrument Settings (continued)

Setting	Value	Batch Id	Install date	
E0 calibration reading stir speed	0%			
Spectrometer calibration stir duration	5 s			
Spectrometer calibration stir speed	30%			
Spectrometer calibration wash pump volume	20.0 mL			
Spectrometer calibration wash stir duration	5 s			
Spectrometer calibration wash stir speed	30%			
Overhead dispense height	10000			

# Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00