

Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

18C-24002 Instrument ID: **T312060**

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH +) 0.82 ±0.06 (n=50) logP (neutral X) 4.04 ±0.02 (n=50)

RMSD 0.362

18C-24002 Points 1 to 24

M04_octanol concentration factor 0.923
Carbonate 0.1572 mM
Acidity error -0.23457 mM

18C-24002 Points 25 to 51

M04_octanol concentration factor 0.920
Carbonate 0.1065 mM
Acidity error -0.29812 mM

18C-24002 Points 52 to 75

M04_octanol concentration factor 0.952
Carbonate 0.1039 mM
Acidity error -0.10230 mM

Warnings and errors

Errors None Warnings None

Sample logD and percent species

рн	MU4_octanol	MU4_octanol	MU4_octanol	MU4_octanol	MU4_octanol	Comment
-	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.83	12.93 %	0.00 %	85.56 %	1.52 %	
1.200	0.83	12.81 %	0.00 %	84.80 %	2.38 %	Stomach pH
2.000	0.89	11.37 %	0.00 %	75.28 %	13.34 %	
3.000	1.26	5.17 %	0.01 %	34.21 %	60.62 %	
4.000	2.09	0.80 %	0.01 %	5.30 %	93.89 %	
5.000	3.03	0.08 %	0.01 %	0.56 %	99.35 %	
6.000	3.75	0.01 %	0.01 %	0.06 %	99.93 %	
6.500	3.93	0.00 %	0.01 %	0.02 %	99.97 %	
7.000	4.00	0.00 %	0.01 %	0.01 %	99.98 %	
7.400	4.02	0.00 %	0.01 %	0.00 %	99.99 %	Blood pH
8.000	4.04	0.00 %	0.01 %	0.00 %	99.99 %	•
9.000	4.04	0.00 %	0.01 %	0.00 %	99.99 %	
10.000	4.04	0.00 %	0.01 %	0.00 %	99.99 %	
11.000	4.04	0.00 %	0.01 %	0.00 %	99.99 %	
12.000	4.04	0.00 %	0.01 %	0.00 %	99.99 %	



Sample name: Assay name: Assay ID:

Filename:


M04_octanol pH-metric high logP

18C-24002

Experiment start time: 3/24/2018 1:34:06 AM Analyst: **Dorothy Levorse**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r






Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM Analyst: Assay name: pH-metric high logP **Dorothy Levorse** Assay ID:

18C-24002 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

Graphs (continued)





Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM pH-metric high logP Analyst: Assay name: **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-24002

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-24002 Points 1 to 24

Overall results

RMSD 0.282 Average ionic strength 0.156 M 24.9°C Average temperature Partition ratio 0.0186:1

Analyte concentration range 2528.6 µM to 2603.3 µM

Total points considered 19 of 24

Warnings and errors

Errors None Warnings None

Four-Plus parameters

Alpha 0.119 3/24/2018 1:34:06 AM C:\Sirius_T3\HCl18C23.t3r S 0.9972 3/24/2018 1:34:06 AM C:\Sirius T3\HCl18C23.t3r jΗ 0.9 3/24/2018 1:34:06 AM C:\Sirius_T3\HCl18C23.t3r jOH -0.33/24/2018 1:34:06 AM C:\Sirius_T3\HCl18C23.t3r


Titrants


0.50 M KOH 1.003190 3/24/2018 1:34:06 AM C:\Sirius_T3\KOH18C23.t3r

Sample

M04_octanol concentration factor 0.923 Base pKa 1 5.97 logP(XH +)0.19 logP (neutral X) 3.95

Sample graphs








Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

18C-24002 Instrument ID: **T312060**

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species


рН	M04_octanol	M04_octanol	M04_octanol	M04_octanol	M04_octanol	Comment
	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.22	97.04 %	0.00 %	2.79 %	0.17 %	
1.200	0.23	96.94 %	0.00 %	2.79 %	0.27 %	Stomach pH
2.000	0.40	95.55 %	0.01 %	2.75 %	1.70 %	
3.000	1.05	82.83 %	0.09 %	2.38 %	14.70 %	
4.000	1.98	35.54 %	0.38 %	1.02 %	63.06 %	
5.000	2.94	5.30 %	0.57 %	0.15 %	93.98 %	
6.000	3.66	0.56 %	0.60 %	0.02 %	98.83 %	
6.500	3.84	0.18 %	0.60 %	0.01 %	99.22 %	
7.000	3.91	0.06 %	0.60 %	0.00 %	99.34 %	
7.400	3.93	0.02 %	0.60 %	0.00 %	99.38 %	Blood pH
8.000	3.95	0.01 %	0.60 %	0.00 %	99.39 %	
9.000	3.95	0.00 %	0.60 %	0.00 %	99.40 %	
10.000	3.95	0.00 %	0.60 %	0.00 %	99.40 %	
11.000	3.95	0.00 %	0.60 %	0.00 %	99.40 %	
12.000	3.95	0.00 %	0.60 %	0.00 %	99.40 %	

Carbonate and acidity



Other graphs







Sample name: M04_octanol Assay name:

pH-metric high logP

Assay ID: Filename:

18C-24002


Experiment start time: 3/24/2018 1:34:06 AM

Analyst: **Dorothy Levorse**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

Other graphs (continued)





Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM Analyst: Assay name: pH-metric high logP **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-24002

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-24002 Points 25 to 51

Overall results

RMSD 0.334 Average ionic strength 0.162 M Average temperature 25.0°C Partition ratio 0.0644:1

Analyte concentration range 2267.2 µM to 2336.0 µM

Total points considered 24 of 27

Warnings and errors

Errors None Warnings None

Four-Plus parameters


Alpha 0.119 3/24/2018 1:34:06 AM C:\Sirius_T3\HCl18C23.t3r S 0.9972 3/24/2018 1:34:06 AM C:\Sirius T3\HCl18C23.t3r jΗ 0.9 3/24/2018 1:34:06 AM C:\Sirius_T3\HCl18C23.t3r jOH -0.33/24/2018 1:34:06 AM C:\Sirius_T3\HCl18C23.t3r


Titrants

Sample

M04_octanol concentration factor 0.920 Base pKa 1 5.97 logP(XH +)0.19 logP (neutral X) 4.01

Sample graphs








Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM Assay name: pH-metric high logP Analyst: **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-24002

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r


Sample graphs (continued)



Sample logD and percent species


рН	M04_octanol	M04_octanol	M04_octanol	M04_octanol	M04_octanol	Comment
	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.22	90.35 %	0.00 %	9.01 %	0.64 %	
1.200	0.24	90.02 %	0.00 %	8.97 %	1.01 %	Stomach pH
2.000	0.42	85.44 %	0.01 %	8.52 %	6.03 %	
3.000	1.10	55.36 %	0.06 %	5.52 %	39.07 %	
4.000	2.04	12.24 %	0.13 %	1.22 %	86.40 %	
5.000	3.00	1.39 %	0.15 %	0.14 %	98.32 %	
6.000	3.72	0.14 %	0.15 %	0.01 %	99.69 %	
6.500	3.90	0.04 %	0.15 %	0.00 %	99.80 %	
7.000	3.97	0.01 %	0.15 %	0.00 %	99.83 %	
7.400	3.99	0.01 %	0.15 %	0.00 %	99.84 %	Blood pH
8.000	4.01	0.00 %	0.15 %	0.00 %	99.85 %	
9.000	4.01	0.00 %	0.15 %	0.00 %	99.85 %	
10.000	4.01	0.00 %	0.15 %	0.00 %	99.85 %	
11.000	4.01	0.00 %	0.15 %	0.00 %	99.85 %	
12.000	4.01	0.00 %	0.15 %	0.00 %	99.85 %	


Carbonate and acidity



0.106 mM Acidity error -0.298 mM

Other graphs







Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP

18C-24002


Experiment start time: 3/24/2018 1:34:06 AM

Analyst: **Dorothy Levorse**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

Other graphs (continued)





Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM pH-metric high logP Analyst: Assay name: **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-24002

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-24002 Points 52 to 75

Overall results

RMSD 0.452 Average ionic strength 0.168 M 25.0°C Average temperature Partition ratio 0.2805:1

Analyte concentration range 1766.3 μM to 1811.0 μM

Total points considered 19 of 24

Warnings and errors

Errors None Warnings None

Four-Plus parameters


Alpha 0.119 3/24/2018 1:34:06 AM C:\Sirius_T3\HCl18C23.t3r S 0.9972 3/24/2018 1:34:06 AM C:\Sirius T3\HCl18C23.t3r jΗ 0.9 3/24/2018 1:34:06 AM C:\Sirius_T3\HCl18C23.t3r jOH -0.33/24/2018 1:34:06 AM C:\Sirius_T3\HCl18C23.t3r


Titrants

Sample

M04_octanol concentration factor 0.952 Base pKa 1 5.97 logP(XH +)0.19 logP (neutral X) 3.70

Sample graphs







Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP

18C-24002


Experiment start time: 3/24/2018 1:34:06 AM

Analyst: **Dorothy Levorse**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

Sample graphs (continued)



Sample logD and percent species


рН	M04_octanol	M04_octanol	M04_octanol	M04_octanol	M04_octanol	Comment
	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.20	68.98 %	0.00 %	29.97 %	1.05 %	
1.200	0.21	68.56 %	0.00 %	29.78 %	1.65 %	Stomach pH
2.000	0.32	63.02 %	0.01 %	27.38 %	9.60 %	
3.000	0.84	33.81 %	0.04 %	14.69 %	51.47 %	
4.000	1.74	6.00 %	0.06 %	2.61 %	91.33 %	
5.000	2.69	0.65 %	0.07 %	0.28 %	99.00 %	
6.000	3.42	0.07 %	0.07 %	0.03 %	99.84 %	
6.500	3.59	0.02 %	0.07 %	0.01 %	99.90 %	
7.000	3.67	0.01 %	0.07 %	0.00 %	99.92 %	
7.400	3.69	0.00 %	0.07 %	0.00 %	99.93 %	Blood pH
8.000	3.70	0.00 %	0.07 %	0.00 %	99.93 %	
9.000	3.70	0.00 %	0.07 %	0.00 %	99.93 %	
10.000	3.70	0.00 %	0.07 %	0.00 %	99.93 %	
11.000	3.70	0.00 %	0.07 %	0.00 %	99.93 %	
12.000	3.70	0.00 %	0.07 %	0.00 %	99.93 %	

Carbonate and acidity



Carbonate 0.104 mM Acidity error -0.102 mM

Other graphs





Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP


18C-24002


Experiment start time: 3/24/2018 1:34:06 AM

Analyst: **Dorothy Levorse**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r







Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM Analyst: Dorothy Levorse

18C-24002 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

Assay Model

Assay ID:

Settings	Value	Date/Time changed	Imported from
Sample name	M04_octanol	3/9/2018 4:34:21 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.001110 g	3/23/2018 5:00:32 PM	User entered value
Formula weight	269.73 g/mol	3/9/2018 4:34:21 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	269.73	3/9/2018 4:34:21 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	1	3/9/2018 4:34:21 PM	User entered value
Sample is a	Base	3/9/2018 4:34:21 PM	User entered value
pKa 1	5.97	3/9/2018 4:34:21 PM	User entered value
logp (XH +)	0.19	3/9/2018 4:34:33 PM	User entered value
logP (neutral X)	3.50	3/23/2018 2:32:46 PM	User entered value

Events

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
5:00.2	Initial pH = 9.63	4 50000	0.04000 !	0.00400 1	0.000041	0.044	0.00400	0.00040	0.00044	40.5
7:59.7	Data point 1		0.04892 mL						0.00014	
8:46.4	Data point 2		0.04892 mL						0.00010	
9:21.9	Data point 3		0.04892 mL					0.05785	0.00060	
9:57.5	Data point 4		0.04892 mL						0.00035	
10:33.1			0.04892 mL						0.00027	
11:08.6	•		0.04892 mL						0.00087	
11:44.0	Data point 7		0.04892 mL						0.00017	
12:34.9			0.04892 mL						0.00073	
	Data point 9		0.04892 mL						0.00092	
14:11.6			0.04892 mL						0.00090	
	Data point 11		0.04892 mL						0.00050	
15:22.5	Data point 12		0.04892 mL					0.81314	0.00048	
	Data point 13		0.04892 mL						0.00077	
	Data point 14		0.04892 mL						0.00075	
17:22.1	Data point 15	1.50000 mL	0.04892 mL	0.04711 mL	0.03001 mL	4.865	-0.01667	0.86899	0.00088	14.5 s
18:12.3	Data point 16	1.50000 mL	0.04892 mL	0.04751 mL	0.03001 mL	5.110	-0.01586	0.79314	0.00088	16.0 s
18:53.6	Data point 17	1.50000 mL	0.04892 mL	0.04770 mL	0.03001 mL	5.348	-0.01986	0.98346	0.00099	20.5 s
19:44.8	Data point 18	1.50000 mL	0.04892 mL	0.04795 mL	0.03001 mL	5.739	-0.01915	0.96564	0.00096	26.5 s
20:36.6	Data point 19	1.50000 mL	0.04892 mL	0.04807 mL	0.03001 mL	5.946	-0.01903	0.92495	0.00098	30.5 s
21:37.5	Data point 20	1.50000 mL	0.04892 mL	0.04819 mL	0.03001 mL	6.643	-0.02497	0.94842	0.00127	Timed out at
	·									59.5 s
23:08.1	Data point 21	1.50000 mL	0.04892 mL	0.04833 mL	0.03001 mL	8.221	-0.07094	0.99161	0.00352	Timed out at
	'									59.5 s
24:43.7	Data point 22	1.50000 mL	0.04892 mL	0.04840 mL	0.03001 mL	8.555	-0.03178	0.97425	0.00159	Timed out at
										59.5 s
26:19.4	Data point 23	1.50000 mL	0.04892 mL	0.04847 mL	0.03001 mL	8.794	-0.01850	0.95661	0.00093	
	Data point 24		0.04892 mL						0.00100	
	Data point 25		0.10310 mL						0.00088	
	Data point 26		0.10310 mL						0.00096	
	Data point 27		0.10310 mL						0.00024	
	Data point 28		0.10310 mL					0.08086	0.00098	
32:05.1	•		0.10310 mL						0.00017	
	Data point 30		0.10310 mL						0.00029	
	Data point 31		0.10310 mL						0.00024	
	Data point 32		0.10310 mL						0.00024	
	Data point 33		0.10310 mL						0.00063	
35:44.5			0.10310 mL						0.00037	
	Data point 35		0.10310 mL				-0.00653		0.00037	
ა≀.∪ხ. I	Data point 36	TITI DUUUC.1	0.10310 mL	U. 1003 I IIIL	U. 11002 IIIL	4.1/8	-0.00915	0.01492	0.00058	10.0 8

1.50000 mL 0.10310 mL 0.10071 mL 0.11002 mL 4.385 -0.00988 0.57323

37:41.5 Data point 37

0.00064 10.0 s



Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM Assay name: pH-metric high logP Analyst: **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-24002

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

Events (continued)

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
38:16.9	Data point 38	1.50000 mL	0.10310 mL	0.10101 mL	0.11002 mL	4.612	-0.00861	0.68608	0.00051	10.5 s
38:52.8	Data point 39		0.10310 mL					0.66939	0.00074	
39:28.7	Data point 40		0.10310 mL						0.00080	
40:05.1	Data point 41		0.10310 mL					0.59589	0.00074	
40:47.6	Data point 42		0.10310 mL						0.00092	
41:31.6	Data point 43		0.10310 mL						0.00094	
42:33.2	Data point 44		0.10310 mL						0.00097	
43:51.7	Data point 45		0.10310 mL					0.98601	0.00209	
10.01.7	Bata point 10	1.00000 1112	0.100101112	0.10100 IIIL	0.11002 IIIE	1.001	0.01100	0.00001	0.00200	at 59.5 s
45:22.1	Data point 46	1 50000 ml	0.10310 mL	0 10174 ml	0 11002 ml	7 209	-0.04636	0 98987	0.00230	
TO.ZZ. 1	Data point 40	1.50000 1112	0.10010 IIIL	0.1017 + IIIL	0.11002 IIIL	1.200	-0.04000	0.30307	0.00200	at 59.5 s
46:52.7	Data point 47	1 50000 ml	0.10310 mL	0 10170 ml	0 11002 ml	7 9/17	0.06750	0.00607	0.00335	Timed out
40.32.7	Data point 47	1.50000 IIIL	0.10310111L	0.10179111L	0.11002 IIIL	1.041	-0.00739	0.99007	0.00333	at 59.5 s
48:23.2	Data point 48	1 50000 ml	0.10310 mL	0 10193 ml	0 11002 ml	8 503	0.03401	0.00546	0.00169	Timed out
40.23.2	Data point 46	1.50000 IIIL	0.10310111L	0.10163 IIIL	0.11002 IIIL	0.503	-0.03401	0.99540	0.00100	at 59.5 s
40·E2 7	Data point 40	1 50000 ml	0.10210 ml	0 10100 ml	0 11002 ml	0 600	0.01520	0.60040	0.00007	
49:53.7	Data point 49		0.10310 mL					0.60040	0.00097	
51:11.7	Data point 50		0.10310 mL				-0.01850	0.93390	0.00095	
52:34.3	Data point 51		0.10310 mL					0.96553	0.00064	
54:34.0	Data point 52		0.16030 mL					0.52317		
55:20.3	Data point 53		0.16030 mL					0.84480	0.00080	
55:56.0	Data point 54		0.16030 mL					0.49233	0.00021	
56:32.0	Data point 55		0.16030 mL					0.59238	0.00090	
57:07.6	Data point 56		0.16030 mL					0.08785	0.00054	
57:43.1	Data point 57		0.16030 mL					0.10471	0.00040	
58:18.5	Data point 58		0.16030 mL					0.09026	0.00037	10.0 s
59:04.3	Data point 59	1.50000 mL	0.16030 mL	0.15433 mL	0.51002 mL	3.351	-0.00880	0.28962	0.00081	10.0 s
59:39.8	Data point 60	1.50000 mL	0.16030 mL	0.15600 mL	0.51002 mL	3.582	-0.01386	0.91128	0.00072	10.5 s
1:00:15.8	Data point 61	1.50000 mL	0.16030 mL	0.15722 mL	0.51002 mL	3.861	-0.01252	0.59534	0.00080	10.0 s
1:00:56.3	Data point 62	1.50000 mL	0.16030 mL	0.15779 mL	0.51002 mL	4.058	-0.00480	0.17045	0.00057	10.0 s
1:01:42.1	Data point 63	1.50000 mL	0.16030 mL	0.15818 mL	0.51002 mL	4.256	-0.00431	0.26492	0.00041	10.0 s
1:02:32.9	Data point 64	1.50000 mL	0.16030 mL	0.15849 mL	0.51002 mL	4.464	0.00183	0.00838	0.00099	18.5 s
1:03:16.7	Data point 65	1.50000 mL	0.16030 mL	0.15873 mL	0.51002 mL	4.850	-0.01054	0.48317	0.00075	10.5 s
1:03:57.8	Data point 66		0.16030 mL					0.55670	0.00071	12.0 s
1:04:40.3			0.16030 mL						0.00098	
1:05:25.7	Data point 68	1.50000 mL		0.15910 mL				0.97545	0.00092	
1:06:50.3	Data point 69	1.50000 mL		0.15924 mL			-0.10642		0.00528	
	Bata point oo		0.10000 1112	0.100211112	0.010021112		000.2	0.00201	0.00020	at 59.5 s
1:08:25.9	Data point 70	1 50000 ml	0.16030 mL	0 15931 ml	0.51002 ml	8 203	-0 07202	0.98050	0.00359	Timed out
1.00.20.0	Bata point 10	1.00000 1112	0.10000 IIIL	0.100011112	0.010021112	0.200	0.07202	0.00000	0.00000	at 59.5 s
1.00.56 /	Data point 71	1 50000 ml	0.16030 mL	0 15036 ml	0.51002 ml	8 406	-0.04547	0 07023	0.00227	
1.03.50.4	Data point 7 i	1.50000 IIIL	0.10030 IIIL	0.15550 IIIL	0.51002 IIIL	0.700	-0.0-0-1	0.91925	0.00221	at 59.5 s
1.11.26 0	Data point 72	1 50000 ml	0.16030 mL	0.150/11 ml	0.51002 ml	8 600	0.01914	0.91212	0.00099	
			0.16030 mL							Timed out
1.12.39.4	Data point 73	1.50000 IIIL	0. 10030 IIIL	0.10945 IIIL	0.5 1002 IIIL	0.714	-0.02109	0.00709	0.00110	
1.11.15 1	Data point 74	1 50000	0.16020!	0.15052	0.51000 m	0.000	0.04740	0.72252	0.00000	at 59.5 s
	Data point 74		0.16030 mL						0.00099	
	Data point 75		0.16030 mL			ყ.სპნ	-0.01441	0.58124	0.00093	19.5 8

1:15:40.3 Assay volumes 1.50000 mL 0.16030 mL 0.15957 mL 0.51002 mL



Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM Assay name: pH-metric high logP **Dorothy Levorse**

Instrument ID: Assay ID: 18C-24002 T312060

Filename: C:\Sirius_T3\Meh	tap\20180323_exp3			312060 I_pH-metric high logF
Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Dorothy Levorse			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	9.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings	, ,			
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium				
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.030 mL			
Partition solvent added	Automatic			
After partition addition, stir for	1 seconds			
Sample Sonication				
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	60 seconds			
After sonication stir for	5 seconds			
Sample Dissolution				
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution				
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge				
Perform a carbonate purge	No			
Tomporature Control				

Temperature Control

Wait for temperature Yes Required start temperature 25.0°C Acceptable deviation 0.5°C Time to wait 60 seconds 50%

Stir speed of Titration 1

Titrate from Low to high pH

Adjust to start pH Yes

After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds Stirrer speed for partitioning 50%

Titration 2

Titrate from Low to high pH Add additional water 0.00 mL Additional partition solvent volume 0.080 mL Additional partition solvent added Automatic After pH adjust stir for 30 seconds

Stir to allow partitioning for 15 seconds Stirrer speed for partitioning 55%

Report by: Dorothy Levorse 3/26/2018 11:48:40 AM



Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM Analyst: Dorothy Levorse

Assay ID: 18C-24002 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

Assay Settings (continued)

Value	Original Value	Date/Time changed	Imported from
	_	_	•
Low to high pH			
0.00 mL			
0.400 mL			
Automatic			
30 seconds			
15 seconds			
60%			
No			
0 seconds			
20 points			
0.50 seconds			
0.00100 dpH/dt			
60 seconds			
	Low to high pH 0.00 mL 0.400 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt	Low to high pH 0.00 mL 0.400 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt	Low to high pH 0.00 mL 0.400 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.119	3/24/2018 1:34:06 AM	C:\Sirius_T3\HCl18C23.t3r
Four-Plus S	0.9972	3/24/2018 1:34:06 AM	C:\Sirius_T3\HCl18C23.t3r
Four-Plus jH	0.9	3/24/2018 1:34:06 AM	C:\Sirius_T3\HCl18C23.t3r
Four-Plus jOH	-0.3	3/24/2018 1:34:06 AM	C:\Sirius_T3\HCl18C23.t3r
Base concentration factor	1.003	3/24/2018 1:34:06 AM	C:\Sirius_T3\KOH18C23.t3r
Acid concentration factor	0.997	3/24/2018 1:34:06 AM	C:\Sirius_T3\HCl18C23.t3r

Instrument Settings

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module Dispenser 0 Syringe volume Firmware version	Water 2.5 mL 1.2.1(r2)	T3DM1200361	3/31/2009 6:24:52 AM 3/31/2009 6:25:05 AM
Titrant Dispenser 2 Syringe volume Firmware version	Water (0.15 M KCI) Acid 0.5 mL 1.2.1(r2)	02-06-2018	3/16/2018 11:09:18 AM 3/31/2009 6:25:11 AM
Titrant Dispenser 1 Syringe volume Firmware version	Acid (0.5 M HCI) Base 0.5 mL 1.2.1(r2)	03-16-2018	3/16/2018 10:56:23 AM 3/31/2009 6:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2)	3/22/2018	3/23/2018 9:34:17 AM 3/31/2009 6:26:24 AM
Distribution valve 5 Firmware version Port A	Distribution Valve 1.1.3	02.09.2019	3/31/2009 6:28:19 AM
Port A Port B Dispenser 3 Syringe volume Firmware version	Methanol (80%, 0.15 M KCI) Cyclohexane Buffer 0.5 mL 1.2.1(r2)	02-08-2018 11-01-17	3/6/2018 10:28:59 AM 2/27/2018 11:37:57 AM 8/3/2010 6:05:16 AM
Titrant Dispenser 6	Dodecane Octanol	2018/01/31	2/28/2018 11:18:04 AM 10/22/2010 11:52:43 AM




Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM Analyst: Dorothy Levorse

Assay ID: 18C-24002 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting Syringe volume	Value 0.5 mL	Batch Id	Install date
Firmware version Titrant Titrator	1.2.1(r2) Octanol	01-31-2018 T3TM1200161	2/27/2018 10:59:35 AM 3/31/2009 6:24:17 AM
Horizontal axis firmware version Vertical axis firmware version Chassis I/O firmware version Probe I/O firmware version	1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2 1.11 Al1Dl0DO4 Norgren I/O 1.1.1	131W1230101	0/01/2000 0.24. 17 AW
Electrode E0 calibration Filling solution	T3 Electrode +4.89 mV 3M KCI	T3E0923 KCL097	1/23/2018 3:01:00 PM 3/24/2018 1:34:34 AM 3/23/2018 9:29:07 AM
Liquids Wash 1 Wash 2 Buffer position 1	50% IPA:50% Water 0.5% Trition X-100 in H20 pH7 Wash		3/23/2018 9:29:12 AM 3/23/2018 9:29:15 AM 3/23/2018 9:29:19 AM
Buffer position 2 Storage position Wash water Waste	pH 7 7.7e+003 mL 2.5e+003 mL	03-12-2018	3/23/2018 9:29:21 AM 3/23/2018 9:30:23 AM 3/12/2018 9:25:04 AM 3/12/2018 9:24:49 AM
Temperature controller Turbidity detector Spectrometer		074811	8/5/2010 7:35:13 AM 3/31/2009 6:24:45 AM 11/23/2010 12:22:28 PM
Dip probe Wavelength coefficient A0 Wavelength coefficient A1 Wavelength coefficient A2	183.333 2.21568 -0.000289308	10196	
Total lamp lit time Calibrated on Integration time	162:53:01 2/27/2018 11:40:38 AM 40		11/23/2010 12:22:28 PM
Scans averaged Autoloader Left-right axis firmware version Front-back axis firmware version	10 1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2	T3AL1200345	11/10/2015 10:34:13 AM
Vertical axis firmware version Chassis I/O firmware version Configuration	1.17 Al1Dl2DO2 Stepper 2 1.11 Al1Dl0DO4 Norgren I/O		
Alternate titration position Alternate reference position Maximum standard vial volume	Titration position Reference position 3.50 mL		
Maximum alternate vial volume Automatic action idle period Titrant tube volume Syringe flush count	25.00 mL 5 minute(s) 1.3 mL 3.50		
Flowing wash pump volume Flowing wash stir duration Flowing wash stir speed Solvent wash stir duration	20.0 mL 5 s 30% 5 s		
Solvent wash stir speed Surfactant wash stir duration Surfactant wash stir speed	30% 5 s 30%		
E0 calibration minimum number of points E0 calibration maximum standard deviation E0 calibration timeout period E0 calibration stir duration	10 0.01500 60 s 5 s		
E0 calibration preparation stir speed E0 calibration buffer wash stir duration E0 calibration buffer wash stir speed	30% 5 s 30%		
E0 calibration reading stir speed	0%		





Sample name: M04_octanol Experiment start time: 3/24/2018 1:34:06 AM Analyst: Dorothy Levorse

18C-24002 Instrument ID: **T312060**

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24002_M04_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	e 20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		
,			

Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00