

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH2 +) -1.34 logP (neutral XH) 3.04 logP (X -) -1.62

18C-01008 Points 1 to 32

M16_octanol concentration factor 0.795
Carbonate 0.0308 mM
Acidity error 3.10957 mM

18C-01008 Points 33 to 60

M16_octanol concentration factor 0.841
Carbonate 0.2347 mM
Acidity error 2.72359 mM

18C-01008 Points 61 to 95

M16_octanol concentration factor 0.890
Carbonate 0.3111 mM
Acidity error 2.55032 mM

Warnings and errors

Errors None

Warnings User stopped refinement before final result convergence occured

Sample logD and percent species

рн	M16_octanoi	M16_octanoi	M16_octanoi	M16_octanoi	M16_octanoi	M16_octanoi	M16_octanoi	Comment
	logD	M16_octanolH2	M16_octanolH	M16_octanol	M16_octanolH2*	M16_octanolH*	M16_octanol*	
1.000	-0.58	79.04 %	0.02 %	0.00 %	3.57 %	17.37 %	0.00 %	
1.200	-0.41	71.75 %	0.02 %	0.00 %	3.24 %	24.99 %	0.00 %	Stomach pH
2.000	0.35	30.82 %	0.06 %	0.00 %	1.39 %	67.73 %	0.00 %	P · · ·
3.000	1.33	4.34 %	0.09 %	0.00 %	0.20 %	95.38 %	0.00 %	
4.000	2.26	0.45 %	0.09 %	0.00 %	0.02 %	99.44 %	0.00 %	
5.000	2.87	0.05 %	0.09 %	0.00 %	0.00 %	99.86 %	0.00 %	
6.000	3.02	0.00 %	0.09 %	0.00 %	0.00 %	99.90 %	0.00 %	
6.500	3.03	0.00 %	0.09 %	0.00 %	0.00 %	99.91 %	0.00 %	
7.000	3.03	0.00 %	0.09 %	0.00 %	0.00 %	99.91 %	0.00 %	
7.400	3.03	0.00 %	0.09 %	0.00 %	0.00 %	99.91 %	0.00 %	Blood pH
8.000	2.99	0.00 %	0.09 %	0.01 %	0.00 %	99.90 %	0.00 %	-
9.000	2.71	0.00 %	0.09 %	0.10 %	0.00 %	99.80 %	0.00 %	
10.000	1.95	0.00 %	0.09 %	1.03 %	0.00 %	98.85 %	0.02 %	
11.000	0.98	0.00 %	0.08 %	9.41 %	0.00 %	90.28 %	0.23 %	
12.000	-0.01	0.00 %	0.04 %	50.40 %	0.00 %	48.35 %	1.21 %	



Assay ID:

Sample name: M16_octanol Assay name:


pH-metric high logP

18C-01008 Filename:

Experiment start time: 3/1/2018 10:41:05 AM

Pion Analyst: Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r






Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-01008 Points 1 to 32

Overall results

RMSD 0.203
Average ionic strength 0.159 M
Average temperature 25.0°C
Partition ratio 0.0247 : 1

Analyte concentration range 4379.7 µM to 4702.9 µM

Total points considered 23 of 32

Warnings and errors

Errors None

Warnings Excessive acidity error present

Four-Plus parameters

à	Alpha	0.130	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
à	S	0.9970	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
à	jΗ	8.0	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
à	jОН	-0.4	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r

Titrants

0.50 M HCI 0.993513 3/1/2018 10:41:05 AM C:\Sirius_T3\HCl18B27.t3r 0.50 M KOH 0.999845 3/1/2018 10:41:05 AM C:\Sirius_T3\KOH18B27.t3r

Sample

 ▼ M16_octanol concentration factor
 0.795

 ♣ Base pKa 1
 4.70


 ♣ Acid pKa 2
 8.94


 ♣ logP (XH2 +)
 -0.95

 logP (neutral XH)
 3.01

 ♣ logP (X -)
 -1.62

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M16_octanol		M16_octanol	M16_octanol		M16_octanol	M16_octanol	
	logD	M16_octanolH2	M16_octanolH	M16_octanol	M16_octanolH2*	M16_octanolH*	M16_octanol*	
1.000	-0.50	99.20 %	0.02 %	0.00 %	0.28 %	0.50 %	0.00 %	
1.200	-0.36	98.90 %	0.03 %	0.00 %	0.27 %	0.79 %	0.00 %	Stomach pH
2.000	0.33	94.77 %	0.19 %	0.00 %	0.26 %	4.78 %	0.00 %	
3.000	1.30	65.49 %	1.31 %	0.00 %	0.18 %	33.02 %	0.00 %	
4.000	2.23	16.01 %	3.19 %	0.00 %	0.04 %	80.75 %	0.00 %	
5.000	2.83	1.87 %	3.73 %	0.00 %	0.01 %	94.39 %	0.00 %	
6.000	2.99	0.19 %	3.80 %	0.00 %	0.00 %	96.01 %	0.00 %	
6.500	3.00	0.06 %	3.80 %	0.01 %	0.00 %	96.12 %	0.00 %	
7.000	3.00	0.02 %	3.80 %	0.04 %	0.00 %	96.13 %	0.00 %	
7.400	3.00	0.01 %	3.80 %	0.11 %	0.00 %	96.08 %	0.00 %	Blood pH
8.000	2.96	0.00 %	3.79 %	0.44 %	0.00 %	95.77 %	0.00 %	•
9.000	2.68	0.00 %	3.65 %	4.19 %	0.00 %	92.16 %	0.00 %	ļ
10.000	1.91	0.00 %	2.65 %	30.40 %	0.00 %	66.93 %	0.02 %	ļ
11.000	0.95	0.00 %	0.71 %	81.34 %	0.00 %	17.91 %	0.05 %	ļ
12.000	-0.04	0.00 %	0.09 %	97.71 %	0.00 %	2.15 %	0.06 %	

Carbonate and acidity



Other graphs






pH-metric high logP Analyst: **Pion** Assay name: 18C-01008 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-01008 Points 33 to 60

Overall results

RMSD 0.036
Average ionic strength 0.171 M
Average temperature 25.0°C
Partition ratio 0.0751 : 1

Analyte concentration range 3612.7 µM to 3864.7 µM

Total points considered 19 of 28

Warnings and errors

Errors None

Warnings Excessive acidity error present

Four-Plus parameters


à	Alpha	0.130	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
à	S	0.9970	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
à	jΗ	8.0	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
à	jОН	-0.4	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r


Titrants

Sample

Φ	M16_octanol concentration factor	0.841
	Base pKa 1	4.70
	Acid pKa 2	8.94
	logP (XH2 +)	0.31
₩	logP (neutral XH)	3.10
	logP (X -)	-1.62

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	
	logD	M16_octanolH2	M16_octanolH	M16_octanol	M16_octanolH2*	M16_octanolH*	M16_octanol*	
1.000	0.36	85.40 %	0.02 %	0.00 %	12.96 %	1.62 %	0.00 %	
1.200	0.38	84.59 %	0.03 %	0.00 %	12.83 %	2.55 %	0.00 %	Stomach pH
2.000	0.66	74.43 %	0.15 %	0.00 %	11.29 %	14.13 %	0.00 %	•
3.000	1.43	32.57 %	0.65 %	0.00 %	4.94 %	61.84 %	0.00 %	
4.000	2.33	4.92 %	0.98 %	0.00 %	0.75 %	93.36 %	0.00 %	
5.000	2.93	0.52 %	1.03 %	0.00 %	0.08 %	98.37 %	0.00 %	
6.000	3.08	0.05 %	1.04 %	0.00 %	0.01 %	98.90 %	0.00 %	
6.500	3.09	0.02 %	1.04 %	0.00 %	0.00 %	98.94 %	0.00 %	
7.000	3.10	0.01 %	1.04 %	0.01 %	0.00 %	98.94 %	0.00 %	
7.400	3.09	0.00 %	1.04 %	0.03 %	0.00 %	98.93 %	0.00 %	Blood pH
8.000	3.06	0.00 %	1.04 %	0.12 %	0.00 %	98.84 %	0.00 %	
9.000	2.77	0.00 %	1.03 %	1.18 %	0.00 %	97.79 %	0.00 %	
10.000	2.01	0.00 %	0.93 %	10.66 %	0.00 %	88.39 %	0.02 %	
11.000	1.04	0.00 %	0.47 %	54.37 %	0.00 %	45.06 %	0.10 %	
12.000	0.05	0.00 %	0.08 %	92.12 %	0.00 %	7.64 %	0.17 %	

Carbonate and acidity



Other graphs






Assay name: pH-metric high logP Analyst: **Pion** 18C-01008 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-01008 Points 61 to 95

Overall results

RMSD 0.274
Average ionic strength 0.182 M
Average temperature 25.0°C
Partition ratio 0.2049 : 1

Analyte concentration range 2812.2 µM to 2977.4 µM

Total points considered 26 of 35

Warnings and errors

Errors None

Warnings One or more logP values out of range

Excessive acidity error present

Four-Plus parameters

Alpha	0.130	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
S	0.9970	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCI18B27.t3r
jΗ	8.0	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
jОН	-0.4	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r


Titrants


0.50 M HCI 0.993513 3/1/2018 10:41:05 AM C:\Sirius_T3\HCl18B27.t3r 0.50 M KOH 0.999845 3/1/2018 10:41:05 AM C:\Sirius_T3\KOH18B27.t3r

Sample

•	M16_octanol concentration factor	0.890
	Base pKa 1	4.70
	Acid pKa 2	8.94
	logP (XH2 +)	-3.82
	logP (neutral XH)	3.06
	logP (X -)	-1.62

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	
	logD	M16_octanolH2	M16_octanolH	M16_octanol	M16_octanolH2*	M16_octanolH*	M16_octanol*	
1.000	-0.64	95.51 %	0.02 %	0.00 %	0.00 %	4.47 %	0.00 %	
1.200	-0.44	93.06 %	0.03 %	0.00 %	0.00 %	6.91 %	0.00 %	Stomach pH
2.000	0.36	68.01 %	0.14 %	0.00 %	0.00 %	31.85 %	0.00 %	•
3.000	1.35	17.53 %	0.35 %	0.00 %	0.00 %	82.12 %	0.00 %	
4.000	2.28	2.08 %	0.42 %	0.00 %	0.00 %	97.50 %	0.00 %	
5.000	2.88	0.21 %	0.42 %	0.00 %	0.00 %	99.36 %	0.00 %	
6.000	3.04	0.02 %	0.42 %	0.00 %	0.00 %	99.55 %	0.00 %	
6.500	3.05	0.01 %	0.42 %	0.00 %	0.00 %	99.57 %	0.00 %	
7.000	3.05	0.00 %	0.42 %	0.00 %	0.00 %	99.57 %	0.00 %	
7.400	3.05	0.00 %	0.42 %	0.01 %	0.00 %	99.56 %	0.00 %	Blood pH
8.000	3.01	0.00 %	0.42 %	0.05 %	0.00 %	99.53 %	0.00 %	·
9.000	2.73	0.00 %	0.42 %	0.48 %	0.00 %	99.09 %	0.00 %	
10.000	1.96	0.00 %	0.40 %	4.64 %	0.00 %	94.93 %	0.02 %	
11.000	1.00	0.00 %	0.28 %	32.70 %	0.00 %	66.85 %	0.16 %	
12.000	0.01	0.00 %	0.07 %	82.63 %	0.00 %	16.89 %	0.41 %	

Carbonate and acidity



Other graphs






Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name	M16_octanol	2/27/2018 6:14:13 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.001570 g	2/28/2018 4:25:47 PM	User entered value
Formula weight	210.23 g/mol	2/27/2018 5:08:55 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	210.23	2/27/2018 5:08:55 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	2	2/27/2018 5:08:55 PM	User entered value
Sample is a	Ampholyte	2/27/2018 5:08:55 PM	User entered value
pKa 1	4.70	2/27/2018 5:08:55 PM	User entered value
Туре	Base	2/27/2018 5:08:55 PM	
pKa 2	8.94	2/27/2018 5:08:55 PM	User entered value
Туре	Acid	2/27/2018 5:08:55 PM	
logp (XH2 +)	-0.97	2/27/2018 5:09:30 PM	User entered value
logP (neutral XH)	1.47	2/27/2018 5:09:44 PM	User entered value
logP (X -)	-1.62	2/27/2018 5:09:50 PM	User entered value

Events

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
9:20.9	Initial pH = 3.61									•
	Data point 1			0.00275 mL			-0.00290		0.00017	
	Data point 2			0.01799 mL			0.00158	0.11489	0.00023	
	Data point 3			0.02752 mL			0.00139	0.25998	0.00013	
	Data point 4			0.03368 mL			-0.00059		0.00009	
15:09.4	Data point 5			0.03805 mL			-0.00565		0.00086	
15:44.9	Data point 6			0.04080 mL			-0.00330		0.00021	10.0 s
	Data point 7			0.04353 mL			-0.00387		0.00023	
	Data point 8			0.04567 mL			-0.00110		0.00033	
	Data point 9			0.04793 mL			-0.01147		0.00099	
19:19.1	Data point 10			0.05000 mL			-0.00492		0.00028	
	Data point 11			0.05148 mL			-0.01175		0.00081	
	Data point 12			0.05247 mL			-0.01315		0.00079	
21:26.7	Data point 13			0.05322 mL			-0.00889		0.00053	
22:02.7	Data point 14			0.05407 mL			-0.01660		0.00092	
22:39.6	Data point 15	1.50000 mL	0.04518 mL	0.05489 mL	0.04000 mL	7.954	-0.05119	0.98315	0.00255	Timed out at 59.5 s
24:04.9	Data point 16	1.50000 mL	0.04518 mL	0.05506 mL	0.04000 mL	8.206	-0.02013	0.98810	0.00100	38.0 s
	Data point 17	1.50000 mL	0.04518 mL	0.05534 mL	0.04000 mL	8.859	-0.01483	0.91186	0.00077	10.5 s
	Data point 18	1.50000 mL	0.04518 mL	0.05611 mL	0.04000 mL	9.410	-0.01450	0.98052	0.00072	10.5 s
26:20.3	Data point 19	1.50000 mL	0.04518 mL	0.05696 mL	0.04000 mL	9.691	-0.01710	0.91406	0.00088	10.0 s
27:06.1	Data point 20	1.50000 mL	0.04518 mL	0.05804 mL	0.04000 mL	9.876	-0.01009	0.92731	0.00052	10.0 s
27:57.0	Data point 21	1.50000 mL	0.04518 mL	0.05915 mL	0.04000 mL	10.060	-0.00874	0.33867	0.00074	10.0 s
	Data point 22	1.50000 mL	0.04518 mL	0.06018 mL	0.04000 mL	10.223	-0.01285	0.40934	0.00099	10.0 s
29:23.5	Data point 23	1.50000 mL	0.04518 mL	0.06159 mL	0.04000 mL	10.406	-0.00727	0.93231	0.00037	10.5 s
	Data point 24	1.50000 mL	0.04518 mL	0.06329 mL	0.04000 mL	10.593	-0.00688	0.41094	0.00053	10.0 s
31:06.0	Data point 25	1.50000 mL	0.04518 mL	0.06543 mL	0.04000 mL	10.776	-0.01364	0.90456	0.00071	10.5 s
31:57.5	Data point 26	1.50000 mL	0.04518 mL	0.06839 mL	0.04000 mL	10.961	-0.01272	0.56692	0.00083	10.0 s
32:43.3	Data point 27	1.50000 mL	0.04518 mL	0.07150 mL	0.04000 mL	11.147	0.00312	0.33680	0.00027	10.0 s
33:34.4	Data point 28	1.50000 mL	0.04518 mL	0.07596 mL	0.04000 mL	11.337	-0.00151	0.26906	0.00014	10.0 s
34:20.4	Data point 29	1.50000 mL	0.04518 mL	0.08177 mL	0.04000 mL	11.533	-0.00131	0.26464	0.00013	10.0 s
34:56.1	Data point 30	1.50000 mL	0.04518 mL	0.09099 mL	0.04000 mL	11.712	0.01582	0.69550	0.00094	10.0 s
35:31.9	Data point 31	1.50000 mL	0.04518 mL	0.10525 mL	0.04000 mL	11.903	-0.01334	0.54570	0.00089	10.0 s
	Data point 32	1.50000 mL	0.04518 mL	0.11997 mL	0.04000 mL	12.036	-0.00001	0.00000	0.00060	10.0 s
	Data point 33			0.11997 mL			-0.01011	0.31324	0.00089	10.0 s
38:11.4	Data point 34			0.13967 mL			-0.01017	0.51336	0.00070	10.0 s
38:47.1	Data point 35			0.15202 mL			0.01692	0.77326	0.00095	



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Events (continued)

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	nH SD	dpH/dt time
39:27.2	Data point 36						-0.00131			•
40:02.8	Data point 37						0.01355	0.83353		
40:53.9	Data point 38						-0.00343			
40.55.9 41:39.9	Data point 39						0.00068	0.03293	0.00074	
							0.00068	0.0270 4 0.52106		
42:25.7	Data point 40									
43:12.1	Data point 41						0.00068	0.03379		
43:57.8	Data point 42							0.15775		
44:39.8			0.17239 mL				-0.00521			
45:15.3	Data point 44							0.23890		
45:51.2	Data point 45	1.50000 mL	0.17239 mL	0.18203 mL	0.14000 mL	6.710	-0.05359	0.98654	0.00266	Timed out at 59.5 s
47:21.8	Data point 46	1.50000 mL	0.17239 mL	0.18234 mL	0.14000 mL	8.402	-0.02008	0.98880	0.00100	
48:30.8	Data point 47						-0.01599		0.00096	
49:15.2	Data point 48						-0.01391			
49:52.1	Data point 49						-0.01192		0.00063	
50:28.1	Data point 50									
51:19.0	Data point 51								0.00097	
52:10.0	Data point 52									
53:01.4	Data point 52								0.00040	
53:52.4	Data point 54							0.52199		
53.52.4 54:43.5	Data point 55						-0.00005			
54:43.5 55:29.4	Data point 56							0.14476		
56:15.3	Data point 57									
56:51.0	Data point 58							0.83806		
57:26.7	Data point 59							0.93309		
58:02.7	Data point 60							0.92215	0.00040	
59:24.9	Data point 61						-0.00463		0.00036	
1:00:11.3	•						0.01276	0.89016	0.00067	
1:00:54.6							0.01126	0.59387	0.00072	
1:01:30.1							0.01630	0.71883		
1:02:16.1							-0.00161			
1:03:02.1							0.01037	0.37883		
1:03:37.6							0.00768	0.16641	0.00093	
1:04:20.5							0.00085	0.00280	0.00079	10.0 s
1:04:56.0	Data point 69	1.50000 mL	0.31343 mL	0.32025 mL	0.44000 mL	3.661	0.00260	0.02567		
	•						-0.00269			
1:06:06.9							-0.00164			
							-0.01091		0.00097	
	Data point 73						0.01056	0.34237	0.00089	
	Data point 74						-0.01091		0.00090	
	Data point 75						-0.01377		0.00095	
	Data point 76						-0.01377		0.00033	
	Data point 77						-0.01473		0.00100	
	Data point 78						-0.01933			Timed out at
1.12.00.0	Dάια μοπιτο	1.50000 1112	U.3 1343 IIIL	U.32333 IIIL	U.44000 IIIL	1.101	-0.00000	0.80302	0.00200	59.5 s
1:13:31.0	Data point 79	1.50000 mL	0.31343 mL	0.32340 mL	0.44000 mL	7.330	-0.03553	0.97995	0.00177	Timed out at
1.15.01 5	Data point 80	4 50000 ml	^ 212/12 ml	0 22254 ml	0 44000 ml	୦ ୨୦୫	-0.03259	0.06414	0.00164	59.5 s Timed out at
1.10.01.0	Data ропт оо	1.50000 m∟	U.3 I343 III∟	U.32304 IIIL	U.44000 III∟	ბ. ∠უ∪	-U.U3233	U.904 14	U.UU 1U -1	59.5 s
1:16:26.9	Data point 81	1.50000 mL	0.31343 mL	0.32366 mL	0.44000 mL	8.546	-0.01671	0.68295	0.00100	
	Data point 82						-0.01490		0.00094	
	Data point 83						0.00173	0.00806	0.00095	
	Data point 84						0.00175		0.00099	
	Data point 85								0.00098	
	Data point 86								0.00098	
	Data point 87								0.00094	
	Data point 88								0.00054	
1.41.10.4	Data μυπτ συ	1.50000 IIIL	0.31343 IIIL	0.5500 <i>i</i> iiiL	0.44000 IIIL	10.730	0.01100	0.70090	0.00000	10.0 5
1										,





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Events (continued)

Time	Event	Water	Acid	Base	Octanol	pН	dpH/dt	pH R-squared	pH SD	dpH/dt time
1:21:51.7	Data point 89	1.50000 mL	0.31343 mL	0.33321 mL	0.44000 mL	10.914	0.00908	0.75807	0.00051	10.0 s
1:22:27.2	Data point 90	1.50000 mL	0.31343 mL	0.33655 mL	0.44000 mL	11.100	0.00124	0.00385	0.00099	10.0 s
1:23:02.7	Data point 91	1.50000 mL	0.31343 mL	0.34132 mL	0.44000 mL	11.258	0.01768	0.78273	0.00099	16.0 s
1:23:60.0	Data point 92	1.50000 mL	0.31343 mL	0.34927 mL	0.44000 mL	11.450	0.01237	0.90868	0.00064	10.0 s
1:24:46.2	Data point 93	1.50000 mL	0.31343 mL	0.36051 mL	0.44000 mL	11.643	0.01058	0.60009	0.00067	10.0 s
1:25:21.9	Data point 94	1.50000 mL	0.31343 mL	0.37681 mL	0.44000 mL	11.828	0.01552	0.68214	0.00093	10.0 s
1:25:58.0	Data point 95	1.50000 mL	0.31343 mL	0.40214 mL	0.44000 mL	12.016	0.00685	0.34934	0.00057	10.0 s
1:26:17.0	Assay volumes	1.50000 mL	0.31343 mL	0.40214 mL	0.44000 mL					



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	12.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings				
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium	4 = 0 1			
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.040 mL			
Partition solvent added	Automatic			
After partition addition, stir for	1 seconds			
Sample Sonication				
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	300 seconds			
After sonication stir for	5 seconds			
Sample Dissolution	Vaa			
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution	To start pH			
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge	No			
Perform a carbonate purge Temperature Control	No			
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			
Stir speed of	50%			
Titration 1	30 70			
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	50%			
Titration 2	J J / U			
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume				
Additional partition solvent added	Automatic			
After pU adjust stir for	30 seconds			

Reported at: 3/2/2018 2:47:55 PM

30 seconds

15 seconds

55%

After pH adjust stir for

Stir to allow partitioning for

Stirrer speed for partitioning



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Assay Settings (continued)

Value	Original Value	Date/Time changed	Imported from
	•	_	-
Low to high pH			
0.00 mL			
0.300 mL			
Automatic			
30 seconds			
15 seconds			
60%			
No			
0 seconds			
20 points			
0.50 seconds			
0.00100 dpH/dt			
60 seconds			
	Low to high pH 0.00 mL 0.300 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt	Low to high pH 0.00 mL 0.300 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt	Low to high pH 0.00 mL 0.300 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.130	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
Four-Plus S	0.9970	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
Four-Plus jH	8.0	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
Four-Plus jOH	-0.4	3/1/2018 10:41:04 AM	C:\Sirius_T3\HCl18B27.t3r
Base concentration factor	1.000	3/1/2018 10:41:05 AM	C:\Sirius_T3\KOH18B27.t3r
Acid concentration factor	0.994	3/1/2018 10:41:05 AM	C:\Sirius_T3\HCl18B27.t3r

Instrument Settings

Setting Instrument owner Instrument ID	Value Merck T312060	Batch Id	Install date
Instrument type Software version Dispenser module Dispenser 0 Syringe volume	T3 Simulator 1.1.3.0 Water 2.5 mL	T3DM1200361	3/31/2009 5:24:52 AM 3/31/2009 5:25:05 AM
Firmware version Titrant Dispenser 2 Syringe volume	1.2.1(r2) Water (0.15 M KCI) Acid 0.5 mL	02-06-2018	2/27/2018 10:05:59 AM 3/31/2009 5:25:11 AM
Firmware version Titrant Dispenser 1 Syringe volume Firmware version	1.2.1(r2) Acid (0.5 M HCI) Base 0.5 mL 1.2.1(r2)	02-27-2018	2/27/2018 10:27:22 AM 3/31/2009 5:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2)	9/22/2017	2/27/2018 10:21:22 AM 3/31/2009 5:26:24 AM
Distribution valve 5 Firmware version Port A	Distribution Valve 1.1.3 Methanol (80%, 0.15 M KCI)	09-26-17	3/31/2009 5:28:19 AM 2/7/2018 9:42:01 AM
Port B Dispenser 3 Syringe volume Firmware version	Cyclohexane Buffer 0.5 mL 1.2.1(r2)	11-01-17	2/27/2018 10:37:57 AM 8/3/2010 5:05:16 AM
Titrant Dispenser 6	Dodecane Octanol	2018/01/31	2/28/2018 10:18:04 AM 10/22/2010 10:52:43 AM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting Suringa valuma	Value 0.5 mL	Batch Id	Install date
Syringe volume Firmware version	1.2.1(r2)		
Titrant	Octanol	01-31-2018	2/27/2018 9:59:35 AM
Titrator			3/31/2009 5:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1	T0=000	1/00/00/10 0 0 1 00 514
Electrode	T3 Electrode	T3E0923	1/23/2018 2:01:00 PM
E0 calibration	+4.68 mV	KCL097	3/1/2018 10:41:53 AM
Filling solution Liquids	3M KCI	KCL097	2/27/2018 9:49:43 AM
Wash 1	50% IPA:50% Water		2/28/2018 10:23:32 AM
Wash 2	0.5% Trition X-100 in H20		2/28/2018 10:23:34 AM
Buffer position 1	pH7 Wash		2/28/2018 10:24:06 AM
Buffer position 2	pH 7		2/28/2018 10:24:08 AM
Storage position	•		2/28/2018 10:21:14 AM
Wash water	8.2e+003 mL	02-27-2018	2/27/2018 9:54:39 AM
Waste	7.2e+003 mL		11/28/2017 10:36:29 AM
Temperature controller			8/5/2010 6:35:13 AM
Turbidity detector		074044	3/31/2009 5:24:45 AM
Spectrometer Dip probe		074811 10196	11/23/2010 11:22:28 AM
Dip probe Wavelength coefficient A0	183.333	10196	
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	112:08:55		11/23/2010 11:22:28 AM
Calibrated on	2/27/2018 10:40:38 AM		
Integration time	40		
Scans averaged	10		
Autoloader		T3AL1200345	11/10/2015 9:34:13 AM
Left-right axis firmware version	1.17 Al1DI2DO2 Stepper 2		
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version Chassis I/O firmware version	1.17 Al1Dl2DO2 Stepper 2 1.11 Al1Dl0DO4 Norgren I/O		
Configuration	1.11 Al IDIODO4 Noigien I/O		
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed Solvent wash stir duration	30% 5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed	30% 0%		
E0 calibration reading stir speed	U /0		



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting	Value	Batch Id	d Install date
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		

Refinement Settings

• 441		
Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00
• •		

Experiment Log

- [2:58] Air gap created for Water (0.15 M KCI)
- [2:58] Air gap created for Acid (0.5 M HCI)
- [2:58] Air gap created for Base (0.5 M KOH)
- [2:59] Air gap released for Water (0.15 M KCI)
- [3:03] Titrator arm moved over Titration position
- [3:03] Titration 1 of 3
- [3:03] Adding initial titrants
- [3:03] Automatically add 1.50000 mL of water
- [3:28] Dispensed 1.500000 mL of Water (0.15 M KCI)
- [3:32] Titrator arm moved over Drain
- [9:13] Titrator arm moved to Titration position
- [9:13] Argon flow rate set to 100
- [9:13] Stirrer speed set to 10
- [9:18] Automatically add 0.04000 mL of Octanol
- [9:19] Dispensed 0.040005 mL of Octanol
- [9:21] Initial pH = 3.61
- [9:21] Iterative adjust 3.61 -> 2.00
- [9:21] pH 3.61 -> 2.00
- [9:22] Air gap released for Acid (0.5 M HCI)
- [9:23] Dispensed 0.041627 mL of Acid (0.5 M HCI)
- [9:28] pH 2.04 -> 2.00
- [9:28] Dispensed 0.003551 mL of Acid (0.5 M HCl)
- [9:33] Holding pH 2.00
- [11:33] Stirrer speed set to 0
- [11:33] Stirrer speed set to 50
- [11:33] Iterative adjust 1.97 -> 2.00
- [11:33] pH 1.97 -> 2.00
- [11:34] Air gap released for Base (0.5 M KOH)
- [11:35] Dispensed 0.002752 mL of Base (0.5 M KOH)
- [12:25] Stirrer speed set to 0
- [12:36] Datapoint id 1 collected
- [12:36] Stirrer speed set to 50
- [12:41] pH 2.02 -> 2.22
- [12:41] Using cautious pH adjust
- [12:41] Dispensed 0.007667 mL of Base (0.5 M KOH)
- [12:46] Stepping pH = 2.10
- [12:46] Dispensed 0.005974 mL of Base (0.5 M KOH)
- [12:52] Stepping pH = 2.19
- [12:52] Dispensed 0.001599 mL of Base (0.5 M KOH)
- [12:57] Stepping pH = 2.22
- [13:12] Stirrer speed set to 0
- [13:22] Datapoint id 2 collected



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [13:22] Charge balance equation is out by 0.5%
- [13:22] Stirrer speed set to 50
- [13:27] pH 2.22 -> 2.42
- [13:27] Using charge balance adjust
- [13:27] Dispensed 0.009525 mL of Base (0.5 M KOH)
- [13:48] Stirrer speed set to 0
- [13:58] Datapoint id 3 collected
- [13:58] Charge balance equation is out by -8.2%
- [13:58] Stirrer speed set to 50
- [14:03] pH 2.41 -> 2.61
- [14:03] Using charge balance adjust
- [14:03] Dispensed 0.006162 mL of Base (0.5 M KOH)
- [14:23] Stirrer speed set to 0
- [14:33] Datapoint id 4 collected
- [14:33] Charge balance equation is out by -16.7%
- [14:33] Stirrer speed set to 50
- [14:38] pH 2.58 -> 2.78
- [14:38] Using cautious pH adjust
- [14:38] Dispensed 0.002093 mL of Base (0.5 M KOH)
- [14:43] Stepping pH = 2.67
- [14:44] Dispensed 0.001623 mL of Base (0.5 M KOH)
- [14:49] Stepping pH = 2.74
- [14:49] Dispensed 0.000659 mL of Base (0.5 M KOH)
- [14:54] Stepping pH = 2.78
 - 4.54] Stepping pri 2.76
- [15:09] Stirrer speed set to 0
- [15:19] Datapoint id 5 collected
- [15:19] Charge balance equation is out by -4.9%
- [15:19] Stirrer speed set to 50
- [15:24] pH 2.78 -> 2.98
- [15:24] Using charge balance adjust
- [15:24] Dispensed 0.002752 mL of Base (0.5 M KOH)
- [15:45] Stirrer speed set to 0
- [15:55] Datapoint id 6 collected
- [15:55] Charge balance equation is out by -31.1%
- [15:55] Stirrer speed set to 50
- [16:00] pH 2.92 -> 3.12
- [16:00] Using cautious pH adjust
- [16:00] Dispensed 0.001035 mL of Base (0.5 M KOH)
- [16:05] Stepping pH = 2.99
- [16:05] Dispensed 0.001129 mL of Base (0.5 M KOH)
- [16:10] Stepping pH = 3.07
- [16:10] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [16:15] Stepping pH = 3.11
- [16:15] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [16:20] Stepping pH = 3.12
- [16:36] Stirrer speed set to 0
- [16:46] Datapoint id 7 collected
- [16:46] Charge balance equation is out by -32.9%
- [16:46] Stirrer speed set to 50
- [16:51] pH 3.11 -> 3.31
- [16:51] Using cautious pH adjust
- [16:51] Dispensed 0.000729 mL of Base (0.5 M KOH)
- [16:56] Stepping pH = 3.18
- [16:56] Dispensed 0.000776 mL of Base (0.5 M KOH)
- [17:01] Stepping pH = 3.25
- [17:01] Dispensed 0.000494 mL of Base (0.5 M KOH)
- [17:06] Stepping pH = 3.30
- [17:06] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [17:11] Stepping pH = 3.31



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [17:26] Stirrer speed set to 0
- [17:37] Datapoint id 8 collected
- [17:37] Charge balance equation is out by -46.4%
- [17:37] Stirrer speed set to 50
- [17:42] pH 3.30 -> 3.50
- [17:42] Using cautious pH adjust
- [17:42] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [17:47] Stepping pH = 3.36
- [17:47] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [17:52] Stepping pH = 3.43
- [17:52] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [17:57] Stepping pH = 3.48
- [17:57] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [18:02] Stepping pH = 3.49
- [18:02] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [18:08] Stepping pH = 3.50
- [18:23] Stirrer speed set to 0
- [18:33] Datapoint id 9 collected
- [18:33] Charge balance equation is out by -90.2%
- [18:33] Stirrer speed set to 50
- [18:38] pH 3.49 -> 3.69
- [18:38] Using cautious pH adjust
- [18:38] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [18:43] Stepping pH = 3.55
- [18:43] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [18:48] Stepping pH = 3.63
- [18:48] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [18:53] Stepping pH = 3.67
- [18:53] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [18:58] Stepping pH = 3.68
- [18:59] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [19:04] Stepping pH = 3.69 [19:19] Stirrer speed set to 0
- [19:29] Datapoint id 10 collected
- [19:29] Charge balance equation is out by -90.4%
- [19:29] Stirrer speed set to 50
- [19:34] pH 3.68 -> 3.88
- [19:34] Using cautious pH adjust
- [19:34] Dispensed 0.000564 mL of Base (0.5 M KOH)
- [19:40] Stepping pH = 3.76
- [19:40] Dispensed 0.000564 mL of Base (0.5 M KOH)
- [19:45] Stepping pH = 3.85
- [19:45] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [19:50] Stepping pH = 3.87
- [19:50] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [19:55] Stepping pH = 3.89
- [20:10] Stirrer speed set to 0
- [20:20] Datapoint id 11 collected
- [20:20] Charge balance equation is out by -31.4%
- [20:20] Stirrer speed set to 50
- [20:25] pH 3.88 -> 4.08
- [20:25] Using cautious pH adjust
- [20:26] Dispensed 0.000659 mL of Base (0.5 M KOH)
- [20:31] Stepping pH = 4.01
- [20:31] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [20:36] Stepping pH = 4.07
- [20:51] Stirrer speed set to 0
- [21:01] Datapoint id 12 collected
- [21:01] Charge balance equation is out by 25.7%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [21:01] Stirrer speed set to 50
- [21:06] pH 4.07 -> 4.27
- [21:06] Using cautious pH adjust
- [21:06] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [21:11] Stepping pH = 4.29
- [21:26] Stirrer speed set to 0
- [21:37] Datapoint id 13 collected
- [21:37] Charge balance equation is out by 50.0%
- [21:37] Stirrer speed set to 50
- [21:42] pH 4.28 -> 4.48
- [21:42] Using cautious pH adjust
- [21:42] Dispensed 0.000847 mL of Base (0.5 M KOH)
- [21:47] Stepping pH = 4.69
- [22:02] Stirrer speed set to 0
- [22:14] Datapoint id 14 collected
- [22:14] Charge balance equation is out by 50.0%
- 20.441 Officer and and to 50
- [22:14] Stirrer speed set to 50
- [22:19] pH 4.68 -> 4.88
- [22:19] Using cautious pH adjust
- [22:19] Dispensed 0.000823 mL of Base (0.5 M KOH)
- [22:24] Stepping pH = 7.46
- [22:39] Stirrer speed set to 0
- [23:39] Datapoint id 15 collected
- [23:39] Charge balance equation is out by 50.0%
- [23:39] Stirrer speed set to 50
- [23:44] pH 7.86 -> 8.06
- [23:44] Using cautious pH adjust
- [23:44] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [23:49] Stepping pH = 8.28
- [24:04] Stirrer speed set to 0
- [24:43] Datapoint id 16 collected
- [24:43] Charge balance equation is out by 50.0%
- [24:43] Stirrer speed set to 50
- [24:48] pH 8.16 -> 8.36
- [24:48] Using cautious pH adjust
- [24:48] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [24:53] Stepping pH = 8.87
- [25:08] Stirrer speed set to 0
- [25:18] Datapoint id 17 collected
- [25:18] Charge balance equation is out by 50.0%
- [25:18] Stirrer speed set to 50
- [25:24] pH 8.85 -> 9.05
- [25:24] Using cautious pH adjust
- [25:24] Dispensed 0.000776 mL of Base (0.5 M KOH)
- [25:29] Stepping pH = 9.43
- [25:44] Stirrer speed set to 0
- [25:54] Datapoint id 18 collected
- [25:54] Charge balance equation is out by 50.0%
- [25:54] Stirrer speed set to 50
- [26:00] pH 9.41 -> 9.61
- [26:00] Using cautious pH adjust
- [26:00] Dispensed 0.000847 mL of Base (0.5 M KOH)
- [26:05] Stepping pH = 9.70
- [26:20] Stirrer speed set to 0
- [26:30] Datapoint id 19 collected
- [26:30] Charge balance equation is out by 50.0%
- [26:30] Stirrer speed set to 50
- [26:35] pH 9.69 -> 9.89
- [26:35] Using cautious pH adjust



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [26:35] Dispensed 0.000706 mL of Base (0.5 M KOH)
- [26:40] Stepping pH = 9.85
- [26:40] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [26:45] Stepping pH = 9.86
- [26:45] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [26:51] Stepping pH = 9.89
- [27:06] Stirrer speed set to 0
- [27:16] Datapoint id 20 collected
- [27:16] Charge balance equation is out by 23.1%
- [27:16] Stirrer speed set to 50
- [27:21] pH 9.87 -> 10.07
- [27:21] Using cautious pH adjust
- [27:21] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [27:26] Stepping pH = 10.00
- [27:26] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [27:31] Stepping pH = 10.04
- [27:31] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [27:36] Stepping pH = 10.06
- [27:36] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [27:42] Stepping pH = 10.07
- [27:57] Stirrer speed set to 0
- [28:07] Datapoint id 21 collected
- [28:07] Charge balance equation is out by 5.9%
- [28:07] Stirrer speed set to 50
- [28:12] pH 10.06 -> 10.26
- [28:12] Using charge balance adjust
- [28:12] Dispensed 0.001035 mL of Base (0.5 M KOH)
- [28:32] Stirrer speed set to 0
- [28:42] Datapoint id 22 collected
- [28:42] Charge balance equation is out by -19.8%
- [28:42] Stirrer speed set to 50
- [28:47] pH 10.23 -> 10.43
- [28:47] Using cautious pH adjust
- [28:47] Dispensed 0.000494 mL of Base (0.5 M KOH)
- [28:52] Stepping pH = 10.30
- [28:53] Dispensed 0.000517 mL of Base (0.5 M KOH)
- [28:58] Stepping pH = 10.37
- [28:58] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [29:03] Stepping pH = 10.41
- [29:03] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [29:08] Stepping pH = 10.42
- [29:23] Stirrer speed set to 0
- [29:34] Datapoint id 23 collected
- [29:34] Charge balance equation is out by -43.7%
- [29:34] Stirrer speed set to 50
- [29:39] pH 10.41 -> 10.61
- [29:39] Using cautious pH adjust
- [29:39] Dispensed 0.000517 mL of Base (0.5 M KOH)
- [29:44] Stepping pH = 10.48
- [29:44] Dispensed 0.000611 mL of Base (0.5 M KOH)
- [29:49] Stepping pH = 10.55
- [29:49] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [29:54] Stepping pH = 10.59
- [29:54] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [29:59] Stepping pH = 10.60
- [30:15] Stirrer speed set to 0
- [30:25] Datapoint id 24 collected
- [30:25] Charge balance equation is out by -63.3%
- [30:25] Stirrer speed set to 50



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [30:30] pH 10.60 -> 10.80
- [30:30] Using cautious pH adjust
- [30:30] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [30:35] Stepping pH = 10.66
- [30:35] Dispensed 0.000776 mL of Base (0.5 M KOH)
- [30:40] Stepping pH = 10.73
- [30:40] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [30:45] Stepping pH = 10.78
- [30:45] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [30:51] Stepping pH = 10.79
- [31:06] Stirrer speed set to 0
- [31:16] Datapoint id 25 collected
- [31:16] Charge balance equation is out by -66.0%
- [31:16] Stirrer speed set to 50
- [31:21] pH 10.78 -> 10.98
- [31:21] Using cautious pH adjust
- [31:21] Dispensed 0.000870 mL of Base (0.5 M KOH)
- [31:26] Stepping pH = 10.84
- [31:27] Dispensed 0.001082 mL of Base (0.5 M KOH)
- [31:32] Stepping pH = 10.91
- [31:32] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [31:37] Stepping pH = 10.96
- [31:37] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [31:42] Stepping pH = 10.97
- [31:57] Stirrer speed set to 0
- [32:07] Datapoint id 26 collected
- [32:07] Charge balance equation is out by -71.2%
- [32:07] Stirrer speed set to 50
- [32:12] pH 10.96 -> 11.16
- [32:12] Using cautious pH adjust
- [32:12] Dispensed 0.001246 mL of Base (0.5 M KOH)
- [32:17] Stepping pH = 11.04
- [32:17] Dispensed 0.001270 mL of Base (0.5 M KOH)
- [32:23] Stepping pH = 11.12
- [32:23] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [32:28] Stepping pH = 11.15
- [32:43] Stirrer speed set to 0
- [32:53] Datapoint id 27 collected
- [32:53] Charge balance equation is out by -24.0%
- [32:53] Stirrer speed set to 50
- [32:58] pH 11.15 -> 11.35
- [32:58] Using cautious pH adjust
- [32:58] Dispensed 0.001881 mL of Base (0.5 M KOH)
- [33:03] Stepping pH = 11.25
- [33:03] Dispensed 0.001505 mL of Base (0.5 M KOH)
- [33:09] Stepping pH = 11.31
- [33:09] Dispensed 0.000729 mL of Base (0.5 M KOH)
- [33:14] Stepping pH = 11.34
- [33:14] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [33:19] Stepping pH = 11.34
- [33:34] Stirrer speed set to 0
- [33:44] Datapoint id 28 collected
- [33:44] Charge balance equation is out by -18.0%
- [33:44] Stirrer speed set to 50
- [33:49] pH 11.34 -> 11.54
- [33:49] Using cautious pH adjust
- [33:49] Dispensed 0.002916 mL of Base (0.5 M KOH)
- [33:54] Stepping pH = 11.44
- [33:55] Dispensed 0.002328 mL of Base (0.5 M KOH)



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [34:00] Stepping pH = 11.52
- [34:00] Dispensed 0.000564 mL of Base (0.5 M KOH)
- [34:05] Stepping pH = 11.54
- [34:20] Stirrer speed set to 0
- [34:30] Datapoint id 29 collected
- [34:30] Charge balance equation is out by 0.7%
- [34:30] Stirrer speed set to 50
- [34:35] pH 11.54 -> 11.74
- [34:35] Using charge balance adjust
- [34:36] Dispensed 0.009219 mL of Base (0.5 M KOH)
- [34:56] Stirrer speed set to 0
- [35:06] Datapoint id 30 collected
- [35:06] Charge balance equation is out by -12.5%
- [35:06] Stirrer speed set to 50
- [35:11] pH 11.72 -> 11.92
- [35:11] Using charge balance adjust [35:11] Dispensed 0.014252 mL of Base (0.5 M KOH)
- [35:31] Stirrer speed set to 0
- [35:42] Datapoint id 31 collected
- [35:42] Charge balance equation is out by -7.9%
- [35:42] Stirrer speed set to 50
- [35:47] pH 11.91 -> 12.05
- [35:47] Using charge balance adjust
- [35:47] Dispensed 0.014722 mL of Base (0.5 M KOH)
- [36:07] Stirrer speed set to 0
- [36:17] Datapoint id 32 collected
- [36:17] Charge balance equation is out by -36.5% [36:17] Titration 2 of 3
- [36:17] Adding initial titrants [36:17] Automatically add 0.10000 mL of Octanol
- [36:20] Dispensed 0.100000 mL of Octanol
- [36:20] Stirrer speed set to 10
- [36:21] Stirrer speed set to 55
- [36:21] Iterative adjust 12.04 -> 2.00
- [36:21] pH 12.04 -> 2.00
- [36:23] Dispensed 0.100000 mL of Acid (0.5 M HCI)
- [36:28] pH 2.29 -> 2.00
- [36:29] Dispensed 0.024106 mL of Acid (0.5 M HCI)
- [36:34] pH 2.03 -> 2.00
- [36:34] Dispensed 0.003104 mL of Acid (0.5 M HCI)
- [37:25] Stirrer speed set to 0
- [37:35] Datapoint id 33 collected
- [37:35] Stirrer speed set to 55
- [37:40] pH 1.97 -> 2.17
- [37:40] Using cautious pH adjust
- [37:40] Dispensed 0.009995 mL of Base (0.5 M KOH)
- [37:45] Stepping pH = 2.06
- [37:46] Dispensed 0.007291 mL of Base (0.5 M KOH)
- [37:51] Stepping pH = 2.13
- [37:51] Dispensed 0.002422 mL of Base (0.5 M KOH)
- [37:56] Stepping pH = 2.17
- [38:11] Stirrer speed set to 0
- [38:21] Datapoint id 34 collected
- [38:21] Charge balance equation is out by 1.4%
- [38:21] Stirrer speed set to 55
- [38:26] pH 2.17 -> 2.37
- [38:26] Using charge balance adjust
- [38:26] Dispensed 0.012347 mL of Base (0.5 M KOH)
- [38:47] Stirrer speed set to 0



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [39:01] Datapoint id 35 collected
- [39:01] Charge balance equation is out by -4.4%
- [39:01] Stirrer speed set to 55
- [39:06] pH 2.37 -> 2.57
- [39:06] Using charge balance adjust
- [39:07] Dispensed 0.007855 mL of Base (0.5 M KOH)
- [39:27] Stirrer speed set to 0
- [39:37] Datapoint id 36 collected
- [39:37] Charge balance equation is out by -10.3%
- [39:37] Stirrer speed set to 55
- [39:42] pH 2.56 -> 2.76
- [39:42] Using charge balance adjust
- [39:42] Dispensed 0.005198 mL of Base (0.5 M KOH)
- [40:02] Stirrer speed set to 0
- [40:12] Datapoint id 37 collected
- [40:12] Charge balance equation is out by -18.8%
- [40:12] Stirrer speed set to 55
- [40:18] pH 2.73 -> 2.93
- [40:18] Using cautious pH adjust
- [40:18] Dispensed 0.001811 mL of Base (0.5 M KOH)
- [40:23] Stepping pH = 2.81
- [40:23] Dispensed 0.001552 mL of Base (0.5 M KOH)
- [40:28] Stepping pH = 2.88
- [40:28] Dispensed 0.000729 mL of Base (0.5 M KOH)
- [40:33] Stepping pH = 2.91
- [40:33] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [40:38] Stepping pH = 2.93
- [40:54] Stirrer speed set to 0
- [41:04] Datapoint id 38 collected
- [41:04] Charge balance equation is out by -19.6%
- [41:04] Stirrer speed set to 55
- [41:09] pH 2.93 -> 3.13
- [41:09] Using cautious pH adjust
- [41:09] Dispensed 0.001223 mL of Base (0.5 M KOH)
- [41:14] Stepping pH = 2.98
- [41:14] Dispensed 0.001576 mL of Base (0.5 M KOH)
- [41:19] Stepping pH = 3.07
- [41:19] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [41:24] Stepping pH = 3.12
- [41:40] Stirrer speed set to 0
- [41:50] Datapoint id 39 collected
- [41:50] Charge balance equation is out by -46.4%
- [41:50] Stirrer speed set to 55
- [41:55] pH 3.12 -> 3.32
- [41:55] Using cautious pH adjust
- [41:55] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [42:00] Stepping pH = 3.18
- [42:00] Dispensed 0.001105 mL of Base (0.5 M KOH)
- [42:05] Stepping pH = 3.28
- [42:05] Dispensed 0.000400 mL of Base (0.5 M KOH)
- [42:10] Stepping pH = 3.32
- [42:25] Stirrer speed set to 0
- [42:36] Datapoint id 40 collected
- [42:36] Charge balance equation is out by -32.0%
- [42:36] Stirrer speed set to 55
- [42:41] pH 3.32 -> 3.52
- [42:41] Using cautious pH adjust
- [42:41] Dispensed 0.000776 mL of Base (0.5 M KOH)
- [42:46] Stepping pH = 3.40



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [42:46] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [42:51] Stepping pH = 3.48
- [42:51] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [42:57] Stepping pH = 3.51
- [43:12] Stirrer speed set to 0
- [43:22] Datapoint id 41 collected
- [43:22] Charge balance equation is out by -16.3%
- [43:22] Stirrer speed set to 55
- [43:27] pH 3.52 -> 3.72
- [43:27] Using cautious pH adjust
- [43:27] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [43:32] Stepping pH = 3.61
- [43:32] Dispensed 0.000564 mL of Base (0.5 M KOH)
- [43:37] Stepping pH = 3.69
- [43:37] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [43:42] Stepping pH = 3.71
- [43:57] Stirrer speed set to 0
- [44:14] Datapoint id 42 collected
- [44:14] Charge balance equation is out by 4.1%
- [44:14] Stirrer speed set to 55
- [44:19] pH 3.71 -> 3.91
- [44:19] Using charge balance adjust
- [44:19] Dispensed 0.001599 mL of Base (0.5 M KOH)
- [44:39] Stirrer speed set to 0
- [44:49] Datapoint id 43 collected
- [44:49] Charge balance equation is out by 96.0%
- [44:49] Stirrer speed set to 55
- [44:55] pH 4.11 -> 4.31
- [44:55] Using cautious pH adjust
- [44:55] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [45:00] Stepping pH = 4.50
- [45:15] Stirrer speed set to 0
- [45:25] Datapoint id 44 collected
- [45:25] Charge balance equation is out by 50.0%
- [45:25] Stirrer speed set to 55
- [45:31] pH 4.51 -> 4.71
- [45:31] Using cautious pH adjust
- [45:31] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [45:36] Stepping pH = 6.63
- [45:51] Stirrer speed set to 0
- [46:51] Datapoint id 45 collected
- [46:51] Charge balance equation is out by 50.0%
- [46:51] Stirrer speed set to 55
- [46:56] pH 6.57 -> 6.77
- [46:56] Using cautious pH adjust
- [46:56] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [47:01] Stepping pH = 6.51
- [47:01] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [47:06] Stepping pH = 8.44
- [47:21] Stirrer speed set to 0
- [48:05] Datapoint id 46 collected
- [48:05] Charge balance equation is out by -289.9%
- [48:05] Stirrer speed set to 55
- [48:10] pH 8.37 -> 8.57
- [48:10] Using cautious pH adjust
- [48:10] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [48:15] Stepping pH = 9.01
- [48:30] Stirrer speed set to 0
- [48:49] Datapoint id 47 collected



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [48:49] Charge balance equation is out by 50.0%
- [48:49] Stirrer speed set to 55
- [48:55] pH 8.98 -> 9.18
- [48:55] Using cautious pH adjust
- [48:55] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [49:00] Stepping pH = 9.60
- [49:15] Stirrer speed set to 0
- [49:26] Datapoint id 48 collected
- [49:26] Charge balance equation is out by 50.0%
- [49:26] Stirrer speed set to 55
- [49:31] pH 9.59 -> 9.79
- [49:31] Using cautious pH adjust
- [49:32] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [49:37] Stepping pH = 9.96
- [49:52] Stirrer speed set to 0
- [50:02] Datapoint id 49 collected
- [50:02] Charge balance equation is out by 50.0%
- [50:02] Stirrer speed set to 55
- [50:07] pH 9.95 -> 10.15
- [50:07] Using cautious pH adjust
- [50:07] Dispensed 0.000776 mL of Base (0.5 M KOH)
- [50:13] Stepping pH = 10.16
- [50:28] Stirrer speed set to 0
- [50:38] Datapoint id 50 collected
- [50:38] Charge balance equation is out by 50.0%
- [50:38] Stirrer speed set to 55
- [50:43] pH 10.16 -> 10.36
- [50:43] Using cautious pH adjust
- [50:43] Dispensed 0.000706 mL of Base (0.5 M KOH)
- [50:48] Stepping pH = 10.30
- [50:48] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [50:53] Stepping pH = 10.34
- [50:53] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [50:58] Stepping pH = 10.35
- [50:58] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [51:03] Stepping pH = 10.36
- [51:19] Stirrer speed set to 0
 - (1:10) Determine in E4 collector
- [51:29] Datapoint id 51 collected
- [51:29] Charge balance equation is out by 19.2%
- [51:29] Stirrer speed set to 55
- [51:34] pH 10.35 -> 10.55
- [51:34] Using cautious pH adjust
- [51:34] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [51:39] Stepping pH = 10.45
- [51:39] Dispensed 0.000517 mL of Base (0.5 M KOH)
- [51:44] Stepping pH = 10.51
- [51:44] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [51:49] Stepping pH = 10.54
- [51:49] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [51:55] Stepping pH = 10.55
- [52:10] Stirrer speed set to 0
- [52:20] Datapoint id 52 collected
- [52:20] Charge balance equation is out by -16.9%
- [52:20] Stirrer speed set to 55
- [52:25] pH 10.54 -> 10.74
- [52:25] Using cautious pH adjust
- [52:25] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [52:30] Stepping pH = 10.61
- [52:30] Dispensed 0.000823 mL of Base (0.5 M KOH)



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [52:36] Stepping pH = 10.68
- [52:36] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [52:41] Stepping pH = 10.72
- [52:41] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [52:46] Stepping pH = 10.73
- [53:01] Stirrer speed set to 0
- [53:11] Datapoint id 53 collected
- [53:11] Charge balance equation is out by -57.5%
- [53:11] Stirrer speed set to 55
- [53:16] pH 10.73 -> 10.93
- [53:16] Using cautious pH adjust
- [53:16] Dispensed 0.000964 mL of Base (0.5 M KOH)
- [53:21] Stepping pH = 10.79
- [53:21] Dispensed 0.001176 mL of Base (0.5 M KOH)
- [53:27] Stepping pH = 10.87
- [53:27] Dispensed 0.000659 mL of Base (0.5 M KOH)
- [53:32] Stepping pH = 10.92
- [53:32] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [53:37] Stepping pH = 10.92
- [53:52] Stirrer speed set to 0
- [54:02] Datapoint id 54 collected
- [54:02] Charge balance equation is out by -53.2%
- [54:02] Stirrer speed set to 55
- [54:07] pH 10.92 -> 11.12
- [54:07] Using cautious pH adjust
- [54:07] Dispensed 0.001364 mL of Base (0.5 M KOH)
- [54:12] Stepping pH = 11.00
- [54:13] Dispensed 0.001317 mL of Base (0.5 M KOH)
- [54:18] Stepping pH = 11.07
- [54:18] Dispensed 0.000706 mL of Base (0.5 M KOH)
- [54:23] Stepping pH = 11.11
- [54:23] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [54:28] Stepping pH = 11.12
- [54:43] Stirrer speed set to 0
- [54:53] Datapoint id 55 collected
- [54:53] Charge balance equation is out by -31.5%
- [54:53] Stirrer speed set to 55
- [54:58] pH 11.11 -> 11.31
- [54:58] Using cautious pH adjust
- [54:58] Dispensed 0.002023 mL of Base (0.5 M KOH)
- [55:03] Stepping pH = 11.19
- [55:04] Dispensed 0.002093 mL of Base (0.5 M KOH)
- [55:09] Stepping pH = 11.27
- [55:09] Dispensed 0.001058 mL of Base (0.5 M KOH)
- [55:14] Stepping pH = 11.30 [55:29] Stirrer speed set to 0
- [55:39] Datapoint id 56 collected
- [55:39] Charge balance equation is out by -27.4%
- [55:39] Stirrer speed set to 55
- [55:44] pH 11.30 -> 11.50
- [55:44] Using cautious pH adjust
- [55:44] Dispensed 0.003081 mL of Base (0.5 M KOH)
- [55:49] Stepping pH = 11.40
- [55:50] Dispensed 0.002446 mL of Base (0.5 M KOH)
- [55:55] Stepping pH = 11.46
- [55:55] Dispensed 0.001129 mL of Base (0.5 M KOH)
- [56:00] Stepping pH = 11.49
- [56:15] Stirrer speed set to 0
- [56:25] Datapoint id 57 collected



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [56:25] Charge balance equation is out by -7.6%
- [56:25] Stirrer speed set to 55
- [56:30] pH 11.49 -> 11.69
- [56:30] Using charge balance adjust
- [56:30] Dispensed 0.009572 mL of Base (0.5 M KOH)
- [56:51] Stirrer speed set to 0
- [57:01] Datapoint id 58 collected
- [57:01] Charge balance equation is out by -13.0%
- [57:01] Stirrer speed set to 55
- [57:06] pH 11.67 -> 11.87
- [57:06] Using charge balance adjust
- [57:06] Dispensed 0.014652 mL of Base (0.5 M KOH)
- [57:26] Stirrer speed set to 0
- [57:36] Datapoint id 59 collected
- [57:36] Charge balance equation is out by -6.6%
- [57:36] Stirrer speed set to 55
- [57:41] pH 11.86 -> 12.05
- [57:41] Using charge balance adjust
- [57:42] Dispensed 0.022013 mL of Base (0.5 M KOH)
- [58:02] Stirrer speed set to 0
- [58:12] Datapoint id 60 collected
- [58:12] Charge balance equation is out by -13.7%
- [58:12] Titration 3 of 3
- [58:12] Adding initial titrants
- [58:12] Automatically add 0.30000 mL of Octanol
- [58:19] Dispensed 0.300000 mL of Octanol
- [58:19] Stirrer speed set to 10
- [58:20] Stirrer speed set to 60
- [58:20] Iterative adjust 12.04 -> 2.00
- [58:20] pH 12.04 -> 2.00
- [58:23] Dispensed 0.100000 mL of Acid (0.5 M HCI)
- [58:28] pH 2.44 -> 2.00
- [58:29] Dispensed 0.036007 mL of Acid (0.5 M HCI)
- [58:34] pH 2.04 -> 2.00
- [58:34] Dispensed 0.005033 mL of Acid (0.5 M HCI)
- [59:25] Stirrer speed set to 0
- [59:35] Datapoint id 61 collected
- [59:35] Stirrer speed set to 60
- [59:40] pH 1.97 -> 2.17
- [59:40] Using cautious pH adjust
- [59:40] Dispensed 0.011359 mL of Base (0.5 M KOH)
- [59:45] Stepping pH = 2.05
- [59:45] Dispensed 0.009595 mL of Base (0.5 M KOH)
- [59:51] Stepping pH = 2.15
- [59:51] Dispensed 0.002140 mL of Base (0.5 M KOH)
- [59:56] Stepping pH = 2.17
- [1:00:11] Stirrer speed set to 0
- [1:00:28] Datapoint id 62 collected
- [1:00:28] Charge balance equation is out by -1.7%
- [1:00:28] Stirrer speed set to 60
- [1:00:34] pH 2.18 -> 2.38
- [1:00:34] Using charge balance adjust
- [1:00:34] Dispensed 0.013993 mL of Base (0.5 M KOH)
- [1:00:54] Stirrer speed set to 0
- [1:01:04] Datapoint id 63 collected
- [1:01:04] Charge balance equation is out by -6.3%
- [1:01:04] Stirrer speed set to 60
- [1:01:09] pH 2.38 -> 2.58
- [1:01:09] Using charge balance adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [1:01:10] Dispensed 0.009008 mL of Base (0.5 M KOH)
- [1:01:30] Stirrer speed set to 0
- [1:01:40] Datapoint id 64 collected
- [1:01:40] Charge balance equation is out by -15.0%
- [1:01:40] Stirrer speed set to 60
- [1:01:45] pH 2.56 -> 2.76
- [1:01:45] Using cautious pH adjust
- [1:01:45] Dispensed 0.003057 mL of Base (0.5 M KOH)
- [1:01:50] Stepping pH = 2.62
- [1:01:50] Dispensed 0.003198 mL of Base (0.5 M KOH)
- [1:01:55] Stepping pH = 2.72
- [1:01:56] Dispensed 0.001035 mL of Base (0.5 M KOH)
- 1:02:01 Stepping pH = 2.75
- [1:02:16] Stirrer speed set to 0
- 1:02:76] Detencint id 65 collec
- [1:02:26] Datapoint id 65 collected
- [1:02:26] Charge balance equation is out by -19.0%
- [1:02:26] Stirrer speed set to 60
- [1:02:31] pH 2.76 -> 2.96
- [1:02:31] Using cautious pH adjust
- [1:02:31] Dispensed 0.002046 mL of Base (0.5 M KOH)
- [1:02:36] Stepping pH = 2.84
- [1:02:36] Dispensed 0.001834 mL of Base (0.5 M KOH)
- [1:02:41] Stepping pH = 2.92
- [1:02:41] Dispensed 0.000776 mL of Base (0.5 M KOH)
- [1:02:47] Stepping pH = 2.95
- [1:03:02] Stirrer speed set to 0
- [1:03:12] Datapoint id 66 collected
- [1:03:12] Charge balance equation is out by -13.9%
- [1:03:12] Stirrer speed set to 60
- [1:03:17] pH 2.96 -> 3.16
- [1:03:17] Using charge balance adjust
- [1:03:17] Dispensed 0.002940 mL of Base (0.5 M KOH)
- [1:03:37] Stirrer speed set to 0
- [1:03:55] Datapoint id 67 collected
- [1:03:55] Charge balance equation is out by -11.3%
- [1:03:55] Stirrer speed set to 60
- [1:04:00] pH 3.14 -> 3.34
- [1:04:00] Using charge balance adjust
- [1:04:00] Dispensed 0.002399 mL of Base (0.5 M KOH)
- [1:04:20] Stirrer speed set to 0
- [1:04:30] Datapoint id 68 collected
- [1:04:30] Charge balance equation is out by 13.3%
- [1:04:30] Stirrer speed set to 60
- [1:04:35] pH 3.37 -> 3.57
- [1:04:35] Using charge balance adjust
- [1:04:35] Dispensed 0.002117 mL of Base (0.5 M KOH)
- [1:04:56] Stirrer speed set to 0
- [1:05:06] Datapoint id 69 collected
- [1:05:06] Charge balance equation is out by 44.0%
- [1:05:06] Stirrer speed set to 60
- [1:05:11] pH 3.67 -> 3.87
- [1:05:11] Using cautious pH adjust
- [1:05:11] Dispensed 0.001035 mL of Base (0.5 M KOH)
- [1:05:16] Stepping pH = 3.90
- [1:05:31] Stirrer speed set to 0
- [1:05:41] Datapoint id 70 collected
- [1:05:41] Charge balance equation is out by 50.0%
- [1:05:41] Stirrer speed set to 60
- [1:05:46] pH 3.90 -> 4.10



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [1:05:46] Using cautious pH adjust
- [1:05:46] Dispensed 0.000988 mL of Base (0.5 M KOH)
- [1:05:51] Stepping pH = 4.36
- [1:06:06] Stirrer speed set to 0
- [1:06:17] Datapoint id 71 collected
- [1:06:17] Charge balance equation is out by 50.0%
- [1:06:17] Stirrer speed set to 60
- [1:06:22] pH 4.36 -> 4.56
- [1:06:22] Using cautious pH adjust
- [1:06:22] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [1:06:27] Stepping pH = 5.66
- [1:06:42] Stirrer speed set to 0
- [1:07:02] Datapoint id 72 collected
- [1:07:02] Charge balance equation is out by 50.0%
- [1:07:02] Stirrer speed set to 60
- [1:07:08] pH 5.46 -> 5.66
- [1:07:08] Using cautious pH adjust
- [1:07:08] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:07:13] Stepping pH = 5.60
- [1:07:13] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:07:18] Stepping pH = 5.79
- [1:07:33] Stirrer speed set to 0
- [1:07:53] Datapoint id 73 collected
- [1:07:53] Charge balance equation is out by 32.1%
- [1:07:53] Stirrer speed set to 60
- [1:07:58] pH 5.73 -> 5.93
- [1:07:58] Using cautious pH adjust
- [1:07:58] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:08:03] Stepping pH = 6.04
- [1:08:18] Stirrer speed set to 0
- [1:08:36] Datapoint id 74 collected
- [1:08:36] Charge balance equation is out by 50.0%
- [1:08:36] Stirrer speed set to 60
- [1:08:41] pH 5.99 -> 6.19
- [1:08:41] Using cautious pH adjust
- [1:08:41] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:08:46] Stepping pH = 6.27
- [1:09:01] Stirrer speed set to 0
- [1:09:29] Datapoint id 75 collected
- [1:09:29] Charge balance equation is out by 50.0%
- [1:09:29] Stirrer speed set to 60
- [1:09:34] pH 6.09 -> 6.29
- [1:09:34] Using cautious pH adjust
- [1:09:35] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:09:40] Stepping pH = 6.48
- [1:09:55] Stirrer speed set to 0
- [1:10:29] Datapoint id 76 collected
- [1:10:29] Charge balance equation is out by 50.0%
- [1:10:29] Stirrer speed set to 60
- [1:10:34] pH 6.42 -> 6.62
- [1:10:34] Using cautious pH adjust
- [1:10:34] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:10:39] Stepping pH = 6.76
- [1:10:54] Stirrer speed set to 0
- [1:11:35] Datapoint id 77 collected
- [1:11:35] Charge balance equation is out by 50.0%
- [1:11:35] Stirrer speed set to 60
- [1:11:40] pH 6.79 -> 6.99
- [1:11:40] Using cautious pH adjust



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [1:11:40] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:11:45] Stepping pH = 6.86
- [1:11:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:11:50] Stepping pH = 7.15
- [1:12:05] Stirrer speed set to 0
- [1:13:05] Datapoint id 78 collected
- [1:13:05] Charge balance equation is out by -4.3%
- [1:13:05] Stirrer speed set to 60
- [1:13:10] pH 6.87 -> 7.07
- [1:13:10] Using charge balance adjust
- [1:13:10] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:13:31] Stirrer speed set to 0
- [1:14:31] Datapoint id 79 collected
- [1:14:31] Charge balance equation is out by 129.5%
- [1:14:31] Stirrer speed set to 60
- [1:14:36] pH 7.26 -> 7.46
- [1:14:36] Using cautious pH adjust
- [1:14:36] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:14:41] Stepping pH = 7.23
- [1:14:41] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [1:14:46] Stepping pH = 8.40
- [1:15:01] Stirrer speed set to 0
- [1:16:01] Datapoint id 80 collected
- [1:16:01] Charge balance equation is out by -238.7%
- [1:16:01] Stirrer speed set to 60
- [1:16:06] pH 8.30 -> 8.50
- [1:16:06] Using cautious pH adjust
- [1:16:06] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [1:16:11] Stepping pH = 8.60
- [1:16:27] Stirrer speed set to 0
- [1:16:46] Datapoint id 81 collected
- [1:16:46] Charge balance equation is out by 50.0%
- [1:16:46] Stirrer speed set to 60
- [1:16:51] pH 8.54 -> 8.74
- [1:16:51] Using cautious pH adjust
- [1:16:51] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:16:56] Stepping pH = 8.96
- [1:17:11] Stirrer speed set to 0
- [1:17:25] Datapoint id 82 collected
- [1:17:25] Charge balance equation is out by 50.0%
- [1:17:25] Stirrer speed set to 60
- [1:17:30] pH 8.94 -> 9.14
- [1:17:30] Using cautious pH adjust
- [1:17:30] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [1:17:35] Stepping pH = 9.41
- [1:17:50] Stirrer speed set to 0
- [1:18:02] Datapoint id 83 collected
- [1:18:02] Charge balance equation is out by 50.0%
- [1:18:02] Stirrer speed set to 60
- [1:18:07] pH 9.39 -> 9.59
- [1:18:07] Using cautious pH adjust
- [1:18:08] Dispensed 0.000776 mL of Base (0.5 M KOH)
- [1:18:13] Stepping pH = 9.83
- [1:18:28] Stirrer speed set to 0
- [1:18:38] Datapoint id 84 collected
- [1:18:38] Charge balance equation is out by 50.0%
- [1:18:38] Stirrer speed set to 60
- [1:18:43] pH 9.82 -> 10.02
- [1:18:43] Using cautious pH adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [1:18:43] Dispensed 0.000988 mL of Base (0.5 M KOH)
- [1:18:48] Stepping pH = 10.15
- [1:19:04] Stirrer speed set to 0
- [1:19:14] Datapoint id 85 collected
- [1:19:14] Charge balance equation is out by 50.0%
- [1:19:14] Stirrer speed set to 60
- [1:19:19] pH 10.14 -> 10.34
- [1:19:19] Using cautious pH adjust
- [1:19:19] Dispensed 0.000988 mL of Base (0.5 M KOH)
- [1:19:24] Stepping pH = 10.37
- [1:19:39] Stirrer speed set to 0
- [1:19:49] Datapoint id 86 collected
- [1:19:49] Charge balance equation is out by 50.1%
- [1:19:49] Stirrer speed set to 60
- [1:19:54] pH 10.36 -> 10.56
- [1:19:54] Using cautious pH adjust
- [1:19:54] Dispensed 0.000964 mL of Base (0.5 M KOH)
- [1:19:59] Stepping pH = 10.51
- [1:19:59] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [1:20:05] Stepping pH = 10.54
- [1:20:05] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:20:10] Stepping pH = 10.55
- [1:20:25] Stirrer speed set to 0
- [1:20:35] Datapoint id 87 collected
- [1:20:35] Charge balance equation is out by 27.0%
- [1:20:35] Stirrer speed set to 60
- [1:20:40] pH 10.55 -> 10.75
- [1:20:40] Using cautious pH adjust
- [1:20:40] Dispensed 0.001058 mL of Base (0.5 M KOH)
- [1:20:45] Stepping pH = 10.66
- [1:20:45] Dispensed 0.000659 mL of Base (0.5 M KOH)
- [1:20:50] Stepping pH = 10.71
- [1:20:50] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [1:20:56] Stepping pH = 10.73
- [1:20:56] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:21:01] Stepping pH = 10.74
- [1:21:16] Stirrer speed set to 0
- [1:21:26] Datapoint id 88 collected
- [1:21:26] Charge balance equation is out by -5.7%
- [1:21:26] Stirrer speed set to 60
- [1:21:31] pH 10.74 -> 10.94
- [1:21:31] Using charge balance adjust
- [1:21:31] Dispensed 0.002540 mL of Base (0.5 M KOH)
- [1:21:51] Stirrer speed set to 0
- [1:22:01] Datapoint id 89 collected
- [1:22:01] Charge balance equation is out by -13.9%
- [1:22:01] Stirrer speed set to 60
- [1:22:06] pH 10.92 -> 11.12
- [1:22:06] Using charge balance adjust
- [1:22:07] Dispensed 0.003340 mL of Base (0.5 M KOH)
- [1:22:27] Stirrer speed set to 0
- [1:22:37] Datapoint id 90 collected
- [1:22:37] Charge balance equation is out by -9.0%
- [1:22:37] Stirrer speed set to 60
- [1:22:42] pH 11.10 -> 11.30
- [1:22:42] Using charge balance adjust
- [1:22:42] Dispensed 0.004774 mL of Base (0.5 M KOH)
- [1:23:02] Stirrer speed set to 0
- [1:23:18] Datapoint id 91 collected



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01008 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01008_M16_octanol_pH-metric high logP.t3r

- [1:23:18] Charge balance equation is out by -23.1%
- [1:23:18] Stirrer speed set to 60
- [1:23:23] pH 11.26 -> 11.46
- [1:23:23] Using cautious pH adjust
- [1:23:24] Dispensed 0.003340 mL of Base (0.5 M KOH)
- [1:23:29] Stepping pH = 11.35
- [1:23:29] Dispensed 0.002893 mL of Base (0.5 M KOH)
- [1:23:34] Stepping pH = 11.43
- [1:23:34] Dispensed 0.001270 mL of Base (0.5 M KOH)
- [1:23:39] Stepping pH = 11.45
- [1:23:39] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [1:23:45] Stepping pH = 11.46
- [1:24:00] Stirrer speed set to 0
- [1:24:10] Datapoint id 92 collected
- [1:24:10] Charge balance equation is out by -19.0%
- [1:24:10] Stirrer speed set to 60
- [1:24:15] pH 11.46 -> 11.66
- [1:24:15] Using cautious pH adjust
- [1:24:15] Dispensed 0.005174 mL of Base (0.5 M KOH)
- [1:24:20] Stepping pH = 11.56
- [1:24:20] Dispensed 0.004092 mL of Base (0.5 M KOH)
- [1:24:25] Stepping pH = 11.62
- [1:24:26] Dispensed 0.001976 mL of Base (0.5 M KOH)
- [1:24:31] Stepping pH = 11.65
- [1:24:46] Stirrer speed set to 0
- [1:24:56] Datapoint id 93 collected
- [1:24:56] Charge balance equation is out by -7.9%
- [1:24:56] Stirrer speed set to 60
- [1:25:01] pH 11.65 -> 11.85
- [1:25:01] Using charge balance adjust
- [1:25:01] Dispensed 0.016298 mL of Base (0.5 M KOH)
- [1:25:22] Stirrer speed set to 0
- [1:25:32] Datapoint id 94 collected
- [1:25:32] Charge balance equation is out by -11.2%
- [1:25:32] Stirrer speed set to 60
- [1:25:37] pH 11.84 -> 12.04
- [1:25:37] Using charge balance adjust
- [1:25:37] Dispensed 0.025329 mL of Base (0.5 M KOH)
- [1:25:58] Stirrer speed set to 0
- [1:26:08] Datapoint id 95 collected
- [1:26:08] Charge balance equation is out by -9.5%
- [1:26:08] Argon flow rate set to 0
- [1:26:12] Titrator arm moved over Titration position