

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH2 +) -0.21 ±0.12 (n=50) logP (neutral XH) 2.63 ±0.02 (n=50) logP (X -) 0.46

RMSD 0.583

18C-09016 Points 2 to 28

M18_octanol concentration factor 1.091
Carbonate 0.1016 mM
Acidity error -0.00215 mM

18C-09016 Points 29 to 56

M18_octanol concentration factor 1.051
Carbonate 0.0884 mM
Acidity error -0.10411 mM

18C-09016 Points 57 to 86

M18_octanol concentration factor 1.063
Carbonate 0.1449 mM
Acidity error 0.00187 mM

Warnings and errors

Errors None Warnings None

Sample logD and percent species

M18_octanol	M18_octanol	M18_octanol	M18_octanoi	M18_octanoi	M18_octanol	M18_octanol	Comment
logD	M18_octanolH2	M18_octanolH	M18_octanol	M18_octanolH2*	M18_octanolH*	M18_octanol*	
-0.19	61.03 %	0.00 %	0.00 %	37.85 %	1.12 %	0.00 %	
-0.19	60.63 %	0.00 %	0.00 %	37.60 %	1.76 %	0.00 %	Stomach pH
-0.10	55.44 %	0.02 %	0.00 %	34.39 %	10.15 %	0.00 %	
0.39	28.95 %	0.12 %	0.00 %	17.95 %	52.98 %	0.00 %	
1.26	5.01 %	0.21 %	0.00 %	3.11 %	91.67 %	0.00 %	
2.11	0.54 %	0.23 %	0.00 %	0.34 %	98.89 %	0.00 %	
2.54	0.05 %	0.23 %	0.00 %	0.03 %	99.68 %	0.00 %	
2.60	0.02 %	0.23 %	0.00 %	0.01 %	99.74 %	0.00 %	
2.62	0.01 %	0.23 %	0.00 %	0.00 %	99.76 %	0.00 %	
2.63	0.00 %	0.23 %	0.00 %	0.00 %	99.76 %	0.00 %	Blood pH
2.63	0.00 %	0.23 %	0.00 %	0.00 %	99.76 %	0.00 %	- 1
2.62	0.00 %	0.23 %	0.01 %	0.00 %	99.75 %	0.02 %	ļ
2.55	0.00 %	0.23 %	0.05 %	0.00 %	99.57 %	0.15 %	ļ
2.13	0.00 %	0.23 %	0.51 %	0.00 %	97.79 %	1.47 %	
1.32	0.00 %	0.19 %	4.33 %	0.00 %	82.99 %	12.49 %	
	-0.19 -0.19 -0.10 0.39 1.26 2.11 2.54 2.60 2.62 2.63 2.63 2.62 2.55 2.13	logD M18_octanolH2 -0.19 61.03 % -0.19 60.63 % -0.10 55.44 % 0.39 28.95 % 1.26 5.01 % 2.11 0.54 % 2.54 0.05 % 2.60 0.02 % 2.62 0.01 % 2.63 0.00 % 2.62 0.00 % 2.55 0.00 % 2.13 0.00 %	logD M18_octanolH2 M18_octanolH -0.19 61.03 % 0.00 % -0.19 60.63 % 0.00 % -0.10 55.44 % 0.02 % 0.39 28.95 % 0.12 % 1.26 5.01 % 0.21 % 2.11 0.54 % 0.23 % 2.54 0.05 % 0.23 % 2.60 0.02 % 0.23 % 2.62 0.01 % 0.23 % 2.63 0.00 % 0.23 % 2.62 0.00 % 0.23 % 2.62 0.00 % 0.23 % 2.55 0.00 % 0.23 % 2.13 0.00 % 0.23 %	logD M18_octanolH2 M18_octanolH M18_octanolH M18_octanolH -0.19 61.03 % 0.00 % 0.00 % -0.19 60.63 % 0.00 % 0.00 % -0.10 55.44 % 0.02 % 0.00 % 0.39 28.95 % 0.12 % 0.00 % 1.26 5.01 % 0.21 % 0.00 % 2.11 0.54 % 0.23 % 0.00 % 2.54 0.05 % 0.23 % 0.00 % 2.60 0.02 % 0.23 % 0.00 % 2.62 0.01 % 0.23 % 0.00 % 2.63 0.00 % 0.23 % 0.00 % 2.63 0.00 % 0.23 % 0.01 % 2.62 0.00 % 0.23 % 0.01 % 2.63 0.00 % 0.23 % 0.01 % 2.62 0.00 % 0.23 % 0.01 % 2.63 0.00 % 0.23 % 0.01 % 2.62 0.00 % 0.23 % 0.05 % 2.55 0.00 %	logD M18_octanolH2 M18_octanolH M18_octanol M18_octanolH2* -0.19 61.03 % 0.00 % 0.00 % 37.85 % -0.19 60.63 % 0.00 % 0.00 % 37.60 % -0.10 55.44 % 0.02 % 0.00 % 34.39 % 0.39 28.95 % 0.12 % 0.00 % 17.95 % 1.26 5.01 % 0.21 % 0.00 % 3.11 % 2.11 0.54 % 0.23 % 0.00 % 0.34 % 2.54 0.05 % 0.23 % 0.00 % 0.01 % 2.60 0.02 % 0.23 % 0.00 % 0.01 % 2.62 0.01 % 0.23 % 0.00 % 0.00 % 2.63 0.00 % 0.23 % 0.00 % 0.00 % 2.63 0.00 % 0.23 % 0.00 % 0.00 % 2.62 0.00 % 0.23 % 0.00 % 0.00 % 2.63 0.00 % 0.23 % 0.00 % 0.00 % 2.63 0.00 % 0.23 % <t< td=""><td>logD M18_octanolH2 M18_octanolH M18_octanol M18_octanolH2* M18_octanolH* -0.19 61.03 % 0.00 % 0.00 % 37.85 % 1.12 % -0.19 60.63 % 0.00 % 0.00 % 37.60 % 1.76 % -0.10 55.44 % 0.02 % 0.00 % 34.39 % 10.15 % 0.39 28.95 % 0.12 % 0.00 % 17.95 % 52.98 % 1.26 5.01 % 0.21 % 0.00 % 3.11 % 91.67 % 2.11 0.54 % 0.23 % 0.00 % 0.34 % 98.89 % 2.54 0.05 % 0.23 % 0.00 % 0.03 % 99.68 % 2.60 0.02 % 0.23 % 0.00 % 0.01 % 99.74 % 2.62 0.01 % 0.23 % 0.00 % 0.00 % 99.76 % 2.63 0.00 % 0.23 % 0.00 % 0.00 % 99.76 % 2.62 0.00 % 0.23 % 0.00 % 0.00 % 99.75 % 2.62 0.00 %</td><td>logD M18_octanolH2 M18_octanolH M18_octanol M18_octanolH2* M18_octanolH* M18_octanol* -0.19 61.03 % 0.00 % 0.00 % 37.85 % 1.12 % 0.00 % -0.19 60.63 % 0.00 % 0.00 % 37.60 % 1.76 % 0.00 % -0.10 55.44 % 0.02 % 0.00 % 34.39 % 10.15 % 0.00 % 0.39 28.95 % 0.12 % 0.00 % 17.95 % 52.98 % 0.00 % 1.26 5.01 % 0.21 % 0.00 % 3.11 % 91.67 % 0.00 % 2.11 0.54 % 0.23 % 0.00 % 0.34 % 98.89 % 0.00 % 2.54 0.05 % 0.23 % 0.00 % 0.03 % 99.68 % 0.00 % 2.60 0.02 % 0.23 % 0.00 % 0.01 % 99.74 % 0.00 % 2.62 0.01 % 0.23 % 0.00 % 0.00 % 99.76 % 0.00 % 2.63 0.00 % 0.23 % 0.00 % 0.00 % <t< td=""></t<></td></t<>	logD M18_octanolH2 M18_octanolH M18_octanol M18_octanolH2* M18_octanolH* -0.19 61.03 % 0.00 % 0.00 % 37.85 % 1.12 % -0.19 60.63 % 0.00 % 0.00 % 37.60 % 1.76 % -0.10 55.44 % 0.02 % 0.00 % 34.39 % 10.15 % 0.39 28.95 % 0.12 % 0.00 % 17.95 % 52.98 % 1.26 5.01 % 0.21 % 0.00 % 3.11 % 91.67 % 2.11 0.54 % 0.23 % 0.00 % 0.34 % 98.89 % 2.54 0.05 % 0.23 % 0.00 % 0.03 % 99.68 % 2.60 0.02 % 0.23 % 0.00 % 0.01 % 99.74 % 2.62 0.01 % 0.23 % 0.00 % 0.00 % 99.76 % 2.63 0.00 % 0.23 % 0.00 % 0.00 % 99.76 % 2.62 0.00 % 0.23 % 0.00 % 0.00 % 99.75 % 2.62 0.00 %	logD M18_octanolH2 M18_octanolH M18_octanol M18_octanolH2* M18_octanolH* M18_octanol* -0.19 61.03 % 0.00 % 0.00 % 37.85 % 1.12 % 0.00 % -0.19 60.63 % 0.00 % 0.00 % 37.60 % 1.76 % 0.00 % -0.10 55.44 % 0.02 % 0.00 % 34.39 % 10.15 % 0.00 % 0.39 28.95 % 0.12 % 0.00 % 17.95 % 52.98 % 0.00 % 1.26 5.01 % 0.21 % 0.00 % 3.11 % 91.67 % 0.00 % 2.11 0.54 % 0.23 % 0.00 % 0.34 % 98.89 % 0.00 % 2.54 0.05 % 0.23 % 0.00 % 0.03 % 99.68 % 0.00 % 2.60 0.02 % 0.23 % 0.00 % 0.01 % 99.74 % 0.00 % 2.62 0.01 % 0.23 % 0.00 % 0.00 % 99.76 % 0.00 % 2.63 0.00 % 0.23 % 0.00 % 0.00 % <t< td=""></t<>



Assay ID:

Filename:

Sample name: M18_octanol Assay name:


pH-metric high logP

18C-09016

Experiment start time: 3/9/2018 11:02:44 PM

Pion Analyst: Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r






Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

Graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-09016 Points 2 to 28

Overall results

RMSD 0.731
Average ionic strength 0.156 M
Average temperature 24.9°C
Partition ratio 0.0626 : 1

Analyte concentration range 2691.8 µM to 2763.3 µM

Total points considered 25 of 27

Warnings and errors

Errors None Warnings None

Four-Plus parameters

à	Alpha	0.102	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r
à	S	0.9967	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r
à	jΗ	1.2	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r
à	jОН	0.0	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r


Titrants


0.50 M HCI 0.999843 3/9/2018 11:02:44 PM C:\Sirius_T3\HCl18C09.t3r 0.50 M KOH 0.999845 3/9/2018 11:02:44 PM C:\Sirius_T3\KOH18B27.t3r

Sample

7	M18_octanol concentration factor	1.091
à	Base pKa 1	5.37
à	Acid pKa 2	10.65
7	logP (XH2 +)	-1.02
7	logP (neutral XH)	2.58
à	logP (X -)	0.46

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M18_octanol	M18_octanol	M18_octanol	M18_octanol		M18_octanol	M18_octanol	Comment
	logD	M18_octanolH2	M18_octanolH	M18_octanol	M18_octanolH2*	M18_octanolH*	M18_octanol*	
1.000	-0.95	99.30 %	0.00 %	0.00 %	0.60 %	0.10 %	0.00 %	
1.200	-0.91	99.24 %	0.01 %	0.00 %	0.60 %	0.16 %	0.00 %	Stomach pH
2.000	-0.59	98.37 %	0.04 %	0.00 %	0.59 %	1.00 %	0.00 %	
3.000	0.23	89.97 %	0.38 %	0.00 %	0.54 %	9.11 %	0.00 %	
4.000	1.19	48.52 %	2.07 %	0.00 %	0.29 %	49.12 %	0.00 %	ľ
5.000	2.06	8.65 %	3.69 %	0.00 %	0.05 %	87.60 %	0.00 %	
6.000	2.49	0.94 %	4.00 %	0.00 %	0.01 %	95.05 %	0.00 %	
6.500	2.55	0.30 %	4.03 %	0.00 %	0.00 %	95.67 %	0.00 %	
7.000	2.57	0.09 %	4.04 %	0.00 %	0.00 %	95.86 %	0.00 %	
7.400	2.57	0.04 %	4.04 %	0.00 %	0.00 %	95.92 %	0.00 %	Blood pH
8.000	2.58	0.01 %	4.04 %	0.01 %	0.00 %	95.94 %	0.00 %	•
9.000	2.57	0.00 %	4.04 %	0.09 %	0.00 %	95.85 %	0.02 %	ļ
10.000	2.49	0.00 %	4.00 %	0.90 %	0.00 %	94.94 %	0.16 %	ļ
11.000	2.08	0.00 %	3.65 %	8.18 %	0.00 %	86.70 %	1.48 %	ļ
12.000	1.28	0.00 %	1.95 %	43.76 %	0.00 %	46.39 %	7.90 %	ļ

Carbonate and acidity



Other graphs






Assay ID: Filename:

Sample name: M18_octanol Experiment start time: 3/9/2018 11:02:44 PM

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-09016 Points 29 to 56

Overall results

RMSD 0.550
Average ionic strength 0.162 M
Average temperature 25.0°C
Partition ratio 0.1763 : 1

Analyte concentration range 2278.7 µM to 2341.3 µM

Total points considered 27 of 28

Warnings and errors

Errors None

Warnings One or more logP values out of range

Four-Plus parameters


Alpha	0.102	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r
S	0.9967	3/9/2018 11:02:43 PM	C:\Sirius T3\HCl18C09.t3r
jΗ	1.2	3/9/2018 11:02:43 PM	C:\Sirius T3\HCl18C09.t3r
jОН	0.0	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r


Titrants

Sample

₩	M18_octanol concentration factor	1.051
	Base pKa 1	5.37
	Acid pKa 2	10.65
₩	logP (XH2 +)	7.44
₩	logP (neutral XH)	9.35
À	logP (X -)	0.46

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M18_octanol logD	M18_octanol M18_octanolH2	M18_octanol M18_octanolH	M18_octanol M18_octanol		M18_octanol M18_octanolH*	M18_octanol M18_octanol*	
1.000	7.44	0.00 %	0.00 %	0.00 %	99.66 %	0.34 %	0.00 %	
1.200	7.44	0.00 %	0.00 %	0.00 %	99.46 %	0.54 %		Stomach pH
2.000	7.45	0.00 %	0.00 %	0.00 %	96.66 %	3.34 %	0.00 %	
3.000	7.57	0.00 %	0.00 %	0.00 %	74.32 %	25.68 %	0.00 %	
4.000	8.07	0.00 %	0.00 %	0.00 %	22.45 %	77.55 %	0.00 %	
5.000	8.84	0.00 %	0.00 %	0.00 %	2.81 %	97.19 %	0.00 %	
6.000	9.26	0.00 %	0.00 %	0.00 %	0.29 %	99.71 %	0.00 %	
6.500	9.32	0.00 %	0.00 %	0.00 %	0.09 %	99.91 %	0.00 %	
7.000	9.34	0.00 %	0.00 %	0.00 %	0.03 %	99.97 %	0.00 %	
7.400	9.34	0.00 %	0.00 %	0.00 %	0.01 %	99.99 %	0.00 %	Blood pH
8.000	9.35	0.00 %	0.00 %	0.00 %	0.00 %	100.00 %	0.00 %	•
9.000	9.34	0.00 %	0.00 %	0.00 %	0.00 %	100.00 %	0.00 %	
10.000	9.26	0.00 %	0.00 %	0.00 %	0.00 %	100.00 %	0.00 %	I
11.000	8.84	0.00 %	0.00 %	0.00 %	0.00 %	100.00 %	0.00 %	
12.000	7.98	0.00 %	0.00 %	0.00 %	0.00 %	100.00 %	0.00 %	

Carbonate and acidity



Other graphs





Assay ID: Filename:

Sample name: M18_octanol Assay name:

pH-metric high logP

18C-09016


Experiment start time: 3/9/2018 11:02:44 PM

Analyst: **Pion**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-09016 Points 57 to 86

Overall results

RMSD 0.447
Average ionic strength 0.168 M
Average temperature 25.0°C
Partition ratio 0.4400 : 1

Analyte concentration range 1742.6 µM to 1782.5 µM

Total points considered 28 of 30

Warnings and errors

Errors None

Warnings One or more logP values out of range

Four-Plus parameters


Alpha	0.102	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r
S	0.9967	3/9/2018 11:02:43 PM	C:\Sirius T3\HCl18C09.t3r
jΗ	1.2	3/9/2018 11:02:43 PM	C:\Sirius T3\HCl18C09.t3r
jОН	0.0	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r


Titrants

Sample

Φ	M18_octanol concentration factor	1.063
	Base pKa 1	5.37
	Acid pKa 2	10.65
₩	logP (XH2 +)	-3.10
₩	logP (neutral XH)	2.50
	logP (X -)	0.46

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M18_octanol	M18_octanol	M18_octanol	M18_octanol		M18_octanol	M18_octanol	Comment
	logD	M18_octanolH2	M18_octanolH	M18_octanol	M18_octanolH2*	M18_octanolH*	M18_octanol*	ļ
1.000	-1.84	99.37 %	0.00 %	0.00 %	0.04 %	0.59 %	0.00 %	ļ
1.200	-1.65	99.02 %	0.01 %	0.00 %	0.03 %	0.94 %	0.00 %	Stomach pH
2.000	-0.87	94.30 %	0.04 %	0.00 %	0.03 %	5.63 %	0.00 %	·
3.000	0.13	62.45 %	0.27 %	0.00 %	0.02 %	37.26 %	0.00 %	,
4.000	1.11	14.26 %	0.61 %	0.00 %	0.01 %	85.12 %	0.00 %	1
5.000	1.98	1.64 %	0.70 %	0.00 %	0.00 %	97.66 %	0.00 %	ľ
6.000	2.41	0.17 %	0.71 %	0.00 %	0.00 %	99.13 %	0.00 %	ļ
6.500	2.47	0.05 %	0.71 %	0.00 %	0.00 %	99.24 %	0.00 %	ļ
7.000	2.49	0.02 %	0.71 %	0.00 %	0.00 %	99.27 %	0.00 %	ľ
7.400	2.50	0.01 %	0.71 %	0.00 %	0.00 %	99.28 %	0.00 %	Blood pH
8.000	2.50	0.00 %	0.71 %	0.00 %	0.00 %	99.28 %	0.00 %	- 1
9.000	2.49	0.00 %	0.71 %	0.02 %	0.00 %	99.25 %	0.02 %	ļ
10.000	2.42	0.00 %	0.71 %	0.16 %	0.00 %	98.93 %	0.20 %	ļ
11.000	2.00	0.00 %	0.69 %	1.53 %	0.00 %	95.84 %	1.95 %	ļ
12.000	1.21	0.00 %	0.52 %	11.68 %	0.00 %	72.98 %	14.82 %	ļ

Carbonate and acidity



Other graphs






Assay name: pH-metric high logP Analyst: **Pion** 18C-09016 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

Other graphs (continued)





pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-09016 Instrument ID: T312060

Filename: $C:\Sirius_T3\\Mehtap\\20180309_exp31_logP_T3-2\\18C-09016_M18_octanol_pH-metric\ high\ logP.t3r$

Assay Model

Settings Sample name		Date/Time changed 2/27/2018 7:08:39 PM	
Sample by	Weight	2/0/2040 2:22:E2 DM	Default value
Sample weight Formula weight		3/9/2018 2:22:52 PM 2/27/2018 7:08:39 PM	
Solubility	Unknown	2/2//2010 / .00.59 F W	Default value
Molecular weight	267.11	2/27/2018 7:08:39 PM	
Individual pKa ionic environments	No		Default value
Number of pKas	2	2/27/2018 7:08:39 PM	User entered value
Sample is a	Ampholyte	2/27/2018 7:08:39 PM	User entered value
pKa 1	5.37	2/27/2018 7:08:39 PM	User entered value
Туре	Base	2/27/2018 7:08:39 PM	User entered value
pKa 2	10.65	2/27/2018 7:08:39 PM	User entered value
Туре	Acid	2/27/2018 7:08:39 PM	User entered value
logp (XH2 +)	-0.35	2/28/2018 3:20:28 PM	User entered value
logP (neutral XH)	2.57	3/2/2018 4:34:50 PM	User entered value
logP (X -)	0.46	2/27/2018 7:09:34 PM	User entered value

Events	;									
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/d time
4:58.7 4:59.8	Manual volume addition Initial pH = 7.99				0.10000 mL					unie
7:59.3	Data point 2	1.50000 mL	0.04812 mL	0.00477 mL	0.10000 mL	2.062	-0.00681	0.27502	0.00064	
8:45.5	Data point 3	1.50000 mL	0.04812 mL	0.01698 mL	0.10000 mL	2.268	0.00821	0.44883	0.00061	
9:21.1	Data point 4	1.50000 mL	0.04812 mL	0.02495 mL	0.10000 mL	2.471	-0.00977	0.70424	0.00058	
9:56.6	Data point 5	1.50000 mL	0.04812 mL	0.03001 mL	0.10000 mL	2.662	-0.00069	0.07464	0.00012	
10:32.1	Data point 6	1.50000 mL	0.04812 mL	0.03335 mL	0.10000 mL	2.850	-0.00288	0.33859	0.00024	_
11:07.6	Data point 7	1.50000 mL	0.04812 mL	0.03565 mL	0.10000 mL	3.085	-0.01383	0.83922	0.00074	s 10.0
11:43.3	Data point 8	1.50000 mL	0.04812 mL	0.03725 mL	0.10000 mL	3.283	-0.00418	0.74326	0.00024	10.0
12:18.8	Data point 9	1.50000 mL	0.04812 mL	0.03859 mL	0.10000 mL	3.456	-0.01450	0.85812	0.00077	10.0
13:09.7	Data point 10	1.50000 mL	0.04812 mL	0.04017 mL	0.10000 mL	3.653	-0.01030	0.40559	0.00080	
14:00.8	Data point 11	1.50000 mL	0.04812 mL	0.04156 mL	0.10000 mL	3.848	-0.00945	0.70466	0.00056	
14:36.7	Data point 12	1.50000 mL	0.04812 mL	0.04278 mL	0.10000 mL	4.049	-0.00976	0.86426	0.00052	
15:12.1	Data point 13	1.50000 mL	0.04812 mL	0.04396 mL	0.10000 mL	4.256	-0.00996	0.79638	0.00055	
15:47.4	Data point 14	1.50000 mL	0.04812 mL	0.04499 mL	0.10000 mL	4.460	-0.01429	0.84980	0.00077	
16:22.8	Data point 15	1.50000 mL	0.04812 mL	0.04584 mL	0.10000 mL	4.667	-0.01676	0.72668	0.00097	
16:58.2	Data point 16	1.50000 mL	0.04812 mL	0.04647 mL	0.10000 mL	4.852	-0.01624	0.82288	0.00088	s 10.5
17:34.1	Data point 17	1.50000 mL	0.04812 mL	0.04694 mL	0.10000 mL	4.977	-0.06194	0.99731	0.00306	s Time out at



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

Events (continued)

Time 19:04.6	Event Data point 18	Water 1.50000 mL	Acid 0.04812 mL	Base 0.04744 mL	Octanol 0.10000 mL	pH 5.138	dpH/dt -0.02805	pH R-squared 0.99638		dpH/dt time Timed out at 59.5 s
20:40.3	Data point 19	1.50000 mL	0.04812 mL	0.04774 mL	0.10000 mL	5.264	-0.01845	0.96248	0.00093	
21:34.7	Data point 20	1.50000 mL	0.04812 mL	0.04798 mL	0.10000 mL	5.641	-0.01992	0.97402	0.00100	21.0 s
22:26.3	Data point 21	1.50000 mL	0.04812 mL	0.04807 mL	0.10000 mL	5.855	-0.01912	0.93098	0.00098	26.5 s
23:18.2	Data point 22								0.00093	
24:11.1	Data point 23								0.00099	
25:39.1	Data point 24									Timed out at 59.5 s
27:14.7	Data point 25									Timed out at 59.5 s
28:45.1	Data point 26								0.00100	
30:08.9	Data point 27								0.00097	
31:15.5	Data point 28								0.00098	
32:47.3	Data point 29							0.00545	0.00033	
33:34.0	Data point 30							0.77392	0.00078	
34:09.7	Data point 31							0.41139	0.00011	
34:45.3	Data point 32								0.00012	
35:20.8	Data point 33								0.00072	
35:56.2	Data point 34								0.00045	
36:42.0	Data point 35								0.00067	
37:18.0	Data point 36								0.00033	
37:53.4	Data point 37								0.00032	
38:39.1	Data point 38								0.00039	
39:30.1	Data point 39								0.00072	
40:06.1	Data point 40								0.00036	
40:51.7	Data point 41								0.00039	
41:38.0	Data point 42								0.00064	
42:29.5	Data point 43								0.00082	
43:05.4	Data point 44								0.00088	
43:47.3	Data point 45								0.00074	
44:30.0 45:36.0	Data point 46 Data point 47								0.00099 0.00098	
45.36.0 47:02.6	Data point 48									Timed out at
48:28.0	Data point 49									59.5 s Timed out at
49:58.5	Data point 50									59.5 s Timed out at
51:29.0	Data point 51									59.5 s Timed out at
52:59.4	Data point 52									59.5 s Timed out at
	•									59.5 s
54:30.0	Data point 53									Timed out at 59.5 s
56:00.4	Data point 54								0.00098	
57:24.5	Data point 55								0.00090	
58:39.9	Data point 56								0.00099	
1:00:24.9									0.00059	
1:01:11.2								0.04813	0.00018	
1:01:46.9	Data point 59 Data point 60							0.09971 0.68746	0.00021 0.00072	
	Data point 61							0.56892	0.00072	
	Data point 62								0.00047	
	Data point 63							0.08135	0.00031	
	Data point 64								0.00039	
	Data point 65								0.00037	
1										

Assay Events



Sample name: M18_octanol Experiment start time: 3/9/2018 11:02:44 PM

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

Events (continued)

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
1:06:16.0	Data point 66	1.50000 mL	0.16011 mL	0.15691 mL	0.80000 mL	3.758	-0.00478	0.09544	0.00076	10.0 s
1:06:51.4	Data point 67	1.50000 mL	0.16011 mL	0.15793 mL	0.80000 mL	3.973	0.00088	0.00414	0.00068	10.0 s
1:07:26.8	Data point 68	1.50000 mL	0.16011 mL	0.15863 mL	0.80000 mL	4.186	-0.01639	0.67438	0.00099	10.0 s
1:08:02.3	Data point 69	1.50000 mL	0.16011 mL	0.15912 mL	0.80000 mL	4.378	-0.01336	0.45806	0.00098	10.0 s
1:08:37.7	Data point 70	1.50000 mL	0.16011 mL	0.15945 mL	0.80000 mL	4.589	-0.00985	0.60220	0.00063	11.0 s
1:09:14.1	Data point 71	1.50000 mL	0.16011 mL	0.15967 mL	0.80000 mL	4.736	-0.00475	0.10626	0.00072	11.5 s
1:09:56.1	Data point 72				0.80000 mL			0.37022	0.00096	12.0 s
1:10:43.8	Data point 73				0.80000 mL			0.93142	0.00098	34.5 s
1:11:48.9	Data point 74	1.50000 mL	0.16011 mL	0.16018 mL	0.80000 mL	5.755	-0.01516	0.77806	0.00085	45.5 s
1:13:05.0	Data point 75	1.50000 mL	0.16011 mL	0.16025 mL	0.80000 mL	6.047	-0.01855	0.94284	0.00094	50.5 s
1:14:26.2	Data point 76				0.80000 mL			0.96730	0.00099	58.0 s
1:15:49.6	Data point 77	1.50000 mL	0.16011 mL	0.16035 mL	0.80000 mL	6.474	-0.02760	0.99228	0.00137	Timed out at 59.5 s
1:17:20.1	Data point 78	1.50000 mL	0.16011 mL	0.16040 mL	0.80000 mL	6.897	-0.06395	0.99250	0.00317	Timed out at 59.5 s
1:18:45.5	Data point 79	1.50000 mL	0.16011 mL	0.16044 mL	0.80000 mL	7.185	-0.06453	0.99687	0.00319	Timed out at 59.5 s
1:20:16.1	Data point 80	1.50000 mL	0.16011 mL	0.16049 mL	0.80000 mL	7.570	-0.08020	0.98798	0.00398	Timed out at 59.5 s
1:21:46.6	Data point 81	1.50000 mL	0.16011 mL	0.16054 mL	0.80000 mL	7.910	-0.05115	0.99413	0.00253	Timed out at 59.5 s
1:23:17.1	Data point 82	1.50000 mL	0.16011 mL	0.16058 mL	0.80000 mL	8.105	-0.04177	0.93249	0.00214	Timed out at 59.5 s
1:24:52.7	Data point 83	1.50000 mL	0.16011 mL	0.16065 mL	0.80000 mL	8.398	-0.02637	0.89962	0.00137	Timed out at 59.5 s
1:26:33.5	Data point 84	1.50000 mL	0.16011 mL	0.16075 mL	0.80000 mL	8.655	-0.01834	0.88375	0.00096	28.0 s
1:27:37.2		1.50000 mL	0.16011 mL	0.16084 mL	0.80000 mL	8.817	-0.01433	0.54002	0.00096	16.5 s
1:28:29.4	Data point 86	1.50000 mL	0.16011 mL	0.16096 mL	0.80000 mL	9.030	-0.01224	0.37247	0.00099	14.0 s
1:28:52.4	Assay volumes	1.50000 mL	0.16011 mL	0.16096 mL	0.80000 mL					



Assay name: Analyst: pH-metric high logP Pion Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings	Diam			
Analyst name	Pion			
Standard Experiment Settings	0			
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	9.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL 0.10000 mL			
Maximum titrant addition	100%			
Argon flow rate				
Start titration using	Cautious pH adjust			
Advanced General Settings	Mana			
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose	Nama			
Titrant pre-dose	None			
Assay Medium	1 50 ml			
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.100 mL			
Partition solvent added	Manual in advance			
After partition addition, stir for	1 seconds			
Sample Sonication	Vaa			
Sonicate	Yes			
Adjust pH for sonication	No co			
Sonicate for	60 seconds			
After sonication stir for	5 seconds			
Sample Dissolution	V			
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution	To start pH			
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge	NI.			
Perform a carbonate purge	No			
Temperature Control	Vaa			
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			
Stir speed of	50%			
Titration 1	المطامئة ما منا			
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	50%			
Titration 2	Lance to John John all			
Titrate from	Low to high pH			
Add additional water	0.00 mL			
•	0.200 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			

Report by: Dorothy Levorse 3/16/2018 1:54:28 PM

Stir to allow partitioning for

Stirrer speed for partitioning

15 seconds

55%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3			_	
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.500 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00100 dpH/dt			
Stability timeout after	60 seconds			

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.102	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r
Four-Plus S	0.9967	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r
Four-Plus jH	1.2	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r
Four-Plus jOH	0.0	3/9/2018 11:02:43 PM	C:\Sirius_T3\HCl18C09.t3r
Base concentration factor	1.000	3/9/2018 11:02:44 PM	C:\Sirius_T3\KOH18B27.t3r
Acid concentration factor	1.000	3/9/2018 11:02:44 PM	C:\Sirius_T3\HCl18C09.t3r

Instrument Settings

Setting	Value	Batch Id	Install date
Instrument owner	Merck		
Instrument ID	T312060		
Instrument type	T3 Simulator		
Software version	1.1.3.0		
Dispenser module		T3DM1200361	3/31/2009 6:24:52 AM
Dispenser 0	Water		3/31/2009 6:25:05 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Water (0.15 M KCI)	02-06-2018	2/27/2018 11:05:59 AM
Dispenser 2	Acid		3/31/2009 6:25:11 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Acid (0.5 M HCI)	02-27-2018	2/27/2018 11:27:22 AM
Dispenser 1	Base		3/31/2009 6:25:21 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Base (0.5 M KOH)	9/22/2017	2/27/2018 11:21:22 AM
Dispenser 5	Cosolvent		3/31/2009 6:26:24 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Distribution valve 5	Distribution Valve		3/31/2009 6:28:19 AM
Firmware version	1.1.3		
Port A	Methanol (80%, 0.15 M KCI)		3/6/2018 10:28:59 AM
Port B	Cyclohexane	11-01-17	2/27/2018 11:37:57 AM
Dispenser 3	Buffer		8/3/2010 6:05:16 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Dodecane	2018/01/31	2/28/2018 11:18:04 AM
Dispenser 6	Octanol		10/22/2010 11:52:43 AM




Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting Syringe volume	Value 0.5 mL	Batch Id	Install date
Firmware version	1.2.1(r2)		
Titrant	Octanol	01-31-2018	2/27/2018 10:59:35 AM
Titrator	Cotanor		3/31/2009 6:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2	101W1200101	0/01/2000 0.24.17 /tivi
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0923	1/23/2018 3:01:00 PM
E0 calibration	+4.01 mV	1020020	3/9/2018 11:03:12 PM
Filling solution	3M KCI	KCL097	3/9/2018 11:05:42 AM
Liquids		NOL007	0/0/2010 11:00:12740
Wash 1	50% IPA:50% Water		3/9/2018 11:04:22 AM
Wash 2	0.5% Trition X-100 in H20		3/9/2018 11:04:25 AM
Buffer position 1	pH7 Wash		3/9/2018 11:04:27 AM
Buffer position 2	pH 7		3/9/2018 11:04:30 AM
Storage position	p		3/9/2018 11:05:04 AM
Wash water	5e+003 mL	02-27-2018	2/27/2018 10:54:39 AM
Waste	1.1e+004 mL	02 27 2010	11/28/2017 11:36:29 AM
Temperature controller	1.10.0011112		8/5/2010 7:35:13 AM
Turbidity detector			3/31/2009 6:24:45 AM
Spectrometer		074811	11/23/2010 12:22:28 PM
Dip probe		10196	11/20/2010 12:22:201 101
Wavelength coefficient A0	183.333	10100	
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	123:16:41		11/23/2010 12:22:28 PM
Calibrated on	2/27/2018 11:40:38 AM		11/20/2010 12:22:201 101
Integration time	40		
Scans averaged	10		
Autoloader		T3AL1200345	11/10/2015 10:34:13 AM
Left-right axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Configuration			
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		
5 ,			





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09016 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09016_M18_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		

Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00
• •		