

18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

## pH-metric Result

logP (XH +) 1.14 ±0.04 (n=50) logP (neutral X) 4.16 ±0.01 (n=50)

#### 18C-06015 Points 1 to 20

M02 octanol concentration factor 0.918 Carbonate Acidity error -0.42538 mM

#### 18C-06015 Points 21 to 34

M02\_octanol concentration factor 0.898 Carbonate 0.0001 mM Acidity error -0.16274 mM

#### 18C-06015 Points 35 to 48

M02 octanol concentration factor 0.730 Carbonate 0.1987 mM Acidity error -0.19134 mM

## Warnings and errors

Errors None Warnings None

рΗ

## Sample logD and percent species

M02\_octanol M02\_octanol M02\_octanol

	logD	M02_octanolH	M02_octanol	M02_octanolH*	M02_octanol*	
1.000	1.18	<del>6</del> .22 %	0.00 %	85.42 %	8.36 %	
1.200	1.20	5.93 %	0.00 %	81.44 %	12.64 %	Stomach pH
2.000	1.43	3.55 %	0.00 %	48.74 %	47.71 %	
3.000	2.17	0.67 %	0.01 %	9.21 %	90.12 %	
4.000	3.09	0.07 %	0.01 %	1.01 %	98.91 %	
5.000	3.84	0.01 %	0.01 %	0.10 %	99.88 %	
6.000	4.11	0.00 %	0.01 %	0.01 %	99.98 %	
6.500	4.14	0.00 %	0.01 %	0.00 %	99.99 %	
7.000	4.15	0.00 %	0.01 %	0.00 %	99.99 %	
7.400	4.16	0.00 %	0.01 %	0.00 %	99.99 %	Blood pH
8.000	4.16	0.00 %	0.01 %	0.00 %	99.99 %	
9.000	4.16	0.00 %	0.01 %	0.00 %	99.99 %	
10.000	4.16	0.00 %	0.01 %	0.00 %	99.99 %	
11.000	4.16	0.00 %	0.01 %	0.00 %	99.99 %	
12.000	4.16	0.00 %	0.01 %	0.00 %	99.99 %	

M02\_octanol Comment



Assay ID:

Filename:


Sample name: M02\_octanol Assay name: pH-metric high logP

18C-06015

Experiment start time: 3/6/2018 5:41:59 PM Analyst: **Dorothy Levorse** 

Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r






Filename:

Sample name: M02\_octanol Experiment start time: 3/6/2018 5:41:59 PM pH-metric high logP Analyst: **Dorothy Levorse** Assay name: Assay ID:

18C-06015 Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# **Graphs** (continued)





Sample name: M02\_octanol Experiment start time: 3/6/2018 5:41:59 PM Analyst: Assay name: pH-metric high logP **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-06015

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 1 of 3 18C-06015 Points 1 to 20

#### Overall results

RMSD 2.836 Average ionic strength 0.157 M Average temperature 24.9°C Partition ratio 0.0122:1

Analyte concentration range 4481.9 μM to 4615.4 μM

Total points considered 16 of 20

## Warnings and errors

Errors None Warnings None

## Four-Plus parameters

Alpha 0.124 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r S 0.9973 jΗ 0.9 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r jΟH -0.7

## Titrants


0.50 M HCI 0.989131 3/6/2018 5:41:59 PM C:\Sirius T3\18C-06006 Blank standardisation.t3r


0.50 M KOH 0.999845 3/6/2018 5:41:59 PM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M02\_octanol concentration factor 0.918 Base pKa 1 5.03 logP(XH +)1.32 logP (neutral X) 4.20

#### Sample graphs








Instrument ID: T312060 18C-06015

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# Sample graphs (continued)



## Sample logD and percent species


рН	M02_octanol	M02_octanol	M02_octanol	_	M02_octanol	Comment
	logD	M02_octanolH	MU2_octanoi	M02_octanolH*	M02_octanol*	
1.000	1.35	78.52 %	0.01 %	20.04 %	1.43 %	
1.200	1.37	77.87 %	0.01 %	19.87 %	2.25 %	Stomach pH
2.000	1.55	69.52 %	0.06 %	17.74 %	12.68 %	•
3.000	2.23	32.38 %	0.30 %	8.26 %	59.05 %	
4.000	3.14	5.11 %	0.48 %	1.30 %	93.11 %	
5.000	3.89	0.54 %	0.51 %	0.14 %	98.81 %	
6.000	4.16	0.05 %	0.51 %	0.01 %	99.42 %	
6.500	4.19	0.02 %	0.51 %	0.00 %	99.47 %	
7.000	4.20	0.01 %	0.51 %	0.00 %	99.48 %	
7.400	4.20	0.00 %	0.51 %	0.00 %	99.49 %	Blood pH
8.000	4.20	0.00 %	0.51 %	0.00 %	99.49 %	·
9.000	4.20	0.00 %	0.51 %	0.00 %	99.49 %	
10.000	4.20	0.00 %	0.51 %	0.00 %	99.49 %	
11.000	4.20	0.00 %	0.51 %	0.00 %	99.49 %	
12.000	4.20	0.00 %	0.51 %	0.00 %	99.49 %	

# Carbonate and acidity



Carbonate 0.000 mM Acidity error -0.425 mM

# Other graphs





Assay ID: Filename:

Sample name: M02\_octanol Assay name:

pH-metric high logP


18C-06015

Experiment start time: 3/6/2018 5:41:59 PM Analyst: **Dorothy Levorse** 

Instrument ID: T312060

 $C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric\ high\ logP.t3r$ 

# Other graphs (continued)





Sample name: M02\_octanol Experiment start time: 3/6/2018 5:41:59 PM Analyst: Assay name: pH-metric high logP **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-06015

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 2 of 3 18C-06015 Points 21 to 34

## Overall results

RMSD 4.143 Average ionic strength 0.163 M Average temperature 25.0°C Partition ratio 0.0289:1

Analyte concentration range 4116.8 µM to 4254.1 µM

Total points considered 11 of 14

## Warnings and errors

Errors None Warnings None

## Four-Plus parameters

Alpha 0.124 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r S 0.9973 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r jΗ 0.9 jOH 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r -0.7

## Titrants


0.50 M HCI 0.989131 3/6/2018 5:41:59 PM C:\Sirius T3\18C-06006 Blank standardisation.t3r


0.50 M KOH 0.999845 3/6/2018 5:41:59 PM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M02\_octanol concentration factor 0.898 Base pKa 1 5.03 logP(XH +)1.32 logP (neutral X) 4.19

#### Sample graphs





Analyst:

Experiment start time: 3/6/2018 5:41:59 PM

**Dorothy Levorse** 



Assay ID:

Filename:


Sample name: M02\_octanol Assay name: pH-metric high logP

18C-06015

C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

Instrument ID: T312060

# Sample graphs (continued)



## Sample logD and percent species


рН	M02_octanol	M02_octanol	M02_octanol	M02_octanol	M02_octanol	Comment
	logD	M02_octanolH	M02_octanol	M02_octanolH*	M02_octanol*	
1.000	1.35	60.82 %	0.01 %	36.66 %	2.52 %	
1.200	1.36	59.93 %	0.01 %	36.13 %	3.93 %	Stomach pH
2.000	1.55	49.56 %	0.05 %	29.87 %	20.52 %	
3.000	2.21	17.38 %	0.16 %	10.48 %	71.97 %	
4.000	3.12	2.32 %	0.22 %	1.40 %	96.06 %	
5.000	3.87	0.24 %	0.22 %	0.14 %	99.39 %	
6.000	4.14	0.02 %	0.22 %	0.01 %	99.74 %	
6.500	4.17	0.01 %	0.22 %	0.00 %	99.76 %	
7.000	4.18	0.00 %	0.22 %	0.00 %	99.77 %	
7.400	4.19	0.00 %	0.22 %	0.00 %	99.77 %	Blood pH
8.000	4.19	0.00 %	0.22 %	0.00 %	99.77 %	
9.000	4.19	0.00 %	0.22 %	0.00 %	99.78 %	
10.000	4.19	0.00 %	0.22 %	0.00 %	99.78 %	
11.000	4.19	0.00 %	0.22 %	0.00 %	99.78 %	
12.000	4.19	0.00 %	0.22 %	0.00 %	99.78 %	


# Carbonate and acidity



Carbonate 0.000 mM Acidity error -0.163 mM

# Other graphs







Assay ID:

Filename:

Sample name: M02\_octanol Assay name:

pH-metric high logP


18C-06015

Experiment start time: 3/6/2018 5:41:59 PM Analyst: **Dorothy Levorse** 

Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Sample name: M02\_octanol Experiment start time: 3/6/2018 5:41:59 PM Analyst: Assay name: pH-metric high logP **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-06015

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 3 of 3 18C-06015 Points 35 to 48

#### Overall results

RMSD 0.729 Average ionic strength 0.169 M Average temperature 25.0°C Partition ratio 0.1637:1

Analyte concentration range 3405.3 µM to 3503.8 µM

Total points considered 14 of 14

## Warnings and errors

Errors None Warnings None

## Four-Plus parameters

Alpha 0.124 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r S 0.9973 jΗ 0.9 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r jOH 3/6/2018 5:41:59 PM C:\Sirius\_T3\18C-06006\_Blank standardisation.t3r -0.7

#### Titrants


0.50 M HCI 0.989131 3/6/2018 5:41:59 PM C:\Sirius T3\18C-06006 Blank standardisation.t3r


0.50 M KOH 0.999845 3/6/2018 5:41:59 PM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M02 octanol concentration factor 0.730 Base pKa 1 5.03 logP(XH +)1.32 logP (neutral X) 4.31

#### Sample graphs





Analyst:

Experiment start time: 3/6/2018 5:41:59 PM

**Dorothy Levorse** 




Assay ID:

Sample name: M02\_octanol Assay name: pH-metric high logP

Instrument ID: T312060 18C-06015

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# Sample graphs (continued)



## Sample logD and percent species


рН	M02_octanol	M02_octanol	M02_octanol	M02_octanol	M02_octanol	Comment
	logD	M02_octanolH	M02_octanol	M02_octanolH*	M02_octanol*	
1.000	1.36	21.15 %	0.00 %	72.32 %	6.53 %	
1.200	1.38	20.37 %	0.00 %	69.66 %	9.97 %	Stomach pH
2.000	1.60	13.32 %	0.01 %	45.54 %	41.14 %	
3.000	2.32	2.83 %	0.03 %	9.68 %	87.46 %	
4.000	3.24	0.32 %	0.03 %	1.09 %	98.56 %	
5.000	3.99	0.03 %	0.03 %	0.11 %	99.83 %	
6.000	4.26	0.00 %	0.03 %	0.01 %	99.96 %	
6.500	4.29	0.00 %	0.03 %	0.00 %	99.97 %	
7.000	4.30	0.00 %	0.03 %	0.00 %	99.97 %	
7.400	4.30	0.00 %	0.03 %	0.00 %	99.97 %	Blood pH
8.000	4.31	0.00 %	0.03 %	0.00 %	99.97 %	
9.000	4.31	0.00 %	0.03 %	0.00 %	99.97 %	
10.000	4.31	0.00 %	0.03 %	0.00 %	99.97 %	
11.000	4.31	0.00 %	0.03 %	0.00 %	99.97 %	
12.000	4.31	0.00 %	0.03 %	0.00 %	99.97 %	


# Carbonate and acidity



Carbonate 0.199 mM Acidity error -0.191 mM

# Other graphs





Reported at: 3/9/2018 11:30:33 AM



Assay ID: Filename:

Sample name: M02\_octanol Assay name:

pH-metric high logP


18C-06015

Experiment start time: 3/6/2018 5:41:59 PM **Dorothy Levorse** Analyst:

Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# **Assay Model**

Value	Date/Time changed	Imported from
M02_octanol	12/6/2017 4:20:03 PM	User entered value
Weight		Default value
0.002110 g	3/6/2018 3:40:26 PM	User entered value
289.26 g/mol	12/6/2017 4:20:03 PM	User entered value
Unknown		Default value
289.26	12/6/2017 4:20:03 PM	User entered value
No		Default value
1	12/6/2017 4:20:03 PM	User entered value
Base	12/6/2017 4:20:03 PM	User entered value
5.03	12/6/2017 4:20:03 PM	User entered value
1.32	3/2/2018 3:38:13 PM	User entered value
4.10	3/2/2018 3:38:07 PM	User entered value
	M02_octanol Weight 0.002110 g 289.26 g/mol Unknown 289.26 No 1 Base 5.03 1.32	M02_octanol

#### **Events**

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
5:08.1	Initial pH = 5.88					1	Ja		F	
8:07.7	Data point 1	1.50000 mL	0.05466 mL	0.00583 mL	0.01999 mL	2.051	0.00164	0.09342	0.00026	10.0 s
8:53.9	Data point 2		0.05466 mL						0.00029	
9:30.0	Data point 3		0.05466 mL						0.00015	
10:06.1	Data point 4		0.05466 mL						0.00062	
10:41.7	Data point 5	1.50000 mL	0.05466 mL	0.03996 mL	0.01999 mL	2.883	-0.00877	0.93188	0.00045	10.0 s
11:17.2	Data point 6	1.50000 mL	0.05466 mL	0.04363 mL	0.01999 mL	3.088	-0.01598	0.80563	0.00088	10.0 s
11:52.7	Data point 7	1.50000 mL	0.05466 mL	0.04652 mL	0.01999 mL	3.307	-0.01840	0.83749	0.00099	12.0 s
12:30.2	Data point 8	1.50000 mL	0.05466 mL	0.04871 mL	0.01999 mL	3.520	-0.01804	0.93551	0.00092	16.0 s
13:11.8	Data point 9	1.50000 mL	0.05466 mL	0.05031 mL	0.01999 mL	3.755	-0.01934	0.98232	0.00096	20.0 s
14:07.6	Data point 10		0.05466 mL						0.00095	33.5 s
15:06.5	Data point 11	1.50000 mL	0.05466 mL	0.05151 mL	0.01999 mL	4.071	-0.01887	0.88004	0.00099	34.0 s
16:06.0			0.05466 mL						0.00089	37.5 s
	Data point 13		0.05466 mL						0.00099	40.0 s
18:29.8			0.05466 mL						0.00099	
19:38.8	Data point 15		0.05466 mL						0.00093	46.5 s
20:55.9	Data point 16		0.05466 mL						0.00096	
22:16.0	Data point 17		0.05466 mL						0.00095	
23:53.8	Data point 18		0.05466 mL						0.00099	
	Data point 19		0.05466 mL						0.00095	
26:33.4			0.05466 mL						0.00095	
27:57.2			0.11178 mL						0.00043	
	Data point 22		0.11178 mL						0.00066	
29:19.2			0.11178 mL						0.00044	
29:54.8	Data point 24		0.11178 mL						0.00056	
30:30.9			0.11178 mL						0.00063	
31:06.8			0.11178 mL						0.00065	
	Data point 27		0.11178 mL						0.00073	
	Data point 28		0.11178 mL						0.00038	
32:53.6	•		0.11178 mL						0.00082	
33:45.0	•		0.11178 mL						0.00059	
34:20.9	Data point 31		0.11178 mL						0.00100	
35:10.1	Data point 32	1.50000 mL	0.11178 mL	0.10992 mL	0.05000 mL	7.825	-0.06475	0.99712	0.00320	Timed out at 59.5 s
36:45.7	Data point 33	1.50000 mL	0.11178 mL	0.11011 mL	0.05000 mL	8.913	-0.01987	0.97328	0.00099	
37:37.5	Data point 34		0.11178 mL						0.00098	
39:14.6	Data point 35		0.17180 mL					0.00822	0.00082	
40:01.6	Data point 36		0.17180 mL						0.00094	
40:38.3	Data point 37		0.17180 mL						0.00083	
41:13.9	Data point 38		0.17180 mL						0.00082	
41:50.1	Data point 39		0.17180 mL					0.62249	0.00079	
	Data point 40		0.17180 mL					0.63323	0.00080	
	,									

Reported at: 3/9/2018 11:30:33 AM





18C-06015 Instrument ID: T312060

 $C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric\ high\ logP.t3r$ Filename:

# **Events (continued)**

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
43:11.9	Data point 41	1.50000 mL	0.17180 mL	0.16597 mL	0.30000 mL	3.316	-0.01737	0.82670	0.00094	14.5 s
43:57.2	Data point 42	1.50000 mL	0.17180 mL	0.16710 mL	0.30000 mL	3.518	-0.01034	0.83937	0.00056	10.5 s
44:48.5	Data point 43	1.50000 mL	0.17180 mL	0.16799 mL	0.30000 mL	3.731	0.01219	0.49410	0.00086	10.0 s
45:23.9	Data point 44	1.50000 mL	0.17180 mL	0.16858 mL	0.30000 mL	4.002	-0.01494	0.80259	0.00082	10.0 s
46:04.5	Data point 45	1.50000 mL	0.17180 mL	0.16900 mL	0.30000 mL	4.312	-0.01556	0.81871	0.00085	10.5 s
46:45.7	Data point 46	1.50000 mL	0.17180 mL	0.16928 mL	0.30000 mL	4.694	-0.01928	0.95047	0.00098	29.0 s
47:50.4	Data point 47	1.50000 mL	0.17180 mL	0.16980 mL	0.30000 mL	8.255	-0.04492	0.97778	0.00224	Timed out at
										59.5 s
	Data point 48					9.190	-0.01473	0.61216	0.00093	15.0 s
49:55.4	Assay volumes	1.50000 mL	0.17180 mL	0.17034 mL	0.30000 mL					



Assav ID: 18C-06015 Instrument ID: T312060

Assay ID: 18C-06015 Filename: C:\Sirius_T3\Meh	ntap\20180306_exp3			312060 l_pH-metric high logP.
Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				-
Analyst name	Dorothy Levorse			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	9.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings	N.I.			
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose	Nana			
Titrant pre-dose	None			
Assay Medium	1 F0 mal			
ISA water volume	1.50 mL			
Water added	Automatic Octanol			
Partition solvent type Partition volume	0.020 mL			
Partition volume Partition solvent added	Automatic			
After partition addition, stir for	1 seconds			
Sample Sonication	i seconus			
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	120 seconds			
After sonication stir for	5 seconds			
Sample Dissolution	0 00001100			
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution	To start pH			
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge	2 / <del>-</del>			
Perform a carbonate purge	No			

#### Temperature Control

Yes Wait for temperature Required start temperature 25.0°C Acceptable deviation 0.5°C Time to wait 60 seconds Stir speed of 50%

Titration 1

Titrate from Low to high pH

Adjust to start pH Yes After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds Stirrer speed for partitioning 50%

Titration 2

Stirrer speed for partitioning

Titrate from Low to high pH Add additional water 0.00 mL Additional partition solvent volume 0.030 mL Additional partition solvent added Automatic After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds

Reported at: 3/9/2018 11:30:33 AM

55%



Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3				
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.250 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00100 dpH/dt			
Stability timeout after	60 seconds			

## **Calibration Settings**

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.124	3/6/2018 5:41:59 PM	C:\Sirius_T3\18C-06006_Blank standardisation.t3r
Four-Plus S	0.9973	3/6/2018 5:41:59 PM	C:\Sirius_T3\18C-06006_Blank standardisation.t3r
Four-Plus jH	0.9	3/6/2018 5:41:59 PM	C:\Sirius_T3\18C-06006_Blank standardisation.t3r
Four-Plus jOH	-0.7	3/6/2018 5:41:59 PM	C:\Sirius_T3\18C-06006_Blank standardisation.t3r
Base concentration factor	1.000	3/6/2018 5:41:59 PM	C:\Sirius_T3\KOH18B27.t3r
Acid concentration factor	0.989	3/6/2018 5:41:59 PM	C:\Sirius_T3\18C-06006_Blank standardisation.t3r

## **Instrument Settings**

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module Dispenser 0 Syringe volume Firmware version	Water 2.5 mL 1.2.1(r2)	T3DM1200361	3/31/2009 5:24:52 AM 3/31/2009 5:25:05 AM
Titrant Dispenser 2 Syringe volume Firmware version	Water (0.15 M KCI) Acid 0.5 mL 1.2.1(r2)	02-06-2018	2/27/2018 10:05:59 AM 3/31/2009 5:25:11 AM
Titrant Dispenser 1 Syringe volume Firmware version	Acid (0.5 M HCI) Base 0.5 mL 1.2.1(r2)	02-27-2018	2/27/2018 10:27:22 AM 3/31/2009 5:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2)	9/22/2017	2/27/2018 10:21:22 AM 3/31/2009 5:26:24 AM
Distribution valve 5 Firmware version Port A	Distribution Valve 1.1.3 Methanol (80%, 0.15 M KCI)	02-08-2018	3/31/2009 5:28:19 AM 3/6/2018 9:28:59 AM
Port B Dispenser 3 Syringe volume Firmware version	Cyclohexane Buffer 0.5 mL 1.2.1(r2)	11-01-17	2/27/2018 10:37:57 AM 8/3/2010 5:05:16 AM
Titrant Dispenser 6	Dodecane Octanol	2018/01/31	2/28/2018 10:18:04 AM 10/22/2010 10:52:43 AM

Reported at: 3/9/2018 11:30:33 AM



Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

# Instrument Settings (continued)

Setting Syringe volume	<b>Value</b> 0.5 mL	Batch Id	Install date
Firmware version Titrant	1.2.1(r2) Octanol	01-31-2018 T3TM1200161	2/27/2018 9:59:35 AM
Titrator Horizontal axis firmware version Vertical axis firmware version Chassis I/O firmware version	1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2 1.11 Al1Dl0DO4 Norgren I/O 1.1.1	1311111200161	3/31/2009 5:24:17 AM
Probe I/O firmware version Electrode E0 calibration	T3 Electrode +7.02 mV	T3E0923	1/23/2018 2:01:00 PM 3/6/2018 5:42:27 PM
Filling solution Liquids	3M KCI	KCL097	3/6/2018 9:23:20 AM
Wash 1 Wash 2 Buffer position 1 Buffer position 2 Storage position Wash water	50% IPA:50% Water 0.5% Trition X-100 in H20 pH7 Wash pH 7 5.8e+003 mL	02-27-2018	3/6/2018 9:24:32 AM 3/6/2018 9:24:35 AM 3/6/2018 9:24:38 AM 3/6/2018 9:24:40 AM 3/6/2018 9:24:07 AM 2/27/2018 9:54:39 AM
Waste Temperature controller Turbidity detector	9.8e+003 mL		11/28/2017 10:36:29 AM 8/5/2010 6:35:13 AM 3/31/2009 5:24:45 AM
Spectrometer Dip probe	400.000	074811 10196	11/23/2010 11:22:28 AM
Wavelength coefficient A0 Wavelength coefficient A1 Wavelength coefficient A2	183.333 2.21568 -0.000289308		
Total lamp lit time Calibrated on Integration time	123:01:40 2/27/2018 10:40:38 AM 40		11/23/2010 11:22:28 AM
Scans averaged	10		
Autoloader Left-right axis firmware version Front-back axis firmware version Vertical axis firmware version Chassis I/O firmware version	1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2 1.11 Al1Dl0DO4 Norgren I/O	T3AL1200345	11/10/2015 9:34:13 AM
Configuration	1.11741Blobb 1 Holgion #6		
Alternate titration position Alternate reference position	Titration position Reference position		
Maximum standard vial volume Maximum alternate vial volume Automatic action idle period Titrant tube volume	3.50 mL 25.00 mL 5 minute(s) 1.3 mL		
Syringe flush count Flowing wash pump volume	3.50 20.0 mL		
Flowing wash stir duration Flowing wash stir speed Solvent wash stir duration	5 s 30% 5 s		
Solvent wash stir speed Surfactant wash stir duration Surfactant wash stir speed E0 calibration minimum number of points	30% 5 s 30% 10		
E0 calibration maximum standard deviation E0 calibration timeout period E0 calibration stir duration	0.01500 60 s 5 s		
E0 calibration preparation stir speed E0 calibration buffer wash stir duration E0 calibration buffer wash stir speed	30% 5 s 30%		
E0 calibration reading stir speed	0%		

Reported at: 3/9/2018 11:30:33 AM



Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

## Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		

## Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00

#### Experiment Log

- [49] Air gap released for Acid (0.5 M HCI) [49] Air gap released for Base (0.5 M KOH)
- [1:45] Air gap created for Water (0.15 M KCl)
- [1:46] Air gap created for Acid (0.5 M HCl)
- [1:46] Air gap created for Base (0.5 M KOH)
- [1:46] Air gap released for Water (0.15 M KCI)
  - 1.40] All gap released for water (0.15 W KCI)
- [1:50] Titrator arm moved over Titration position
- [1:50] Titration 1 of 3
- [1:50] Adding initial titrants
- [1:50] Automatically add 1.50000 mL of water
- [2:15] Dispensed 1.500000 mL of Water (0.15 M KCI)
- [2:20] Titrator arm moved over Drain
- [5:01] Titrator arm moved to Titration position
- [5:01] Argon flow rate set to 100
- [5:01] Stirrer speed set to 10
- [5:06] Automatically add 0.02000 mL of Octanol
- [5:07] Dispensed 0.019991 mL of Octanol
- [5:08] Initial pH = 5.88
- [5:08] Iterative adjust 5.88 -> 2.00
- [5:08] pH 5.88 -> 2.00
- [5:10] Air gap released for Acid (0.5 M HCl)
- [5:10] Dispensed 0.054657 mL of Acid (0.5 M HCI)
- [5:15] Holding pH 2.00
- [7:15] Stirrer speed set to 0
- [7:15] Stirrer speed set to 50
- [7:15] Iterative adjust 1.95 -> 2.00
- [7:15] pH 1.95 -> 2.00
- [7:16] Air gap released for Base (0.5 M KOH)
- [7:17] Dispensed 0.005833 mL of Base (0.5 M KOH)
- [8:07] Stirrer speed set to 0
- [8:17] Datapoint id 1 collected
- [8:17] Stirrer speed set to 50
- [8:22] pH 2.06 -> 2.26
- [8:22] Using cautious pH adjust
- [8:23] Dispensed 0.007197 mL of Base (0.5 M KOH)
- [8:28] Stepping pH = 2.15
- [8:28] Dispensed 0.005127 mL of Base (0.5 M KOH)
- [8:33] Stepping pH = 2.23
- [8:33] Dispensed 0.001317 mL of Base (0.5 M KOH)
- [8:38] Stepping pH = 2.26
- [8:54] Stirrer speed set to 0
- [9:04] Datapoint id 2 collected
- Reported at: 3/9/2018 11:30:33 AM Page 18 of 26



Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

- [9:04] Charge balance equation is out by 5.2%
- [9:04] Stirrer speed set to 50
- [9:09] pH 2.26 -> 2.46
- [9:09] Using charge balance adjust
- [9:09] Dispensed 0.009501 mL of Base (0.5 M KOH)
- [9:30] Stirrer speed set to 0
- [9:40] Datapoint id 3 collected
- [9:40] Charge balance equation is out by 5.9%
- 9:40 Stirrer speed set to 50
- [9:45] pH 2.48 -> 2.68
- [9:45] Using charge balance adjust
- [9:46] Dispensed 0.006397 mL of Base (0.5 M KOH)
- [10:06] Stirrer speed set to 0
- [10:16] Datapoint id 4 collected
- [10:16] Charge balance equation is out by 10.6%
- [10:16] Stirrer speed set to 50
- [10:21] pH 2.71 -> 2.91
- [10:21] Using charge balance adjust
- [10:21] Dispensed 0.004586 mL of Base (0.5 M KOH)
- [10:41] Stirrer speed set to 0
- [10:51] Datapoint id 5 collected
- [10:51] Charge balance equation is out by -13.3%
- [10:51] Stirrer speed set to 50
- [10:56] pH 2.89 -> 3.09
- [10:56] Using charge balance adjust
- [10:57] Dispensed 0.003669 mL of Base (0.5 M KOH)
- [11:17] Stirrer speed set to 0
- [11:27] Datapoint id 6 collected
- [11:27] Charge balance equation is out by -1.0%
- [11:27] Stirrer speed set to 50
- [11:32] pH 3.09 -> 3.29
- [11:32] Using charge balance adjust
- [11:32] Dispensed 0.002893 mL of Base (0.5 M KOH)
- [11:52] Stirrer speed set to 0
- [12:04] Datapoint id 7 collected
- [12:04] Charge balance equation is out by 7.5%
- [12:04] Stirrer speed set to 50
- [12:09] pH 3.31 -> 3.51
- [12:09] Using charge balance adjust
- [12:10] Dispensed 0.002187 mL of Base (0.5 M KOH)
- [12:30] Stirrer speed set to 0
- [12:46] Datapoint id 8 collected
- [12:46] Charge balance equation is out by 5.8%
- [12:46] Stirrer speed set to 50
- [12:51] pH 3.52 -> 3.72
- [12:51] Using charge balance adjust
- [12:51] Dispensed 0.001599 mL of Base (0.5 M KOH)
- [13:11] Stirrer speed set to 0
- [13:32] Datapoint id 9 collected
- [13:32] Charge balance equation is out by 17.6%
- [13:32] Stirrer speed set to 50
- [13:37] pH 3.76 -> 3.96
- [13:37] Using cautious pH adjust
- [13:37] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [13:42] Stepping pH = 3.90
- [13:42] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [13:47] Stepping pH = 3.93
- [13:47] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [13:52] Stepping pH = 3.95



Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

- [14:07] Stirrer speed set to 0
- [14:41] Datapoint id 10 collected
- [14:41] Charge balance equation is out by 24.5%
- [14:41] Stirrer speed set to 50
- [14:46] pH 3.92 -> 4.12
- [14:46] Using cautious pH adjust
- [14:46] Dispensed 0.000400 mL of Base (0.5 M KOH)
- [14:51] Stepping pH = 4.14
- [15:06] Stirrer speed set to 0
- [15:40] Datapoint id 11 collected
- [15:40] Charge balance equation is out by 50.0%
- [15:40] Stirrer speed set to 50
- [15:45] pH 4.08 -> 4.28
- [15:45] Using cautious pH adjust
- [15:45] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [15:51] Stepping pH = 4.29
- [16:06] Stirrer speed set to 0
- [16:43] Datapoint id 12 collected
- [16:43] Charge balance equation is out by 50.0%
- [16:43] Stirrer speed set to 50
- [16:48] pH 4.23 -> 4.43
- [16:48] Using cautious pH adjust
- [16:48] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [16:53] Stepping pH = 4.38
- [16:54] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [16:59] Stepping pH = 4.39
- [16:59] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [17:04] Stepping pH = 4.48
- [17:19] Stirrer speed set to 0
- [17:59] Datapoint id 13 collected
- [17:59] Charge balance equation is out by 17.5%
- [17:59] Stirrer speed set to 50
- [18:04] pH 4.49 -> 4.69
- [18:04] Using cautious pH adjust
- [18:04] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [18:09] Stepping pH = 4.58
- [18:09] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [18:14] Stepping pH = 4.81
- [18:29] Stirrer speed set to 0
- [19:13] Datapoint id 14 collected
- [19:13] Charge balance equation is out by 11.3%
- [19:13] Stirrer speed set to 50
- [19:18] pH 4.89 -> 5.09
- [19:18] Using charge balance adjust
- [19:18] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [19:38] Stirrer speed set to 0
- [20:25] Datapoint id 15 collected
- [20:25] Charge balance equation is out by -19.7%
- [20:25] Stirrer speed set to 50
- [20:30] pH 5.10 -> 5.30
- [20:30] Using cautious pH adjust
- [20:30] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [20:35] Stepping pH = 5.11
- [20:35] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [20:40] Stepping pH = 5.54
- [20:56] Stirrer speed set to 0
- [21:45] Datapoint id 16 collected
- [21:45] Charge balance equation is out by -87.8%
- [21:45] Stirrer speed set to 50



Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

- [21:50] pH 6.10 -> 6.30
- [21:50] Using cautious pH adjust
- [21:50] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [21:55] Stepping pH = 6.11
- [21:55] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [22:01] Stepping pH = 6.60
- [22:16] Stirrer speed set to 0
- [23:13] Datapoint id 17 collected
- [23:13] Charge balance equation is out by -96.1%
- [23:13] Stirrer speed set to 50
- [23:18] pH 7.61 -> 7.81
- [23:18] Using cautious pH adjust
- [23:18] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [23:23] Stepping pH = 7.64
- [23:23] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [23:28] Stepping pH = 7.64
- [23:28] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [23:33] Stepping pH = 7.71
- [23:33] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [23:38] Stepping pH = 7.95
- [23:53] Stirrer speed set to 0
- [24:46] Datapoint id 18 collected
- [24:46] Charge balance equation is out by -1,094.0%
- [24:46] Stirrer speed set to 50
- [24:51] pH 8.41 -> 8.61 [24:51] Using cautious pH adjust
- [24:51] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [24:56] Stepping pH = 8.43
- [24:56] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [25:01] Stepping pH = 8.43
- [25:01] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [25:06] Stepping pH = 8.72
- [25:21] Stirrer speed set to 0
- [25:57] Datapoint id 19 collected [25:57] Charge balance equation is out by -649.6%
- [25:57] Stirrer speed set to 50
- [26:02] pH 8.94 -> 9.05
- [26:02] Using cautious pH adjust
- [26:03] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [26:08] Stepping pH = 8.95
- [26:08] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [26:13] Stepping pH = 8.97
- [26:13] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [26:18] Stepping pH = 9.10
- [26:33] Stirrer speed set to 0
- [27:03] Datapoint id 20 collected
- [27:03] Charge balance equation is out by -476.6%
- [27:03] Titration 2 of 3
- [27:03] Adding initial titrants
- [27:03] Automatically add 0.03000 mL of Octanol
- [27:04] Dispensed 0.030009 mL of Octanol
- [27:04] Stirrer speed set to 10
- [27:05] Stirrer speed set to 55
- [27:05] Iterative adjust 9.21 -> 2.00
- [27:05] pH 9.21 -> 2.00
- [27:06] Dispensed 0.057126 mL of Acid (0.5 M HCI)
- [27:57] Stirrer speed set to 0
- [28:07] Datapoint id 21 collected
- [28:07] Stirrer speed set to 55
- Reported at: 3/9/2018 11:30:33 AM



Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

- [28:12] pH 1.97 -> 2.17
- [28:12] Using cautious pH adjust
- [28:12] Dispensed 0.009666 mL of Base (0.5 M KOH)
- [28:17] Stepping pH = 2.05
- [28:18] Dispensed 0.007361 mL of Base (0.5 M KOH)
- [28:23] Stepping pH = 2.14
- [28:23] Dispensed 0.001623 mL of Base (0.5 M KOH)
- [28:28] Stepping pH = 2.17
- [28:43] Stirrer speed set to 0
- [28:53] Datapoint id 22 collected
- [28:53] Charge balance equation is out by 3.5%
- [28:53] Stirrer speed set to 55
- [28:58] pH 2.18 -> 2.38
- [28:58] Using charge balance adjust
- [28:59] Dispensed 0.012300 mL of Base (0.5 M KOH)
- [29:19] Stirrer speed set to 0
- [29:29] Datapoint id 23 collected
- [29:29] Charge balance equation is out by 2.4%
- [29:29] Stirrer speed set to 55
- [29:34] pH 2.40 -> 2.60
- [29:34] Using charge balance adjust
- [29:34] Dispensed 0.008255 mL of Base (0.5 M KOH)
- [29:54] Stirrer speed set to 0
- [30:05] Datapoint id 24 collected
- [30:05] Charge balance equation is out by 12.9%
- [30:05] Stirrer speed set to 55
- [30:10] pH 2.63 -> 2.83
- [30:10] Using charge balance adjust
- [30:10] Dispensed 0.005597 mL of Base (0.5 M KOH)
- [30:31] Stirrer speed set to 0
- [30:41] Datapoint id 25 collected
- [30:41] Charge balance equation is out by 0.3%
- [30:41] Stirrer speed set to 55
- [30:46] pH 2.83 -> 3.03
- [30:46] Using charge balance adjust
- [30:46] Dispensed 0.004069 mL of Base (0.5 M KOH)
- [31:06] Stirrer speed set to 0
- [31:16] Datapoint id 26 collected
- [31:16] Charge balance equation is out by 8.9%
- [31:16] Stirrer speed set to 55
- [31:22] pH 3.06 -> 3.26
- [31:22] Using charge balance adjust
- [31:22] Dispensed 0.002869 mL of Base (0.5 M KOH)
- [31:42] Stirrer speed set to 0
- [31:52] Datapoint id 27 collected
- [31:52] Charge balance equation is out by 0.5%
- [31:52] Stirrer speed set to 55
- [31:57] pH 3.26 -> 3.46
- [31:57] Using charge balance adjust
- [31:58] Dispensed 0.002023 mL of Base (0.5 M KOH)
- [32:18] Stirrer speed set to 0
- [32:28] Datapoint id 28 collected
- [32:28] Charge balance equation is out by -2.6%
- [32:28] Stirrer speed set to 55
- [32:33] pH 3.47 -> 3.67
- [32:33] Using charge balance adjust
- [32:33] Dispensed 0.001411 mL of Base (0.5 M KOH)
- [32:53] Stirrer speed set to 0
- [33:04] Datapoint id 29 collected



Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

- [33:04] Charge balance equation is out by 27.7%
- [33:04] Stirrer speed set to 55
- [33:09] pH 3.73 -> 3.93
- [33:09] Using cautious pH adjust
- [33:09] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [33:14] Stepping pH = 3.84
- [33:14] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [33:19] Stepping pH = 3.92
- [33:19] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [33:24] Stepping pH = 3.92
- [33:24] Dispensed 0.000071 mL of Base (0.5 M KOH)
- 00.24] Disperised 0.00007
- [33:30] Stepping pH = 3.94
- [33:45] Stirrer speed set to 0
- 55.45] Stiller speed set to 0
- [33:55] Datapoint id 30 collected
- [33:55] Charge balance equation is out by 7.6%
- [33:55] Stirrer speed set to 55
- [34:00] pH 3.95 -> 4.15
- [34:00] Using charge balance adjust
- [34:00] Dispensed 0.000517 mL of Base (0.5 M KOH)
- [34:21] Stirrer speed set to 0
- [34:39] Datapoint id 31 collected
- [34:39] Charge balance equation is out by 139.7%
- [34:39] Stirrer speed set to 55
- [34:44] pH 4.47 -> 4.67
- [34:44] Using cautious pH adjust
- [34:44] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [34:49] Stepping pH = 4.47
- [34:50] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [34:55] Stepping pH = 7.49
- [35:10] Stirrer speed set to 0
- [36:10] Datapoint id 32 collected
- [36:10] Charge balance equation is out by -201.6%
- [36:10] Stirrer speed set to 55
- [36:15] pH 8.18 -> 8.38
- [36:15] Using cautious pH adjust
- [36:15] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [36:20] Stepping pH = 8.14
- [36:20] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [36:25] Stepping pH = 8.10
- [36:25] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [36:30] Stepping pH = 8.92
- [36:45] Stirrer speed set to 0
- [37:17] Datapoint id 33 collected
- [37:17] Charge balance equation is out by -2,089.2%
- [37:17] Stirrer speed set to 55
- [37:22] pH 9.03 -> 9.05
- [37:22] Using cautious pH adjust
- [37:37] Stirrer speed set to 0
- [38:10] Datapoint id 34 collected
- [38:10] Charge balance equation is out by 100.0%
- [38:10] Titration 3 of 3
- [38:10] Adding initial titrants
- [38:10] Automatically add 0.25000 mL of Octanol
- [38:16] Dispensed 0.250000 mL of Octanol
- [38:16] Stirrer speed set to 10
- [38:17] Stirrer speed set to 60
- [38:17] Iterative adjust 9.00 -> 2.00
- [38:17] pH 9.00 -> 2.00
- [38:19] Dispensed 0.058725 mL of Acid (0.5 M HCI)



Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

- [38:24] pH 2.01 -> 2.00
- [38:24] Dispensed 0.001294 mL of Acid (0.5 M HCl)
- [39:14] Stirrer speed set to 0
- [39:25] Datapoint id 35 collected
- [39:25] Stirrer speed set to 60
- [39:30] pH 1.98 -> 2.18
- [39:30] Using cautious pH adjust
- [39:30] Dispensed 0.010254 mL of Base (0.5 M KOH)
- [39:35] Stepping pH = 2.06
- [39:36] Dispensed 0.008043 mL of Base (0.5 M KOH)
- [39:41] Stepping pH = 2.15
- [39:41] Dispensed 0.002493 mL of Base (0.5 M KOH)
- [39:46] Stepping pH = 2.18
- [40:01] Stirrer speed set to 0
- [40:01] Stiffer speed set to 0 [40:12] Datapoint id 36 collected
- [40:12] Charge balance equation is out by -1.3%
- [40:12] Stirrer speed set to 60
- [40:17] pH 2.18 -> 2.38
- [40:17] Using charge balance adjust
- [40:18] Dispensed 0.013594 mL of Base (0.5 M KOH)
- [40:38] Stirrer speed set to 0
- [40:48] Datapoint id 37 collected
- [40:48] Charge balance equation is out by -1.2%
- [40:48] Stirrer speed set to 60
- [40:53] pH 2.38 -> 2.58
- [40:53] Using charge balance adjust
- [40:53] Dispensed 0.009196 mL of Base (0.5 M KOH)
- [41:14] Stirrer speed set to 0
- [41:24] Datapoint id 38 collected
- [41:24] Charge balance equation is out by 9.4%
- [41:24] Stirrer speed set to 60
- [41:29] pH 2.60 -> 2.80
- [41:29] Using charge balance adjust
- [41:30] Dispensed 0.006044 mL of Base (0.5 M KOH)
- [41:50] Stirrer speed set to 0
- [42:00] Datapoint id 39 collected
- [42:00] Charge balance equation is out by 16.8%
- [42:00] Stirrer speed set to 60
- [42:05] pH 2.85 -> 3.05
- [42:05] Using cautious pH adjust
- [42:05] Dispensed 0.001929 mL of Base (0.5 M KOH)
- [42:10] Stepping pH = 2.94
- [42:10] Dispensed 0.001435 mL of Base (0.5 M KOH)
- [42:15] Stepping pH = 3.03
- [42:15] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [42:20] Stepping pH = 3.04
- [42:36] Stirrer speed set to 0
- [42:46] Datapoint id 40 collected
- [42:46] Charge balance equation is out by 5.9%
- [42:46] Stirrer speed set to 60
- [42:51] pH 3.05 -> 3.25
- [42:51] Using charge balance adjust
- [42:51] Dispensed 0.002611 mL of Base (0.5 M KOH)
- [43:12] Stirrer speed set to 0
- [43:26] Datapoint id 41 collected
- [43:26] Charge balance equation is out by 35.1%
- [43:26] Stirrer speed set to 60
- [43:31] pH 3.32 -> 3.52
- [43:31] Using cautious pH adjust



Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

- [43:31] Dispensed 0.000729 mL of Base (0.5 M KOH)
- [43:36] Stepping pH = 3.43
- [43:37] Dispensed 0.000400 mL of Base (0.5 M KOH)
- [43:42] Stepping pH = 3.52
- [43:57] Stirrer speed set to 0
- [44:07] Datapoint id 42 collected
- [44:07] Charge balance equation is out by 22.8%
- [44:07] Stirrer speed set to 60
- [44:12] pH 3.53 -> 3.73
- [44:12] Using cautious pH adjust
- [44:12] Dispensed 0.000470 mL of Base (0.5 M KOH)
- [44:18] Stepping pH = 3.63
- [44:18] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [44:23] Stepping pH = 3.71
- [44:23] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [44:28] Stepping pH = 3.72
- [44:28] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [44:33] Stepping pH = 3.73
- [44:48] Stirrer speed set to 0
- [44:58] Datapoint id 43 collected
- [44:58] Charge balance equation is out by 5.5%
- [44:58] Stirrer speed set to 60
- [45:03] pH 3.75 -> 3.95
- [45:03] Using charge balance adjust
- [45:03] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [45:24] Stirrer speed set to 0
- [45:34] Datapoint id 44 collected
- [45:34] Charge balance equation is out by 26.0%
- [45:34] Stirrer speed set to 60
- [45:39] pH 4.03 -> 4.23
- [45:39] Using cautious pH adjust
- [45:39] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [45:44] Stepping pH = 4.07
- [45:44] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [45:49] Stepping pH = 4.31
- [46:04] Stirrer speed set to 0
- [46:15] Datapoint id 45 collected
- [46:15] Charge balance equation is out by -34.4%
- [46:15] Stirrer speed set to 60
- [46:20] pH 4.37 -> 4.57
- [46:20] Using cautious pH adjust
- [46:20] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [46:25] Stepping pH = 4.37
- [46:25] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [46:30] Stepping pH = 4.68 [46:45] Stirrer speed set to 0
- [47:14] Datapoint id 46 collected
- [47:14] Charge balance equation is out by -97.4%
- [47:14] Stirrer speed set to 60
- [47:19] pH 4.84 -> 5.04
- [47:19] Using cautious pH adjust
- [47:19] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [47:25] Stepping pH = 4.84
- [47:25] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [47:30] Stepping pH = 4.86
- [47:30] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [47:35] Stepping pH = 8.18
- [47:50] Stirrer speed set to 0
- [48:50] Datapoint id 47 collected
- Reported at: 3/9/2018 11:30:33 AM

## **Experiment Log**



Sample name: M02\_octanol Experiment start time: 3/6/2018 5:41:59 PM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

Assay ID: 18C-06015 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180306\_exp30\_logP\_T3-2\18C-06015\_M02\_octanol\_pH-metric high logP.t3r

## Experiment Log (continued)

[48:50] Charge balance equation is out by -866.5%

[48:50] Stirrer speed set to 60

[48:55] pH 8.41 -> 8.61

[48:55] Using cautious pH adjust

[48:55] Dispensed 0.000024 mL of Base (0.5 M KOH)

[49:00] Stepping pH = 8.38

[49:00] Dispensed 0.000047 mL of Base (0.5 M KOH)

[49:05] Stepping pH = 8.35

[49:06] Dispensed 0.000235 mL of Base (0.5 M KOH)

[49:11] Stepping pH = 8.46

[49:11] Dispensed 0.000235 mL of Base (0.5 M KOH)

[49:16] Stepping pH = 9.22

[49:31] Stirrer speed set to 0

[49:46] Datapoint id 48 collected

[49:46] Charge balance equation is out by -3,486.6%

[49:46] Argon flow rate set to 0

[49:50] Titrator arm moved over Titration position

[50:12] The autoloader failed to pick at location "Sample position"