

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH2 +) -0.24 ±0.09 (n=50) logP (neutral XH) 2.63 ±0.01 (n=50) logP (X -) 0.46

RMSD 0.266

18C-09014 Points 2 to 24

M18_octanol concentration factor 0.988
Carbonate 0.1577 mM
Acidity error -0.03698 mM

18C-09014 Points 25 to 48

M18_octanol concentration factor 0.944
Carbonate 0.1161 mM
Acidity error -0.54687 mM

18C-09014 Points 49 to 73

M18_octanol concentration factor 0.993
Carbonate 0.1503 mM
Acidity error -1.10219 mM

Warnings and errors

Errors None Warnings None

Sample logD and percent species

рН	M18_octanol logD	M18_octanol M18_octanolH2					M18_octanol*	
1.000	-0.23	63.00 %	0.00 %	0.00 %	35.85 %	1.14 %	0.00 %	
1.200	-0.22	62.58 %	0.00 %	0.00 %	35.61 %	1.80 %	0.00 %	Stomach pH
2.000	-0.12	57.11 %	0.02 %	0.00 %	32.50 %	10.37 %	0.00 %	-
3.000	0.38	29.51 %	0.13 %	0.00 %	16.79 %	53.57 %	0.00 %	
4.000	1.25	5.06 %	0.22 %	0.00 %	2.88 %	91.85 %	0.00 %	
5.000	2.11	0.54 %	0.23 %	0.00 %	0.31 %	98.91 %	0.00 %	
6.000	2.54	0.05 %	0.23 %	0.00 %	0.03 %	99.68 %	0.00 %	
6.500	2.60	0.02 %	0.23 %	0.00 %	0.01 %	99.74 %	0.00 %	
7.000	2.62	0.01 %	0.23 %	0.00 %	0.00 %	99.76 %	0.00 %	
7.400	2.62	0.00 %	0.23 %	0.00 %	0.00 %	99.76 %	0.00 %	Blood pH
8.000	2.63	0.00 %	0.23 %	0.00 %	0.00 %	99.76 %	0.00 %	·
9.000	2.62	0.00 %	0.23 %	0.01 %	0.00 %	99.75 %	0.02 %	
10.000	2.54	0.00 %	0.23 %	0.05 %	0.00 %	99.56 %	0.15 %	
11.000	2.13	0.00 %	0.23 %	0.51 %	0.00 %	97.77 %	1.48 %	
12.000	1.32	0.00 %	0.19 %	4.36 %	0.00 %	82.87 %	12.57 %	



Sample name: M18_octanol Assay name:


pH-metric high logP

18C-09014 Assay ID:

Experiment start time: 3/9/2018 8:37:21 PM

Pion Analyst: Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r



Log R


pH (Concentration scale)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

Graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-09014 Points 2 to 24

Overall results

RMSD 0.227
Average ionic strength 0.157 M
Average temperature 24.9°C
Partition ratio 0.0624 : 1

Analyte concentration range 3258.0 µM to 3348.4 µM

Total points considered 14 of 23

Warnings and errors

Errors None Warnings None

Four-Plus parameters

 Alpha
 0.102
 3/9/2018 8:37:21 PM C:\Sirius_T3\HCl18C09.t3r

 S
 0.9967
 3/9/2018 8:37:21 PM C:\Sirius_T3\HCl18C09.t3r

 jH
 1.2
 3/9/2018 8:37:21 PM C:\Sirius_T3\HCl18C09.t3r


 jOH
 0.0
 3/9/2018 8:37:21 PM C:\Sirius_T3\HCl18C09.t3r


Titrants

Sample

M18_octanol concentration factor 0.988
Base pKa 1 5.37
Acid pKa 2 10.65
logP (XH2 +) -1.08
logP (neutral XH) 2.57
logP (X -) 0.46

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M18_octanol	M18_octanol	M18_octanol	M18_octanol		M18_octanol	M18_octanol	Comment
	logD	M18_octanolH2	M18_octanolH	M18_octanol	M18_octanolH2*	M18_octanolH*	M18_octanol*	
1.000	-1.01	99.39 %	0.00 %	0.00 %	0.51 %	0.10 %	0.00 %	
1.200	-0.97	99.33 %	0.01 %	0.00 %	0.51 %	0.16 %	0.00 %	Stomach pH
2.000	-0.62	98.48 %	0.04 %	0.00 %	0.51 %	0.97 %	0.00 %	
3.000	0.22	90.24 %	0.38 %	0.00 %	0.46 %	8.91 %	0.00 %	
4.000	1.18	49.14 %	2.10 %	0.00 %	0.25 %	48.51 %	0.00 %	
5.000	2.05	8.85 %	3.77 %	0.00 %	0.05 %	87.33 %	0.00 %	
6.000	2.48	0.96 %	4.10 %	0.00 %	0.00 %	94.93 %	0.00 %	
6.500	2.54	0.31 %	4.13 %	0.00 %	0.00 %	95.56 %	0.00 %	
7.000	2.56	0.10 %	4.14 %	0.00 %	0.00 %	95.76 %	0.00 %	
7.400	2.56	0.04 %	4.14 %	0.00 %	0.00 %	95.82 %	0.00 %	Blood pH
8.000	2.57	0.01 %	4.14 %	0.01 %	0.00 %	95.84 %	0.00 %	- 1
9.000	2.56	0.00 %	4.14 %	0.09 %	0.00 %	95.75 %	0.02 %	ļ
10.000	2.48	0.00 %	4.10 %	0.92 %	0.00 %	94.82 %	0.17 %	ļ
11.000	2.07	0.00 %	3.73 %	8.36 %	0.00 %	86.40 %	1.50 %	ļ
12.000	1.27	0.00 %	1.98 %	44.28 %	0.00 %	45.77 %	7.97 %	ļ

Carbonate and acidity



Other graphs





Assay ID: Filename:

Sample name: M18_octanol Assay name:

pH-metric high logP


18C-09014

Experiment start time: 3/9/2018 8:37:21 PM

Analyst: Pion Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-09014 Points 25 to 48

Overall results

RMSD 0.131
Average ionic strength 0.163 M
Average temperature 25.0°C
Partition ratio 0.1756 : 1

Analyte concentration range 2756.3 µM to 2832.3 µM

Total points considered 16 of 24

Warnings and errors

Errors None Warnings None

Four-Plus parameters

 Alpha
 0.102
 3/9/2018 8:37:21 PM C:\Sirius_T3\HCl18C09.t3r

 S
 0.9967
 3/9/2018 8:37:21 PM C:\Sirius_T3\HCl18C09.t3r

 jH
 1.2
 3/9/2018 8:37:21 PM C:\Sirius_T3\HCl18C09.t3r


 jOH
 0.0
 3/9/2018 8:37:21 PM C:\Sirius_T3\HCl18C09.t3r


Titrants

Sample

M18_octanol concentration factor 0.944
Base pKa 1 5.37
Acid pKa 2 10.65
logP (XH2 +) -1.92
logP (neutral XH) 2.67
logP (X -) 0.46

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M18_octanol	M18_octanol	M18_octanol	M18_octanol	M18_octanol	M18_octanol	M18_octanol	
	logD	M18_octanolH2	—	M18_octanol	M18_octanolH2*	_	M18_octanol*	
1.000	-1.50	99.44 %	0.00 %	0.00 %	0.21 %	0.35 %	0.00 %	
1.200	-1.36	99.23 %	0.01 %	0.00 %	0.21 %	0.55 %	0.00 %	Stomach pH
2.000	-0.68	96.38 %	0.04 %	0.00 %	0.20 %	3.37 %	0.00 %	
3.000	0.30	73.74 %	0.31 %	0.00 %	0.16 %	25.79 %	0.00 %	
4.000	1.28	22.02 %	0.94 %	0.00 %	0.05 %	76.99 %	0.00 %	
5.000	2.14	2.75 %	1.17 %	0.00 %	0.01 %	96.07 %	0.00 %	
6.000	2.58	0.28 %	1.20 %	0.00 %	0.00 %	98.52 %	0.00 %	
6.500	2.64	0.09 %	1.20 %	0.00 %	0.00 %	98.71 %	0.00 %	
7.000	2.66	0.03 %	1.20 %	0.00 %	0.00 %	98.77 %	0.00 %	
7.400	2.66	0.01 %	1.21 %	0.00 %	0.00 %	98.78 %	0.00 %	Blood pH
8.000	2.67	0.00 %	1.21 %	0.00 %	0.00 %	98.79 %	0.00 %	·
9.000	2.66	0.00 %	1.20 %	0.03 %	0.00 %	98.75 %	0.01 %	
10.000	2.58	0.00 %	1.20 %	0.27 %	0.00 %	98.39 %	0.14 %	
11.000	2.16	0.00 %	1.16 %	2.59 %	0.00 %	94.94 %	1.31 %	
12.000	1.36	0.00 %	0.86 %	19.18 %	0.00 %	70.24 %	9.72 %	

Carbonate and acidity



Other graphs





Assay ID: Filename:

Sample name: M18_octanol Assay name:

pH-metric high logP


18C-09014

Experiment start time: 3/9/2018 8:37:21 PM

Analyst: **Pion** Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-09014 Points 49 to 73

Overall results

RMSD 0.342
Average ionic strength 0.169 M
Average temperature 25.0°C
Partition ratio 0.4375 : 1

Analyte concentration range 2105.7 µM to 2154.4 µM

Total points considered 22 of 25

Warnings and errors

Errors None

Warnings One or more logP values out of range

Excessive acidity error present

Four-Plus parameters


Alpha	0.102	3/9/2018 8:37:21 PM	C:\Sirius_T3\HCl18C09.t3r
S	0.9967	3/9/2018 8:37:21 PM	C:\Sirius_T3\HCl18C09.t3r
jΗ	1.2	3/9/2018 8:37:21 PM	C:\Sirius_T3\HCl18C09.t3r
jOH	0.0	3/9/2018 8:37:21 PM	C:\Sirius_T3\HCl18C09.t3r


Titrants

Sample

>	M18 octanol concentration factor	0.993
	Base pKa 1	5.37
	Acid pKa 2	10.65
>	logP (XH2 +)	-2.03
>	logP (neutral XH)	2.50
ı	logP (X -)	0.46

Sample graphs








Assay ID:

Sample name: M18_octanol Experiment start time: 3/9/2018 8:37:21 PM
Assay name: pH-metric high logP Analyst: Pion

pH-metric high logP Analyst: Pion 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M18_octanol	M18_octanol	M18_octanol	M18_octanol	M18_octanol	M18_octanol	M18_octanol	
	logD	M18_octanolH2	—	M18_octanol	M18_octanolH2*	_	M18_octanol*	
1.000	-1.64	99.00 %	0.00 %	0.00 %	0.40 %	0.59 %	0.00 %	
1.200	-1.51	98.66 %	0.01 %	0.00 %	0.40 %	0.93 %	0.00 %	Stomach pH
2.000	-0.84	93.97 %	0.04 %	0.00 %	0.38 %	5.61 %	0.00 %	
3.000	0.14	62.30 %	0.27 %	0.00 %	0.25 %	37.18 %	0.00 %	
4.000	1.12	14.25 %	0.61 %	0.00 %	0.06 %	85.08 %	0.00 %	
5.000	1.98	1.64 %	0.70 %	0.00 %	0.01 %	97.66 %	0.00 %	
6.000	2.41	0.17 %	0.71 %	0.00 %	0.00 %	99.12 %	0.00 %	
6.500	2.47	0.05 %	0.71 %	0.00 %	0.00 %	99.24 %	0.00 %	
7.000	2.49	0.02 %	0.71 %	0.00 %	0.00 %	99.27 %	0.00 %	
7.400	2.50	0.01 %	0.71 %	0.00 %	0.00 %	99.28 %	0.00 %	Blood pH
8.000	2.50	0.00 %	0.71 %	0.00 %	0.00 %	99.29 %	0.00 %	·
9.000	2.50	0.00 %	0.71 %	0.02 %	0.00 %	99.25 %	0.02 %	
10.000	2.42	0.00 %	0.71 %	0.16 %	0.00 %	98.93 %	0.20 %	
11.000	2.00	0.00 %	0.69 %	1.53 %	0.00 %	95.85 %	1.93 %	
12.000	1.22	0.00 %	0.52 %	11.69 %	0.00 %	73.05 %	14.75 %	

Carbonate and acidity



Other graphs





Assay ID: Filename:

Sample name: M18_octanol Assay name:

pH-metric high logP


18C-09014

Experiment start time: 3/9/2018 8:37:21 PM

Analyst: **Pion** Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name	M18_octanol	2/27/2018 7:08:39 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.001480 g	3/9/2018 2:22:36 PM	User entered value
Formula weight	267.11 g/mol	2/27/2018 7:08:39 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	267.11	2/27/2018 7:08:39 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	2	2/27/2018 7:08:39 PM	User entered value
Sample is a	Ampholyte	2/27/2018 7:08:39 PM	User entered value
pKa 1	5.37	2/27/2018 7:08:39 PM	User entered value
Туре	Base	2/27/2018 7:08:39 PM	User entered value
pKa 2	10.65	2/27/2018 7:08:39 PM	User entered value
Туре	Acid	2/27/2018 7:08:39 PM	User entered value
logp (XH2 +)	-0.35	2/28/2018 3:20:28 PM	User entered value
logP (neutral XH)	2.57	3/2/2018 4:34:50 PM	User entered value
logP (X -)	0.46	2/27/2018 7:09:34 PM	User entered value

Events

⊏vents										
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/c
4:57.9 4:59.0	Manual volume addition Initial pH = 7.17				0.10000 mL					tille
7:58.6	Data point 2	1.50000 mL	0.05002 mL	0.00475 mL	0.10000 mL	2.041	-0.00537	0.75188	0.00031	10.5 s
8:45.2	Data point 3	1.50000 mL	0.05002 mL	0.01719 mL	0.10000 mL	2.248	-0.00630	0.13591	0.00084	-
9:20.9	Data point 4	1.50000 mL	0.05002 mL	0.02556 mL	0.10000 mL	2.455	-0.00491	0.10647	0.00074	-
9:56.5	Data point 5	1.50000 mL	0.05002 mL	0.03079 mL	0.10000 mL	2.645	-0.00008	0.00076	0.00014	-
10:31.9	Data point 6	1.50000 mL	0.05002 mL	0.03429 mL	0.10000 mL	2.849	-0.00867	0.76331	0.00049	-
11:07.9	Data point 7	1.50000 mL	0.05002 mL	0.03667 mL	0.10000 mL	3.088	-0.00680	0.26943	0.00065	-
11:58.8	Data point 8	1.50000 mL	0.05002 mL	0.03845 mL	0.10000 mL	3.287	-0.00858	0.64011	0.00053	-
12:34.8	Data point 9	1.50000 mL	0.05002 mL	0.03991 mL	0.10000 mL	3.481	-0.00517	0.71303	0.00030	•
13:10.2	Data point 10	1.50000 mL	0.05002 mL	0.04132 mL	0.10000 mL	3.675	-0.01100	0.68415	0.00066	-
13:45.7	Data point 11	1.50000 mL	0.05002 mL	0.04276 mL	0.10000 mL	3.880	-0.01228	0.70505	0.00072	10.5 s
14:21.7	Data point 12	1.50000 mL	0.05002 mL	0.04421 mL	0.10000 mL	4.108	-0.01324	0.78765	0.00074	•
14:57.1	Data point 13	1.50000 mL	0.05002 mL	0.04558 mL	0.10000 mL	4.320	-0.00760	0.32722	0.00066	•
15:32.6	Data point 14	1.50000 mL	0.05002 mL	0.04675 mL	0.10000 mL	4.509	-0.01875	0.88887	0.00098	10.5 s
16:08.5	Data point 15	1.50000 mL	0.05002 mL	0.04770 mL	0.10000 mL	4.638	-0.09531	0.99824	0.00471	Timed out
17:39.0	Data point 16	1.50000 mL	0.05002 mL	0.04857 mL	0.10000 mL	4.731	-0.03485	0.99545	0.00173	at Timed out at
19:04.3	Data point 17	1.50000 mL	0.05002 mL	0.04934 mL	0.10000 mL	4.989	-0.01963	0.98594	0.00098	40.5 s



Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-09014 Instrument ID: T312060

Filename:

 $C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric\ high\ logP.t3r$

Events ((continued)
_,0	Jointinada

FACIII2 (Continueu)									
Time 20:15.5	Event Data point 18	Water 1.50000 mL	Acid 0.05002 mL	Base 0.05012 mL	Octanol 0.10000 mL	pH 6.544	dpH/dt -0.03976	pH R-squared 0.99389		dpH/dt time Timed out at
21:46.0	Data point 19	1.50000 mL	0.05002 mL	0.05028 mL	0.10000 mL	6.855	-0.03344	0.99646	0.00165	59.5 s Timed out at 59.5 s
23:16.5	Data point 20	1.50000 mL	0.05002 mL	0.05038 mL	0.10000 mL	7.476	-0.05450	0.99713	0.00270	Timed out at 59.5 s
24:47.0	Data point 21	1.50000 mL	0.05002 mL	0.05042 mL	0.10000 mL	8.023	-0.06285	0.99743	0.00311	Timed out at 59.5 s
26:17.5	Data point 22								0.00206	Timed out at 59.5 s
27:53.2	Data point 23								0.00099	
29:14.8	Data point 24								0.00099	
30:56.2	Data point 25								0.00045	
31:42.4	Data point 26								0.00052	
32:18.0	Data point 27								0.00066	
32:53.6	Data point 28							0.00001	0.00050	
33:39.6	Data point 29								0.00025	
34:15.1	Data point 30								0.00065	
34:51.1	Data point 31								0.00020	
35:26.5	Data point 32								0.00049	
36:01.9	Data point 33								0.00090	
36:37.4 37:18.0	Data point 34								0.00076	
38:08.9	Data point 35 Data point 36								0.00051 0.00029	
38:44.4	Data point 37								0.00029	
39:30.1	Data point 38								0.00040	
40:06.0	Data point 39								0.00079	
40:47.1	Data point 40								0.00073	
41:24.5	Data point 41								0.00086	
42:06.6	Data point 42								0.00096	
42:52.8	Data point 43									Timed out at 59.5 s
44:28.4	Data point 44								0.00326	Timed out at 59.5 s
46:04.1	Data point 45									Timed out at 59.5 s
47:39.9	Data point 46								0.00098	
49:12.6	Data point 47								0.00100	
50:23.1	Data point 48								0.00092	
52:23.2	Data point 49								0.00069	
53:10.0 53:45.8	Data point 50 Data point 51							0.50246 0.76262	0.00026 0.00083	
53.45.6 54:21.4	Data point 52								0.00085	
54:56.9	Data point 52								0.00043	
55:32.9	Data point 54								0.00072	
56:09.0	Data point 55								0.00091	
56:44.4	Data point 56								0.00030	
57:20.0	Data point 57								0.00063	
57:55.8	Data point 58								0.00044	
58:31.3	Data point 59							0.02519	0.00090	
59:07.2	Data point 60								0.00038	
59:43.0	Data point 61	1.50000 mL	0.16729 mL	0.16223 mL	0.80000 mL	4.381	-0.00427	0.08186	0.00074	
	Data point 62								0.00075	11.0 s
	Data point 63								0.00096	
	Data point 64								0.00097	
1:03:10.9	Data point 65	1.50000 mL	0.16729 mL	0.16345 mL	0.80000 mL	6.233	-0.02014	0.98026	0.00100	Timed out at 59.5 s





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

Events (continued)

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
1:04:41.4	Data point 66	1.50000 mL	0.16729 mL	0.16352 mL	0.80000 mL	6.668	-0.04896	0.99508	0.00242	Timed out at 59.5 s
1:06:12.0	Data point 67	1.50000 mL	0.16729 mL	0.16359 mL	0.80000 mL	7.159	-0.06899	0.99765	0.00341	Timed out at 59.5 s
1:07:42.4	Data point 68	1.50000 mL	0.16729 mL	0.16366 mL	0.80000 mL	7.675	-0.06553	0.99682	0.00324	Timed out at 59.5 s
1:09:18.0	Data point 69	1.50000 mL	0.16729 mL	0.16373 mL	0.80000 mL	8.043	-0.04078	0.98943	0.00202	Timed out at 59.5 s
1:10:48.5	Data point 70	1.50000 mL	0.16729 mL	0.16378 mL	0.80000 mL	8.282	-0.02444	0.97943	0.00122	Timed out at 59.5 s
1:13:30.4	Data point 71 Data point 72 Data point 73 Assay volumes	1.50000 mL 1.50000 mL	0.16729 mL 0.16729 mL 0.16729 mL 0.16729 mL	0.16397 mL 0.16406 mL	0.80000 mL 0.80000 mL	8.923		0.67890 0.25217 0.95656	0.00094 0.00097 0.00098	14.0 s



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	9.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings	' '			
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose	1070			
Titrant pre-dose	None			
Assay Medium	None			
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.100 mL			
Partition solvent added	Manual in advance			
After partition addition, stir for	1 seconds			
Sample Sonication	Vaa			
Sonicate	Yes			
Adjust pH for sonication	No co			
Sonicate for	60 seconds			
After sonication stir for	5 seconds			
Sample Dissolution				
Perform a dissolution stage	Yes			
	To start pH			
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge				
Perform a carbonate purge	No			
Temperature Control				
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			
Stir speed of	50%			
Titration 1				
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	50%			
Titration 2				
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume				
Additional partition solvent added	Automatic			
After pU adjust stir for	20 cocondo			

30 seconds

15 seconds

55%

After pH adjust stir for

Stir to allow partitioning for

Stirrer speed for partitioning



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3		_	_	
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.500 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00100 dpH/dt			
Stability timeout after	60 seconds			

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.102	3/9/2018 8:37:21 PM	C:\Sirius_T3\HCl18C09.t3r
Four-Plus S	0.9967	3/9/2018 8:37:21 PM	C:\Sirius_T3\HCl18C09.t3r
Four-Plus jH	1.2	3/9/2018 8:37:21 PM	C:\Sirius_T3\HCl18C09.t3r
Four-Plus jOH	0.0	3/9/2018 8:37:21 PM	C:\Sirius_T3\HCl18C09.t3r
Base concentration factor	1.000	3/9/2018 8:37:21 PM	C:\Sirius_T3\KOH18B27.t3r
Acid concentration factor	1.000	3/9/2018 8:37:21 PM	C:\Sirius_T3\HCl18C09.t3r

Instrument Settings

Setting	Value	Batch Id	Install date
Instrument owner	Merck		
Instrument ID	T312060		
Instrument type	T3 Simulator		
Software version	1.1.3.0		
Dispenser module		T3DM1200361	3/31/2009 6:24:52 AM
Dispenser 0	Water		3/31/2009 6:25:05 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Water (0.15 M KCI)	02-06-2018	2/27/2018 11:05:59 AM
Dispenser 2	Acid		3/31/2009 6:25:11 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Acid (0.5 M HCI)	02-27-2018	2/27/2018 11:27:22 AM
Dispenser 1	Base		3/31/2009 6:25:21 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Base (0.5 M KOH)	9/22/2017	2/27/2018 11:21:22 AM
Dispenser 5	Cosolvent		3/31/2009 6:26:24 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Distribution valve 5	Distribution Valve		3/31/2009 6:28:19 AM
Firmware version	1.1.3		
Port A	Methanol (80%, 0.15 M KCI)		3/6/2018 10:28:59 AM
Port B	Cyclohexane	11-01-17	2/27/2018 11:37:57 AM
Dispenser 3	Buffer		8/3/2010 6:05:16 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Dodecane	2018/01/31	2/28/2018 11:18:04 AM
Dispenser 6	Octanol		10/22/2010 11:52:43 AM




Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting Syringe volume	Value 0.5 mL	Batch Id	Install date
Firmware version Titrant	1.2.1(r2) Octanol	01-31-2018	2/27/2018 10:59:35 AM
Titrator		T3TM1200161	3/31/2009 6:24:17 AM
Horizontal axis firmware version Vertical axis firmware version Chassis I/O firmware version Probe I/O firmware version	1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2 1.11 Al1Dl0DO4 Norgren I/O 1.1.1		
Electrode E0 calibration	T3 Electrode +5.00 mV	T3E0923	1/23/2018 3:01:00 PM 3/9/2018 8:37:49 PM
Filling solution	3M KCI	KCL097	3/9/2018 11:05:42 AM
Liquids			
Wash 1 Wash 2 Buffer position 1 Buffer position 2	50% IPA:50% Water 0.5% Trition X-100 in H20 pH7 Wash pH 7		3/9/2018 11:04:22 AM 3/9/2018 11:04:25 AM 3/9/2018 11:04:27 AM 3/9/2018 11:04:30 AM
Storage position			3/9/2018 11:05:04 AM
Wash water Waste	5.1e+003 mL 1e+004 mL	02-27-2018	2/27/2018 10:54:39 AM 11/28/2017 11:36:29 AM
Temperature controller Turbidity detector	re+004 IIIL		8/5/2010 7:35:13 AM 3/31/2009 6:24:45 AM
Spectrometer		074811	11/23/2010 12:22:28 PM
Dip probe		10196	17/20/2010 12:22:201 11
Wavelength coefficient A0	183.333		
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		44/00/0040 40 00 00 00
Total lamp lit time	123:16:41		11/23/2010 12:22:28 PM
Calibrated on	2/27/2018 11:40:38 AM		
Integration time Scans averaged	40 10		
Autoloader	10	T3AL1200345	11/10/2015 10:34:13 AM
Left-right axis firmware version Front-back axis firmware version Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2	10/1E1200040	11/10/2010 10:04:10/NW
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Configuration Alternate titration position	Titration position		
Alternate unation position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration Flowing wash stir speed	5 s 30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration E0 calibration preparation stir speed	5 s 30%		
E0 calibration preparation still speed E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09014 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09014_M18_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		

Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00
• •		