

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH2 +) -0.17 ±0.20 (n=50) logP (neutral XH) 3.04 ±0.01 (n=50) logP (X -) -1.62

-g. (/ / / ...-

18C-01009 Points 1 to 31

M16_octanol concentration factor 0.767
Carbonate 0.0840 mM
Acidity error 3.05344 mM

18C-01009 Points 32 to 62

M16_octanol concentration factor 0.794
Carbonate 0.1938 mM
Acidity error 2.68082 mM

18C-01009 Points 63 to 92

M16_octanol concentration factor 0.836
Carbonate 0.2810 mM
Acidity error 2.47693 mM

Warnings and errors

Errors None Warnings None

Sample logD and percent species

рН	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	Comment
	logD	M16_octanolH2	M16_octanolH	M16_octanol	M16_octanolH2*	M16_octanolH*	M16_octanol*	
1.000	-0.05	52.62 %	0.01 %	0.00 %	35.72 %	11.64 %	0.00 %	
1.200	0.01	49.26 %	0.02 %	0.00 %	33.44 %	17.28 %	0.00 %	Stomach pH
2.000	0.46	25.68 %	0.05 %	0.00 %	17.44 %	56.83 %	0.00 %	•
3.000	1.35	4.20 %	0.08 %	0.00 %	2.85 %	92.87 %	0.00 %	
4.000	2.27	0.45 %	0.09 %	0.00 %	0.30 %	99.16 %	0.00 %	
5.000	2.87	0.05 %	0.09 %	0.00 %	0.03 %	99.83 %	0.00 %	
6.000	3.02	0.00 %	0.09 %	0.00 %	0.00 %	99.90 %	0.00 %	
6.500	3.04	0.00 %	0.09 %	0.00 %	0.00 %	99.91 %	0.00 %	
7.000	3.04	0.00 %	0.09 %	0.00 %	0.00 %	99.91 %	0.00 %	
7.400	3.03	0.00 %	0.09 %	0.00 %	0.00 %	99.91 %	0.00 %	Blood pH
8.000	3.00	0.00 %	0.09 %	0.01 %	0.00 %	99.90 %	0.00 %	-
9.000	2.71	0.00 %	0.09 %	0.10 %	0.00 %	99.80 %	0.00 %	
10.000	1.95	0.00 %	0.09 %	1.02 %	0.00 %	98.86 %	0.02 %	
11.000	0.98	0.00 %	0.08 %	9.35 %	0.00 %	90.34 %	0.22 %	
12.000	-0.00	0.00 %	0.04 %	50.23 %	0.00 %	48.52 %	1.20 %	



Assay ID:


Sample name: M16_octanol Assay name:

pH-metric high logP 18C-01009

Experiment start time: 3/1/2018 12:08:09 PM

Analyst: Pion Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r



Log R


pH (Concentration scale)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

Graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-01009 Points 1 to 31

Overall results

RMSD 0.138
Average ionic strength 0.159 M
Average temperature 25.0°C
Partition ratio 0.0247 : 1

Analyte concentration range 4991.2 µM to 5365.1 µM

Total points considered 23 of 31

Warnings and errors

Errors None

Warnings Excessive acidity error present

Four-Plus parameters

1	Alpha	0.130	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
	s ·	0.9970	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
j	Н	0.8	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
j	OH	-0.4	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r


Titrants


0.50 M HCI 0.993513 3/1/2018 12:08:08 PM C:\Sirius_T3\HCl18B27.t3r 0.50 M KOH 0.999845 3/1/2018 12:08:09 PM C:\Sirius_T3\KOH18B27.t3r

Sample

₩	M16_octanol concentration factor	0.767
	Base pKa 1	4.70
	Acid pKa 2	8.94
₩	logP (XH2 +)	0.85
₩	logP (neutral XH)	3.08
۱	logP (X -)	-1.62

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M16_octanol logD	M16_octanol M16_octanolH2	M16_octanol M16_octanolH	M16_octanol M16_octanol		M16_octanol M16_octanolH*	M16_octanol M16_octanol*	
1.000	0.86	84.71 %	0.02 %	0.00 %	14.77 %	0.50 %	0.00 %	
1.200	0.87	84.45 %	0.03 %	0.00 %	14.73 %	0.79 %	0.00 %	Stomach pH
2.000	0.97	80.93 %	0.16 %	0.00 %	14.12 %	4.80 %	0.00 %	P · · ·
3.000	1.48	55.95 %	1.12 %	0.00 %	9.76 %	33.17 %	0.00 %	
4.000	2.31	13.69 %	2.73 %	0.00 %	2.39 %	81.18 %	0.00 %	
5.000	2.90	1.60 %	3.19 %	0.00 %	0.28 %	94.92 %	0.00 %	
6.000	3.06	0.16 %	3.25 %	0.00 %	0.03 %	96.56 %	0.00 %	
6.500	3.07	0.05 %	3.25 %	0.01 %	0.01 %	96.67 %	0.00 %	
7.000	3.07	0.02 %	3.25 %	0.04 %	0.00 %	96.69 %	0.00 %	
7.400	3.07	0.01 %	3.25 %	0.09 %	0.00 %	96.65 %	0.00 %	Blood pH
8.000	3.03	0.00 %	3.24 %	0.37 %	0.00 %	96.38 %	0.00 %	·
9.000	2.75	0.00 %	3.14 %	3.60 %	0.00 %	93.26 %	0.00 %	
10.000	1.98	0.00 %	2.37 %	27.21 %	0.00 %	70.41 %	0.02 %	
11.000	1.02	0.00 %	0.69 %	78.86 %	0.00 %	20.41 %	0.05 %	
12.000	0.03	0.00 %	0.08 %	97.34 %	0.00 %	2.52 %	0.06 %	

Carbonate and acidity



Other graphs






pH-metric high logP Analyst: Pion Assay name: Assay ID: 18C-01009 Instrument ID: T312060 Filename:

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-01009 Points 32 to 62

Overall results

RMSD 0.147
Average ionic strength 0.171 M
Average temperature 25.0°C
Partition ratio 0.0751 : 1

Analyte concentration range 4115.3 µM to 4404.0 µM

Total points considered 23 of 31

Warnings and errors

Errors None

Warnings Excessive acidity error present

Four-Plus parameters


Alpha	0.130	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
S	0.9970	3/1/2018 12:08:08 PM	C:\Sirius T3\HCl18B27.t3r
jΗ	0.8	3/1/2018 12:08:08 PM	C:\Sirius T3\HCl18B27.t3r
jОН	-0.4	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r


Titrants

Sample

Φ	M16_octanol concentration factor	0.794
	Base pKa 1	4.70
	Acid pKa 2	8.94
	logP (XH2 +)	0.03
₩	logP (neutral XH)	3.06
	logP (X -)	-1.62

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M16_octanol logD	M16_octanol M16_octanolH2	M16_octanol M16_octanolH	M16_octanol M16_octanol	M16_octanol M16_octanolH2*	M16_octanol M16_octanolH*	M16_octanol M16_octanol*	
1.000	0.11	91.17 %	0.02 %	0.00 %	7.26 %	1.55 %	0.00 %	
1.200	0.15	90.33 %	0.03 %	0.00 %	7.20 %	2.44 %		Stomach pH
2.000	0.52	79.86 %	0.16 %	0.00 %	6.36 %	13.62 %	0.00 %	•
3.000	1.37	35.65 %	0.71 %	0.00 %	2.84 %	60.80 %	0.00 %	
4.000	2.28	5.45 %	1.09 %	0.00 %	0.43 %	93.02 %	0.00 %	
5.000	2.88	0.58 %	1.15 %	0.00 %	0.05 %	98.23 %	0.00 %	
6.000	3.03	0.06 %	1.16 %	0.00 %	0.00 %	98.78 %	0.00 %	
6.500	3.05	0.02 %	1.16 %	0.00 %	0.00 %	98.82 %	0.00 %	
7.000	3.05	0.01 %	1.16 %	0.01 %	0.00 %	98.82 %	0.00 %	
7.400	3.04	0.00 %	1.16 %	0.03 %	0.00 %	98.81 %	0.00 %	Blood pH
8.000	3.01	0.00 %	1.15 %	0.13 %	0.00 %	98.71 %	0.00 %	·
9.000	2.72	0.00 %	1.14 %	1.31 %	0.00 %	97.55 %	0.00 %	
10.000	1.96	0.00 %	1.02 %	11.72 %	0.00 %	87.24 %	0.02 %	
11.000	0.99	0.00 %	0.50 %	56.98 %	0.00 %	42.42 %	0.10 %	
12.000	0.01	0.00 %	0.08 %	92.84 %	0.00 %	6.91 %	0.17 %	

Carbonate and acidity



Other graphs






pH-metric high logP Pion Assay name: Analyst: 18C-01009 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-01009 Points 63 to 92

Overall results

RMSD 0.627
Average ionic strength 0.182 M
Average temperature 25.0°C
Partition ratio 0.2052 : 1

Analyte concentration range 3198.3 µM to 3391.8 µM

Total points considered 22 of 30

Warnings and errors

Errors None

Warnings One or more logP values out of range

Excessive acidity error present

Four-Plus parameters


Alpha	0.130	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
S	0.9970	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
jΗ	8.0	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
jΟH	-0.4	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r


Titrants

Sample

>	M16_octanol concentration factor	0.836
	Base pKa 1	4.70
	Acid pKa 2	8.94
	logP (XH2 +)	-5.10
>	logP (neutral XH)	3.01
ı	logP (X -)	-1.62

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	M16_octanol	
	logD	M16_octanolH2	M16_octanolH	M16_octanol	M16_octanolH2*	M16_octanolH*	M16_octanol*	
1.000	-0.69	95.95 %	0.02 %	0.00 %	0.00 %	4.03 %	0.00 %	
1.200	-0.49	93.73 %	0.03 %	0.00 %	0.00 %	6.24 %	0.00 %	Stomach pH
2.000	0.31	70.32 %	0.14 %	0.00 %	0.00 %	29.54 %	0.00 %	•
3.000	1.30	19.16 %	0.38 %	0.00 %	0.00 %	80.46 %	0.00 %	
4.000	2.23	2.31 %	0.46 %	0.00 %	0.00 %	97.22 %	0.00 %	
5.000	2.83	0.24 %	0.47 %	0.00 %	0.00 %	99.29 %	0.00 %	
6.000	2.99	0.02 %	0.47 %	0.00 %	0.00 %	99.50 %	0.00 %	
6.500	3.00	0.01 %	0.47 %	0.00 %	0.00 %	99.52 %	0.00 %	
7.000	3.00	0.00 %	0.47 %	0.01 %	0.00 %	99.52 %	0.00 %	
7.400	3.00	0.00 %	0.47 %	0.01 %	0.00 %	99.51 %	0.00 %	Blood pH
8.000	2.96	0.00 %	0.47 %	0.05 %	0.00 %	99.47 %	0.00 %	·
9.000	2.68	0.00 %	0.47 %	0.54 %	0.00 %	98.99 %	0.00 %	
10.000	1.92	0.00 %	0.45 %	5.15 %	0.00 %	94.38 %	0.03 %	
11.000	0.95	0.00 %	0.31 %	35.12 %	0.00 %	64.40 %	0.17 %	
12.000	-0.04	0.00 %	0.07 %	84.09 %	0.00 %	15.42 %	0.41 %	

Carbonate and acidity



Other graphs






Experiment start time: 3/1/2018 12:08:09 PM Sample name: M16_octanol

pH-metric high logP Analyst: Pion Assay name: 18C-01009 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name	_	2/27/2018 6:14:13 PM	
Sample by	Weight		Default value
Sample weight		2/28/2018 4:25:54 PM	
Formula weight	210.23 g/mol	2/27/2018 5:08:55 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	210.23	2/27/2018 5:08:55 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	2	2/27/2018 5:08:55 PM	User entered value
Sample is a	Ampholyte	2/27/2018 5:08:55 PM	User entered value
pKa 1	4.70	2/27/2018 5:08:55 PM	User entered value
Туре	Base	2/27/2018 5:08:55 PM	User entered value
pKa 2	8.94	2/27/2018 5:08:55 PM	User entered value
Туре	Acid	2/27/2018 5:08:55 PM	User entered value
logp (XH2 +)	-0.97	2/27/2018 5:09:30 PM	User entered value
logP (neutral XH)	1.47	2/27/2018 5:09:44 PM	User entered value
logP (X -)	-1.62	2/27/2018 5:09:50 PM	User entered value

Events

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
8:39.5	Initial pH = 3.74									
11:44.1	Data point 1	1.50000 mL	0.04499 mL	0.00202 mL	0.04000 mL	2.009	-0.00786	0.25633	0.00077	10.0 s
12:30.3	Data point 2	1.50000 mL	0.04499 mL	0.01726 mL	0.04000 mL	2.218	-0.00029	0.00763	0.00016	10.0 s
13:06.0	Data point 3		0.04499 mL		0.04000 mL		-0.00594	0.13958	0.00079	10.0 s
13:41.5	Data point 4		0.04499 mL		0.04000 mL			0.06954	0.00016	10.0 s
14:27.4	Data point 5	1.50000 mL	0.04499 mL	0.03702 mL	0.04000 mL	2.780	-0.00378	0.06388	0.00074	10.0 s
15:02.9	Data point 6	1.50000 mL	0.04499 mL		0.04000 mL		-0.00866	0.21191	0.00093	10.0 s
15:48.8	Data point 7	1.50000 mL	0.04499 mL	0.04259 mL	0.04000 mL	3.115	-0.00483	0.80380	0.00027	10.0 s
16:39.8	Data point 8	1.50000 mL			0.04000 mL		-0.00755	0.21309	0.00081	10.0 s
17:30.8	Data point 9	1.50000 mL			0.04000 mL		-0.00612		0.00034	10.0 s
18:21.9	Data point 10		0.04499 mL		0.04000 mL		-0.00337		0.00091	10.0 s
19:12.8	Data point 11		0.04499 mL		0.04000 mL		-0.00649	0.12208	0.00092	
20:03.7	Data point 12	1.50000 mL			0.04000 mL		-0.01178	0.65123	0.00072	10.0 s
20:39.2	Data point 13	1.50000 mL			0.04000 mL		-0.01852		0.00096	
21:15.1	Data point 14	1.50000 mL			0.04000 mL			0.79081	0.00088	13.0 s
21:53.5	Data point 15	1.50000 mL	0.04499 mL		0.04000 mL				0.00098	
22:53.5	Data point 16	1.50000 mL			0.04000 mL		-0.01692		0.00086	10.0 s
23:28.9	Data point 17		0.04499 mL		0.04000 mL		-0.01207		0.00061	10.0 s
24:04.3	Data point 18		0.04499 mL		0.04000 mL			0.54950	0.00080	10.0 s
24:50.2	Data point 19		0.04499 mL		0.04000 mL		-0.00092		0.00028	10.5 s
25:26.2	Data point 20		0.04499 mL		0.04000 mL			0.82899	0.00044	
26:01.6	Data point 21	1.50000 mL			0.04000 mL			0.62163	0.00087	
26:52.5	Data point 22	1.50000 mL	0.04499 mL		0.04000 mL		-0.00623	0.21568	0.00066	10.0 s
27:43.4	Data point 23		0.04499 mL		0.04000 mL		-0.01087	0.60836	0.00069	10.0 s
28:34.5	Data point 24		0.04499 mL		0.04000 mL		0.00564	0.16265	0.00069	10.5 s
29:26.1	Data point 25		0.04499 mL		0.04000 mL		-0.00670	0.84026	0.00036	10.0 s
30:12.0	Data point 26	1.50000 mL			0.04000 mL		-0.00717	0.22192	0.00075	10.0 s
30:47.5	Data point 27		0.04499 mL		0.04000 mL		0.00185	0.09507	0.00030	10.0 s
31:33.5	Data point 28	1.50000 mL			0.04000 mL		-0.00034	0.00072	0.00063	10.0 s
32:09.1	Data point 29	1.50000 mL	0.04499 mL		0.04000 mL		-0.00357	0.09765	0.00056	10.0 s
32:45.0	Data point 30		0.04499 mL	0.11117 mL	0.04000 mL		0.00119	0.00998	0.00059	10.0 s
33:20.6	Data point 31		0.04499 mL		0.04000 mL		0.00213	0.35467	0.00018	10.0 s
34:38.0	Data point 32		0.17241 mL		0.14000 mL		-0.00814	0.79043	0.00045	10.0 s
35:24.3	Data point 33		0.17241 mL		0.14000 mL		-0.00761	0.40848	0.00059	10.0 s
35:59.9	Data point 34		0.17241 mL		0.14000 mL		0.01002	0.94066	0.00051	10.0 s
36:35.6	Data point 35	1.50000 mL	0.17241 mL		0.14000 mL		0.00480	0.68312	0.00029	10.0 s
37:11.1	Data point 36	1.50000 mL	0.17241 mL	0.16472 mL	0.14000 mL	2.737	0.00482	0.88092	0.00025	10.0 s



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

Events (continued)

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
38:02.2	Data point 37	1.50000 mL	0.17241 mL	0.16898 mL	0.14000 mL	2.939	-0.00481	0.19222	0.00054	•
38:53.3	Data point 38	1.50000 mL	0.17241 mL	0.17272 mL	0.14000 mL	3.135	0.00012	0.00403	0.00010	10.0 s
39:39.2	Data point 39	1.50000 mL	0.17241 mL	0.17512 mL	0.14000 mL	3.328	0.00518	0.27140	0.00049	10.0 s
40:25.0	Data point 40	1.50000 mL	0.17241 mL	0.17693 mL	0.14000 mL	3.530	-0.00192	0.46039	0.00014	10.0 s
41:00.4	Data point 41	1.50000 mL	0.17241 mL	0.17858 mL	0.14000 mL	3.764	-0.00045	0.04173	0.00011	10.0 s
41:41.0	Data point 42						0.01048	0.42558	0.00079	10.0 s
42:16.4	Data point 43						0.00274	0.05649	0.00057	10.5 s
42:52.4	Data point 44	1.50000 mL	0.17241 mL	0.18156 mL	0.14000 mL	5.845	-0.01422	0.55499	0.00094	24.0 s
43:47.0	Data point 45	1.50000 mL	0.17241 mL	0.18175 mL	0.14000 mL	6.336	-0.01876	0.92228	0.00096	33.0 s
44:50.5	Data point 46	1.50000 mL	0.17241 mL	0.18184 mL	0.14000 mL	6.906	-0.02960	0.95500	0.00150	Timed out at
										59.5 s
46:15.9	Data point 47	1.50000 mL	0.17241 mL	0.18191 mL	0.14000 mL	7.300	-0.02968	0.97977	0.00148	Timed out at
										59.5 s
47:46.5	Data point 48						-0.01552		0.00078	
48:53.6	Data point 49						-0.01766		0.00095	
49:34.5	Data point 50						-0.00134		0.00097	
50:13.0	Data point 51						-0.00738			10.0 s
50:58.7	Data point 52						0.00095	0.02065	0.00033	10.0 s
51:49.7	Data point 53						-0.00379		0.00022	
52:25.6	Data point 54						-0.00945		0.00082	
53:16.7	Data point 55						-0.00222		0.00018	10.0 s
54:07.8	Data point 56						-0.00583		0.00082	10.0 s
54:58.8	Data point 57						-0.00232			10.0 s
55:44.8	Data point 58						-0.00054		0.00013	10.0 s
56:30.8	Data point 59						0.00672	0.53103	0.00046	
57:06.5	Data point 60						0.00008	0.00005		
57:42.4	Data point 61						0.00523	0.08361		10.0 s
58:18.3	Data point 62							0.20969	0.00091	10.0 s
59:40.6	Data point 63						-0.00419			10.0 s
1:00:27.0							-0.01035		0.00065	
1:01:12.7							-0.00423	0.05683	0.00088	10.0 s
1:01:48.4							0.00526	0.10194	0.00081	10.0 s
1:02:24.0	Data point 67						0.00672	0.26884	0.00064	10.0 s
1:03:09.9	Data point 68						0.01126	0.42408	0.00085	18.0 s
1:03:53.4							0.00666	0.70178	0.00039	
1:04:29.5							0.00922	0.44168	0.00069	
1:05:05.5							0.00860	0.24150	0.00086	10.0 s
	Data point 72						0.00769	0.25705		
	Data point 73						0.00226	0.01769	0.00084	
1:06:52.3	Data point 74	1.50000 mL	0.31378 mL	0.32298 mL	0.44000 mL	6.389	-0.01886		0.00097	
1:08:18.4	Data point 75	1.50000 mL	0.31378 ML	0.32326 ML	0.44000 ML	7.922	-0.05162	0.99435	0.00256	Timed out at
4.00.40.0	Data ==:=1.70	4 50000!	0.04070	0.000401	0.440001	0.400	0.04040	0.00000	0.00000	59.5 s
	Data point 76						-0.01649		0.00099	
	Data point 77						-0.01726		0.00098	
	Data point 78						-0.01698		0.00097	
	Data point 79						-0.01727		0.00099	
	Data point 80						-0.01163		0.00078	
	Data point 81						-0.00037		0.00056	
	Data point 82								0.00018	
	Data point 83							0.22886	0.00067	
	Data point 84							0.45250	0.00058	
	Data point 85							0.47694	0.00061	
	Data point 86							0.16909	0.00087	
	Data point 87							0.13436	0.00087	
	Data point 88							0.16576	0.00074	
1:19:12.1								0.61053	0.00078	
1.20.00.5	Data point 90	1.50000 ML	U.313/8 ML	U.3/0/6 ML	0.44000 ML	11.810	0.00711	0.92311	0.00037	10.0 8





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

Events (continued)

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
1:21:12.1	Data point 91 Data point 92 Assay volumes	1.50000 mL	0.31378 mL	0.40837 mL					0.00069 0.00077	



Assay name: pH-metric high logP Pion 18C-01009 Assay ID: Instrument ID: T312060

Assay ID: 18C-01009 Filename: C:\Sirius_T3\Meh	tap\20180228_exp2			312060 I_pH-metric high log <mark>f</mark>
Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings		J	· ·	•
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	12.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings	, , , , , , , , , , , , , , , , , , ,			
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium				
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.040 mL			
Partition solvent added	Automatic			
After partition addition, stir for	1 seconds			
Sample Sonication				
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	300 seconds			
After sonication stir for	5 seconds			
Sample Dissolution				
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution	To start pH			
Stir to dissolve for	120 seconds			
For discolution of at	100/			

10% For dissolution, stir at

Carbonate purge

Perform a carbonate purge No

Temperature Control

Yes Wait for temperature 25.0°C Required start temperature Acceptable deviation 0.5°C Time to wait 60 seconds 50%

Stir speed of

Titration 1

Titrate from Low to high pH

Adjust to start pH Yes After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds

Stirrer speed for partitioning 50%

Titration 2

Titrate from Low to high pH Add additional water 0.00 mL Additional partition solvent volume 0.100 mL Additional partition solvent added Automatic

After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds Stirrer speed for partitioning 55%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3		_	_	
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.300 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00100 dpH/dt			
Stability timeout after	60 seconds			

Calibration Settings

Value	Date/Time changed	Imported from
0.130	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
0.9970	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
8.0	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
-0.4	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
1.000	3/1/2018 12:08:09 PM	C:\Sirius_T3\KOH18B27.t3r
0.994	3/1/2018 12:08:08 PM	C:\Sirius_T3\HCl18B27.t3r
	0.130 0.9970 0.8 -0.4 1.000	1.000 3/1/2018 12:08:09 PM

Instrument Settings

Setting	Value	Batch Id	Install date
Instrument owner	Merck		
Instrument ID	T312060		
Instrument type	T3 Simulator		
Software version	1.1.3.0		
Dispenser module		T3DM1200361	3/31/2009 5:24:52 AM
Dispenser 0	Water		3/31/2009 5:25:05 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Water (0.15 M KCI)	02-06-2018	2/27/2018 10:05:59 AM
Dispenser 2	Acid		3/31/2009 5:25:11 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Acid (0.5 M HCI)	02-27-2018	2/27/2018 10:27:22 AM
Dispenser 1	Base		3/31/2009 5:25:21 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Base (0.5 M KOH)	9/22/2017	2/27/2018 10:21:22 AM
Dispenser 5	Cosolvent		3/31/2009 5:26:24 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Distribution valve 5	Distribution Valve		3/31/2009 5:28:19 AM
Firmware version	1.1.3		
Port A	Methanol (80%, 0.15 M KCl)	09-26-17	2/7/2018 9:42:01 AM
Port B	Cyclohexane	11-01-17	2/27/2018 10:37:57 AM
Dispenser 3	Buffer		8/3/2010 5:05:16 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Dodecane	2018/01/31	2/28/2018 10:18:04 AM
Dispenser 6	Octanol		10/22/2010 10:52:43 AM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Syringe volume Firmware version	0.5 mL 1.2.1(r2)		
Titrant	Octanol	01-31-2018	2/27/2018 9:59:35 AM
Titrator	Octanol		3/31/2009 5:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2	1011111200101	0/01/2000 0.24.17 / WI
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1DI0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0923	1/23/2018 2:01:00 PM
E0 calibration	+4.26 mV		3/1/2018 12:08:36 PM
Filling solution	3M KCI	KCL097	2/27/2018 9:49:43 AM
Liquids			
Wash 1	50% IPA:50% Water		2/28/2018 10:23:32 AM
Wash 2	0.5% Trition X-100 in H20		2/28/2018 10:23:34 AM
Buffer position 1	pH7 Wash		2/28/2018 10:24:06 AM
Buffer position 2	pH 7		2/28/2018 10:24:08 AM
Storage position	9.201002 ml	02 27 2010	2/28/2018 10:21:14 AM
Wash water Waste	8.2e+003 mL 7.3e+003 mL	02-27-2018	2/27/2018 9:54:39 AM 11/28/2017 10:36:29 AM
Temperature controller	7.30+003 IIIL		8/5/2010 6:35:13 AM
Turbidity detector			3/31/2009 5:24:45 AM
Spectrometer		074811	11/23/2010 11:22:28 AM
Dip probe		10196	11/20/2010 11.22.20 / tivi
Wavelength coefficient A0	183.333	10100	
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	112:08:55		11/23/2010 11:22:28 AM
Calibrated on	2/27/2018 10:40:38 AM		
Integration time	40		
Scans averaged	10		
Autoloader		T3AL1200345	11/10/2015 9:34:13 AM
Left-right axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1DI0DO4 Norgren I/O		
Configuration	Titratian position		
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume Maximum alternate vial volume	3.50 mL 25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting Spectrometer calibration stir duration	Value 5 s	Batch Id	Install date
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume Spectrometer calibration wash stir duration	20.0 mL 5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		

Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity detection method Turbidity wavelength to assess		
Turbidity maximum absorbance		0.100
Turbidity probe threshold	50.00	50.00
Experiment Log		

- [2:16] Air gap created for Water (0.15 M KCI)
- [2:17] Air gap created for Acid (0.5 M HCI)
- [2:17] Air gap created for Base (0.5 M KOH)
- [2:17] Air gap released for Water (0.15 M KCI)
- [2:21] Titrator arm moved over Titration position
- [2:21] Titration 1 of 3
- [2:21] Adding initial titrants
- [2:21] Automatically add 1.50000 mL of water
- [2:46] Dispensed 1.500000 mL of Water (0.15 M KCI)
- [2:51] Titrator arm moved over Drain
- [8:32] Titrator arm moved to Titration position
- [8:32] Argon flow rate set to 100
- [8:32] Stirrer speed set to 10
- [8:37] Automatically add 0.04000 mL of Octanol
- [8:38] Dispensed 0.040005 mL of Octanol
- [8:39] Initial pH = 3.74
- [8:39] Iterative adjust 3.74 -> 2.00
- [8:39] pH 3.74 -> 2.00
- [8:41] Air gap released for Acid (0.5 M HCl)
- [8:41] Dispensed 0.042615 mL of Acid (0.5 M HCI)
- [8:46] pH 2.02 -> 2.00
- [8:47] Dispensed 0.002375 mL of Acid (0.5 M HCl)
- [8:52] Holding pH 2.00
- [10:52] Stirrer speed set to 0
- [10:52] Stirrer speed set to 50
- [10:52] Iterative adjust 1.98 -> 2.00
- [10:52] pH 1.98 -> 2.00
- [10:53] Air gap released for Base (0.5 M KOH)
- [10:53] Dispensed 0.002023 mL of Base (0.5 M KOH)
- [11:44] Stirrer speed set to 0
- [11:54] Datapoint id 1 collected
- [11:54] Stirrer speed set to 50
- [11:59] pH 2.02 -> 2.22
- [11:59] Using cautious pH adjust
- [11:59] Dispensed 0.007596 mL of Base (0.5 M KOH)
- [12:04] Stepping pH = 2.10
- [12:04] Dispensed 0.005997 mL of Base (0.5 M KOH)
- [12:10] Stepping pH = 2.19
- [12:10] Dispensed 0.001646 mL of Base (0.5 M KOH)
- [12:15] Stepping pH = 2.22
- [12:30] Stirrer speed set to 0
- [12:40] Datapoint id 2 collected
- Reported at: 3/2/2018 2:54:07 PM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [12:40] Charge balance equation is out by -0.2%
- [12:40] Stirrer speed set to 50
- [12:45] pH 2.23 -> 2.43
- [12:45] Using charge balance adjust
- [12:45] Dispensed 0.009407 mL of Base (0.5 M KOH)
- [13:06] Stirrer speed set to 0
- [13:16] Datapoint id 3 collected
- [13:16] Charge balance equation is out by -6.0%
- [13:16] Stirrer speed set to 50
- [13:21] pH 2.42 -> 2.62
- [13:21] Using charge balance adjust
- [13:21] Dispensed 0.006021 mL of Base (0.5 M KOH)
- [13:41] Stirrer speed set to 0
- [13:51] Datapoint id 4 collected
- [13:51] Charge balance equation is out by -18.1%
- [13:51] Stirrer speed set to 50
- [13:56] pH 2.59 -> 2.79
- [13:56] Using cautious pH adjust
- [13:56] Dispensed 0.002046 mL of Base (0.5 M KOH)
- [14:01] Stepping pH = 2.67
- [14:02] Dispensed 0.001787 mL of Base (0.5 M KOH)
- [14:07] Stepping pH = 2.76
- [14:07] Dispensed 0.000494 mL of Base (0.5 M KOH)
- [14:12] Stepping pH = 2.78
- 14.071 Otimor and add a
- [14:27] Stirrer speed set to 0
- [14:37] Datapoint id 5 collected
- [14:37] Charge balance equation is out by -5.4%
- [14:37] Stirrer speed set to 50
- [14:42] pH 2.79 -> 2.99
- [14:42] Using charge balance adjust
- [14:42] Dispensed 0.002705 mL of Base (0.5 M KOH)
- [15:02] Stirrer speed set to 0
- [15:12] Datapoint id 6 collected
- [15:12] Charge balance equation is out by -31.5%
- [15:12] Stirrer speed set to 50
- [15:18] pH 2.93 -> 3.13
- [15:18] Using cautious pH adjust
- [15:18] Dispensed 0.001035 mL of Base (0.5 M KOH)
- [15:23] Stepping pH = 2.99
- [15:23] Dispensed 0.001294 mL of Base (0.5 M KOH)
- [15:28] Stepping pH = 3.08
- [15:28] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [15:33] Stepping pH = 3.12
- [15:48] Stirrer speed set to 0
- [15:58] Datapoint id 7 collected
- [15:58] Charge balance equation is out by -39.0%
- [15:58] Stirrer speed set to 50
- [16:03] pH 3.12 -> 3.32
- [16:03] Using cautious pH adjust
- [16:04] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [16:09] Stepping pH = 3.18
- [16:09] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [16:14] Stepping pH = 3.27
- [16:14] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [16:19] Stepping pH = 3.30
- [16:19] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [16:24] Stepping pH = 3.32
- [16:39] Stirrer speed set to 0
- [16:49] Datapoint id 8 collected



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [16:49] Charge balance equation is out by -50.3%
- [16:49] Stirrer speed set to 50
- [16:55] pH 3.32 -> 3.52
- [16:55] Using cautious pH adjust
- [16:55] Dispensed 0.000611 mL of Base (0.5 M KOH)
- [17:00] Stepping pH = 3.37
- [17:00] Dispensed 0.000823 mL of Base (0.5 M KOH)
- [17:05] Stepping pH = 3.45
- [17:05] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [17:10] Stepping pH = 3.49
- [17:10] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [17:15] Stepping pH = 3.51
- [17:30] Stirrer speed set to 0
- [17:40] Datapoint id 9 collected
- [17:40] Charge balance equation is out by -79.1%
- [17:40] Stirrer speed set to 50
- [17:45] pH 3.51 -> 3.71
- [17:45] Using cautious pH adjust
- [17:46] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [17:51] Stepping pH = 3.56
- [17:51] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [17:56] Stepping pH = 3.63
- [17:56] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [18:01] Stepping pH = 3.69
- [18:01] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [18:06] Stepping pH = 3.70
- [18:21] Stirrer speed set to 0
- 10.21] Stiller Speed Set to 0
- [18:31] Datapoint id 10 collected
- [18:31] Charge balance equation is out by -81.1%
- [18:31] Stirrer speed set to 50
- [18:37] pH 3.70 -> 3.90
- [18:37] Using cautious pH adjust
- [18:37] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [18:42] Stepping pH = 3.77
- [18:42] Dispensed 0.000729 mL of Base (0.5 M KOH)
- [18:47] Stepping pH = 3.86
- [18:47] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [18:52] Stepping pH = 3.89
- [18:52] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [18:57] Stepping pH = 3.89
- [19:12] Stirrer speed set to 0
- [19:22] Datapoint id 11 collected
- [19:22] Charge balance equation is out by -36.1%
- [19:22] Stirrer speed set to 50
- [19:27] pH 3.89 -> 4.09
- [19:27] Using cautious pH adjust
- [19:28] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [19:33] Stepping pH = 4.04
- [19:33] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [19:38] Stepping pH = 4.08
- [19:38] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [19:43] Stepping pH = 4.08
- [19:43] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [19:48] Stepping pH = 4.11
- [20:03] Stirrer speed set to 0
- [20:13] Datapoint id 12 collected
- [20:13] Charge balance equation is out by 17.6%
- [20:13] Stirrer speed set to 50
- [20:18] pH 4.11 -> 4.31



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [20:18] Using cautious pH adjust
- [20:18] Dispensed 0.000894 mL of Base (0.5 M KOH)
- [20:24] Stepping pH = 4.40
- [20:39] Stirrer speed set to 0
- [20:49] Datapoint id 13 collected
- [20:49] Charge balance equation is out by 50.0%
- [20:49] Stirrer speed set to 50
- [20:54] pH 4.39 -> 4.59
- [20:54] Using cautious pH adjust
- [20:54] Dispensed 0.001011 mL of Base (0.5 M KOH)
- [21:00] Stepping pH = 5.11
- [21:15] Stirrer speed set to 0
- [21:28] Datapoint id 14 collected
- [21:28] Charge balance equation is out by 50.0%
- [21:28] Stirrer speed set to 50
- [21:33] pH 5.12 -> 5.32
- [21:33] Using cautious pH adjust
- [21:33] Dispensed 0.000611 mL of Base (0.5 M KOH)
- [21:38] Stepping pH = 8.25
- [21:53] Stirrer speed set to 0
- [22:28] Datapoint id 15 collected
- [22:28] Charge balance equation is out by 50.0%
- [22:28] Stirrer speed set to 50
- [22:33] pH 8.53 -> 8.73
- [22:33] Using cautious pH adjust
- [22:33] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [22:38] Stepping pH = 9.19
- [22:53] Stirrer speed set to 0
- [23:03] Datapoint id 16 collected
- [23:03] Charge balance equation is out by 50.0%
- [23:03] Stirrer speed set to 50
- [23:08] pH 9.18 -> 9.38
- [23:08] Using cautious pH adjust
- [23:08] Dispensed 0.001011 mL of Base (0.5 M KOH)
- [23:13] Stepping pH = 9.60
- [23:28] Stirrer speed set to 0
- [23:38] Datapoint id 17 collected
- [23:38] Charge balance equation is out by 50.0%
- [23:38] Stirrer speed set to 50
- [23:44] pH 9.59 -> 9.79
- [23:44] Using cautious pH adjust
- [23:44] Dispensed 0.000847 mL of Base (0.5 M KOH)
- [23:49] Stepping pH = 9.79
- [24:04] Stirrer speed set to 0
- [24:14] Datapoint id 18 collected
- [24:14] Charge balance equation is out by 50.0%
- [24:14] Stirrer speed set to 50
- [24:19] pH 9.78 -> 9.98
- [24:19] Using cautious pH adjust
- [24:19] Dispensed 0.000729 mL of Base (0.5 M KOH)
- [24:24] Stepping pH = 9.91
- [24:24] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [24:29] Stepping pH = 9.94
- [24:30] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [24:35] Stepping pH = 9.98
- [24:50] Stirrer speed set to 0
- [25:00] Datapoint id 19 collected
- [25:00] Charge balance equation is out by 10.9%
- [25:00] Stirrer speed set to 50



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [25:05] pH 9.98 -> 10.18
- [25:05] Using charge balance adjust
- [25:05] Dispensed 0.001199 mL of Base (0.5 M KOH)
- [25:26] Stirrer speed set to 0
- [25:36] Datapoint id 20 collected
- [25:36] Charge balance equation is out by -11.3%
- [25:36] Stirrer speed set to 50
- [25:41] pH 10.16 -> 10.36
- [25:41] Using charge balance adjust
- [25:41] Dispensed 0.001082 mL of Base (0.5 M KOH)
- [26:01] Stirrer speed set to 0
- [26:11] Datapoint id 21 collected
- [26:11] Charge balance equation is out by -32.3%
- [26:11] Stirrer speed set to 50
- [26:16] pH 10.29 -> 10.49
- [26:16] Using cautious pH adjust
- [26:16] Dispensed 0.000517 mL of Base (0.5 M KOH)
- [26:21] Stepping pH = 10.36
- [26:21] Dispensed 0.000611 mL of Base (0.5 M KOH)
- [26:27] Stepping pH = 10.43
- [26:27] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [26:32] Stepping pH = 10.47
- [26:32] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [26:37] Stepping pH = 10.49
- [26:52] Stirrer speed set to 0
- [27:02] Datapoint id 22 collected
- [27:02] Charge balance equation is out by -58.4%
- [27:02] Stirrer speed set to 50
- [27:07] pH 10.48 -> 10.68
- [27:07] Using cautious pH adjust
- [27:07] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [27:12] Stepping pH = 10.55
- [27:12] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [27:18] Stepping pH = 10.61
- [27:18] Dispensed 0.000494 mL of Base (0.5 M KOH)
- [27:23] Stepping pH = 10.66
- [27:23] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [27:28] Stepping pH = 10.68
- [27:43] Stirrer speed set to 0
- [27:53] Datapoint id 23 collected
- [27:53] Charge balance equation is out by -67.1%
- [27:53] Stirrer speed set to 50
- [27:58] pH 10.67 -> 10.87
- [27:58] Using cautious pH adjust
- [27:58] Dispensed 0.000729 mL of Base (0.5 M KOH)
- [28:03] Stepping pH = 10.73
- [28:03] Dispensed 0.000964 mL of Base (0.5 M KOH)
- [28:08] Stepping pH = 10.80
- [28:09] Dispensed 0.000706 mL of Base (0.5 M KOH)
- [28:14] Stepping pH = 10.84
- [28:14] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [28:19] Stepping pH = 10.86
- [28:34] Stirrer speed set to 0
- [28:45] Datapoint id 24 collected
- [28:45] Charge balance equation is out by -85.3%
- [28:45] Stirrer speed set to 50
- [28:50] pH 10.85 -> 11.05
- [28:50] Using cautious pH adjust
- [28:50] Dispensed 0.001011 mL of Base (0.5 M KOH)
- Reported at: 3/2/2018 2:54:07 PM



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [28:55] Stepping pH = 10.91
- [28:55] Dispensed 0.001270 mL of Base (0.5 M KOH)
- [29:00] Stepping pH = 11.00
- [29:00] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [29:05] Stepping pH = 11.04
- [29:05] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [29:11] Stepping pH = 11.05
- [29:26] Stirrer speed set to 0
- [29:36] Datapoint id 25 collected
- [29:36] Charge balance equation is out by -55.8%
- [29:36] Stirrer speed set to 50
- [29:41] pH 11.04 -> 11.24
- [29:41] Using cautious pH adjust
- [29:41] Dispensed 0.001505 mL of Base (0.5 M KOH)
- [29:46] Stepping pH = 11.13
- [29:46] Dispensed 0.001364 mL of Base (0.5 M KOH)
- [29:51] Stepping pH = 11.21
- [29:51] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [29:56] Stepping pH = 11.23
- [30:11] Stirrer speed set to 0
- [30:22] Datapoint id 26 collected
- [30:22] Charge balance equation is out by -13.7%
- [30:22] Stirrer speed set to 50
- [30:27] pH 11.23 -> 11.43
- [30:27] Using charge balance adjust
- [30:27] Dispensed 0.004516 mL of Base (0.5 M KOH)
- [30:47] Stirrer speed set to 0
- [30:57] Datapoint id 27 collected
- [30:57] Charge balance equation is out by -16.5%
- [30:57] Stirrer speed set to 50
- [31:02] pH 11.40 -> 11.60
- [31:02] Using cautious pH adjust
- [31:02] Dispensed 0.003340 mL of Base (0.5 M KOH)
- [31:07] Stepping pH = 11.51
- [31:08] Dispensed 0.002187 mL of Base (0.5 M KOH)
- [31:13] Stepping pH = 11.57
- [31:13] Dispensed 0.000847 mL of Base (0.5 M KOH)
- [31:18] Stepping pH = 11.59
- [31:33] Stirrer speed set to 0
- [31:43] Datapoint id 28 collected
- [31:43] Charge balance equation is out by 4.5%
- [31:43] Stirrer speed set to 50
- [31:48] pH 11.59 -> 11.79
- [31:48] Using charge balance adjust
- [31:48] Dispensed 0.010583 mL of Base (0.5 M KOH)
- [32:09] Stirrer speed set to 0
- [32:19] Datapoint id 29 collected
- [32:19] Charge balance equation is out by -12.6%
- [32:19] Stirrer speed set to 50
- [32:24] pH 11.77 -> 11.97
- [32:24] Using charge balance adjust
- [32:24] Dispensed 0.016157 mL of Base (0.5 M KOH)
- [32:44] Stirrer speed set to 0
- [32:55] Datapoint id 30 collected
- [32:55] Charge balance equation is out by -7.8%
- [32:55] Stirrer speed set to 50
- [33:00] pH 11.96 -> 12.05
- [33:00] Using charge balance adjust
- [33:00] Dispensed 0.009737 mL of Base (0.5 M KOH)



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [33:20] Stirrer speed set to 0
- [33:30] Datapoint id 31 collected
- [33:30] Charge balance equation is out by -63.2%
- [33:30] Titration 2 of 3
- [33:30] Adding initial titrants
- [33:30] Automatically add 0.10000 mL of Octanol
- [33:33] Dispensed 0.100000 mL of Octanol
- [33:33] Stirrer speed set to 10
- [33:34] Stirrer speed set to 55
- [33:34] Iterative adjust 12.04 -> 2.00
- [33:34] pH 12.04 -> 2.00
- [33:36] Dispensed 0.100000 mL of Acid (0.5 M HCI)
- [33:41] pH 2.31 -> 2.00
- [33:42] Dispensed 0.025141 mL of Acid (0.5 M HCI)
- [33:47] pH 2.02 -> 2.00
- [33:47] Dispensed 0.002281 mL of Acid (0.5 M HCl)
- [34:38] Stirrer speed set to 0
- [34:48] Datapoint id 32 collected
- [34:48] Stirrer speed set to 55
- [34:53] pH 1.98 -> 2.18
- [34:53] Using cautious pH adjust
- [34:53] Dispensed 0.009784 mL of Base (0.5 M KOH)
- [34:58] Stepping pH = 2.06
- [34:58] Dispensed 0.007502 mL of Base (0.5 M KOH)
- [35:03] Stepping pH = 2.15
- [35:04] Dispensed 0.001999 mL of Base (0.5 M KOH)
- [35:09] Stepping pH = 2.18
- [35:24] Stirrer speed set to 0
- [35:34] Datapoint id 33 collected
- [35:34] Charge balance equation is out by 1.4%
- [35:34] Stirrer speed set to 55
- [35:39] pH 2.19 -> 2.39
- [35:39] Using charge balance adjust
- [35:39] Dispensed 0.011994 mL of Base (0.5 M KOH)
- [35:59] Stirrer speed set to 0
- [36:10] Datapoint id 34 collected
- [36:10] Charge balance equation is out by -3.2%
- [36:10] Stirrer speed set to 55
- [36:15] pH 2.39 -> 2.59
- [36:15] Using charge balance adjust
- [36:15] Dispensed 0.007549 mL of Base (0.5 M KOH)
- [36:35] Stirrer speed set to 0
- [36:45] Datapoint id 35 collected
- [36:45] Charge balance equation is out by -9.6%
- [36:45] Stirrer speed set to 55
- [36:50] pH 2.58 -> 2.78
- [36:50] Using charge balance adjust
- [36:50] Dispensed 0.004986 mL of Base (0.5 M KOH)
- [37:11] Stirrer speed set to 0
- [37:21] Datapoint id 36 collected
- [37:21] Charge balance equation is out by -20.0%
- [37:21] Stirrer speed set to 55
- [37:26] pH 2.74 -> 2.94
- [37:26] Using cautious pH adjust
- [37:26] Dispensed 0.001764 mL of Base (0.5 M KOH)
- [37:31] Stepping pH = 2.82
- [37:31] Dispensed 0.001623 mL of Base (0.5 M KOH)
- [37:36] Stepping pH = 2.90
- [37:36] Dispensed 0.000659 mL of Base (0.5 M KOH)



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [37:41] Stepping pH = 2.93
- [37:42] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [37:47] Stepping pH = 2.94
- [38:02] Stirrer speed set to 0
- [38:12] Datapoint id 37 collected
- [38:12] Charge balance equation is out by -20.5%
- [38:12] Stirrer speed set to 55
- [38:17] pH 2.94 -> 3.14
- [38:17] Using cautious pH adjust
- [38:17] Dispensed 0.001223 mL of Base (0.5 M KOH)
- [38:22] Stepping pH = 3.00
- [38:22] Dispensed 0.001529 mL of Base (0.5 M KOH)
- [38:27] Stepping pH = 3.08
- [38:27] Dispensed 0.000823 mL of Base (0.5 M KOH)
- [38:33] Stepping pH = 3.13
- [38:33] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [38:38] Stepping pH = 3.14
- [38:53] Stirrer speed set to 0
- [39:03] Datapoint id 38 collected
- [39:03] Charge balance equation is out by -53.0%
- [39:03] Stirrer speed set to 55
- [39:08] pH 3.14 -> 3.34
- [39:08] Using cautious pH adjust
- [39:08] Dispensed 0.000941 mL of Base (0.5 M KOH)
- [39:13] Stepping pH = 3.21
- [39:13] Dispensed 0.001035 mL of Base (0.5 M KOH)
- [39:18] Stepping pH = 3.30
- [39:18] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [39:24] Stepping pH = 3.33
- [39:39] Stirrer speed set to 0
- [39:49] Datapoint id 39 collected
- [39:49] Charge balance equation is out by -27.4%
- [39:49] Stirrer speed set to 55
- [39:54] pH 3.34 -> 3.54
- [39:54] Using cautious pH adjust
- [39:54] Dispensed 0.000823 mL of Base (0.5 M KOH)
- [39:59] Stepping pH = 3.42
- [39:59] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [40:04] Stepping pH = 3.51
- [40:04] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [40:09] Stepping pH = 3.53
- [40:25] Stirrer speed set to 0
- [40:35] Datapoint id 40 collected
- [40:35] Charge balance equation is out by -8.9%
- [40:35] Stirrer speed set to 55
- [40:40] pH 3.53 -> 3.73
- [40:40] Using charge balance adjust
- [40:40] Dispensed 0.001646 mL of Base (0.5 M KOH)
- [41:00] Stirrer speed set to 0
- [41:10] Datapoint id 41 collected
- [41:10] Charge balance equation is out by 16.2%
- [41:10] Stirrer speed set to 55
- [41:15] pH 3.77 -> 3.97
- [41:15] Using cautious pH adjust
- [41:15] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [41:20] Stepping pH = 3.96
- [41:20] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [41:25] Stepping pH = 3.96
- [41:41] Stirrer speed set to 0



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [41:51] Datapoint id 42 collected
- [41:51] Charge balance equation is out by 47.0%
- [41:51] Stirrer speed set to 55
- [41:56] pH 3.97 -> 4.17
- [41:56] Using cautious pH adjust
- [41:56] Dispensed 0.001011 mL of Base (0.5 M KOH)
- [42:01] Stepping pH = 4.29
- [42:16] Stirrer speed set to 0
- [42:27] Datapoint id 43 collected
- [42:27] Charge balance equation is out by 50.0%
- [42:27] Stirrer speed set to 55
- [42:32] pH 4.29 -> 4.49
- [42:32] Using cautious pH adjust
- [42:32] Dispensed 0.001011 mL of Base (0.5 M KOH)
- [42:37] Stepping pH = 5.79
- [42:52] Stirrer speed set to 0
- [43:16] Datapoint id 44 collected
- [43:16] Charge balance equation is out by 50.0%
- [43:16] Stirrer speed set to 55
- [43:21] pH 5.84 -> 6.04
- [43:21] Using cautious pH adjust
- [43:21] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [43:26] Stepping pH = 5.94
- [43:26] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [43:31] Stepping pH = 6.18
- [43:46] Stirrer speed set to 0
- [44:20] Datapoint id 45 collected
- [44:20] Charge balance equation is out by 18.2%
- [44:20] Stirrer speed set to 55
- [44:25] pH 6.29 -> 6.49
- [44:25] Using cautious pH adjust
- [44:25] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [44:30] Stepping pH = 6.38
- [44:30] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [44:35] Stepping pH = 6.89
- [44:50] Stirrer speed set to 0
- [45:50] Datapoint id 46 collected
- [45:50] Charge balance equation is out by 7.4%
- [45:50] Stirrer speed set to 55
- [45:55] pH 6.83 -> 7.03
- [45:55] Using charge balance adjust
- [45:55] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [46:15] Stirrer speed set to 0
- [47:15] Datapoint id 47 collected
- [47:15] Charge balance equation is out by 136.7%
- [47:15] Stirrer speed set to 55
- [47:21] pH 7.23 -> 7.43
- [47:21] Using cautious pH adjust
- [47:21] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [47:26] Stepping pH = 7.19
- [47:26] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [47:31] Stepping pH = 8.55
- [47:46] Stirrer speed set to 0
- [48:28] Datapoint id 48 collected
- [48:28] Charge balance equation is out by -241.0%
- [48:28] Stirrer speed set to 55
- [48:33] pH 8.44 -> 8.64
- [48:33] Using cautious pH adjust
- [48:33] Dispensed 0.000329 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [48:38] Stepping pH = 9.15
- [48:53] Stirrer speed set to 0
- [49:09] Datapoint id 49 collected
- [49:09] Charge balance equation is out by 50.0%
- [49:09] Stirrer speed set to 55
- [49:14] pH 9.13 -> 9.33
- [49:14] Using cautious pH adjust
- [49:14] Dispensed 0.000894 mL of Base (0.5 M KOH)
- [49:19] Stepping pH = 9.72
- [49:34] Stirrer speed set to 0
- [49:47] Datapoint id 50 collected
- [49:47] Charge balance equation is out by 50.0%
- [49:47] Stirrer speed set to 55
- [49:52] pH 9.72 -> 9.92
- [49:52] Using cautious pH adjust
- [49:52] Dispensed 0.000988 mL of Base (0.5 M KOH)
- [49:57] Stepping pH = 10.04
- [50:12] Stirrer speed set to 0
- [50:23] Datapoint id 51 collected
- [50:23] Charge balance equation is out by 50.0%
- [50:23] Stirrer speed set to 55
- [50:28] pH 10.04 -> 10.24
- [50:28] Using cautious pH adjust
- [50:28] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [50:33] Stepping pH = 10.22
- [50:33] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [50:38] Stepping pH = 10.23
- [50:38] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [50:43] Stepping pH = 10.24
- [50:58] Stirrer speed set to 0
- [51:08] Datapoint id 52 collected
- [51:08] Charge balance equation is out by 39.5%
- [51:08] Stirrer speed set to 55
- [51:13] pH 10.24 -> 10.44
- [51:13] Using cautious pH adjust
- [51:13] Dispensed 0.000729 mL of Base (0.5 M KOH)
- [51:19] Stepping pH = 10.36
- [51:19] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [51:24] Stepping pH = 10.41
- [51:24] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [51:29] Stepping pH = 10.43
- [51:29] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [51:34] Stepping pH = 10.43
- [51:49] Stirrer speed set to 0
- [52:00] Datapoint id 53 collected
- [52:00] Charge balance equation is out by 9.5%
- [52:00] Stirrer speed set to 55
- [52:05] pH 10.42 -> 10.62
- [52:05] Using charge balance adjust
- [52:05] Dispensed 0.001482 mL of Base (0.5 M KOH)
- [52:25] Stirrer speed set to 0
- [52:35] Datapoint id 54 collected
- [52:35] Charge balance equation is out by -24.8%
- [52:35] Stirrer speed set to 55
- [52:40] pH 10.58 -> 10.78
- [52:40] Using cautious pH adjust
- [52:40] Dispensed 0.000823 mL of Base (0.5 M KOH)
- [52:45] Stepping pH = 10.64
- [52:46] Dispensed 0.001011 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [52:51] Stepping pH = 10.71
- [52:51] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [52:56] Stepping pH = 10.76
- [52:56] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [53:01] Stepping pH = 10.77
- [53:16] Stirrer speed set to 0
- [53:26] Datapoint id 55 collected
- [53:26] Charge balance equation is out by -67.6%
- [53:26] Stirrer speed set to 55
- [53:31] pH 10.77 -> 10.97
- [53:31] Using cautious pH adjust
- [53:32] Dispensed 0.001058 mL of Base (0.5 M KOH)
- [53:37] Stepping pH = 10.84
- [53:37] Dispensed 0.001223 mL of Base (0.5 M KOH)
- [53:42] Stepping pH = 10.91
- [53:42] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [53:47] Stepping pH = 10.96
- [53:47] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [53:52] Stepping pH = 10.96
- [54:07] Stirrer speed set to 0
- 54.07] Suiter speed set to 0
- [54:17] Datapoint id 56 collected
- [54:17] Charge balance equation is out by -47.4%
- [54:17] Stirrer speed set to 55
- [54:22] pH 10.96 -> 11.16
- [54:22] Using cautious pH adjust
- [54:23] Dispensed 0.001482 mL of Base (0.5 M KOH)
- [54:28] Stepping pH = 11.04
- [54:28] Dispensed 0.001458 mL of Base (0.5 M KOH)
- [54:33] Stepping pH = 11.11
- [54:33] Dispensed 0.000847 mL of Base (0.5 M KOH)
- [54:38] Stepping pH = 11.14
- [54:38] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [54:43] Stepping pH = 11.15
- [54:58] Stirrer speed set to 0
- [55:08] Datapoint id 57 collected
- [55:08] Charge balance equation is out by -38.3%
- [55:08] Stirrer speed set to 55
- [55:14] pH 11.15 -> 11.35
- [55:14] Using cautious pH adjust
- [55:14] Dispensed 0.002234 mL of Base (0.5 M KOH)
- [55:19] Stepping pH = 11.23
- [55:19] Dispensed 0.002258 mL of Base (0.5 M KOH)
- [55:24] Stepping pH = 11.31
- [55:24] Dispensed 0.001058 mL of Base (0.5 M KOH)
- [55:29] Stepping pH = 11.34
- [55:44] Stirrer speed set to 0
- [55:54] Datapoint id 58 collected
- [55:54] Charge balance equation is out by -24.3%
- [55:54] Stirrer speed set to 55
- [55:59] pH 11.34 -> 11.54
- [55:59] Using cautious pH adjust
- [56:00] Dispensed 0.003434 mL of Base (0.5 M KOH)
- [56:05] Stepping pH = 11.44
- [56:05] Dispensed 0.002775 mL of Base (0.5 M KOH)
- [56:10] Stepping pH = 11.50
- [56:10] Dispensed 0.001364 mL of Base (0.5 M KOH)
- [56:15] Stepping pH = 11.53
- [56:30] Stirrer speed set to 0
- [56:40] Datapoint id 59 collected



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [56:40] Charge balance equation is out by -10.4%
- [56:40] Stirrer speed set to 55
- [56:45] pH 11.54 -> 11.74
- [56:45] Using charge balance adjust
- [56:46] Dispensed 0.010748 mL of Base (0.5 M KOH)
- [57:06] Stirrer speed set to 0
- [57:16] Datapoint id 60 collected
- [57:16] Charge balance equation is out by -12.6%
- [57:16] Stirrer speed set to 55
- [57:21] pH 11.71 -> 11.91
- [57:21] Using charge balance adjust
- [57:22] Dispensed 0.016345 mL of Base (0.5 M KOH)
- [57:42] Stirrer speed set to 0
- [57:52] Datapoint id 61 collected
- [57:52] Charge balance equation is out by -9.7%
- [57:52] Stirrer speed set to 55
- [57:57] pH 11.90 -> 12.05
- [57:57] Using charge balance adjust
- [57:58] Dispensed 0.018156 mL of Base (0.5 M KOH)
- [58:18] Stirrer speed set to 0
- [58:28] Datapoint id 62 collected
- [58:28] Charge balance equation is out by -33.1%
- [58:28] Titration 3 of 3
- [58:28] Adding initial titrants
- [58:28] Automatically add 0.30000 mL of Octanol
- [58:35] Dispensed 0.300000 mL of Octanol
- [58:35] Stirrer speed set to 10
- [58:36] Stirrer speed set to 60
- [58:36] Iterative adjust 12.04 -> 2.00
- [58:36] pH 12.04 -> 2.00
- [58:38] Dispensed 0.100000 mL of Acid (0.5 M HCI)
- [58:44] pH 2.47 -> 2.00
- [58:45] Dispensed 0.037465 mL of Acid (0.5 M HCI)
- [58:50] pH 2.03 -> 2.00
- [58:50] Dispensed 0.003904 mL of Acid (0.5 M HCI)
- [59:40] Stirrer speed set to 0
- [59:50] Datapoint id 63 collected
- [59:50] Stirrer speed set to 60
- [59:55] pH 1.98 -> 2.18
- [59:55] Using cautious pH adjust
- [59:56] Dispensed 0.011265 mL of Base (0.5 M KOH)
- [1:00:01] Stepping pH = 2.06
- [1:00:01] Dispensed 0.008796 mL of Base (0.5 M KOH)
- [1:00:06] Stepping pH = 2.15
- [1:00:06] Dispensed 0.001881 mL of Base (0.5 M KOH)
- [1:00:11] Stepping pH = 2.18
- [1:00:26] Stirrer speed set to 0
- [1:00:47] Datapoint id 64 collected
- [1:00:47] Charge balance equation is out by 2.6%
- [1:00:47] Stirrer speed set to 60
- [1:00:52] pH 2.19 -> 2.39
- [1:00:52] Using charge balance adjust
- [1:00:52] Dispensed 0.013805 mL of Base (0.5 M KOH)
- [1:01:12] Stirrer speed set to 0
- [1:01:22] Datapoint id 65 collected
- [1:01:22] Charge balance equation is out by -4.3%
- [1:01:22] Stirrer speed set to 60
- [1:01:27] pH 2.39 -> 2.59
- [1:01:27] Using charge balance adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [1:01:28] Dispensed 0.008866 mL of Base (0.5 M KOH)
- [1:01:48] Stirrer speed set to 0
- [1:01:58] Datapoint id 66 collected
- [1:01:58] Charge balance equation is out by -13.7%
- [1:01:58] Stirrer speed set to 60
- [1:02:03] pH 2.57 -> 2.77
- [1:02:03] Using charge balance adjust
- [1:02:03] Dispensed 0.006021 mL of Base (0.5 M KOH)
- [1:02:23] Stirrer speed set to 0
- [1:02:34] Datapoint id 67 collected
- [1:02:34] Charge balance equation is out by -21.7%
- [1:02:34] Stirrer speed set to 60
- [1:02:39] pH 2.73 -> 2.93
- [1:02:39] Using cautious pH adjust
- [1:02:39] Dispensed 0.002187 mL of Base (0.5 M KOH)
- [1:02:44] Stepping pH = 2.82
- [1:02:44] Dispensed 0.001834 mL of Base (0.5 M KOH)
- [1:02:49] Stepping pH = 2.89
- [1:02:49] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [1:02:54] Stepping pH = 2.92
- [1:03:09] Stirrer speed set to 0
- [1:03:27] Datapoint id 68 collected
- [1:03:27] Charge balance equation is out by -9.5%
- [1:03:27] Stirrer speed set to 60
- [1:03:33] pH 2.93 -> 3.13
- [1:03:33] Using charge balance adjust
- [1:03:33] Dispensed 0.003175 mL of Base (0.5 M KOH)
- [1:03:53] Stirrer speed set to 0
- [1:04:03] Datapoint id 69 collected
- [1:04:03] Charge balance equation is out by -12.8%
- [1:04:03] Stirrer speed set to 60
- [1:04:09] pH 3.11 -> 3.31
- [1:04:09] Using charge balance adjust
- [1:04:09] Dispensed 0.002587 mL of Base (0.5 M KOH)
- [1:04:29] Stirrer speed set to 0
- [1:04:40] Datapoint id 70 collected
- [1:04:40] Charge balance equation is out by 10.9%
- [1:04:40] Stirrer speed set to 60
- [1:04:45] pH 3.34 -> 3.54
- [1:04:45] Using charge balance adjust
- [1:04:45] Dispensed 0.002305 mL of Base (0.5 M KOH)
- [1:05:05] Stirrer speed set to 0
- [1:05:15] Datapoint id 71 collected
- [1:05:15] Charge balance equation is out by 48.3%
- [1:05:15] Stirrer speed set to 60
- [1:05:20] pH 3.64 -> 3.84
- [1:05:20] Using cautious pH adjust
- [1:05:20] Dispensed 0.001152 mL of Base (0.5 M KOH)
- [1:05:25] Stepping pH = 3.89
- [1:05:40] Stirrer speed set to 0
- [1:05:51] Datapoint id 72 collected
- [1:05:51] Charge balance equation is out by 50.0%
- [1:05:51] Stirrer speed set to 60
- [1:05:56] pH 3.89 -> 4.09
- [1:05:56] Using cautious pH adjust
- [1:05:56] Dispensed 0.001105 mL of Base (0.5 M KOH)
- [1:06:01] Stepping pH = 4.44
- [1:06:16] Stirrer speed set to 0
- [1:06:26] Datapoint id 73 collected



Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [1:06:26] Charge balance equation is out by 50.0%
- [1:06:26] Stirrer speed set to 60
- [1:06:32] pH 4.44 -> 4.64
- [1:06:32] Using cautious pH adjust
- [1:06:32] Dispensed 0.000706 mL of Base (0.5 M KOH)
- [1:06:37] Stepping pH = 6.57
- [1:06:52] Stirrer speed set to 0
- [1:07:47] Datapoint id 74 collected
- [1:07:47] Charge balance equation is out by 50.0%
- [1:07:47] Stirrer speed set to 60
- [1:07:53] pH 6.31 -> 6.51
- [1:07:53] Using cautious pH adjust
- [1:07:53] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:07:58] Stepping pH = 6.28
- [1:07:58] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [1:08:03] Stepping pH = 8.21
- [1:08:18] Stirrer speed set to 0
- [1:09:18] Datapoint id 75 collected
- [1:09:18] Charge balance equation is out by -228.6%
- [1:09:18] Stirrer speed set to 60
- [1:09:23] pH 7.81 -> 8.01
- [1:09:23] Using cautious pH adjust
- [1:09:23] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:09:28] Stepping pH = 7.85
- [1:09:28] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:09:33] Stepping pH = 8.33
- [1:09:48] Stirrer speed set to 0
- [1:10:33] Datapoint id 76 collected
- [1:10:33] Charge balance equation is out by -40.8%
- [1:10:33] Stirrer speed set to 60
- [1:10:38] pH 8.20 -> 8.40
- [1:10:38] Using cautious pH adjust
- [1:10:38] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [1:10:43] Stepping pH = 8.56
- [1:10:58] Stirrer speed set to 0
- [1:11:21] Datapoint id 77 collected
- [1:11:21] Charge balance equation is out by 50.0%
- [1:11:21] Stirrer speed set to 60
- [1:11:26] pH 8.49 -> 8.69
- [1:11:26] Using cautious pH adjust
- [1:11:26] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:11:31] Stepping pH = 8.91
- [1:11:46] Stirrer speed set to 0
- [1:11:58] Datapoint id 78 collected
- [1:11:58] Charge balance equation is out by 50.0%
- [1:11:58] Stirrer speed set to 60
- [1:12:03] pH 8.87 -> 9.07
- [1:12:03] Using cautious pH adjust
- [1:12:03] Dispensed 0.000400 mL of Base (0.5 M KOH)
- [1:12:09] Stepping pH = 9.33
- [1:12:24] Stirrer speed set to 0
- [1:12:35] Datapoint id 79 collected
- [1:12:35] Charge balance equation is out by 50.0%
- [1:12:35] Stirrer speed set to 60
- [1:12:40] pH 9.30 -> 9.50
- [1:12:40] Using cautious pH adjust [1:12:40] Dispensed 0.000800 mL of Base (0.5 M KOH)
- [1:12:45] Stepping pH = 9.79
- [1:13:00] Stirrer speed set to 0



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [1:13:10] Datapoint id 80 collected
- [1:13:10] Charge balance equation is out by 50.0%
- [1:13:10] Stirrer speed set to 60
- [1:13:15] pH 9.77 -> 9.97
- [1:13:15] Using cautious pH adjust
- [1:13:15] Dispensed 0.001105 mL of Base (0.5 M KOH)
- [1:13:20] Stepping pH = 10.15
- [1:13:36] Stirrer speed set to 0
- [1:13:46] Datapoint id 81 collected
- [1:13:46] Charge balance equation is out by 50.0%
- [1:13:46] Stirrer speed set to 60
- [1:13:51] pH 10.14 -> 10.34
- [1:13:51] Using cautious pH adjust
- [1:13:51] Dispensed 0.001082 mL of Base (0.5 M KOH)
- 1:13:56 Stepping pH = 10.37
- [1:14:11] Stirrer speed set to 0
- 1.4 4.041 Determint id 00 celler
- [1:14:21] Datapoint id 82 collected
- [1:14:21] Charge balance equation is out by 50.0%
- [1:14:21] Stirrer speed set to 60
- [1:14:26] pH 10.36 -> 10.56
- [1:14:26] Using cautious pH adjust
- [1:14:26] Dispensed 0.001058 mL of Base (0.5 M KOH)
- [1:14:31] Stepping pH = 10.51
- [1:14:31] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [1:14:37] Stepping pH = 10.55
- [1:14:37] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [1:14:42] Stepping pH = 10.55
- [1:14:57] Stirrer speed set to 0
- [1:15:07] Datapoint id 83 collected
- [1:15:07] Charge balance equation is out by 29.7%
- [1:15:07] Stirrer speed set to 60
- [1:15:12] pH 10.55 -> 10.75
- [1:15:12] Using cautious pH adjust
- [1:15:12] Dispensed 0.001105 mL of Base (0.5 M KOH)
- [1:15:17] Stepping pH = 10.66
- [1:15:17] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [1:15:22] Stepping pH = 10.71
- [1:15:22] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [1:15:27] Stepping pH = 10.73
- [1:15:28] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:15:33] Stepping pH = 10.74
- [1:15:48] Stirrer speed set to 0
- [1:15:58] Datapoint id 84 collected
- [1:15:58] Charge balance equation is out by -4.8%
- [1:15:58] Stirrer speed set to 60
- [1:16:03] pH 10.74 -> 10.94
- [1:16:03] Using charge balance adjust
- [1:16:03] Dispensed 0.002611 mL of Base (0.5 M KOH)
- [1:16:23] Stirrer speed set to 0
- [1:16:33] Datapoint id 85 collected
- [1:16:33] Charge balance equation is out by -12.8%
- [1:16:33] Stirrer speed set to 60
- [1:16:38] pH 10.92 -> 11.12
- [1:16:38] Using charge balance adjust
- [1:16:39] Dispensed 0.003387 mL of Base (0.5 M KOH)
- [1:16:59] Stirrer speed set to 0
- [1:17:09] Datapoint id 86 collected
- [1:17:09] Charge balance equation is out by -11.7%
- [1:17:09] Stirrer speed set to 60



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01009 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180228_exp28_logP_T3-2\18C-01009_M16_octanol_pH-metric high logP.t3r

- [1:17:14] pH 11.10 -> 11.30
- [1:17:14] Using charge balance adjust
- [1:17:14] Dispensed 0.004751 mL of Base (0.5 M KOH)
- [1:17:34] Stirrer speed set to 0
- [1:17:44] Datapoint id 87 collected
- [1:17:44] Charge balance equation is out by -27.1%
- [1:17:44] Stirrer speed set to 60
- [1:17:50] pH 11.25 -> 11.45
- [1:17:50] Using cautious pH adjust
- [1:17:50] Dispensed 0.003269 mL of Base (0.5 M KOH)
- [1:17:55] Stepping pH = 11.34
- [1:17:55] Dispensed 0.002893 mL of Base (0.5 M KOH)
- [1:18:00] Stepping pH = 11.41
- [1:18:00] Dispensed 0.001364 mL of Base (0.5 M KOH)
- [1:18:05] Stepping pH = 11.44
- [1:18:05] Dispensed 0.000470 mL of Base (0.5 M KOH)
- [1:18:10] Stepping pH = 11.45
- [1:18:26] Stirrer speed set to 0
- [1:18:36] Datapoint id 88 collected
- [1:18:36] Charge balance equation is out by -22.0%
- [1:18:36] Stirrer speed set to 60
- [1:18:41] pH 11.45 -> 11.65
- [1:18:41] Using cautious pH adjust
- [1:18:41] Dispensed 0.005103 mL of Base (0.5 M KOH)
- [1:18:46] Stepping pH = 11.54
- [1:18:46] Dispensed 0.004163 mL of Base (0.5 M KOH)
- [1:18:51] Stepping pH = 11.61
- [1:18:51] Dispensed 0.002023 mL of Base (0.5 M KOH)
- [1:18:57] Stepping pH = 11.64
- [1:19:12] Stirrer speed set to 0
- [1:19:34] Datapoint id 89 collected
- [1:19:34] Charge balance equation is out by -10.3%
- [1:19:34] Stirrer speed set to 60
- [1:19:39] pH 11.64 -> 11.84
- [1:19:39] Using charge balance adjust
- [1:19:40] Dispensed 0.015851 mL of Base (0.5 M KOH)
- [1:20:00] Stirrer speed set to 0
- [1:20:10] Datapoint id 90 collected
- [1:20:10] Charge balance equation is out by -12.5%
- [1:20:10] Stirrer speed set to 60
- [1:20:15] pH 11.82 -> 12.02
- [1:20:15] Using charge balance adjust
- [1:20:16] Dispensed 0.024506 mL of Base (0.5 M KOH)
- [1:20:36] Stirrer speed set to 0
- [1:20:46] Datapoint id 91 collected
- [1:20:46] Charge balance equation is out by -12.2%
- [1:20:46] Stirrer speed set to 60
- [1:20:51] pH 12.01 -> 12.05
- [1:20:51] Using charge balance adjust
- [1:20:51] Dispensed 0.007103 mL of Base (0.5 M KOH)
- [1:21:12] Stirrer speed set to 0
- [1:21:22] Datapoint id 92 collected
- [1:21:22] Charge balance equation is out by -85.4%
- [1:21:22] Argon flow rate set to 0
- [1:21:26] Titrator arm moved over Titration position