

Sample name: **M15**
 Assay name: **UV-metric psKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Levorse**
 Instrument ID: **T311053**




Yasuda-Shedlovsky result

Extrapolation type	pKa 0%	SD	Intercept	Slope	R ²	Ionic strength	Temperature
Yasuda-Shedlovsky	2.56	±0.06	7.65	-262.9812	0.9959	0.165 M	24.9°C
Yasuda-Shedlovsky	5.29	±0.06	9.46	-190.6287	0.9935	0.165 M	24.9°C

Component assay results

Titration	Methanol weight%	Direction	Result type	Dielectric constant	[H2O]	Ionic strength	Temperature	psKa 1	psKa 2
17J-12008 Points 4 to 39	49.51 %	Up	UV-metric pKa	56.6	24.7 M	0.157 M	24.9°C	✓	✓
17J-12008 Points 41 to 77	40.06 %	Up	UV-metric pKa	61.0	30.0 M	0.166 M	24.9°C	✓	✓
17J-12008 Points 79 to 116	30.24 %	Up	UV-metric pKa	65.4	35.7 M	0.172 M	24.9°C	✓	✓

Graphs



UV-metric psKa Titration 1 of 3 17J-12008 Points 4 to 39

Results

pKa 1 **1.62**
 pKa 2 **4.71**
 RMSD **0.003 0.003 0.005**
 Chi squared **0.0250**
 PCA calculated number of pKas **3**
 Average ionic strength **0.157 M**
 Average temperature **24.9°C**
 Analyte concentration range **28.5 µM to 26.8 µM**
 Methanol weight % **49.5 %**
 Dielectric constant **56.6**
 Water concentration **24.7 M**

Sample name: **M15**
 Assay name: **UV-metric psKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Levorse**
 Instrument ID: **T311053**

Results (continued)

Number of pKas source **Predicted**
 Wavelength clipping **230.0 nm to 450.0 nm**
 pH clipping **1.467 to 12.528**


Warnings and errors

Errors None
 Warnings PCA calculation disagrees with predicted number of pKas

Assay Settings

Setting	Value	Original Value	Date/Time changed	Imported from
Buffer in use	Yes			
Buffer type	Phosphate Buffer			
Assay Medium				
Volume of buffer introduced	0.025000 mL			
Add buffer manually	Manual			


Graphs



Sample name: **M15**
 Assay name: **UV-metric psKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Levorse**
 Instrument ID: **T311053**

Graphs (continued)



UV-metric psKa Titration 2 of 3 17J-12008 Points 41 to 77

Results

pKa 1 **1.84**
 pKa 2 **4.83**
 RMSD **0.004 0.005 0.003**
 Chi squared **0.0275**
 PCA calculated number of pKas **3**
 Average ionic strength **0.166 M**
 Average temperature **24.9°C**
 Analyte concentration range **23.4 µM to 22.1 µM**
 Methanol weight % **40.1 %**
 Dielectric constant **61.0**
 Water concentration **30.0 M**

Sample name: **M15**
 Assay name: **UV-metric pKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric pKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Leverage**
 Instrument ID: **T311053**

Results (continued)

Number of pKas source **Predicted**
 Wavelength clipping **230.0 nm to 450.0 nm**
 pH clipping **1.510 to 12.508**


Warnings and errors

Errors **None**
 Warnings **PCA calculation disagrees with predicted number of pKas**

Assay Settings

Setting	Value	Original Value	Date/Time changed	Imported from
Buffer in use	Yes			
Buffer type	Phosphate Buffer			
Assay Medium				
Volume of buffer introduced	0.025000 mL			
Add buffer manually	Manual			


Graphs



Sample name: **M15**
 Assay name: **UV-metric psKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Levorse**
 Instrument ID: **T311053**

Graphs (continued)



UV-metric psKa Titration 3 of 3 17J-12008 Points 79 to 116

Results

pKa 1 **2.09**
 pKa 2 **5.00**
 RMSD **0.006 0.009 0.004**
 Chi squared **0.0538**
 PCA calculated number of pKas **3**
 Average ionic strength **0.172 M**
 Average temperature **24.9°C**
 Analyte concentration range **18.0 µM to 17.1 µM**
 Methanol weight % **30.2 %**
 Dielectric constant **65.4**
 Water concentration **35.7 M**

Sample name: **M15**
 Assay name: **UV-metric pKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric pKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Levorse**
 Instrument ID: **T311053**

Results (continued)

Number of pKas source **Predicted**
 Wavelength clipping **230.0 nm to 450.0 nm**
 pH clipping **1.516 to 12.527**


Warnings and errors

Errors None
 Warnings PCA calculation disagrees with predicted number of pKas

Assay Settings

Setting	Value	Original Value	Date/Time changed	Imported from
Buffer in use	Yes			
Buffer type	Phosphate Buffer			
Assay Medium				
Volume of buffer introduced	0.025000 mL			
Add buffer manually	Manual			


Graphs



Sample name: **M15**
 Assay name: **UV-metric psKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Levorse**
 Instrument ID: **T311053**

Graphs (continued)



Assay Model

Settings

Sample name	M15	10/11/2017 4:19:59 PM	User entered value
Sample by	Volume		Default value
Sample volume	0.0010 mL	10/11/2017 4:19:59 PM	User entered value
Solvent	DMSO		Default value
Sample concentration	0.045400 M	10/11/2017 4:19:59 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	209.25	10/11/2017 4:20:11 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	2	10/11/2017 4:19:59 PM	User entered value
Sample is a	Base	10/11/2017 4:19:59 PM	User entered value
pKa 1	2.94	10/11/2017 4:19:59 PM	User entered value
pKa 2	5.25	10/11/2017 4:19:59 PM	User entered value
logP (XH2 2+)	-10.00		Default value
logp (XH +)	-10.00		Default value

Sample name: **M15**
 Assay name: **UV-metric psKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Levorse**
 Instrument ID: **T311053**

Assay Model (continued)

Settings	Value	Date/Time changed	Imported from
logP (neutral X)	-10.00	10/11/2017 4:19:59 PM	User entered value

Events

Time	Event	Water	Acid	Base	Methanol	Buffer	pH	dpH/dt	pH R-squared
3:31.4	Dark spectrum								
3:32.8	Reference spectrum								
4:00.4	Volume reset due to vial change								
4:44.6	Initial pH = 8.34								
5:50.1	Data point 4	0.34995 mL	0.06874 mL	0.00000 mL	1.15005 mL	0.02500 mL	1.967	-0.00525	0.44790
6:19.0	Data point 5	0.34995 mL	0.06874 mL	0.02491 mL	1.15005 mL	0.02500 mL	2.169	-0.00954	0.41836
6:35.8	Data point 6	0.34995 mL	0.06874 mL	0.03996 mL	1.15005 mL	0.02500 mL	2.352	0.01979	0.90627
6:52.7	Data point 7	0.34995 mL	0.06874 mL	0.04976 mL	1.15005 mL	0.02500 mL	2.545	0.02536	0.91550
7:09.4	Data point 8	0.34995 mL	0.06874 mL	0.05614 mL	1.15005 mL	0.02500 mL	2.755	0.01329	0.83619
7:26.1	Data point 9	0.34995 mL	0.06874 mL	0.06011 mL	1.15005 mL	0.02500 mL	2.944	0.01099	0.84718
7:42.7	Data point 10	0.34995 mL	0.06874 mL	0.06270 mL	1.15005 mL	0.02500 mL	3.150	0.00700	0.76597
7:59.4	Data point 11	0.34995 mL	0.06874 mL	0.06432 mL	1.15005 mL	0.02500 mL	3.342	0.00970	0.82302
8:15.9	Data point 12	0.34995 mL	0.06874 mL	0.06536 mL	1.15005 mL	0.02500 mL	3.510	0.01419	0.84907
8:47.9	Data point 13	0.34995 mL	0.06874 mL	0.06637 mL	1.15005 mL	0.02500 mL	3.713	0.02970	0.98099
9:25.1	Data point 14	0.34995 mL	0.06874 mL	0.06729 mL	1.15005 mL	0.02500 mL	3.910	0.03689	0.97433
9:46.8	Data point 15	0.34995 mL	0.06874 mL	0.06780 mL	1.15005 mL	0.02500 mL	4.300	0.07909	0.98985
10:13.5	Data point 16	0.34995 mL	0.06874 mL	0.06823 mL	1.15005 mL	0.02500 mL	5.353	0.09957	0.99392
11:21.8	Data point 17	0.34995 mL	0.06874 mL	0.06841 mL	1.15005 mL	0.02500 mL	6.034	0.13850	0.99424
12:38.5	Data point 18	0.34995 mL	0.06874 mL	0.06851 mL	1.15005 mL	0.02500 mL	6.608	0.09912	0.99017
13:50.4	Data point 19	0.34995 mL	0.06874 mL	0.06863 mL	1.15005 mL	0.02500 mL	7.042	0.10030	0.99113
14:43.6	Data point 20	0.34995 mL	0.06874 mL	0.06874 mL	1.15005 mL	0.02500 mL	7.340	0.09933	0.99340
15:31.6	Data point 21	0.34995 mL	0.06874 mL	0.06886 mL	1.15005 mL	0.02500 mL	7.598	0.09999	0.99204
16:15.0	Data point 22	0.34995 mL	0.06874 mL	0.06900 mL	1.15005 mL	0.02500 mL	7.863	0.09906	0.99051
17:05.4	Data point 23	0.34995 mL	0.06874 mL	0.06914 mL	1.15005 mL	0.02500 mL	8.174	0.09949	0.99270
17:55.3	Data point 24	0.34995 mL	0.06874 mL	0.06926 mL	1.15005 mL	0.02500 mL	8.445	0.09598	0.90762
18:31.8	Data point 25	0.34995 mL	0.06874 mL	0.06936 mL	1.15005 mL	0.02500 mL	8.783	0.09877	0.96284
19:21.0	Data point 26	0.34995 mL	0.06874 mL	0.06945 mL	1.15005 mL	0.02500 mL	9.120	0.09548	0.96493
20:12.5	Data point 27	0.34995 mL	0.06874 mL	0.06957 mL	1.15005 mL	0.02500 mL	9.500	0.09442	0.94739
20:53.2	Data point 28	0.34995 mL	0.06874 mL	0.06968 mL	1.15005 mL	0.02500 mL	9.802	0.09795	0.96406
21:26.0	Data point 29	0.34995 mL	0.06874 mL	0.06983 mL	1.15005 mL	0.02500 mL	10.047	0.09929	0.96654
21:48.5	Data point 30	0.34995 mL	0.06874 mL	0.07001 mL	1.15005 mL	0.02500 mL	10.298	0.02890	0.82114
22:20.3	Data point 31	0.34995 mL	0.06874 mL	0.07030 mL	1.15005 mL	0.02500 mL	10.497	0.01960	0.93055
22:37.0	Data point 32	0.34995 mL	0.06874 mL	0.07072 mL	1.15005 mL	0.02500 mL	10.777	0.00296	0.24851
23:08.8	Data point 33	0.34995 mL	0.06874 mL	0.07166 mL	1.15005 mL	0.02500 mL	10.976	-0.00336	0.40552
23:41.0	Data point 34	0.34995 mL	0.06874 mL	0.07314 mL	1.15005 mL	0.02500 mL	11.169	-0.01019	0.82632
24:13.2	Data point 35	0.34995 mL	0.06874 mL	0.07563 mL	1.15005 mL	0.02500 mL	11.359	-0.00614	0.72743
24:40.3	Data point 36	0.34995 mL	0.06874 mL	0.07949 mL	1.15005 mL	0.02500 mL	11.558	-0.01212	0.86864
25:12.5	Data point 37	0.34995 mL	0.06874 mL	0.08544 mL	1.15005 mL	0.02500 mL	11.754	-0.00675	0.73695
25:45.0	Data point 38	0.34995 mL	0.06874 mL	0.09544 mL	1.15005 mL	0.02500 mL	11.945	-0.00910	0.76255
26:06.9	Data point 39	0.34995 mL	0.06874 mL	0.10085 mL	1.15005 mL	0.02500 mL	12.028	-0.00626	0.60483
27:48.4	Reference spectrum								
28:52.4	Data point 41	0.50000 mL	0.16634 mL	0.10087 mL	1.15005 mL	0.02500 mL	2.010	-0.05068	0.89643
29:19.9	Data point 42	0.50000 mL	0.16634 mL	0.12429 mL	1.15005 mL	0.02500 mL	2.208	0.00950	0.75308
29:36.9	Data point 43	0.50000 mL	0.16634 mL	0.13930 mL	1.15005 mL	0.02500 mL	2.410	0.01380	0.69492
29:53.8	Data point 44	0.50000 mL	0.16634 mL	0.14880 mL	1.15005 mL	0.02500 mL	2.617	0.00757	0.36152
30:10.6	Data point 45	0.50000 mL	0.16634 mL	0.15475 mL	1.15005 mL	0.02500 mL	2.816	0.02783	0.86319
30:27.4	Data point 46	0.50000 mL	0.16634 mL	0.15849 mL	1.15005 mL	0.02500 mL	3.023	0.00277	0.13850
30:44.1	Data point 47	0.50000 mL	0.16634 mL	0.16082 mL	1.15005 mL	0.02500 mL	3.227	0.01696	0.91503
31:00.8	Data point 48	0.50000 mL	0.16634 mL	0.16228 mL	1.15005 mL	0.02500 mL	3.402	0.01558	0.90902
31:17.4	Data point 49	0.50000 mL	0.16634 mL	0.16326 mL	1.15005 mL	0.02500 mL	3.572	0.02207	0.96804
31:33.9	Data point 50	0.50000 mL	0.16634 mL	0.16392 mL	1.15005 mL	0.02500 mL	3.664	0.02780	0.95210
31:55.6	Data point 51	0.50000 mL	0.16634 mL	0.16477 mL	1.15005 mL	0.02500 mL	3.975	0.03761	0.97378

Sample name: **M15**
 Assay name: **UV-metric psKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Levorse**
 Instrument ID: **T311053**

Events (continued)

Time	Event	Water	Acid	Base	Methanol	Buffer	pH	dpH/dt	pH R-squared	pH Slope
32:17.4	Data point 52	0.50000 mL	0.16634 mL	0.16526 mL	1.15005 mL	0.02500 mL	4.386	0.08947	0.99208	0.99208
32:44.2	Data point 53	0.50000 mL	0.16634 mL	0.16562 mL	1.15005 mL	0.02500 mL	5.005	0.10041	0.99004	0.99004
33:34.0	Data point 54	0.50000 mL	0.16634 mL	0.16583 mL	1.15005 mL	0.02500 mL	6.003	0.09781	0.99468	0.99468
34:37.8	Data point 55	0.50000 mL	0.16634 mL	0.16597 mL	1.15005 mL	0.02500 mL	6.472	0.09955	0.99445	0.99445
35:26.0	Data point 56	0.50000 mL	0.16634 mL	0.16606 mL	1.15005 mL	0.02500 mL	6.735	0.09893	0.98708	0.98708
36:04.1	Data point 57	0.50000 mL	0.16634 mL	0.16616 mL	1.15005 mL	0.02500 mL	6.978	0.09618	0.98868	0.98868
36:40.5	Data point 58	0.50000 mL	0.16634 mL	0.16628 mL	1.15005 mL	0.02500 mL	7.260	0.10005	0.98767	0.98767
37:14.8	Data point 59	0.50000 mL	0.16634 mL	0.16639 mL	1.15005 mL	0.02500 mL	7.530	0.09972	0.98386	0.98386
37:55.7	Data point 60	0.50000 mL	0.16634 mL	0.16651 mL	1.15005 mL	0.02500 mL	7.824	0.10067	0.99226	0.99226
38:36.6	Data point 61	0.50000 mL	0.16634 mL	0.16660 mL	1.15005 mL	0.02500 mL	8.083	0.09648	0.97074	0.97074
39:20.7	Data point 62	0.50000 mL	0.16634 mL	0.16670 mL	1.15005 mL	0.02500 mL	8.412	0.09956	0.98382	0.98382
40:10.9	Data point 63	0.50000 mL	0.16634 mL	0.16679 mL	1.15005 mL	0.02500 mL	8.863	0.09764	0.98633	0.98633
40:57.2	Data point 64	0.50000 mL	0.16634 mL	0.16686 mL	1.15005 mL	0.02500 mL	9.224	0.09288	0.94299	0.94299
41:37.4	Data point 65	0.50000 mL	0.16634 mL	0.16693 mL	1.15005 mL	0.02500 mL	9.534	0.10038	0.98616	0.98616
42:07.1	Data point 66	0.50000 mL	0.16634 mL	0.16700 mL	1.15005 mL	0.02500 mL	9.770	0.09351	0.96224	0.96224
42:32.7	Data point 67	0.50000 mL	0.16634 mL	0.16712 mL	1.15005 mL	0.02500 mL	10.020	0.06284	0.97651	0.97651
42:49.2	Data point 68	0.50000 mL	0.16634 mL	0.16731 mL	1.15005 mL	0.02500 mL	10.275	0.01695	0.81940	0.81940
43:20.9	Data point 69	0.50000 mL	0.16634 mL	0.16773 mL	1.15005 mL	0.02500 mL	10.470	0.00881	0.64616	0.64616
43:53.0	Data point 70	0.50000 mL	0.16634 mL	0.16834 mL	1.15005 mL	0.02500 mL	10.666	0.00308	0.32860	0.32860
44:09.6	Data point 71	0.50000 mL	0.16634 mL	0.16919 mL	1.15005 mL	0.02500 mL	10.865	-0.00524	0.65148	0.65148
44:26.2	Data point 72	0.50000 mL	0.16634 mL	0.17053 mL	1.15005 mL	0.02500 mL	11.033	-0.00879	0.79184	0.79184
44:58.4	Data point 73	0.50000 mL	0.16634 mL	0.17317 mL	1.15005 mL	0.02500 mL	11.225	-0.00979	0.78876	0.78876
45:25.5	Data point 74	0.50000 mL	0.16634 mL	0.17681 mL	1.15005 mL	0.02500 mL	11.418	-0.01504	0.90841	0.90841
46:03.1	Data point 75	0.50000 mL	0.16634 mL	0.18380 mL	1.15005 mL	0.02500 mL	11.626	-0.01108	0.62860	0.62860
46:35.2	Data point 76	0.50000 mL	0.16634 mL	0.19365 mL	1.15005 mL	0.02500 mL	11.822	-0.01119	0.87568	0.87568
47:02.5	Data point 77	0.50000 mL	0.16634 mL	0.20884 mL	1.15005 mL	0.02500 mL	12.008	-0.00839	0.82304	0.82304
48:47.3	Reference spectrum									
50:10.6	Data point 79	0.83996 mL	0.29523 mL	0.20887 mL	1.15005 mL	0.02500 mL	2.016	-0.01620	0.88855	0.88855
50:38.2	Data point 80	0.83996 mL	0.29523 mL	0.23544 mL	1.15005 mL	0.02500 mL	2.210	0.01120	0.73558	0.73558
50:55.3	Data point 81	0.83996 mL	0.29523 mL	0.25249 mL	1.15005 mL	0.02500 mL	2.406	0.00019	0.00022	0.00022
51:12.1	Data point 82	0.83996 mL	0.29523 mL	0.26352 mL	1.15005 mL	0.02500 mL	2.600	-0.01948	0.79039	0.79039
51:28.9	Data point 83	0.83996 mL	0.29523 mL	0.27067 mL	1.15005 mL	0.02500 mL	2.802	-0.00148	0.04503	0.04503
51:45.6	Data point 84	0.83996 mL	0.29523 mL	0.27514 mL	1.15005 mL	0.02500 mL	2.990	-0.01144	0.54674	0.54674
52:02.2	Data point 85	0.83996 mL	0.29523 mL	0.27803 mL	1.15005 mL	0.02500 mL	3.209	-0.01542	0.69490	0.69490
52:18.8	Data point 86	0.83996 mL	0.29523 mL	0.27980 mL	1.15005 mL	0.02500 mL	3.383	-0.00720	0.38977	0.38977
52:35.4	Data point 87	0.83996 mL	0.29523 mL	0.28097 mL	1.15005 mL	0.02500 mL	3.538	-0.00290	0.08907	0.08907
53:07.5	Data point 88	0.83996 mL	0.29523 mL	0.28140 mL	1.15005 mL	0.02500 mL	3.830	-0.01232	0.83931	0.83931
53:34.5	Data point 89	0.83996 mL	0.29523 mL	0.28429 mL	1.15005 mL	0.02500 mL	4.972	-0.08188	0.86107	0.86107
54:03.1	Data point 90	0.83996 mL	0.29523 mL	0.28457 mL	1.15005 mL	0.02500 mL	5.781	0.01061	0.01674	0.01674
54:25.4	Data point 91	0.83996 mL	0.29523 mL	0.28471 mL	1.15005 mL	0.02500 mL	6.080	-0.05291	0.40901	0.40901
54:47.5	Data point 92	0.83996 mL	0.29523 mL	0.28481 mL	1.15005 mL	0.02500 mL	6.283	-0.06378	0.75328	0.75328
55:09.2	Data point 93	0.83996 mL	0.29523 mL	0.28492 mL	1.15005 mL	0.02500 mL	6.535	-0.07363	0.74563	0.74563
55:35.7	Data point 94	0.83996 mL	0.29523 mL	0.28504 mL	1.15005 mL	0.02500 mL	6.755	-0.09453	0.90047	0.90047
56:06.4	Data point 95	0.83996 mL	0.29523 mL	0.28516 mL	1.15005 mL	0.02500 mL	6.966	0.01849	0.44080	0.44080
56:33.3	Data point 96	0.83996 mL	0.29523 mL	0.28528 mL	1.15005 mL	0.02500 mL	7.165	0.05052	0.90661	0.90661
57:05.2	Data point 97	0.83996 mL	0.29523 mL	0.28542 mL	1.15005 mL	0.02500 mL	7.406	0.09582	0.94132	0.94132
57:37.1	Data point 98	0.83996 mL	0.29523 mL	0.28554 mL	1.15005 mL	0.02500 mL	7.646	0.09448	0.95830	0.95830
58:11.7	Data point 99	0.83996 mL	0.29523 mL	0.28563 mL	1.15005 mL	0.02500 mL	7.893	0.09904	0.97936	0.97936
58:51.1	Data point 100	0.83996 mL	0.29523 mL	0.28572 mL	1.15005 mL	0.02500 mL	8.200	0.09751	0.97855	0.97855
59:35.9	Data point 101	0.83996 mL	0.29523 mL	0.28582 mL	1.15005 mL	0.02500 mL	8.648	0.09789	0.95408	0.95408
1:00:19.7	Data point 102	0.83996 mL	0.29523 mL	0.28589 mL	1.15005 mL	0.02500 mL	8.948	0.09695	0.94482	0.94482
1:01:00.0	Data point 103	0.83996 mL	0.29523 mL	0.28598 mL	1.15005 mL	0.02500 mL	9.292	0.09665	0.97234	0.97234
1:01:30.1	Data point 104	0.83996 mL	0.29523 mL	0.28608 mL	1.15005 mL	0.02500 mL	9.557	0.09671	0.97597	0.97597
1:01:55.9	Data point 105	0.83996 mL	0.29523 mL	0.28619 mL	1.15005 mL	0.02500 mL	9.823	0.05116	0.98408	0.98408
1:02:12.4	Data point 106	0.83996 mL	0.29523 mL	0.28638 mL	1.15005 mL	0.02500 mL	10.102	-0.00513	0.37093	0.37093
1:02:44.3	Data point 107	0.83996 mL	0.29523 mL	0.28676 mL	1.15005 mL	0.02500 mL	10.299	0.00942	0.46416	0.46416
1:03:00.9	Data point 108	0.83996 mL	0.29523 mL	0.28735 mL	1.15005 mL	0.02500 mL	10.517	-0.02188	0.89061	0.89061

Sample name: **M15**
 Assay name: **UV-metric psKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Levorse**
 Instrument ID: **T311053**

Events (continued)

Time	Event	Water	Acid	Base	Methanol	Buffer	pH	dpH/dt	pH R-squared	pH SD
1:03:17.3	Data point 109	0.83996 mL	0.29523 mL	0.28829 mL	1.15005 mL	0.02500 mL	10.700	-0.02234	0.91518	0.0011
1:03:34.0	Data point 110	0.83996 mL	0.29523 mL	0.28970 mL	1.15005 mL	0.02500 mL	10.875	-0.02240	0.88731	0.0011
1:03:50.7	Data point 111	0.83996 mL	0.29523 mL	0.29177 mL	1.15005 mL	0.02500 mL	11.040	-0.02425	0.93339	0.0012
1:04:17.7	Data point 112	0.83996 mL	0.29523 mL	0.29565 mL	1.15005 mL	0.02500 mL	11.240	-0.01758	0.89791	0.0009
1:04:55.2	Data point 113	0.83996 mL	0.29523 mL	0.30287 mL	1.15005 mL	0.02500 mL	11.461	-0.01669	0.89086	0.0008
1:05:27.6	Data point 114	0.83996 mL	0.29523 mL	0.31244 mL	1.15005 mL	0.02500 mL	11.653	-0.01628	0.88701	0.0008
1:05:54.9	Data point 115	0.83996 mL	0.29523 mL	0.32733 mL	1.15005 mL	0.02500 mL	11.844	-0.01423	0.90910	0.0007
1:06:22.6	Data point 116	0.83996 mL	0.29523 mL	0.35089 mL	1.15005 mL	0.02500 mL	12.027	-0.01514	0.88771	0.0007
1:08:22.0	Assay volumes	1.08996 mL	0.43822 mL	0.35089 mL	1.15005 mL	0.02500 mL				

Assay Settings

Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Dorothy Levorse			
Separate reference vial	Yes			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	12.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings				
Detect turbidity using	Spectrometer			
Monitor at a wavelength of	500.0 nm			
Absorbance threshold of	0.100			
Collect turbidity sensor data	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	15%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium				
Cosolvent in use	Yes			
Cosolvent type	Methanol			
Cosolvent volume	1.15 mL			
Cosolvent added	Automatic			
ISA water volume	0.35 mL			
Water added	Automatic			
After water addition, stir for	5 seconds			
At a speed of	15%			
Buffer in use	Yes			
Buffer type	Phosphate Buffer			
Volume of buffer introduced	0.025000 mL			
Add buffer manually	Manual			
After medium addition, stir for	5 seconds			
Sample Sonication				
Sonicate	No			
Sample Dissolution				
Perform a dissolution stage	No			
Carbonate purge				
Perform a carbonate purge	No			
Temperature Control				
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			

Sample name: **M15**
 Assay name: **UV-metric psKa**
 Assay ID: **17J-12008**
 Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
 Analyst: **Dorothy Levorse**
 Instrument ID: **T311053**

Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Time to wait	60 seconds			
Stir speed of	15%			
Titration 1				
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	10 seconds			
Titration 2				
Titrate from	Low to high pH			
Additional cosolvent volume	0.00 mL			
Add additional water	0.15 mL			
Additional water added	Automatic			
After pH adjust stir for	10 seconds			
Titration 3				
Titrate from	Low to high pH			
Additional cosolvent volume	0.00 mL			
Add additional water	0.34 mL			
Additional water added	Automatic			
After pH adjust stir for	10 seconds			
Data Point Stability				
Stir during data point collection	Yes			
For point collection, stir at	15%			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00500 dpH/dt			
Stability timeout after	60 seconds			
Experiment cleanup				
Adjust pH to cleanup	To start pH			
And then stir for	60 seconds			
For cleaning, stir at	20%			
Then add water volume	0.25 mL			
And then stir for	30 seconds			

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.109	10/12/2017 9:40:50 AM	C:\Sirius_T3\17J-11006_Blank standardisation.t3r
Four-Plus S	1.0007	10/12/2017 9:40:50 AM	C:\Sirius_T3\17J-11006_Blank standardisation.t3r
Four-Plus jH	0.3	10/12/2017 9:40:50 AM	C:\Sirius_T3\17J-11006_Blank standardisation.t3r
Four-Plus jOH	-0.2	10/12/2017 9:40:50 AM	C:\Sirius_T3\17J-11006_Blank standardisation.t3r
Base concentration factor	1.011	10/12/2017 9:40:50 AM	C:\Sirius_T3\KOH17122.t3r
Acid concentration factor	0.995	10/12/2017 9:40:50 AM	C:\Sirius_T3\17J-11005_Blank standardisation.t3r

Instrument Settings

Setting	Value	Batch Id	Install date
Instrument owner	Merck		
Instrument ID	T311053		
Instrument type	T3 Simulator		
Software version	1.1.3.0		
Dispenser module		T3DM1100253	3/31/2009 6:24:52 AM
Dispenser 0	Water		3/31/2009 6:25:05 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Water (0.15 M KCl)	10-10-2017	10/10/2017 10:48:53 AM
Dispenser 2	Acid		3/31/2009 6:25:11 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		

Sample name: **M15**
Assay name: **UV-metric psKa**
Assay ID: **17J-12008**
Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
Analyst: **Dorothy Levorse**
Instrument ID: **T311053**

Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Titration	Acid (0.5 M HCl)	166940 and 172875	10/6/2017 2:55:40 PM
Dispenser 1	Base		3/31/2009 6:25:21 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titration	Base (0.5 M KOH)	9-22-17	9/22/2017 4:02:42 PM
Dispenser 5	Cosolvent		3/31/2009 6:26:24 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Distribution valve 5	Distribution Valve		3/31/2009 6:28:19 AM
Firmware version	1.1.3		
Port A	Methanol (80%, 0.15 M KCl)	9-26-17	10/5/2017 5:02:03 PM
Port B	Cyclohexane		9/19/2017 2:15:02 PM
Port C	MeCN (50%, 0.15 M KCl)	10-2-17	10/2/2017 11:28:55 AM
Dispenser 3	Buffer		8/3/2010 6:05:16 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titration	Phosphate Buffer		10/10/2017 9:57:33 AM
Dispenser 6	Octanol		10/22/2010 11:52:43 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titration	Octanol	9-14-17	9/14/2017 10:30:38 AM
Titration		T3TM1100153	3/31/2009 6:24:17 AM
Horizontal axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Vertical axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Chassis I/O firmware version	1.11 AI1DI0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0769	8/15/2017 10:21:54 AM
E0 calibration	-9.19 mV		10/12/2017 9:41:14 AM
Filling solution	3M KCl	KCL095	10/10/2017 9:58:43 AM
Liquids			
Wash 1	50% IPA:50% Water		10/11/2017 8:31:15 AM
Wash 2	0.5% Triton X-100 in H2O		10/11/2017 8:31:17 AM
Buffer position 1	pH7 Wash		10/11/2017 8:31:21 AM
Buffer position 2	pH 7		10/11/2017 8:31:23 AM
Storage position			10/11/2017 8:31:26 AM
Wash water	4.2e+003 mL	10-6-17	10/6/2017 3:04:25 PM
Waste	5.9e+003 mL		10/6/2017 3:04:33 PM
Temperature controller			8/5/2010 7:35:13 AM
Turbidity detector			3/31/2009 6:24:45 AM
Spectrometer		072390	11/23/2010 12:22:28 PM
Dip probe		11086	
Wavelength coefficient A0	185.563		
Wavelength coefficient A1	2.17439		
Wavelength coefficient A2	-0.000285622		
Total lamp lit time	419:28:33		11/23/2010 12:22:28 PM
Calibrated on	10/11/2017 8:30:19 AM		
Integration time	10		
Scans averaged	10		
Autoloader		T3AL1100237	11/10/2015 10:34:13 AM
Left-right axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Front-back axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Vertical axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Chassis I/O firmware version	1.11 AI1DI0DO4 Norgren I/O		
Configuration			
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		



Sample name: **M15**
Assay name: **UV-metric psKa**
Assay ID: **17J-12008**
Filename: **C:\Sirius_T3\Mehtap\20171011_exp15_pKa\17J-12008_M15_UV-metric psKa.t3r**

Experiment start time: **10/12/2017 9:40:50 AM**
Analyst: **Dorothy Levorse**
Instrument ID: **T311053**

Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		

Refinement Settings

Setting	Value	Default value
Turbidity detection method	Spectrometer	Spectrometer
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00
Exclude turbid points	Yes	Yes
Low intensity warning threshold	100	100
Minimum absorbance change threshold	0.100	0.100
Eigenvector autocorrelation threshold	0.80	0.80
Maximum RMSD severe warning	0.250	0.250
Maximum RMSD warning	0.050	0.050

Tray Information

Title
Location E1