

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH +) -0.33 ±0.82 (n=50) logP (neutral X) 2.12 ±0.01 (n=50)

RMSD 0.537

18C-09010 Points 1 to 24

M11_octanol concentration factor 0.962
Carbonate 0.0600 mM
Acidity error -0.46392 mM

18C-09010 Points 25 to 45

M11_octanol concentration factor 0.930
Carbonate 0.0890 mM
Acidity error -0.75734 mM

18C-09010 Points 46 to 69

M11_octanol concentration factor 0.907
Carbonate 0.1290 mM
Acidity error -0.81079 mM

Warnings and errors

Errors None Warnings None

Sample logD and percent species

рН	M11_octanol logD	M11_octanol M11_octanolH	M11_octanol M11_octanol	M11_octanol M11_octanolH*	M11_octanol M11_octanol*	
1.000	-0.20	61.06 %	0.08 %	28.50 %	1 0 .36 %	
1.200	-0.13	57.55 %	0.12 %	26.86 %	15.48 %	Stomach pH
2.000	0.33	31.48 %	0.41 %	14.69 %	53.42 %	
3.000	1.19	5.39 %	0.69 %	2.51 %	91.40 %	
4.000	1.87	0.58 %	0.75 %	0.27 %	98.40 %	
5.000	2.09	0.06 %	0.75 %	0.03 %	99.16 %	
6.000	2.12	0.01 %	0.75 %	0.00 %	99.24 %	
6.500	2.12	0.00 %	0.75 %	0.00 %	99.24 %	
7.000	2.12	0.00 %	0.75 %	0.00 %	99.25 %	
7.400	2.12	0.00 %	0.75 %	0.00 %	99.25 %	Blood pH
8.000	2.12	0.00 %	0.75 %	0.00 %	99.25 %	·
9.000	2.12	0.00 %	0.75 %	0.00 %	99.25 %	
10.000	2.12	0.00 %	0.75 %	0.00 %	99.25 %	
11.000	2.12	0.00 %	0.75 %	0.00 %	99.25 %	
12.000	2.12	0.00 %	0.75 %	0.00 %	99.25 %	



Sample name: M11_octanol Assay name:

pH-metric high logP

18C-09010

Experiment start time: 3/9/2018 4:01:24 PM


Analyst: **Pion** Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r



Assay ID:

Filename:



Log R


pH (Concentration scale)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-09010 Points 1 to 24

Overall results

RMSD 0.632
Average ionic strength 0.157 M
Average temperature 24.9°C
Partition ratio 0.0123 : 1

Analyte concentration range 4427.7 µM to 4566.6 µM

Total points considered 21 of 24

Warnings and errors

Errors None Warnings None

Four-Plus parameters

Alpha 0.102 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r S 0.9967 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r jH 1.2 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r jOH 0.0 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r


Titrants


0.50 M HCI 0.999843 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r 0.50 M KOH 0.999845 3/9/2018 4:01:24 PM C:\Sirius_T3\KOH18B27.t3r

Sample

M11_octanol concentration factor 0.962
Base pKa 1 3.89
logP (XH +) 0.52
logP (neutral X) 2.10

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M11_octanol	M11_octanol	M11_octanol	M11_octanol	M11_octanol	Comment
	logD	M11_octanolH	M11_octanol	M11_octanolH*	M11_octanol*	
1.000	0.54	95.82 %	0.12 %	3.86 %	0.19 %	
1.200	0.55	95.65 %	0.20 %	3.86 %	0.30 %	Stomach pH
2.000	0.68	93.20 %	1.20 %	3.76 %	1.84 %	
3.000	1.23	73.19 %	9.43 %	2.95 %	14.43 %	
4.000	1.86	23.25 %	29.95 %	0.94 %	45.86 %	
5.000	2.07	2.97 %	38.29 %	0.12 %	58.62 %	
6.000	2.09	0.31 %	39.39 %	0.01 %	60.30 %	
6.500	2.10	0.10 %	39.47 %	0.00 %	60.43 %	
7.000	2.10	0.03 %	39.50 %	0.00 %	60.47 %	
7.400	2.10	0.01 %	39.51 %	0.00 %	60.48 %	Blood pH
8.000	2.10	0.00 %	39.51 %	0.00 %	60.49 %	
9.000	2.10	0.00 %	39.51 %	0.00 %	60.49 %	
10.000	2.10	0.00 %	39.51 %	0.00 %	60.49 %	
11.000	2.10	0.00 %	39.51 %	0.00 %	60.49 %	
12.000	2.10	0.00 %	39.51 %	0.00 %	60.49 %	

Carbonate and acidity



Other graphs





Assay ID:

Sample name: M11_octanol Assay name:

pH-metric high logP


18C-09010 Filename:

Experiment start time: 3/9/2018 4:01:24 PM

Analyst: Pion Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-09010 Points 25 to 45

Overall results

RMSD 0.373
Average ionic strength 0.164 M
Average temperature 25.0°C
Partition ratio 0.0405 : 1

Analyte concentration range 4011.7 µM to 4147.6 µM

Total points considered 18 of 21

Warnings and errors

Errors None Warnings None

Four-Plus parameters

 Alpha
 0.102
 3/9/2018 4:01:24 PM
 C:\Sirius_T3\HCl18C09.t3r

 S
 0.9967
 3/9/2018 4:01:24 PM
 C:\Sirius_T3\HCl18C09.t3r

 jH
 1.2
 3/9/2018 4:01:24 PM
 C:\Sirius_T3\HCl18C09.t3r

 jOH
 0.0
 3/9/2018 4:01:24 PM
 C:\Sirius_T3\HCl18C09.t3r


Titrants


0.50 M HCI 0.999843 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r 0.50 M KOH 0.999845 3/9/2018 4:01:24 PM C:\Sirius_T3\KOH18B27.t3r

Sample

M11_octanol concentration factor 0.930
Base pKa 1 3.89
logP (XH +) 0.55
logP (neutral X) 2.22

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

M11_octanol	M11_octanol				Comment
logD	M11_octanolH	M11_octanol	M11_octanolH*	M11_octanol*	
0.57	86.67 %	0.11 %	12.47 %	0.75 %	
0.59	86.23 %	0.18 %	12.40 %	1.19 %	Stomach pH
0.75	80.41 %	1.04 %	11.57 %	6.98 %	
1.34	46.71 %	6.02 %	6.72 %	40.55 %	
1.98	9.00 %	11.59 %	1.29 %	78.12 %	
2.19	0.99 %	12.77 %	0.14 %	86.09 %	
2.22	0.10 %	12.91 %	0.01 %	86.98 %	
2.22	0.03 %	12.92 %	0.00 %	87.05 %	
2.22	0.01 %	12.92 %	0.00 %	87.07 %	
2.22	0.00 %	12.92 %	0.00 %	87.08 %	Blood pH
2.22	0.00 %	12.92 %	0.00 %	87.08 %	
2.22	0.00 %	12.92 %	0.00 %	87.08 %	
2.22	0.00 %	12.92 %	0.00 %	87.08 %	
2.22	0.00 %	12.92 %	0.00 %	87.08 %	
2.22	0.00 %	12.92 %	0.00 %	87.08 %	
	1.98 2.22 2.22 2.22 2.22 2.22 2.22 2.22 2	logD M11_octanolH 0.57 86.67 % 0.59 86.23 % 0.75 80.41 % 1.34 46.71 % 1.98 9.00 % 2.19 0.99 % 2.22 0.10 % 2.22 0.01 % 2.22 0.00 % 2.22 0.00 % 2.22 0.00 % 2.22 0.00 % 2.22 0.00 % 2.22 0.00 % 2.22 0.00 % 2.22 0.00 % 2.22 0.00 %	logD M11_octanolH M11_octanol 0.57 86.67 % 0.11 % 0.59 86.23 % 0.18 % 0.75 80.41 % 1.04 % 1.34 46.71 % 6.02 % 1.98 9.00 % 11.59 % 2.19 0.99 % 12.77 % 2.22 0.10 % 12.91 % 2.22 0.03 % 12.92 % 2.22 0.01 % 12.92 % 2.22 0.00 % 12.92 % 2.22 0.00 % 12.92 % 2.22 0.00 % 12.92 % 2.22 0.00 % 12.92 % 2.22 0.00 % 12.92 % 2.22 0.00 % 12.92 % 2.22 0.00 % 12.92 % 2.22 0.00 % 12.92 %	logD M11_octanolH M11_octanol M11_octanolH* 0.57 86.67 % 0.11 % 12.47 % 0.59 86.23 % 0.18 % 12.40 % 0.75 80.41 % 1.04 % 11.57 % 1.34 46.71 % 6.02 % 6.72 % 1.98 9.00 % 11.59 % 1.29 % 2.19 0.99 % 12.77 % 0.14 % 2.22 0.10 % 12.91 % 0.01 % 2.22 0.03 % 12.92 % 0.00 % 2.22 0.00 % 12.92 % 0.00 % 2.22 0.00 % 12.92 % 0.00 % 2.22 0.00 % 12.92 % 0.00 % 2.22 0.00 % 12.92 % 0.00 % 2.22 0.00 % 12.92 % 0.00 % 2.22 0.00 % 12.92 % 0.00 % 2.22 0.00 % 12.92 % 0.00 % 2.22 0.00 % 12.92 % 0.00 %	logD M11_octanolH M11_octanol M11_octanolH* M11_octanol* 0.57 86.67 % 0.11 % 12.47 % 0.75 % 0.59 86.23 % 0.18 % 12.40 % 1.19 % 0.75 80.41 % 1.04 % 11.57 % 6.98 % 1.34 46.71 % 6.02 % 6.72 % 40.55 % 1.98 9.00 % 11.59 % 1.29 % 78.12 % 2.19 0.99 % 12.77 % 0.14 % 86.09 % 2.22 0.10 % 12.91 % 0.01 % 86.98 % 2.22 0.03 % 12.92 % 0.00 % 87.05 % 2.22 0.01 % 12.92 % 0.00 % 87.08 % 2.22 0.00 % 12.92 % 0.00 % 87.08 % 2.22 0.00 % 12.92 % 0.00 % 87.08 % 2.22 0.00 % 12.92 % 0.00 % 87.08 % 2.22 0.00 % 12.92 % 0.00 % 87.08 % 2.22 0.00 % 12.92 % </td

Carbonate and acidity



Other graphs






Assay name: pH-metric high logP Analyst: **Pion** 18C-09010 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-09010 Points 46 to 69

Overall results

RMSD 0.548
Average ionic strength 0.170 M
Average temperature 25.0°C
Partition ratio 0.1193 : 1

Analyte concentration range 3477.2 µM to 3586.0 µM

Total points considered 22 of 24

Warnings and errors

Errors None Warnings None

Four-Plus parameters

Alpha 0.102 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r S 0.9967 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r jH 1.2 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r jOH 0.0 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r


Titrants


0.50 M HCI 0.999843 3/9/2018 4:01:24 PM C:\Sirius_T3\HCl18C09.t3r 0.50 M KOH 0.999845 3/9/2018 4:01:24 PM C:\Sirius_T3\KOH18B27.t3r

Sample

M11_octanol concentration factor 0.907
Base pKa 1 3.89
logP (XH +) 0.55
logP (neutral X) 2.25

Sample graphs








Assay ID:

Sample name: M11_octanol Experiment start time: 3/9/2018 4:01:24 PM Assay name:

pH-metric high logP Analyst: Pion Instrument ID: T312060 18C-09010

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Sample graphs (continued)




Sample logD and percent species

рН	M11_octanol	M11_octanol	M11_octanol	M11_octanol	M11_octanol	Comment
	logD	M11_octanolH	M11_octanol	M11_octanolH*	M11_octanol*	
1.000	0.58	68.74 %	0.09 %	29.31 %	1.86 %	
1.200	0.59	67.97 %	0.14 %	28.98 %	2.92 %	Stomach pH
2.000	0.76	58.48 %	0.75 %	24.93 %	15.84 %	
3.000	1.37	23.45 %	3.02 %	10.00 %	63.53 %	
4.000	2.00	3.36 %	4.32 %	1.43 %	90.89 %	
5.000	2.21	0.35 %	4.52 %	0.15 %	94.98 %	
6.000	2.24	0.04 %	4.54 %	0.02 %	95.41 %	
6.500	2.25	0.01 %	4.54 %	0.00 %	95.44 %	
7.000	2.25	0.00 %	4.54 %	0.00 %	95.46 %	
7.400	2.25	0.00 %	4.54 %	0.00 %	95.46 %	Blood pH
8.000	2.25	0.00 %	4.54 %	0.00 %	95.46 %	
9.000	2.25	0.00 %	4.54 %	0.00 %	95.46 %	
10.000	2.25	0.00 %	4.54 %	0.00 %	95.46 %	
11.000	2.25	0.00 %	4.54 %	0.00 %	95.46 %	
12.000	2.25	0.00 %	4.54 %	0.00 %	95.46 %	

Carbonate and acidity



Other graphs



0.16



Assay ID: Filename:

Sample name: M11_octanol Assay name:

pH-metric high logP


18C-09010

Experiment start time: 3/9/2018 4:01:24 PM

Pion Analyst: Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name	M11_octanol	2/27/2018 5:54:30 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.001520 g	3/9/2018 2:22:09 PM	User entered value
Formula weight	211.22 g/mol	2/27/2018 5:54:30 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	211.22	2/27/2018 5:54:30 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	1	2/27/2018 5:54:30 PM	User entered value
Sample is a	Base	2/27/2018 5:54:30 PM	User entered value
pKa 1	3.89	2/27/2018 5:54:30 PM	User entered value
logp (XH +)	0.55	3/2/2018 4:29:35 PM	User entered value
logP (neutral X)	2.19	3/2/2018 4:29:30 PM	User entered value

Events

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
5:00.2	Initial pH = 9.10									
7:59.8	Data point 1		0.05315 mL					0.30733	0.00046	
8:46.1	Data point 2		0.05315 mL					0.69530	0.00017	10.5 s
9:22.2	Data point 3		0.05315 mL					0.32148	0.00033	10.0 s
9:57.8	Data point 4	1.50000 mL	0.05315 mL	0.03048 mL	0.01999 mL	2.655	0.00078	0.10371	0.00012	10.5 s
10:33.9	Data point 5	1.50000 mL	0.05315 mL	0.03452 mL	0.01999 mL	2.854	-0.00337	0.19592	0.00038	10.0 s
11:09.3	Data point 6	1.50000 mL	0.05315 mL	0.03765 mL	0.01999 mL	3.078	0.00146	0.01563	0.00058	10.5 s
11:45.3	Data point 7	1.50000 mL	0.05315 mL	0.04026 mL	0.01999 mL	3.271	-0.00477	0.53594	0.00032	10.0 s
12:20.8			0.05315 mL						0.00038	10.0 s
12:56.3	Data point 9	1.50000 mL	0.05315 mL	0.04478 mL	0.01999 mL	3.707	-0.01214	0.77245	0.00068	10.5 s
13:32.3	Data point 10	1.50000 mL	0.05315 mL	0.04657 mL	0.01999 mL	3.912	-0.01020	0.53326	0.00069	10.0 s
14:07.7	Data point 11	1.50000 mL	0.05315 mL	0.04798 mL	0.01999 mL	4.094	-0.00754	0.73914	0.00043	10.0 s
14:43.1	Data point 12	1.50000 mL	0.05315 mL	0.04908 mL	0.01999 mL	4.274	-0.01378	0.67663	0.00083	10.0 s
15:18.5	Data point 13	1.50000 mL	0.05315 mL	0.04988 mL	0.01999 mL	4.454	-0.01794	0.90958	0.00093	16.5 s
16:00.3	Data point 14	1.50000 mL	0.05315 mL	0.05047 mL	0.01999 mL	4.635	-0.01742	0.86808	0.00092	15.0 s
16:40.7	Data point 15	1.50000 mL	0.05315 mL	0.05087 mL	0.01999 mL	4.824	-0.01787	0.84622	0.00096	24.5 s
17:30.7	Data point 16	1.50000 mL	0.05315 mL	0.05115 mL	0.01999 mL	5.023	-0.01827	0.90373	0.00095	21.0 s
18:17.2	Data point 17	1.50000 mL	0.05315 mL	0.05134 mL	0.01999 mL	5.193	-0.01878	0.91929	0.00097	23.0 s
19:10.8	Data point 18	1.50000 mL	0.05315 mL	0.05155 mL	0.01999 mL	5.715	-0.01901	0.96541	0.00096	37.0 s
20:18.6	Data point 19	1.50000 mL	0.05315 mL	0.05174 mL	0.01999 mL	6.941	-0.08865	0.99621	0.00438	Timed out at 59.5 s
21:54.2	Data point 20	1.50000 mL	0.05315 mL	0.05181 mL	0.01999 mL	7.267	-0.07065	0.99539	0.00350	Timed out at 59.5 s
23:29.8	Data point 21	1.50000 mL	0.05315 mL	0.05188 mL	0.01999 mL	7.738	-0.08731	0.99760	0.00432	Timed out at 59.5 s
25:00.3	Data point 22	1.50000 mL	0.05315 mL	0.05193 mL	0.01999 mL	8.184	-0.04618	0.99543	0.00229	Timed out at 59.5 s
26:36.0	Data point 23	1.50000 mL	0.05315 mL	0.05200 mL	0.01999 mL	8.628	-0.02509	0.99264	0.00124	Timed out at 59.5 s
28:11.7	Data point 24	1.50000 mL	0.05315 mL	0.05212 mL	0.01999 mL	9.017	-0.01906	0.93147	0.00097	40.5 s
29:51.6	Data point 25	1.50000 mL	0.11291 mL	0.05212 mL	0.06999 mL	1.966	0.00082	0.01226	0.00037	10.5 s
30:38.4	Data point 26	1.50000 mL	0.11291 mL	0.06950 mL	0.06999 mL	2.173	-0.00810	0.69034	0.00048	10.5 s
31:14.6	Data point 27	1.50000 mL	0.11291 mL	0.08090 mL	0.06999 mL	2.379	0.00330	0.80719	0.00018	10.0 s
31:50.1		1.50000 mL	0.11291 mL	0.08841 mL	0.06999 mL	2.580	-0.01357	0.47607	0.00097	14.0 s
32:29.7	Data point 29	1.50000 mL	0.11291 mL	0.09374 mL	0.06999 mL	2.779	-0.00359	0.57678	0.00023	10.0 s
33:05.2		1.50000 mL	0.11291 mL	0.09781 mL	0.06999 mL	2.960	-0.00538	0.43414	0.00040	
33:40.7			0.11291 mL						0.00086	
34:27.1			0.11291 mL						0.00074	10.0 s
35:02.5			0.11291 mL						0.00098	10.0 s
35:38.0	Data point 34		0.11291 mL						0.00029	10.0 s
00.40.4	D-1 1-1-05	4 50000	0.44004	0.40040 1	0.00000 1	4 00 4	0.00755	0.00007	0 00000	400 -

1.50000 mL 0.11291 mL 0.10842 mL 0.06999 mL 4.034 -0.00755 0.89227

36:13.4 Data point 35

0.00039 10.0 s



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Events (continued)

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
36:48.8	Data point 36	1 50000 ml	0 11291 ml	0.10915 mL	0 06999 ml	4 284	-0.01862	0.91237	0.00096	
37:29.3	Data point 37			0.10960 mL					0.00092	
38:19.8	Data point 38			0.10988 mL					0.00095	
39:10.7	Data point 39			0.11018 mL					0.00097	
40:16.4	Data point 40			0.11035 mL					0.00097	
41:22.4	Data point 41			0.11042 mL						Timed out
	Bata point 11	1.00000 1112	0.112011112	0.110121112	0.00000 1112	0.010	0.00700	0.00100	0.00100	at 59.5 s
42:53.0	Data point 42	1.50000 mL	0.11291 mL	0.11056 mL	0.06999 mL	7.814	-0.06409	0.99645	0.00317	Timed out
										at 59.5 s
44:28.6	Data point 43	1.50000 mL	0.11291 mL	0.11070 mL	0.06999 mL	8.652	-0.02006	0.98249	0.00100	46.0 s
45:50.3	Data point 44			0.11082 mL					0.00098	38.5 s
46:59.3	Data point 45	1.50000 mL	0.11291 mL	0.11091 mL	0.06999 mL	9.073	-0.01775	0.86902	0.00094	27.5 s
48:28.6	Data point 46	1.50000 mL	0.17582 mL	0.11091 mL	0.21999 mL	1.960	-0.00614	0.61144	0.00039	10.0 s
49:14.9	Data point 47	1.50000 mL	0.17582 mL	0.13074 mL	0.21999 mL	2.169	0.01032	0.30311	0.00093	10.0 s
49:50.6	Data point 48	1.50000 mL	0.17582 mL	0.14341 mL	0.21999 mL	2.367	-0.00532	0.40508	0.00041	10.0 s
50:26.1	Data point 49	1.50000 mL	0.17582 mL	0.15214 mL	0.21999 mL	2.563	-0.00517	0.08291	0.00089	10.5 s
51:02.2	Data point 50	1.50000 mL	0.17582 mL	0.15837 mL	0.21999 mL	2.768	-0.00887	0.21333	0.00095	10.0 s
51:37.8	Data point 51	1.50000 mL	0.17582 mL	0.16289 mL	0.21999 mL	2.967	0.00295	0.02140	0.00100	10.0 s
52:13.3	Data point 52			0.16625 mL				0.17172	0.00041	
52:48.8	Data point 53			0.16860 mL					0.00100	
53:29.5	Data point 54			0.16992 mL				0.85818	0.00077	
54:43.8	Data point 55			0.17110 mL				0.01708	0.00035	10.5 s
55:19.9	Data point 56			0.17180 mL				0.25159	0.00087	
56:00.9	Data point 57			0.17222 mL				0.95151	0.00097	11.5 s
56:37.7	Data point 58			0.17248 mL					0.00092	
57:29.8	Data point 59			0.17277 mL					0.00071	
58:19.3	Data point 60			0.17288 mL					0.00095	
59:10.9	Data point 61			0.17295 mL					0.00080	
1:00:11.1	Data point 62			0.17307 mL						Timed out
	5	4 =0000	==	0.4=0.40	0.04000			0.00440		at 59.5 s
1:01:41.6	Data point 63	1.50000 mL	0.1/582 mL	0.17319 mL	0.21999 mL	7.526	-0.12365	0.99119	0.00614	Timed out at 59.5 s
1:03:12.0	Data point 64	1.50000 mL	0.17582 mL	0.17324 mL	0.21999 mL	7.651	-0.08733	0.99550	0.00433	Timed out
1.04.42 5	Data point 65	1 50000 ml	0 17582 ml	0.17328 mL	0 21000 ml	7 837	_∩ 13623	0.43667	0.01018	at 59.5 s Timed out
1.04.42.0	Data point 05	1.50000 IIIL	0.17302 IIIL	0.17 320 IIIL	0.2 1999 IIIL	1.001	-0.13023	0.43007	0.01010	at 59.5 s
1:06:18.2	Data point 66	1.50000 mL	0.17582 mL	0.17335 mL	0.21999 mL	8.021	-0.06575	0.99623	0.00325	Timed out at 59.5 s
1:07:53.8	Data point 67	1.50000 mL	0.17582 mL	0.17345 mL	0.21999 mL	8.434	-0.11867	0.41261	0.00912	Timed out
1:09:29.4	Data point 68	1.50000 ml	0.17582 ml	0.17361 mL	0.21999 ml	8.870	-0.01721	0.79761	0.00095	at 59.5 s 21.0 s
	Data point 69			0.17373 mL					0.00096	
	A cook values co			0.17070 mL					, .	

1:10:50.6 Assay volumes 1.50000 mL 0.17582 mL 0.17373 mL 0.21999 mL



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	9.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings	' '			
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium	110110			
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.020 mL			
Partition solvent added	Automatic			
After partition addition, stir for	1 seconds			
Sample Sonication	i seconds			
Sonicate	Yes			
	No			
Adjust pH for sonication Sonicate for	60 seconds			
After sonication stir for				
	5 seconds			
Sample Dissolution	Vaa			
Perform a dissolution stage	Yes			
	To start pH			
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge				
Perform a carbonate purge	No			
Temperature Control				
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			
Stir speed of	50%			
Titration 1				
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	50%			
Titration 2				
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume				
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			

30 seconds

15 seconds

55%

After pH adjust stir for

Stir to allow partitioning for

Stirrer speed for partitioning



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3		•	•	•
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.150 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00100 dpH/dt			
Stability timeout after	60 seconds			

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.102	3/9/2018 4:01:24 PM	C:\Sirius_T3\HCl18C09.t3r
Four-Plus S	0.9967	3/9/2018 4:01:24 PM	C:\Sirius_T3\HCl18C09.t3r
Four-Plus jH	1.2	3/9/2018 4:01:24 PM	C:\Sirius_T3\HCl18C09.t3r
Four-Plus jOH	0.0	3/9/2018 4:01:24 PM	C:\Sirius_T3\HCl18C09.t3r
Base concentration factor	1.000	3/9/2018 4:01:24 PM	C:\Sirius_T3\KOH18B27.t3r
Acid concentration factor	1.000	3/9/2018 4:01:24 PM	C:\Sirius_T3\HCl18C09.t3r

Instrument Settings

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module Dispenser 0 Syringe volume Firmware version	Water 2.5 mL 1.2.1(r2)	T3DM1200361	3/31/2009 6:24:52 AM 3/31/2009 6:25:05 AM
Titrant Dispenser 2 Syringe volume Firmware version	Water (0.15 M KCI) Acid 0.5 mL 1.2.1(r2)	02-06-2018	2/27/2018 11:05:59 AM 3/31/2009 6:25:11 AM
Titrant Dispenser 1 Syringe volume Firmware version	Acid (0.5 M HCI) Base 0.5 mL 1.2.1(r2)	02-27-2018	2/27/2018 11:27:22 AM 3/31/2009 6:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2)	9/22/2017	2/27/2018 11:21:22 AM 3/31/2009 6:26:24 AM
Distribution valve 5 Firmware version Port A Port B Dispenser 3	Distribution Valve 1.1.3 Methanol (80%, 0.15 M KCI) Cyclohexane Buffer	02-08-2018 11-01-17	3/31/2009 6:28:19 AM 3/6/2018 10:28:59 AM 2/27/2018 11:37:57 AM 8/3/2010 6:05:16 AM
Syringe volume Firmware version Titrant Dispenser 6	0.5 mL 1.2.1(r2) Dodecane Octanol	2018/01/31	2/28/2018 11:18:04 AM 10/22/2010 11:52:43 AM




Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting Suringa valuma	Value	Batch Id	Install date
Syringe volume Firmware version	0.5 mL 1.2.1(r2)		
Titrant	Octanol	01-31-2018	2/27/2018 10:59:35 AM
Titrator	Cotanor		3/31/2009 6:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2	1011111200101	0,0 1,2000 0.2 1.11 7
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0923	1/23/2018 3:01:00 PM
E0 calibration	+6.40 mV	1401.007	3/9/2018 4:01:52 PM
Filling solution	3M KCI	KCL097	3/9/2018 11:05:42 AM
Liquids Wash 1	50% IPA:50% Water		3/9/2018 11:04:22 AM
Wash 2	0.5% Trition X-100 in H20		3/9/2018 11:04:25 AM
Buffer position 1	pH7 Wash		3/9/2018 11:04:27 AM
Buffer position 2	pH 7		3/9/2018 11:04:30 AM
Storage position	•		3/9/2018 11:05:04 AM
Wash water	5.3e+003 mL	02-27-2018	2/27/2018 10:54:39 AM
Waste	1e+004 mL		11/28/2017 11:36:29 AM
Temperature controller			8/5/2010 7:35:13 AM
Turbidity detector		074044	3/31/2009 6:24:45 AM
Spectrometer		074811	11/23/2010 12:22:28 PM
Dip probe	183.333	10196	
Wavelength coefficient A0 Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	123:16:41		11/23/2010 12:22:28 PM
Calibrated on	2/27/2018 11:40:38 AM		
Integration time	40		
Scans averaged	10		
Autoloader		T3AL1200345	11/10/2015 10:34:13 AM
Left-right axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version Chassis I/O firmware version	1.17 AI1DI2DO2 Stepper 2 1.11 AI1DI0DO4 Norgren I/O		
Configuration	1.11 Al IDIODO4 Noigieil I/O		
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed Solvent wash stir duration	30% 5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed E0 calibration reading stir speed	30% 0%		
Lo campiation reading still speed	O /0		





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-09010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180309_exp31_logP_T3-2\18C-09010_M11_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting
Spectrometer calibration stir duration
Spectrometer calibration stir speed
Spectrometer calibration wash pump volume
Spectrometer calibration wash stir duration
Spectrometer calibration wash stir speed
Spectrometer calibration wash stir speed
Overhead dispense height

Satch Id Install date
Install date

5 s

30%

30%

10000

Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00
• •		