

Assay name: pH-metric high logP Analyst: Pion 18C-01002 Instrument ID: T312060 Assay ID:

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

## pH-metric Result

logP (XH +) -3.66 ±2.11 (n=50) logP (neutral X) 2.01 ±0.02 (n=50)

### 18C-01002 Points 1 to 28

M11 octanol concentration factor 0.940 Carbonate 0.0379 mM Acidity error 0.00456 mM

### 18C-01002 Points 29 to 49

M11\_octanol concentration factor 0.777 Carbonate 0.1154 mM -0.46703 mM Acidity error

### 18C-01002 Points 50 to 75

M11 octanol concentration factor 0.713 Carbonate 0.1025 mM Acidity error -1.33190 mM

## Warnings and errors

Errors None

Warnings One or more logP values out of range

## Sample logD and percent species

рН	M11_octanol logD	M11_octanol M11_octanolH	M11_octanol M11_octanol	M11_octanol M11_octanolH*	M11_octanol M11_octanol*	
1.000	-0.88	88.32 %	0.11 %	0.02 %	11.55 %	
1.200	-0.68	82.68 %	0.17 %	0.02 %	17.13 %	Stomach pH
2.000	0.11	43.10 %	0.56 %	0.01 %	56.34 %	
3.000	1.06	7.04 %	0.91 %	0.00 %	92.05 %	
4.000	1.76	0.75 %	0.97 %	0.00 %	98.28 %	
5.000	1.97	0.08 %	0.98 %	0.00 %	98.95 %	
6.000	2.00	0.01 %	0.98 %	0.00 %	99.02 %	
6.500	2.01	0.00 %	0.98 %	0.00 %	99.02 %	
7.000	2.01	0.00 %	0.98 %	0.00 %	99.02 %	
7.400	2.01	0.00 %	0.98 %	0.00 %	99.02 %	Blood pH
8.000	2.01	0.00 %	0.98 %	0.00 %	99.02 %	
9.000	2.01	0.00 %	0.98 %	0.00 %	99.02 %	
10.000	2.01	0.00 %	0.98 %	0.00 %	99.02 %	
11.000	2.01	0.00 %	0.98 %	0.00 %	99.02 %	
12.000	2.01	0.00 %	0.98 %	0.00 %	99.02 %	



Assay ID:

Sample name: M11\_octanol Assay name:


pH-metric high logP

18C-01002

Experiment start time: 3/1/2018 1:21:08 AM

Analyst: Pion Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r



2

3

5

7

pH (Concentration scale)

-2

Log R

0

0

-4


11



pH-metric high logP Analyst: Assay name: Pion Assay ID: 18C-01002 Instrument ID: T312060 Filename:

C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

# **Graphs** (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 1 of 3 18C-01002 Points 1 to 28

#### Overall results

RMSD 0.953
Average ionic strength 0.157 M
Average temperature 24.9°C
Partition ratio 0.0122 : 1

Analyte concentration range 3149.9 µM to 3249.9 µM

Total points considered 22 of 28

## Warnings and errors

Errors None

Warnings One or more logP values out of range

## Four-Plus parameters

Alpha 0.130 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r S 0.9970 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r jH 0.8 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r jOH -0.4 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r


# Titrants


0.50 M HCI 0.993513 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r 0.50 M KOH 0.999845 3/1/2018 1:21:08 AM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M11\_octanol concentration factor 0.940
Base pKa 1 3.89
logP (XH +) -5.60
logP (neutral X) 1.94

### Sample graphs







Assay ID:


Filename:

Sample name: M11\_octanol Experiment start time: 3/1/2018 1:21:08 AM Assay name:

pH-metric high logP Analyst: Pion Instrument ID: 18C-01002 T312060

C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r


## Sample graphs (continued)




## Sample logD and percent species


рН	M11_octanol	M11_octanol	M11_octanol		M11_octanol	Comment
	logD	M11_octanolH	M11_octanol	M11_octanolH*	M11_octanol*	
1.000	-0.96	99.74 %	0.13 %	0.00 %	0.14 %	
1.200	-0.76	99.58 %	0.20 %	0.00 %	0.21 %	Stomach pH
2.000	0.04	97.42 %	1.26 %	0.00 %	1.32 %	
3.000	0.99	79.07 %	10.19 %	0.00 %	10.74 %	
4.000	1.69	27.42 %	35.32 %	0.00 %	37.26 %	
5.000	1.90	3.64 %	46.90 %	0.00 %	49.46 %	
6.000	1.93	0.38 %	48.48 %	0.00 %	51.14 %	
6.500	1.93	0.12 %	48.61 %	0.00 %	51.27 %	
7.000	1.93	0.04 %	48.65 %	0.00 %	51.31 %	
7.400	1.93	0.02 %	48.66 %	0.00 %	51.33 %	Blood pH
8.000	1.94	0.00 %	48.67 %	0.00 %	51.33 %	
9.000	1.94	0.00 %	48.67 %	0.00 %	51.33 %	
10.000	1.94	0.00 %	48.67 %	0.00 %	51.33 %	
11.000	1.94	0.00 %	48.67 %	0.00 %	51.33 %	
12.000	1.94	0.00 %	48.67 %	0.00 %	51.33 %	

# Carbonate and acidity



## Other graphs








Analyst: Assay name: pH-metric high logP Pion 18C-01002 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 2 of 3 18C-01002 Points 29 to 49

### Overall results

RMSD 0.536
Average ionic strength 0.163 M
Average temperature 25.0°C
Partition ratio 0.0407 : 1

Analyte concentration range 2865.3 µM to 2956.0 µM

Total points considered 17 of 21

## Warnings and errors

Errors None

Warnings One or more logP values out of range

## Four-Plus parameters

Alpha 0.130 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r S 0.9970 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r jH 0.8 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r jOH -0.4 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r


# Titrants


0.50 M HCI 0.993513 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r 0.50 M KOH 0.999845 3/1/2018 1:21:08 AM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M11\_octanol concentration factor 0.777
Base pKa 1 3.89
logP (XH +) -5.10
logP (neutral X) 2.08

### Sample graphs








Assay ID:

Sample name: M11\_octanol Experiment start time: 3/1/2018 1:21:08 AM

Assay name: pH-metric high logP Analyst: Pion Instrument ID: 18C-01002 T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r


## Sample graphs (continued)




## Sample logD and percent species


рН	M11_octanol		M11_octanol		M11_octanol	Comment
	logD	W11_octanoiH	M11_octanoi	M11_octanolH*	M11_octanol*	
1.000	-0.81	99.24 %	0.13 %	0.00 %	0.63 %	
1.200	-0.61	98.81 %	0.20 %	0.00 %	0.99 %	Stomach pH
2.000	0.19	92.92 %	1.20 %	0.00 %	5.88 %	
3.000	1.14	56.77 %	7.31 %	0.00 %	35.91 %	
4.000	1.83	11.61 %	14.96 %	0.00 %	73.44 %	
5.000	2.05	1.30 %	16.70 %	0.00 %	82.00 %	
6.000	2.08	0.13 %	16.90 %	0.00 %	82.97 %	
6.500	2.08	0.04 %	16.91 %	0.00 %	83.05 %	
7.000	2.08	0.01 %	16.92 %	0.00 %	83.07 %	
7.400	2.08	0.01 %	16.92 %	0.00 %	83.08 %	Blood pH
8.000	2.08	0.00 %	16.92 %	0.00 %	83.08 %	·
9.000	2.08	0.00 %	16.92 %	0.00 %	83.08 %	
10.000	2.08	0.00 %	16.92 %	0.00 %	83.08 %	
11.000	2.08	0.00 %	16.92 %	0.00 %	83.08 %	
12.000	2.08	0.00 %	16.92 %	0.00 %	83.08 %	

## Carbonate and acidity



# Other graphs







Assay name:

Assay ID: Filename:

Sample name: M11\_octanol

pH-metric high logP

18C-01002


C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

Experiment start time: 3/1/2018 1:21:08 AM

Analyst: Pion Instrument ID: T312060

pH (Concentration scale)

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 3 of 3 18C-01002 Points 50 to 75

### Overall results

RMSD 0.550
Average ionic strength 0.169 M
Average temperature 25.0°C
Partition ratio 0.1199 : 1

Analyte concentration range 2489.4 µM to 2561.6 µM

Total points considered 21 of 26

## Warnings and errors

Errors None

Warnings One or more logP values out of range

Excessive acidity error present

## Four-Plus parameters

Alpha 0.130 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r S 0.9970 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r jH 0.8 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r jOH -0.4 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r


#### Titrants


0.50 M HCI 0.993513 3/1/2018 1:21:08 AM C:\Sirius\_T3\HCl18B27.t3r 0.50 M KOH 0.999845 3/1/2018 1:21:08 AM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M11\_octanol concentration factor 0.713
Base pKa 1 3.89
logP (XH +) -5.10
logP (neutral X) 1.99

#### Sample graphs








Assay ID:

Sample name: M11\_octanol Experiment start time: 3/1/2018 1:21:08 AM

Assay name: pH-metric high logP Analyst: Pion Instrument ID: 18C-01002 T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r


## Sample graphs (continued)




## Sample logD and percent species


рН	M11_octanol	M11_octanol	M11_octanol		M11_octanol	
	logD	M11_octanolH	M11_octanol	M11_octanolH*	M11_octanol*	
1.000	-0.90	98.39 %	0.13 %	0.00 %	1.49 %	
1.200	-0.70	97.47 %	0.20 %	0.00 %	2.33 %	Stomach pH
2.000	0.09	85.92 %	1.11 %	0.00 %	12.97 %	
3.000	1.05	37.90 %	4.88 %	0.00 %	57.21 %	
4.000	1.74	5.75 %	7.41 %	0.00 %	86.84 %	
5.000	1.96	0.61 %	7.82 %	0.00 %	91.58 %	
6.000	1.99	0.06 %	7.86 %	0.00 %	92.08 %	
6.500	1.99	0.02 %	7.86 %	0.00 %	92.12 %	
7.000	1.99	0.01 %	7.86 %	0.00 %	92.13 %	
7.400	1.99	0.00 %	7.86 %	0.00 %	92.13 %	Blood pH
8.000	1.99	0.00 %	7.86 %	0.00 %	92.14 %	
9.000	1.99	0.00 %	7.86 %	0.00 %	92.14 %	
10.000	1.99	0.00 %	7.86 %	0.00 %	92.14 %	
11.000	1.99	0.00 %	7.86 %	0.00 %	92.14 %	
12.000	1.99	0.00 %	7.86 %	0.00 %	92.14 %	

## Carbonate and acidity



# Other graphs







Sample name: M11\_octanol Assay name:

pH-metric high logP


18C-01002 Assay ID: Filename:

Experiment start time: 3/1/2018 1:21:08 AM

Analyst: Pion Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

## Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name	M11_octanol	2/27/2018 4:54:30 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.001080 g	2/28/2018 4:24:59 PM	User entered value
Formula weight	211.22 g/mol	2/27/2018 4:54:30 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	211.22	2/27/2018 4:54:30 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	1	2/27/2018 4:54:30 PM	User entered value
Sample is a	Base	2/27/2018 4:54:30 PM	User entered value
pKa 1	3.89	2/27/2018 4:54:30 PM	User entered value
logp (XH +)	-5.60	2/28/2018 1:53:21 PM	User entered value
logP (neutral X)	2.09	2/28/2018 1:53:11 PM	User entered value

#### **Events**

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
8:59.7	Initial pH = 8.06									-
11:59.3	Data point 1	1.50000 mL	0.05153 mL	0.00179 mL	0.01999 mL	2.005	-0.00768	0.87550	0.00040	10.0 s
12:45.4	Data point 2	1.50000 mL	0.05153 mL	0.01700 mL	0.01999 mL	2.215	-0.00104	0.14554	0.00013	10.0 s
13:21.1	Data point 3	1.50000 mL	0.05153 mL	0.02672 mL	0.01999 mL	2.436	-0.00042	0.01137	0.00020	10.0 s
13:56.6	Data point 4	1.50000 mL	0.05153 mL	0.03274 mL	0.01999 mL	2.635	-0.00726	0.18423	0.00084	10.0 s
14:32.2	Data point 5	1.50000 mL	0.05153 mL	0.03681 mL	0.01999 mL	2.861	-0.00222	0.27923	0.00021	10.0 s
15:07.7	Data point 6	1.50000 mL	0.05153 mL	0.03965 mL	0.01999 mL	3.060	-0.00837	0.77297	0.00047	10.0 s
15:43.1	Data point 7	1.50000 mL	0.05153 mL	0.04191 mL	0.01999 mL	3.262	-0.00454	0.61543	0.00029	10.0 s
16:18.6	Data point 8				0.01999 mL				0.00029	10.0 s
16:54.1	Data point 9	1.50000 mL	0.05153 mL	0.04553 mL	0.01999 mL	3.705	-0.01127	0.51278	0.00078	10.0 s
17:29.6	Data point 10	1.50000 mL	0.05153 mL	0.04697 mL	0.01999 mL	3.894	-0.00094	0.05467	0.00020	10.5 s
18:05.4	Data point 11	1.50000 mL	0.05153 mL	0.04814 mL	0.01999 mL	4.080	-0.00544	0.86041	0.00029	10.0 s
	Data point 12	1.50000 mL	0.05153 mL	0.04906 mL	0.01999 mL	4.274	-0.01020	0.91768	0.00053	10.5 s
	Data point 13	1.50000 mL	0.05153 mL	0.04974 mL	0.01999 mL	4.463	-0.01138	0.78057	0.00064	
19:52.2	Data point 14				0.01999 mL				0.00075	10.0 s
	Data point 15	1.50000 mL	0.05153 mL	0.05054 mL	0.01999 mL	4.848	-0.01852	0.92972	0.00095	13.5 s
21:06.5	Data point 16				0.01999 mL				0.00094	12.0 s
	Data point 17				0.01999 mL				0.00099	22.5 s
22:47.2	Data point 18				0.01999 mL				0.00094	40.0 s
23:57.7	Data point 19				0.01999 mL				0.00091	44.5 s
25:07.7					0.01999 mL				0.00099	44.5 s
	Data point 21				0.01999 mL				0.00092	
27:39.8	Data point 22	1.50000 mL	0.05153 mL	0.05134 mL	0.01999 mL	7.465	-0.05629	0.99724	0.00278	Timed out at
										59.5 s
29:10.2	Data point 23	1.50000 mL	0.05153 mL	0.05139 mL	0.01999 mL	7.984	-0.04354	0.97478	0.00218	Timed out at
										59.5 s
30:40.7	Data point 24	1.50000 mL	0.05153 mL	0.05143 mL	0.01999 mL	8.234	-0.03056	0.98451	0.00152	Timed out at
										59.5 s
	Data point 25				0.01999 mL				0.00097	
33:46.5	Data point 26				0.01999 mL				0.00098	
	Data point 27				0.01999 mL				0.00096	35.0 s
	Data point 28				0.01999 mL				0.00096	
	Data point 29				0.06999 mL				0.00053	
	Data point 30				0.06999 mL				0.00085	
	Data point 31				0.06999 mL			0.00320	0.00029	
39:26.6					0.06999 mL				0.00013	
40:02.1	Data point 33				0.06999 mL				0.00018	
40:37.6	Data point 34				0.06999 mL				0.00039	
41:13.0					0.06999 mL				0.00016	
41:48.5	Data point 36	1.50000 mL	0.10802 mL	0.10146 mL	0.06999 mL	3.509	-0.00677	0.28436	0.00063	10.0 s

1.50000 mL 0.10802 mL 0.10278 mL 0.06999 mL 3.709 -0.00300 0.43965

Reported at: 3/2/2018 1:46:50 PM

42:34.2 Data point 37

0.00022 10.0 s



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

## **Events (continued)**

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
43:09.6	Data point 38	1.50000 mL	0 10802 ml	0.10388 mL	0 06999 ml	3 957	-0 00406	0 29082	0.00037	
43:55.3	Data point 39		0.10802 mL						0.00028	
44:31.3	Data point 40		0.10802 mL						0.00083	
45:11.8	Data point 41		0.10802 mL					0.82458	0.00095	
45:55.4	Data point 42		0.10802 mL					0.55308	0.00087	
46:41.4	Data point 43		0.10802 mL						0.00094	
48:03.1	Data point 44		0.10802 mL							Timed out at 59.5 s
49:33.5	Data point 45	1.50000 mL	0.10802 mL	0.10607 mL	0.06999 mL	7.455	-0.05872	0.98494	0.00292	Timed out at 59.5 s
51:04.0	Data point 46		0.10802 mL						0.00298	at 59.5 s
52:39.7	Data point 47		0.10802 mL							Timed out at 59.5 s
54:15.4	Data point 48		0.10802 mL						0.00095	28.0 s
55:19.1	Data point 49		0.10802 mL					0.93380	0.00095	
56:45.5	Data point 50		0.16959 mL					0.05617	0.00064	
57:31.7	Data point 51		0.16959 mL						0.00069	
58:07.5	Data point 52		0.16959 mL					0.23781	0.00023	
58:43.1	Data point 53		0.16959 mL					0.26467	0.00077	
59:18.7	Data point 54		0.16959 mL					0.60816	0.00015	
59:54.8	Data point 55		0.16959 mL					0.56070	0.00015	
1:00:30.3	Data point 56		0.16959 mL				0.00108	0.00354	0.00090	
1:01:16.1	Data point 57		0.16959 mL					0.78494	0.00027	
	Data point 58		0.16959 mL					0.00046	0.00096	
			0.16959 mL					0.00106	0.00095	
	Data point 60		0.16959 mL					0.02472	0.00073	
	Data point 61		0.16959 mL					0.00712	0.00084	
	Data point 62		0.16959 mL					0.07278	0.00086	
	Data point 63		0.16959 mL					0.04516	0.00094	
	Data point 64		0.16959 mL 0.16959 mL				-0.01776	0.89415 0.97728	0.00093	
1:07:17.1 1:08:28.0	Data point 65		0.16959 mL					0.97726	0.00089	
	Data point 66  Data point 67		0.16959 mL						0.00144	Timed out at 59.5 s 55.0 s
1:11:13.8	Data point 68		0.16959 mL						0.00088	
	·									at 59.5 s
	Data point 69		0.16959 mL							Timed out at 59.5 s
	Data point 70		0.16959 mL							Timed out at 59.5 s
1:15:45.3	Data point 71	1.50000 mL	0.16959 mL	0.16402 mL	0.21999 mL	8.104	-0.04184	0.98451	0.00208	Timed out at 59.5 s
1:17:31.3	Data point 72	1.50000 mL	0.16959 mL	0.16414 mL	0.21999 mL	8.306	-0.03456	0.94740	0.00175	Timed out at 59.5 s
1:19:06.9	Data point 73	1.50000 mL	0.16959 mL	0.16420 mL	0.21999 mL	8.544	-0.01489	0.84815	0.00080	
1:20:16.2	Data point 74		0.16959 mL						0.00099	44.0 s
1:21:36.0	Data point 75 Assay volumes		0.16959 mL 0.16959 mL			9.012	-0.01523	0.58113	0.00099	13.0 s



Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	9.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings	o a a a o o o o o o o o o o o o o o o o			
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose	10 70			
Titrant pre-dose	None			
Assay Medium	NONE			
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type Partition volume	Octanol 0.020 mL			
Partition solvent added				
	Automatic			
After partition addition, stir for	1 seconds			
Sample Sonication	Voo			
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	300 seconds			
After sonication stir for	5 seconds			
Sample Dissolution	V			
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution				
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge				
Perform a carbonate purge	No			
Temperature Control	.,			
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			
Stir speed of	50%			
Titration 1				
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	50%			
Titration 2				
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.050 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			

Reported at: 3/2/2018 1:46:50 PM

Stir to allow partitioning for

Stirrer speed for partitioning

55%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

# Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3				
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.150 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00100 dpH/dt			
Stability timeout after	60 seconds			

# Calibration Settings

Value	Date/Time changed	Imported from
0.130	3/1/2018 1:21:08 AM	C:\Sirius_T3\HCl18B27.t3r
0.9970	3/1/2018 1:21:08 AM	C:\Sirius_T3\HCl18B27.t3r
8.0	3/1/2018 1:21:08 AM	C:\Sirius_T3\HCl18B27.t3r
-0.4	3/1/2018 1:21:08 AM	C:\Sirius_T3\HCl18B27.t3r
1.000	3/1/2018 1:21:08 AM	C:\Sirius_T3\KOH18B27.t3r
0.994	3/1/2018 1:21:08 AM	C:\Sirius_T3\HCl18B27.t3r
	0.130 0.9970 0.8 -0.4 1.000	-0.4 3/1/2018 1:21:08 AM 1.000 3/1/2018 1:21:08 AM

# Instrument Settings

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module Dispenser 0 Syringe volume Firmware version	Water 2.5 mL 1.2.1(r2)	T3DM1200361	3/31/2009 5:24:52 AM 3/31/2009 5:25:05 AM
Titrant Dispenser 2 Syringe volume Firmware version	Water (0.15 M KCI) Acid 0.5 mL 1.2.1(r2)	02-06-2018	2/27/2018 10:05:59 AM 3/31/2009 5:25:11 AM
Titrant Dispenser 1 Syringe volume Firmware version	Acid (0.5 M HCI) Base 0.5 mL 1.2.1(r2)	02-27-2018	2/27/2018 10:27:22 AM 3/31/2009 5:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2)	9/22/2017	2/27/2018 10:21:22 AM 3/31/2009 5:26:24 AM
Distribution valve 5 Firmware version Port A Port B Dispenser 3	Distribution Valve 1.1.3 Methanol (80%, 0.15 M KCI) Cyclohexane Buffer	09-26-17 11-01-17	3/31/2009 5:28:19 AM 2/7/2018 9:42:01 AM 2/27/2018 10:37:57 AM 8/3/2010 5:05:16 AM
Syringe volume Firmware version Titrant Dispenser 6	0.5 mL 1.2.1(r2) Dodecane Octanol	2018/01/31	2/28/2018 10:18:04 AM 10/22/2010 10:52:43 AM

Reported at: 3/2/2018 1:46:50 PM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

# Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Syringe volume Firmware version	0.5 mL 1.2.1(r2)		
Titrant	Octanol	01-31-2018	2/27/2018 9:59:35 AM
Titrator	Octarior		3/31/2009 5:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2	1011111200101	0,0 1,2000 0.2 1.11 7
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0923	1/23/2018 2:01:00 PM
E0 calibration	+3.07 mV	1401.007	3/1/2018 1:21:36 AM
Filling solution	3M KCI	KCL097	2/27/2018 9:49:43 AM
Liquids Wash 1	50% IPA:50% Water		2/28/2018 10:23:32 AM
Wash 2	0.5% Trition X-100 in H20		2/28/2018 10:23:34 AM
Buffer position 1	pH7 Wash		2/28/2018 10:24:06 AM
Buffer position 2	pH 7		2/28/2018 10:24:08 AM
Storage position	•		2/28/2018 10:21:14 AM
Wash water	8.5e+003 mL	02-27-2018	2/27/2018 9:54:39 AM
Waste	7e+003 mL		11/28/2017 10:36:29 AM
Temperature controller			8/5/2010 6:35:13 AM
Turbidity detector		0=1011	3/31/2009 5:24:45 AM
Spectrometer		074811	11/23/2010 11:22:28 AM
Dip probe	183.333	10196	
Wavelength coefficient A0 Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	112:08:55		11/23/2010 11:22:28 AM
Calibrated on	2/27/2018 10:40:38 AM		
Integration time	40		
Scans averaged	10		
Autoloader		T3AL1200345	11/10/2015 9:34:13 AM
Left-right axis firmware version	1.17 Al1DI2DO2 Stepper 2		
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version Chassis I/O firmware version	1.17 Al1Dl2DO2 Stepper 2 1.11 Al1Dl0DO4 Norgren I/O		
Configuration	1.11 Al IDIODO4 Noigien I/O		
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration Flowing wash stir speed	5 s 30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation			
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed E0 calibration buffer wash stir duration	30% 5 s		
E0 calibration buffer wash stir duration E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		
candidater reading our opood			

Reported at: 3/2/2018 1:46:50 PM



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

## Instrument Settings (continued)

Spectrometer calibration stir duration 5 s
opeon official danishation of a
Spectrometer calibration stir speed 30%
Spectrometer calibration wash pump volume 20.0 mL
Spectrometer calibration wash stir duration 5 s
Spectrometer calibration wash stir speed 30%
Overhead dispense height 10000

## Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00
, ·		

### Experiment Log

[2:37] Air gap created for Water (0.15 M KCl)	
[2:37] Air gap created for Acid (0.5 M HCI)	

- [2:38] Air gap created for Base (0.5 M KOH) [2:38] Air gap released for Water (0.15 M KCI)
- [2:42] Titrator arm moved over Titration position
- [2:42] Titration 1 of 3
- [2:42] Adding initial titrants
- [2:42] Automatically add 1.50000 mL of water
- [3:07] Dispensed 1.500000 mL of Water (0.15 M KCI)
- [3:11] Titrator arm moved over Drain
- [8:53] Titrator arm moved to Titration position
- [8:53] Argon flow rate set to 100
- [8:53] Stirrer speed set to 10
- [8:58] Automatically add 0.02000 mL of Octanol
- [8:58] Dispensed 0.019991 mL of Octanol
- [8:59] Initial pH = 8.06
- [8:59] Iterative adjust 8.06 -> 2.00
- [8:59] pH 8.06 -> 2.00
- [9:01] Air gap released for Acid (0.5 M HCI)
- [9:02] Dispensed 0.051529 mL of Acid (0.5 M HCI)
- [9:07] Holding pH 2.00
- [11:07] Stirrer speed set to 0
- [11:07] Stirrer speed set to 50
- [11:07] Iterative adjust 1.98 -> 2.00
- [11:07] pH 1.98 -> 2.00
- [11:08] Air gap released for Base (0.5 M KOH)
- [11:09] Dispensed 0.001787 mL of Base (0.5 M KOH)
- [11:59] Stirrer speed set to 0
- [12:09] Datapoint id 1 collected
- [12:09] Stirrer speed set to 50
- [12:14] pH 2.01 -> 2.21
- [12:14] Using cautious pH adjust
- [12:14] Dispensed 0.007855 mL of Base (0.5 M KOH)
- [12:19] Stepping pH = 2.09
- [12:20] Dispensed 0.006232 mL of Base (0.5 M KOH)
- [12:25] Stepping pH = 2.19
- [12:25] Dispensed 0.001129 mL of Base (0.5 M KOH)
- [12:30] Stepping pH = 2.21
- [12:45] Stirrer speed set to 0
- [12:55] Datapoint id 2 collected
- [12:55] Charge balance equation is out by 3.1%
- [12:55] Stirrer speed set to 50



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [13:00] pH 2.22 -> 2.42
- [13:00] Using charge balance adjust
- [13:01] Dispensed 0.009713 mL of Base (0.5 M KOH)
- [13:21] Stirrer speed set to 0
- [13:31] Datapoint id 3 collected
- [13:31] Charge balance equation is out by 7.5%
- [13:31] Stirrer speed set to 50
- [13:36] pH 2.44 -> 2.64
- [13:36] Using charge balance adjust
- [13:36] Dispensed 0.006021 mL of Base (0.5 M KOH)
- [13:56] Stirrer speed set to 0
- [14:06] Datapoint id 4 collected
- [14:06] Charge balance equation is out by -3.5%
- [14:06] Stirrer speed set to 50
- [14:11] pH 2.64 -> 2.84
- [14:11] Using charge balance adjust
- [14:12] Dispensed 0.004069 mL of Base (0.5 M KOH)
- [14:32] Stirrer speed set to 0
- [14:42] Datapoint id 5 collected
- [14:42] Charge balance equation is out by 10.1%
- [14:42] Stirrer speed set to 50
- [14:47] pH 2.86 -> 3.06
- [14:47] Using charge balance adjust
- [14:47] Dispensed 0.002846 mL of Base (0.5 M KOH)
- [15:07] Stirrer speed set to 0
- [15:17] Datapoint id 6 collected
- [15:17] Charge balance equation is out by -2.6%
- [15:17] Stirrer speed set to 50
- [15:22] pH 3.07 -> 3.27
- [15:22] Using charge balance adjust
- [15:23] Dispensed 0.002258 mL of Base (0.5 M KOH)
- [15:43] Stirrer speed set to 0
- [15:53] Datapoint id 7 collected
- [15:53] Charge balance equation is out by -2.0%
- [15:53] Stirrer speed set to 50
- [15:58] pH 3.27 -> 3.47
- [15:58] Using charge balance adjust
- [15:58] Dispensed 0.001929 mL of Base (0.5 M KOH)
- [16:18] Stirrer speed set to 0
- [16:28] Datapoint id 8 collected
- [16:28] Charge balance equation is out by 8.1%
- [16:28] Stirrer speed set to 50
- [16:33] pH 3.49 -> 3.69
- [16:33] Using charge balance adjust
- [16:34] Dispensed 0.001693 mL of Base (0.5 M KOH)
- [16:54] Stirrer speed set to 0
- [17:04] Datapoint id 9 collected
- [17:04] Charge balance equation is out by 8.4%
- [17:04] Stirrer speed set to 50
- [17:09] pH 3.71 -> 3.91
- [17:09] Using charge balance adjust
- [17:09] Dispensed 0.001435 mL of Base (0.5 M KOH)
- [17:29] Stirrer speed set to 0
- [17:40] Datapoint id 10 collected
- [17:40] Charge balance equation is out by -7.0%
- [17:40] Stirrer speed set to 50
- [17:45] pH 3.90 -> 4.10
- [17:45] Using charge balance adjust
- [17:45] Dispensed 0.001176 mL of Base (0.5 M KOH)



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [18:05] Stirrer speed set to 0
- [18:15] Datapoint id 11 collected
- [18:15] Charge balance equation is out by -8.8%
- [18:15] Stirrer speed set to 50
- [18:20] pH 4.09 -> 4.29
- [18:20] Using charge balance adjust
- [18:20] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [18:41] Stirrer speed set to 0
- [18:51] Datapoint id 12 collected
- [18:51] Charge balance equation is out by -6.3%
- [18:51] Stirrer speed set to 50
- [18:56] pH 4.28 -> 4.48
- [18:56] Using charge balance adjust
- [18:56] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [19:16] Stirrer speed set to 0
- [19:26] Datapoint id 13 collected
- [19:26] Charge balance equation is out by -8.5%
- [19:26] Stirrer speed set to 50
- [19:32] pH 4.47 -> 4.67
- [19:32] Using charge balance adjust
- [19:32] Dispensed 0.000470 mL of Base (0.5 M KOH)
- [19:52] Stirrer speed set to 0
- [20:02] Datapoint id 14 collected
- [20:02] Charge balance equation is out by -10.3%
- [20:02] Stirrer speed set to 50
- [20:07] pH 4.65 -> 4.85
- [20:07] Using charge balance adjust
- [20:07] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [20:27] Stirrer speed set to 0
- [20:41] Datapoint id 15 collected
- [20:41] Charge balance equation is out by -2.0%
- [20:41] Stirrer speed set to 50
- [20:46] pH 4.86 -> 5.06
- [20:46] Using charge balance adjust
- [20:46] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [21:06] Stirrer speed set to 0
- [21:18] Datapoint id 16 collected
- [21:18] Charge balance equation is out by -17.2%
- [21:18] Stirrer speed set to 50
- [21:23] pH 5.03 -> 5.23
- [21:23] Using cautious pH adjust
- [21:23] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [21:28] Stepping pH = 5.04
- [21:29] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [21:34] Stepping pH = 5.67 [21:49] Stirrer speed set to 0
- [22:11] Datapoint id 17 collected
- [22:11] Charge balance equation is out by -93.1%
- [22:11] Stirrer speed set to 50
- [22:16] pH 5.82 -> 6.02
- [22:16] Using cautious pH adjust
- [22:16] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [22:21] Stepping pH = 5.84
- [22:22] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [22:27] Stepping pH = 5.97
- [22:27] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [22:32] Stepping pH = 6.20
- [22:47] Stirrer speed set to 0
- [23:27] Datapoint id 18 collected



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [23:27] Charge balance equation is out by -130.2%
- [23:27] Stirrer speed set to 50
- [23:32] pH 6.32 -> 6.52
- [23:32] Using cautious pH adjust
- [23:32] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [23:37] Stepping pH = 6.38
- [23:37] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [23:42] Stepping pH = 6.54
- [23:57] Stirrer speed set to 0
- [24:42] Datapoint id 19 collected
- [24:42] Charge balance equation is out by -11.0%
- [24:42] Stirrer speed set to 50
- [24:47] pH 6.49 -> 6.69
- [24:47] Using charge balance adjust
- [24:47] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [25:07] Stirrer speed set to 0
- [25:52] Datapoint id 20 collected
- [25:52] Charge balance equation is out by 7.8%
- [25:52] Stirrer speed set to 50
- [25:57] pH 6.69 -> 6.89
- [25:57] Using charge balance adjust
- [25:57] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [26:17] Stirrer speed set to 0
- [27:09] Datapoint id 21 collected
- [27:09] Charge balance equation is out by 15.7%
- [27:09] Stirrer speed set to 50
- [27:14] pH 6.96 -> 7.16
- [27:14] Using cautious pH adjust
- [27:14] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [27:19] Stepping pH = 7.08
- [27:19] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [27:24] Stepping pH = 7.29
- [27:40] Stirrer speed set to 0
- [28:40] Datapoint id 22 collected
- [28:40] Charge balance equation is out by -110.3%
- [28:40] Stirrer speed set to 50
- [28:45] pH 7.43 -> 7.63
- [28:45] Using cautious pH adjust
- [28:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [28:50] Stepping pH = 7.47
- [28:50] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [28:55] Stepping pH = 7.66
- [29:10] Stirrer speed set to 0
- [30:10] Datapoint id 23 collected
- [30:10] Charge balance equation is out by -359.3%
- [30:10] Stirrer speed set to 50
- [30:15] pH 7.88 -> 8.08
- [30:15] Using cautious pH adjust
- [30:15] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [30:20] Stepping pH = 7.89
- [30:20] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [30.25] Stepping pH = 8.10
- [30:40] Stirrer speed set to 0
- [31:40] Datapoint id 24 collected
- [31:40] Charge balance equation is out by -547.1%
- [31:40] Stirrer speed set to 50
- [31:46] pH 8.15 -> 8.35
- [31:46] Using cautious pH adjust
- [31:46] Dispensed 0.000024 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [31:51] Stepping pH = 8.14
- [31:51] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [31:56] Stepping pH = 8.26
- [31:56] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [32:01] Stepping pH = 8.46
- [32:16] Stirrer speed set to 0
- [33:11] Datapoint id 25 collected
- [33:11] Charge balance equation is out by -675.1%
- [33:11] Stirrer speed set to 50
- [33:16] pH 8.52 -> 8.72
- [33:16] Using cautious pH adjust
- [33:16] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [33:21] Stepping pH = 8.57
- [33:21] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [33:26] Stepping pH = 8.67
  - 0.20] Stepping pri = 0.07
- [33:26] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [33:31] Stepping pH = 8.76
- [33:46] Stirrer speed set to 0
- [34:27] Datapoint id 26 collected
- 74.27 Datapoint la 20 concette
- [34:27] Charge balance equation is out by -301.4%
- [34:27] Stirrer speed set to 50
- [34:32] pH 8.73 -> 8.93
- [34:32] Using cautious pH adjust
- [34:32] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [34:37] Stepping pH = 8.74
- [34:37] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [34:42] Stepping pH = 8.81
- [34:42] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [34:47] Stepping pH = 8.96
- [35:02] Stirrer speed set to 0
- [35:37] Datapoint id 27 collected
- [35:37] Charge balance equation is out by -310.7%
- [35:37] Stirrer speed set to 50
- [35:43] pH 8.95 -> 9.05
- [35:43] Using cautious pH adjust
- [35:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [35:48] Stepping pH = 8.95
- [35:48] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [35:53] Stepping pH = 9.03
- [36:08] Stirrer speed set to 0
- [36:29] Datapoint id 28 collected
- [36:29] Charge balance equation is out by -289.1%
- [36:29] Titration 2 of 3
- [36:29] Adding initial titrants
- [36:29] Automatically add 0.05000 mL of Octanol
- [36:31] Dispensed 0.050000 mL of Octanol
- [36:31] Stirrer speed set to 10
- [36:32] Stirrer speed set to 55
- [36:32] Iterative adjust 9.08 -> 2.00
- [36:32] pH 9.08 -> 2.00
- [36:33] Dispensed 0.054139 mL of Acid (0.5 M HCI)
- 36:38 pH 2.02 -> 2.00
- [36:38] Dispensed 0.002352 mL of Acid (0.5 M HCl)
- [37:29] Stirrer speed set to 0
- [37:39] Datapoint id 29 collected
- [37:39] Stirrer speed set to 55
- [37:44] pH 1.97 -> 2.17
- [37:44] Using cautious pH adjust
- [37:44] Dispensed 0.009313 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [37:49] Stepping pH = 2.06
- [37:50] Dispensed 0.006538 mL of Base (0.5 M KOH)
- [37:55] Stepping pH = 2.15
- [37:55] Dispensed 0.001458 mL of Base (0.5 M KOH)
- [38:00] Stepping pH = 2.17
- [38:15] Stirrer speed set to 0
- [38:25] Datapoint id 30 collected
- [38:25] Charge balance equation is out by 7.0%
- [38:25] Stirrer speed set to 55
- [38:30] pH 2.17 -> 2.37
- [38:30] Using charge balance adjust
- [38:30] Dispensed 0.011783 mL of Base (0.5 M KOH)
- [38:51] Stirrer speed set to 0
- [39:01] Datapoint id 31 collected
- [39:01] Charge balance equation is out by 12.3%
- [39:01] Stirrer speed set to 55
- [39:06] pH 2.40 -> 2.60
- [39:06] Using charge balance adjust
- [39:06] Dispensed 0.007291 mL of Base (0.5 M KOH)
- [39:26] Stirrer speed set to 0
- [39:36] Datapoint id 32 collected
- [39:36] Charge balance equation is out by 5.2%
- [39:36] Stirrer speed set to 55
- [39:41] pH 2.62 -> 2.82
- [39:41] Using charge balance adjust
- [39:42] Dispensed 0.004892 mL of Base (0.5 M KOH)
- [40:02] Stirrer speed set to 0
- [40:12] Datapoint id 33 collected
- [40:12] Charge balance equation is out by 2.4%
- [40:12] Stirrer speed set to 55
- [40:17] pH 2.83 -> 3.03
- [40:17] Using charge balance adjust
- [40:17] Dispensed 0.003551 mL of Base (0.5 M KOH)
- [40:37] Stirrer speed set to 0
- [40:47] Datapoint id 34 collected
- [40:47] Charge balance equation is out by -4.4%
- [40:47] Stirrer speed set to 55
- [40:52] pH 3.02 -> 3.22
- [40:52] Using charge balance adjust
- [40:52] Dispensed 0.002752 mL of Base (0.5 M KOH)
- [41:13] Stirrer speed set to 0
- [41:23] Datapoint id 35 collected
- [41:23] Charge balance equation is out by 5.0%
- [41:23] Stirrer speed set to 55
- [41:28] pH 3.24 -> 3.44
- [41:28] Using charge balance adjust
- [41:28] Dispensed 0.002140 mL of Base (0.5 M KOH)
- [41:48] Stirrer speed set to 0
- [41:58] Datapoint id 36 collected
- [41:58] Charge balance equation is out by 37.0%
- [41:58] Stirrer speed set to 55
- [42:03] pH 3.51 -> 3.71
- [42:03] Using cautious pH adjust
- [42:03] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [42:08] Stepping pH = 3.61
- [42:09] Dispensed 0.000517 mL of Base (0.5 M KOH)
- [42:14] Stepping pH = 3.70
- [42:14] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [42:19] Stepping pH = 3.71



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

## Experiment Log (continued)

- [42:34] Stirrer speed set to 0
- [42:44] Datapoint id 37 collected
- [42:44] Charge balance equation is out by 11.6%
- [42:44] Stirrer speed set to 55
- [42:49] pH 3.71 -> 3.91
- [42:49] Using charge balance adjust
- [42:49] Dispensed 0.001105 mL of Base (0.5 M KOH)
- [43:09] Stirrer speed set to 0
- [43:19] Datapoint id 38 collected
- [43:19] Charge balance equation is out by 21.3%
- [43:19] Stirrer speed set to 55
- [43:24] pH 3.96 -> 4.16
- [43:24] Using cautious pH adjust
- [43:24] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [43:30] Stepping pH = 4.04
- [43:30] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [43:35] Stepping pH = 4.14
- [43:35] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [43:40] Stepping pH = 4.16
- [43:55] Stirrer speed set to 0
- [44:06] Datapoint id 39 collected
- [44:06] Charge balance equation is out by -1.4%
- [44:06] Stirrer speed set to 55
- [44:11] pH 4.17 -> 4.37
- [44:11] Using charge balance adjust
- [44:11] Dispensed 0.000494 mL of Base (0.5 M KOH)
- [44:31] Stirrer speed set to 0
- [44:41] Datapoint id 40 collected
- [44:41] Charge balance equation is out by 26.0%
- [44:41] Stirrer speed set to 55
- [44:46] pH 4.42 -> 4.62
- [44:46] Using cautious pH adjust
- [44:46] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [44:51] Stepping pH = 4.46
- [44:51] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [44:56] Stepping pH = 4.72
- [45:12] Stirrer speed set to 0
- [45:25] Datapoint id 41 collected
- [45:25] Charge balance equation is out by -38.6%
- [45:25] Stirrer speed set to 55
- [45:30] pH 4.78 -> 4.98
- [45:30] Using cautious pH adjust
- [45:30] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [45:35] Stepping pH = 4.79
- [45:35] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [45:40] Stepping pH = 5.22 [45:55] Stirrer speed set to 0
- [46:11] Datapoint id 42 collected
- [46:11] Charge balance equation is out by -92.6%
- [46:11] Stirrer speed set to 55
- [46:16] pH 5.29 -> 5.49
- [46:16] Using cautious pH adjust
- [46:16] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [46:21] Stepping pH = 5.29
- [46:21] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [46:26] Stepping pH = 5.90
- [46:41] Stirrer speed set to 0
- [47:27] Datapoint id 43 collected
- [47:27] Charge balance equation is out by -206.7%

Reported at: 3/2/2018 1:46:50 PM



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [47:27] Stirrer speed set to 55
- [47:32] pH 6.12 -> 6.32
- [47:32] Using cautious pH adjust
- [47:32] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [47:37] Stepping pH = 6.14
- [47:37] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [47:43] Stepping pH = 6.30
- [47:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [47:48] Stepping pH = 6.62
- [48:03] Stirrer speed set to 0 [49:03] Datapoint id 44 collected
- [49:03] Charge balance equation is out by -118.6%
- [49:03] Stirrer speed set to 55
- [49:08] pH 6.74 -> 6.94
- [49:08] Using cautious pH adjust
- [49:08] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [49:13] Stepping pH = 6.76
- [49:13] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [49:18] Stepping pH = 7.01
- [49:33] Stirrer speed set to 0
- [50:33] Datapoint id 45 collected
- [50:33] Charge balance equation is out by -103.4%
- [50:33] Stirrer speed set to 55
- [50:38] pH 7.33 -> 7.53
- [50:38] Using cautious pH adjust
- [50:38] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [50:43] Stepping pH = 7.29
- [50:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [50:49] Stepping pH = 7.56
- [51:04] Stirrer speed set to 0
- [52:04] Datapoint id 46 collected
- [52:04] Charge balance equation is out by -392.4%
- [52:04] Stirrer speed set to 55
- [52:09] pH 7.78 -> 7.98
- [52:09] Using cautious pH adjust
- [52:09] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [52:14] Stepping pH = 7.73
- [52:14] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [52:19] Stepping pH = 7.78
- [52:19] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [52:24] Stepping pH = 8.11
- [52:39] Stirrer speed set to 0
- [53:39] Datapoint id 47 collected
- [53:39] Charge balance equation is out by -1,044.0%
- [53:39] Stirrer speed set to 55
- [53:45] pH 7.92 -> 8.12
- [53:45] Using cautious pH adjust
- [53:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [53:50] Stepping pH = 7.88
- [53:50] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [53:55] Stepping pH = 7.86
- [53:55] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [54:00] Stepping pH = 8.83
- [54:15] Stirrer speed set to 0
- [54:43] Datapoint id 48 collected
- [54:43] Charge balance equation is out by -2,120.7%
- [54:43] Stirrer speed set to 55
- [54:48] pH 8.82 -> 9.02
- [54:48] Using cautious pH adjust



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [54:48] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [54:53] Stepping pH = 8.83
- [54:53] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [54:58] Stepping pH = 8.90
- [54:59] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [55:04] Stepping pH = 9.03
- [55:19] Stirrer speed set to 0
- [55:43] Datapoint id 49 collected
- [55:43] Charge balance equation is out by -286.1%
- [55:43] Titration 3 of 3
- [55:43] Adding initial titrants
- [55:43] Automatically add 0.15000 mL of Octanol
- [55:47] Dispensed 0.150000 mL of Octanol
- [55:47] Stirrer speed set to 10
- [55:48] Stirrer speed set to 60
- [55:48] Iterative adjust 9.05 -> 2.00
- [55:48] pH 9.05 -> 2.00
- [55:50] Dispensed 0.056397 mL of Acid (0.5 M HCI)
- [55:55] pH 2.05 -> 2.00
- [55:55] Dispensed 0.005174 mL of Acid (0.5 M HCI)
- [56:45] Stirrer speed set to 0
- [56:55] Datapoint id 50 collected
- [56:55] Stirrer speed set to 60
- [57:00] pH 1.97 -> 2.17
- 7.00] pri 1.07 2.17
- [57:00] Using cautious pH adjust [57:01] Dispensed 0.010113 mL of Base (0.5 M KOH)
- [57:06] Stepping pH = 2.05
- 77.00] Stepping pri = 2.03
- [57:06] Dispensed 0.007949 mL of Base (0.5 M KOH)
- [57:11] Stepping pH = 2.15
- [57:11] Dispensed 0.001129 mL of Base (0.5 M KOH)
- [57:16] Stepping pH = 2.17
- [57:31] Stirrer speed set to 0
- [57:41] Datapoint id 51 collected
- [57:41] Charge balance equation is out by 5.1%
- [57:41] Stirrer speed set to 60
- [57:47] pH 2.17 -> 2.37
- [57:47] Using charge balance adjust
- [57:47] Dispensed 0.012959 mL of Base (0.5 M KOH)
- [58:07] Stirrer speed set to 0
- [58:17] Datapoint id 52 collected
- [58:17] Charge balance equation is out by 9.6%
- [58:17] Stirrer speed set to 60
- [58:22] pH 2.40 -> 2.60
- [58:22] Using charge balance adjust
- [58:23] Dispensed 0.008278 mL of Base (0.5 M KOH)
- [58:43] Stirrer speed set to 0
- [58:53] Datapoint id 53 collected
- [58:53] Charge balance equation is out by 10.7%
- [58:53] Stirrer speed set to 60
- [58:58] pH 2.62 -> 2.82
- [58:58] Using charge balance adjust
- [58:58] Dispensed 0.005433 mL of Base (0.5 M KOH)
- [59:18] Stirrer speed set to 0
- [59:29] Datapoint id 54 collected
- [59:29] Charge balance equation is out by -3.5%
- [59:29] Stirrer speed set to 60
- [59:34] pH 2.82 -> 3.02
- [59:34] Using charge balance adjust
- [59:34] Dispensed 0.003857 mL of Base (0.5 M KOH)



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [59:55] Stirrer speed set to 0
- [1:00:05] Datapoint id 55 collected
- [1:00:05] Charge balance equation is out by 4.7%
- [1:00:05] Stirrer speed set to 60
- [1:00:10] pH 3.04 -> 3.24
- [1:00:10] Using charge balance adjust
- [1:00:10] Dispensed 0.002681 mL of Base (0.5 M KOH)
- [1:00:30] Stirrer speed set to 0
- [1:00:40] Datapoint id 56 collected
- [1:00:40] Charge balance equation is out by 20.3%
- [1:00:40] Stirrer speed set to 60
- [1:00:45] pH 3.28 -> 3.48
- [1:00:45] Using cautious pH adjust
- [1:00:45] Dispensed 0.000870 mL of Base (0.5 M KOH)
- [1:00:50] Stepping pH = 3.37
- [1:00:50] Dispensed 0.000706 mL of Base (0.5 M KOH)
- [1:00:56] Stepping pH = 3.47
- [1:00:56] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:01:01] Stepping pH = 3.49
- [1:01:16] Stirrer speed set to 0
- [1:01:26] Datapoint id 57 collected
- [1:01:26] Charge balance equation is out by 5.0%
- [1:01:26] Stirrer speed set to 60
- [1:01:31] pH 3.49 -> 3.69
- [1:01:31] Using charge balance adjust
- [1:01:31] Dispensed 0.001152 mL of Base (0.5 M KOH)
- [1:01:51] Stirrer speed set to 0
- [1:02:01] Datapoint id 58 collected
- [1:02:01] Charge balance equation is out by 16.7%
- [1:02:01] Stirrer speed set to 60
- [1:02:06] pH 3.73 -> 3.93
- [1:02:06] Using cautious pH adjust
- [1:02:07] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [1:02:12] Stepping pH = 3.78
- [1:02:12] Dispensed 0.000564 mL of Base (0.5 M KOH)
- [1:02:17] Stepping pH = 3.98
- [1:02:32] Stirrer speed set to 0
- [1:02:42] Datapoint id 59 collected
- [1:02:42] Charge balance equation is out by -28.5%
- [1:02:42] Stirrer speed set to 60
- [1:02:47] pH 4.01 -> 4.21
- [1:02:47] Using cautious pH adjust
- [1:02:47] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [1:02:52] Stepping pH = 4.05
- [1:02:52] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [1:02:57] Stepping pH = 4.32
- [1:03:12] Stirrer speed set to 0
- [1:03:22] Datapoint id 60 collected
- [1:03:22] Charge balance equation is out by -40.2%
- [1:03:22] Stirrer speed set to 60
- [1:03:28] pH 4.36 -> 4.56
- [1:03:28] Using cautious pH adjust
- [1:03:28] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:03:33] Stepping pH = 4.38
- [1:03:33] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [1:03:38] Stepping pH = 4.73
- [1:03:53] Stirrer speed set to 0
- [1:04:05] Datapoint id 61 collected
- [1:04:05] Charge balance equation is out by -65.1%



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [1:04:05] Stirrer speed set to 60
- [1:04:10] pH 4.80 -> 5.00
- [1:04:10] Using cautious pH adjust
- [1:04:10] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:04:15] Stepping pH = 4.83
- [1:04:15] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:04:20] Stepping pH = 4.94
- [1:04:21] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:04:26] Stepping pH = 5.04
- [1:04:41] Stirrer speed set to 0
- [1:04:54] Datapoint id 62 collected
- [1:04:54] Charge balance equation is out by -98.4%
- [1:04:54] Stirrer speed set to 60
- [1:04:59] pH 5.08 -> 5.28
- [1:04:59] Using cautious pH adjust
- [1:04:59] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:05:04] Stepping pH = 5.13
- [1:05:05] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:05:10] Stepping pH = 5.22
- [1:05:10] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:05:15] Stepping pH = 5.32
- [1:05:30] Stirrer speed set to 0
- [1:05:44] Datapoint id 63 collected
- [1:05:44] Charge balance equation is out by -64.2%
- [1:05:44] Stirrer speed set to 60
- [1:05:50] pH 5.43 -> 5.63
- [1:05:50] Using cautious pH adjust
- [1:05:50] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:05:55] Stepping pH = 5.51
- [1:05:55] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:06:00] Stepping pH = 5.65
- [1:06:15] Stirrer speed set to 0
- [1:06:51] Datapoint id 64 collected
- [1:06:51] Charge balance equation is out by -14.9%
- [1:06:51] Stirrer speed set to 60
- [1:06:57] pH 5.73 -> 5.93
- [1:06:57] Using charge balance adjust
- [1:06:57] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:07:17] Stirrer speed set to 0
- [1:07:57] Datapoint id 65 collected
- [1:07:57] Charge balance equation is out by 69.9%
- [1:07:57] Stirrer speed set to 60
- [1:08:02] pH 6.06 -> 6.26
- [1:08:02] Using cautious pH adjust
- [1:08:02] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:08:08] Stepping pH = 6.22
- [1:08:08] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:08:13] Stepping pH = 6.45
- [1:08:28] Stirrer speed set to 0
- [1:09:28] Datapoint id 66 collected
- [1:09:28] Charge balance equation is out by 4.3%
- [1:09:28] Stirrer speed set to 60
- [1:09:33] pH 6.45 -> 6.65
- [1:09:33] Using charge balance adjust
- [1:09:33] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:09:53] Stirrer speed set to 0
- [1:10:48] Datapoint id 67 collected
- [1:10:48] Charge balance equation is out by -24.7%
- [1:10:48] Stirrer speed set to 60



Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [1:10:53] pH 6.58 -> 6.78
- [1:10:53] Using cautious pH adjust
- [1:10:53] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:10:58] Stepping pH = 6.77
- [1:11:13] Stirrer speed set to 0 [1:12:13] Datapoint id 68 collected
- [1:12:13] Charge balance equation is out by 47.7%
- [1:12:13] Stirrer speed set to 60
- [1:12:19] pH 6.87 -> 7.07
- [1:12:19] Using cautious pH adjust
- [1:12:19] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:12:24] Stepping pH = 7.05
- [1:12:24] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:12:29] Stepping pH = 7.33
- [1:12:44] Stirrer speed set to 0
- [1:13:44] Datapoint id 69 collected
- [1:13:44] Charge balance equation is out by -44.9%
- [1:13:44] Stirrer speed set to 60
- [1:13:49] pH 7.36 -> 7.56
- [1:13:49] Using cautious pH adjust
- [1:13:49] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:13:54] Stepping pH = 7.45
- [1:13:54] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:13:59] Stepping pH = 7.73
- [1:14:14] Stirrer speed set to 0
- [1:15:14] Datapoint id 70 collected
- [1:15:14] Charge balance equation is out by -212.1%
- [1:15:14] Stirrer speed set to 60
- [1:15:20] pH 7.77 -> 7.97
- [1:15:20] Using cautious pH adjust
- [1:15:20] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:15:25] Stepping pH = 7.85
- [1:15:25] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:15:30] Stepping pH = 8.04
- [1:15:45] Stirrer speed set to 0
- [1:16:45] Datapoint id 71 collected
- [1:16:45] Charge balance equation is out by -401.0%
- [1:16:45] Stirrer speed set to 60
- [1:16:50] pH 8.05 -> 8.25
- [1:16:50] Using cautious pH adjust
- [1:16:50] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:16:55] Stepping pH = 8.15
- [1:16:55] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:17:00] Stepping pH = 8.21
- [1:17:00] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:17:06] Stepping pH = 8.19
- [1:17:06] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:17:11] Stepping pH = 8.20
- [1:17:11] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:17:16] Stepping pH = 8.32
- [1:17:31] Stirrer speed set to 0
- [1:18:31] Datapoint id 72 collected
- [1:18:31] Charge balance equation is out by -1,251.5%
- [1:18:31] Stirrer speed set to 60
- [1:18:36] pH 8.28 -> 8.48
- [1:18:36] Using cautious pH adjust
- [1:18:36] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:18:41] Stepping pH = 8.27
- [1:18:41] Dispensed 0.000024 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-01002 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180228\_exp28\_logP\_T3-2\18C-01002\_M11\_octanol\_pH-metric high logP.t3r

- [1:18:46] Stepping pH = 8.39
- [1:18:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:18:51] Stepping pH = 8.53
- 1:19:07 Stirrer speed set to 0
- [1:19:40] Datapoint id 73 collected
- [1:19:40] Charge balance equation is out by -521.2%
- [1:19:40] Stirrer speed set to 60
- [1:19:45] pH 8.59 -> 8.79
- [1:19:45] Using cautious pH adjust
- [1:19:45] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:19:51] Stepping pH = 8.65
- [1:19:51] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:19:56] Stepping pH = 8.73
- [1:19:56] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:20:01] Stepping pH = 8.81
- [1:20:16] Stirrer speed set to 0
- [1:21:00] Datapoint id 74 collected
- [1:21:00] Charge balance equation is out by -196.5%
- [1:21:00] Stirrer speed set to 60
- [1:21:05] pH 8.79 -> 8.99
- [1:21:05] Using cautious pH adjust
- [1:21:05] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:21:10] Stepping pH = 8.78
- [1:21:10] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [1:21:15] Stepping pH = 8.93
- [1:21:15] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:21:21] Stepping pH = 9.00
- [1:21:36] Stirrer speed set to 0
- [1:21:49] Datapoint id 75 collected
- [1:21:49] Charge balance equation is out by -327.4%
- [1:21:49] Argon flow rate set to 0
- [1:21:53] Titrator arm moved over Titration position