

Sample name: M04_octanol Experiment start time: 3/24/2018 2:50:31 AM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

18C-24003 Instrument ID: **T312060**

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH +) 0.87 ±0.03 (n=50) logP (neutral X) 3.95 ±0.01 (n=50)

RMSD 0.342

18C-24003 Points 1 to 26

M04_octanol concentration factor 0.985
Carbonate 0.1126 mM
Acidity error 0.07154 mM

18C-24003 Points 27 to 52

M04_octanol concentration factor 0.875
Carbonate 0.1064 mM
Acidity error -0.19519 mM

18C-24003 Points 53 to 77

M04_octanol concentration factor 0.960
Carbonate 0.1188 mM
Acidity error 0.08198 mM

Warnings and errors

Errors None Warnings None

Sample logD and percent species

рH	M04_octanol	M04_octanol	M04_octanol	M04_octanol	M04_octanol	Comment
	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.88	11.71 %	0.00 %	87.17 %	1.12 %	
1.200	0.88	11.63 %	0.00 %	86.60 %	1.77 %	Stomach pH
2.000	0.92	10.63 %	0.00 %	79.16 %	10.20 %	
3.000	1.23	5.54 %	0.01 %	41.26 %	53.19 %	
4.000	2.01	0.96 %	0.01 %	7.13 %	91.90 %	
5.000	2.94	0.10 %	0.01 %	0.77 %	99.12 %	
6.000	3.67	0.01 %	0.01 %	0.08 %	99.90 %	
6.500	3.84	0.00 %	0.01 %	0.02 %	99.96 %	
7.000	3.91	0.00 %	0.01 %	0.01 %	99.98 %	
7.400	3.94	0.00 %	0.01 %	0.00 %	99.99 %	Blood pH
8.000	3.95	0.00 %	0.01 %	0.00 %	99.99 %	
9.000	3.95	0.00 %	0.01 %	0.00 %	99.99 %	
10.000	3.95	0.00 %	0.01 %	0.00 %	99.99 %	
11.000	3.95	0.00 %	0.01 %	0.00 %	99.99 %	
12.000	3.95	0.00 %	0.01 %	0.00 %	99.99 %	



Sample name: M04_octanol Assay name: Assay ID:


pH-metric high logP

18C-24003

Experiment start time: 3/24/2018 2:50:31 AM Analyst: **Dorothy Levorse**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r






Assay ID: Filename:

Sample name: M04_octanol Experiment start time: 3/24/2018 2:50:31 AM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

18C-24003 Instrument ID: **T312060**

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

Graphs (continued)





Sample name: M04_octanol Experiment start time: 3/24/2018 2:50:31 AM pH-metric high logP Analyst: Assay name: **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-24003

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-24003 Points 1 to 26

Overall results

RMSD 0.279 Average ionic strength 0.157 M Average temperature 24.9°C Partition ratio 0.0186:1

Analyte concentration range 2638.8 μM to 2719.9 μM

Total points considered 24 of 26

Warnings and errors

Errors None Warnings None

Four-Plus parameters

Alpha 0.119 3/24/2018 2:50:31 AM C:\Sirius_T3\HCl18C23.t3r S 0.9972 3/24/2018 2:50:31 AM C:\Sirius T3\HCl18C23.t3r jΗ 0.9 3/24/2018 2:50:31 AM C:\Sirius_T3\HCl18C23.t3r jOH -0.33/24/2018 2:50:31 AM C:\Sirius_T3\HCl18C23.t3r


Titrants


0.50 M HCI 0.997124 3/24/2018 2:50:31 AM C:\Sirius T3\HCl18C23.t3r 0.50 M KOH 1.003190 3/24/2018 2:50:31 AM C:\Sirius_T3\KOH18C23.t3r

Sample

M04_octanol concentration factor 0.985 Base pKa 1 5.97 logP(XH +)0.19 logP (neutral X) 3.88

Sample graphs





Analyst:

Experiment start time: 3/24/2018 2:50:31 AM

Dorothy Levorse



Assay ID:

Filename:


Sample name: M04_octanol Assay name: pH-metric high logP

18C-24003

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

Instrument ID: T312060


Sample graphs (continued)




Sample logD and percent species

рН	M04_octanol	M04_octanol	M04_octanol	_	M04_octanol	Comment
	logD	M04_octanolH	_	_	M04_octanol*	
1.000	0.21	97.06 %	0.00 %	2.79 %	0.15 %	
1.200	0.22	96.98 %	0.00 %	2.79 %	0.23 %	Stomach pH
2.000	0.37	95.80 %	0.01 %	2.75 %	1.44 %	
3.000	0.98	84.73 %	0.09 %	2.44 %	12.75 %	
4.000	1.91	39.31 %	0.42 %	1.13 %	59.14 %	
5.000	2.87	6.18 %	0.66 %	0.18 %	92.98 %	
6.000	3.59	0.66 %	0.70 %	0.02 %	98.62 %	
6.500	3.77	0.21 %	0.71 %	0.01 %	99.08 %	
7.000	3.84	0.07 %	0.71 %	0.00 %	99.23 %	
7.400	3.86	0.03 %	0.71 %	0.00 %	99.27 %	Blood pH
8.000	3.87	0.01 %	0.71 %	0.00 %	99.29 %	
9.000	3.88	0.00 %	0.71 %	0.00 %	99.29 %	
10.000	3.88	0.00 %	0.71 %	0.00 %	99.29 %	
11.000	3.88	0.00 %	0.71 %	0.00 %	99.29 %	
12.000	3.88	0.00 %	0.71 %	0.00 %	99.29 %	

Carbonate and acidity



Other graphs





Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP

18C-24003


Experiment start time: 3/24/2018 2:50:31 AM

Analyst: **Dorothy Levorse**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

Other graphs (continued)





Sample name: M04_octanol Experiment start time: 3/24/2018 2:50:31 AM Analyst: Assay name: pH-metric high logP **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-24003

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-24003 Points 27 to 52

Overall results

RMSD 0.447 Average ionic strength 0.163 M 25.0°C Average temperature Partition ratio 0.0643:1

Analyte concentration range 2363.9 µM to 2436.4 µM

Total points considered 25 of 26

Warnings and errors

Errors None Warnings None

Four-Plus parameters

Alpha 0.119 3/24/2018 2:50:31 AM C:\Sirius_T3\HCl18C23.t3r S 0.9972 3/24/2018 2:50:31 AM C:\Sirius T3\HCl18C23.t3r jΗ 0.9 3/24/2018 2:50:31 AM C:\Sirius_T3\HCl18C23.t3r jOH -0.33/24/2018 2:50:31 AM C:\Sirius_T3\HCl18C23.t3r


Titrants


0.50 M HCI 0.997124 3/24/2018 2:50:31 AM C:\Sirius T3\HCl18C23.t3r

Sample

M04_octanol concentration factor 0.875 Base pKa 1 5.97 logP(XH +)0.19 logP (neutral X) 3.90

Sample graphs








Sample name: M04_octanol Experiment start time: 3/24/2018 2:50:31 AM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

18C-24003 Instrument ID: **T312060**

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species


рН	M04_octanol	M04_octanol	M04_octanol	M04_octanol	M04_octanol	Comment
	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.21	90.50 %	0.00 %	9.01 %	0.49 %	
1.200	0.23	90.24 %	0.00 %	8.99 %	0.77 %	Stomach pH
2.000	0.38	86.67 %	0.01 %	8.63 %	4.69 %	
3.000	1.00	60.90 %	0.07 %	6.07 %	32.97 %	
4.000	1.93	15.33 %	0.16 %	1.53 %	82.98 %	
5.000	2.88	1.81 %	0.19 %	0.18 %	97.82 %	
6.000	3.61	0.18 %	0.20 %	0.02 %	99.60 %	
6.500	3.78	0.06 %	0.20 %	0.01 %	99.74 %	
7.000	3.86	0.02 %	0.20 %	0.00 %	99.78 %	
7.400	3.88	0.01 %	0.20 %	0.00 %	99.79 %	Blood pH
8.000	3.89	0.00 %	0.20 %	0.00 %	99.80 %	
9.000	3.89	0.00 %	0.20 %	0.00 %	99.80 %	
10.000	3.90	0.00 %	0.20 %	0.00 %	99.80 %	
11.000	3.90	0.00 %	0.20 %	0.00 %	99.80 %	
12.000	3.90	0.00 %	0.20 %	0.00 %	99.80 %	

Carbonate and acidity



Other graphs







Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP

18C-24003


Experiment start time: 3/24/2018 2:50:31 AM

Analyst: **Dorothy Levorse**

Instrument ID: T312060

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

Other graphs (continued)





Sample name: M04_octanol Experiment start time: 3/24/2018 2:50:31 AM pH-metric high logP Analyst: Assay name: **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-24003

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-24003 Points 53 to 77

Overall results

RMSD 0.260 Average ionic strength 0.168 M 24.9°C Average temperature Partition ratio 0.2796:1

Analyte concentration range 1840.5 μM to 1888.2 μM

Total points considered 23 of 25

Warnings and errors

Errors None Warnings None

Four-Plus parameters

Alpha 0.119 3/24/2018 2:50:31 AM C:\Sirius_T3\HCl18C23.t3r S 0.9972 3/24/2018 2:50:31 AM C:\Sirius T3\HCl18C23.t3r jΗ 0.9 3/24/2018 2:50:31 AM C:\Sirius_T3\HCl18C23.t3r jOH -0.33/24/2018 2:50:31 AM C:\Sirius_T3\HCl18C23.t3r


Titrants


0.50 M HCI 0.997124 3/24/2018 2:50:31 AM C:\Sirius T3\HCl18C23.t3r

Sample

M04 octanol concentration factor 0.960 Base pKa 1 5.97 logP(XH +)0.19 logP (neutral X) 3.59

Sample graphs





Analyst:

Experiment start time: 3/24/2018 2:50:31 AM

Dorothy Levorse



Assay ID:


Sample name: M04_octanol Assay name: pH-metric high logP

18C-24003

Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

Sample graphs (continued)



Sample logD and percent species


рН	M04_octanol	M04_octanol	M04_octanol	M04_octanol	M04_octanol	Comment
	logD	M04_octanolH	M04_octanol	M04_octanolH*	M04_octanol*	
1.000	0.20	69.22 %	0.00 %	29.97 %	0.81 %	
1.200	0.21	68.89 %	0.00 %	29.83 %	1.28 %	Stomach pH
2.000	0.29	64.50 %	0.01 %	27.93 %	7.57 %	
3.000	0.76	38.35 %	0.04 %	16.61 %	45.00 %	
4.000	1.63	7.59 %	0.08 %	3.29 %	89.04 %	
5.000	2.58	0.84 %	0.09 %	0.36 %	98.70 %	
6.000	3.31	0.09 %	0.09 %	0.04 %	99.79 %	
6.500	3.48	0.03 %	0.09 %	0.01 %	99.87 %	
7.000	3.55	0.01 %	0.09 %	0.00 %	99.90 %	
7.400	3.58	0.00 %	0.09 %	0.00 %	99.90 %	Blood pH
8.000	3.59	0.00 %	0.09 %	0.00 %	99.91 %	
9.000	3.59	0.00 %	0.09 %	0.00 %	99.91 %	
10.000	3.59	0.00 %	0.09 %	0.00 %	99.91 %	
11.000	3.59	0.00 %	0.09 %	0.00 %	99.91 %	
12.000	3.59	0.00 %	0.09 %	0.00 %	99.91 %	


Carbonate and acidity



Carbonate 0.119 mM Acidity error 0.082 mM

Other graphs







Filename:

Sample name: M04_octanol Assay name:

pH-metric high logP

18C-24003

C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r


Experiment start time: 3/24/2018 2:50:31 AM Analyst:

Dorothy Levorse

Instrument ID:

T312060

Other graphs (continued)





18C-24003 Instrument ID: **T312060**

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

Assay Model

Assay ID:

Settings	Value	Date/Time changed	Imported from
Sample name	M04_octanol	3/9/2018 4:34:21 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.001160 g	3/23/2018 5:00:38 PM	User entered value
Formula weight	269.73 g/mol	3/9/2018 4:34:21 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	269.73	3/9/2018 4:34:21 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	1	3/9/2018 4:34:21 PM	User entered value
Sample is a	Base	3/9/2018 4:34:21 PM	User entered value
pKa 1	5.97	3/9/2018 4:34:21 PM	User entered value
logp (XH +)	0.19	3/9/2018 4:34:33 PM	User entered value
logP (neutral X)	3.50	3/23/2018 2:32:46 PM	User entered value

Events

Time	Event	Water	Acid	Base	Octanol	рΗ	dpH/dt	pH R-squared	pH SD	dpH/dt time
4:59.5	Initial pH = 5.97									
8:04.2	Data point 1		0.04965 mL					0.35010	0.00036	10.0 s
8:50.5	Data point 2	1.50000 mL	0.04965 mL	0.01651 mL	0.03001 mL	2.218	-0.00672	0.12678	0.00093	10.0 s
9:26.1	Data point 3		0.04965 mL					0.38859	0.00034	10.0 s
10:01.6	Data point 4	1.50000 mL	0.04965 mL	0.03151 mL	0.03001 mL	2.606	-0.00797	0.19529	0.00089	10.0 s
10:37.2	Data point 5		0.04965 mL					0.44711	0.00014	10.0 s
11:12.6	Data point 6	1.50000 mL	0.04965 mL	0.03775 mL	0.03001 mL	3.003	0.00280	0.09562	0.00045	10.0 s
11:48.1	Data point 7		0.04965 mL					0.14586	0.00048	10.0 s
	Data point 8		0.04965 mL						0.00020	
13:24.7	Data point 9		0.04965 mL						0.00055	10.0 s
14:00.2	Data point 10		0.04965 mL						0.00050	
14:46.0			0.04965 mL					0.86864	0.00063	10.0 s
	Data point 12		0.04965 mL						0.00090	10.0 s
15:56.8	Data point 13	1.50000 mL	0.04965 mL	0.04659 mL	0.03001 mL	4.268	-0.01706	0.91582	0.00088	10.5 s
16:32.8	Data point 14		0.04965 mL						0.00095	11.5 s
17:14.9			0.04965 mL						0.00089	13.0 s
	Data point 16		0.04965 mL						0.00092	13.0 s
18:31.7	Data point 17		0.04965 mL						0.00091	
19:11.1	Data point 18		0.04965 mL						0.00095	
20:02.3			0.04965 mL						0.00098	
20:54.3	Data point 20		0.04965 mL						0.00093	
21:55.7	Data point 21		0.04965 mL						0.00099	
23:19.4	Data point 22	1.50000 mL	0.04965 mL	0.04981 mL	0.03001 mL	6.965	-0.05661	0.98918	0.00281	Timed out at
										59.5 s
24:55.0	Data point 23	1.50000 mL	0.04965 mL	0.04988 mL	0.03001 mL	7.732	-0.10014	0.99537	0.00496	Timed out at
										59.5 s
26:25.5	Data point 24	1.50000 mL	0.04965 mL	0.04993 mL	0.03001 mL	8.300	-0.05511	0.99507	0.00273	Timed out at
										59.5 s
28:01.1	Data point 25		0.04965 mL						0.00099	
	Data point 26		0.04965 mL						0.00098	43.0 s
	Data point 27		0.10503 mL						0.00031	
	Data point 28		0.10503 mL						0.00013	10.0 s
	Data point 29		0.10503 mL						0.00018	
33:13.7			0.10503 mL					0.01014	0.00034	
	Data point 31		0.10503 mL						0.00030	
34:35.2			0.10503 mL						0.00017	
35:11.1	Data point 33		0.10503 mL						0.00055	
	Data point 34		0.10503 mL						0.00041	
36:53.1	Data point 35		0.10503 mL						0.00030	
37:28.6	Data point 36	1.50000 mL	0.10503 mL	0.10080 mL	0.11002 mL	3.806	-0.00526	0.75008	0.00030	10.0 s

1.50000 mL 0.10503 mL 0.10162 mL 0.11002 mL 3.999 -0.01309 0.90328

38:09.2 Data point 37

0.00068 10.0 s



18C-24003 Instrument ID: T312060

 $C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric\ high\ logP.t3r$ Filename:

Events (continued)

	(oontinaoa)									
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
38:49.7	Data point 38	1.50000 mL	0.10503 mL	0.10228 mL	0.11002 mL	4.196	-0.00633	0.53369	0.00043	10.0 s
39:35.4	Data point 39	1.50000 mL	0.10503 mL	0.10285 mL	0.11002 mL	4.425	-0.01381	0.73265	0.00080	10.0 s
40:10.8	Data point 40		0.10503 mL					0.00129	0.00069	
40:57.0	Data point 41		0.10503 mL					0.77724	0.00088	
41:38.5	Data point 42		0.10503 mL					0.38707	0.00067	
42:25.7	Data point 43		0.10503 mL						0.00089	
43:08.2	Data point 44		0.10503 mL					0.85957	0.00094	
43:52.2	Data point 45		0.10503 mL						0.00094	
45:04.4	Data point 46		0.10503 mL							Timed out
	Data point 40									at 59.5 s
46:29.9	Data point 47	1.50000 mL	0.10503 mL	0.10398 mL	0.11002 mL	7.265	-0.09351	0.99809	0.00462	Timed out at 59.5 s
48:00.3	Data point 48	1.50000 mL	0.10503 mL	0.10402 mL	0.11002 mL	7.780	-0.09755	0.99711	0.00483	Timed out at 59.5 s
49:35.9	Data point 49	1.50000 mL	0.10503 mL	0.10409 mL	0.11002 mL	8.206	-0.05783	0.98965	0.00287	Timed out
51:06.4	Data point 50	1.50000 mL	0.10503 mL	0.10414 mL	0.11002 mL	8.404	-0.03670	0.98260	0.00183	at 59.5 s Timed out
										at 59.5 s
52:36.9	Data point 51	1.50000 mL	0.10503 mL	0.10419 mL	0.11002 mL	8.772	-0.02623	0.91998	0.00135	Timed out at 59.5 s
54:07.4	Data point 52	1.50000 mL	0.10503 mL	0.10423 mL	0.11002 mL	9.008	-0.01967	0.95537	0.00099	
56:07.6	Data point 53		0.16334 mL					0.69174	0.00030	
56:53.8	Data point 54		0.16334 mL					0.42530	0.00086	
57:29.6	Data point 55		0.16334 mL					0.37291	0.00089	
58:05.2	Data point 56		0.16334 mL					0.57469	0.00065	
58:40.8	Data point 57		0.16334 mL					0.84718	0.00039	
59:16.3	Data point 57		0.16334 mL					0.17701	0.00053	
1:00:02.6			0.16334 mL					0.17701	0.00033	
1:00:02.0			0.16334 mL					0.12433	0.00086	
	Data point 60									
1:01:20.1	Data point 61		0.16334 mL					0.90591	0.00043	
1:02:11.1	Data point 62		0.16334 mL					0.62944	0.00063	
1:02:46.4	•		0.16334 mL					0.72392	0.00082	
	Data point 64		0.16334 mL						0.00070	
	Data point 65		0.16334 mL						0.00096	
1:05:11.2			0.16334 mL					0.79185	0.00087	
	Data point 67		0.16334 mL					0.62099	0.00099	
	Data point 68		0.16334 mL						0.00091	
1:07:35.2	Data point 69		0.16334 mL						0.00097	
1:08:50.9	Data point 70	1.50000 mL	0.16334 mL	0.16287 mL	0.51002 mL	6.378	-0.03153	0.95619	0.00159	Timed out at 59.5 s
1:10:21.3	Data point 71	1.50000 mL	0.16334 mL	0.16291 mL	0.51002 mL	6.803	-0.06207	0.98528	0.00309	Timed out at 59.5 s
1:11:51.8	Data point 72	1.50000 mL	0.16334 mL	0.16296 mL	0.51002 mL	7.279	-0.09237	0.99346	0.00458	Timed out at 59.5 s
1:13:22.3	Data point 73	1.50000 mL	0.16334 mL	0.16301 mL	0.51002 mL	7.684	-0.09598	0.99125	0.00476	Timed out at 59.5 s
1:14:52.8	Data point 74	1.50000 mL	0.16334 mL	0.16305 mL	0.51002 mL	8.114	-0.08212	0.99372	0.00407	Timed out at 59.5 s
1:16:23.3	Data point 75	1.50000 mL	0.16334 mL	0.16310 mL	0.51002 mL	8.426	-0.04721	0.97560	0.00236	Timed out at 59.5 s
1·17·50 0	Data point 76	1 50000 ml	0.16334 mL	0 16317 ml	0.51002 ml	8 764	_0 01910	U 03800	0.00092	
1:19:28.4	Data point 77 Assay volumes	1.50000 mL	0.16334 mL	0.16324 mL	0.51002 mL				0.00092	
1	, ,									



htap\20180323_exp3	3_logP_T3-2\180	C-24003_M04_octano	ol_pH-metric high logP		
Value	Original Value	Date/Time changed	Imported from		
Dorothy Levorse					
3					
2.000					
9.000					
0.200					
0.00002 mL					
0.10000 mL					
100%					
Cautious pH adjust					
None					
No					
No					
5 seconds					
10%					
None					
1.50 mL					
Automatic					
Octanol					
0.030 mL					
Automatic					
1 seconds					
Yes					
No					
60 seconds					
5 seconds					
Yes					
n To start pH					
	Value Dorothy Levorse 3 2.000 9.000 0.200 0.00002 mL 0.10000 mL 100% Cautious pH adjust None No No 5 seconds 10% None 1.50 mL Automatic Octanol 0.030 mL Automatic 1 seconds Yes No 60 seconds 5 seconds	Value Original Value Dorothy Levorse 3 2.000 9.000 0.200 0.00002 mL 0.10000 mL 100% Cautious pH adjust None No No S seconds 10% None 1.50 mL Automatic Octanol 0.030 mL Automatic 1 seconds Yes No 60 seconds 5 seconds 7 yes	Value Original Value Date/Time changed Dorothy Levorse 3 2.000 9.000 0.200 0.00002 mL 0.10000 mL 100% Cautious pH adjust None No No S seconds 10% None 1.50 mL Automatic Octanol 0.030 mL Automatic 1 seconds Yes No 60 seconds 5 seconds 5 seconds 7 yes		

Sample Dissolution	
Perform a dissolution stage	Yes
Adjust and hold pH for dissolution	To start pH
Stir to dissolve for	120 seconds
For dissolution, stir at	10%

Carbonate purge

Perform	a carbonate	purae	No
CITCITI	a carbonate	pargo	110

Temperature Control

Wait for temperature	Yes
Required start temperature	25.0°C
Acceptable deviation	0.5°C
Time to wait	60 seconds
Stir speed of	50%
Titration 1	

110 40011 1	
Titrate from	Low to high pH
Adjust to start pH	Yes
After pH adjust stir for	30 seconds
Stir to allow partitioning for	15 seconds
Stirrer speed for partitioning	50%
Titration 2	
Titrate from	Low to high pH
Add additional water	0.00 ml

Additional partition solvent volume 0.080 mL Additional partition solvent added Automatic After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds Stirrer speed for partitioning 55%



Assay ID: 18C-24003 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3		•	•	•
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.400 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation				
Stability timeout after	60 seconds			

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.119	3/24/2018 2:50:31 AM	C:\Sirius_T3\HCl18C23.t3r
Four-Plus S	0.9972	3/24/2018 2:50:31 AM	C:\Sirius_T3\HCl18C23.t3r
Four-Plus jH	0.9	3/24/2018 2:50:31 AM	C:\Sirius_T3\HCl18C23.t3r
Four-Plus jOH	-0.3	3/24/2018 2:50:31 AM	C:\Sirius_T3\HCl18C23.t3r
Base concentration factor	1.003	3/24/2018 2:50:31 AM	C:\Sirius_T3\KOH18C23.t3r
Acid concentration factor	0.997	3/24/2018 2:50:31 AM	C:\Sirius_T3\HCl18C23.t3r

Instrument Settings

Setting Instrument owner Instrument ID Instrument type Software version	Value Merck T312060 T3 Simulator 1.1.3.0	Batch Id	Install date
Dispenser module Dispenser 0 Syringe volume Firmware version	Water 2.5 mL 1.2.1(r2)	T3DM1200361	3/31/2009 6:24:52 AM 3/31/2009 6:25:05 AM
Titrant Dispenser 2 Syringe volume Firmware version	Water (0.15 M KCI) Acid 0.5 mL 1.2.1(r2)	02-06-2018	3/16/2018 11:09:18 AM 3/31/2009 6:25:11 AM
Titrant Dispenser 1 Syringe volume Firmware version	Acid (0.5 M HCI) Base 0.5 mL 1.2.1(r2)	03-16-2018	3/16/2018 10:56:23 AM 3/31/2009 6:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2)	3/22/2018	3/23/2018 9:34:17 AM 3/31/2009 6:26:24 AM
Distribution valve 5 Firmware version Port A	Distribution Valve 1.1.3 Methanol (80%, 0.15 M KCI)	02-08-2018	3/31/2009 6:28:19 AM 3/6/2018 10:28:59 AM
Port B Dispenser 3 Syringe volume Firmware version	Cyclohexane Buffer 0.5 mL 1.2.1(r2)	11-01-17	2/27/2018 11:37:57 AM 8/3/2010 6:05:16 AM
Titrant Dispenser 6	Dodecane Octanol	2018/01/31	2/28/2018 11:18:04 AM 10/22/2010 11:52:43 AM




Assay ID: 18C-24003 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting Syringe volume	Value 0.5 mL	Batch Id	Install date
Firmware version Titrant	1.2.1(r2) Octanol	01-31-2018	2/27/2018 10:59:35 AM
Titrator Horizontal axis firmware version Vertical axis firmware version Chassis I/O firmware version	1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2 1.11 Al1Dl0DO4 Norgren I/O	T3TM1200161	3/31/2009 6:24:17 AM
Probe I/O firmware version Electrode E0 calibration Filling solution	1.1.1 T3 Electrode +5.05 mV 3M KCI	T3E0923 KCL097	1/23/2018 3:01:00 PM 3/24/2018 2:50:59 AM 3/23/2018 9:29:07 AM
Liquids Wash 1 Wash 2 Buffer position 1 Buffer position 2 Storage position Wash water Waste Temperature controller Turbidity detector	50% IPA:50% Water 0.5% Trition X-100 in H20 pH7 Wash pH 7 7.6e+003 mL 2.6e+003 mL	03-12-2018	3/23/2018 9:29:12 AM 3/23/2018 9:29:15 AM 3/23/2018 9:29:19 AM 3/23/2018 9:29:21 AM 3/23/2018 9:30:23 AM 3/12/2018 9:25:04 AM 3/12/2018 9:24:49 AM 8/5/2010 7:35:13 AM 3/31/2009 6:24:45 AM
Spectrometer Dip probe Wavelength coefficient A0	183.333	074811 10196	11/23/2010 12:22:28 PM
Wavelength coefficient A1 Wavelength coefficient A2 Total lamp lit time Calibrated on Integration time	2.21568 -0.000289308 162:53:01 2/27/2018 11:40:38 AM 40		11/23/2010 12:22:28 PM
Scans averaged Autoloader Left-right axis firmware version Front-back axis firmware version Vertical axis firmware version Chassis I/O firmware version	1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2 1.17 Al1Dl2DO2 Stepper 2 1.11 Al1Dl2DO4 Normen I/O	T3AL1200345	11/10/2015 10:34:13 AM
Chassis I/O firmware version Configuration Alternate titration position Alternate reference position Maximum standard vial volume Maximum alternate vial volume Automatic action idle period Titrant tube volume Syringe flush count Flowing wash pump volume Flowing wash stir duration Flowing wash stir speed Solvent wash stir duration Solvent wash stir speed Surfactant wash stir speed Surfactant wash stir speed E0 calibration minimum number of points E0 calibration maximum standard deviation E0 calibration timeout period E0 calibration stir duration E0 calibration preparation stir speed E0 calibration buffer wash stir duration E0 calibration buffer wash stir speed E0 calibration reading stir speed	1.11 Al1DI0DO4 Norgren I/O Titration position Reference position 3.50 mL 25.00 mL 5 minute(s) 1.3 mL 3.50 20.0 mL 5 s 30% 5 s 30% 5 s 30% 10 0.01500 60 s 5 s 30% 5 s 30% 5 s 30% 5 s 30% 60 s		





Assay ID: 18C-24003 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180323_exp33_logP_T3-2\18C-24003_M04_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting	value	Batch Id	install date
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		

Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00