

18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

## pH-metric Result

logP (XH +) 0.45 ±0.04 (n=50) logP (neutral X) 4.08 ±0.01 (n=50)

#### 18C-03011 Points 1 to 20

M02 octanol concentration factor 0.967 Carbonate 0.1632 mM Acidity error -1.49856 mM

#### 18C-03011 Points 21 to 48

M02\_octanol concentration factor 1.157 Carbonate 0.1227 mM -1.56956 mM Acidity error

#### 18C-03011 Points 49 to 70

M02 octanol concentration factor 1.597 Carbonate 0.1477 mM Acidity error -1.25180 mM

## Warnings and errors

Errors None Warnings None

рΗ

# Sample logD and percent species

M02\_octanol M02\_octanol M02\_octanol

•	_	_	_	<b>—</b>	_	
	logD	M02_octanolH	M02_octanol	M02_octanolH*	M02_octanol*	
1.000	0.59	20.31 %	0.00 %	56.73 %	22.95 %	
1.200	0.66	17.91 %	0.00 %	50.02 %	32.07 %	Stomach pH
2.000	1.15	6.63 %	0.01 %	18.50 %	74.87 %	
3.000	2.06	0.86 %	0.01 %	2.39 %	96.74 %	
4.000	3.02	0.09 %	0.01 %	0.25 %	99.66 %	
5.000	3.77	0.01 %	0.01 %	0.02 %	99.96 %	
6.000	4.04	0.00 %	0.01 %	0.00 %	99.99 %	
6.500	4.07	0.00 %	0.01 %	0.00 %	99.99 %	
7.000	4.08	0.00 %	0.01 %	0.00 %	99.99 %	
7.400	4.08	0.00 %	0.01 %	0.00 %	99.99 %	Blood pH
8.000	4.08	0.00 %	0.01 %	0.00 %	99.99 %	
9.000	4.08	0.00 %	0.01 %	0.00 %	99.99 %	
10.000	4.08	0.00 %	0.01 %	0.00 %	99.99 %	
11.000	4.08	0.00 %	0.01 %	0.00 %	99.99 %	
12.000	4.08	0.00 %	0.01 %	0.00 %	99.99 %	

M02\_octanol

M02\_octanol Comment



Sample name: M02\_octanol Assay name: Assay ID:

Filename:

pH-metric high logP


18C-03011

Experiment start time: 3/3/2018 1:46:01 PM

Analyst: **Dorothy Levorse** 

Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r






Filename:

Sample name: M02\_octanol Experiment start time: 3/3/2018 1:46:01 PM pH-metric high logP Analyst: **Dorothy Levorse** Assay name: Assay ID:

18C-03011 Instrument ID: T312060

 $C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric\ high\ logP.t3r$ 

# **Graphs** (continued)





Sample name: M02\_octanol Experiment start time: 3/3/2018 1:46:01 PM Analyst: Assay name: pH-metric high logP **Dorothy Levorse** Assay ID:

Instrument ID: T312060 18C-03011

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 1 of 3 18C-03011 Points 1 to 20

#### Overall results

RMSD 0.135 Average ionic strength 0.157 M Average temperature 24.9°C Partition ratio 0.0123:1

Analyte concentration range 4180.2 μM to 4304.6 μM

Total points considered 16 of 20

## Warnings and errors

Errors None

Warnings Excessive acidity error present

#### Four-Plus parameters

Alpha 0.111 3/3/2018 1:46:01 PM C:\Sirius\_T3\HCl18C02.t3r 0.9988 3/3/2018 1:46:01 PM C:\Sirius\_T3\HCl18C02.t3r S 3/3/2018 1:46:01 PM C:\Sirius\_T3\HCl18C02.t3r jΗ 1.0 jOH -0.83/3/2018 1:46:01 PM C:\Sirius\_T3\HCl18C02.t3r


# Titrants


0.50 M HCI 0.999058 3/3/2018 1:46:01 PM C:\Sirius T3\HCI18C02.t3r 0.50 M KOH 0.999845 3/3/2018 1:46:01 PM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M02 octanol concentration factor 0.967 Base pKa 1 5.03 logP(XH +)1.31 logP (neutral X) 4.13

#### Sample graphs





Analyst:

Experiment start time: 3/3/2018 1:46:01 PM

**Dorothy Levorse** 



Assay ID:

Filename:


Sample name: M02\_octanol Assay name: pH-metric high logP

18C-03011

C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

Instrument ID: T312060


# Sample graphs (continued)




## Sample logD and percent species


рН	M02_octanol logD	M02_octanol M02_octanolH	M02_octanol M02_octanol	M02_octanol M02_octanolH*	M02_octanol M02_octanol*	Comment
1.000	1.33	79.02 %	0.01 %		1.22 %	
1.200	1.35	78.46 %	0.01 %	19.60 %	1.93 %	Stomach pH
2.000	1.52	71.14 %	0.07 %	17.77 %	11.02 %	•
3.000	2.16	35.60 %	0.33 %	8.90 %	55.17 %	
4.000	3.07	5.94 %	0.55 %	1.48 %	92.02 %	
5.000	3.82	0.64 %	0.59 %	0.16 %	98.61 %	
6.000	4.09	0.06 %	0.60 %	0.02 %	99.32 %	
6.500	4.12	0.02 %	0.60 %	0.01 %	99.38 %	
7.000	4.13	0.01 %	0.60 %	0.00 %	99.39 %	
7.400	4.13	0.00 %	0.60 %	0.00 %	99.40 %	Blood pH
8.000	4.13	0.00 %	0.60 %	0.00 %	99.40 %	•
9.000	4.13	0.00 %	0.60 %	0.00 %	99.40 %	
10.000	4.13	0.00 %	0.60 %	0.00 %	99.40 %	
11.000	4.13	0.00 %	0.60 %	0.00 %	99.40 %	
12.000	4.13	0.00 %	0.60 %	0.00 %	99.40 %	

# Carbonate and acidity



# Other graphs





Analyst:

Experiment start time: 3/3/2018 1:46:01 PM

**Dorothy Levorse** 

T312060



Assay ID:

Filename:


Sample name: M02\_octanol
Assay name: pH-metric hig

pH-metric high logP

18C-03011 Instrument ID:

C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

# Other graphs (continued)





Instrument ID: T312060 18C-03011

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 2 of 3 18C-03011 Points 21 to 48

## Overall results

RMSD 0.183 Average ionic strength 0.163 M Average temperature 25.0°C Partition ratio 0.0290:1

Analyte concentration range 3843.8 µM to 3969.0 µM

Total points considered 20 of 28

## Warnings and errors

None Errors

Warnings Excessive acidity error present

## Four-Plus parameters

Alpha	0.111	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r
S	0.9988	3/3/2018 1:46:01 PM	C:\Sirius T3\HCl18C02.t3r
jΗ	1.0	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r
jОН	-0.8	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r


# Titrants


0.50 M HCI 0.999058 3/3/2018 1:46:01 PM C:\Sirius T3\HCI18C02.t3r 0.50 M KOH 0.999845 3/3/2018 1:46:01 PM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M02\_octanol concentration factor 1.157 Base pKa 1 5.03 logP (XH +) 1.35 logP (neutral X) 4.30

## Sample graphs





Analyst:

Experiment start time: 3/3/2018 1:46:01 PM

**Dorothy Levorse** 

T312060




Assay ID:

Sample name: M02\_octanol Assay name: pH-metric high logP

Instrument ID: 18C-03011

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r


# Sample graphs (continued)



## Sample logD and percent species


рН	M02_octanol logD	M02_octanol M02_octanolH	M02_octanol M02_octanol	M02_octanol M02_octanolH*	M02_octanol M02_octanol*	Comment
1.000	1.39	58.55 %	0.01 %	38.26 %	3.18 %	
1.200	1.41	57.48 %	0.01 %	37.56 %	4.95 %	Stomach pH
2.000	1.61	45.50 %	0.04 %	29.73 %	24.72 %	
3.000	2.32	14.10 %	0.13 %	9.21 %	76.56 %	
4.000	3.24	1.78 %	0.17 %	1.17 %	96.88 %	
5.000	3.99	0.18 %	0.17 %	0.12 %	99.53 %	
6.000	4.26	0.02 %	0.17 %	0.01 %	99.80 %	
6.500	4.29	0.01 %	0.17 %	0.00 %	99.82 %	
7.000	4.30	0.00 %	0.17 %	0.00 %	99.83 %	
7.400	4.30	0.00 %	0.17 %	0.00 %	99.83 %	Blood pH
8.000	4.30	0.00 %	0.17 %	0.00 %	99.83 %	·
9.000	4.30	0.00 %	0.17 %	0.00 %	99.83 %	
10.000	4.30	0.00 %	0.17 %	0.00 %	99.83 %	
11.000	4.30	0.00 %	0.17 %	0.00 %	99.83 %	
12.000	4.30	0.00 %	0.17 %	0.00 %	99.83 %	


# Carbonate and acidity



0.123 mM

# Other graphs







Assay ID: Filename:

Sample name: M02\_octanol Assay name:

pH-metric high logP


18C-03011

C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

Experiment start time: 3/3/2018 1:46:01 PM Analyst: **Dorothy Levorse** 

Instrument ID: T312060

# Other graphs (continued)





Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

# pH-metric high logP Titration 3 of 3 18C-03011 Points 49 to 70

#### Overall results

RMSD 0.256
Average ionic strength 0.169 M
Average temperature 25.0°C
Partition ratio 0.1643 : 1

Analyte concentration range 3175.3 µM to 3268.6 µM

Total points considered 16 of 22

## Warnings and errors

Errors None

Warnings Sample concentration factor out of range

Excessive acidity error present

## Four-Plus parameters

Alpha	0.111	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r
S	0.9988	3/3/2018 1:46:01 PM	C:\Sirius T3\HCl18C02.t3r
jΗ	1.0	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r
jOH	-0.8	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r


#### **Titrants**


0.50 M HCI 0.999058 3/3/2018 1:46:01 PM C:\Sirius\_T3\HCl18C02.t3r 0.50 M KOH 0.999845 3/3/2018 1:46:01 PM C:\Sirius\_T3\KOH18B27.t3r

#### Sample

M02\_octanol concentration factor 1.597
Base pKa 1 5.03
logP (XH +) 1.32
logP (neutral X) 4.56

#### Sample graphs





Analyst:

Experiment start time: 3/3/2018 1:46:01 PM

**Dorothy Levorse** 

T312060




Assay ID:

Sample name: M02\_octanol
Assay name: pH-metric high logP

18C-03011 Instrument ID:

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r


# Sample graphs (continued)




## Sample logD and percent species

рН	M02_octanol	M02_octanol	M02_octanol	M02_octanol	M02_octanol	Comment
	logD	M02_octanolH	M02_octanol	M02_octanolH*	M02_octanol*	
1.000	1.39	20.02 %	0.00 %	68.72 %	11.26 %	
1.200	1.42	18.78 %	0.00 %	64.47 %	16.75 %	Stomach pH
2.000	1.74	9.94 %	0.01 %	34.13 %	55.92 %	
3.000	2.56	1.65 %	0.02 %	5.66 %	92.68 %	
4.000	3.50	0.18 %	0.02 %	0.61 %	99.20 %	
5.000	4.25	0.02 %	0.02 %	0.06 %	99.90 %	
6.000	4.52	0.00 %	0.02 %	0.01 %	99.98 %	
6.500	4.55	0.00 %	0.02 %	0.00 %	99.98 %	
7.000	4.56	0.00 %	0.02 %	0.00 %	99.98 %	
7.400	4.56	0.00 %	0.02 %	0.00 %	99.98 %	Blood pH
8.000	4.56	0.00 %	0.02 %	0.00 %	99.98 %	
9.000	4.56	0.00 %	0.02 %	0.00 %	99.98 %	
10.000	4.56	0.00 %	0.02 %	0.00 %	99.98 %	
11.000	4.56	0.00 %	0.02 %	0.00 %	99.98 %	
12.000	4.56	0.00 %	0.02 %	0.00 %	99.98 %	

# Carbonate and acidity



# Other graphs





Assay ID: Filename:

Sample name: M02\_octanol Assay name:


pH-metric high logP


18C-03011

Experiment start time: 3/3/2018 1:46:01 PM Analyst: **Dorothy Levorse** 

Instrument ID: T312060

C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r







Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

## Assay Model

Settings	Value	Date/Time changed	Imported from
Sample name	M02_octanol	12/6/2017 4:20:03 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.001960 g	3/2/2018 5:09:57 PM	User entered value
Formula weight	289.26 g/mol	12/6/2017 4:20:03 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	289.26	12/6/2017 4:20:03 PM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	1	12/6/2017 4:20:03 PM	User entered value
Sample is a	Base	12/6/2017 4:20:03 PM	User entered value
pKa 1	5.03	12/6/2017 4:20:03 PM	User entered value
logp (XH +)	1.32	3/2/2018 3:38:13 PM	User entered value
logP (neutral X)	4.10	3/2/2018 3:38:07 PM	User entered value

#### **Events**

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
6:00.1 8:59.7	Initial pH = 8.70	1 50000 ml	0.05252 mL	0.00162 ml	0.01000	2.079	0.01626	0.75364	0.00093	15.00
6.59.7 9:50.7	Data point 1 Data point 2		0.05252 IIIL 0.05252 mL						0.00093	
10:26.8			0.05252 mL						0.00003	
11:02.9			0.05252 mL						0.00011	
11:38.4	Data point 5		0.05252 mL						0.00021	
12:13.9			0.05252 mL						0.00078	
12:49.9			0.05252 mL						0.00084	
13:25.4			0.05252 mL						0.00084	
	Data point 9		0.05252 mL						0.00098	
15:11.8			0.05252 mL						0.00096	
	Data point 10		0.05252 mL						0.00096	
	Data point 12		0.05252 mL						0.00096	
	Data point 12		0.05252 mL						0.00098	
19:56.1			0.05252 mL							Timed out at
	·									59.5 s
21:31.8	Data point 15	1.50000 mL	0.05252 mL	0.04800 mL	0.01999 mL	6.482	-0.04707	0.98476	0.00234	Timed out at 59.5 s
23:07.4	Data point 16	1.50000 mL	0.05252 mL	0.04809 mL	0.01999 mL	6.961	-0.05553	0.98035	0.00277	Timed out at
24:38.0	Data point 17	1.50000 mL	0.05252 mL	0.04817 mL	0.01999 mL	7.521	-0.08676	0.99361	0.00430	59.5 s Timed out at
	•									59.5 s
26:13.5	Data point 18	1.50000 mL	0.05252 mL	0.04824 mL	0.01999 mL	8.125	-0.07978	0.98996	0.00396	Timed out at 59.5 s
27:49.3	Data point 19	1.50000 mL	0.05252 mL	0.04840 mL	0.01999 mL	8.917	-0.01848	0.88522	0.00097	
29:15.8			0.05252 mL						0.00091	
30:51.7	•		0.10875 mL						0.00059	
	Data point 22		0.10875 mL						0.00045	10.0 s
32:13.7		1.50000 mL	0.10875 mL	0.07660 mL	0.05000 mL	2.396	-0.00146	0.31714	0.00013	
32:49.9			0.10875 mL						0.00026	10.0 s
33:25.5	Data point 25	1.50000 mL	0.10875 mL	0.08998 mL	0.05000 mL	2.785	-0.00991	0.29494	0.00090	10.0 s
34:01.0	Data point 26		0.10875 mL						0.00035	10.0 s
34:36.5	Data point 27	1.50000 mL	0.10875 mL	0.09722 mL	0.05000 mL	3.160	-0.00652	0.78690	0.00036	10.0 s
35:11.9	Data point 28	1.50000 mL	0.10875 mL	0.09951 mL	0.05000 mL	3.395	-0.01414	0.94922	0.00072	10.5 s
36:03.3	Data point 29	1.50000 mL	0.10875 mL	0.10073 mL	0.05000 mL	3.592	-0.01172	0.70688	0.00069	10.0 s
36:54.3	Data point 30	1.50000 mL	0.10875 mL	0.10158 mL	0.05000 mL	3.795	-0.01307	0.64357	0.00080	10.0 s
37:45.1	Data point 31	1.50000 mL	0.10875 mL	0.10219 mL	0.05000 mL	4.005	-0.01735	0.75875	0.00098	10.5 s
38:21.1		1.50000 mL	0.10875 mL	0.10261 mL	0.05000 mL	4.238	-0.01660	0.84161	0.00089	12.5 s
38:59.1	Data point 33	1.50000 mL	0.10875 mL	0.10287 mL	0.05000 mL	4.447	-0.01934	0.93729	0.00099	16.0 s
39:40.4	Data point 34	1.50000 mL	0.10875 mL	0.10303 mL	0.05000 mL	4.596	-0.01901	0.88676	0.00100	30.5 s
40.44 =	D-1 1-1 0F	4 50000	0.40075	0.40000	0.05000	4 0 40	0.04004	0.77007	0.0000	400

1.50000 mL 0.10875 mL 0.10322 mL 0.05000 mL 4.948 -0.01691 0.77907

Reported at: 3/6/2018 2:10:02 PM

40:41.5 Data point 35

0.00095 16.0 s



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

# **Events (continued)**

r vents (	(continueu)									
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/dt time
41:28.0	Data point 36	1.50000 mL	0.10875 mL	0.10334 mL	0.05000 mL	5.244	-0.01809	0.82366	0.00098	32.5 s
42:31.1	Data point 37	1.50000 mL	0.10875 mL	0.10341 mL	0.05000 mL	5.578	-0.01894	0.89362	0.00099	47.0 s
43:48.7	Data point 38	1.50000 mL	0.10875 mL	0.10348 mL	0.05000 mL	6.124	-0.02064	0.93374	0.00106	Timed out
										at 59.5 s
45:19.2	Data point 39	1.50000 mL	0.10875 mL	0.10358 mL	0.05000 mL	6.796	-0.05867	0.98597	0.00292	Timed out
										at 59.5 s
46:54.9	Data point 40	1.50000 mL	0.10875 mL	0.10365 mL	0.05000 mL	7.247	-0.06483	0.98904	0.00322	Timed out
	5	4 =0000						0.00400		at 59.5 s
48:30.7	Data point 41	1.50000 mL	0.10875 mL	0.103/2 mL	0.05000 mL	7.610	-0.08577	0.99139	0.00426	Timed out
50.04.0	Data ==:=1.40	4 500001	0.40075!	0.400701	0.05000!	7 7 7 4	0.00000	0.00400	0.00044	at 59.5 s
50:01.2	Data point 42	1.50000 ML	0.10875 ML	0.10376 ML	0.05000 mL	7.754	-0.06230	0.98106	0.00311	Timed out
51:31.7	Data point 42	1 50000 ml	0 10075 ml	0 10201 ml	0.05000 mL	0 020	0.05002	0.00074	0 00202	at 59.5 s Timed out
31.31.7	Data point 43	1.50000 IIIL	0.10075 IIIL	0.10361 IIIL	0.03000 IIIL	0.030	-0.05665	0.90974	0.00292	at 59.5 s
53:02.2	Data point 44	1 50000 ml	0 10875 ml	0 10386 ml	0.05000 mL	8 260	-0.04245	0.97500	0.00212	Timed out
33.02.2	Data point 44	1.50000 IIIL	0.10075 IIIL	0.10300 IIIL	0.03000 IIIL	0.209	-0.04243	0.97500	0.00212	at 59.5 s
54:32.7	Data point 45	1 50000 ml	0 10875 ml	0.10390 ml	0.05000 mL	8 461	-0.02925	0 94322	0 00149	Timed out
J4.52.7	Data point 40	1.50000 1112	0.10075 IIIL	0.10000 IIIL	0.03000 IIIL	0.401	-0.02020	0.04022	0.00143	at 59.5 s
56:08.3	Data point 46	1 50000 ml	0 10875 ml	0 10398 ml	0.05000 mL	8 742	-0 01524	0 61499	0.00096	
57:14.5	Data point 47				0.05000 mL			0.86242	0.00099	
58:21.1	Data point 48				0.05000 mL			0.74240	0.00096	
59:51.2	Data point 49				0.30000 mL			0.07579	0.00065	
	Data point 50				0.30000 mL			0.52474	0.00091	
	Data point 51				0.30000 mL			0.50607	0.00087	10.0 s
	Data point 52	1.50000 mL	0.16893 mL	0.14699 mL	0.30000 mL	2.577	-0.00872	0.34259	0.00074	10.0 s
1:02:24.5	Data point 53	1.50000 mL	0.16893 mL	0.15303 mL	0.30000 mL	2.768	-0.00055	0.00118	0.00079	10.0 s
1:02:59.9	Data point 54				0.30000 mL			0.04207	0.00096	14.5 s
	Data point 55				0.30000 mL			0.03061	0.00083	10.0 s
	Data point 56				0.30000 mL			0.28200	0.00096	
1:05:06.3					0.30000 mL			0.32693	0.00090	
	Data point 58				0.30000 mL			0.52606	0.00032	
	Data point 59				0.30000 mL			0.91405	0.00064	
1:07:09.1					0.30000 mL				0.00097	
	Data point 61				0.30000 mL			0.55571	0.00087	
	Data point 62				0.30000 mL				0.00094	
1:10:35.1	Data point 63	1.50000 ML	0.16893 IIIL	0.16446 IIIL	0.30000 mL	0.140	-0.03968	0.91076	0.00205	Timed out
1.12.10 7	Data point 64	1 50000 ml	0 16903 ml	0.16459 ml	0.30000 mL	6 094	0.00034	0.08030	0.00405	at 59.5 s Timed out
1.12.10.7	Data point 04	1.50000 IIIL	0.10093 IIIL	0.10430 IIIL	0.50000 IIIL	0.304	-0.03324	0.90030	0.00433	at 59.5 s
1.13.113	Data point 65	1 50000 ml	0 16803 ml	0 16/63 ml	0.30000 mL	7 173	_0 07212	0 07043	0 00363	Timed out
1.10.41.0	Data point 00	1.50000 IIIL	0.10095 IIIL	0.10 <del>4</del> 05 IIIL	0.50000 IIIL	7.175	-0.07212	0.97043	0.00302	at 59.5 s
1.15.11 8	Data point 66	1 50000 ml	0 16893 ml	0 16467 ml	0.30000 mL	7 339	-0.06600	0 98427	0.00328	Timed out
1.10.11.0	Bata point oo	1.00000 1112	0.10000 IIIL	0.10407 IIIL	0.00000 IIIL	7.000	0.00000	0.00421	0.00020	at 59.5 s
1.16.42.3	Data point 67	1 50000 ml	0 16893 ml	0 16472 ml	0.30000 mL	7 729	-0.07810	0 99269	0.00387	Timed out
	Data point of		0.10000 1112	0022	0.000002	0	0.07010	0.00200	0.00001	at 59.5 s
1:18:12.8	Data point 68	1.50000 mL	0.16893 mL	0.16477 mL	0.30000 mL	8.153	-0.05661	0.98395	0.00282	Timed out
	5.15. p 5.1.1.					21700			<b></b>	at 59.5 s
1:19:48.4	Data point 69	1.50000 mL	0.16893 mL	0.16484 mL	0.30000 mL	8.561	-0.03628	0.93717	0.00185	Timed out
	,	. –	· —	_	· <del>-</del>	-				at 59.5 s
1:21:24.0	Data point 70	1.50000 mL	0.16893 mL	0.16505 mL	0.30000 mL	9.190	-0.01034	0.34587	0.00087	15.0 s
1:21:48.1	Assay volumes	1.50000 mL	0.16893 mL	0.16505 mL	0.30000 mL					



Assay ID: Instrument ID: 18C-03011 T312060

Filename: C:\Sirius_T3\Meh	tap\20180302_exp2	9_logP_T3-2\180	C-03011_M02_octano	I_pH-metric high logP.t
Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Dorothy Levorse			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	9.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
Advanced General Settings				
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium				
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.020 mL			
Partition solvent added	Automatic			
After partition addition, stir for	1 seconds			
Sample Sonication				
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	120 seconds			
After sonication stir for	5 seconds			
Sample Dissolution				
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution	To start pH			
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge				
Perform a carbonate purge	No			
Temperature Control				
Wait for temperature	Yes			
Required start temperature	25.0°C			
A scentable deviation	0 F°C			

Acceptable deviation 0.5°C Time to wait 60 seconds

Stir speed of Titration 1

Titrate from Low to high pH

50%

Adjust to start pH Yes

After pH adjust stir for 30 seconds Stir to allow partitioning for 15 seconds Stirrer speed for partitioning 50%

Titration 2

Titrate from Low to high pH 0.00 mL

Add additional water Additional partition solvent volume 0.030 mL Additional partition solvent added Automatic After pH adjust stir for 30 seconds

Stir to allow partitioning for 15 seconds Stirrer speed for partitioning 55%



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

# Assay Settings (continued)

Value	Original Value	Date/Time changed	Imported from
	_	_	
Low to high pH			
0.00 mL			
0.250 mL			
Automatic			
30 seconds			
15 seconds			
60%			
No			
0 seconds			
20 points			
0.50 seconds			
0.00100 dpH/dt			
60 seconds			
	Low to high pH 0.00 mL 0.250 mL Automatic 30 seconds 15 seconds 60% No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt	Low to high pH 0.00 mL 0.250 mL Automatic 30 seconds 15 seconds 60%  No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt	Low to high pH 0.00 mL 0.250 mL Automatic 30 seconds 15 seconds 60%  No 0 seconds 20 points 0.50 seconds 0.00100 dpH/dt

# Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.111	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r
Four-Plus S	0.9988	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r
Four-Plus jH	1.0	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r
Four-Plus jOH	-0.8	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r
Base concentration factor	1.000	3/3/2018 1:46:01 PM	C:\Sirius_T3\KOH18B27.t3r
Acid concentration factor	0.999	3/3/2018 1:46:01 PM	C:\Sirius_T3\HCl18C02.t3r

# Instrument Settings

Setting	Value	Batch Id	Install date
Instrument owner	Merck		
Instrument ID	T312060		
Instrument type	T3 Simulator		
Software version	1.1.3.0		
Dispenser module		T3DM1200361	3/31/2009 5:24:52 AM
Dispenser 0	Water		3/31/2009 5:25:05 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Water (0.15 M KCI)	02-06-2018	2/27/2018 10:05:59 AM
Dispenser 2	Acid		3/31/2009 5:25:11 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Acid (0.5 M HCI)	02-27-2018	2/27/2018 10:27:22 AM
Dispenser 1	Base		3/31/2009 5:25:21 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Base (0.5 M KOH)	9/22/2017	2/27/2018 10:21:22 AM
Dispenser 5	Cosolvent		3/31/2009 5:26:24 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Distribution valve 5	Distribution Valve		3/31/2009 5:28:19 AM
Firmware version	1.1.3		
Port A	Methanol (80%, 0.15 M KCl)	09-26-17	2/7/2018 9:42:01 AM
Port B	Cyclohexane	11-01-17	2/27/2018 10:37:57 AM
Dispenser 3	Buffer		8/3/2010 5:05:16 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Dodecane	2018/01/31	2/28/2018 10:18:04 AM
Dispenser 6	Octanol		10/22/2010 10:52:43 AM



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

# Instrument Settings (continued)

Setting Syringe volume	Value 0.5 mL	Batch Id	Install date
Firmware version	1.2.1(r2)		
Titrant	Octanol	01-31-2018	2/27/2018 9:59:35 AM
Titrator		T3TM1200161	3/31/2009 5:24:17 AM
Horizontal axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1	T05000	4/00/0040 0 04 00 DM
Electrode	T3 Electrode	T3E0923	1/23/2018 2:01:00 PM
E0 calibration	+6.59 mV	KCI 007	3/3/2018 1:46:29 PM
Filling solution Liquids	3M KCI	KCL097	3/2/2018 9:43:24 AM
Wash 1	50% IPA:50% Water		3/2/2018 9:45:12 AM
Wash 2	0.5% Trition X-100 in H20		3/2/2018 9:45:15 AM
Buffer position 1	pH7 Wash		3/2/2018 9:45:18 AM
Buffer position 2	pH 7		3/2/2018 9:45:21 AM
Storage position	•		3/2/2018 9:44:44 AM
Wash water	6.7e+003 mL	02-27-2018	2/27/2018 9:54:39 AM
Waste	8.7e+003 mL		11/28/2017 10:36:29 AM
Temperature controller			8/5/2010 6:35:13 AM
Turbidity detector			3/31/2009 5:24:45 AM
Spectrometer		074811	11/23/2010 11:22:28 AM
Dip probe	100.000	10196	
Wavelength coefficient A0	183.333		
Wavelength coefficient A1 Wavelength coefficient A2	2.21568 -0.000289308		
Total lamp lit time	120:41:49		11/23/2010 11:22:28 AM
Calibrated on	2/27/2018 10:40:38 AM		11/23/2010 11.22.20 AW
Integration time	40		
Scans averaged	10		
Autoloader		T3AL1200345	11/10/2015 9:34:13 AM
Left-right axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Configuration	Tituation position		
Alternate titration position	Titration position		
Alternate reference position Maximum standard vial volume	Reference position 3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points E0 calibration maximum standard deviation	10 0.01500		
E0 calibration timeout period	60 s		
E0 calibration timeout period E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

## Instrument Settings (continued)

Setting	Value	Batch Id	l Install date
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		

### Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00
, ,		

### Experiment Log

[2:38] Air gap created for Water	(0.15 M KCI)
[2:38] Air gap created for Acid (	

- [2:38] Air gap created for Base (0.5 M KOH)
- [2:39] Air gap released for Water (0.15 M KCI)
- [2:42] Titrator arm moved over Titration position
- [2:42] Titration 1 of 3
- [2:42] Adding initial titrants
- [2:42] Automatically add 1.50000 mL of water
- [3:08] Dispensed 1.500000 mL of Water (0.15 M KCI)
- [3:12] Titrator arm moved over Drain
- [5:54] Titrator arm moved to Titration position
- [5:54] Argon flow rate set to 100 [5:54] Stirrer speed set to 10
- [5:59] Automatically add 0.02000 mL of Octanol
- [5:59] Dispensed 0.019991 mL of Octanol
- [6:00] Initial pH = 8.70
- [6:00] Iterative adjust 8.70 -> 2.00
- [6:00] pH 8.70 -> 2.00
- [6:02] Air gap released for Acid (0.5 M HCI)
- [6:03] Dispensed 0.052516 mL of Acid (0.5 M HCI)
- [6:08] Holding pH 2.00
- [8:08] Stirrer speed set to 0
- [8:08] Stirrer speed set to 50
- [8:08] Iterative adjust 1.98 -> 2.00
- [8:08] pH 1.98 -> 2.00
- [8:09] Air gap released for Base (0.5 M KOH)
- [8:10] Dispensed 0.001623 mL of Base (0.5 M KOH)
- [9:00] Stirrer speed set to 0
- [9:15] Datapoint id 1 collected
- [9:15] Stirrer speed set to 50
- [9:20] pH 2.08 -> 2.28
- [9:20] Using cautious pH adjust
- [9:20] Dispensed 0.006656 mL of Base (0.5 M KOH)
- [9:25] Stepping pH = 2.17
- [9:26] Dispensed 0.005456 mL of Base (0.5 M KOH)
- [9:31] Stepping pH = 2.25
- [9:31] Dispensed 0.001787 mL of Base (0.5 M KOH)
- [9:36] Stepping pH = 2.29
- [9:51] Stirrer speed set to 0
- [10:02] Datapoint id 2 collected
- [10:02] Charge balance equation is out by -4.4%
- [10:02] Stirrer speed set to 50



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [10:07] pH 2.29 -> 2.49
- [10:07] Using charge balance adjust
- [10:07] Dispensed 0.008749 mL of Base (0.5 M KOH)
- [10:27] Stirrer speed set to 0
- [10:38] Datapoint id 3 collected
- [10:38] Charge balance equation is out by -0.3%
- [10:38] Stirrer speed set to 50
- [10:43] pH 2.50 -> 2.70
- [10:43] Using charge balance adjust
- [10:43] Dispensed 0.006044 mL of Base (0.5 M KOH)
- [11:03] Stirrer speed set to 0
- [11:13] Datapoint id 4 collected
- [11:13] Charge balance equation is out by -8.9%
- [11:13] Stirrer speed set to 50
- [11:18] pH 2.68 -> 2.88
- [11:18] Using charge balance adjust
- [11:19] Dispensed 0.004563 mL of Base (0.5 M KOH)
- [11:39] Stirrer speed set to 0
- [11:49] Datapoint id 5 collected
- [11:49] Charge balance equation is out by -3.6%
- [11:49] Stirrer speed set to 50
- [11:54] pH 2.88 -> 3.08
- [11:54] Using charge balance adjust
- [11:54] Dispensed 0.003551 mL of Base (0.5 M KOH)
- [12:14] Stirrer speed set to 0
- [12:25] Datapoint id 6 collected
- [12:25] Charge balance equation is out by 8.3%
- [12:25] Stirrer speed set to 50
- [12:30] pH 3.10 -> 3.30
- [12:30] Using charge balance adjust
- [12:30] Dispensed 0.002728 mL of Base (0.5 M KOH)
- [12:50] Stirrer speed set to 0
- [13:00] Datapoint id 7 collected
- [13:00] Charge balance equation is out by -4.3%
- [13:00] Stirrer speed set to 50
- [13:05] pH 3.29 -> 3.49
- [13:05] Using charge balance adjust
- [13:06] Dispensed 0.002117 mL of Base (0.5 M KOH)
- [13:26] Stirrer speed set to 0
- [13:49] Datapoint id 8 collected
- [13:49] Charge balance equation is out by 12.1%
- [13:49] Stirrer speed set to 50
- [13:54] pH 3.51 -> 3.71
- [13:54] Using charge balance adjust
- [13:54] Dispensed 0.001505 mL of Base (0.5 M KOH)
- [14:14] Stirrer speed set to 0
- [14:31] Datapoint id 9 collected
- [14:31] Charge balance equation is out by 23.9%
- [14:31] Stirrer speed set to 50
- [14:36] pH 3.76 -> 3.96
- [14:36] Using cautious pH adjust
- [14:36] Dispensed 0.000494 mL of Base (0.5 M KOH)
- [14:42] Stepping pH = 3.86
- [14:42] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [14:47] Stepping pH = 3.94
- [14:47] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [14:52] Stepping pH = 3.92
- [14:52] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [14:57] Stepping pH = 4.16



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [15:12] Stirrer speed set to 0
- [15:32] Datapoint id 10 collected
- [15:32] Charge balance equation is out by -51.7%
- [15:32] Stirrer speed set to 50
- [15:37] pH 4.11 -> 4.31
- [15:37] Using cautious pH adjust
- [15:37] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [15:42] Stepping pH = 4.20
- [15:43] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [15:48] Stepping pH = 4.28
- [15:48] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [15:53] Stepping pH = 4.28
- [15:53] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [15:58] Stepping pH = 4.50
- [16:13] Stirrer speed set to 0
- [16:51] Datapoint id 11 collected
- [16:51] Charge balance equation is out by -54.0%
- [16:51] Stirrer speed set to 50
- [16:56] pH 4.45 -> 4.65
- [16:56] Using cautious pH adjust
- [16:56] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [17:01] Stepping pH = 4.46
- [17:01] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [17:06] Stepping pH = 4.97
- [17:21] Stirrer speed set to 0
- [18:04] Datapoint id 12 collected
- [18:04] Charge balance equation is out by -97.0%
- [18:04] Stirrer speed set to 50
- [18:09] pH 4.93 -> 5.13
- [18:09] Using cautious pH adjust
- [18:09] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [18:15] Stepping pH = 4.94
- [18:15] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [18:20] Stepping pH = 5.22
- [18:35] Stirrer speed set to 0
- [19:26] Datapoint id 13 collected [19:26] Charge balance equation is out by -95.2%
- [19:26] Stirrer speed set to 50
- [19:31] pH 5.23 -> 5.43
- [19:31] Using cautious pH adjust
- [19:31] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [19:36] Stepping pH = 5.23
- [19:36] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [19:41] Stepping pH = 5.64
- [19:56] Stirrer speed set to 0
- [20:56] Datapoint id 14 collected
- [20:56] Charge balance equation is out by -204.5%
- [20:56] Stirrer speed set to 50
- [21:02] pH 5.71 -> 5.91
- [21:02] Using cautious pH adjust
- [21:02] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [21:07] Stepping pH = 5.75
- [21:07] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [21:12] Stepping pH = 5.78
- [21:12] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [21:17] Stepping pH = 6.41
- [21:32] Stirrer speed set to 0
- [22:32] Datapoint id 15 collected
- [22:32] Charge balance equation is out by -213.0%



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [22:32] Stirrer speed set to 50
- [22:37] pH 6.50 -> 6.70
- [22:37] Using cautious pH adjust
- [22:37] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [22:42] Stepping pH = 6.50
- [22:42] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [22:47] Stepping pH = 6.64
- [22:48] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [22:53] Stepping pH = 6.88
- [23:08] Stirrer speed set to 0
- [24:08] Datapoint id 16 collected
- [24:08] Charge balance equation is out by -164.9%
- [24:08] Stirrer speed set to 50
- [24:13] pH 6.94 -> 7.14
- [24:13] Using cautious pH adjust
- [24:13] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [24:18] Stepping pH = 6.91
- [24:18] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [24:23] Stepping pH = 7.15
- [24:38] Stirrer speed set to 0
- [25:38] Datapoint id 17 collected
- [25:38] Charge balance equation is out by -314.2%
- [25:38] Stirrer speed set to 50
- [25:43] pH 7.53 -> 7.73
- [25:43] Using cautious pH adjust
- [25:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [25:49] Stepping pH = 7.51
- [25:49] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [25:54] Stepping pH = 7.55
- [25:54] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [25:59] Stepping pH = 7.84
- [26:14] Stirrer speed set to 0
- [27:14] Datapoint id 18 collected
- [27:14] Charge balance equation is out by -916.6%
- [27:14] Stirrer speed set to 50
- [27:19] pH 8.10 -> 8.30
- [27:19] Using cautious pH adjust
- [27:19] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [27:24] Stepping pH = 8.05
- [27:24] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [27:29] Stepping pH = 8.01
- [27:29] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [27:35] Stepping pH = 8.98
- [27:50] Stirrer speed set to 0
- [28:46] Datapoint id 19 collected
- [28:46] Charge balance equation is out by -2,236.6%
- [28:46] Stirrer speed set to 50
- [28:51] pH 8.92 -> 9.05
- [28:51] Using cautious pH adjust
- [28:51] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [28:56] Stepping pH = 8.92
- [28:56] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [29:01] Stepping pH = 9.11
- [29:16] Stirrer speed set to 0
- [29:53] Datapoint id 20 collected
- [29:53] Charge balance equation is out by -259.1%
- [29:53] Titration 2 of 3
- [29:53] Adding initial titrants
- [29:53] Automatically add 0.03000 mL of Octanol



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [29:54] Dispensed 0.030009 mL of Octanol
- [29:54] Stirrer speed set to 10
- [29:55] Stirrer speed set to 55
- [29:55] Iterative adjust 9.14 -> 2.00
- [29:55] pH 9.14 -> 2.00
- [29:57] Dispensed 0.054492 mL of Acid (0.5 M HCI)
- [30:02] pH 2.02 -> 2.00
- [30:02] Dispensed 0.001740 mL of Acid (0.5 M HCI)
- [30:52] Stirrer speed set to 0
- [31:02] Datapoint id 21 collected
- [31:02] Stirrer speed set to 55
- [31:07] pH 2.02 -> 2.22
- [31:07] Using cautious pH adjust
- [31:08] Dispensed 0.008396 mL of Base (0.5 M KOH)
- [31:13] Stepping pH = 2.10
- [31:13] Dispensed 0.007079 mL of Base (0.5 M KOH)
- [31:18] Stepping pH = 2.19
- [31:18] Dispensed 0.001576 mL of Base (0.5 M KOH)
- [31:23] Stepping pH = 2.21
- [31:38] Stirrer speed set to 0
- 01.30 Stiller Speed Set to 0
- [31:48] Datapoint id 22 collected
- [31:48] Charge balance equation is out by -1.5%
- [31:48] Stirrer speed set to 55
- [31:54] pH 2.22 -> 2.42
- [31:54] Using charge balance adjust
- [31:54] Dispensed 0.011077 mL of Base (0.5 M KOH)
- [32:14] Stirrer speed set to 0
- [32:25] Datapoint id 23 collected
- [32:25] Charge balance equation is out by -11.5%
- [32:25] Stirrer speed set to 55
- [32:30] pH 2.40 -> 2.60
- [32:30] Using charge balance adjust
- [32:30] Dispensed 0.007879 mL of Base (0.5 M KOH)
- [32:50] Stirrer speed set to 0
- [33:00] Datapoint id 24 collected
- [33:00] Charge balance equation is out by 4.0%
- [33:00] Stirrer speed set to 55
- [33:05] pH 2.61 -> 2.81
- [33:05] Using charge balance adjust
- [33:06] Dispensed 0.005503 mL of Base (0.5 M KOH)
- [33:26] Stirrer speed set to 0
- [33:36] Datapoint id 25 collected
- [33:36] Charge balance equation is out by -14.7%
- [33:36] Stirrer speed set to 55
- [33:41] pH 2.79 -> 2.99
- [33:41] Using charge balance adjust
- [33:41] Dispensed 0.004163 mL of Base (0.5 M KOH)
- [34:01] Stirrer speed set to 0
- [34:11] Datapoint id 26 collected
- [34:11] Charge balance equation is out by -5.7%
- [34:11] Stirrer speed set to 55
- [34:16] pH 2.98 -> 3.18
- [34:16] Using charge balance adjust
- [34:17] Dispensed 0.003081 mL of Base (0.5 M KOH)
- [34:37] Stirrer speed set to 0
- [34:47] Datapoint id 27 collected
- [34:47] Charge balance equation is out by -10.6%
- [34:47] Stirrer speed set to 55
- [34:52] pH 3.16 -> 3.36



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [34:52] Using charge balance adjust
- [34:52] Dispensed 0.002281 mL of Base (0.5 M KOH)
- [35:12] Stirrer speed set to 0
- [35:23] Datapoint id 28 collected
- [35:23] Charge balance equation is out by 15.1%
- [35:23] Stirrer speed set to 55
- [35:28] pH 3.40 -> 3.60
- [35:28] Using cautious pH adjust
- [35:28] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [35:33] Stepping pH = 3.52
- [35:33] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [35:38] Stepping pH = 3.58
- [35:38] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [35:43] Stepping pH = 3.59
- [35:43] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [35:49] Stepping pH = 3.61
- [36:04] Stirrer speed set to 0
- [36:14] Datapoint id 29 collected
- [36:14] Charge balance equation is out by 17.7%
- [36:14] Stirrer speed set to 55
- 26:101 pl 2 60 > 2 00
- [36:19] pH 3.60 -> 3.80
- [36:19] Using cautious pH adjust
- [36:19] Dispensed 0.000517 mL of Base (0.5 M KOH)
- [36:24] Stepping pH = 3.73
- [36:24] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [36:29] Stepping pH = 3.77
- [36:29] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [36:34] Stepping pH = 3.78
- [36:34] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [36:40] Stepping pH = 3.80
- [36:55] Stirrer speed set to 0
- [37:05] Datapoint id 30 collected
- [37:05] Charge balance equation is out by 16.3%
- [37:05] Stirrer speed set to 55
- [37:10] pH 3.80 -> 4.00
- [37:10] Using cautious pH adjust
- [37:10] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [37:15] Stepping pH = 3.91
- [37:15] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [37:20] Stepping pH = 3.99
- [37:20] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [37:25] Stepping pH = 3.99
- [37:25] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [37:30] Stepping pH = 4.02
- [37:45] Stirrer speed set to 0
- [37:56] Datapoint id 31 collected
- [37:56] Charge balance equation is out by 9.2%
- [37:56] Stirrer speed set to 55
- [38:01] pH 4.02 -> 4.22
- [38:01] Using charge balance adjust
- [38:01] Dispensed 0.000423 mL of Base (0.5 M KOH)
- [38:21] Stirrer speed set to 0
- [38:34] Datapoint id 32 collected
- [38:34] Charge balance equation is out by 10.4%
- [38:34] Stirrer speed set to 55
- [38:39] pH 4.26 -> 4.46
- [38:39] Using charge balance adjust
- [38:39] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [38:59] Stirrer speed set to 0



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [39:15] Datapoint id 33 collected
- [39:15] Charge balance equation is out by -4.0%
- [39:15] Stirrer speed set to 55
- [39:21] pH 4.47 -> 4.67
- [39:21] Using charge balance adjust
- [39:21] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [39:41] Stirrer speed set to 0
- [40:11] Datapoint id 34 collected
- [40:11] Charge balance equation is out by -37.4%
- [40:11] Stirrer speed set to 55
- [40:16] pH 4.64 -> 4.84
- [40:16] Using cautious pH adjust
- [40:17] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [40:22] Stepping pH = 4.65
- [40:22] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [40:27] Stepping pH = 4.98
- [40:42] Stirrer speed set to 0
- [40:58] Datapoint id 35 collected
- [40:58] Charge balance equation is out by -87.0%
- [40:58] Stirrer speed set to 55
- [41:03] pH 5.00 -> 5.20
- [41:03] Using cautious pH adjust
- [41:03] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [41:08] Stepping pH = 5.00
- [41:08] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [41:13] Stepping pH = 5.27
- [41:28] Stirrer speed set to 0
- [42:01] Datapoint id 36 collected
- [42:01] Charge balance equation is out by -98.9%
- [42:01] Stirrer speed set to 55
- [42:06] pH 5.33 -> 5.53
- [42:06] Using cautious pH adjust
- [42:06] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [42:11] Stepping pH = 5.34
- [42:11] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [42:16] Stepping pH = 5.59 [42:31] Stirrer speed set to 0
- [43:19] Datapoint id 37 collected
- [43:19] Charge balance equation is out by -92.9%
- [43:19] Stirrer speed set to 55
- [43:24] pH 5.70 -> 5.90
- [43:24] Using cautious pH adjust
- [43:24] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [43:29] Stepping pH = 5.72
- [43:29] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [43:34] Stepping pH = 6.12
- [43:49] Stirrer speed set to 0
- [44:49] Datapoint id 38 collected
- [44:49] Charge balance equation is out by -87.7%
- [44:49] Stirrer speed set to 55
- [44:54] pH 6.16 -> 6.36
- [44:54] Using cautious pH adjust
- [44:54] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [44:59] Stepping pH = 6.16
- [44:59] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [45:04] Stepping pH = 6.69
- [45:20] Stirrer speed set to 0
- [46:20] Datapoint id 39 collected
- [46:20] Charge balance equation is out by -96.1%



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [46:20] Stirrer speed set to 55
- [46:25] pH 6.83 -> 7.03
- [46:25] Using cautious pH adjust
- [46:25] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [46:30] Stepping pH = 6.85
- [46:30] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [46:35] Stepping pH = 6.96
- [46:35] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [46:40] Stepping pH = 7.20
- [46:55] Stirrer speed set to 0
- [47:55] Datapoint id 40 collected
- [47:55] Charge balance equation is out by -198.9%
- [47:55] Stirrer speed set to 55
- [48:00] pH 7.31 -> 7.51
- [48:00] Using cautious pH adjust
- [48:00] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [48:05] Stepping pH = 7.34
- [48:06] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [48:11] Stepping pH = 7.46
- [48:11] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [48:16] Stepping pH = 7.68
- [48:31] Stirrer speed set to 0
- [49:31] Datapoint id 41 collected
- [49:31] Charge balance equation is out by -452.4%
- [49:31] Stirrer speed set to 55
- [49:36] pH 7.54 -> 7.74
- [49:36] Using cautious pH adjust
- [49:36] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [49:41] Stepping pH = 7.55
- [49:41] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [49:46] Stepping pH = 7.73
- [50:02] Stirrer speed set to 0
- [51:02] Datapoint id 42 collected
- [51:02] Charge balance equation is out by -414.6%
- [51:02] Stirrer speed set to 55
- [51:07] pH 7.69 -> 7.89
- [51:07] Using cautious pH adjust
- [51:07] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [51:12] Stepping pH = 7.66
- [51:12] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [51:17] Stepping pH = 7.90
- [51:32] Stirrer speed set to 0
- [52:32] Datapoint id 43 collected
- [52:32] Charge balance equation is out by -498.4%
- [52:32] Stirrer speed set to 55
- [52:37] pH 8.03 -> 8.23
- [52:37] Using cautious pH adjust
- [52:37] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [52:42] Stepping pH = 8.05
- [52:42] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [52:47] Stepping pH = 8.24
- [53:02] Stirrer speed set to 0
- [54:03] Datapoint id 44 collected
- 54:03] Charge balance equation is out by -454.8%
- [54:03] Stirrer speed set to 55
- [54:08] pH 8.27 -> 8.47
- [54:08] Using cautious pH adjust
- [54:08] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [54:13] Stepping pH = 8.30



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [54:13] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [54:18] Stepping pH = 8.46
- [54:33] Stirrer speed set to 0
- [55:33] Datapoint id 45 collected
- [55:33] Charge balance equation is out by -298.4%
- [55:33] Stirrer speed set to 55
- [55:38] pH 8.50 -> 8.70
- [55:38] Using cautious pH adjust
- [55:38] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [55:43] Stepping pH = 8.53
- [55:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [55:48] Stepping pH = 8.65
- [55:48] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [55:54] Stepping pH = 8.76
- [56:09] Stirrer speed set to 0
- [56:39] Datapoint id 46 collected [56:39] Charge balance equation is out by -284.7%
- [56:39] Stirrer speed set to 55
- [56:44] pH 8.74 -> 8.94
- [56:44] Using cautious pH adjust
- [56:44] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [56:49] Stepping pH = 8.74
- [56:50] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [56:55] Stepping pH = 8.90
- [56:55] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [57:00] Stepping pH = 8.96
- [57:15] Stirrer speed set to 0
- [57:51] Datapoint id 47 collected
- [57:51] Charge balance equation is out by -204.3%
- [57:51] Stirrer speed set to 55
- [57:56] pH 8.98 -> 9.05
- [57:56] Using cautious pH adjust
- [57:56] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [58:01] Stepping pH = 8.98
- [58:01] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [58:06] Stepping pH = 9.02
- [58:21] Stirrer speed set to 0
- [58:48] Datapoint id 48 collected
- [58:48] Charge balance equation is out by -216.1%
- [58:48] Titration 3 of 3
- [58:48] Adding initial titrants
- [58:48] Automatically add 0.25000 mL of Octanol
- [58:53] Dispensed 0.250000 mL of Octanol
- [58:53] Stirrer speed set to 10
- [58:55] Stirrer speed set to 60
- [58:55] Iterative adjust 9.03 -> 2.00
- [58:55] pH 9.03 -> 2.00
- [58:56] Dispensed 0.056044 mL of Acid (0.5 M HCI)
- [59:01] pH 2.04 -> 2.00
- [59:01] Dispensed 0.004139 mL of Acid (0.5 M HCI)
- [59:52] Stirrer speed set to 0
- [1:00:02] Datapoint id 49 collected
- [1:00:02] Stirrer speed set to 60
- [1:00:07] pH 1.99 -> 2.19
- [1:00:07] Using cautious pH adjust
- [1:00:07] Dispensed 0.009666 mL of Base (0.5 M KOH)
- [1:00:12] Stepping pH = 2.07
- [1:00:12] Dispensed 0.009055 mL of Base (0.5 M KOH)
- [1:00:17] Stepping pH = 2.16



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [1:00:18] Dispensed 0.002258 mL of Base (0.5 M KOH)
- [1:00:23] Stepping pH = 2.19
- [1:00:38] Stirrer speed set to 0
- [1:00:48] Datapoint id 50 collected
- [1:00:48] Charge balance equation is out by -8.6%
- [1:00:48] Stirrer speed set to 60
- [1:00:53] pH 2.19 -> 2.39
- [1:00:53] Using charge balance adjust
- [1:00:53] Dispensed 0.012912 mL of Base (0.5 M KOH)
- [1:01:14] Stirrer speed set to 0
- [1:01:24] Datapoint id 51 collected
- [1:01:24] Charge balance equation is out by -8.4%
- [1:01:24] Stirrer speed set to 60
- [1:01:29] pH 2.37 -> 2.57
- [1:01:29] Using charge balance adjust
- [1:01:29] Dispensed 0.008984 mL of Base (0.5 M KOH)
- [1:01:49] Stirrer speed set to 0
- [1:01:59] Datapoint id 52 collected
- [1:01:59] Charge balance equation is out by 2.2%
- [1:01:59] Stirrer speed set to 60
- [1:02:04] pH 2.58 -> 2.78
- [1:02:04] Using charge balance adjust
- [1:02:05] Dispensed 0.006044 mL of Base (0.5 M KOH)
- [1:02:25] Stirrer speed set to 0
- [1:02:35] Datapoint id 53 collected
- [1:02:35] Charge balance equation is out by -7.6%
- [1:02:35] Stirrer speed set to 60
- [1:02:40] pH 2.77 -> 2.97
- [1:02:40] Using charge balance adjust
- [1:02:40] Dispensed 0.004233 mL of Base (0.5 M KOH)
- [1:03:00] Stirrer speed set to 0
- [1:03:15] Datapoint id 54 collected
- [1:03:15] Charge balance equation is out by 7.6%
- [1:03:15] Stirrer speed set to 60
- [1:03:20] pH 2.99 -> 3.19
- [1:03:20] Using charge balance adjust
- [1:03:20] Dispensed 0.002752 mL of Base (0.5 M KOH)
- [1:03:40] Stirrer speed set to 0
- [1:03:50] Datapoint id 55 collected
- [1:03:50] Charge balance equation is out by 19.2%
- [1:03:50] Stirrer speed set to 60
- [1:03:55] pH 3.23 -> 3.43
- [1:03:55] Using cautious pH adjust
- [1:03:55] Dispensed 0.000847 mL of Base (0.5 M KOH)
- [1:04:01] Stepping pH = 3.33
- [1:04:01] Dispensed 0.000541 mL of Base (0.5 M KOH)
- [1:04:06] Stepping pH = 3.41
- [1:04:06] Dispensed 0.000118 mL of Base (0.5 M KOH)
  - .04.00] Dispensed 0.000 110
- [1:04:11] Stepping pH = 3.42 [1:04:11] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [1:04:16] Stepping pH = 3.43
- [1:04:31] Stirrer speed set to 0
- [1:04:41] Datapoint id 56 collected
- [1:04:41] Charge balance equation is out by 1.6%
- [1:04:41] Stirrer speed set to 60
- [1:04:46] pH 3.44 -> 3.64
- [1:04:46] Using charge balance adjust
- [1:04:46] Dispensed 0.001105 mL of Base (0.5 M KOH)
- [1:05:07] Stirrer speed set to 0



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [1:05:22] Datapoint id 57 collected
- [1:05:22] Charge balance equation is out by 27.7%
- [1:05:22] Stirrer speed set to 60
- [1:05:27] pH 3.69 -> 3.89
- [1:05:27] Using cautious pH adjust
- [1:05:27] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [1:05:32] Stepping pH = 3.76
- [1:05:32] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [1:05:37] Stepping pH = 3.92
- [1:05:52] Stirrer speed set to 0
- [1:06:03] Datapoint id 58 collected
- [1:06:03] Charge balance equation is out by -0.4%
- [1:06:03] Stirrer speed set to 60
- [1:06:08] pH 3.93 -> 4.13
- [1:06:08] Using charge balance adjust
- [1:06:08] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [1:06:28] Stirrer speed set to 0
- [1:06:39] Datapoint id 59 collected
- [1:06:39] Charge balance equation is out by 41.3%
- [1:06:39] Stirrer speed set to 60
- [1:06:44] pH 4.23 -> 4.43
- [1:06:44] Using cautious pH adjust
- [1:06:44] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [1:06:49] Stepping pH = 4.24
- [1:06:49] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [1:06:54] Stepping pH = 4.79
- [1:07:09] Stirrer speed set to 0
- 1.07.09] Stiller speed set to
- [1:07:54] Datapoint id 60 collected
- [1:07:54] Charge balance equation is out by -83.8%
- [1:07:54] Stirrer speed set to 60
- [1:08:00] pH 4.78 -> 4.98
- [1:08:00] Using cautious pH adjust
- [1:08:00] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:08:05] Stepping pH = 4.77
- [1:08:05] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [1:08:10] Stepping pH = 5.22
- [1:08:25] Stirrer speed set to 0
- [1:08:41] Datapoint id 61 collected
- [1:08:41] Charge balance equation is out by -208.6%
- [1:08:41] Stirrer speed set to 60
- [1:08:46] pH 5.30 -> 5.50
- [1:08:46] Using cautious pH adjust
- [1:08:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:08:51] Stepping pH = 5.31
- [1:08:51] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:08:56] Stepping pH = 5.45
- [1:08:56] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:09:01] Stepping pH = 5.56
- [1:09:16] Stirrer speed set to 0
- [1:10:05] Datapoint id 62 collected
- [1:10:05] Charge balance equation is out by -166.1%
- [1:10:05] Stirrer speed set to 60
- [1:10:10] pH 5.72 -> 5.92
- [1:10:10] Using cautious pH adjust
- [1:10:10] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:10:15] Stepping pH = 5.75
- [1:10:15] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:10:20] Stepping pH = 6.09
- [1:10:35] Stirrer speed set to 0



Assay ID: 18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

- [1:11:35] Datapoint id 63 collected
- [1:11:35] Charge balance equation is out by -74.3%
- [1:11:35] Stirrer speed set to 60
- [1:11:41] pH 6.22 -> 6.42
- [1:11:41] Using cautious pH adjust
- [1:11:41] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:11:46] Stepping pH = 6.24
- [1:11:46] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:11:51] Stepping pH = 6.31
- [1:11:51] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [1:11:56] Stepping pH = 6.87
- [1:12:11] Stirrer speed set to 0
- [1:13:11] Datapoint id 64 collected
- [1:13:11] Charge balance equation is out by -176.8%
- [1:13:11] Stirrer speed set to 60
- [1:13:16] pH 7.12 -> 7.32
- [1:13:16] Using cautious pH adjust
- [1:13:16] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:13:21] Stepping pH = 7.20
- [1:13:21] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:13:27] Stepping pH = 7.33
- [1:13:42] Stirrer speed set to 0
- [1:14:42] Datapoint id 65 collected
- [1:14:42] Charge balance equation is out by -102.3%
- [1:14:42] Stirrer speed set to 60
- [1:14:47] pH 7.16 -> 7.36
- [1:14:47] Using cautious pH adjust
- [1:14:47] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:14:52] Stepping pH = 7.24
- [1:14:52] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:14:57] Stepping pH = 7.41
- [1:15:12] Stirrer speed set to 0
- [1:16:12] Datapoint id 66 collected
- [1:16:12] Charge balance equation is out by -115.3%
- [1:16:12] Stirrer speed set to 60
- [1:16:17] pH 7.58 -> 7.78
- [1:16:17] Using cautious pH adjust
- [1:16:17] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:16:22] Stepping pH = 7.73
- [1:16:22] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:16:28] Stepping pH = 7.93
- [1:16:43] Stirrer speed set to 0
- [1:17:43] Datapoint id 67 collected
- [1:17:43] Charge balance equation is out by -302.3%
- [1:17:43] Stirrer speed set to 60
- [1:17:48] pH 7.75 -> 7.95
- [1:17:48] Using cautious pH adjust
- [1:17:48] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:17:53] Stepping pH = 7.87
- [1:17:53] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:17:58] Stepping pH = 8.15
- [1:18:13] Stirrer speed set to 0
- [1:19:13] Datapoint id 68 collected
- [1:19:13] Charge balance equation is out by -373.8%
- [1:19:13] Stirrer speed set to 60
- [1:19:18] pH 8.24 -> 8.44
- [1:19:18] Using cautious pH adjust
- [1:19:18] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:19:23] Stepping pH = 8.26

#### **Experiment Log**



Assay ID:

Sample name: M02\_octanol Experiment start time: 3/3/2018 1:46:01 PM
Assay name: pH-metric high logP Analyst: Dorothy Levorse

18C-03011 Instrument ID: T312060

Filename: C:\Sirius\_T3\Mehtap\20180302\_exp29\_logP\_T3-2\18C-03011\_M02\_octanol\_pH-metric high logP.t3r

## Experiment Log (continued)

[1:19:23] Dispensed 0.000024 mL of Base (0.5 M KOH)

[1:19:28] Stepping pH = 8.42

[1:19:29] Dispensed 0.000024 mL of Base (0.5 M KOH)

[1:19:34] Stepping pH = 8.58

[1:19:49] Stirrer speed set to 0

[1:20:49] Datapoint id 69 collected

[1:20:49] Charge balance equation is out by -452.3%

[1:20:49] Stirrer speed set to 60

[1:20:54] pH 8.61 -> 8.81

[1:20:54] Using cautious pH adjust

[1:20:54] Dispensed 0.000024 mL of Base (0.5 M KOH)

[1:20:59] Stepping pH = 8.61

[1:20:59] Dispensed 0.000024 mL of Base (0.5 M KOH)

[1:21:04] Stepping pH = 8.60

[1:21:04] Dispensed 0.000165 mL of Base (0.5 M KOH)

[1:21:09] Stepping pH = 9.23

[1:21:24] Stirrer speed set to 0

[1:21:39] Datapoint id 70 collected

[1:21:39] Charge balance equation is out by -904.7%

[1:21:39] Argon flow rate set to 0

[1:21:43] Titrator arm moved over Titration position