

Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

pH-metric Result

logP (XH +) 0.09 ±0.08 (n=50) logP (neutral X) 2.99 ±0.03 (n=50)

18C-26010 Points 2 to 16

M13_octanol concentration factor 1.006
Carbonate 0.0000 mM
Acidity error 2.38959 mM

18C-26010 Points 17 to 40

M13_octanol concentration factor 0.639
Carbonate 0.1054 mM
Acidity error 2.15628 mM

18C-26010 Points 41 to 68

M13_octanol concentration factor 0.745
Carbonate 0.1125 mM
Acidity error 2.15554 mM

Warnings and errors

Errors None Warnings None

рΗ

Sample logD and percent species

M13_octanol M13_octanol M13_octanol

•	_	_	_	_	_	
	logD	M13_octanolH	M13_octanol	M13_octanolH*	M13_octanol*	
1.000	0.09	44.67 %	0.00 %	54.60 %	0.74 %	
1.200	0.10	44.47 %	0.00 %	54.36 %	1.17 %	Stomach pH
2.000	0.14	41.88 %	0.01 %	51.19 %	6.93 %	·
3.000	0.46	25.78 %	0.04 %	31.52 %	42.65 %	
4.000	1.24	5.32 %	0.09 %	6.51 %	88.08 %	
5.000	2.15	0.60 %	0.10 %	0.73 %	98.57 %	
6.000	2.79	0.06 %	0.10 %	0.07 %	99.76 %	
6.500	2.91	0.02 %	0.10 %	0.02 %	99.86 %	
7.000	2.96	0.01 %	0.10 %	0.01 %	99.88 %	
7.400	2.98	0.00 %	0.10 %	0.00 %	99.89 %	Blood pH
8.000	2.99	0.00 %	0.10 %	0.00 %	99.90 %	
9.000	2.99	0.00 %	0.10 %	0.00 %	99.90 %	
10.000	2.99	0.00 %	0.10 %	0.00 %	99.90 %	
11.000	2.99	0.00 %	0.10 %	0.00 %	99.90 %	
12.000	2.99	0.00 %	0.10 %	0.00 %	99.90 %	

M13_octanol

M13_octanol Comment

Experiment start time: 3/26/2018 4:16:35 PM

Pion



Assay ID:


Filename:

Sample name: M13_octanol Assay name:

pH-metric high logP 18C-26010

Analyst:

Instrument ID: T312060 C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r






Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

Graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 1 of 3 18C-26010 Points 2 to 16

Overall results

RMSD 0.079
Average ionic strength 0.158 M
Average temperature 24.9°C
Partition ratio 0.1253 : 1

Analyte concentration range 2374.9 µM to 2437.5 µM

Total points considered 8 of 15

Warnings and errors

Errors None

Warnings One or more logP values out of range

Excessive acidity error present

Four-Plus parameters

Alpha 0.122 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r S 0.9974 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r jH 1.1 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r jOH -0.5 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r

Titrants


0.50 M HCI 0.994478 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r


0.50 M KOH 1.003190 3/26/2018 4:16:35 PM C:\Sirius_T3\KOH18C23.t3r

Sample

M13_octanol concentration factor 1.006
Base pKa 1 5.77
logP (XH +) -4.56
logP (neutral X) 2.93

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species


рН	M13_octanol	M13_octanol	M13_octanol		M13_octanol	
	logD	W13_octanoiH	M13_octanoi	M13_octanolH*	M13_octanol*	
1.000	-1.84	99.82 %	0.00 %	0.00 %	0.18 %	
1.200	-1.64	99.71 %	0.00 %	0.00 %	0.28 %	Stomach pH
2.000	-0.84	98.21 %	0.02 %	0.00 %	1.77 %	•
3.000	0.16	84.62 %	0.14 %	0.00 %	15.23 %	
4.000	1.15	35.50 %	0.60 %	0.00 %	63.90 %	
5.000	2.09	5.22 %	0.89 %	0.00 %	93.90 %	
6.000	2.73	0.55 %	0.93 %	0.00 %	98.52 %	
6.500	2.85	0.17 %	0.93 %	0.00 %	98.89 %	
7.000	2.90	0.05 %	0.93 %	0.00 %	99.01 %	
7.400	2.92	0.02 %	0.93 %	0.00 %	99.04 %	Blood pH
8.000	2.92	0.01 %	0.93 %	0.00 %	99.06 %	·
9.000	2.93	0.00 %	0.93 %	0.00 %	99.06 %	
10.000	2.93	0.00 %	0.93 %	0.00 %	99.07 %	
11.000	2.93	0.00 %	0.93 %	0.00 %	99.07 %	
12.000	2.93	0.00 %	0.93 %	0.00 %	99.07 %	

Carbonate and acidity



Other graphs








Assay name: pH-metric high logP Analyst: **Pion** Instrument ID: T312060 Assay ID: 18C-26010 Filename:

C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 2 of 3 18C-26010 Points 17 to 40

Overall results

RMSD 0.567
Average ionic strength 0.163 M
Average temperature 25.0°C
Partition ratio 0.2937 : 1

Analyte concentration range 1935.6 µM to 1984.1 µM

Total points considered 18 of 24

Warnings and errors

Errors None

Warnings Sample concentration factor out of range

One or more logP values out of range

Excessive acidity error present

Four-Plus parameters

Alpha 0.122 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r
S 0.9974 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r
jH 1.1 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r
jOH -0.5 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r


Titrants


0.50 M KOH 1.003190 3/26/2018 4:16:35 PM C:\Sirius T3\KOH18C23.t3r

Sample

M13_octanol concentration factor 0.639
Base pKa 1 5.77
logP (XH +) -4.56
logP (neutral X) 2.85

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M13_octanol	M13_octanol	M13_octanol	M13_octanol	M13_octanol	Comment
	logD	M13_octanolH	M13_octanol	M13_octanolH*	M13_octanol*	
1.000	-1.92	99.65 %	0.00 %	0.00 %	0.35 %	
1.200	-1.72	99.44 %	0.00 %	0.00 %	0.55 %	Stomach pH
2.000	-0.92	96.60 %	0.02 %	0.00 %	3.39 %	
3.000	0.08	73.95 %	0.13 %	0.00 %	25.93 %	
4.000	1.07	22.11 %	0.38 %	0.00 %	77.52 %	
5.000	2.01	2.76 %	0.47 %	0.00 %	96.77 %	
6.000	2.65	0.28 %	0.48 %	0.00 %	99.24 %	
6.500	2.77	0.09 %	0.48 %	0.00 %	99.43 %	
7.000	2.82	0.03 %	0.48 %	0.00 %	99.49 %	
7.400	2.84	0.01 %	0.48 %	0.00 %	99.51 %	Blood pH
8.000	2.84	0.00 %	0.48 %	0.00 %	99.52 %	
9.000	2.85	0.00 %	0.48 %	0.00 %	99.52 %	
10.000	2.85	0.00 %	0.48 %	0.00 %	99.52 %	
11.000	2.85	0.00 %	0.48 %	0.00 %	99.52 %	
12.000	2.85	0.00 %	0.48 %	0.00 %	99.52 %	

Carbonate and acidity



Other graphs






Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

pH-metric high logP Titration 3 of 3 18C-26010 Points 41 to 68

Overall results

RMSD 0.220
Average ionic strength 0.169 M
Average temperature 25.0°C
Partition ratio 0.6600 : 1

Analyte concentration range 1412.7 µM to 1440.3 µM

Total points considered 23 of 28

Warnings and errors

Errors None

Warnings One or more logP values out of range

Excessive acidity error present

Four-Plus parameters

Alpha 0.122 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r S 0.9974 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r jH 1.1 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r jOH -0.5 3/26/2018 4:16:34 PM C:\Sirius_T3\18C-26006_Blank standardisation.t3r


Titrants


0.50 M KOH 1.003190 3/26/2018 4:16:35 PM C:\Sirius_T3\KOH18C23.t3r

Sample

M13_octanol concentration factor 0.745
Base pKa 1 5.77
logP (XH +) -4.56
logP (neutral X) 2.73

Sample graphs








Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r


Sample graphs (continued)




Sample logD and percent species

рН	M13_octanol	M13_octanol	M13_octanol	M13_octanol	M13_octanol	Comment
	logD	M13_octanolH	M13_octanol	M13_octanolH*	M13_octanol*	
1.000	-2.04	99.39 %	0.00 %	0.00 %	0.60 %	
1.200	-1.84	99.04 %	0.00 %	0.00 %	0.95 %	Stomach pH
2.000	-1.04	94.26 %	0.02 %	0.00 %	5.72 %	
3.000	-0.04	62.18 %	0.11 %	0.00 %	37.72 %	
4.000	0.96	14.12 %	0.24 %	0.00 %	85.64 %	
5.000	1.90	1.62 %	0.27 %	0.00 %	98.11 %	
6.000	2.53	0.16 %	0.28 %	0.00 %	99.56 %	
6.500	2.66	0.05 %	0.28 %	0.00 %	99.67 %	
7.000	2.71	0.02 %	0.28 %	0.00 %	99.70 %	
7.400	2.72	0.01 %	0.28 %	0.00 %	99.71 %	Blood pH
8.000	2.73	0.00 %	0.28 %	0.00 %	99.72 %	
9.000	2.73	0.00 %	0.28 %	0.00 %	99.72 %	
10.000	2.73	0.00 %	0.28 %	0.00 %	99.72 %	
11.000	2.73	0.00 %	0.28 %	0.00 %	99.72 %	
12.000	2.73	0.00 %	0.28 %	0.00 %	99.72 %	

Carbonate and acidity



Other graphs






Pion Assay name: pH-metric high logP Analyst: 18C-26010 Instrument ID: T312060 Assay ID: Filename:

C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

Other graphs (continued)





Assay name: pH-metric high logP Analyst: Pion Assay ID: 18C-26010 Instrument ID: T312060

Filename: $C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric\ high\ logP.t3r$

Assay Model

Assay	Model												
Setting	S		Value		Date/	Time change	d I	Imported	from				
Sample	name		M13_o	ctanol	2/27/2	018 5:57:49 F	PM (User ente	ered valu	ıe			
Sample	by		Weight					Default va	alue				
Sample	weight		0.00126	30 g	3/26/2	018 4:15:03 F	PM (User ente	ered valu	ıe			
Formula	a weight		295.34	g/mol	2/27/2	:018 5:57:49 F	PM (User ente	ered valu	ıe			
Solubilit	ty		Unknov	vn			I	Default va	alue				
Molecul	ar weight		295.34		2/27/2	018 5:57:49 F	PM (User ente	ered valu	ıe			
Individu	al pKa ionic environments		No				- 1	Default va	alue				
Numbei	r of pKas		1		2/27/2	:018 5:57:49 F	PM (User ente	ered valu	ue			
Sample	is a		Base		2/27/2	018 5:57:49 F	PM (User ente	ered valu	ue			
pKa 1			5.77		2/27/2	018 5:57:49 F	PM (User ente	ered valu	ue			
logp (XI	H +)		-4.56		3/2/20	18 4:30:48 PM	M (User ente	ered valu	ue			
logP (ne	eutral X)		2.99		3/2/20	18 4:30:43 PI	M (User ente	ered valu	ne			
Event	S												
Time	Event	Wa	ter	Acid		Base	Oc	tanol	рН	dpH/dt	pH R-squared	pH SD	dpH/ time
5:17.0 5:18.1	Manual volume addition Initial pH = 3.49						0.2	0000 mL					unie
8.22.9	Data point 2	1.5	0000 ml	0 044	73 ml	0.00553 ml	0.2	0000 ml	2 003	-0 00959	0.73625	0.00055	10.5

Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	pH SD	dpH/e
5:17.0 5:18.1	Manual volume addition Initial pH = 3.49				0.20000 mL					unie
8:22.9	Data point 2	1.50000 mL	0.04473 mL	0.00553 mL	0.20000 mL	2.003	-0.00959	0.73625	0.00055	10.5 s
9:09.4	Data point 3	1.50000 mL	0.04473 mL	0.02016 mL	0.20000 mL	2.209	0.00564	0.12526	0.00079	10.0 S
9:45.0	Data point 4	1.50000 mL	0.04473 mL	0.02970 mL	0.20000 mL	2.418	0.01231	0.87847	0.00065	10.5
10:21.0	Data point 5	1.50000 mL	0.04473 mL	0.03568 mL	0.20000 mL	2.638	0.00365	0.22847	0.00038	10.5
10:57.1	Data point 6	1.50000 mL	0.04473 mL	0.03944 mL	0.20000 mL	2.855	-0.00035	0.00242	0.00035	10.0
11:32.6	Data point 7	1.50000 mL	0.04473 mL	0.04198 mL	0.20000 mL	3.087	0.00584	0.40770	0.00045	10.5

10:21.0	Data point 5	1.50000 mL	0.04473 mL	0.03568 mL	0.20000 mL	2.638	0.00365	0.22847	0.00038	10.5 s
10:57.1	Data point 6	1.50000 mL	0.04473 mL	0.03944 mL	0.20000 mL	2.855	-0.00035	0.00242	0.00035	10.0
11:32.6	Data point 7	1.50000 mL	0.04473 mL	0.04198 mL	0.20000 mL	3.087	0.00584	0.40770	0.00045	10.5
12:08.5	Data point 8	1.50000 mL	0.04473 mL	0.04384 mL	0.20000 mL	3.358	-0.00076	0.00844	0.00041	10.0 s
12:59.5	Data point 9	1.50000 mL	0.04473 mL	0.04501 mL	0.20000 mL	3.555	-0.00079	0.00582	0.00051	10.5 s
13:45.7	Data point 10	1.50000 mL	0.04473 mL	0.04603 mL	0.20000 mL	3.736	-0.01231	0.48567	0.00087	10.5 s
14:32.0	Data point 11	1.50000 mL	0.04473 mL	0.04734 mL	0.20000 mL	3.933	-0.00938	0.30311	0.00084	11.0 s
15:08.5	Data point 12	1.50000 mL	0.04473 mL	0.04842 mL	0.20000 mL	4.073	-0.01778	0.77327	0.00100	13.5 s
16:08.2	Data point 13	1.50000 mL	0.04473 mL	0.04972 mL	0.20000 mL	4.259	-0.01176	0.47760	0.00084	14.5 s
16:58.3	Data point 14	1.50000 mL	0.04473 mL	0.05045 mL	0.20000 mL	4.506	0.01843	0.89339	0.00096	21.0 s
17:44.7	Data point 15	1.50000 mL	0.04473 mL	0.05096 mL	0.20000 mL	4.931	0.01593	0.72074	0.00093	22.5 s
18:37.6	Data point 16	1.50000 mL	0.04473 mL	0.05172 mL	0.20000 mL	9.043	-0.08412	0.99442	0.00416	Time

Reported	d at: 3/30/2018 5:50:25 PM	М							Page 13	of 29
22:49.6	Data point 20	1.50000 mL	0.09854 mL	0.08826 mL	0.50000 mL	2.578	-0.01356	0.80667	0.00075	10.0 s
22:13.6	Data point 19	1.50000 mL	0.09854 mL	0.08095 mL	0.50000 mL	2.366	0.00068	0.00245	0.00068	10.5 s
21:38.0	Data point 18	1.50000 mL	0.09854 mL	0.06919 mL	0.50000 mL	2.151	-0.00296	0.10195	0.00046	10.0 s
20:42.7	Data point 17	1.50000 mL	0.09854 mL	0.05172 mL	0.50000 mL	1.945	0.00587	0.10684	0.00089	at 19.0 s
18:37.6	Data point 16	1.50000 mL	0.04473 mL	0.05172 mL	0.20000 mL	9.043	-0.08412	0.99442	0.00416	s Time out
17:44.7	Data point 15	1.50000 mL	0.04473 mL	0.05096 mL	0.20000 mL	4.931	0.01593	0.72074	0.00093	22.5
16:58.3	Data point 14	1.50000 mL	0.04473 mL	0.05045 mL	0.20000 mL	4.506	0.01843	0.89339	0.00096	21.0
16:08.2	Data point 13	1.50000 mL	0.04473 mL	0.04972 mL	0.20000 mL	4.259	-0.01176	0.47760	0.00084	14.5
15:08.5	Data point 12	1.50000 mL	0.04473 mL	0.04842 mL	0.20000 mL	4.073	-0.01778	0.77327	0.00100	13.5
14:32.0	Data point 11	1.50000 mL	0.04473 mL	0.04734 mL	0.20000 mL	3.933	-0.00938	0.30311	0.00084	11.0
13:45.7	Data point 10	1.50000 mL	0.04473 mL	0.04603 mL	0.20000 mL	3.736	-0.01231	0.48567	0.00087	10.5
12:59.5	Data point 9	1.50000 mL	0.04473 mL	0.04501 mL	0.20000 mL	3.555	-0.00079	0.00582	0.00051	10.5
12:08.5	Data point 8	1.50000 mL	0.04473 mL	0.04384 mL	0.20000 mL	3.358	-0.00076	0.00844	0.00041	10.0
11:32.6	Data point 7	1.50000 mL	0.04473 mL	0.04198 mL	0.20000 mL	3.087	0.00584	0.40770	0.00045	•
10.07.1	Bata point o	1.00000 1112	0.011101112	0.000111112	0.200002		0.00000	0.002.12	0.00000	S



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

Events (continued)

Lvenis (continu c a)									
Time	Event	Water	Acid	Base	Octanol	рН	dpH/dt	pH R-squared	-	dpH/dt time
23:25.1	Data point 21	1.50000 mL	0.09854 mL	0.09304 mL	0.50000 mL	2.792	0.00291	0.03033	0.00083	10.5 s
24:01.1	Data point 22	1.50000 mL	0.09854 mL	0.09631 mL	0.50000 mL	2.996	-0.00625	0.16494	0.00076	10.0 s
24:36.5	Data point 23	1.50000 mL	0.09854 mL	0.09878 mL	0.50000 mL	3.202	-0.00833	0.32865	0.00072	10.0 s
25:11.7	Data point 24	1.50000 mL	0.09854 mL	0.10073 mL	0.50000 mL	3.524	-0.01622	0.91829	0.00084	10.5 s
25:52.8	Data point 25	1.50000 mL	0.09854 mL	0.10167 mL	0.50000 mL	3.706	-0.01693	0.76864	0.00095	11.0 s
26:44.5	Data point 26	1.50000 mL	0.09854 mL	0.10261 mL	0.50000 mL	3.881	-0.01256	0.44869	0.00093	12.0 s
27:32.0	Data point 27	1.50000 mL	0.09854 mL	0.10332 mL	0.50000 mL	4.048	-0.01675	0.86339	0.00089	13.0 s
28:20.7	Data point 28	1.50000 mL	0.09854 mL	0.10383 mL	0.50000 mL	4.234	-0.01972	0.95699	0.00100	10.0 s
29:06.4	Data point 29	1.50000 mL	0.09854 mL	0.10430 mL	0.50000 mL	4.495	-0.01298	0.53508	0.00088	12.0 s
29:43.8	Data point 30	1.50000 mL	0.09854 mL	0.10459 mL	0.50000 mL	4.751	0.01337	0.62718	0.00083	23.5 s
30:37.9	Data point 31		0.09854 mL					0.50324	0.00083	
31:21.0	Data point 32							0.75428	0.00095	
32:05.0	Data point 33		0.09854 mL					0.99413		Timed out at 59.5 s
33:35.6	Data point 34	1.50000 mL	0.09854 mL	0.10529 mL	0.50000 mL	7.882	-0.07582	0.99362	0.00376	Timed out at 59.5 s
35:06.1	Data point 35	1.50000 mL	0.09854 mL	0.10536 mL	0.50000 mL	7.874	-0.03303	0.99176	0.00164	
36:36.6	Data point 36	1.50000 mL	0.09854 mL	0.10541 mL	0.50000 mL	7.975	-0.02481	0.97940	0.00124	Timed out at 59.5 s
38:07.1	Data point 37	1.50000 mL	0.09854 mL	0.10546 mL	0.50000 mL	8.262	-0.01947	0.94137	0.00099	50.5 s
39:28.1	Data point 38		0.09854 mL					0.77622	0.00096	
40:37.7	Data point 39		0.09854 mL					0.52717	0.00093	
41:37.8	Data point 40		0.09854 mL					0.70869	0.00091	
43:42.7	Data point 41		0.15644 mL					0.57402	0.00087	
44:28.9	Data point 42		0.15644 mL					0.01426	0.00040	
45:05.1	Data point 43		0.15644 mL					0.34363	0.00092	
45:48.7	Data point 44		0.15644 mL					0.09759	0.00069	
46:24.3	Data point 45		0.15644 mL					0.02683	0.00042	
47:00.3	Data point 46		0.15644 mL						0.00097	
47:35.8	Data point 47		0.15644 mL					0.33990	0.00047	
48:11.2	Data point 48		0.15644 mL					0.44316	0.00036	
48:57.0	Data point 49		0.15644 mL					0.96977	0.00078	
49:32.4	Data point 50		0.15644 mL					0.02822	0.00093	
50:16.3	Data point 51		0.15644 mL						0.00098	
50:58.8	Data point 52		0.15644 mL					0.75358	0.00099	
51:39.4	Data point 53		0.15644 mL						0.00098	
52:19.8	Data point 54		0.15644 mL						0.00083	
53:05.6	Data point 55		0.15644 mL						0.00045	
53:51.2	Data point 56		0.15644 mL						0.00043	
54:36.8			0.15644 mL						0.00004	
54.30.8 55:17.8	Data point 57		0.15644 mL						0.00089	
56:08.7	Data point 58		0.15644 mL							
50.06.7 57:12.7	Data point 59		0.15644 mL						0.00094	
57.12.7 58:37.2	Data point 60		0.15644 mL						0.00100	
59:49.7	Data point 61 Data point 62		0.15644 mL						0.00097 0.00137	Timed out
1:01:20.3	Data point 63	1.50000 mL	0.15644 mL	0.16326 mL	1.20000 mL	7.854	-0.03775	0.96659	0.00190	at 59.5 s Timed out at 59.5 s
1·02·50 a	Data point 64	1.50000 ml	0.15644 mL	0.16331 ml	1 20000 ml	8 232	-0 01636	0.69169	0.00097	
	Data point 65		0.15644 mL						0.00097	
	Data point 66		0.15644 mL						0.00098	
	Data point 67		0.15644 mL						0.00098	
	Data point 68		0.15644 mL						0.00094	
	Assay volumes					9.001	0.00042	0.00270	0.00034	20.0 3



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

o.tomas_rotwen	tapizo 1000zi _cxpo-	+_logi _10 ±1100	3 200 TO_IN TO_OCIANO	_prr metre mgn regr
Assay Settings				
Setting	Value	Original Value	Date/Time changed	Imported from
General Settings				
Analyst name	Pion			
Standard Experiment Settings				
Number of titrations	3			
Minimum pH	2.000			
Maximum pH	9.000			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.100002 mL			
Argon flow rate	100%			
Start titration using	Cautious pH adjust			
	Cautious pri aujust			
Advanced General Settings	Mana			
Detect turbidity using	None			
Collect turbidity sensor data	No			
Collect UV spectra	No .			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	10%			
Titrant Pre-Dose				
Titrant pre-dose	None			
Assay Medium				
ISA water volume	1.50 mL			
Water added	Automatic			
Partition solvent type	Octanol			
Partition volume	0.200 mL			
Partition solvent added	Manual in advance			
After partition addition, stir for	1 seconds			
Sample Sonication				
Sonicate	Yes			
Adjust pH for sonication	No			
Sonicate for	60 seconds			
After sonication stir for	20 seconds			
Sample Dissolution	20 00001100			
Perform a dissolution stage	Yes			
	To start pH			
Adjust and hold pH for dissolution	•			
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
Carbonate purge				
Perform a carbonate purge	No			
Temperature Control	.,			
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			
Stir speed of	50%			
Titration 1				
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	50%			
Titration 2				
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.300 mL			
Additional partition solvent added	Automatic			
enerne sould diring	ALL SELECTIONS			

Reported at: 3/30/2018 5:50:25 PM

30 seconds

15 seconds

55%

After pH adjust stir for

Stir to allow partitioning for

Stirrer speed for partitioning



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Titration 3		•	•	•
Titrate from	Low to high pH			
Add additional water	0.00 mL			
Additional partition solvent volume	0.700 mL			
Additional partition solvent added	Automatic			
After pH adjust stir for	30 seconds			
Stir to allow partitioning for	15 seconds			
Stirrer speed for partitioning	60%			
Data Point Stability				
Stir during data point collection	No			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00100 dpH/dt			
Stability timeout after	60 seconds			

Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.122	3/26/2018 4:16:34 PM	C:\Sirius_T3\18C-26006_Blank standardisation.t3r
Four-Plus S	0.9974	3/26/2018 4:16:34 PM	C:\Sirius_T3\18C-26006_Blank standardisation.t3r
Four-Plus jH			C:\Sirius_T3\18C-26006_Blank standardisation.t3r
Four-Plus jOH	-0.5	3/26/2018 4:16:34 PM	C:\Sirius_T3\18C-26006_Blank standardisation.t3r
Base concentration factor	1.003	3/26/2018 4:16:35 PM	C:\Sirius_T3\KOH18C23.t3r
Acid concentration factor	0.994	3/26/2018 4:16:34 PM	C:\Sirius_T3\18C-26006_Blank standardisation.t3r

Instrument Settings

Setting Instrument owner Instrument ID Instrument type	Value Merck T312060 T3 Simulator	Batch Id	Install date
Software version Dispenser module Dispenser 0 Syringe volume Firmware version	1.1.3.0 Water 2.5 mL 1.2.1(r2)	T3DM1200361	3/31/2009 6:24:52 AM 3/31/2009 6:25:05 AM
Titrant Dispenser 2 Syringe volume Firmware version	Water (0.15 M KCI) Acid 0.5 mL 1.2.1(r2)	02-06-2018	3/16/2018 11:09:18 AM 3/31/2009 6:25:11 AM
Titrant Dispenser 1 Syringe volume Firmware version	Acid (0.5 M HCl) Base 0.5 mL 1.2.1(r2)	03-16-2018	3/16/2018 10:56:23 AM 3/31/2009 6:25:21 AM
Titrant Dispenser 5 Syringe volume Firmware version	Base (0.5 M KOH) Cosolvent 2.5 mL 1.2.1(r2)	3/22/2018	3/23/2018 9:34:17 AM 3/31/2009 6:26:24 AM
Distribution valve 5 Firmware version	Distribution Valve 1.1.3	00 00 0040	3/31/2009 6:28:19 AM
Port A Port B Dispenser 3 Syringe volume Firmware version	Methanol (80%, 0.15 M KCI) Cyclohexane Buffer 0.5 mL 1.2.1(r2)	02-08-2018 11-01-17	3/6/2018 10:28:59 AM 2/27/2018 11:37:57 AM 8/3/2010 6:05:16 AM
Titrant Dispenser 6	Dodecane Octanol	2018/01/31	2/28/2018 11:18:04 AM 10/22/2010 11:52:43 AM

Reported at: 3/30/2018 5:50:25 PM



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Setting Syringe volume	Value 0.5 mL	Batch Id	Install date
Firmware version	1.2.1(r2)		
Titrant	Octanol	01-31-2018	2/27/2018 10:59:35 AM
Titrator	4.47.414.010.00.01	T3TM1200161	3/31/2009 6:24:17 AM
Horizontal axis firmware version	1.17 Al1DI2DO2 Stepper 2		
Vertical axis firmware version Chassis I/O firmware version	1.17 Al1Dl2DO2 Stepper 2 1.11 Al1Dl0DO4 Norgren I/O		
Probe I/O firmware version	1.11 AT DIODO4 Not great #0		
Electrode	T3 Electrode	T3E0923	1/23/2018 3:01:00 PM
E0 calibration	+5.08 mV	1020020	3/26/2018 4:17:19 PM
Filling solution	3M KCI	KCL097	3/26/2018 9:20:42 AM
Liquids			
Wash 1	50% IPA:50% Water		3/26/2018 9:21:48 AM
Wash 2	0.5% Trition X-100 in H20		3/26/2018 9:21:51 AM
Buffer position 1	pH7 Wash		3/26/2018 9:21:54 AM
Buffer position 2	pH 7		3/26/2018 9:21:57 AM
Storage position	6.651.0021	02 42 2040	3/26/2018 9:21:20 AM
Wash water Waste	6.6e+003 mL 3.7e+003 mL	03-12-2018	3/12/2018 9:25:04 AM 3/12/2018 9:24:49 AM
Temperature controller	3.7 6 1 003 IIIL		8/5/2010 7:35:13 AM
Turbidity detector			3/31/2009 6:24:45 AM
Spectrometer		074811	11/23/2010 12:22:28 PM
Dip probe		10196	
Wavelength coefficient A0	183.333		
Wavelength coefficient A1	2.21568		
Wavelength coefficient A2	-0.000289308		
Total lamp lit time	168:19:11		11/23/2010 12:22:28 PM
Calibrated on	2/27/2018 11:40:38 AM		
Integration time	40 10		
Scans averaged Autoloader	10	T3AL1200345	11/10/2015 10:34:13 AM
Left-right axis firmware version	1.17 Al1Dl2DO2 Stepper 2	13AL1200343	11/10/2015 10:54:15 AM
Front-back axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Vertical axis firmware version	1.17 Al1Dl2DO2 Stepper 2		
Chassis I/O firmware version	1.11 Al1Dl0DO4 Norgren I/O		
Configuration	_		
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period Titrant tube volume	5 minute(s) 1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points E0 calibration maximum standard deviation	10 0.01500		
E0 calibration maximum standard deviation E0 calibration timeout period	60 s		
E0 calibration timeout period E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		

Reported at: 3/30/2018 5:50:25 PM



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

Instrument Settings (continued)

Spectrometer calibration stir duration 5 s Spectrometer calibration stir speed 30% Spectrometer calibration wash pump volume 20.0 mL Spectrometer calibration wash stir duration 5 s Spectrometer calibration wash stir speed 30% Overhead dispense height 10000	Spectrometer calibration stir speed Spectrometer calibration wash pump volume Spectrometer calibration wash stir duration Spectrometer calibration wash stir speed	30% 20.0 mL 5 s 30%	Batch Id	d Install date	
--	---	------------------------------	----------	----------------	--

Refinement Settings

Setting	Value	Default value
Turbidity detection method	None	None
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance		0.100
Turbidity probe threshold	50.00	50.00
• •		

Experiment Log

- [2:41] Air gap created for Water (0.15 M KCI)
- [2:41] Air gap created for Acid (0.5 M HCI)
- [2:41] Air gap created for Base (0.5 M KOH)
- [2:42] Air gap released for Water (0.15 M KCI)
- [2:46] Titrator arm moved over Titration position
- [2:46] Titration 1 of 3
- [2:46] Adding initial titrants
- [2:46] Automatically add 1.50000 mL of water
- [3:11] Dispensed 1.500000 mL of Water (0.15 M KCI)
- [3:15] Titrator arm moved over Drain
- [4:56] Titrator arm moved to Titration position
- [4:56] Argon flow rate set to 100
- [4:56] Stirrer speed set to 10
- [5:18] Initial pH = 3.49
- [5:18] Iterative adjust 3.49 -> 2.00
- [5:18] pH 3.49 -> 2.00
- [5:19] Air gap released for Acid (0.5 M HCI)
- [5:20] Dispensed 0.042404 mL of Acid (0.5 M HCI)
- [5:25] pH 2.02 -> 2.00
- [5:25] Dispensed 0.002328 mL of Acid (0.5 M HCI)

- [5:30] Holding pH 2.00
- [7:30] Stirrer speed set to 0
- [7:30] Stirrer speed set to 50
- [7:30] Iterative adjust 1.95 -> 2.00
- [7:30] pH 1.95 -> 2.00
- [7:31] Air gap released for Base (0.5 M KOH)
- [7:32] Dispensed 0.005527 mL of Base (0.5 M KOH)
- [8:22] Stirrer speed set to 0
- [8:33] Datapoint id 2 collected
- [8:33] Stirrer speed set to 50
- [8:38] pH 2.01 -> 2.21
- [8:38] Using cautious pH adjust
- [8:38] Dispensed 0.007808 mL of Base (0.5 M KOH)
- [8:43] Stepping pH = 2.09
- [8:44] Dispensed 0.006068 mL of Base (0.5 M KOH)
- [8:49] Stepping pH = 2.19
- [8:49] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [8:54] Stepping pH = 2.21 [9:09] Stirrer speed set to 0
- [9:19] Datapoint id 3 collected
- [9:19] Charge balance equation is out by 6.4%
- [9:19] Stirrer speed set to 50



pH-metric high logP **Pion** Assay name: Analyst: Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [9:24] pH 2.22 -> 2.42
- [9:24] Using charge balance adjust
- [9:24] Dispensed 0.009548 mL of Base (0.5 M KOH)
- [9:45] Stirrer speed set to 0
- [9:55] Datapoint id 4 collected
- [9:55] Charge balance equation is out by -0.4%
- [9:55] Stirrer speed set to 50
- [10:00] pH 2.43 -> 2.63
- [10:00] Using charge balance adjust
- [10:00] Dispensed 0.005974 mL of Base (0.5 M KOH)
- [10:21] Stirrer speed set to 0
- [10:31] Datapoint id 5 collected
- [10:31] Charge balance equation is out by 5.5%
- [10:31] Stirrer speed set to 50
- [10:36] pH 2.65 -> 2.85
- [10:36] Using charge balance adjust
- [10:36] Dispensed 0.003763 mL of Base (0.5 M KOH)
- [10:57] Stirrer speed set to 0
- [11:07] Datapoint id 6 collected
- [11:07] Charge balance equation is out by 4.1%
- [11:07] Stirrer speed set to 50
- [11:12] pH 2.87 -> 3.07
- [11:12] Using charge balance adjust
- [11:12] Dispensed 0.002540 mL of Base (0.5 M KOH)
- [11:32] Stirrer speed set to 0
- [11:43] Datapoint id 7 collected
- [11:43] Charge balance equation is out by 10.5%
- [11:43] Stirrer speed set to 50
- [11:48] pH 3.10 -> 3.30
- [11:48] Using charge balance adjust
- [11:48] Dispensed 0.001858 mL of Base (0.5 M KOH)
- [12:08] Stirrer speed set to 0
- [12:18] Datapoint id 8 collected
- [12:18] Charge balance equation is out by 30.9%
- [12:18] Stirrer speed set to 50
- [12:23] pH 3.36 -> 3.56
- [12:23] Using cautious pH adjust
- [12:23] Dispensed 0.000753 mL of Base (0.5 M KOH)
- [12:28] Stepping pH = 3.52
- [12:29] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [12:34] Stepping pH = 3.55
- [12:34] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [12:39] Stepping pH = 3.55
- [12:39] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [12:44] Stepping pH = 3.57 [12:59] Stirrer speed set to 0
- [13:10] Datapoint id 9 collected
- [13:10] Charge balance equation is out by 23.3%
- [13:10] Stirrer speed set to 50
- [13:15] pH 3.56 -> 3.76
- [13:15] Using cautious pH adjust
- [13:15] Dispensed 0.000682 mL of Base (0.5 M KOH)
- [13:20] Stepping pH = 3.73
- [13:20] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [13:25] Stepping pH = 3.74
- [13:25] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [13:30] Stepping pH = 3.77
- [13:45] Stirrer speed set to 0
- [13:56] Datapoint id 10 collected



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [13:56] Charge balance equation is out by 25.9%
- [13:56] Stirrer speed set to 50
- [14:01] pH 3.74 -> 3.94
- [14:01] Using cautious pH adjust
- [14:01] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [14:06] Stepping pH = 3.89
- [14:06] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [14:11] Stepping pH = 3.90
- [14:11] Dispensed 0.000494 mL of Base (0.5 M KOH)
- [14:16] Stepping pH = 4.00
- [14:32] Stirrer speed set to 0
- [14:43] Datapoint id 11 collected
- [14:43] Charge balance equation is out by -4.2%
- [14:43] Stirrer speed set to 50
- [14:48] pH 3.94 -> 4.14
- [14:48] Using charge balance adjust
- [14:48] Dispensed 0.001082 mL of Base (0.5 M KOH)
- [15:08] Stirrer speed set to 0
- [15:22] Datapoint id 12 collected
- [15:22] Charge balance equation is out by -32.8%
- [15:22] Stirrer speed set to 50
- [15:27] pH 4.08 -> 4.28
- [15:27] Using cautious pH adjust
- [15:27] Dispensed 0.000470 mL of Base (0.5 M KOH)
- [15:32] Stepping pH = 4.19
- [15:32] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [15:37] Stepping pH = 4.23
- [15:37] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [15:42] Stepping pH = 4.26
- [15:42] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [15:47] Stepping pH = 4.26
- [15:48] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [15:53] Stepping pH = 4.34
- [16:08] Stirrer speed set to 0
- [16:22] Datapoint id 13 collected
- [16:22] Charge balance equation is out by -36.9%
- [16:22] Stirrer speed set to 50
- [16:27] pH 4.26 -> 4.46
- [16:27] Using cautious pH adjust
- [16:27] Dispensed 0.000376 mL of Base (0.5 M KOH)
- [16:33] Stepping pH = 4.41
- [16:33] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [16:38] Stepping pH = 4.41
- [16:38] Dispensed 0.000259 mL of Base (0.5 M KOH)
- [16:43] Stepping pH = 4.58 [16:58] Stirrer speed set to 0
- [17:19] Datapoint id 14 collected
- [17:19] Charge balance equation is out by -0.5%
- [17:19] Stirrer speed set to 50
- [17:24] pH 4.50 -> 4.70
- [17:24] Using charge balance adjust
- [17:24] Dispensed 0.000517 mL of Base (0.5 M KOH)
- [17:44] Stirrer speed set to 0
- [18:07] Datapoint id 15 collected
- [18:07] Charge balance equation is out by 114.9%
- [18:07] Stirrer speed set to 50
- [18:12] pH 4.88 -> 5.08
- [18:12] Using cautious pH adjust
- [18:12] Dispensed 0.000118 mL of Base (0.5 M KOH)



pH-metric high logP **Pion** Assay name: Analyst: Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [18:17] Stepping pH = 4.87
- [18:17] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [18:22] Stepping pH = 9.43
- [18:37] Stirrer speed set to 0
- [19:37] Datapoint id 16 collected
- [19:37] Charge balance equation is out by -207.2%
- [19:37] Titration 2 of 3
- [19:37] Adding initial titrants
- [19:37] Automatically add 0.30000 mL of Octanol
- [19:44] Dispensed 0.300000 mL of Octanol
- [19:44] Stirrer speed set to 10
- [19:45] Stirrer speed set to 55
- [19:45] Iterative adjust 9.03 -> 2.00
- [19:45] pH 9.03 -> 2.00
- [19:47] Dispensed 0.052140 mL of Acid (0.5 M HCI)
- [19:52] pH 2.02 -> 2.00
- [19:52] Dispensed 0.001670 mL of Acid (0.5 M HCl)
- [20:42] Stirrer speed set to 0
- [21:01] Datapoint id 17 collected
- [21:01] Stirrer speed set to 55
- [21:06] pH 1.95 -> 2.15
- [21:06] Using cautious pH adjust [21:07] Dispensed 0.009525 mL of Base (0.5 M KOH)
- [21:12] Stepping pH = 2.05
- [21:12] Dispensed 0.006279 mL of Base (0.5 M KOH)
- [21:17] Stepping pH = 2.13
- [21:17] Dispensed 0.001670 mL of Base (0.5 M KOH)
- [21:22] Stepping pH = 2.15
- [21:38] Stirrer speed set to 0
- [21:48] Datapoint id 18 collected
- [21:48] Charge balance equation is out by 8.3%
- [21:48] Stirrer speed set to 55
- [21:53] pH 2.16 -> 2.36
- [21:53] Using charge balance adjust
- [21:53] Dispensed 0.011759 mL of Base (0.5 M KOH)
- [22:13] Stirrer speed set to 0
- [22:24] Datapoint id 19 collected
- [22:24] Charge balance equation is out by 2.3%
- [22:24] Stirrer speed set to 55
- [22:29] pH 2.38 -> 2.58
- [22:29] Using charge balance adjust
- [22:29] Dispensed 0.007314 mL of Base (0.5 M KOH)
- [22:49] Stirrer speed set to 0
- [22:59] Datapoint id 20 collected
- [22:59] Charge balance equation is out by 0.9%
- [22:59] Stirrer speed set to 55
- [23:04] pH 2.59 -> 2.79
- [23:04] Using charge balance adjust
- [23:05] Dispensed 0.004774 mL of Base (0.5 M KOH)
- [23:25] Stirrer speed set to 0
- [23:35] Datapoint id 21 collected
- [23:35] Charge balance equation is out by 2.4%
- [23:35] Stirrer speed set to 55
- [23:40] pH 2.81 -> 3.01
- [23:40] Using charge balance adjust
- [23:40] Dispensed 0.003269 mL of Base (0.5 M KOH)
- [24:01] Stirrer speed set to 0
- [24:11] Datapoint id 22 collected
- [24:11] Charge balance equation is out by -4.9%



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [24:11] Stirrer speed set to 55
- [24:16] pH 3.01 -> 3.21
- [24:16] Using charge balance adjust
- [24:16] Dispensed 0.002469 mL of Base (0.5 M KOH)
- [24:36] Stirrer speed set to 0
- [24:46] Datapoint id 23 collected
- [24:46] Charge balance equation is out by -3.0%
- [24:46] Stirrer speed set to 55
- [24:51] pH 3.21 -> 3.41
- [24:51] Using charge balance adjust
- [24:51] Dispensed 0.001952 mL of Base (0.5 M KOH)
- [25:11] Stirrer speed set to 0
- [25:22] Datapoint id 24 collected
- [25:22] Charge balance equation is out by 54.8%
- [25:22] Stirrer speed set to 55
- [25:27] pH 3.53 -> 3.73
- [25:27] Using cautious pH adjust
- [25:27] Dispensed 0.000706 mL of Base (0.5 M KOH)
- [25:32] Stepping pH = 3.68
- [25:32] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [25:37] Stepping pH = 3.73
- [25:52] Stirrer speed set to 0
- [26:03] Datapoint id 25 collected
- [26:03] Charge balance equation is out by 34.0%
- [26:03] Stirrer speed set to 55
- [26:08] pH 3.72 -> 3.92
- [26:08] Using cautious pH adjust
- [26:09] Dispensed 0.000564 mL of Base (0.5 M KOH)
- [26:14] Stepping pH = 3.85
- [26:14] Dispensed 0.000188 mL of Base (0.5 M KOH)
- [26:19] Stepping pH = 3.89
- [26:19] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [26:24] Stepping pH = 3.90
- [26:24] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [26:29] Stepping pH = 3.92
- [26:44] Stirrer speed set to 0
- [26:56] Datapoint id 26 collected
- [26:56] Charge balance equation is out by 15.5%
- [26:56] Stirrer speed set to 55
- [27:01] pH 3.89 -> 4.09
- [27:01] Using cautious pH adjust
- [27:01] Dispensed 0.000447 mL of Base (0.5 M KOH)
- [27:06] Stepping pH = 4.03
- [27:06] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [27:11] Stepping pH = 4.06
- [27:11] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [27:16] Stepping pH = 4.09
- [27:32] Stirrer speed set to 0
- [27:45] Datapoint id 27 collected
- [27:45] Charge balance equation is out by 19.9%
- [27:45] Stirrer speed set to 55
- [27:50] pH 4.06 -> 4.26
- [27:50] Using cautious pH adjust
- [27:50] Dispensed 0.000329 mL of Base (0.5 M KOH)
- [27:55] Stepping pH = 4.19
- [27:55] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [28:00] Stepping pH = 4.23
- [28:00] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [28:05] Stepping pH = 4.25



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [28:20] Stirrer speed set to 0
- [28:30] Datapoint id 28 collected
- [28:30] Charge balance equation is out by 21.2%
- [28:30] Stirrer speed set to 55
- [28:35] pH 4.24 -> 4.44
- [28:35] Using cautious pH adjust
- [28:35] Dispensed 0.000235 mL of Base (0.5 M KOH)
- [28:41] Stepping pH = 4.39
- [28:41] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [28:46] Stepping pH = 4.39
- [28:46] Dispensed 0.000165 mL of Base (0.5 M KOH)
- [28:51] Stepping pH = 4.52
- [29:06] Stirrer speed set to 0
- [29:18] Datapoint id 29 collected
- [29:18] Charge balance equation is out by -0.7%
- [29:18] Stirrer speed set to 55
- [29:23] pH 4.50 -> 4.70
- [29:23] Using charge balance adjust [29:23] Dispensed 0.000282 mL of Base (0.5 M KOH)
- [29:43] Stirrer speed set to 0
- [30:07] Datapoint id 30 collected
- [30:07] Charge balance equation is out by 24.6%
- [30:07] Stirrer speed set to 55
- [30:12] pH 4.75 -> 4.95
- [30:12] Using cautious pH adjust
- [30:12] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [30:17] Stepping pH = 4.80
- [30:17] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [30:22] Stepping pH = 5.09
- [30:37] Stirrer speed set to 0
- [30:50] Datapoint id 31 collected
- [30:50] Charge balance equation is out by -26.9%
- [30:50] Stirrer speed set to 55
- [30:55] pH 5.05 -> 5.25
- [30:55] Using cautious pH adjust
- [30:55] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [31:00] Stepping pH = 5.07
- [31:00] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [31:05] Stepping pH = 5.37
- [31:21] Stirrer speed set to 0
- [31:34] Datapoint id 32 collected
- [31:34] Charge balance equation is out by -86.8%
- [31:34] Stirrer speed set to 55
- [31:39] pH 5.41 -> 5.61
- [31:39] Using cautious pH adjust
- [31:39] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [31:44] Stepping pH = 5.40
- [31:44] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [31:49] Stepping pH = 6.97
- [32:05] Stirrer speed set to 0
- [33:05] Datapoint id 33 collected
- [33:05] Charge balance equation is out by -201.9%
- [33:05] Stirrer speed set to 55
- [33:10] pH 6.45 -> 6.65
- [33:10] Using cautious pH adjust
- [33:10] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [33:15] Stepping pH = 6.42
- [33:15] Dispensed 0.000118 mL of Base (0.5 M KOH)
- [33:20] Stepping pH = 8.04



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [33:35] Stirrer speed set to 0
- [34:35] Datapoint id 34 collected
- [34:35] Charge balance equation is out by -241.6%
- [34:35] Stirrer speed set to 55
- [34:40] pH 7.57 -> 7.77
- [34:40] Using cautious pH adjust
- [34:40] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [34:45] Stepping pH = 7.53
- [34:45] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [34:51] Stepping pH = 7.84
- [35:06] Stirrer speed set to 0
- [36:06] Datapoint id 35 collected
- [36:06] Charge balance equation is out by -393.8%
- [36:06] Stirrer speed set to 55
- [36:11] pH 7.75 -> 7.95
- [36:11] Using cautious pH adjust
- [36:11] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [36:16] Stepping pH = 7.80
- [36:16] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [36:21] Stepping pH = 8.02
- [36:36] Stirrer speed set to 0
- [37:36] Datapoint id 36 collected
- [37:36] Charge balance equation is out by -332.8%
- [37:36] Stirrer speed set to 55
- [37:41] pH 8.03 -> 8.23
- [37:41] Using cautious pH adjust
- [37:41] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [37:46] Stepping pH = 8.20
- [37:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [37:52] Stepping pH = 8.36
- [38:07] Stirrer speed set to 0
- [38:57] Datapoint id 37 collected
- [38:57] Charge balance equation is out by -361.1%
- [38:57] Stirrer speed set to 55
- [39:02] pH 8.24 -> 8.44
- [39:02] Using cautious pH adjust
- [39:02] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [39:07] Stepping pH = 8.35
- [39:07] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [39:13] Stepping pH = 8.60
- [39:28] Stirrer speed set to 0
- [40:02] Datapoint id 38 collected
- [40:02] Charge balance equation is out by -280.6%
- [40:02] Stirrer speed set to 55
- [40:07] pH 8.66 -> 8.86
- [40:07] Using cautious pH adjust
- [40:07] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [40:12] Stepping pH = 8.70
- [40:12] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [40:17] Stepping pH = 8.80
- [40:17] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [40:22] Stepping pH = 8.89
- [40:37] Stirrer speed set to 0
- [41:02] Datapoint id 39 collected
- [41:02] Charge balance equation is out by -172.2%
- [41:02] Stirrer speed set to 55
- [41:07] pH 8.88 -> 9.05
- [41:07] Using cautious pH adjust
- [41:07] Dispensed 0.000024 mL of Base (0.5 M KOH)



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [41:12] Stepping pH = 8.90
- [41:12] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [41:17] Stepping pH = 8.98
- [41:17] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [41:22] Stepping pH = 9.07
- [41:37] Stirrer speed set to 0
- [41:56] Datapoint id 40 collected
- [41:56] Charge balance equation is out by -170.5%
- [41:56] Titration 3 of 3
- [41:56] Adding initial titrants
- [41:56] Automatically add 0.70000 mL of Octanol
- [42:44] Dispensed 0.700000 mL of Octanol
- [42:44] Stirrer speed set to 10
- [42:45] Stirrer speed set to 60
- [42:45] Iterative adjust 9.02 -> 2.00
- [42:45] pH 9.02 -> 2.00
- [42:47] Dispensed 0.054727 mL of Acid (0.5 M HCl)
- [42:52] pH 2.03 -> 2.00
- [42:52] Dispensed 0.003175 mL of Acid (0.5 M HCI)
- [43:42] Stirrer speed set to 0
- [43:52] Datapoint id 41 collected
- [43:52] Stirrer speed set to 60
- [43:57] pH 1.96 -> 2.16
- 10.67] pri 1.66 2.16
- [43:57] Using cautious pH adjust
- [43:58] Dispensed 0.009972 mL of Base (0.5 M KOH)
- [44:03] Stepping pH = 2.05
- [44:03] Dispensed 0.007596 mL of Base (0.5 M KOH)
- [44:08] Stepping pH = 2.14
- [44:08] Dispensed 0.001317 mL of Base (0.5 M KOH)
- [44:13] Stepping pH = 2.16
- [44:28] Stirrer speed set to 0
- [44:39] Datapoint id 42 collected
- [44:39] Charge balance equation is out by 5.2%
- [44:39] Stirrer speed set to 60
- [44:44] pH 2.16 -> 2.36
- [44:44] Using charge balance adjust
- [44:44] Dispensed 0.012817 mL of Base (0.5 M KOH)
- [45:05] Stirrer speed set to 0
- [45:23] Datapoint id 43 collected
- [45:23] Charge balance equation is out by 5.4%
- [45:23] Stirrer speed set to 60
- [45:28] pH 2.38 -> 2.58
- [45:28] Using charge balance adjust
- [45:28] Dispensed 0.007973 mL of Base (0.5 M KOH)
- [45:48] Stirrer speed set to 0
- [45:58] Datapoint id 44 collected
- [45:58] Charge balance equation is out by 3.7%
- [45:58] Stirrer speed set to 60
- [46:03] pH 2.59 -> 2.79
- [46:03] Using charge balance adjust
- [46:04] Dispensed 0.005292 mL of Base (0.5 M KOH)
- [46:24] Stirrer speed set to 0
- [46:34] Datapoint id 45 collected
- [46:34] Charge balance equation is out by 12.0%
- [46:34] Stirrer speed set to 60
- [46:39] pH 2.82 -> 3.02
- [46:39] Using charge balance adjust
- [46:40] Dispensed 0.003575 mL of Base (0.5 M KOH)
- [47:00] Stirrer speed set to 0



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [47:10] Datapoint id 46 collected
- [47:10] Charge balance equation is out by -4.8%
- [47:10] Stirrer speed set to 60
- [47:15] pH 3.02 -> 3.22
- [47:15] Using charge balance adjust
- [47:15] Dispensed 0.002634 mL of Base (0.5 M KOH)
- [47:35] Stirrer speed set to 0
- [47:45] Datapoint id 47 collected
- [47:45] Charge balance equation is out by 7.7%
- [47:45] Stirrer speed set to 60
- [47:50] pH 3.25 -> 3.45
- [47:50] Using charge balance adjust
- [47:51] Dispensed 0.001881 mL of Base (0.5 M KOH)
- [48:11] Stirrer speed set to 0
- [48:21] Datapoint id 48 collected
- [48:21] Charge balance equation is out by 17.5%
- [48:21] Stirrer speed set to 60
- [48:26] pH 3.49 -> 3.69
- [48:26] Using cautious pH adjust
- [48:26] Dispensed 0.000635 mL of Base (0.5 M KOH)
- [48:31] Stepping pH = 3.59
- [48:31] Dispensed 0.000400 mL of Base (0.5 M KOH)
- [48:36] Stepping pH = 3.66
- [48:36] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [48:41] Stepping pH = 3.68
- [48:57] Stirrer speed set to 0
- [49:07] Datapoint id 49 collected
- 49.07] Datapoint id 49 collected
- [49:07] Charge balance equation is out by 7.9%
- [49:07] Stirrer speed set to 60
- [49:12] pH 3.68 -> 3.88
- [49:12] Using charge balance adjust
- [49:12] Dispensed 0.000917 mL of Base (0.5 M KOH)
- [49:32] Stirrer speed set to 0
- [49:51] Datapoint id 50 collected
- [49:51] Charge balance equation is out by 2.2%
- [49:51] Stirrer speed set to 60
- [49:56] pH 3.90 -> 4.10
- [49:56] Using charge balance adjust
- [49:56] Dispensed 0.000588 mL of Base (0.5 M KOH)
- [50:16] Stirrer speed set to 0
- [50:28] Datapoint id 51 collected
- [50:28] Charge balance equation is out by -30.8%
- [50:28] Stirrer speed set to 60
- [50:33] pH 4.05 -> 4.25
- [50:33] Using cautious pH adjust
- [50:33] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [50:38] Stepping pH = 4.09
- [50:38] Dispensed 0.000353 mL of Base (0.5 M KOH)
- [50:43] Stepping pH = 4.29
- [50:58] Stirrer speed set to 0
- [51:08] Datapoint id 52 collected
- [51:08] Charge balance equation is out by -28.1%
- [51:08] Stirrer speed set to 60
- [51:14] pH 4.27 -> 4.47
- [51:14] Using cautious pH adjust
- [51:14] Dispensed 0.000141 mL of Base (0.5 M KOH)
- [51:19] Stepping pH = 4.30
- [51:19] Dispensed 0.000306 mL of Base (0.5 M KOH)
- [51:24] Stepping pH = 4.59



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [51:39] Stirrer speed set to 0
- [51:49] Datapoint id 53 collected
- [51:49] Charge balance equation is out by -65.4%
- [51:49] Stirrer speed set to 60
- [51:54] pH 4.58 -> 4.78
- [51:54] Using cautious pH adjust
- [51:54] Dispensed 0.000071 mL of Base (0.5 M KOH)
- [51:59] Stepping pH = 4.59
- [51:59] Dispensed 0.000212 mL of Base (0.5 M KOH)
- [52:04] Stepping pH = 4.99
- [52:19] Stirrer speed set to 0
- [52:29] Datapoint id 54 collected
- [52:29] Charge balance equation is out by -96.2%
- [52:29] Stirrer speed set to 60
- [52:35] pH 4.97 -> 5.17
- [52:35] Using cautious pH adjust
- [52:35] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [52:40] Stepping pH = 4.97
- [52:40] Dispensed 0.000094 mL of Base (0.5 M KOH)
- [52:45] Stepping pH = 5.07
- [52:45] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [52:50] Stepping pH = 5.19
- 52.30] Stepping pri = 3.19
- [53:05] Stirrer speed set to 0
- [53:15] Datapoint id 55 collected
- [53:15] Charge balance equation is out by -183.3%
- [53:15] Stirrer speed set to 60
- [53:20] pH 5.19 -> 5.39
- [53:20] Using cautious pH adjust
- [53:20] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [53:25] Stepping pH = 5.20
- [53:25] Dispensed 0.000047 mL of Base (0.5 M KOH)
- [53:31] Stepping pH = 5.35
- [53:31] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [53:36] Stepping pH = 5.46
- [53:51] Stirrer speed set to 0
- [54:01] Datapoint id 56 collected
- [54:01] Charge balance equation is out by -136.9%
- [54:01] Stirrer speed set to 60
- [54:06] pH 5.47 -> 5.67
- [54:06] Using cautious pH adjust
- [54:06] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [54:11] Stepping pH = 5.50
- [54:11] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [54:16] Stepping pH = 5.63
- [54:16] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [54:21] Stepping pH = 5.79
- [54:36] Stirrer speed set to 0
- [54:47] Datapoint id 57 collected
- [54:47] Charge balance equation is out by -142.7%
- [54:47] Stirrer speed set to 60
- [54:52] pH 5.85 -> 6.05
- [54:52] Using cautious pH adjust
- [54:52] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [54:57] Stepping pH = 5.91
- [54:57] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [55:02] Stepping pH = 6.14
- [55:17] Stirrer speed set to 0
- [55:38] Datapoint id 58 collected
- [55:38] Charge balance equation is out by -37.8%
- Reported at: 3/30/2018 5:50:25 PM



pH-metric high logP Assay name: Analyst: Pion Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [55:38] Stirrer speed set to 60
- [55:43] pH 6.09 -> 6.29
- [55:43] Using cautious pH adjust
- [55:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [55:48] Stepping pH = 6.20
- [55:48] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [55:53] Stepping pH = 6.54
- [56:08] Stirrer speed set to 0
- [56:42] Datapoint id 59 collected
- [56:42] Charge balance equation is out by -15.1%
- [56:42] Stirrer speed set to 60
- [56:47] pH 6.56 -> 6.76
- [56:47] Using cautious pH adjust
- [56:47] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [56:52] Stepping pH = 6.65
- [56:52] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [56:57] Stepping pH = 6.97
- [57:12] Stirrer speed set to 0
- [58:11] Datapoint id 60 collected
- [58:11] Charge balance equation is out by 3.8%
- [58:11] Stirrer speed set to 60
- [58:16] pH 6.99 -> 7.19
- [58:16] Using charge balance adjust
- [58:17] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [58:37] Stirrer speed set to 0
- [59:19] Datapoint id 61 collected
- [59:19] Charge balance equation is out by -18.2%
- [59:19] Stirrer speed set to 60
- [59:24] pH 6.94 -> 7.14
- [59:24] Using cautious pH adjust
- [59:24] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [59:29] Stepping pH = 6.99
- [59:29] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [59:34] Stepping pH = 7.34
- [59:49] Stirrer speed set to 0
- [1:00:49] Datapoint id 62 collected
- [1:00:49] Charge balance equation is out by -30.8%
- [1:00:49] Stirrer speed set to 60
- [1:00:54] pH 7.34 -> 7.54
- [1:00:54] Using cautious pH adjust
- [1:00:55] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:01:00] Stepping pH = 7.38
- [1:01:00] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:01:05] Stepping pH = 7.84
- [1:01:20] Stirrer speed set to 0
- [1:02:20] Datapoint id 63 collected
- [1:02:20] Charge balance equation is out by -73.9%
- [1:02:20] Stirrer speed set to 60
- [1:02:25] pH 7.88 -> 8.08
- [1:02:25] Using cautious pH adjust
- [1:02:25] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:02:30] Stepping pH = 7.97
- [1:02:30] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:02:35] Stepping pH = 8.26 [1:02:50] Stirrer speed set to 0
- [1:03:33] Datapoint id 64 collected
- [1:03:33] Charge balance equation is out by -237.4%
- [1:03:33] Stirrer speed set to 60
- [1:03:38] pH 8.20 -> 8.40



Assay name: pH-metric high logP Analyst: Pion
Assay ID: 18C-26010 Instrument ID: T312060

Filename: C:\Sirius_T3\Mehtap\20180327_exp34_logP_T3-2\18C-26010_M13_octanol_pH-metric high logP.t3r

- [1:03:38] Using cautious pH adjust
- [1:03:38] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:03:43] Stepping pH = 8.28
- [1:03:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- 1:03:48 Stepping pH = 8.53
- [1:04:03] Stirrer speed set to 0
- [1:04:38] Datapoint id 65 collected
- [1:04:38] Charge balance equation is out by -230.0%
- [1:04:38] Stirrer speed set to 60
- [1:04:43] pH 8.49 -> 8.69
- [1:04:43] Using cautious pH adjust
- [1:04:43] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:04:48] Stepping pH = 8.55
- [1:04:48] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:04:53] Stepping pH = 8.71
- [1:05:08] Stirrer speed set to 0
- [1:05:37] Datapoint id 66 collected
- [1:05:37] Charge balance equation is out by -127.2%
- [1:05:37] Stirrer speed set to 60
- [1:05:42] pH 8.71 -> 8.91
- [1:05:42] Using cautious pH adjust
- [1:05:42] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:05:47] Stepping pH = 8.74
- [1:05:47] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:05:53] Stepping pH = 8.87
- [1:05:53] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:05:58] Stepping pH = 8.98
- [1:06:13] Stirrer speed set to 0
- [1:06:41] Datapoint id 67 collected
- [1:06:41] Charge balance equation is out by -152.4%
- [1:06:41] Stirrer speed set to 60
- [1:06:46] pH 8.97 -> 9.05
- [1:06:46] Using cautious pH adjust
- [1:06:46] Dispensed 0.000024 mL of Base (0.5 M KOH)
- [1:06:51] Stepping pH = 9.00
- [1:07:06] Stirrer speed set to 0
- [1:07:33] Datapoint id 68 collected
- [1:07:33] Charge balance equation is out by -27.6%
- [1:07:33] Argon flow rate set to 0
- [1:07:37] Titrator arm moved over Titration position
- [1:07:58] The autoloader failed to pick at location "Sample position"